

EASY SETUP

EN

IN THE BOX:



Nursery Unit



Sensor Pad



1 AC Adapter



4 AAA Alkaline Batteries

NOT INCLUDED:



Prepare one hardboard
Thickness: **6 mm (1/4")**
Minimum dimensions:
33 cm (13") x 33 cm (13")
Maximum dimensions:
equal to the mattress size

IMPORTANT

The hardboard
must be larger
than Sensor Pad
surface area.



Miniature
Philips
Screwdriver

STEP 1 - INSTALL APPLICATION

- 1) Scan **QR Code** with your smartphone
OR Search for **Angelcare baby monitor**
in your App Store.
- 2) Open application and accept
PUSH Notifications.
- 3) Sign up and create your account.

.....
Shortly after signing up, Angelcare
will send you a automatically generated
password by email.

1.1



OR




1.3




STEP 2 - INSTALL NURSERY UNIT

- 1) Plug in your Nursery Unit.
- 2) Press Nursery Unit Power Button
and wait for activation sound (30 seconds).
.....
Power light will turn from white to red
and from red to white until a twinkling sound
confirms complete activation.
- 3) Rotate Camera Button to top position.


2.1



2.2



2.3



STEP 3 - SET UP CAMERA


VIA WIRELESS SETUP

- 1) Select arrow to add camera.
- 2) Type in your network name and password. Use your wireless network name and password to successfully set up your camera(s).
- 3) **QR Code** will be displayed on your device screen.



STEP 4 - CONNECT

- 1) Press and hold **Neck Button** until you hear a twinkling sound (3 seconds).
- 2) Place **Nursery Unit** and device approximately 6 inch apart until you hear a twinkling sound. (Align QR code with camera lens)
- 3) Wait for **30 seconds** for camera to connect. Another twinkling sound will confirm camera connection: **Success Message** will appear on your screen.




NOTE:

If the connection is not successful, you will hear a failed sound. Please try to connect again.

Power light will turn red during the connecting waiting period.

STEP 5 - EXPLORE THE APPLICATION

- 1) Edit your camera name.
- 2) Select your camera.
- 3) Modify your settings.



OR (IF YOU DON'T HAVE AN INTERNET CONNECTION)

VIA ACCESS POINT SETUP (No Internet connection needed)

- 1) Set camera to AP mode (switch underneath camera)
- 2) Select Angelcare camera via your WiFi settings device using password: 12345678
- 3) Restart the app (if it does not work immediately)



FEDERAL COMMUNICATIONS COMMISSION INTERFERENCE STATEMENT

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

CAUTION:

Any changes or modifications not expressly approved by the grantee of this device could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF exposure warning

This equipment must be installed and operated in accordance with provided instructions and the antenna(s) used for this transmitter must be installed to provide a separation distance of at least 20 cm from all persons and must not be co-located or operating in conjunction with any other antenna or transmitter. End-users and installers must be provide with antenna installation instructions and transmitter operating conditions for satisfying RF exposure compliance.

Canada, Industry Canada (IC) Notices

This Class B digital apparatus complies with Canadian ICES-003 and RSS-210.

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Canada, avis d'Industry Canada (IC)

Cet appareil numérique de classe B est conforme aux normes canadiennes ICES-003 et RSS-210.

Son fonctionnement est soumis aux deux conditions suivantes : (1) cet appareil ne doit pas causer d'interférence et (2) cet appareil doit accepter toute interférence, notamment les interférences qui peuvent affecter son fonctionnement.

Radio Frequency (RF) Exposure Information

The radiated output power of the Wireless Device is below the Industry Canada (IC) radio frequency exposure limits. The Wireless Device should be used in such a manner such that the potential for human contact during normal operation is minimized.

This device has also been evaluated and shown compliant with the IC RF Exposure limits under mobile exposure conditions. (antennas are greater than 20cm from a person's body).



Informations concernant l'exposition aux fréquences radio (RF)

La puissance de sortie émise par l'appareil de sans fil est inférieure à la limite d'exposition aux fréquences radio d'Industry Canada (IC). Utilisez l'appareil de sans fil de façon à minimiser les contacts humains lors du fonctionnement normal.

Ce périphérique a également été évalué et démontré conforme aux limites d'exposition aux RF d'IC dans des conditions d'exposition à des appareils mobiles (les antennes se situent à moins de 20 cm du corps d'une personne).

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and**
- (2) this device must accept any interference received, including interference that may cause undesired operation**



NURSERY UNIT
Rotate camera to top position

Sensor Pad Introduction

Spec:

1. Dimension: 240*240*10mm
2. 2M length cable with 3.5mm pin
3. Piezoelectricity to detect vital signs

Function:

1. Connect to Babycam by 3.5mm pin.
2. Baby lies down on the sensor pad, once baby without vital signs, pad will send signal to Babycam for alarm.