97050096  
rev. 000  
00/2007

ITALIANO ..... 2

ENGLISH..... 23

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## Indice

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# OPERATING INSTRUCTIONS

## Contents

1. General precautions .....	24
1.1. Description of the symbols used .....	24
1.2. Purpose and use of the equipment .....	25
1.2.1. Classification .....	25
1.2.2. Environmental conditions .....	25
1.2.3. Guarantee .....	25
2. Description of the equipment .....	26
2.1. Description of the x-ray unit .....	26
3. Operation .....	26
3.1. Turning The X-ray equipment on and off .....	26
3.1.1. Turning on the basic X-ray unit .....	26
3.1.2. Turning on the handheld .....	27
3.1.3. Control panel .....	28
3.1.4. Factory settings .....	28
3.1.5. Automatic handheld shut off .....	29
3.1.6. Hand-held stand-by time .....	29
3.2. Checking the set parameters .....	29
3.3. Batteries and charge level status .....	30
3.4. Position of the patient .....	30
3.5. Putting the x-ray unit cone into the required position .....	30
3.6. Position of the X-ray plate or sensor .....	31
3.7. Checking the exposure time on the display .....	32
3.8. Procedure to be followed when taking the x-ray .....	32
3.9. Setting the exposure time .....	33
4. Advanced options .....	33
5. Routine maintenance .....	35
6. Cleaning and disinfection .....	35
7. Dismantling the equipment when no longer used .....	35
8. Specifications .....	35
9. Specifications .....	36
9.1. Dimensions .....	36
9.2. Identification nameplates .....	38
9.3 Fault table .....	40

## 1. General precautions

These instructions describe how to use the DC x-ray unit correctly. Please read this manual carefully before using the equipment.

 **NOTE:** This manual does not contain all the recommendations and rules on holding a source of ionogenic radiation as these vary from place to place. Only the most common ones shall be mentioned and this means that it is the user's responsibility to check local standards and observe the relevant laws.

No part of this manual is to be reproduced, stored in a retrieval system or transmitted in any form or by any means, i.e. electronic, mechanical, photocopying, translation or otherwise, without the prior written permission of the manufacturer.

CEFLA SC - CEFLA DENTAL GROUP has a company policy of continual development. Although every effort is made to keep technical documentation up-to-date at all times some of the instructions, specifications and figures given in this manual may slightly differ from the purchased product. The manufacturer reserves the right to make changes without prior notice.

The original version of this manual is written in Italian.

### 1.1. Description of the symbols used

The following symbols are used:

	Type of protection against direct and indirect contact: <b>Class I</b> . Level of protection against direct and indirect contact: <b>Type B</b> .
	<b>WARNING!</b> Indicates a situation where failure to observe the instructions could cause damage to the equipment or injury to the user and/or patient.
	<b>NOTE:</b> Indicates important information for the user and/or technical support personnel.
	Earth connection.
	Alternating current.
	On.
	Off.
	Non-ionizing radiation.
	ionizing radiation.
	Equipment that meets the essential requirements established by Directive EC 99/05 art.12 (Class I devices).
	The equipment conforms to the requirements stated in EEC directive 93/42.
	Equipment abides by the essential requirements established by USA and Canada.
	I.M.Q. mark (Italian Institute for the Quality Mark)
	Symbol indicating separate collection for electrical and electronic equipment in compliance with Directives 2002/95/ EC and 2003/108/ EC.
FCC ID	F.C.C. mark (Federal Communication Commission) *

\* This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## 1.2. Purpose and use of the equipment

This x-ray unit is designed for use in the dental surgery to make endo-oral x-rays for diagnostic purposes. This equipment can be used to produce traditional x-rays developed using chemicals or, alternatively, it can be used with digital x-ray sensors.

### 1.2.1. Classification

#### MEDICAL EQUIPMENT Classification

Classification of the equipment according to the rules indicated in Annex IX of EEC Directive 93/42: **Class IIb**.

#### ELECTRO-MEDICAL EQUIPMENT Classification

Classification of the equipment according to standard I.E.C. 60601-1 on the safety of medical equipment: **Class I - Type B**.

#### RADIO EQUIPMENT AND TELECOMMUNICATIONS TERMINAL EQUIPMENT Classification

Equipment classification in accordance with directive 99/05/EEC art.12: **Class I**.

#### EMC classification

Classification of the equipment according to standard CEI EN 55011: **Group 1 - Type B**.

### 1.2.2. Environmental conditions

The equipment should be installed in an area where the following conditions are present:

- Temperature from +10° to +40° C;
- Storage temperature: -40°C to +70°C;
- Temperatura di immagazzinamento: -40°C a +70°C;
- relative humidity from 25 to 75% without condensate;
- atmospheric pressure from 700 to 1060 hPa;
- The electrical wiring in the room where the equipment is installed must conform to the I.E.C. 60364-7-710;V2 specification (i.e. the regulations concerning the electrical wiring to be used in surgeries) or equivalent standards in force in the country where the equipment is installed.

### 1.2.3. Guarantee

CEFLA SC - CEFLA DENTAL GROUP guarantees the safety of the operator as well as the reliability and high performance of the equipment. The guarantee shall be voided if the following precautions are not taken:

- Observe the conditions specified in the guarantee certificate itself.
- The equipment must only be used by following the instructions contained herein.
- Equipment installation, expansion and technical support must be performed exclusively by personnel authorized by the manufacturer to carry out these operations.
- Never open the equipment casing. Installation, repairs and, in general, any other operations requiring the casing to be opened are to be performed exclusively by personnel authorized by the manufacturer to carry out these operations.
- The equipment must be installed in areas which comply with that stated above at point 1.2.2 "Ambient Conditions".
- The area where the x-ray unit is installed must comply with official regulations regarding protection against radiation in the country where the equipment is used.

#### **SAFETY PRECAUTIONS**

- If any person who is not an authorized technician changes the product in any way by replacing parts or components with other ones not used by the manufacturer they shall assume responsibility for the product.
- Do not forget to turn off the main switch on the equipment before leaving the surgery.
- This equipment is not water-proof (risk of electrocution).
- This equipment is not suitable for use in the presence of a mix of inflammable anaesthetic gas with oxygen or nitrous oxide.
- This equipment must be stored properly so that it is kept in top working order at all times. The manufacturer shall not be held responsible for misuse, carelessness or improper use of the equipment.
- This equipment is to be used exclusively by qualified personnel with the proper medical training.
- The user must be present at all times when the equipment is turned on or ready for start-up. In particular, never leave the equipment unattended in the presence of children/the mentally disabled or other unauthorised personnel in general.
- If the x-ray equipment is damaged or oil leaks, do not use the equipment and contact customer service immediately.

#### **PROTECTION AGAINST RADIATION**


**X-rays are dangerous and it is therefore essential to take suitable safety precautions before using them. Areas where it is possible to be exposed to the x-rays shall be clearly indicated by using this symbol which should remind personnel to observe the safety rules laid down by the laws in force in the country where the equipment is used.**

- Control the emission of x-rays from the greatest distance possible (several metres) from the focal spot and the X-ray irradiation beam in the opposite direction to where the rays are emitted.
- Only the authorized personnel and the patient can remain in the area when x-rays are being emitted.
- Always protect the patient's thyroid and gonads under all circumstances.
- **Whenever the patient is a child or disabled person requiring the presence of the dentist to keep the film in position it is advisable to use a pair of tongs and a special glove to protect the hand against x-rays. Use a suitable overall to protect the rest of the body against exposure to x-rays.**

## 2. Description of the equipment

### 2.1. Description of the x-ray unit

Description of the equipment:

- **Double pantograph arm:** This is a double-jointed arm to which a further extension arm can be added. The extension arm may be: 40 cm, 60 cm or 90 cm.

- **Generator with collimator ( a + b ).** The generator is of the constant power type 60kV / 7-3,5mA. The generator can turn endlessly on the horizontal plane. On the other hand, rotation on the vertical plane is restricted upwards by a mechanical end-stop.

The X-ray generator has two different beam collimations recognized in automatic mode:

- **8" CYLINDRICAL COLLIMATOR:**  
minimum focus/skin distance 20 cm and size of x-ray beam leaving collimator 60 mm (with rectangular collimator shut off).
- **12" RECTANGULAR COLLIMATOR:**  
minimum focus/skin distance 30 cm and size of x-ray beam leaving collimator 45x35 mm (with rectangular collimator turned on).

- **Handheld:** the handheld can be placed near the control unit itself or in a remote location. This allows the dentist to move around the office freely and go out of the radiation area.


## 3. Operation

### 3.1. Turning The X-ray equipment on and off

#### 3.1.1. Turning on the basic X-ray unit


The control unit is turned on and off with the main switch (A), as illustrated in the figure below. The switch lights up when the control unit is energized. Whenever turned on, the equipment performs an operational test that takes a few seconds. Once the test has been completed, a buzzer rings and the indicator light (B) on the X-ray generator lights up at the same time.

**NOTE:** The exposure time and the parameters displayed when the unit is turned on are the last ones set before the central control unit was turned off. If the central control unit is left untouched for a few minutes it will go into stand-by mode. Simply press any key on the control panel to reactivate it.


### 3.1.2. Turning on the handheld

The handheld is turned on by pressing any key, except for the one for x-ray emission. A buzzer rings to confirm the apparatus has been turned on. The unit will be in the standard configuration described in detail in paragraph 3.1.3 and then search for the base it works with. If the base is off, the handheld will not indicate the field nor status "ready". If the base is latter turned on, the handheld will detect it within thirty seconds or by pressing any function key on the pushbutton panel.

### 3.1.3. Control panel

As illustrated in the figure below, the handheld has four function keys and the exposure key.

The main functions of the keys on the handheld vary according to how they are pressed:

KEY	BRIEFLY PRESSED (less than 3 sec.)	PRESSED LONGER (more than 3 sec.)
	Changes over from ADULT to CHILD and vice versa (takes place when key is released).	Saves, if permitted, the sensitivity of the new time selected. The memo icon () lights up when the data item can be saved.
	Changes amongst the various types of teeth to select the area to be examined.	Displays the exposition time of the tooth in mGy and, if the key is held down a few more seconds, in mGy*cm2.
	Increases the exposure times in steps, according to the set scale.	Increases the scroll speed of the values in increasing order.
	Increases the exposure times in steps, according to the set scale.	Increases the scroll speed of the values in decreasing order.
	NO EFFECTS ARE OBTAINED IF THE KEY IS PRESSED LESS THAN A SECOND.	STARTS X-RAY EXPOSURE (the button has to be held down while the x-rays are being emitted, "dead man" function).

**NOTE: "Dead man" function:** the system that starts x-ray exposure with the dedicated key on the wireless handheld allows x-rays to be emitted only when the user presses and holds down the exposure key. X-ray emission will stop if the key is released ahead of time.

**NOTE:** The function related to pressing the key briefly is performed by pressing the key which will activate the function assigned to it. On the other hand, to perform the function carried out when the key is held down longer, press the key until the relative function is started. The buzzer will ring shortly to signal the function has started.

### 3.1.4. Factory settings

The x-ray unit is supplied with the following factory settings:

- Anode current: 7 mA (NORM mode).
- Sensitivity: level 19.
- Handheld stand by: 5 minutes
- Patient's built: adult (ADULT symbol selected).
- Collimator presence signaled on the display (collimator symbol off if the rectangular collimator is not turned on in the head).
- Exposure times according to standard R20: 0,010-0,011-0,012-0,014-0,016-0,018-0,020-0,022-0,025-0,028-0,032-0,036-0,040-0,045-0,050-0,056-0,063-0,071-0,080-0,090-0,100-0,112-0,125-0,140-0,160-0,180-0,200-0,224-0,250-0,280-0,315-0,355-0,400-0,500-0,560-0,630-0,710-0,800-0,900-1,000.

**NOTE:**

These times comply with current standards I.E.C. 60601-2-7 (1999) and the ISO 497 series R'20 recommendations and CANNOT BE MODIFIED.

### 3.1.5. Automatic handheld shut off

About one minute after the base is turned off the handheld automatically shuts off. Similarly, the handheld will automatically shut off if it is far away or in any case outside the operating range for interfacing with the base.

### 3.1.6. Hand-held stand-by time

The entire x-ray unit will switch over to stand-by (even if the base is on) and the handheld will automatically shut off after approximately five minutes of non-use to save battery power. Press any key, except for the “exposure” key, to turn the handheld back on showing the last selection made by the user. To change the stand by time, refer to chapter 4 regarding the handheld’s “Advanced options”.

## 3.2. Checking the set parameters

Before actually taking an exposure, make sure the exposure parameters for the examination in progress are correctly set.

### - Checking the collimator used.

The icon on the handheld’s screen should be on or off, depending on the operating mode selected:

**ICON**  **ON:** indicates that the cylindrical collimator is activated.


**ICON**  **OFF:** indicates the rectangular collimator is activated.

**NOTE:** After turning the collimator on or off, the preset exposure times and icon will automatically be modified within a few seconds.

### - Checking the set patient’s build.

- **“Child” selected:** indicates the x-ray unit is set for patients with small builds.

- **“Adult” selected:** Indicates the x-ray unit is set for patients with average-large builds.

To change the selection, press the relative button.

**NOTE:** After the change has been made, the preset exposure times will automatically be modified.

### - Checking the set intraoral examination

Examination of upper molars

Examination of upper premolars/canine or rear “bite wing”

Examination of upper incisors or front “bite wing”

Examination of lower incisors

Examination of lower premolars/canine

Examination of lower molars

### 3.3. Batteries and charge level status

The handheld runs on two widely available AA alkaline batteries to assure sufficient stand-alone operation. The charge level of the batteries is given on the screen as follows:

Battery fully charged (a symbol does not appear in the area that shows the battery charge level).

Battery half-charged.

Battery charge level low or almost dead (causing the handheld to automatically shut off).

**The batteries should be removed from the handheld if it is not going to be used for an extended period.**


### 3.4. Position of the patient

A positioner or alignment device specific for the selected image receiver should always be used to assure the x-rays are correctly aligned regardless of the position the patient's head is in.

### 3.5. Putting the x-ray unit cone into the required position

The x-ray unit cone is initially locked in place by an electric brake. To release the brake and turn the cone on the pivot point, use only the touch-sensitive unlock areas (see figure below).

Hold down the unlock buttons to tilt the cone the extent desired by the dentist to take the exposure and then simply release the unlock buttons to lock it in place.


**NOTE:** The cone should be firmly held in place with both hands while pressing the release buttons and turning the cone on the pivot point.

### 3.6. Position of the X-ray plate or sensor

The parallel technique, where applicable, provides more accurate images in terms of size compared to the bisecting technique. A rectangular collimator, with 30 cm focus-skin distance, is always preferable to obtain better quality pictures. To avoid exposing the image receiver only partly (whether it is a sensor or photostimulable phosphorus plate system) an alignment device that gives rectangular collimators guidelines should be used. These lines are usually given on the alignment ring.

#### Parallel technique


- The x ray emission axis is perpendicular to the image receiver (for example a sensor or photostimulated phosphor plate) which in turn is parallel with the tooth's long axis.
- As a result, the picture of the tooth will only be deformed by the divergence of the x rays in relation to the focus spot.
- Radiographic enlargement may reach up to 15%.
- For some "special" pictures, for example occluded ones, it may be necessary to remove the rectangular collimator and use the round one if a positioner is not present.

### 3.7. Checking the exposure time on the display

Before starting exposure, check the time setting on the handheld's screen (see the tables with the original exposure times, paragraph 3.1.4). To change the value, use keys "+" and "-".

**NOTE:**

Changing the exposure time is only temporary: if the new time is not stored in the memory it will be lost (see paragraph 3.8).

### 3.8. Procedure to be followed when taking the x-ray

- Pick up the handheld and go a safe distance away (at least 2 meters) maintaining visual contact with the patient and x-ray unit during the exposure. Make sure "ready" is indicated.

**READY**

- Tell the patient to stay still.
- Press and hold down the "Exposure" key on the handheld until the audible warning sound (beep) stops and the yellow light goes out.


"Exposure" key

Light on panel illuminated during exposure.

**NOTE:** If the "EMIT X-RAY" key is released at any time, exposure will be interrupted and error code E01 will appear on the display.

- Once exposure has been completed, it is possible to proceed with the next exposure unless the x-ray unit has reached the maximum allowable temperature. The percentage the cone exceeds the maximum allowable temperature is always shown on the screen (see icon below).

Once the maximum temperature is reached, it is necessary to wait a certain time to allow the instrument to cool down. This condition is indicated by this symbol.

- At this point the exposure function will be disabled until the screen shows "ready" again.

**READY**

- As soon as "READY" appears on the handheld, another exposure can be taken.

### 3.9. Setting the exposure time

The factory set exposure time limits cannot be changed in this x-ray unit.

**WARNING!**

Once this customisation has been made, the values indicated in the “Chart indicating original exposure times” no longer apply.

Proceed as follows to set the new exposure times:

- Set the exposure parameters according to the x-ray to be taken (see paragraph 3.2.).

- Press keys “INCREASE” or “DECREASE” to change the exposure time.

**NOTE:** When the desired button is held down, function “repeat” is activated and the displayed time scrolls faster.

- To make the new entry operative, make sure the symbol for button “MEMO” on the handheld’s screen is illuminated.

- **MEMO on:** indicates that the new setting can be stored in the memory.

Keep the “MEMO” button pressed for approx. 3 seconds. A BEEP will be given to indicate that the setting has been stored in the memory.

- **MEMO off:** indicates that the new setting cannot be stored in the memory.


**NOTE:** It is not possible to set the time when it exceeds the factory set exposure time limits.

### 4. Advanced options

The handheld allows the user to view, edit and set some operating parameters by simply combining the keys provided. Follow the steps given below to access:

Key combination	Description of command
	Press these two keys to adjust the sensitivity levels (set according to the table given below and type of sensor/receiver being used) and change the value from the minimum to maximum allowable ones (on a scale from 1 to 25) with keys “+” and “-”; to confirm the desired level and go back to the main screen, press key “adult”.
	Hold down these two keys to go to a setup menu (from P 01 to P 03). Press key “adult” to make the selection. Once within the individual configurations, they can be scrolled with keys “+” and “-” and selected by pressing key “adult” again. The configurations are given in detail below: P 01: Sets the stand by time (from a minimum of 5 to a maximum of 30 minutes). P 02: Assigns an identification tag to the x-ray unit’s base (from 1 to 5 or none, using “+” and “-”). P 03: Shows the list of software versions (from 1 to 3, using “+” and “-”).
	Hold down these two keys at the same time to set the rated current for x-ray exposure. The normal configuration, identified by “NORM” corresponding to 7mA, should always be used unless otherwise indicated by technical support.
	Hold down these two keys at the same time to immediately restore the “factory settings” (see paragraph 3.1.4) NOTE: This option can be accessed only if the advanced options menu has already been opened, i.e. after pressing the “tooth” and “adult” keys at the same time.

- Sensitivity levels

- Key:**  
 I = Incisors  
 CP = Canines / Premolars  
 M = Molars  
 Bp = Rear bite wing posteriore  
 Ba = Front bite wing

**NOTE:** For further information regarding any fault messages that may be displayed, see the table in paragraph 9.2.

## 5. Routine maintenance

### ATTENTION!

**Any technical maintenance work required must be carried out by qualified personnel or by a specialised technician authorised by the manufacturer.**

It is the user's responsibility to check that routine maintenance is carried out by an authorised technician at least every 2 years. The maintenance methods are specified in the Technical Service Manual possessed by the Authorised Technicians.

## 6. Cleaning and disinfection

The x-ray unit may be a source of infection passed from one patient to another if it is not cleaned properly. For this reason it should be disinfected on the outside every day after use.

If digital X-ray sensors are used make sure they are always used with disposable hygienic covers. When disinfecting the x-ray unit use soft disposable paper and avoid using corrosive substances or immersing the unit in liquids.

To avoid damaging the plastic materials use products containing:

- **Ethanol 96%**

Concentration: maximum 30 g per 100 g of disinfectant.

- **Propanole**


Concentration: maximum 20 g per 100 g of disinfectant.

- **Combination of ethanol and propanole**

Concentration: the combination of the two should be maximum 40 g per 100 g of disinfectant.

Compatibility tests between plastics and the following products have been carried out with no negative consequences:

- Incidin Spezial (Henkel Ecolab).
- FD 312 (DURR)
- Omnidid (Omnident)
- Plastisept (ALPRO) (not tuberculocide as not an alcohol-based disinfectant)
- RelyOn Virkosept (DuPont).

### ATTENTION!

- **Do not use products containing isopropyl alcohol (2-propanol, iso-propanol).**
- **Do not use products that contain sodium hypochlorite (bleach).**
- **Do not use cleaners that contain phenol.**
- **Do not spray the selected products directly on the surfaces.**
- **Never combine products with each other or with liquids other than the products listed above.**
- **All products must be used as directed by the manufacturer.**

### Cleaning and disinfecting instructions.

Use soft disposable paper towels or sterile gauze to clean and disinfect.

Do not use sponges or in any case, any material that can be reused.

### ATTENTION!

- **Shut off the unit prior to clean and disinfecting the external parts.**
- **Never lubricate the pivot point of the x-ray cone as proper operation of the locking system may be compromised.**
- **All material used to clean and disinfect must be thrown away.**

## 7. Dismantling the equipment when no longer used

As set out in Directives 2002/95/ EC, 2002/96/ EC and 2003/108/ EC, on the restrictions of the use of certain hazardous substances in electrical and electronic equipment along with collection, treatment, recycling and disposal of waste electrical and electronic equipment the latter must be treated as municipal waste, therefore sorted and collected separately. When new equipment of equivalent type is purchased the waste equipment should be returned to the distributor on a one-to-one basis for disposal.

As far as reuse, recycling and other forms of waste recovery mentioned above are concerned, the manufacturer is responsible for the actions specified by individual local laws.

Efficient collection of sorted waste separately to recycle and treat waste electrical and electronic equipment aids in preventing negative environmental impacts while protecting human health. In addition it facilitates recycling of the materials used to construct the equipment.

The crossed out wheeled bin placed on the equipment indicates that the waste equipment must be collected separately from other waste.

### ATTENTION!

**Illegal waste disposal carries heavy fines defined by local laws.**

## 8. Specifications

- Rated voltage: 230VAC/115VAC (according to the model).
- Max. mains voltage fluctuation:  $\pm 15\%$ .
- Rated current: 6A for 230 VAC model; 10A for 115VAC model.
- Frequency: 50 / 60 Hz.
- Maximum absorbed power: 1,4kVA.
- Apparent line resistance: 0,5 $\Omega$ .
- Fuses for electronic parts: 6,3A-230VAC; 10A-115VAC
- Generator: constant potential type.
- Rated high voltage: 60kV.
- Rated current: 3,5mA - 7mA.
- Rated power at 0.1 sec: 420W (60kV 7mA), 210W (60kV 3,5mA).
- Reference time current product: 0,7 mAs (7mA - 0,1s) - 0,35 mAs (3,5mA - 0,1s).
- Total filtration: 2,5mm.
- Hemivalent layer (HLV) at 60kV: >2mm Al.
- Leaked radiation: <0,25mGy / h at 1 metre from focusing point.


- Ability to be reproduced: 0,05.
- Electrical classification: Class I - Type B, intermittent service.
- Set exposure time:  $\pm 0,010$  a 1,000 se.
- Accuracy of times indicated:  $\pm 10\%$  (maximum).

#### WEIGHTS

- Weight of the unit with packaging: 25 Kg.
- Weight of the hand-held control panel: 0,3 Kg.
- Weight of the x-ray unit: 23 Kg (with 90 cm extension).

#### RADIOGENIC TUBE

- Radiogenic tube: TOSHIBA D-041.
- Focal spot: 0,4 mm in compliance with IEC 336 / 1993.
- Tolerance for position of the focal spot along the reference axis:  $\pm 2\%$ .
- Rated high voltage and maximum current: (60kV 7mA)  $\pm 10\%$ .
- Anode construction material: Tungsten (W).
- Anode inclination:  $12,5^\circ$  (see figure to the side).
- Anode thermal load: 4,3 KJ (6 KUT).
- Maximum continuous heat dissipation: 7 W.
- Operating cycle: 1:60 (1 second exposure - 60 seconds pause time).


#### COLLIMATOR SPECIFICATIONS

- Long rectangular collimator: SSD = 30 cm, (12"), x-ray patterns less than or equal to 45 x 35 mm,
- Short round collimator: SSD = 20 cm, (8"), x-ray patterns less than or equal to 60 mm diameter.

#### BATTERIES

Type: AA alkaline.




A - REFERENCE AXIS

B - FOCUSING POINT IDENTIFICATION

## 9. Specifications

### 9.1. Dimensions

## 9.2. Identification nameplates

### WARNING!

Never remove the identification nameplates provided on the generator, central control unit and collimator cone.



#### Central control unit (NAMEPLATE)

The nameplate is found next to the main switch.

NOTE: The voltage given on this nameplate will vary depending on the model of x-ray unit used (230 Vac or 115 Vac).

Details shown on the nameplate:

- Name of the manufacturer
- Name of the equipment
- Rated voltage
- Type of current
- Rated frequency
- Maximum absorbed power
- Serial number
- Year of manufacture


#### X-ray unit head

The nameplate is found on the lower cover at the back of the radiogenic unit.

Details shown on the nameplate:

- Name of the manufacturer
- Name of the equipment
- Specifications
- Serial number
- Date of manufacture.



### Collimator

The nameplate for the rectangular collimator is found outside it.

Details shown on the nameplate:

- Name of the manufacturer
- Type of cone
- Serial number
- Date of manufacture.


### Handheld

The nameplate for the handheld is found in the battery compartment.

Details shown on the nameplate:

- Name of the manufacturer
- Name of the equipment
- Rated voltage
- Number and type of batteries.
- Serial number.

### 9.3 Fault table

FAULT	CAUSE	SOLUTION
<b>E01</b>	EXPOSURE KEY RELEASED TOO SOON.	HOLD DOWN THE KEY UNTIL THE EXPOSURE CYCLE HAS FINISHED.
<b>E02</b>	COLLIMATOR SELECTION INCONSISTENT.	AFTER TURNING THE RECTANGULAR COLLIMATOR ON OR OFF WAIT A FEW SECONDS FOR THE ICON ON THE HANDHELD TO BE REFRESHED.

**FOR OTHER FAULT CODES, CONTACT CUSTOMER SERVICE.**

DICHIARAZIONE “CE” DI CONFORMITÀ  
“CE” CONFORMITY DECLARATION  
DECLARATION DE CONFORMITÉ “CE”  
ERKLÄRUNG VON “CE” ZUSTIMMUNG  
DECLARACION DE CONFORMIDAD “CE”  
DECLARAÇÃO DE CONFORMIDADE “CE”  
ΔΗΛΩΣΗ ΠΙΣΤΟΤΗΤΑΣ “CE”  
ДЕКЛАРАЦИЯ СООТВЕТСТВИЯ “CE”  
CERTYFIKATU CE

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