

# Cochlear™

## Nucleus® CP900 Series Sound Processors

User Guide



*Hear now. And always*



**Cochlear™**

# Symbols



**Note**

Important information or advice. Can save inconvenience.

---



**Tip**

Time saving hint or saves inconvenience.

---



**Caution (no harm)**

Special care to be taken to ensure safety and effectiveness.

Could cause damage to equipment.

---



**Warning (harmful)**

Potential safety hazards and serious adverse reactions.

Could cause harm to person.

---

# Contents

|                                     |    |
|-------------------------------------|----|
| About .....                         | 4  |
| Turn on.....                        | 6  |
| Wear.....                           | 10 |
| Store.....                          | 31 |
| Charge.....                         | 32 |
| Care .....                          | 39 |
| Lights, beeps and buttons .....     | 54 |
| General warnings and cautions ..... | 58 |
| Other information.....              | 60 |

## About

The Cochlear Nucleus CP900 Series Sound Processors works with your Cochlear Implant to transfer sound to the ear. Your processor is made up of a processing unit, earhook, coil, coil cable and a battery module.

You can choose either the CP910 or the CP920 processing unit. The CP920 processing unit is smaller than the CP910 and has no accessory port.


You can power your processing unit with a range of battery modules. The range includes a tamper resistant battery module that uses disposable batteries, a standard rechargeable and compact rechargeable battery module, and a LiteWear cable that allows the battery to be worn on the body.

You can control your sound processor as a stand-alone device by pressing the buttons on the processor. All stand-alone functions are covered in this guide.

You can also choose a Cochlear Nucleus CR210 Remote Control or CR230 Remote Assistant to control your processor. The remotes let you control one or two bilateral processors at a time and can provide extra troubleshooting functions. All remote functions are covered in separate user guides.

People with certain types of hearing loss can wear the processor in Hybrid mode by adding an accessory called the acoustic component. The acoustic component sends extra acoustic sound into the ear canal. Hybrid mode must be set up by your clinician.

A number of tools and accessories are supplied with your sound processor. The most important of these are your drying pack, screwdriver, and microphone protector applicator with removal tool.



**Warning**

Do not allow children under three to play with the system. Removable parts of the system can be lost or may be a choking hazard.



Behind the Ear (BTE) parts (standard)



BTE parts (mini)

## Turn On




### Turn on

|                               |   |
|-------------------------------|---|
| Turn on and off .....         | 7 |
| Change programs .....         | 8 |
| Lock and unlock buttons ..... | 9 |

## Turn on and off



Press and hold the lower button for two seconds to turn on or off.

| Indicator light  | What it means  |
|--|--|
| <br>Green flashes                                 | Turning on processor. The number of flashes indicates the number of the current program. |
| <br>Steady orange while pressing the lower button | Turning off processor.   |
| <br>Quick green flashes                           | Microphones are receiving sound  |



### Note

Your clinician may set up your processor to turn off automatically when it has not been used or worn for more than two minutes. If this happens, simply turn your processor back on.


**Turn On**

## Change programs

You can choose between programs to change the way your sound processor deals with sound e.g. in windy, noisy or quiet places. Usually two programs are all you need, but your clinician can give you up to four programs.



Quick-press and release the lower button to switch between programs.

| Indicator light  | What it means   |
|--|---|
| <br>Green flashes | Changing the program. The number of flashes indicates the number of the selected program. |




Ask your clinician to write your programs below

P1 \_\_\_\_\_      P2 \_\_\_\_\_  
P3 \_\_\_\_\_      P4 \_\_\_\_\_

## Lock and unlock buttons



Quick-press and release both buttons at the same time to lock and unlock the processor.

| Indicator lights  | What it means                 |
|---|-------------------------------|
|  | Locking processor buttons.    |
|  | Unlocking processor buttons.  |
|  | Processor buttons are locked. |

# Wear

- Attach to your head .....11
- Talk on the phone..... 12
- Water, sand and dirt ..... 14
- Sport and exercise..... 16
- Wear LiteWear .....22
- Travel.....23
- Audio accessories .....24
- Turn audio accessories on and off .....26
- Personalise your processor .....27
- Wear an acoustic component.....28

## Attach to your head




① Place the BTE on your ear, letting the coil dangle.




② Move the coil sideways and onto your implant.

| Indicator light                   | What it means  |
|-----------------------------------|--|
| ●●●● Flash of orange every second | Coil is off, or processor is not connected to the correct implant. |

 **Tip**

Twist the coil clockwise (right BTE) or counter clockwise (left BTE) once **ONLY** before placing it on the head. The twist will help the coil stay in place.



 **Wear**




## Talk on the phone

Listen to hearing loops and hear voices on the telephone more clearly with the telecoil (t-switch). Your processor can be set up with either manual or auto telecoil. Auto telecoil turns the telecoil on and off whenever you use a phone or move in and out of range of a hearing loop.

Whether you choose auto or manual telecoil, you can control the telecoil with your remote or just turn it on and off with your processor.



Quick-press the top button to turn telecoil on and off.

| Indicator light   | What it means   |
|---|---|
|  Long flash of blue          | Changing between microphones and the telecoil or audio accessory.     |
|  Long flash of green         | Changing between the telecoil or audio accessory and the microphones. |
|  Constant quick blue flashes | Processor is receiving sound  |


 **Tip**

Slowly move the phone's sound outlet over the microphones on your processor to find your best position for talking on the phone.

Phone being held over microphones

 **Note**

Certain electronic devices or machinery may automatically turn on Auto Telecoil. If this happens, move away from the device or machinery and wait a few seconds for Auto Telecoil to turn itself off, or turn it off manually.

 **Wear**

### Water, sand and dirt

Your sound processor is resistant to water, sand and dust. However, it is still a precision electronic device so the following precautions should be taken.



Your processor is water resistant with rechargeable batteries attached. But do not wear it while bathing or swimming.



If your processor ever does get wet, dry it with a soft cloth, change the microphone protectors and place it in a drying kit for 12 hours.



Your processor is dust resistant with rechargeable batteries attached. But do not get sand or dirt into any part of your sound processor.



If sand or dirt ever do enter the system, shake the components carefully to remove it and change the microphone protectors.

Fully assembled as described

Your processor is protected against failure from dust and temporary immersion in water (IP57) when you wear it with:

- a rechargeable battery module
- a coil and coil cable
- a closed accessory port
- no acoustic component


Fully assembled as described

Your processor is protected against failure from splashing water or access of foreign objects 1.0 mm in diameter or larger (IP44) when you wear it with:


- a tamper resistant battery cover
- a coil and coil cable
- a closed accessory port
- an acoustic component worn in the ear.

Aqua Accessory TBA



 **Wear**

### Sport and exercise



Use Cochlear CP900 series accessories such as the Snugfit or Mic Lock to help hold your processor in place when you play sport or exercise.



After exercise, wipe your processor with a soft cloth to remove sweat or grime.



Then check your microphone protectors.




Change them whenever you notice signs of sweat or grime, or any change in performance.




#### Childproofing Tip


Accessories like the Snugfit or Mic Lock can help hold a processor on a child's ear. Children under three years of age should only use tamper resistant accessories such as the tamper resistant earhooks, tamper resistant battery cover and tamper resistant Snugfit. Always ensure these are properly locked together before use.



A thin sweatband or skullcap worn over the sound processor can make wearing a helmet easier.



Rechargeable batteries give you the maximum water and sweat resistance.




## Wear

### Wear a Snugfit


The Cochlear Nucleus CP900 Series Snugfit holds your sound processor more securely in place than an earhook alone. It gives parents and carers easy access to controls and accessories and is available in large, medium and tamper resistant (small).



To attach a Snugfit:



- ① Press firmly on the top of the earhook with your thumb to click it off the sound processor.




- ② Click your Snugfit into place.



- ③ Carefully bend the lower part into shape to suit your ear.

To remove a Snugfit:



- ① Remove your Snugfit by clicking it off the processor.



**Childproofing Tip**

To attach the Tamper Resistant Snugfit:



- ① The smallest Snugfit comes with a Snugfit band to fix it more securely to the processor.



- ② Thread the Snugfit band onto the Snugfit as shown.



- ③ Click the Snugfit into place on the processor and push the Snugfit band onto the battery module.



- ④ Push the band securely onto the battery module.

Wear a Mic Lock


The Cochlear Nucleus Mic Lock helps hold the processor securely on the ear.



To attach the Mic Lock the first time:



- ① Push the Mic Lock band onto the processor with the tubing facing forward.
- ② Place the processor on your ear and bring the tubing up to the earhook.



- ③ Trim the tubing to a length that will hold the processor firmly in place without discomfort. Do not cut it too short.



- ④ Soften the end of the tubing in warm water and push it onto your earhook.
- ⑤ Let the tubing cool before wearing.

## Wear LiteWear

The Cochlear Nucleus CP900 Series LiteWear Cable allows any standard (not compact) battery module to be worn on the body - so only the processing unit is worn on the ear. You can order LiteWear at <TBA> or talk to your clinician.

IMAGE SHOWING  
LITEWEAR BEING WORN  
BY A CHILD

LiteWear is made up of a LiteWear cable, case and a choice of fixing aids.



### Note

You can wear a Snugfit with LiteWear to hold the processing unit more securely on the ear.



### Childproofing Tip

Attach the LiteWear case between the shoulder blades of very young children to help prevent tampering.

## Travel

You can travel with your sound processor with confidence because Cochlear is supported by clinics all over the world.

- Take a printout of your most recent program in case you need help with your processor.
- If you have a back-up sound processor, check that it is programmed correctly and take it with you.
- It's okay to move through metal detectors and full body scanners with your sound processor on. To avoid any possible buzzing sounds in your ear, turn off the telecoil.
- Ask your clinician for a Nucleus 6 ID card. In the unlikely event that your implant sets off a metal detector the ID card will help you explain that you have an implanted medical device.
- Do not place sound processors directly onto conveyor belts as you move through airports. Place them in a bag or plastic tray.
- Your sound processor will not interfere with a plane's navigation system so you won't need to turn it off during takeoff and landing.



### Note

Visit [www.cochlear.com/clinic-finder](http://www.cochlear.com/clinic-finder) to find the nearest clinic in the places you are travelling

## Wear


### Audio accessories




The accessory port on the CP910 processing unit lets you connect a wide range of audio accessories.



For example, a Cochlear Nucleus personal audio cable let you listen to iphones, computers, TV sets etc. Bilateral recipients can listen in stereo.



The Cochlear Nucleus Lapel Microphone can make conversation easier in noisy environments.



Parents and carers can use the Cochlear Nucleus Monitor Earphones to check the processor and audio accessories are working.

You can order audio accessories at <TBA> or talk to your clinician.

### Note

The CP910 processing unit is compatible with the CP800 range of audio accessories. You can order adaptors to make your CP910 processor compatible with Freedom and Euro accessories. Your CP910 processor is NOT compatible with the BAHA range of audio accessories.

### Connect audio accessories

To connect an audio accessory to the processor:



- ① Carefully lift open the accessory cover with the screwdriver. Do not twist the cover.



- ② Hold the side grips on the accessory connector and carefully push it into the accessory port until it clicks into place.



- ③ When you connect an audio accessory, your processor automatically detects it.
- ④ If applicable, insert the other end of the audio accessory cable into the compatible device e.g. Ipod.

### Warning

Always use the Mains Isolation Cable when connecting any Personal Audio Cable (including Bilateral) to:

- a mains powered device e.g. TV.
- a battery powered device that is connected to mains power (e.g. a laptop connected to an electric outlet for charging).

Do not use excessive force, twist or pull on accessories when connecting them to the processor or at any time while they are connected.




**Wear**

### Turn audio accessories on and off

Audio accessories turn on automatically when you plug them into your processor. You can turn them off and back on again as follows.

Quick-press the top button to switch between microphones and the audio accessory.



| Indicator light   | What it means   |
|---|---|
|  Long flash of blue           | Changing between microphones and the telecoil or audio accessory.     |
|  Long flash of green         | Changing between the telecoil or audio accessory and the microphones. |
|  Constant quick blue flashes | Processor is receiving sound from audio accessory                     |

### Personalise your processor

You can order Cochlear Nucleus Sound Processor Covers and Coil Covers to personalise your processor and protect it from scratches.

Attaching and removing processor covers:



- ① To attach, slide the cover onto the processor.
- ② To remove, pull the cover away from the processor.

Attaching and removing coil covers:



- ① To attach, discard the clear plastic protector on the cover and click the cover into place on the coil.
- ② To remove, lift the cover off the coil.

### Bilateral stickers

If you have two Cochlear implants, you must match the correct processor to the correct implant. Ask your clinician to give you bilateral stickers (red for right, blue for left) to make identifying left and right processors easier.

## Wear

### Wear an acoustic component

If your sound processor is set up for use in Hybrid mode, you can wear it with an acoustic component.



- ① Hold the coil in one hand and place the processor on your ear with the other.
- ② Place the coil on your head.
- ③ Place your dome or earmould in your ear as instructed by your clinician.

### Remove your acoustic component

Wear your processor without the acoustic component by replacing it with a standard earhook:




- ① Remove the battery module from the processing unit by twisting the two parts as shown.



- ② Insert your Cochlear screwdriver into the opening under the earhook to unlatch the receiver cable.



- ③ Carefully push and twist the screwdriver until the receiver cable disconnects.



- ④ Re-attach the battery module.



- ⑤ Press firmly on top of the earhook with your thumb to click it off the sound processor.



- ⑥ Click the standard earhook into place.

### Caution

Use only the Hybrid screwdriver – not the CP800 screwdriver – to remove the acoustic component.

## Wear

### Attach your acoustic component




- ① Press firmly on the top of the standard earhook with your thumb to click it off the processor.



- ② Click the Hybrid earhook into place.



- ③ Carefully slide the receiver cable plug into the sound processor.



- ④ It should click securely into place. If necessary, carefully twist it until you feel the click.

## Store



Store your processor at night in a UV protected dry aid kit or drying pack. You can store the processor fully assembled for short periods.



For long term storage, remove the batteries and store them separately from each other to prevent shorting. Storage cases are available from Cochlear.



Storage cases are available from Cochlear. Buy dry aid kits and drying packs from Cochlear or pharmacies and read the instructions for use.

### Warning

All drying chemicals must be kept away from young children. Swallowing this material can cause serious internal injuries.



## Charge

|   |    |
|---|----|
| Batteries .....   | 33 |
| Battery life expectancy .....   | 33 |
| Replace the battery.....  | 34 |
| Lock the battery to the processing unit .....   | 35 |
| Lock and unlock the tamper resistant battery cover.....   | 36 |
| The standard battery has a tamper resistant lock to help prevent children opening the battery cover. .... | 36 |
| Charge rechargeable battery modules.....  | 37 |
| Change disposable batteries .....   | 38 |

## Batteries

You have a choice of three battery types with the CP900 series system:



The Cochlear Nucleus CP900 Series Standard Rechargeable Battery module.



The Cochlear Nucleus CP900 Series Compact Rechargeable Battery module.



The Cochlear Nucleus CP900 Series Standard Battery module, which uses disposable batteries and has a tamper resistant lock.


## Battery life expectancy

Battery life varies according to the programs used each day, the thickness of skin covering your implant, and the type of battery. Use 675 zinc/air batteries in the standard battery module for longest battery life. To get the longest life from the rechargeable batteries, always recharge before use. Rechargeable batteries may take longer to fully recharge with age and should be replaced as needed.

## Charge


### Replace the battery

To remove the battery:



- ① Twist the battery module as shown to remove it from the processing unit.

To attach the battery:



- ① Touch the battery module and processing unit at an angle so their parts fit together.



- ② Twist the battery module as shown to attach the two parts.



## Childproofing Tip

### Lock the battery to the processing unit

#### Caution

Always check this lock is unlocked before attaching or removing the battery module.

### Lock the CP910 Processing Unit:



- ① Carefully lift open the accessory port cover with your Cochlear screwdriver. Do not twist the cover.



- ② Push the lock to the far left to lock the battery module.



- ③ Push the lock to the far right with your Cochlear screwdriver to unlock the battery module.




- ④ Carefully close the cover by pushing in the top hinge first.


**Charge**

**Childproofing Tip**

Lock the CP920 processing unit



- ① Push the lock to the far left with your Cochlear screwdriver to lock the battery module.



- ② Push the lock to the far right to unlock the battery module.

Lock and unlock the tamper resistant battery cover

The standard battery has a tamper resistant lock to help prevent children opening the battery cover.



- ① To lock, turn the locking screw clockwise with the CP800 screwdriver until it is in a horizontal position.
- ② To unlock, turn the locking screw anticlockwise until it is in a vertical position.

Charge rechargeable battery modules



- ① Angle and fit the rechargeable battery into a socket on the Cochlear Nucleus Battery Charger.







- ② Twist clockwise to connect.

This image will show lights.



- ③ Connect the battery charger to the power adaptor, plug it into mains power and switch on.

| Battery wall power light  | What it means   |
|---|---|
|  Steady green    | Battery charger power is on   |
| Does not light up   | Power is not plugged in or turned on at the power socket, or there is no power. |
| Battery module light  | What it means   |
|  Steady orange   | Battery module is recharging.   |
|  Steady green    | Battery module is fully recharged.  |
|  Flashing orange | There is a problem with the battery module.                                     |
| Does not light up   | Battery module is not properly placed, or there is no power.                    |

## Charge

### Change disposable batteries

The Cochlear Nucleus CP900 Series Tamper Resistant Battery Cover uses two high power zinc/air disposable batteries. Cochlear recommends Power One p675 Implant plus zinc/air batteries. Do not use silver oxide or alkaline batteries.



- ① Pull the battery cover away from the battery holder.




- ② Remove the batteries from the battery holder.



- ③ Remove the new batteries from the packet and let them stand for a few seconds.

- ④ Insert the batteries into the holder with the flat side (positive terminal) facing up.



- ⑤ Replace the battery cover by sliding it up towards the processing unit.

## Care

|  |    |
|--|----|
| Care for your sound processor .....    | 40 |
| Replace microphone protectors .....    | 41 |
| Care for your battery charger .....    | 43 |
| Change earhooks .....                  | 44 |
| Change the coil cable .....            | 48 |
| Change the coil magnet .....           | 49 |
| Care for your acoustic component ..... | 50 |

## Care for your sound processor

### Every day

- ❑ Check all parts for dirt and moisture. Wipe the processor, coil cable, earhook, covers, Snugfits and Mic Lock cables with a soft dry cloth. (You can leave the earhooks, Snugfits, Mic Locks or coil cables on the processor while cleaning.)
- ❑ Keep your processor free from moisture by drying it every night in your drying pack.
- ❑ Remove the battery and make sure all the contacts are clean. Carefully tap or blow on them to remove any dirt.
- ❑ Make sure the accessory port and audio accessories are clean. Blow on the accessory port to remove any dust and clean accessories with a soft dry cloth.
- ❑ Check the microphone protectors for signs of dirt or grime and replace if needed.

### Every month

- ❑ Replace the dry brick in your drying pack.
- ❑ Check if earhooks or Snugfits are becoming loose or showing signs of wear. Replace as needed.

### Every three months

- ❑ Replace the microphone protectors.


#### **Note**

If your sound processor is subject to extreme conditions, you should clean and replace parts more often to keep it working at its best.

## Replace microphone protectors

Replace your microphone protectors every three months, whenever they look dirty or whenever you notice any loss in sound quality. Always replace both protectors at the same time.

### To remove microphone protectors:



- ① Remove the protective cap from the Cochlear Nucleus Microphone Protector Removal Tool.




- ② Push the tip of the removal tool under the edge of the microphone protector. Do not push the tool into the middle of the protector as you might damage the microphone.
- ③ Carefully lift out the microphone protector and discard.

#### **Warning**

The removal tool is a sharp object. It is not recommended for use by children under the age of 12 years.

## Care

To insert new microphone protectors:



- ① Pull the Cochlear Nucleus Microphone Protector Applicator out of its sleeve.



- ② The top end of the applicator is marked with an arrow. Hook the top end onto the earhook.


- ③ Hook the bottom end onto the coil cable.



- ④ Align the microphone protectors with the microphones.



- ⑤ Press the microphone protectors into the microphone holes one at a time.



- ⑥ Remove the applicator by carefully peeling it upwards from the coil cable end.

## Care for your battery charger

Every day:

- ❑ Check your battery charger is clean. If you notice any dust or dirt:
  1. Unplug the power adaptor and remove any battery modules.
  2. Hold the charger upside down and tap it to remove any dirt from the battery charger sockets. Carefully blowing on the battery charger sockets may also help remove dirt.
  3. Wipe the battery charger sockets with a soft dry cloth.
- ❑ Cycle batteries in different sockets when recharging for a more even wear of the sockets.

If it gets wet:

- ❑ If the battery charger ever gets splashed with liquid, carefully shake out the liquid and dry the charger for 24 hours. Do not use the charger until it is dry.

### Change earhooks

Cochlear Nucleus CP900 Series Standard Earhooks come in large, medium and small sizes. The smallest earhook can be made tamper resistant and is recommended for children under three. The acoustic component is worn with a Cochlear Hybrid Earhook. (To change the hybrid earhook, see "Wearing/Removing your acoustic component".)

To change a standard earhook:



① Press firmly on the top of the earhook with your thumb and click it off the sound processor.

② Click the standard earhook in place.



**Note**  
Remove your earhook only when necessary - it may become loose if it is removed too often.

### About the tamper resistant earhook:

The tamper resistant earhook can be locked to the processor with a pin to reduce the risk of the earhook becoming a choking hazard for young children. Older children and adults can wear it without the pin.



The tamper resistant earhook comes with a locking tool and fitting pins.

**Warning**  
The locking tool and pin is not recommended for use by children under 12. Do not let children attach or remove the tamper resistant earhook without adult supervision.

 **Childproofing Tip**

Attach the tamper resistant earhook



① Click the tamper resistant earhook onto the processor. Make sure the hole on the earhook is aligned with the hole on the processing unit.



② Insert a fitting pin into the hole on the blunt end of the locking tool.



③ Partly insert the fitting pin into the hole on the earhook.



④ Flip the tool and use the pointed end to push the pin fully into the earhook. The pin should not protrude from either side of the earhook.

 **Childproofing Tip**

Remove the tamper resistant earhook



① Hold the processor as shown and insert the pointed end of the locking tool into the hole on the earhook.



② Push the tool into the hole to partly remove the fitting pin. Be gentle as too much force may damage the tool.



③ Lift out the earhook fitting pin with tweezers.







④ Press firmly on the top of the earhook with your thumb and click it off the sound processor



## Change the coil cable



The Cochlear Nucleus CP900 Series Coil Cable come in four lengths so you can choose the most comfortable fit. Only remove the coil cable when you are changing it.

-  ① Hold the cable grip as shown and firmly pull it straight out of the processor. Do not tug on the flexible part of the cable.
-  ② Use your fingernails to hold the cable grip as shown and pull the cable firmly off the coil. Do not tug on the flexible part of the cable.
-  ③ Push the new coil cable into the processor until it clicks into place.
-  ④ Push the new coil cable into the coil until it clicks into place.



## Change the coil magnet

The strength of the Cochlear Nucleus CP900 Series Coil Magnet can be adjusted by changing the magnet depth or by replacing with a magnet of a different strength. If the magnet is too weak or strong the coil may fall off or cause discomfort.

To change the strength of the magnet:

-  ① With finger grips facing up, turn the magnet clockwise to increase the magnet's strength.
-  ② Turn the magnet anticlockwise to reduce its strength.

To change the coil magnet:

-  ① With finger grips facing up, unscrew the magnet anticlockwise and remove.
-  ② Insert the new magnet into the coil hole and turn it clockwise until it is at a comfortable depth.

## Care for your acoustic component

The acoustic component accessory is used only by people who are set up to use their processor in Hybrid mode.

### Caution

Handle your sound processor and acoustic component with care above flat surfaces. Never drop them or twist the receiver cable.

### If you use Power or Plus domes:

Every day:

- Check the wax filter in the receiver and change it whenever it is clogged with wax or the system does not sound normal.
- Domes have an additional inbuilt filter. Check the inbuilt filter and change the dome whenever it is clogged with wax, or whenever the dome is dirty, discoloured or losing shape.

### If you use Power or Micro moulds:


Every day:

- Clean the earmould and receiver cable with the supplied brush, a soft dry cloth or our recommended cleaning fluids, Hadeo™ and Cedis® ONLY. Wipe any traces of wax away from the opening of the vent.
- Check if the vent is clogged, and use the supplied vent cleaner to remove wax as needed. Do NOT use the vent cleaner to clean any other part of the system.
- Check the wax filter and change whenever it is clogged with wax or the system does not sound normal.

### Note

Micro moulds can use WaxStop filters ONLY. Power moulds can use either NoWax or WaxStop filters.

## Replace Power domes



- ① Pull the dome off using your fingernails if necessary and discard.
- ② Push the new dome directly onto the receiver.
- ③ Click it securely into place.

## Replace Plus domes

Follow the directions for replacing Power domes and add the following step:



- ④ When fitting a Plus dome, the outer sleeve should face down.

## Care


### Change wax filters

#### Caution


Wax filters must ALWAYS be used with the acoustic component and replaced regularly.

#### Replace a NoWax filter


Power moulds (only) can use a NoWax filter in the sound outlet. Change the filter whenever it is clogged with wax or whenever the system does not sound normal.



- ① Remove a new NoWax tool from its shell. One point on the tool has a removal tool and the other point has the new filter.



- ② Insert the removal tool into the existing wax filter and work it out of the Power mould.




- ③ Push the new filter into the sound outlet on the mould and dispose of the tool and used filter.

#### Replace a WaxStop filter


Micro moulds use a WaxStop filter in the sound outlet. Power moulds use either a WaxStop filter or a NoWax filter. Power and Plus domes must be worn with a WaxStop filter in the receiver.




- ① Remove a WaxStop tool from its shell. One end of the tool has a new filter and the other has a removal tool.



- ② Insert the removal tool into the existing wax filter.



- ③ Pull the used filter out slowly, keeping the tool straight.



- ④ Flip the tool and push the new filter into the sound outlet. Carefully pull the tool free and dispose of the used tool and filter.












## Lights, beeps and buttons

### Lights, beeps and buttons

|                                       |    |
|---------------------------------------|----|
| Lights.....                           | 55 |
| Beeps.....                            | 56 |
| Button controls in advanced mode..... | 57 |

## Lights

Your clinician can set up your processor to show some or all of the following light indications.

| Light  | What it means  |
|--|--|
| <br>Short green flashes                         | Turning on processor. The number of flashes indicates the number of the selected program.                              |
| <br>Steady orange while lower button is pressed | Turning off processor.   |
| <br>Short green flashes                         | Changing the program. The number of flashes indicates the number of the selected program e.g. 2 flashes for program 2. |
| <br>Single green flash                          | Changing volume or sensitivity level.  |
| <br>Green followed by orange flash              | Locking processor buttons.   |
| <br>Orange followed by green flash              | Unlocking processor buttons.   |
| <br>Orange flash when pressing buttons          | Processor buttons are locked.  |
| <br>Long flash of blue                          | Changing between microphones to the telecoil or audio accessory  |
| <br>Long flash of green                       | Changing between the telecoil or audio accessory to the microphones.   |
| <br>Constant quick green flashes              | Processor is receiving sound from the microphones  |
| <br>Constant quick blue flashes               | Processor is receiving sound from telecoil or audio accessory  |

## Lights, beeps and buttons

|                                   |  |
|-----------------------------------|--|
| ●●●● Flash of orange every second | Coil is off, or processor is not connected to the correct implant.   |
| ●●●●● Quick flashes of orange     | Processor battery is empty.  |
| ▬ Steady orange                   | General fault: <ul style="list-style-type: none"> <li>• Fault with processor</li> <li>• Corrupt program</li> <li>• Fault and coil off error at the same time</li> <li>• Fault with acoustic component while in Hybrid mode.</li> </ul> |

## Beeps

Your clinician can set up your processor so you can hear the following beeps.

| Beep                                     | What it means   |
|--|---|
| 🔔🔔🔔 Short high beeps                     | Changing the program. The number of beeps indicates the number of the selected program. |
| 🔔 High beep                              | Changing volume or sensitivity level.   |
| 🔔... Long high beep                      | Switching between using the microphones and the telecoil or audio accessory.            |
| 🔔...🔔...🔔 Long low beeps                 | General fault.  |
| 🔔 Short low beep when you press a button | Processor buttons are locked.   |
| 🔔🔔 Short low beeps                       | Processor battery is low.   |
| 🔔🔔🔔 Short low beeps                      | Battery is empty and processor is turning off.  |

## Button controls in advanced mode

Sound processors can be set up in Advanced mode by your clinician to give you extra control over the volume and sensitivity of the sound you hear.

### Turn processor on and off



Press and hold upper and lower buttons at the same time. Pressing and holding the lower button will only turn on the processor.

### Change program



Press and hold the lower button.

### Turn telecoil/audio accessory on and off.



Press and hold the upper button to switch between microphones and telecoil/audio accessory.

## Lights, beeps and buttons

### Lock/unlock processor buttons




Quick-press both the upper and lower buttons at the same time.

### Change volume or sensitivity

In Advanced mode, you can choose to control volume and sensitivity with your remote and the buttons on your processor. See the user guides for the CR210 and CR230 remotes for more information.

Select "Volume" via your remote and use the buttons below to control the loudness of sound.

Select "Sensitivity" via your remote and use the buttons below to control the *range* of sounds picked up by the processor e.g. soft, background, close or far away.



Press the upper button to increase volume or sensitivity.

Press the lower button to decrease volume or sensitivity.

## General warnings and cautions

### Processor and parts

- Do not push the master volume too high for comfort in case a loud noise occurs nearby.
- Do not place your processor or parts in nose or mouth or other body openings.
- Do not use your processor or parts in a way that will restrict your airways (e.g. when using LiteWear or long coil cable).
- Do not let your processor or parts entangle with any jewellery (e.g. earhook and earrings) or machinery.
- Do not apply continued pressure to the coil when in contact with the skin (e.g. sleeping while lying on coil, or using tight-fitting headwear).
- Remove the processor and coil immediately if there is any discomfort or pain (e.g. if device becomes hot, or sound is uncomfortably loud) and inform clinician.
- Do not attach or wear a body-worn battery configuration beneath layers of clothing in direct contact with the skin (e.g. in pocket or under child's vest).
- Make sure all cables used by a child are securely attached to their clothing.
- Do not place processor or parts in any household devices (e.g. microwave oven, dryer).
- Do not modify your processor. Warranty will be void if modified.
- Do not use a drying aid that has an Ultra Violet C (UVC) lamp (e.g. do not use the Freedom Dry and Store).
- Store spare magnets safely and away from cards that may have a magnetic strip (e.g. credit cards, bus tickets, etc).

## General Warnings

### Batteries

- Dispose of used batteries promptly and carefully, in accordance with local regulations. Keep away from children.
- Wash hands after handling disposable batteries.
- Do not recharge disposable batteries.
- Do not disassemble, deform, immerse in water or dispose of batteries in fire.
- Do not mix old and new batteries or batteries of different types or brands.
- Replace disposable batteries with high power Cochlear supplied or recommended zinc air batteries only. We recommend Power One p675 Implant Plus zinc/air batteries. We do not recommend using silver oxide or alkaline batteries.
- Only use rechargeable batteries and battery chargers supplied or recommended by Cochlear. Use of other batteries or battery chargers may result in harm or injury.
- Do not allow children to replace batteries without adult supervision.
- Do not touch the battery charger contacts or allow children to use the battery charger without adult supervision.
- Do not short-circuit batteries (e.g. do not let terminals of batteries contact each other, do not place batteries loose in pockets, etc.).
- Store unused batteries in original packaging, in a cool dry place. When processor is not in use, remove the disposable or rechargeable batteries and store separately in a cool dry place.
- Do not expose batteries to heat (e.g. never leave batteries in sunlight, behind a window or in a car).
- Do not use damaged or deformed batteries. If skin or eyes come into contact with battery fluid or liquid, wash out with water and seek medical attention immediately.
- Never put batteries in mouth. If swallowed, contact your physician or local Poison Information Centre.
- Please contact your clinician for battery service.

### Acoustic Component

- Refit and replace hard earmoulds regularly. Hard earmoulds can cause deformation of the ear canal. Have them refitted every two years for adults and more regularly for children.

## Other information

### Product component dimensions (Typical values)

| Component                            | Length  | Width   | Depth   | External diameter |
|--------------------------------------|---------|---------|---------|-------------------|
| CP910 processing unit                | 51.3 mm | 47.7 mm | 9.0 mm  | NA                |
| CP920 processing unit                | 39.1 mm | 47.7 mm | 9.0 mm  | NA                |
| Standard battery module              | 28.0 mm | 9.0 mm  | 19.0 mm | NA                |
| Compact rechargeable battery module  | 19.0 mm | 9.0 mm  | 19.0 mm | NA                |
| Standard rechargeable battery module | 28.0 mm | 9.0 mm  | 19.0 mm | NA                |
| Coil                                 |         |         | 8.0 mm  | 32.0 mm           |

### Weight (Typical values. All weights are measured with medium earhook)

| Component  | Weight |
|--|--------|
| CP910 processor (no battery)                       | 5.7 g  |
| CP920 processing unit (no battery)                 | 5.0 g  |
| CP910 processor with compact rechargeable battery  | 10.5 g |
| CP910 processor with standard rechargeable battery | 13.0 g |
| CP910 processor with Zinc/Air batteries            | 12.8 g |
| CP920 processor with compact rechargeable battery  | 9.8 g  |
| CP920 processor with standard rechargeable battery | 12.4 g |
| CP920 processor with Zinc/Air batteries            | 12.1 g |
| Coil (without coil magnet)                         | 5.0 g  |

## Other information

### Operating characteristics

#### Processing unit

| Characteristic              | Value/range   |
|-----------------------------|---|
| Sound input frequency range | 100 Hz to 8 kHz   |
| Operating voltage           | 2.00 V to 4.25 V  |
| Power consumption           | 20 mW to 100 mW   |
| Push-button functions       | Turn processor on and off, turn telecoil/audio accessory on and off, change program, lock/unlock buttons, change sensitivity or volume level. |

#### Battery module

| Type                                 | Capacity/voltage range  |
|--------------------------------------|---|
| Standard battery module              | Two PR44 (Zinc/Air) button cell batteries. 1.45 V (nominal) each. Cochlear recommends Power One p675 Implant Plus zinc/air batteries. |
| Compact rechargeable battery module  | 120 mAh / 3.0 V to 4.2 V  |
| Standard rechargeable battery module | 205 mAh / 3.0 V to 4.2 V  |

#### Coil

| Characteristic      | Value/range    |
|---------------------|----------------|
| Operating voltage   | 2.0 V to 2.6 V |
| Operating frequency | 5 MHz          |

### Environmental conditions

#### Processing unit, coil, coil cable and accessories

| Condition                   | Minimum | Maximum  |
|-----------------------------|---------|----------|
| Storage temperature         | -40 °C  | +50 °C   |
| Storage relative humidity   | 0% RH   | 90% RH   |
| Operating pressure          | 700 hPa | 1030 hPa |
| Operating temperature       | +5 °C   | +50 °C   |
| Operating relative humidity | 0% RH   | 90% RH   |

#### Standard battery module

Check the battery manufacturer's recommended operating conditions for disposable batteries used in your processor.

The processing unit comprises:

- Two omni-directional microphones for receiving sound.
- An internal telecoil for receiving magnetic fields radiated by phones, neckloops and roomloops.
- Custom analogue and digital integrated circuits with digital signal processing (DSP) and bi-directional wireless communication capabilities.
- A tri-colour visual indication of processor function or problem.
- Two push-buttons to allow user control of key features.
- CP910 only – Custom 4-pin accessory connector for connection of audio accessories.
- Custom 4-pin coil connector for connection of the coil cable.
- A range of earhooks.

The batteries provide power to the processor.

The coil acts as a transformer coupling that transfers energy and data to the implant.

#### Materials

- Processing unit: copolyester
- Battery modules (all types) are made of copolyester
- Coil is made of polypropylene (PP), thermoplastic elastomer (TPE)
- Magnet casing is made of acrylonitrile butadiene styrene (ABS)
- Coil cable sheath is made of polyvinyl chloride (PVC)
- Coil cable plugs are made of PP & TPE
- LiteWear cable is of polyester, TPE & PVC (same construction as the coil cable) sheath is PVC; plugs/strain relief are PP and TPE
- LiteWear connection boots are made of copolyester
- LiteWear retention case is made of polyamide



## Other information

- LiteWear Fixing Aids are made of polycarbonate
- Processor covers are made of polystyrene and polycarbonate
- Coil covers are made of polycarbonate
- Hybrid materials for outer sheathing of the housing adaptor and the outer cable sheathing is made of polyamide
- Custom moulds are made of Stereo Lithography Acrylic covered with a UV-cured lacquer
- Plus Domes are made of elastosil silicone polymer
- Power Domes are made of TPE

## Certification and applied standards

The CP900 series sound processors fulfil the essential requirements listed in Annex 1 of the EC directive 90/385/EEC on Active Implantable Medical Devices as last amended by EC Directive 2007/47/EEC.

**CE0197**

## Equipment classification

Your processor is internally powered equipment Type B as described in the international standard IEC 60601-1:2005, Medical Electrical Equipment - Part 1: Basic Safety and Essential Performance.

## FCC (Federal Communications Commission) and Canadian IC compliance

This device complies with part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications made to this equipment not expressly approved by Cochlear Limited may void the FCC authorization to operate this equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Consult the dealer or an experienced radio/TV technician for help.

FCC ID: WTOP900

IC ID: 8039A-P900

## Other information

### Labelling symbols

The symbols below are found on your processor components and packaging:



Consult instructions for use



Catalogue number




Batch code



Serial number



Temperature limitations



Relative humidity limitations

**Rx Only** By prescription



Certification Japan



Certification Australia



Certification New Zealand



Type B applied part



Recyclable material



Segregate electronic waste



Do not dispose in fire



CE Registration Mark

### Legal statement

The statements made in this guide are believed to be true and correct as of the date of publication. However, specifications are subject to change without notice.

Nucleus® Cochlear™ implant systems are covered by one or more international patents.

© Cochlear Limited 2012

# Cochlear™



**Cochlear Ltd** (ABN 96 002 618 073) 14 Mars Road, Lane Cove NSW 2066, Australia Tel: 61 2 9428 6555 Fax: 61 2 9428 6352

**Cochlear Americas** 13059 E Peakview Avenue, Centennial, CO 80111, USA Tel: 1 303 790 9010 Fax: 1 303 792 9025

**Cochlear AG** European Headquarters, Peter Merian-Weg 4, CH - 4052 Basel, Switzerland Tel: 41 61 205 0404 Fax: 41 61 205 0405

**ECIREP Cochlear Deutschland GmbH & Co. KG** Karl-Wiechert-Allee 76A, D-30625 Hannover Germany Tel: 49 511 542 770 Fax: 49 511 542 7770

**Cochlear Europe Ltd** 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey KT15 2HJ, United Kingdom Tel: 44 1932 87 1500 Fax: 44 1932 87 152

**Nihon Cochlear Co Ltd** Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 024

**Cochlear (HK) Ltd** Unit 1810, Hopewell Centre, 183 Queens Road East, Wan Chai, Hong Kong SAR Tel: 852 2530 5773 Fax: 852 2530 5182

**Cochlear Medical Device (Beijing) Co Ltd** Unit 2208 Gemdale Tower B, 91 Jianguo Road, Chaoyang District, Beijing 100022 P.R. China Tel: 86 10 5909 7800 Fax: 86 10 5909 7900

**Cochlear Ltd** (Singapore Branch) 6 Sin Ming Road, #01-16 Sin Ming Plaza Tower 2, Singapore 575585 Tel: 65 6553 3814 Fax: 65 6451 4

**Cochlear Korea Ltd** 1st floor, Cheongwon building, 828-5, Yuksam dong, Kangnam gu, Seoul, Korea Tel: 82 2 533 4663 Fax: 82 2 533 8

**Cochlear Benelux NV** Schaliënhoevedreef 20i, B - 2800 Mechelen, Belgium Tel: 32 1579 5511 Fax: 32 1579 5500

**Cochlear Italia S.r.l.** Via Larga 33, 40138 Bologna, Italia Tel: 39 051 601 53 11 Fax: 39 051 39 20 62

**Cochlear France S.A.S.** Route de l'Orme aux Merisiers, Z.I. Les Algorithmes - Bât. Homère, 91190 Saint Aubin, France Tel: 33 811 111 993 Fax: 33 160 19

**Cochlear Nordic AB** Konstruktionsvägen 14, SE - 435 33 Mölnlycke, Sweden Tel: 46 31 335 14 61 Fax: 46 31 335 14 60

**Cochlear Tıbbi Cihazlar ve Sağlık Hizmetleri Ltd. Sti.** Cubuklu Mah. Bogazici Cad., Bogazici Plaza No: 6/1, Kavacik

TR - 34805 Beykoz-Istanbul, Turkey Tel: 90 216 538 5900 Fax: 90 216 538 5919

**Cochlear Canada Inc** 2500-120 Adelaide Street West, Toronto, ON M5H 1T1 Canada Tel: 1 416 972 5082 Fax: 1 416 972 5083

**Cochlear Medical Device Company India Pvt Ltd** Ground floor Platina Bldg, Plot no C 59, G Block, Bandra Kurla Complex, Bandra East, Mumbai 400051, India Tel: 91 22 61121111 Fax: 91 22 61121100

Nucleus is a registered trademark of Cochlear Limited.

Cochlear and the elliptical logo, Snugfit, Mic Lock and Freedom are trademarks of Cochlear Limited.

Dry and Store is a registered trademark of Ear Technology Corporation.

The statements made in this guide are believed to be true and correct in every detail as of the date of publication. However, specifications are subject to change without notice.

*Hear now. And always*

Printed in Australia  
323649 ISS3 APR12