

Cochlear™

Nucleus® CR210 Remote Control

User Guide



Hear now. And always



Symbols

**Note**

Important information or advice. Can save inconvenience.

**Tip**

Time saving hint or saves inconvenience.

**Caution (no harm)**

Special care to be taken to ensure safety and effectiveness.
Could cause damage to equipment.

**Warning (harmful)**

Potential safety hazards and serious adverse reactions.
Could cause harm to person.

Contents

About	4
Turn on.....	6
Turn on and off.....	6
Pair	7
Out of range	8
Use.....	9
Change programs.....	9
Change volume or sensitivity.....	10
Talk on the phone.....	11
Battery.....	13
Care.....	15
General warnings and cautions	16
Other information.....	17

About


About

Cochlear Nucleus CR210 Remote Control


The Cochlear Nucleus CR210 Remote Control is a simple hand-held controller for your Cochlear Nucleus CP900 Series Sound Processors. The remote allows you to:

1. change programs
2. turn volume or sensitivity up and down
3. turn the telecoil on and off
4. turn auto telecoil on and off
5. turn audio accessories on and off

If you are a bilateral recipient, your remote control can communicate with both your processors at once. For more information on your sound processor, refer to the Cochlear Nucleus CP900 Series Sound Processor User Guide.

 **Note**

"Volume" describes the level of the sounds you hear. "Sensitivity" describes the *range* of sounds you hear e.g. soft, background, close or far away. Usually, your clinician will set up your remote control to control volume. But your clinician can set it up to control sensitivity instead.



Cochlear CR210 Basic Remote Control

Turn On

Turn on

Turn on and off



Slide the switch up to turn on.

When the switch is in the up position, the display will "go to sleep" when not in use. Press any button on the keypad to "wake" it again.



Slide the switch down to turn off.



Tip

Turn off your remote to stop it "waking" and draining the battery when not in use e.g. when carried in your pocket.

Pair

Your clinician may have already paired your remote with your sound processor. If not, you will need to pair it the first time you use it.



If your remote is not paired to any processor, the screen appears as above and the display "blinks" when you press keys.

To pair your processor and remote:



- ① Turn on your processor and remote control. They must both be "awake" and on at all times when pairing.



- ② Place the coil on the back of your remote control.

Image of coil over battery will be briefed into next version.

Turn On



- ③ The screen will show pairing in progress.



- ④ The screen then shows successful pairing of left (L) or right (R) processor.

- ⑤ Bilateral users only - repeat step 2 for your second sound processor.

Out of range



The screen above displays when the processor is out of range or experiencing interference.

Use

Change programs

Your clinician may program your sound processor with up to four programs to control with your remote.

Press to move to the next program

Press to move to the previous program



Ask your clinician to write your programs below

P8 _____

P8 _____

P8 _____

P8 _____

Change volume or sensitivity

Usually, your clinician will set up your remote control to control the volume of the sound you hear. But your clinician can instead set it up to control sensitivity. "Sensitivity" describes the range of sounds you hear e.g. soft, background, close or far away.



Press to turn up the volume or change sensitivity.

Press to turn down the volume or change sensitivity.

Talk on the phone

Your sound processor's telecoil (t-switch) lets you listen to hearing loops and hear voices on the telephone more clearly. You can control it with the telecoil button on your sound processor. Your clinician will set up your remote with ONE of the following functions:

1. **No telecoil function** - you can not control the telecoil with your remote. The telecoil button on your remote will not do anything.
2. **Telecoil only function** - you can turn the telecoil on manually with the remote when you want to listen with the telecoil, and turn it off when you are finished. You will not be able to turn on Auto Telecoil.

To use the telecoil only function:



To turn on the telecoil: Press the button to display the telecoil icon.



To turn off the telecoil: Press again to turn the telecoil off.

3. Telecoil AND auto telecoil function - you will be able to turn both the telecoil and the auto telecoil feature on and off with the remote as needed. When auto telecoil is left on, the telecoil in your sound processor turns on automatically whenever a hearing loop or phone is in range, and turns off a few seconds after you move the processor away from the source.

To control the Telecoil and Auto Telecoil feature:



To turn on the telecoil: Press the button to display the telecoil icon.



To turn on auto telecoil: Press again to display the auto telecoil icon.



To turn off the telecoil: Press again to turn the telecoil off.

Battery

Battery

Change the batteries whenever you can see the low battery indicator or the screen stays blank when you push the side switch upwards into the unlock position.



To change the battery:



① Unlock the battery cover by turning the screw anticlockwise with a Phillips head screwdriver.



② Slide off the battery cover and remove the battery.



- ③ Insert the new battery. Make sure the plus sign "+" is facing up.



- ④ Replace the battery cover bottom edge first and turn the screw clockwise to lock the cover in place.

⚠ Caution

Use disposable CR2032/5004LC 3V standard lithium coin cell batteries ONLY in your remote control. Cochlear does not recommend any other battery type. Do not use rechargeable batteries.

Care


Every day

- ❑ To keep your remote control clean, wipe its display and buttons with a soft dry cloth.


Keep your remote dry. If it ever gets wet:



Remove the battery immediately.



Dry it with a soft dry cloth.



Leave it in a well ventilated place to dry for 12 hours. Do not use a drying kit.



Put a fresh battery in place.

General warnings and cautions

- Use a soft dry cloth ONLY to clean the surface of your remote control.
- Do not use alcohol, household cleaners or abrasive materials to clean your remote control.
- Do not use the Zephyr dry aid kit or any household appliance to dry your remote control.
- Do not bend or compress the display on your remote control.
- Do not disassemble, modify, deform or immerse your remote control in water. If your remote control is not working, return it to your clinician. Warranty will be void if opened.
- Dispose of electrical components in accordance with your local regulations.
- Do not allow metal objects e.g. keys near the display screen.
- The coin battery could pose a choking hazard to small children. Keep the battery cover locked at all times. Keep the remote out of reach of very young children.

Other information

The remote control wireless communication link operates in the 2.4 GHz ISM band using GFSK (Gaussian Frequency Shift Keying) on 10 channels. The link uses a proprietary bi-directional communication protocol and operates over a distance of up to 2m from the CP900 series processor. When interference is present, the wireless communication link switches between the 10 channels to find a channel where the interference least affects the operation of the link. The remote control indicates via its display when the processor is not within operating distance, and when the link has been interrupted due to interference.

Specifications

Physical configuration

The remote control comprises:

- Custom analogue and digital integrated circuits based on a micro-processor design, with bi-directional wireless communication capabilities.
- Push-buttons and a slider for control of processor features and its own system features.
- A display for indicating system status and activity.
- In-built coil sensor.
- Lanyard pin.

Materials

The remote control housing and buttons are made of Polycarbonate and Acrylonitrile Butadiene Styrene (ABS). The display is made of ABS.

Other information

Dimensions

Length	Width	Height
66 mm	36 mm	10.5 mm

Weight

The remote control with battery in place weighs approximately 35 g.

Operating characteristics

Characteristic	Range
Wireless technology	Proprietary Low Power Bidirectional Wireless Link
RF frequency	2.4 GHz
Operating battery voltage	2.0 V - 3.0 V
Battery capacity	230 mAh
Push-key functions	Remote control lock and unlock, telecoil/auto telecoil on and off, program change, volume or sensitivity settings change
Wireless transmission range	Up to 2 m
Display	Custom, icon based monochrome LCD

Environmental conditions

Condition	Minimum	Maximum
Storage temperature	-30 °C (-40 °F)	+50 °C (+122 °F)
Storage relative humidity	0% RH	90% RH
Operating temperature	+5 °C (+41 °F)	+50 °C (+122 °F)
Operating relative humidity	0% RH	90% RH

Certification and applied standards

The CR210 Remote Control fulfils the essential requirements listed in Annex 1 of the EC directive 90/385/EEC on Active Implantable Medical Devices as last amended by EC Directive 2007/47/EEC.

CE0197

Equipment classification

Your remote control is internally powered equipment Type B as described in the international standard IEC 60601-1:2005, Medical Electrical Equipment - Part 1 Basic Safety and Essential Performance.

FCC (Federal Communications Commission) and Canadian IC compliance

This device complies with part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications made to this equipment not expressly approved by Cochlear Limited may void the FCC authorization to operate this equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC ID: WTOR200BA

IC ID: 8039A-R200BA

Labelling symbols



Consult instructions for use



Storage temperature limits



Relative humidity limits



Disposal - do not dispose in fire



Disposal – dispose of electrical components in accordance with your local regulations



Certification - Japan



Certification - Australia



Certification - New Zealand



Serial number



Date of manufacture



Type B applied part

Rx Only

This device restricted to sale by or on the order of a physician.



Recycle



CE Registration Mark

Legal statement

The statements made in this guide are believed to be true and correct as of the date of publication. However, specifications are subject to change without notice.

Nucleus® cochlear implant systems are covered by one or more international patents.

© Cochlear Limited 2012

Cochlear™



Cochlear Ltd (ABN 96 002 618 073) 14 Mars Road, Lane Cove NSW 2066, Australia Tel: 61 2 9428 6555 Fax: 61 2 9428 6352

Cochlear Americas 13059 E Peakview Avenue, Centennial, CO 80111, USA Tel: 1 303 790 9010 Fax: 1 303 792 9025

Cochlear AG European Headquarters, Peter Merian-Weg 4, CH - 4052 Basel, Switzerland Tel: 41 61 205 0404 Fax: 41 61 205 0405

CECIREP Cochlear Deutschland GmbH & Co. KG Karl-Wiechert-Allee 76A, D-30625 Hannover

Germany Tel: 49 511 542 770 Fax: 49 511 542 7770

Cochlear Europe Ltd 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey KT15 2HJ, United Kingdom Tel: 44 1932 87 1500 Fax: 44 1932 87 1520

Nihon Cochlear Co Ltd Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 024

Cochlear (HK) Ltd Unit 1810, Hopewell Centre, 183 Queens Road East, Wan Chai, Hong Kong SAR Tel: 852 2530 5773 Fax: 852 2530 5183

Cochlear Medical Device (Beijing) Co Ltd Unit 2208 Gemdale Tower B, 91 Jianguo Road, Chaoyang District, Beijing 100022

P.R. China Tel: 86 10 5909 7800 Fax: 86 10 5909 7900

Cochlear Ltd (Singapore Branch) 6 Sin Ming Road, #01-16 Sin Ming Plaza Tower 2, Singapore 575585 Tel: 65 6553 3814 Fax: 65 6451 4

Cochlear Korea Ltd 1st floor, Cheongwon building, 828-5, Yuksam dong, Kangnam gu, Seoul, Korea Tel: 82 2 533 4663 Fax: 82 2 533 8

Cochlear Benelux NV Schaliënhoedreef 20i, B - 2800 Mechelen, Belgium Tel: 32 1579 5511 Fax: 32 1579 5500

Cochlear Italia S.r.l. Via Larga 33, 40138 Bologna, Italia Tel: 39 051 601 53 11 Fax: 39 051 39 20 62

Cochlear France S.A.S. Route de l'Orme aux Menisiers, Z.I. Les Algorithmes - Bât. Homère, 91190 Saint Aubin, France Tel: 33 811 111 993 Fax: 33 160 191

Cochlear Nordic AB Konstruktionsvägen 14, SE - 435 33 Mölnlycke, Sweden Tel: 46 31 335 14 61 Fax: 46 31 335 14 60

Cochlear Tibbi Cihazlar ve Sağlık Hizmetleri Ltd. Sti. Cubuklu Mah. Bogazici Cad., Bogazici Plaza No: 6/1, Kavacic

TR - 34805 Beykoz-Istanbul, Turkey Tel: 90 216 538 5900 Fax: 90 216 538 5919

Cochlear Canada Inc 2500-120 Adelaide Street West, Toronto, ON M5H 1T1 Canada Tel: 1 416 972 5082 Fax: 1 416 972 5083

Cochlear Medical Device Company India Pvt Ltd Ground floor Platina Bldg, Plot no C 59, G Block, Bandra Kurla Complex,

Bandra East, Mumbai 400051, India Tel: 91 22 61121111 Fax: 91 22 61121100

Nucleus is a registered trademark of Cochlear Limited.

Cochlear and the elliptical logo are trademarks of Cochlear Limited.

The statements made in this guide are believed to be true and correct in every detail as of the date of publication. However, specifications are subject to change without notice.

Hear now. And always

Printed in Australia
323650 ISS2 APR12