

Owner's Booklet

Remote Control

for your SENSOTREND
Insulin Pump



Diagnostics

D[®] DISETRONIC

Contents

1	The remote control for your SENSOTREND ¹ insulin pump	4
2	Handling Keys Indicator light	12
3	Preparation	18
4	Adding the remote control ID number	22
5	Program a Standard bolus	28
6	Cancel a Standard bolus	38
7	Recommendations for daily use	46
8	Troubleshooting	50
9	Technical features	58
10	Symbols	64
11	Index	68
12	Addresses	70
13	Warranty	72

1 The remote control for your SENSOTREND¹ insulin pump

Congratulations on your new remote control for the SENSOTREND insulin pump. The remote control allows for discreet one-handed programming of a Standard bolus without touching or looking at your insulin pump. Your insulin pump beeps and vibrates to assist you with the programming and the indicator light of the remote control confirms each press of a key with a light signal.

This owner's booklet provides you with the necessary information for safe and efficient use of the remote control with your insulin pump. The remote control is an accessory for a medical device and must be trained on

¹ SENSOTREND is a trademark of a Member of the Roche Group.

its operation prior to use. Prior to using your SENSOTREND insulin pump read the Owner's Booklet to become familiar with its safe use and the appropriate warnings and precautions for its operation.

Prior to using your SENSOTREND insulin pump with your remote control be aware of the following warnings:



Warning: Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings. Work with your doctor or healthcare team to determine the timing,

1 The remote control for your SENSOTREND insulin pump

amount and type of bolus you need to deliver. Make sure you know your personal carbohydrate insulin ratio and your correction bolus ratio.



Warning: Do not use the SENSOTREND insulin pump and remote control if you are not able to read the display, hear the beeps or you are not able to feel the vibrations. The display, beeps and vibrations indicate critical alert information and must be addressed immediately by the insulin pump user.

Failure to see, hear or feel the alerts from the insulin pump could result in serious harm or death of the user.



Warning: Keep the remote control out of unauthorized third persons and of children. Activation of the remote control may lead to unintended bolus delivery.



Warning: Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings.



Warning: If your insulin pump “times out” before you finish programming a bolus, no bolus is delivered. Check the bolus history screens on your insulin pump for the bolus amount delivered with time and date and program a new bolus if necessary (see section 4.8.1 of your SENSOTREND owner's booklet).



Warning: Avoid deliberate contact with water. Check daily that your insulin pump and the remote control are not chipped, cracked or damaged in any way. In the presence of chips and cracks, water, dust, insulin, or other foreign substances may enter your insulin pump and your remote control and lead to malfunction.



Warning: You must inspect the remote control and your insulin pump immediately after an accidental immersion in water or after any accidental contact with foreign liquids that could damage or contaminate the insulin pump system.



Warning: Do not perform any servicing or repair on the remote control by yourself or by an unauthorized person. Any changes or modifications to the remote control that are not expressly approved by Disetronic may

void the warranty. For any questions please contact your local Customer Care Area.



Warning: Avoid electromagnetic fields of radar or antenna installations, high-voltage sources, X-Ray sources, MRI, CAT scan or other sources of electrical current. Do not use your insulin pump and the remote control in such areas. Electromagnetic fields may cause malfunction to your insulin pump and remote control. Always stop your insulin pump beforehand. In other cases, the insulin delivery may be stopped immediately and an error E7: ELECTRONIC ERROR occurs.



Warning: Your insulin pump has not been tested with cardiac pacemakers. Do not use the SENSOTREND insulin pump if you have a cardiac pacemaker.



Warning: Any changes or modification to the remote control or insulin pump that are not expressly approved by Disetronic could void

your ability to operate the SENSOTREND insulin pump.



Warning: Be sure that the remote control ID numbers programmed into your insulin pump belong to your remote control. Do not have other remote ID numbers programmed into the insulin pump besides the remote(s) that you are personally using. Keep the remote control out of reach of children. Do not let your remote control be handled by unauthorized persons as a bolus could be accidentally delivered.

Models available

Two different remote control models are available:

- ▶ The **pocket remote control** which looks like an ordinary key chain.
- ▶ The **watch remote control** is a wristwatch of Swiss quality that integrates both the time piece and the remote control.

The pocket and watch remote control are available only from Disetronic. Please consult your SENSOTREND owner's booklet or contact your local Customer Care Area for more information.


For information on your insulin pump, please refer to the SENSOTREND owner's booklet.

2 Handling


Keys and indicator light of the remote control

The use of the remote control requires three keys. With each key press you hear a soft click (pocket remote control only) and the indicator light confirms that a key has been pressed.

i Please note: The remote control does not respond with beeps and vibrations.



Pocket remote control



Watch remote control

⏻ **Activation** to **activate** the communication between the remote control and your SENSOTREND insulin pump. The indicator light will remain on as **long** as you press the ⏻ key.



to **program** and **cancel** a Standard bolus, and to **increase** or **reset** the bolus amount. For each press of the or key the indicator light flashes.

If you accidentally press the and (or) keys simultaneously, the indicator light remains on during the press and your insulin pump tries to establish the communication. If the communication is already active, your insulin pump will beep** and vibrate periodically (a short vibration every second), if is pressed.

Please note the , , and symbols refer to the keys of your SENSOTREND insulin pump.

** Depending on the beep volume setting.

Your SENSOTREND insulin pump beeps** when



you **activate** the remote control communication, or the **bolus delivery starts**.



you **start** the **Standard bolus menu**, or you **reset** the **bolus amount** to zero.



you **increase** to the desired **bolus amount**, or during the **confirmation** of the **bolus amount**.



you **cancel** the **bolus**.

Your insulin pump beeps only if its beep volume has been turned on (see section 4.7.2 in your SENSOTREND owner's booklet).

** Depending on the beep volume setting.

Your **SENSOTREND** insulin pump vibrates when

(•) you **activate** the remote control (short **communication**, or you **increase** vibration) to the desired **bolus amount**.

(••) you **start** the **Standard bolus** (long **menu, during the confirmation** of vibration) the **bolus amount**, you **reset** the **bolus amount** to zero, or you **cancel** the **bolus**.

Pocket remote control

Use your thumb to slide the remote control out of its cover until you see all three keys. After use, close the cover properly to avoid unintentional key operation and to protect the remote control.


Watch remote control

Pull gently at the crown of the watch remote control until you hear a soft click.

Set the time as desired by turning the crown. When finished push the crown in.


The watch remote control has a conventional watchstrap. It may be replaced according to your personal preference by an authorized jeweller.

Position of the remote control to your **SENSOTREND insulin pump**



The remote control works best at an arm's length distance from your insulin pump. Position the remote control in the vicinity of your insulin pump to establish the communication.

3 Preparation



Warning: Be sure that the remote control ID numbers programmed into your insulin pump belong to your remote control. Do not have other remote ID numbers programmed into the pump besides the remote(s) that you are personally using. Keep the remote control out of reach of children. Do not let your remote control be handled by unauthorized persons as a bolus could be accidentally delivered.

Each remote control has an unique ID number to provide secure communication between your insulin pump and the remote control. This ID number is engraved in the rear casing of the remote control. To view the ID number of the pocket remote control use your thumb to slide the remote control out of its cover. The ID number has to be added to the REGISTERED IDs list of your SENSOTREND insulin pump before use.

Add only the ID(s) of your remote control(s) to your personal SENSOTREND insulin pump. If you are using a different SENSOTREND insulin pump, make sure to add the remote control ID(s) prior to use. A remote control, which is not registered on your insulin pump, cannot activate the communication.

3 Preparation



In order to make use of the remote control, the remote control function on your insulin pump must be turned on. Please refer to the sections 5.1 and 5.3 in your SENSOTREND owner's booklet for further information.


On the rear of the casing you will find additional information engraved. These include the SN and the REF number. Please refer to the section 10 in this owner's booklet for further information on the labelling.

4 Adding the remote control ID number



Warning: Be sure that the remote control ID numbers programmed into your insulin pump belong to your remote control. Do not have other remote ID numbers programmed into the insulin pump besides the remote(s) that you are personally using. Keep the remote control out of reach of children. Do not let your remote control be handled by unauthorized persons as a bolus could be accidentally delivered.

Use the remote control and your SENSOTREND insulin pump and check your insulin pump's display during the entire registry process.


1. Press your insulin pump's  key to move to the REMOTE CONTROL menu. Make sure that the REMOTE CONTROL


menu is turned on. Please refer to section 5.1 of your owner's booklet for further information.



2. Select with . Check that the current status is ON. (Press  or  to turn it to ON () , if necessary.)





3. Press  to move to the ADD ID screen. Any currently registered IDs are displayed.




4. If the ID number of your remote control is not displayed, hold the remote control within 75 cm (2.5 feet) of your insulin pump.


4 Adding the remote control ID number

Press . Your insulin pump scans for 30 seconds for a valid remote control ID.



ADD ID
SCANNING

5. Within these 30 seconds press the  key of the remote control to register its ID.
6. After the ID is added, the REGISTERED IDs screen automatically appears. Compare the display with the engraved ID number on the remote control casing to make sure that the ID number is correctly registered.





REGISTERED IDs
12345678

If no (additional) ID number is registered, the ADD ID screen will appear again. Hold your remote control closer to your insulin pump and repeat steps 4 to 6.

i If you want to **register a further remote control now**, repeat steps 4 to 6 with the additional remote control.

7. Save and exit with .

- i**
- You can register up to three different remote controls to your insulin pump. For safety reasons (unintended bolus activation) we recommend to register only the one remote control in use on your insulin pump.
 - You can **review any registered remote control(s)** at the REGISTERED IDs screen.
 - Check that the remote control function is turned on () in the quick info screen. To view the quick info screen press  from the **RUN** or **STOP** screen.

4 Adding the remote control ID number

- If you do not use the remote control(s), turn the REMOTE CONTROL function to OFF (🔒, see step 2).



Warning: Keep the remote control out of unauthorized third persons and of children. Activation of the remote control may lead to unintended bolus delivery.

5 Program a Standard bolus



Warning: Incorrect programming of your insulin pump may cause inappropriate insulin delivery. Your insulin pump must be programmed with your personal settings prior to starting insulin pump therapy. Do not operate your insulin pump without knowing your personal settings. If you are not sure about your personal settings or you are less experienced, have your doctor or healthcare team check your personal settings.




Warning: If your insulin pump “times out” before you finish programming a bolus, no bolus is delivered. Check the bolus history screens on your insulin pump for the bolus amount delivered with time and date and program a new bolus if necessary (see section 4.8.1 of your SENSOTREND owner's booklet).

The remote control allows you to easily and discreetly program a bolus of insulin. Your SENSOTREND insulin pump must be in RUN. Each key press and each bolus increment you program is confirmed with beeps** and vibrations. Pay attention and count these confirming beeps** and vibrations for secure programming.

The remote programming follows a defined sequence which involves key presses and pre-set time intervals.




If you use the **KeyLock** function on your insulin pump, the keys are automatically unlocked as you activate the communication with the  key and remains unlocked until the bolus is delivered.

** Depending on the beep volume setting.

5 Program a Standard bolus

i The **bolus increment** is initially set to 0.5 units per key press. You can change it to 0.1, 0.2, 0.5, 1.0 or 2.0 units on your insulin pump (See section 5.5.3 in your SENSOTREND owner's booklet).

Program a Standard bolus



Use the remote control, ...

1. Press the  key till communication is activated.

... your **SENSOTREND** insulin pump responds with

beeps** vibrations



display



** Depending on the beep volume setting.

5 Program a Standard bolus

Use the remote control, ...

2. Press the  key for 3 seconds (or ***).



... your **SENSOTREND** insulin pump responds with

beeps** vibrations








display






Do not wait longer than 5 seconds to press the  or  key. Otherwise the communication “times out” and stops.

** Depending on the beep volume setting.

*** This programming sequence is “reversible”, meaning you can start with either the  or  key. If you start the  key, program the bolus amount with  and reset with .

Use the remote control, ...

3. Press  +  + ... (or  + ...***) to program the desired bolus amount.

... your **SENSOTREND** insulin pump responds with

beeps**

vibrations


display



+

+

...



+

...








Use the remote control, ...

4. Wait for 5 seconds.

Your **SENSOTREND** insulin pump confirms the bolus amount with

** Depending on the beep volume setting.

*** This programming sequence is “reversible”, meaning you can start with either the  or  key. If you start the  key, program the bolus amount with  and reset with .

5 Program a Standard bolus

beeps** vibrations display




Use the remote control, ...

5. After 5 seconds (start delay), the bolus delivery starts.

... your **SENSOTREND** insulin pump responds with

beeps** vibrations display



If communication does not activate, move the remote control closer to your insulin pump while pressing the  key until you hear a beep** sequence and a vibration of your insulin pump.

** Depending on the beep volume setting.


Once the communication is established, keep the remote control at the same close position to your insulin pump while programming the Standard bolus. If your insulin pump does not provide any feedback, the communication has not been established. Please refer to the section 8 in this owner's booklet for further information.

Start delays

Start delays in the programming sequence allow you to cancel a bolus even before its delivery has started (see section 6 in this owner's booklet).

After programming the bolus amount, the bolus amount value blinks.




Between confirmation of amount and bolus delivery, the  will blink (start delay).



5 Program a Standard bolus

Communication timeouts

If you do not press the keys within 5 seconds during the fixed programming sequence, the communication between the remote control and your insulin pump “times out” and your insulin pump returns to **RUN**. If you want to continue and program a bolus, you may activate the communication again with the .




Please refer to step 1 in this section for further information.

6 Cancel a Standard bolus

You have six options to cancel a Standard bolus:

– **during programming** (the bolus amount blinks on your insulin pump's display)

Use the remote control, ...

1. Press the  key to reset the bolus amount to 0.0 units (or  if you started programming with .

... your SENSOTREND insulin pump responds with

beeps** vibrations





(-)

display

STANDARD BOLUS
L 0.0U

Use the remote control, ...

2. Your insulin pump times out after 5 seconds if  or  are not pressed.

... your SENSOTREND insulin pump responds with

beeps** vibrations



–


display

15:28
10.0 U/h

or

Use your SENSOTREND insulin pump

Press the  or  key for 3 seconds on your insulin pump.

This will result in an alert A8: BOLUS CANCELLED. Please confirm the alert A8 twice with  on your insulin pump.

** Depending on the beep volume setting.

** Depending on the beep volume setting.

6 Cancel a Standard bolus

... **your SENSOTREND insulin pump responds with**


beeps*



vibrations



display



No bolus is delivered!

* Depending on the beep volume and alarm signal settings.

– **during confirmation** (when the beeps** and vibrations occur) **or**
during the start delay (⊥ blinks on your insulin pump's display)


Use the remote control, ...

Press the or key until you hear a melody*.

This will result in an alert A8: BOLUS CANCELLED. Please confirm the alert A8 twice with on your insulin pump.

... **your SENSOTREND insulin pump responds with**


beeps*



vibrations*



display



or


* Depending on the beep volume and alarm signal settings.

** Depending on the beep volume setting.




6 Cancel a Standard bolus

Use your **SENSOTREND** insulin pump

Press the  or  key on your insulin pump. You hear a melody*.

This will result in an alert A8: BOLUS CANCELLED. Please confirm the alert A8 twice with  on your insulin pump.

... your **SENSOTREND** insulin pump responds with

beeps*	vibrations*	display
		


No bolus is delivered!

* Depending on the beep volume and alarm signal settings.




– **during bolus delivery** (countdown of bolus amount on your insulin pump's display)

Use the remote control, ...

Press the  or  key until you hear a melody*.

This will result in an alert A8: BOLUS CANCELLED. Please confirm the alert A8 twice with  on your insulin pump.

... your **SENSOTREND** insulin pump responds with

beeps*	vibrations*	display
		


or

* Depending on the beep volume and alarm signal settings.

6 Cancel a Standard bolus

Use your **SENSOTREND** insulin pump

Press the  or  key for 3 seconds on your insulin pump.

This will result in an alert A8: BOLUS CANCELLED. Please confirm the alert A8 twice with  on your insulin pump.

...your **SENSOTREND** insulin pumps responds with


beeps*



vibrations*



display



Bolus delivery is interrupted!

Ensure that the cancellation was intended and program a new bolus, if necessary.

* Depending on the beep volume and alarm signal settings.

If the cancellation occurred during the bolus delivery, please go to your insulin pump's BOLUS HISTORY in the INFORMATION MENU (see section 4.8.1 in your SENSOTREND owner's booklet) to check for the actual amount of the delivered bolus. The delivered amount is listed in the initial bolus history screen, if no other bolus was delivered.

In case of re-programming a new bolus, please take this delivered amount into account.

7 Recommendations for daily use

Cleaning: Use only a soft cloth to clean the remote control. Never use solvents, strong detergents, corrosive chemicals or similar substances. If the remote control had an accidental contact with such a liquid, clean it with a damp cloth.

Contact with water and other liquids:



Warning: Avoid deliberate contact with water. Check daily that your insulin pump and the remote control are not chipped, cracked or damaged in any way. In the presence of chips and cracks, water, dust, insulin, or other foreign substances may enter your insulin pump and your remote control and lead to malfunction.



Warning: You must inspect the remote control and your insulin pump immediately after an accidental immersion in water or after any accidental contact with foreign liquids that could damage or contaminate the insulin pump system.

Use a soft cloth to clean and/or dry the outside casing of the remote control if necessary.

Please refer to the sections 8.1.3 and 9.2 in your SENSOTREND owner's booklet for further information.

In case of contact with a liquid, dry the remote control gently with a soft cloth. For the pocket remote control, slide it open, and place it leaning in an upright position with the fixing eyelet pointing downwards to let it dry.

7 Recommendations for daily use

Heat and cold: Do not expose the remote control to extremes in temperature and do not exceed the admissible temperature range. Fast temperature changes combined with water or high humidity may harm the remote control. Do not take the remote control e.g. into saunas, baths, hot tubs, Jacuzzis or whirlpools.

Please refer to section 9.1 in this owner's booklet for further information.

Dropping or damaging the remote control:

Avoid shocks and hits to the remote control as this may break the remote control's seal against liquids. In case of suspected damage inspect the remote control for chips and cracks and press a key to check for the indicator light. If the remote control is damaged, please call your local Customer Care Area.

8 Troubleshooting



Warning: Do not perform any servicing or repair on the remote control by yourself or by an unauthorized person. Any changes or modifications to the remote control that are not expressly approved by Disetronic may void the warranty.

For any questions please contact your local Customer Care Area.

Disturbance of communication




Warning: Avoid electromagnetic fields of radar or antenna installations, high-voltage sources, X-Ray sources, MRI, CAT scan or other sources of electrical current. Do not use your insulin pump and the remote control in such areas. Electromagnetic fields may cause malfunction to your insulin pump and remote control. Always stop your

insulin pump beforehand. In other cases, the insulin delivery may be stopped immediately and an error E7: ELECTRONIC ERROR occurs.

Keep the remote control at a secure distance of at least 60 cm (25 inches) from potential sources of electromagnetic interference.



Warning: Your insulin pump has not been tested with cardiac pacemakers. Do not use the SENSOTREND insulin pump if you have a cardiac pacemaker.

If your insulin pump does not provide feedback to the key presses on the remote control, the communication may be interrupted. Move the remote control into the vicinity of your insulin pump while pressing the  key until your insulin pump responds. Once communication is active, keep the remote control in the same vicinity

8 Troubleshooting

of your insulin pump while programming the bolus.

Please note that there may be some potential for interference by metal objects that are placed within the path of your insulin pump and the remote control. This interference may disturb or interrupt communication between the remote control and your insulin pump.

If you are not able to establish communication between the remote control and your insulin pump, check that the ID number of the remote control is listed in your insulin pump REGISTERED IDs list in the REMOTE CONTROL menu. If it is not listed, add the ID number of your remote control. Please refer to the section 5.3 in your SENSOTREND owner's booklet for further information.

In case that all these corrective actions are not successful please contact your local Customer Care Area.

Battery of the ...

... Pocket remote control

A series of 10 quick flashes by the indicator light of the remote control indicates that the battery is low. Please contact your local Customer Care area for information on obtaining another pocket remote control.

... Watch remote control

The watch remote control has two batteries, one for the remote control function and one to power the watch.

– Remote control function:

A series of 10 quick flashes by the remote control indicates that the battery is low and needs replacement.

– **Watch function:**

When the battery for the watch gets low, the second hand will pause for 4 seconds prior to 4 seconds in one step.

When the battery is low contact your local Customer Care Area to replace the battery.

Repair and disposal ...

Do not modify the remote control in any way and do not perform any service or repair on the remote control. This might damage the seal against water as well as the functioning of your remote control, and in any case would void the warranty.

A remote control may be returned to Disetronic when approved by your local Customer Care Area.

... Pocket remote control

When the pocket remote control is faulty, damaged or the indicator light does not flash as described, please contact your local Customer Care Area.

... Watch remote control

When the watch remote control is faulty or a battery is low, please contact your local Customer Care Area.

Please find the corresponding address at the end of these instructions. Remember to include **a note** with your name, address, and a description of the reason for returning the remote control.

Disposal

After the end of its operation time the battery cannot be replaced in the pocket remote control. Please do not dispose of the remote control in your normal garbage. Return the remote control to your local Customer Care Area.

9 Technical features

9.1 Technical data

Pocket remote control

Water resistance

IPX7: Protected against unintentional water contact, like splashes.

Temperature range

Operating: 0 to +55 °C (+32 to +131 °F)
Storage: -20 to +55 °C (-4 to +131 °F)

Memory None

Battery

Lithium 3V, CR2032
Approx. two years operating time

Communication range

Up to 75 cm to 150 cm (up to 2.5 to 5 feet)

Watch remote control

Water resistance

IPX7: Protected against unintentional water contact, like splashes.

Temperature range

Operating: 0 to +55 °C (+32 to +131 °F)
Storage: -20 to +55 °C (-4 to +131 °F)

Memory None

Watch model

Ronda 762, 2.5mm high, 2 hands

Battery

RC function: Lithium 3V, CR1620
Approx. 1 year operating time
Watch Lithium V364 SR 621
function: Approx. 4 years operating time

Communication range

Up to 75 cm to 100 cm (2.5 to 3 feet)

9.2 Environment

The remote control is equipped with coded low power wireless communication technology. It has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and with RS-210 of Industry Canada. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. this device must accept any interference received, including interference that may cause undesired operation.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. Care should be taken for the remote control to be used in such a manner that the potential for human contact during normal operation is minimized.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. There is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate your insulin pump and/or the remote control.
- Increase the separation between your insulin pump/remote control and the source of the receiving/emitting interference.

9 Technical features

- ▶ Avoid electromagnetic fields of monitors, fluorescent lamps, radar or antenna installations, high-voltage sources, X-Ray sources, MRI, CAT scan or similar equipment.
- ▶ Please contact your local Customer Care Area for assistance.



Warning: Any changes or modification to the remote control or insulin pump that are not expressly approved by Disetronic could void your ability to operate the SENSOTREND insulin pump.

10 Symbols

ID Identification number of the remote control, necessary for communication with your insulin pump

IPX 7 Symbol for protection against the effects of temporary immersion in water, according to IEC 60529: Up to 30 minutes and 1 meter (40 inches) depth

SN Serial number



See instructions in this owner's booklet or in your SENSOTREND owner's booklet

REF Item number

swiss made Manufactured in Switzerland

CE 0123 Marking of conformity to the European Medical Device Directive MDD 93/42/EEC with the number of the notified body



Date of manufacture

FCC ID Federal Communications Commission, part 15 on radio frequency devices



Electronic device of type BF according to the standard IEC 60601-1. Protection against electric shock.



Keys of the remote control (see section 2 in this owner's booklet)



Keys of your SENSOTREND insulin pump (see section 2.1.1 in your SENSOTREND owner's booklet)



SENSOTREND's beeps (see section 2 in this owner's booklet)



SENSOTREND's vibrations (see section 2 in this owner's booklet)

10 Symbols



Do not throw away

LOT

Batch number

11 Index

Add ID number	22
Alert A8: BOLUS CANCELLED	39, 44
Battery	53
Beeps	15
Bolus history	44
Bolus increment	30
Cancel a Standard bolus	38
during bolus delivery	43
during confirmation	41
during programming	38
during start delay	41
Cold	48
Communication	34
Communication timeouts	35
Electromagnetic fields	50
Environment	60
Heat	48

ID number	19
Indicator light	12
KeyLock	30
Keys	12
Pocket remote control	12, 58
Position of the remote control	17
Preparation	18
Program a Standard bolus	31
Recommendations	46
Repair	54
Start delays	35
Symbols	64
Technical data	58
Troubleshooting	50
Vibrations	16
Watch remote control	12, 59
Water	46

12 Addresses

For patients in the Netherlands:

Your local Customer Care Area for your assistance and distributed by:

Roche Diagnostics Nederland BV
Consumer Systems
Transistorstraat 41
NL-1322 CK Almere
The Netherlands

Phone 0800 - 022 05 85

Fax +31 (0)36 53 94 231

Manufactured by:

Disetronic Medical Systems AG

Kirchbergstrasse 190

CH-3401 Burgdorf, Switzerland

Phone: +41 34 424 2424

Fax: +41 34 424 2400

www.disetronic.com

13 Warranty

Disetronic offers an extensive, one-year, global warranty, valid from the date of purchase. The warranty covers the repair or replacement according to Disetronic discretion, in case of a material or processing defect. No claim is admitted for regular use, overuse, improper use or care, chemical impact, force majeure and on the condition that all servicing of the remote control is done by Disetronic Customer Service.

In case of a damage caused by a third person or a none-authorized service centre the warranty expires.

Manufacturer:

Disetronic Medical Systems AG
Kirchbergstrasse 190
CH-3401 Burgdorf
Phone +41 34 424 24 24
Fax +41 34 424 24 00
www.disetronic.com

Distributed by:

Roche Diagnostics Nederland BV
Consumer Systems
Transistorstraat 41
NL-1322 CK Almere
Phone 0800 – 022 05 85
Fax +31 (0)36-53 94 874

SENSOTREND is a trademark
of a Member of the Roche Group.

04451309001/V01/03.04

Rx only

Swiss made



CE 0123



Diagnostics

D *DISETRONIC*