

# **Optotrak Application Programmer's Interface Guide**

**Revision 1.0  
December 2003**

**IMPORTANT**  
**Please read this entire document**  
**before attempting to operate**  
**the Measurement System**

**Part Number: IL-1070086**

Copyright 1992 - 2003 Northern Digital Inc. All Rights Reserved.

NDI and Optotrak are registered trademarks of Northern Digital Inc.

Measurement You Can Trust and Certus are trademarks of Northern Digital Inc.

🍁 Printed in Canada.

---

Published by:

Northern Digital Inc.  
103 Randall Dr.  
Waterloo, Ontario, Canada N2V 1C5

Telephone: + (519) 884-5142  
Toll Free: + (877) 634-6340  
Global: + (800) 634-634-00  
Facsimile: + (519) 884-5184  
Website: www.ndigital.com

Copyright 1992 - 2003, Northern Digital Inc.

All rights reserved. No part of this document may be reproduced, transcribed, transmitted, distributed, modified, merged, translated into any language or used in any form by any means - graphic, electronic, or mechanical, including but not limited to photocopying, recording, taping or information storage and retrieval systems - without the prior written consent of Northern Digital Inc. Certain copying of the software included herein is unlawful. Refer to your software license agreement for information respecting permitted copying.

### **Disclaimer of Warranties and Limitation of Liabilities**

Northern Digital Inc. has taken due care in preparing this document and the programs and data on the electronic media accompanying this document including research, development, and testing.

This document describes the state of Northern Digital Inc.'s knowledge respecting the subject matter herein at the time of its publication, and may not reflect its state of knowledge at all times in the future. Northern Digital Inc. has carefully reviewed this document for technical accuracy. If errors are suspected, the user should consult with Northern Digital Inc. prior to proceeding. Northern Digital Inc. makes no expressed or implied warranty of any kind with regard to these programs nor the supplemental documentation in this book.

Northern Digital Inc. makes no representation, condition or warranty to the user or any other party with respect to the adequacy of this document for any particular purpose or with respect to its adequacy to produce a particular result. The user's right to recover damages caused by fault or negligence on the part of Northern Digital Inc. shall be limited to the amount paid by the user to Northern Digital Inc. for the provision of this document. In no event shall Northern Digital Inc. be liable for special, collateral, incidental, direct, indirect or consequential damages, losses, costs, charges, claims, demands, or claim for lost profits, data, fees or expenses of any nature or kind.

NDI and Optotrak are registered trademarks of Northern Digital Inc. Measurement You Can Trust and Certus are trademarks of Northern Digital Inc.

Product names listed are trademarks of their respective manufacturers. Company names listed are trademarks or trade names of their respective companies.



---

# Table of Contents

<b>1 Symbols and Variables</b> .....	<b>1</b>
1.1 Symbols .....	1
1.2 Abbreviations and Acronyms .....	1
1.3 Variables .....	2
<b>2 How to Use the Optotrak API Guide</b> .....	<b>3</b>
2.1 Optotrak Sample Programs .....	4
2.2 Additional Optotrak Manuals .....	5
<b>3 System Overview</b> .....	<b>7</b>
3.1 Optotrak Certus System Components .....	7
3.2 Detection Region and Characterized Measurement Volume .....	11
3.3 Comparison of Optotrak Certus and Optotrak 3020 Systems .....	12
3.4 Mixed System Capability .....	13
3.5 Optotrak Communications .....	14
3.6 Initializing the Optotrak System .....	15
3.7 Data Conversions and Transformations on the Host Computer .....	17
3.8 Camera Parameter Files .....	18
3.9 Connecting Two Host Computers to the Optotrak System .....	19
3.10 API Quick Guide .....	20
<b>4 Optotrak Programmer's Guide</b> .....	<b>23</b>
4.1 Initializing, Retrieving System Status and Exiting from the Optotrak System .....	24
4.2 Retrieving Real-time Optotrak Data .....	29
4.3 Retrieving Buffered Optotrak Data .....	34

---

<b>5 ODAU Programmer's Guide</b> .....	<b>41</b>
5.1 Setting Up Data Collection from the ODAU .....	42
5.2 Retrieving ODAU Real-time Data .....	43
5.3 Retrieving Buffered ODAU Data .....	46
<b>6 Real-time Rigid Body Programmer's Guide</b> .....	<b>51</b>
6.1 Retrieving Real-time Rigid Body Data .....	51
6.2 Changing Rigid Body Settings .....	54
6.3 Changing the Rigid Body Coordinate System .....	57
6.4 Transforming Previously Obtained Data .....	61
6.5 Checking for Undetermined Transforms .....	63
<b>7 Floating Point Programmer's Guide</b> .....	<b>65</b>
7.1 The Northern Digital Floating Point Format .....	66
7.2 Converting Optotrak and ODAU Raw Data Files .....	73
7.3 Processing NDFP Format Files .....	76
<b>8 Retrieving Data With a Secondary Host Computer</b> .....	<b>83</b>
8.1 Retrieving Optotrak System Real-time Data on a Secondary Host Computer .....	84
8.2 Retrieving Buffered Data on a Secondary Host Computer .....	89
<b>9 Optotrak API Routines</b> .....	<b>93</b>
9.1 Overview .....	93
9.2 Optotrak Certus Specific Routines .....	94
9.3 Optotrak API Routines .....	96
9.4 Optotrak Specific Routines .....	101
9.5 Optotrak Device Handle Routines .....	135
9.6 ODAU Specific Routines .....	143

---

9.7	Real-time Data Retrieval Routines . . . . .	157
9.8	Buffered Data Retrieval Routines . . . . .	198
9.9	Rigid Body Specific Routines . . . . .	206
9.10	Rigid Body Related Routines . . . . .	214
9.11	File Processing Routines . . . . .	220
9.12	Registration and Alignment Routines . . . . .	231
<b>10</b>	<b>Real-time Data Types . . . . .</b>	<b>237</b>
10.1	“Missing” Marker Constants . . . . .	237
10.2	Optotrak Raw and Full Raw Data . . . . .	238
10.3	Optotrak 3D Data . . . . .	241
10.4	Optotrak Rigid Body Transformation Data . . . . .	242
10.5	ODAU Raw Data . . . . .	246
<b>Appendix A</b>	<b>Libraries and Sample Application Programs . . . . .</b>	<b>249</b>
A.1	API Installation CD . . . . .	249
A.2	Sample Programs for All Optotrak Systems . . . . .	249
A.3	Sample Programs for Optotrak Certus Systems . . . . .	264
<b>Appendix B</b>	<b>Error Messages and Constants . . . . .</b>	<b>273</b>
B.1	Error Constants . . . . .	273
B.2	Error Messages . . . . .	273
B.3	Message System Related Error Messages . . . . .	274
B.4	Transputer Related Error Messages . . . . .	277
B.5	Optotrak Related Error Messages . . . . .	277
B.6	ODAU Related Error Messages . . . . .	280
B.7	Real-time Related Error Messages . . . . .	282
B.8	Data Buffer Spooling Related Error Messages . . . . .	285
B.9	Rigid Body Related Error Messages . . . . .	287
B.10	File Processing Related Error Messages . . . . .	288

---

<b>Appendix C Flags and Settings Associated with Rigid Bodies</b> .....	<b>295</b>
C.1 Rigid Body Concepts and Terms .....	295
C.2 Accessing the Rigid Body with the API .....	295
C.3 Flags Affecting Rigid Bodies .....	297
C.4 Error Settings for Rigid Bodies .....	298



# 1 Symbols and Variables

## 1.1 Symbols

### Symbol



### Meaning

Follow the information in this paragraph to avoid either crashing the system or overwriting data.

## 1.2 Abbreviations and Acronyms

**Table 1-1: Abbreviations and Acronyms**

<b>Abbreviation or Acronym</b>	<b>Definition</b>
ANSI	American National Standards Institute
API	Application Programming Interface
IEEE	Institute (of) Electrical (and) Electronic Engineers
IREL	InfraRed light Emitting Diodes
ISA	Industry Standard Architecture
LED	Light Emitting Diode
NDFP	Northern Digital Floating Point
ODAU	Optotrak Data Acquisition Unit
PCI	Peripheral Component Interconnect
RISC	Reduced Instruction Set Computer
RMS	Root Mean Square
SCSI	Small Computer Systems Interface
SCU	System Control Unit
SIG	Silicon Graphics Inc.

## 1.3 Variables

**Table 1-2: Variables**

<b>Variable Symbol</b>	<b>Variable Definition</b>
M	number of markers
N	number of frames
S	number of sensors
C	number of channels
D	digital data
B	number of rigid bodies
$T_{x,y,z}$	translation, in mm, in the x, y, or z direction
$R_{x,y,z}$	rotation values, in radians, from a Euler transformation. $R_x = \text{Yaw}$ $R_y = \text{Pitch}$ $R_z = \text{Roll}$
$R_{00,\dots,22}$	values for a 3 x 3 rotation matrix transformation
r	centroid data
a	raw analog data from an ODAU device
v	voltage
q	vector values describing the orientation component of a quaternion transformation.

## 2 How to Use the Optotrak API Guide

This Application Programmer's Interface (API) guide explains the functions and routines you will need to write custom application software that can communicate with the Optotrak System. Users of this guide are expected to be familiar with the C/C++ programming language.

If you are unfamiliar with the configuration and functions of Optotrak Systems, refer to [“System Overview” on page 7](#) for a brief description of the devices and a model of the computer communications.

The next six chapters contain descriptions of the Optotrak System and how to write application programs for the system. The chapters include information on:

- the overview of the Optotrak System and how it communicates with the an application program
- the structure of an application program: the Program Initialization section, the Program Body Code section and the Program Exit Code section
- how to control and retrieve data from the Optotrak Data Acquisition Unit (ODAU)
- using the real-time Rigid Body extension
- manipulating data stored in a Floating Point File
- using a secondary host computer with the Optotrak System

Chapter 9 [“Optotrak API Routines” on page 93](#) is a comprehensive listing of the routines.

Chapter 10 [“Real-time Data Types” on page 237](#) describes each type of real-time data returned by the devices in the Optotrak System.

A description of the organization of the API CD included with this manual is located in Appendix A [“Libraries and Sample Application Programs” on page 249](#). This includes a description of the sample programs.

Error messages are listed in Appendix B [“Error Messages and Constants” on page 273](#).

Appendix D [“Flags and Settings Associated with Rigid Bodies” on page 295](#) contains information about the variables and error messages associated with rigid bodies.

## 2.1 Optotrak Sample Programs

Sample programs showing the concepts introduced in this guide are included on the API CD accompanying this manual. The CD contains sample programs that use the routine-based interface. The CD also contains the necessary NDI include files and Optotrak System interface libraries for use with application programs compiled with:

- Microsoft Windows 32-bit
- SGI, Linux and Sun workstations.

See [“Libraries and Sample Application Programs” on page 249](#) for a list of the sample programs and the operations they perform.

## 2.2 Additional Optotrak Manuals

Other manuals in the Optotrak System family include:

- Optotrak Certus User Guide
- Optotrak Certus System Guide
- Optotrak Certus Marker Strober User Guide
- Optotrak Certus Tool Strober Guide
- Optotrak Certus 3020 Strober Adapter Guide
- Optotrak Certus Axon Strober Guide
- Optotrak SCSI Interface Guide
- Optotrak Accelerated Processing Option Installation Guide \*
- Optotrak Data Acquisition Unit II Guide (ODAU II)
- Optotrak PCI Interface Card Intallation Guide
- Optotrak Tetherless Strober Controller User Guide \*
- Optotrak 16-Channel Terminal Strip Strober Wiring Guide \*
- Optotrak 24-Channel Strober Wiring Guide \*
- NDI ToolBench User Guide
- NDI 6D Architect User Guide
- NDI DataView User Guide
- Optotrak Certus Econo-Stand Guide
- Optotrak Certus Stand Adapter Guide

(\* for Optotrak 3020 Systems only)

To order these manuals, contact our technical support at:



**INTERNATIONAL HEADQUARTERS:  
NORTHERN DIGITAL INC.**

103 Randall Drive  
Waterloo, ON, Canada N2V 1C5

Phone: + (519) 884-5142  
Toll Free: + (877) 634-6340  
Global: + (800) 634-634-00  
Fax: + (519) 884-5184

Email: [optotrak@ndigital.com](mailto:optotrak@ndigital.com)  
Website: [www.ndigital.com](http://www.ndigital.com)

**EUROPEAN OFFICE:  
NDI EUROPE GmbH**

Fritz-Reichle-Ring 2  
D-78315 Radolfzell  
Germany

Phone: + 49 (77 32) 939 19 00  
Global: + (800) 634 634 00  
Fax: + 49 (77 32) 939 19 09

Email: [optotrak@ndieurope.com](mailto:optotrak@ndieurope.com)  
Website: [www.ndieurope.com](http://www.ndieurope.com)

**ASIA PACIFIC OFFICE:  
NDI ASIA PACIFIC**

Room 2603, 26th Floor  
Office Tower, Convention Plaza  
1 Harbour Road  
Wanchai, Hong Kong

Phone: + 852 2802 2205  
Fax: + 852 2802 0060

Email: [optotrak@ndigital.com](mailto:optotrak@ndigital.com)  
Website: [www.ndigital.com](http://www.ndigital.com)



## 3 System Overview

Read this section for an overview of the Optotrak System. Eight general areas are covered:

- devices that can be included in an Optotrak System
- the communication model for the Optotrak System
- initializing the Optotrak System
- comparisons of the Optotrak Certus System with the Optotrak 3020 System
- capability of a mixed system
- converting data on the host computer instead of the Optotrak System
- camera parameter files
- using two host computers with the Optotrak System

### 3.1 Optotrak Certus System Components

An Optotrak Certus System is designed to track diodes that emit infrared light within the measurement volume. The basic Optotrak Certus System, shown in [Figure 3-1 on page 9](#), consists of:

- a Position Sensor (more are optional)
- a System Control Unit
- one or more strobers, possibly connected with strober extension cables
- infrared light emitting diodes incorporated into a marker
- a host cable to connect the System Control Unit to the host computer
- a link cable to connect the System Control Unit to a Position Sensor
- power cables for both the System Control Unit and the Position Sensor
- a stand or bracket to support the Position Sensor, such as the econo-stand or Brunson stand (optional)
- an NDI PCI interface card (not shown)
- NDI ToolBench software for controlling of, and receiving data from, the Optotrak Certus System
- NDI 6D Architect software, for characterizing rigid bodies and tools.

- NDI Register and Align Wizard
- NDI File Convert Wizard
- NDI DataView software, which is used for viewing data files

There are two types of System Control Unit, strobers and some cables. The standard type is designed for systems where there is electrical isolation from the power mains. All of the strober ports in the standard type provide this type of electrical isolation. The second type, called e-type, is designed for applications that require type BF isolation on some strober ports. The e-type System Control Unit includes two strober ports with type BF isolation and one standard strober port. E-type strobers and cables are designed to maintain the type BF isolation. All e-type devices are labelled with the type BF symbol.


A host computer must be supplied by the user. If the host computer uses a Windows based operating system, an NDI PCI interface card can be installed in this computer to communicate with the Optotrak Certus System. Alternatively, communication can be established through either an Ethernet connection (if the user supplies both an Ethernet card and a cable) or a SCSI interface. To be fully compliant with the approvals listed in the “*System Guide*”, the host computer must also have IEC950 certification.

Additional software and accessories that can be used with the Optotrak Certus System include:

- Small Computer System Interface (SCSI) Interface Unit, an interface device alternative to the NDI PCI interface card and Ethernet options. If the host computer is using either Linux or one of the supported Unix operating systems, the System Control Unit must be connected to the host computer through either the NDI SCSI Interface Unit or an Ethernet card.
- Optotrak Data Acquisition Unit II, which can be used to accept and output both analog and digital signals.
- A vertical or multi-axis head adapter, which is used to attach the Position Sensor to the Brunson stand
- An econo-stand
- A Brunson stand
- A wall bracket
- Optotrak Application Programmer’s Interface, for users who want to write their own applications to control and obtain data from the Optotrak Certus System. This is required if you do not use NDI ToolBench software.



- Digitize software, which is used to create imaginary markers and to determine the coordinates of positions within the volume.



**Figure 3-1: Basic Optotrak Certus System**

The Position Sensor measures the location of the markers and sends the raw data to the System Control Unit. It includes a second link port to allow multiple Position Sensors to be daisy chained together. A total of eight Position Sensors can be used with one System Control Unit.

The System Control Unit receives the raw marker location data from the Position Sensor(s). The System Control Unit can calculate the 3D marker location from the raw data and forward the 3D data to the host computer. Alternatively, the System Control Unit can forward the raw data for conversion to 3D on the host computer.

3D marker positions are calculated in the default coordinate system of the Position Sensor. The default coordinate system is characterized at the factory and is specific to

each Position Sensor, as the exact location of the origin varies slightly with each unit. Systems consisting of multiple Position Sensors must first be configured by the user to register all the Position Sensors to a common coordinate system.

The System Control Unit also communicates with and controls the attached strobers as specified by software used on the host computer. The System Control Unit interfaces with external devices through a synchronization port and communicates with the host computer through either the link or Ethernet ports.

The host computer receives data from the System Control Unit and sends collection instructions, including strober setup, to the System Control Unit. The instructions are based on settings made in either NDI ToolBench software or through the Optotrak Application Programmer's Interface.

Strobers activate the markers in a sequence defined by software on the host computer, but controlled by the System Control Unit. They can communicate with the System Control Unit to inform the System Control Unit of their status and configuration.


Markers consist of an infrared light emitting diode fixed to a base. Markers can be tracked individually or as a group. If a group of markers is to be tracked, the positions of the markers relative to each other must be fixed; this is called a rigid body.

## 3.2 Detection Region and Characterized Measurement Volume


The Optotrak Certus System can track markers within one of two detection regions. The detection region is the volume defined by the field of view of the sensors on the Position Sensor. The Position Sensor can be focused to only one detection region. To determine if your Position Sensor is close or far focus, look at the label on the back of the Position Sensor.

The characterized measurement volume is a subset of the detection region. Within the characterized measurement volume, the accuracy of the measurement for a single marker is known. Outside the characterized measurement volume, the accuracy is unknown. See [Figure 3-2 on page 11](#) and [Figure 3-3 on page 12](#) for diagrams of the detection regions and the characterized measurement volumes. The characterized measurement volume and the detection region are the same for the close focus Position Sensor. Due to continuous product improvement, the dimensions in these diagrams are subject to change without notice.

Data is reported in the default coordinate system of the Position Sensor that is developed as part of the characterization process at NDI. The characterization process is used to develop the camera parameter file for the Position Sensor. This places the origin (0, 0, 0) on the middle sensor as shown in [Figure 3-2 on page 11](#) and [Figure 3-3 on page 12](#). The coordinate system can be aligned to place the origin and axes of the coordinate system in a location that is meaningful to you, as described in the “*Optotrak Certus User Guide*.”



**Figure 3-2 Close Focus Detection Region**



**Figure 3-3 Far Focus Detection Region**

### **3.3 Comparison of Optotrak Certus and Optotrak 3020 Systems**

The Optotrak Certus System improves performance, while retaining backwards compatibility to the Optotrak 3020.

Optotrak 3020 components and Optotrak Certus components have different internal processors. Optotrak 3020 components use transputers while Optotrak Certus components (Optotrak Certus Position Sensor and Optotrak Certus System Control Unit) use an improved processor called the SHARC processor. This processor can process information faster than transputers. All existing API routines that include transputer in their names such as `TransputerDetermineSystemCfg` and `TransputerInitializeSystem`, are generalized to handle the Optotrak Certus processor as well as Optotrak 3020 transputers. These routines will work in the same way in both the Optotrak 3020 System and the Optotrak Certus System.

The Optotrak Certus System is compatible with all Optotrak 3020 components and all operating systems supported by the Optotrak 3020. All previous file types are supported.

Optotrak Certus System devices can be used with Optotrak 3020 System devices in a mixed system.

**Table 3-1: Comparison of Optotrak Certus and Optotrak 3020**

Feature	Optotrak Certus	Optotrak 3020
Number of Markers	512	256
Sampling Rates (Marker Frequency)	4600 Hz	3500 Hz
Real-time Per Marker Data Rates (Hz) N = number of markers	$\frac{4600}{N + 2}$	$\frac{3500}{N + 1}$
Strober Ports	3	4
Maximum Number of Rigid Bodies (3 markers/rigid body)	170	85
Hot Connect/Disconnect	yes	no
Number of Position Sensors that can be used simultaneously	8 (6 if ODAU attached)	8
Host Computer Connection	Ethernet (10/100 mbps) SCSI PCI	PCI SCSI ISA

### 3.4 Mixed System Capability

Optotrak Certus and Optotrak 3020 Systems can be mixed. However, the capabilities of the mixed system is limited to the 3020 device included in the system. The Position Sensor and System Control Unit affect different settings:

**Table 3-2: Capabilities of a Mixed System**

Setting	Minimum	Maximum
Duty Cycle (%)	10	85
Voltage (V)	7.0	12.0
Marker Frequency (Hz)	1500	3500

**Table 3-2: Capabilities of a Mixed System**

<b>Setting</b>	<b>Minimum</b>	<b>Maximum</b>
Number of Markers	1	256

## **3.5 Optotrak Communications**

### **3.5.1 Client-Server Model**

The client-server model illustrates the interface between the Optotrak System and a user-written application program. Using this model, the Optotrak System acts as the server, and the host computer application programs are the clients. The services provided to these clients consist of supplying packets of information of a special type, such as 3D marker positions, rigid body positions, or the raw values produced by the Position Sensors. All clients generally communicate with either the Optotrak Administrator (in the Optotrak System) or the Optoscope Administrator (in the ODAU II) when making requests.

The client-server analogy can be carried further:

- To obtain information from the server, the client must first introduce itself to the server.
- After the introduction the client may make requests to the server.
- The server acknowledges that it understands the requests and attempts to satisfy them by completing each request. An example of this would be setting up a collection or returning data.
- Finally, when the client has finished using the provided services, communications with the server stop and the application program is exited.

### **3.5.2 Messages and Addresses**

The Optotrak System is a message-based system. All requests are made to the Optotrak System via messages to one of the processes (the Optotrak Administrator, the Optoscope Administrator, or the Data Proprietor) residing within the Optotrak System.

The message system can be compared to sending a letter through the mail. As with a letter, all messages contain a destination address and a source address. When you mail a letter, there is no need to know the physical location of the destination if you know the address. In the Optotrak System, the client does not need to know where the Optotrak Administrator is located; it could be on the System Control Unit or

somewhere on one of the Position Sensors. All that is required is the address of the Optotrak Administrator or the Data Proprietor and the message will be successfully relayed. The layer in the API responsible for the transmission of messages is called the message-passing layer.

A message identifier is attached to every message. This identifier determines the type of message being received and tells the receiver how to interpret the message.

### 3.5.3 Optotrak System Programming Interface

The programming interface defines how the application programs communicate with the Optotrak System. The application program must transmit information to and from the devices in the Optotrak System (the System Control Unit and the ODAU).

With the routine based interface, the application program can control devices in the Optotrak System by calling the appropriate routines with the proper parameters. Parameters are passed to the routines using C defined types (int, float, etc.), so the application program does not have to pack the structures. The routines handle the packing of structures, the delivery of appropriate message identification to the Optotrak System, and the receiving of the results from the Optotrak System. By examining the return value of a routine, the application program can determine if the routine was completed successfully.

---

**Note** Since the routine based interface frees the application program from having to use C structures, the library routines can easily be called from other programming languages. Other languages handle data structures, such as strings and arrays, differently than C, so it is essential that you understand the intricacies of calling the C library routines from another programming language.

---

## 3.6 Initializing the Optotrak System

When the Optotrak System is initialized, the current Optotrak System configuration is determined. Transputer programs are downloaded into all Optotrak 3020 components that are connected in the system. Predetermined startup code is downloaded into all Optotrak Certus components from the host computer to switch the components into run-time mode.

After the configuration has been determined, the results are recorded in a network information file (nif) and camera parameter file (cam). The default files are system.nif and standard.cam. They are written to the “realtime” subdirectory under the path specified in the ND\_USER\_DIR environment variable. If the environment variable ND\_USER\_DIR does not exist, then the files are written to the “realtime” subdirectory under the path specified in the ND\_DIR environment variable.

Network information files define the system configuration and the transputer programs that need to be downloaded to the Optotrak 3020 components if they are present in the Optotrak System. Camera parameter files define the operational characteristics for each Position Sensor in the system.



---

**You must re-determine the system configuration and re-initialize the Optotrak System if either the cabling connections to the communication ports at the back of the System Control Unit or the Position Sensor are changed, or if the order of the Position Sensors is changed. The cabling connection has changed once a cable is disconnected. It is a good practice to re-initialize your setup and re-determine the system configuration if you are not certain whether the cabling connections have been changed.**

---

Initialize the system with one of the following:

- application program, call API routines as part of its initialization
- command-line utility programs, before running the application program
- NDI ToolBench, select **Optotrak>Build Network Information File** in the Optotrak socket options window.

Calling the API routines in the application program is easier than using the command-line utility programs and is recommended for most applications.

### 3.6.1 Initializing the Optotrak System From an Application Program

Two steps are required to initialize the system from within the application program. The first step calls the routine `TransputerDetermineSystemCfg` and the second step calls `TransputerLoadSystem`.

`TransputerDetermineSystemCfg` generates the default system configuration files `system.nif` and `standard.cam`. The files are written to the “realtime” subdirectory under the path specified in the `ND_USER_DIR` environment variable. If the environment variable `ND_USER_DIR` does not exist, then the files are written to the “realtime” subdirectory under the path specified in the `ND_DIR` environment variable. `TransputerLoadSystem` downloads the appropriate transputer programs and startup code to the system according to the network information stored in `system.nif`.

If writing the system configurations to disk is undesirable, you can choose to use `TransputerDetermineSystemCfg` to store the system configuration information internally in the API. See “[TransputerDetermineSystemCfg](#)” on page 101.



### 3.6.2 Initializing the Optotrak System From the Command Line

Two steps are required to initialize the system from the command line. The first step determines the system configuration and is platform specific; use one of the following utility programs:

- `optset32.exe` for Windows NT/2000/XP
- `buildnif` for Linux, SGI and Sun systems

The second step downloads the system code to the processors. Use one of the following utility programs:

- `dld.exe` (Windows)
- `dld` (Linux, SGI and Sun).

The command line syntax for the DLD program is:

```
dld - <v>#<nif>
```

The optional argument `<v>` specifies the verbosity level that determines the amount of status information displayed by the DLD. The argument is a number ranging from 0 to 9 inclusively. At setting 0, no status information is provided.

The argument `'<nif>'` specifies the network information file. This name of this file is normally `system.nif`.

### 3.6.3 Initializing the Optotrak System During Program Runtime

After the transputer programs and startup code have been downloaded to the system, the API initiates communications with the system of processors with a call to the routine `TransputerInitializeSystem`. This reads the Optotrak System parameter initialization file `optotrak.ini`, located in the 'settings' subdirectory of the standard NDI system directory and initializes the system with the settings specified in the file.

## 3.7 Data Conversions and Transformations on the Host Computer

By default, the Optotrak System performs conversion from raw data to 3D data and all rigid body transformations internally. If your host computer contains a Pentium processor (or faster), floating point operations can be performed faster on the host computer than on the internal processors. The API can arrange to have all the conversions and transformations done on the host computer. This can be done either by setting appropriate flags with the routine `OptotrakSetProcessingFlags`, or by specifying the options in the initialization file `optotrak.ini`.

To specify the data conversion and rigid body transformation options in the initialization file `optotrak.ini`, add these two lines to the `[Optotrak System]` section:

```
bConvertOnHost = TRUE
bRigidOnHost   = TRUE
```

Setting both values to `FALSE` or removing the lines completely will cause all the conversions and transformations to be done on the Optotrak System.

To have the 3D data conversion done on the Optotrak System and the rigid body transformations on the host computer, use:

```
bConvertOnHost = FALSE
bRigidOnHost   = TRUE
```

`bRigidOnHost` can be set to `FALSE` only if `bConvertOnHost` is also set to `FALSE`. If `bConvertOnHost` is set to `TRUE`, the API will implicitly set `bRigidOnHost` to `TRUE`.

## 3.8 Camera Parameter Files

Camera parameters are written to the default camera parameter file called `standard.cam` during initialization and are stored internally in the Optotrak System. An API must load the camera parameters using `OptotrakLoadCameraParameters` if either 3D position data or rigid body data are to be collected or manipulated.

Camera parameter files define the operational characteristics for each Position Sensor in the system. The file contents are determined by calibration and registration procedures. If you have a single Position Sensor, the calibration parameters are specific to your unit and you do not need to re-register your unit. These parameters are in a file called `c3xxxxx.cam` where `xxxxxx` is the serial number of the Position Sensor. If you have multiple Position Sensors in your system you will need to generate a suitable camera parameter file. This can be accomplished by either performing a registration procedure using the routine `nOptotrakRegisterSystem`, or by using the Registration and Alignment wizard in the NDI ToolBench software.

### 3.8.1 Extended Camera Parameter Files

Extended camera parameter files have a different format but the same file extension (`.cam`) as standard camera parameter files. All API routines automatically distinguish between the two file formats.

Extended camera parameter files are comprised of one or more camera parameter sets associated with each Position Sensor. These files always contain the default camera parameter set corresponding to the default marker type, marker wavelength and original lens model. The extended camera parameter files may also contain parameter sets for additional lens models, marker types and marker wavelengths.

When an extended camera parameter file is loaded by the routine `OptotrakLoadCameraParameters`, the default camera parameter set is automatically selected. Select a different set by calling the routine `OptotrakSetCameraParameters`, having set the appropriate identifiers within the routine. For more information, please see [“OptotrakSetCameraParameters” on page 118](#). Call the routine `OptotrakGetCameraParameterStatus` to find a complete listing of all the most recently loaded camera parameter sets.

### 3.9 Connecting Two Host Computers to the Optotrak System

If you have purchased the additional hardware option that allows you to connect two host computers to the Optotrak System, then you can include two host computers in the system and run application programs from either computer. The physical connection of the second host computer depends on the devices connected to the Optotrak System. For complete instructions on connecting the second host computer, refer to the guide *“Installing the Secondary PC Interface Kit”*.

You must designate one computer as the primary host and the other computer as the secondary host. All of the sample programs provided in the API CD, except Sample Program 16, are designed to be used on the primary host. However, all of the sample programs can be modified to be used on the secondary host computer by making a few simple changes in the code (see [“Retrieving Optotrak System Real-time Data on a Secondary Host Computer” on page 84](#) and [“Sample Program 16” on page 263](#)). Sample Programs 16 is designed to be used on a secondary host computer.

---

**Note** The primary host initializes the system and loads the camera parameters. The secondary host must not repeat these procedures or the programs running on the primary host may display unpredictable behaviors.

---

An example application of a secondary host computer is to run an application program to retrieve and manipulate data while the primary host computer runs the NDI ToolBench software.

A different approach is for the primary host computer to run an application program that performs all the functions associated with setting up collections and retrieving data, while the application program on the secondary host is restricted to retrieving real-time and buffered data. You can also call an API routine to manipulate the data on the secondary computer.

---

**Note** Before you start developing your own custom application programs, try running the sample programs provided in the API CD in the primary and secondary modes on both computers. If you experience any problems when attempting to run the applications on both computers simultaneously, try using the NDI ToolBench software to ensure that your connections have been made correctly.

---

## 3.10 API Quick Guide

This section is a summary of the routines and procedures that are discussed in Chapters 4 through 9. You will need to go through the chapters to fully understand each routine.

### 3.10.1 Initializing the Optotrak System

#### Within the API

- TransputerDetermineSystemCfg, then TransputerLoadSystem
- TransputerInitializeSystem

#### From the Command Line

- optset32.exe or buildnif and dld, depending on your system

#### Using NDI ToolBench

- NDI ToolBench if your host computer is running Windows 9x/NT/2000/Me
- In NDI ToolBench, select **Optotrak>Build Network Information File** in the Optotrak socket options window

### 3.10.2 Basic Optotrak Routines

#### Initialize the System

After TransputerDetermineSystemCfg:

- TransputerLoadSystem
- TransputerInitializeSystem

#### Shut Down the System

- TransputerShutdownSystem

#### Retrieve Real-time Data

- OptotrakSetupCollection
- DataGetLatesttype (blocking method)

- RequestLatest*type*/DataIsReady/DataReceiveLatest*type* (non-blocking)

### Retrieve Buffered Data

- define spool mapping with DataBufferInitializeFile or DataBufferInitializeMem
- spool data with DataBufferSpoolData (blocking method)
- DataBufferStart/DataBufferWriteData (non-blocking)

## 3.10.3 Basic Optotrak Certus Specific Routines

### Device Handle Routines

- OptotrakGetNumberDeviceHandles
- OptotrakGetDeviceHandles
- OptotrakDeviceHandleEnable
- OptotrakDeviceHandleGetNumberProperties
- OptotrakDeviceHandleGetProperties

## 3.10.4 Basic ODAU Routines

---

**Note** Always call OptotrakSetupCollection after calling OdauSetupCollection.

---

### Real-time Data Retrieval

- OdauSetupCollection
- OptotrakSetupCollection
- DataGetLatestOdauRaw

### Retrieving Buffered Data

- OdauSetupCollection
- OptotrakSetupCollection
- DataBufferInitializeFile - for Optotrak System

- `DataBufferInitializeFile` - for ODAU
- `DataBufferSpoolData`

### 3.10.5 Rigid Body Routines

#### Retrieve Real-time Data

- `RigidBodyAdd` or `RigidBodyAddFromFile`
- `DataGetLatestTransforms` (blocking)
- `RequestLatestTransforms/DataIsReady/RequestLatestTransforms` (non-blocking)

#### Change the Transformation Settings

- `RigidBodyChangeSettings`

#### Determine Transformations from Previously Obtained 3D Data

- `OptotrakConvertTransforms`

### 3.10.6 File Handling

#### File Conversion

- `FileConvert` for Optotrak raw to 3D or ODAU raw to voltages

#### Processing Files

- `FileOpen/FileOpenAll` for all file types
- `FileRead/FileReadAll` for read only files
- `FileWrite/FileWriteAll` for read/write files
- `FileClose/FileCloseAll` for all file types

### 3.10.7 Using a Secondary Host

- primary host determines the system configuration and loads the processors
- primary host designates which computer is the secondary host with a flag in the routine `TransputerInitializeSystem`

## 4 Optotrak Programmer's Guide

In this section, you will find basic and advanced information on how to:

- initialize the application program, obtain system status and exit the application program
- obtain real-time data with blocking and non-blocking methods
- obtain buffered data using blocking and non-blocking methods

The API CD contains complete sample programs that can be compiled and run on a host computer. This section will refer to code fragments from these sample programs to illustrate the correct usage of the routines and the proper structure of application programs that can communicate with the Optotrak System effectively. The code fragments are for illustrating the concepts discussed in the documentation and may not reflect exactly the actual code in the sample programs.

Follow through the code fragments of the sample programs to learn how API routines function. Samples begin at the basic level and progress to more advanced topics. In most cases, the code is simplified and may not include error handling and storing data to hard disk and other details in order to focus on the major aspects of the Optotrak System interface.

The flow for each sample program is similar. Each sample program contains a Program Initialization Code section, a Program Body Code section, and a Program Exit Code section. The Program Initialization Code section and the Program Exit Code section are essentially the same for each of the sample programs. These sections are introduced in the first sample program. In other sample programs, they are only referred to when necessary and the main focus is on the code fragments comprising the Program Body Code section.

For a complete listing and description of the available routines, look in [“Optotrak API Routines” on page 93](#).

## 4.1 Initializing, Retrieving System Status and Exiting from the Optotrak System

To initialize the Optotrak System, retrieve system status information, and exit the application, three API code sections must be used: Program Initialization Code, Program Body Code and Program Exit Code.

### 4.1.1 Program Initialization Code

The Program Initialization Code section consists of three API routine calls: `TransputerLoadSystem`, `TransputerInitializeSystem`, and `OptotrakLoadCameraParameters`.

The `TransputerLoadSystem` routine loads the processors in the Optotrak System with the transputer program files outlined in the specified network information file (.nif). This NIF file defines the network connections between all of the Optotrak System components, and is generated by the `TransputerDetermineSystemConfiguration` routine.

---

**Note** After a call to the routine `TransputerLoadSystem`, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may be generated.

---

The routine `TransputerInitializeSystem` initializes a message-passing layer, which allows the application program to send information to and receive information from the Optotrak System.

`OptotrakLoadCameraParameters` loads the camera parameters from the specified camera parameter file and sends the camera parameters to the Optotrak System. The default camera parameter file is called `standard.cam`.

---

**Note** All application programs that calculate 3D or 6D data must call `OptotrakLoadCameraParameters`.

---

### 4.1.2 Program Body Code

In the Program Body Code section the application program sets up the collection parameters and retrieves data from the Optotrak System. For example, the code in [Figure 4-1 on page 28](#) obtains the current Optotrak System status from the Optotrak Administrator process. The status information returned describes the collection parameters (e.g. `fFrameFrequency`, `nMarkers`) and other process entities in the current system configuration (e.g. `nNumOdaus`, `nNumSensors`). Use the status information to determine the sizes of real-time data packets that are returned.



### 4.1.3 Program Exit Code

The Program Exit Code section invokes the routine `TransputerShutdownSystem`. This routine removes the host computer from the message-passing layer.

### 4.1.4 Program Sample Showing How to Initialize and Retrieve System Status

Figure 4-1 on page 28 is an example of how to initialize the Optotrak System and retrieve the system status information. This program is similar, but not identical, to Sample Program 1 on the API CD.

In this sample, the network information file `system.nif` is assumed to have been created with the routine `TransputerDetermineSystemCfg`. Alternatively, you can initialize the Optotrak System using the NDI ToolBench software, or by using a command-line utility program. For more information, please see “[Initializing the Optotrak System](#)” on page 15.

---

**Note** If the Optotrak System was previously initialized by another program the transputer program files do not need to be reloaded.

---

In this sample the results obtained in the body code are output to the display screen.

To initialize the Optotrak System, retrieve the system status and exit the program, follow these steps:

1. Load the Optotrak System with the appropriate transputer programs and startup code.
2. Connect the application program to the message-passing layer.
3. Load the Optotrak System with the appropriate camera parameters.
4. Request and display the current Optotrak System status.
5. Disconnect the application program from the message-passing layer.

```

/*****
Name:                SAMPLE1.C

```

```
Description:
```

```
    Optotrak Sample Program #1.
```

- ```

1. Load the system of processors with the appropriate
   transputer programs and startup code.
2. Initiate communications with the system of processors.
3. Load the appropriate camera parameters.
4. Request/receive/display the current Optotrak System status.

```

Pass NULL for those status variables that are not requested.

5. Disconnect the PC application program from the system of processors.

```
*****/
/*****
C Library Files Included
*****/
#include <stdio.h>
#include <string.h>
#include <stdlib.h>

#ifdef _MSC_VER
void sleep( unsigned int uSec );
#elif __BORLANDC__
#include <dos.h>
#elif __WATCOMC__
#include <dos.h>
#endif

/*****
ND Library Files Included
*****/
#include "ndtypes.h"
#include "ndpack.h"
#include "ndopto.h"
/*****
Name:          main

Input Values:
    int
        argc          :Number of command line parameters.
    unsigned char
        *argv[]       :Pointer array to each parameter.

Output Values:
    None.

Return Value:
    None.

Description:
The main program routine performs all steps listed in the above program
description.
*****/

void main( in argc, unsigned char *argv[])
{
    int
        nNumSensors,
        nNumOdaus,
        nMarkers;
```

```
char
    szNDErrorString[MAX_ERROR_STRING_LENGTH + 1];

/*
 * STEP 1
 * Load the system of processors.
 */
if( TransputerLoadSystem( "system" ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 *Wait one second to let system finish loading
 */
sleep(1);

/*
 * STEP 2
 * Initialize the system of processors.
 */
if( TransputerInitializeSystem( OPTO_LOG_ERRORS_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 3
 * Load the standard camera parameters.
 */
if( OptotrakLoadCameraParameters( "standard" ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 4
 * Request and receive the Optotrak status.
 */
if( OptotrakGetStatus(
    &nNumSensors, /* Number of sensors in the Optotrak System. */
    &nNumOdaus,   /* Number of ODAUs in the Optotrak System. */
    NULL,        /* Number of rigid bodies being tracked by the O/T.*/
    &nMarkers,    /* Number of markers in the collection. */
    NULL,        /* Frequency that data frames are being collected.*/
    NULL,        /* Marker firing frequency. */
    NULL,        /* Dynamic or Static Threshold value being used.*/
    NULL,        /* Minimum gain code amplification being used.*/
    NULL,        /* Stream mode indication for the data buffers*/
    NULL,        /* Marker Duty Cycle being used.*/
```

```
        NULL,          /* Voltage being used when turning on markers */
        NULL,          /* Number of seconds data is being collected.*/
        NULL,          /* Number of seconds data is being pre-triggered.*/
        NULL ) )      /* Configuration flags.*/
{
    goto ERROR_EXIT;
} /* if */

/*
 * Display elements of the status received.
 */

    fprintf( stdout, "Sensors in system      :%3d\n", nNumSensors );
    fprintf( stdout, "ODAU's in system      :%3d\n", nNumOdaus );
    fprintf( stdout, "Default Optotrak Markers:%3d\n", nMarkers );

/*
 * STEP 5
 * Shutdown the message-passing system.
 */
if( TransputerShutdownSystem() )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Exit the program.
 */

fprintf(stdout, "\nProgram execution complete.\n");
exit( 0 );

ERROR_EXIT:
    if(Optotrak Get Error String(szNDErrorString,
                                MAX_ERROR_STRING_LENGTH + 1 ) == 0 )
    {
        fprintf( stdout, szNDErrorString );
    } /* if */
    TransputerShutdownSystem();
    exit( 1 );
} /* main */
```

**Figure 4-1: Retrieving Optotrak System Status**

## 4.2 Retrieving Real-time Optotrak Data

The primary function of the Optotrak System is to collect position data and relay it to the host computer. Data can be retrieved by an application program in real-time, or from the data buffer. When real-time data is requested, the Optotrak System returns the most recent frame of data to the application program. Retrieving real-time data is discussed before retrieving buffered data (see [“Retrieving Buffered Optotrak Data”](#) on page 34).

There are three types of data that can be retrieved from an Optotrak System: raw, full raw and 3D data. Full raw data is raw (centroid) data with additional status information including signal strength, amplification, and error code. Full raw data files are twice as large as raw data files and are used primarily to diagnose problems in the equipment setup. 3D data cannot be collected as quickly as raw or full raw data because of the time required to convert the data. You may define the file names for the data files. In general, filenames for raw data files begin with R# while filenames for 3D data files begin with C#.

---

**Note** Centroid data is *raw* data. Routines that access *raw* data have ‘centroid’ in their name (e.g. `DataGetLatestCentroid`). *Full raw* data is centroid data with additional status information. Routines that access *full raw* data have ‘raw’ in their name (e.g. `DataGetNextRaw`).

---

You must configure the Optotrak System before data can be collected. This sets the values for the collection parameters such as the number of markers and frame rate. See [“OptotrakSetupCollection”](#) on page 124 for details on the collection parameters.

Raw and 3D data can be collected in real-time or as buffered data.

### 4.2.1 Basic Real-time Data Retrieval

There are several ways to retrieve real-time data. The sample below illustrates the simplest, single routine invocation for retrieving data in real-time. An advanced method is described in the next section, [“Advanced Real-time Data Retrieval Without Blocking”](#) on page 32. Once the data has been received from the Optotrak System, it can be displayed or processed with the application program. Refer to [“Real-time Data Retrieval Routines”](#) on page 157 for a more detailed explanation of real-time data retrieval.

Before the application program exits, ensure that the IRED markers are de-activated. This operation can be considered to be part of the Program Exit Code section described in the [“Program Exit Code”](#) on page 25. Hence, it does not appear in any further discussions or code fragments.



**You must ensure that enough memory is allocated to accommodate a frame of real-time data to prevent the possibility of invalid memory regions being overwritten and possible system crashes. The routines for data retrieval copy the data directly into the memory allocated by the host application program. Each of the real-time data retrieval routine descriptions describes how to determine the required amount of memory for its associated data type. See [“Real-time Data Retrieval Routines”](#) on page 157 for the routines. Refer to [“Retrieving Optotrak System Real-time Data on a Secondary Host Computer”](#) on page 84 for sample code that includes memory allocation.**

## 4.2.2 Sample Program to Retrieve and Display Real-time 3D Marker Positions

[Figure 4-2 on page 32](#) shows an example of how to retrieve and display real-time 3D marker positions using the simplest method, the routine `DataGetLatest3D`. This routine blocks while the Optotrak System calculates the latest frame of 3D coordinates and returns once the data has been received. This example is similar, but not identical, to Sample Program 2 on the API CD.

This sample program retrieves real-time 3D data via a routine call and displays the data on the screen. In addition to returning the 3D data, the Optotrak System returns three other related data items: `uFrameNumber`, `uElements`, and `uFlags`. For a description of each of these data items, refer to [“DataGetLatestCentroid”](#) on page 157.

To retrieve and display real-time data from the Optotrak System:

1. Configure the Optotrak System collection settings.
2. Activate the IRED markers.
3. Retrieve the latest frame of 3D marker data.
4. Display the 3D data and the status information.
5. De-activate the IRED markers.

```

/*
 * STEP 1
 * Set up a collection for the Optotrak.
 */
if( OptotrakSetupCollection(
    NUM_MARKERS, /* Number of markers in the collection. */
    (float)100.0, /* Frequency to collect data frames at. */
    (float)2500.0, /* Marker frequency for marker maximum on-time. */
    30, /* Dynamic or Static Threshold value to use. */
    160, /* Minimum gain code amplification to use. */
    0, /* Stream mode for the data buffers. */
    (float)0.35, /* Marker Duty Cycle to use. */
    (float)7.0, /* Voltage to use when turning on markers. */

```

```

        (float)1.0,      /* Number of seconds of data to collect. */
        (float)0.0,    /* Number of seconds to pre-trigger data by. */
        OPTOTRAK_BUFFER_RAW_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 2
 * Activate the markers.
 */
if( OptotrakActivateMarkers() )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Get and display ten frames of 3D data.
 */

fprintf( stdout, "\n\n3D Data Display\n" );
for( uFrameCnt = 0; uFrameCnt < 10; ++uFrameCnt )
{
    /*
     * STEP 3
     * Get a frame of data.
     */
    fprintf( stdout, "\n" );
    if( DataGetLatest3D( &uFrameNumber,
                        &uElements,
                        &uFlags,
                        p3dData ) )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
     * STEP 4
     */

    fprintf( stdout, "Frame Number: %8u\n", uFrameNumber );
    fprintf( stdout, "Elements      : %8u\n", uElements );
    fprintf( stdout, "Flags          : 0x%04x\n", uFlags );
    for( uMarkerCnt = 0; uMarkerCnt < NUM_MARKERS; ++uMarkerCnt )
    {
        fprintf( stdout, "Marker %u X %f Y %f Z %f\n",
                uMarkerCnt + 1,
                p3dData[ uMarkerCnt ].x,
                p3dData[ uMarkerCnt ].y,
                p3dData[ uMarkerCnt ].z );
    }
}

```

```
        } /* for */  
    } /* for */  
  
/*  
 * STEP 5  
 * De-activate the markers.  
 */  
  
if( OptotrakDeActivateMarkers() )  
{  
    goto ERROR_EXIT;  
} /* if */
```

**Figure 4-2: Retrieving and Displaying Real-time 3D Marker Positions**

### 4.2.3 Advanced Real-time Data Retrieval Without Blocking

The Optotrak API also provides a method where the host computer does not block while waiting for the Optotrak System to return the data. This allows the host computer to perform other operations while the Optotrak System is determining the latest frame of data. For example, the application program can process input from another device, perform calculations, or control other devices.

The non-blocking method uses three routines: *RequestLatesttype*, *DataIsReady*, and *DataReceiveLatesttype*.

**RequestLatesttype** requests that the Optotrak System compute the latest frame of data and to send it back to the host computer. This routine does not block and wait for the data to be returned. The application program must determine when the real-time data is ready and then receive it accordingly.

**DataIsReady** is invoked to determine if the real-time data is ready to be received. If the real-time data is waiting to be received, *DataIsReady* returns a non-zero value. The application program should then promptly receive the data.

**DataReceiveLatesttype** must be used to retrieve the data that is waiting to be received. *DataReceiveLatesttype* copies the data waiting to be received directly into the pre-allocated memory location on the host computer.





**Do not send a new request for data until the data from the previous request has been received. This is especially true if the application program is requesting two types of data. Suppose an application was requesting both ODAU II raw data and Position Sensor 3D data, and the requests for the data were made immediately following one another. Once the routine `DataIsReady` returns `TRUE`, the application program has no way of determining whether it was receiving ODAU or Position Sensor data.**

## 4.2.4 Sample Program to Collect Real-time Data Without Blocking

This sample program is an example of how to collect real-time 3D data using the non-blocking method:

1. Request the latest 3D data from the Optotrak System.
2. Wait for the data to be returned from the Optotrak System.
3. Receive the data.

```

/*
 * STEP 1
 * Request the latest 3D data.
 */

if( RequestLatest3D() )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 2
 * Loop until the data is ready to be received.
 */

while( !DataIsReady() )
{
    ; /* We could do other processing here while we're waiting. */
}
/*
 * STEP 3
 * Get the latest data that is waiting to be received.
 */
if( DataReceiveLatest3D( &uFrameNumber, &uElements, &uFlags, p3dData ) )
{
    goto ERROR_EXIT;
} /* if */

```

**Figure 4-3: Advanced Real-time Data Retrieval**

## 4.3 Retrieving Buffered Optotrak Data

Once a collection has been configured, the Optotrak System collects data according to the specified collection parameters and stores each frame of data in a circular buffer on the SCU. An application program can retrieve the data stored in this buffer in fixed-size blocks, and either writes the data to disk, or stores it in the memory previously allocated by the application program.

---

**Note** You must configure the collection and activate the markers before buffered data can be retrieved. Several of the collection parameters affect the amount of data that is returned. The program can also specify whether raw (centroid) data, full raw data or 3D data is to be buffered. For complete details, refer to [“OptotrakSetupCollection” on page 124](#).

---

Before any data can be retrieved from the data buffer, a spool mapping must be made between the data buffer on the SCU and the spool destination for the data on the host computer. A spool mapping defines the destination (e.g. a data file or an allocated memory block) for the buffered data from a particular source. The `DataBufferInitializeFile` routine is invoked to define the spool mapping between the data buffer and the specified file. The routine `DataBufferInitializeMem` is invoked to define the spool mapping between the data buffer and the specified memory block. Once this relationship is defined, any buffered data the host computer receives from the Optotrak System is written directly to the file or memory block.

There are two methods to spool the data. The routine `DataBufferSpoolData` does not return until all buffered data has been received from the Optotrak System. Some application programs may need to do other operations while the buffered data is being spooled. To use an advanced, non-blocking procedure for these programs, refer to [“Advanced Buffered Data Retrieval Without Blocking” on page 36](#).

### 4.3.1 Basic Buffered Data Retrieval

The simplest method of data spooling uses the `DataBufferSpoolData` routine. This routine instructs the Optotrak System to start sending buffered data back to the host computer starting at the current frame. The host computer receives buffered data and writes them directly to the appropriate destination; in the sample below, this is the file “C#001.SM2.” After all the buffered data has been written, the file is closed, the appropriate header information is written and the status of the spooling procedure is returned. If the spool status returned is not zero, then a buffering error has occurred during the spooling procedure, and the spooled data may not be valid.

---

**Note** Once the spool procedure has finished, the previously defined spool mappings are no longer valid. All spool mappings must be reset before performing the next spool procedure.

---

### 4.3.2 Basic Sample Program to Collect Buffered Data to Disk

In this sample program, buffered data is written to the file “C#001.SM2”. This code fragment is similar to Sample Program 3 in the API CD. To collect Optotrak System buffered data to disk:

1. Configure the Optotrak System collection.
2. Activate the IRED markers.
3. Define the spool mapping between the data buffer and the file “C#001.SM2.”
4. Spool the data and print the spool status.

```

/*
*STEP 1
* Set up a collection for the Optotrak.
*/
if( OptotrakSetupCollection(
    6,                /* Number of markers in the collection. */
    (float)50.0,      /* Frequency to collect data frames at. */
    (float)2500.0,    /* Marker frequency for marker maximum on-time.*/
    30,               /* Dynamic or Static Threshold value to use. */
    160,              /* Minimum gain code amplification to use. */
    0,                /* Stream mode for the data buffers. */
    (float)6.5,       /* Voltage to use when turning on markers. */
    (float)2.0,       /* Number of seconds of data to collect. */
    (float)0.0,       /* Number of seconds to pre-trigger data by. */
    OPTOTRAK_FULL_DATA_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
* STEP 2
* Activate the markers.
*/

if( OptotrakActivateMarkers() )
{
    goto ERROR_EXIT;
} /* if */

/*
* STEP 3
* Initialize a data file for spooling of the Optotrak data.
*/

if( DataBufferInitializeFile( OPTOTRAK, "C#001.SM2" ) )

```

```
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 4
 * Spool data to the previously initialized file.
 */
if( DataBufferSpoolData( &uSpoolStatus ) )
{
    goto ERROR_EXIT;
} /* if */
fprintf( stdout, "Spool Status: 0x%04x\n", uSpoolStatus );
```

**Figure 4-4: Collecting Optotrak Buffered Data to Disk**

### **4.3.3 Advanced Buffered Data Retrieval Without Blocking**

Some application programs may need to do other operations while the buffered data is being spooled. These programs use a non-blocking spool procedure.

Three API routines are used to implement the non-blocking spooling process while collecting 3D data: `DataBufferStart`, `DataBufferWriteData`, and `DataBufferStop`. Once all spool mappings are initialized, the data buffers on the SCU and the ODAU (if connected) start spooling data back to the host computer. This is done by invoking the routine `DataBufferStart`.

While the Optotrak System is spooling buffered data back to the host computer, the application program must repeatedly call the routine `DataBufferWriteData`. This routine receives any buffered data from the Optotrak System and writes it to the appropriate spool destination. This routine also sets the following status variables to determine the current status of the Optotrak System: `uRealtimeDataReady`, `uSpoolComplete`, `uSpoolStatus`, and `ulFramesBuffered`.

**uRealtimeDataReady** is set to a non-zero value if there is real-time data to be received by the application program.

**uSpoolComplete** is set to a non-zero value when the spooling procedure is finished. Once `uSpoolComplete` is set to non-zero, the spooling procedure is done, and the application program need not make further calls to `DataBufferWriteData`.

**uSpoolStatus** is zero unless a buffering error has occurred, and the spooled data may not be valid.

**ulFramesBuffered** is set to the number of frames of data returned from the Optotrak System in the collection at the point when this routine is called. If this information is not of interest, set this parameter to `NULL`.

Use the `DataBufferStop` routine to stop the spooling of buffered data before all the data has been sent by the Optotrak System. A call to this routine stops the Optotrak

System from sending buffered data at the current frame being collected. However, there may still be a block of buffered data sent by the Optotrak System: the application program must stay in the loop until the `uSpoolComplete` parameter is set to be non-zero by the routine `DataBufferWriteData`.

#### 4.3.4 Advanced Sample Program to Retrieve Buffered Data Without Blocking

Figure 4-5 on page 39 is an example of the non-blocking spooling procedure. This code fragment is similar to Sample Program 7 on the API CD. The Optotrak System is initialized, a collection is configured and IRED markers are activated, as in previous examples. This sample then waits for marker 1 to come into view. At this point, a data file collection is started and data is collected until all the buffered data has either been spooled or marker 1 goes out of view. Program Exit Code is run before returning back to the operating system.

To wait for a marker to come into view, a pause is created by requesting real-time 3D data in a loop. During each iteration through the loop, the data is examined to determine if a 3D coordinate could be calculated for marker 1. If a position for marker 1 cannot be calculated, the X, Y, and Z coordinates of the 3D position are set to a very large negative number (`BAD_FLOAT`) by the Optotrak System and the program compares the returned values with the constant `MAX_NEGATIVE`. In this sample, once the X value becomes greater than `MAX_NEGATIVE`, the program determines that a valid 3D coordinate has been obtained, and continues with the data file collection.

Similarly, if marker 1 is out of view, the program repeatedly retrieves 3D data and compares the X coordinate to the value `MAX_NEGATIVE`. If the X value is less than `MAX_NEGATIVE`, a 3D number cannot be calculated, the marker is assumed to be out of view and the spooling procedure is terminated.

To collect buffered data to disk:

1. Initialize spooling variables that are used later by the routine `DataBufferWriteData`.
2. Wait for marker 1 to come into view.
3. Start the Optotrak System spooling buffered data to the host computer.
4. Get a frame of 3D marker data.
5. If marker 1 is out of view, stop the Optotrak System spooling buffered data.
6. Write data to the appropriate spool destination and check for spooling completion.

```
/*
 * STEP 1
 * Initialize the necessary spooling variables and a file for
```

```
* spooling of the Optotrak data.
*/
uSpoolStatus      =
uSpoolComplete    =
uRealtimeDataReady = 0;

/*
* STEP 2
* Loop until marker 1 comes into view.
*/

fprintf( stdout, "Waiting for Marker 1...\n" );
do
{
    /*
    * Get a frame of 3D data.
    */

    if( DataGetLatest3D( &uFrameNumber,&uElements,&uFlags,p3dData ) )
    {
        goto ERROR_EXIT;
    }
    /* if */
} /* do */
while( p3dData[ 0].x < MAX_NEGATIVE );

/*
* STEP 3
* Start the Optotrak spooling data to us.
*/

if( DataBufferStart() )
{
    goto ERROR_EXIT;
} /* if */

fprintf( stdout, "Collecting Data File\n" );

/*
* Loop around spooling data to the file until there is no data left.
*/
do
{
    /*
    * STEP 4
    * Get a frame of 3D data.
    */
```

```
if( DataGetLatest3D( &uFrameNumber, &uElements, &uFlags, p3dData ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 5
 * Check to see if marker 1 is out of view and stop the Optotrak
 * from spooling data if this is the case.
 */
if( p3dData[ 0].x < MAX_NEGATIVE )
{
    if( DataBufferStop() )
    {
        goto ERROR_EXIT;
    } /* if */
} /* if */

/*
 * STEP 6
 * Write data if there is any to write.
 */

if( DataBufferWriteData( &uRealtimeDataReady,
                        &uSpoolComplete,
                        &uSpoolStatus,
                        NULL ) )
{
    goto ERROR_EXIT;
} /* if */
} /* do */
while( !uSpoolComplete );
fprintf( stdout, "Spool Status: 0x%04x\n", uSpoolStatus );
```

**Figure 4-5: Advanced Buffered Data Retrieval**





## 5 ODAU Programmer's Guide

This section will help you to understand the routines used to control and retrieve data from an Optotrak Data Acquisition Unit (ODAU). You are assumed to have read and understood all the concepts presented in the previous section, the “[Optotrak Programmer's Guide](#)” on page 23. An application program can retrieve both real-time and buffered data from an ODAU device. Use the sample code fragments in this section to learn how to retrieve both of these data types.

There are two types of data acquisition units: the ODAU I has been superseded by the ODAU II. The ODAU II includes all of the ODAU I functions, but supports higher data rates and has additional features. The multiplexer, analog to digital input/output and timer functions have all been enhanced in the ODAU II. Some routines in this section apply only to the ODAU II. The introduction of each routine in “[ODAU Specific Routines](#)” on page 143 states which ODAU unit the routine is applicable to.

When the Optotrak System includes an ODAU, it is capable of collecting synchronized analog and digital data in addition to centroid (raw), full raw, 3D and 6D transformation data. The ODAU II has 16 analog voltage input channels and 8 digital input/output pins, which are collected synchronously with the Optotrak System data. See [Table 9-11 on page 154](#) for an ODAU I digital I/O port pin configuration or [Table 9-13 on page 156](#) for an ODAU II digital I/O port pin configuration. An application program cannot be written for the ODAU alone because the ODAU connects to and is controlled by the Optotrak System Control Unit. ODAU devices must be connected to a complete Optotrak System.

## 5.1 Setting Up Data Collection from the ODAU

You must set up a valid collection before any data can be retrieved from the ODAU; all ODAU devices must have their collections set up before the Optotrak System's. The order in which these collections are set up is important. After the Optotrak System receives its collection parameters, it signals all ODAU devices to start sampling data. This ensures that all devices start collecting data at the same point in time. The application program can begin to request data after all devices have been successfully configured.

Use the routine `OdaSetupCollection` to configure the ODAU collection parameters. Up to four (4) ODAU devices can be connected to the Optotrak System in a series. A separate call to the `OdaSetupCollection` routine must be made for each ODAU connected in the Optotrak System.

The routine `OdaSetupCollection` uses the first parameter to identify which set of collection parameters are for which ODAU device. There are four pre-defined constants used to identify the ODAU devices: `ODAU1`, `ODAU2`, `ODAU3`, and `ODAU4`. The constants identify the ODAU devices in the following manner: `ODAU1` is the first ODAU device connected to the Optotrak System (i.e. connected to the SCU). `ODAU2`, `ODAU3`, and `ODAU4` are the second, third, and fourth ODAU units in the serial connection. For a complete description of the parameters for `OdaSetupCollection`, refer to "[OdaSetupCollection](#)" on page 149.

---

**Note** After a call to the routine `OptotrakSetupCollection` or `OdaSetupCollection`, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may generated.

---

## 5.2 Retrieving ODAU Real-time Data

This section explains how to retrieve real-time data from an ODAU device connected to the Optotrak System.

Use the routine `DataGetLatestOdauraw` to retrieve real-time data from the ODAU device. This routine retrieves the latest complete frame of raw data from the specified ODAU device. Since the ODAU can collect both analog and digital data, the data frame returned from the ODAU may contain analog data, digital data or both, depending on the collection setup. If both analog and digital data are being sampled, the analog data always precedes the digital data in the data frame. See “[Sample Formats](#)” on page 68 and “[ODAU Raw Data](#)” on page 246 for a complete description of a data frame.

The ODAU device returns analog data only in raw format. If the application requires the data to be expressed as a voltage, the conversion must be done on the host computer.

The formula for calculating voltages for an ODAU I, corrected for offset and gain from the given analog raw data is:

$$Voltage = \frac{RawData - 2048}{204.8 \times Gain}$$

where *RawData* is the analog raw data for a particular ODAU channel and *Gain* is the analog gain used in the current collection setup for that ODAU.

For an ODAU II device, the voltage is calculated using:

$$Voltage = \frac{RawData \times 0.000305175}{Gain}$$

where *RawData* is the analog raw data for a particular ODAU channel and *Gain* is the analog gain used in the current collection setup for that ODAU.

In [Figure 5-1 on page 45](#), the formula for an ODAU II device is used to convert the raw analog data to a voltage, which is then displayed. The upper 8 bits contain the digital data.

For an ODAU I device, the code fragment for calculating the voltage is replaced with:

```
((float) ((int) (puRawData[ uChannelCnt] & 0x0FFF) - 2048) / 204.8);
```

The formula is implemented in the code fragment with one modification – a logical “AND” is performed with the raw data and the hexadecimal constant 0FFF. This operation is added to remove the upper 4 bits of the 16-bit raw data value. Since the analog to digital converter on the ODAU has 12-bit resolution, only the lower 12 bits are required. The upper 4 bits can contain status information, which, if used in the formula, would cause the wrong voltage to be computed.

When all the analog data has been displayed, the digital data is output as hexadecimal values. The digital data is in a 16-bit integer format, where each bit represents the status of the corresponding channel on the digital port. If a bit is 1, then the corresponding channel was at a logic high when the frame was sampled. If a bit is 0, the corresponding channel was at a logic low when the frame was sampled.

## Program Sample of Retrieving Real-time Data from ODAU

Figure 5-1 on page 45 contains the code that shows how to retrieve and display ODAU real-time data. The code fragment assumes that a valid network information file has been set up (“[Initializing the Optotrak System](#)” on page 15) and that the standard Program Initialization Code section has been run successfully. This code fragment is similar to Sample Program 12 on the API CD.

To retrieve ODAU real-time data:

1. Configure the ODAU collection.
2. Configure the Optotrak System collection.
3. Activate the IRED markers.
4. Convert the raw data to voltage and display the voltages.
5. Display the ODAU digital data in hexadecimal format.

```

/*
 * STEP 1
 * Set up a collection for the ODAU.
 */
if( OdauSetupCollection(
    ODAU1,                /*Id the ODAU parameters are for.*/
    NUM_CHANNELS,        /*Number of analog channels to collect.*/
    ODAU_GAIN,           /*Gain to use for the analog channels.*/
    ODAU_DIGITAL_INPB_INPA, /*Mode for the Digital I/O port.*/
    (float)100.0,        /*Frequency to collect data frames at.*/
    (float)90000.0,      /*Frequency to scan channels at. */
    0,                   /*Stream mode for the data buffers. */
    (float)1.0,          /* Number of seconds of data to collect.*/
    (float)0.0,          /* Number of seconds to pretrigger data.*/
    0 ) )
    goto ERROR_EXIT;
} /* if */
/*
 * STEP 2
 * Set up a collection for the OPTOTRAK.
 */
if( OptotrakSetupCollection(
    NUM_MARKERS,         /* Number of markers in the collection. */
    (float)100.0,        /* Frequency to collect data frames at. */

```

```

        (float)2500.0, /* Marker frequency for marker maximum on-time. */
        30,           /* Dynamic or Static Threshold value to use. */
        160,         /* Minimum gain code amplification to use. */
        0,           /* Stream mode for the data buffers. */
        (float)0.20, /* Marker Duty Cycle to use. */
        (float)6.0,  /* Voltage to use when turning on markers. */
        (float)1.0,  /* Number of seconds of data to collect. */
        (float)0.0,  /* Number of seconds to pre-trigger data by. */
        OPTOTRAK_BUFFER_RAW_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */
/*
 * Get and display ten frames of ODAU data.
 */
fprintf( stdout, "ODAU Data Display\n" );
for( uFrameCnt = 0; uFrameCnt < 10; ++uFrameCnt )
{
    /*
     * STEP 3
     * Get a frame of ODAU raw data.
     */
    if( DataGetLatestOdauRaw( ODAU1, &uFrameNumber, &uElements,
                             &uFlags,pRawData ) )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
     * STEP 4
     * Print out the data.
     */
    fprintf( stdout, "\n" );
    fprintf( stdout, "Frame Number: %8u\n", uFrameNumber );
    fprintf( stdout, "Elements      : %8u\n", uElements );
    fprintf( stdout, "Flags          : 0x%04x\n", uFlags );
    for( uChannelCnt = 0; uChannelCnt < NUM_CHANNELS; ++uChannelCnt )
    {
        fprintf( stdout, "Channel %u Raw 0x%04x Voltage %f\n",
                uChannelCnt + 1, puRawData[
                uChannelCnt], (float) (
                (int) (puRawData[uChannelCnt])) * 0.000305175 / (float) ODAU_GAIN );
    } /* for */
    /*
     * STEP 5
     * Print out the digital data in hexadecimal.
     */
    fprintf( stdout, "Digital 0x%04x\n", puRawData[ NUM_CHANNELS ] );
} /* for */

```

**Figure 5-1: Retrieving/Displaying ODAU Real-time Data**

## 5.3 Retrieving Buffered ODAU Data

This section explains the steps that must be followed to retrieve buffered data from the ODAU, and write that information to a Northern Digital Floating Point format data file. Please see the [“Floating Point Programmer’s Guide” on page 65](#) and [“Real-time Data Types” on page 237](#) for an explanation of file types.

Each ODAU device in the Optotrak System has its own circular data buffer. Unlike the Optotrak System data buffer, which can store either raw data or 3D data, the ODAU data buffer can store only ODAU raw data. If the data in the file must be in voltage format, the application program can call the FileConvert routine (see [“FileConvert” on page 221](#)).

A spool mapping must be set for each device that spools buffered data to the host computer. If a particular device in the Optotrak System has no spool mapping set, then any buffered data available to the host computer from that device is discarded. To set a spool mapping for a device’s data buffer, use either of the routines `DataBufferInitializeFile` or `DataBufferInitializeMem`. These two routines identify the device (e.g. ODAU1), and either the file or memory block to which the buffered data is written (in [Figure 5-2 on page 49](#)).

Once the actual spooling of buffered data has started, all devices in the system send blocks of buffered data to the application program. The application can use either the blocking or non-blocking method of data spooling.

The blocking method uses the routine `DataBufferSpoolData`. This routine starts all devices in the Optotrak System spooling data back to the host computer. The Optotrak System ensures that all devices start sending data simultaneously at the same frame, and that the data remains synchronized throughout the spooling procedure. As the buffered data arrives at the host computer, it is written to the appropriate destination, based on the previously initialized spool mappings. When all of the buffered data has been received from all devices spooling data, the spooling procedure is completed and the `uSpoolStatus` variable is updated.

The application program should check the `uSpoolStatus` variable to ensure that it was set to a zero value. If it has been set to a non-zero value, then a buffering error has occurred during the spooling procedure and the data that has been received might not be valid. For more details on the `uSpoolStatus` parameter, refer to the routine documentation for [“DataBufferSpoolData” on page 201](#).

---

**Note** After the spool procedure is complete, the spool mappings that were last used are no longer valid. If the application program is to repeat the spooling procedure, it must reset the spool mappings for each device's data buffer. This prevents an application program from inadvertently overwriting the data just spooled. If a new spool mapping is not set for a device's data buffer, then any buffered data received by that device during the next spool procedure is discarded.

---

## Program Sample of Retrieving Buffered Data from ODAU

Figure 5-2 on page 49 is an example of how to retrieve ODAU buffered data. The `DataBufferInitializeFile` routine sets a spool mapping for a device's data buffer and identifies both the device (in this case ODAU1) and the file to which the buffered data is written (in this case O1#001.S13). The buffered data is retrieved from the ODAU using the blocking method invoked by the routine `DataBufferSpoolData`. This code fragment is similar to Sample Program 13 on the API CD.

To retrieve ODAU buffered data:

1. Configure the ODAU collection.
2. Configure the Optotrak System collection.
3. Activate the IRED markers.
4. Set the spool mapping for the Optotrak System buffered data.
5. Set the spool mapping for the ODAU buffered data.
6. Spool the buffered data to the files and display the spool status.

```

/*
 * STEP 1
 * Set up a collection for the ODAU1.
 */

    if( OdauSetupCollection(
        ODAU1,                /* Id the ODAU parameters are for. */
        NUM_CHANNELS,        /* Number of analog channels to collect.*/
        ODAU_GAIN,           /* Gain to use for the analog channels.*/
        ODAU_DIGITAL_INPB_INPA, /* Mode for the Digital I/O port.*/
        (float)100.0,        /* Frequency to collect data frames at.*/
        (float)90000.0,      /* Frequency to scan channels at.*/
        1,                   /* Stream mode for the data buffers.*/
        (float)2.0,          /* Number of seconds of data to collect.*/
        (float)0.0,          /* Number of seconds to pretrigger data.*/
        0 ) )                /* Flags. */
    {
        goto ERROR_EXIT;
    }

```

```
    } /* if */

    /*
    * STEP 2
    * Set up a collection for the OPTOTRAK.
    */
    if( OptotrakSetupCollection(
        NUM_MARKERS,    /* Number of markers in the collection. */
        (float)50.0,    /* Frequency to collect data frames at. */
        (float)2500.0,  /* Marker frequency for marker maximum on-time.*/
        30,             /* Dynamic or Static Threshold value to use. */
        160,            /* Minimum gain code amplification to use. */
        1,              /* Stream mode for the data buffers. */
        (float)0.20,    /* Marker Duty Cycle to use. */
        (float)6.0,     /* Voltage to use when turning on markers. */
        (float)2.0,     /* Number of seconds of data to collect. */
        (float)0.0,     /* Number of seconds to pre-trigger data by. */
        0 ) )
    {
    goto ERROR_EXIT;
    } /* if */

    /*
    * Activate the markers.
    */

    if( OptotrakActivateMarkers() )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
    * STEP 4
    * Initialize a file for spooling of the OPTOTRAK 3D data.
    */

    if( DataBufferInitializeFile( OPTOTRAK, "C#001.S13" ) )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
    * STEP 5
    * Initialize a file for spooling of the ODAU raw data.
    */
    if( DataBufferInitializeFile( ODAU1, "O1#001.S13" ) )
        goto ERROR_EXIT;
    } /* if */

    /*
    * STEP 6
    * Spool data to the previously initialized files.
```



```
*/  
fprintf( stdout, "Collecting Data Files\n" );  
if( DataBufferSpoolData( &uSpoolStatus ) )  
{  
    goto ERROR_EXIT;  
} /* if */  
fprintf( stdout, "Spool Status: 0x%04x\n", uSpoolStatus );
```

**Figure 5-2: : Collecting ODAU Data to Disk**



## 6 Real-time Rigid Body Programmer's Guide

Use this section of the guide to learn how to use the API rigid body routines. You need to have a sound understanding of rigid body concepts, terminology, algorithms, and data formats to understand the concepts discussed in this section. These are described in [“Flags and Settings Associated with Rigid Bodies”](#) on page 295 and in the *Optotrak Certus Rigid Body and Tool Design Guide*. You should also have read and understood all the API concepts presented in the [“Optotrak Programmer's Guide”](#) on page 23.

You may perform rigid body calculations if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform conversions internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to perform rigid body calculations.

This section of the API describes how to:

- retrieve real-time rigid body data
- change rigid body settings
- change the rigid body coordinate system
- transform raw or 3D data to 6D data
- check for undetermined transformations

The rigid body routines are listed in [“Rigid Body Specific Routines”](#) on page 206 and [“Rigid Body Related Routines”](#) on page 214.

### 6.1 Retrieving Real-time Rigid Body Data

Rigid body data is retrieved by informing the Optotrak System about the rigid body and then requesting data.

Use the routine `RigidBodyAddFromFile` to add the rigid body specified in a rigid body file (plate.rig in [Figure 6-1](#) on page 54) to the Optotrak System's rigid body list. `RigidBodyAddFromFile` reads the definition for the rigid body from the specified rigid body file, and adds the rigid body to the tracking list using the specified rigid body identification. Rigid bodies can also be added using the routine `RigidBodyAdd`.

The Optotrak System maintains a list of up to 10 rigid bodies for which it calculates transformations on the system (for transformations done on the host computer, the list can contain up to 85 rigid bodies).

---

**Note** The rigid body ID is used to identify the transformation data for this rigid body later on in the program. Each rigid body inserted in the Optotrak System's tracking list must have a unique rigid body ID (0 – 9 for on-system conversions, 0 – 84 for on-host conversions). The Optotrak System uses the rigid body ID to insert the definition into the tracking list. If an ID is used a second time the previous definition for the rigid body is lost.

---

Retrieve the rigid body transformation data from the Optotrak System with the routine `DataGetLatestTransforms`. This routine causes the Optotrak System to determine transformations for all rigid bodies in the tracking list using the latest complete frame of data. The transformation data is returned along with three variables; `uFrameNumber`, `uElements`, and `uFlags`.

**Table 6-1: Definition of Variables Returned with Transformation Data**

| Variable                  | Definition                                                                        |
|---------------------------|-----------------------------------------------------------------------------------|
| <code>uFrameNumber</code> | The number of frame data used in the calculations.                                |
| <code>uFlags</code>       | Indicates the current status of the Optotrak System data buffer.                  |
| <code>uElement</code>     | Is set to the number of rigid bodies for which transformation data are available. |

The rigid body transformation data returned by the Optotrak System contains the calculated transformation, the rigid body's ID and, the flags associated with the rigid body. The rigid body ID is the same as the ID used when the rigid body was added into the Optotrak System's tracking list. For a description of the flags parameter, refer to [“Checking for Undetermined Transforms” on page 63](#), and the documentation for the routines [“DataGetLatestTransforms” on page 163](#) and [“RigidBodyChangeSettings” on page 211](#).

The angular components of the transformation data are returned in one of three formats; Euler Angle, Quaternion, or Rotation Matrix. The return format of the transformation is specified when the rigid body is added to the tracking list. If no format is specified, the default format is Euler Angle. For a description of how to access the data elements of these three formats, see [“DataGetLatestTransforms” on page 163](#) and [“Optotrak Rigid Body Transformation Data” on page 242](#).

## Sample Program to Retrieve and Display Real-time Rigid Body Data

The code fragments in this section, and all following sections in this guide, assume that a collection with the appropriate parameters has been set up on the Optotrak System, and that the markers are activated. In this sample, the routine

RigidBodyAddFromFile adds the rigid body specified in the file “plate.rig” to the Optotrak System’s rigid body list. This code fragment is similar to Sample Program 9 on the API CD.

To retrieve and display real-time rigid body data from the Optotrak System:

1. Add the rigid body to the Optotrak System's rigid body tracking list.
2. Get the latest frame of transformation data.
3. Display the rigid body transformation data.

```

/*
 * STEP 1
 * Add a rigid body for tracking to the OPTOTRAK System from a.RIG file.
 */
if( RigidBodyAddFromFile(
    RIGID_BODY_ID, /* ID associated with this rigid body.*/
    1,             /* First marker in the rigid body.*/
    "plate",      /* RIG file containing rigid body coordinates.*/
    0 ) )        /* Flags. */
{
    goto ERROR_EXIT;
} /* if */
/*
 * Get and display ten frames of rigid body data.
 */
fprintf( stdout, "Rigid Body Data Display\n" );
for( uFrameCnt = 0; uFrameCnt < 10; ++uFrameCnt )
{
    /*
     * STEP 2
     * Get a frame of data.
     */
    if( DataGetLatestTransforms( &uFrameNumber, &uElements, &uFlags,
        &RigidBodyData ) )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
     * STEP 3
     * Print out the data.
     */
    fprintf( stdout, "\n" );
    fprintf( stdout, "Frame Number: %8u\n", uFrameNumber );
    fprintf( stdout, "Transforms : %8u\n", uElements );
    fprintf( stdout, "Flags      : 0x%04x\n", uFlags );
    for( uRigidCnt = 0; uRigidCnt < uElements; ++uRigidCnt )
    {
        fprintf( stdout, "Rigid Body %u\n",
            RigidBodyData.pRigidData[uRigidCnt].RigidId );
    }
}

```

```
fprintf( stdout, "XT = %8.2f YT = %8.2f ZT = %8.2f\n",
  RigidBodyData.pRigidData[ uRigidCnt].transformation.
  euler.translation.x, RigidBodyData.pRigidData
  [ uRigidCnt].transformation.          euler.translation.y,
  RigidBodyData.pRigidData[ uRigidCnt].transformation.
  euler.translation.z );
fprintf( stdout, "Y = %8.2f P = %8.2f R = %8.2f\n",
  RigidBodyData.pRigidData
  [ uRigidCnt].transformation.          euler.rotation.yaw,
  RigidBodyData.pRigidData[ uRigidCnt].transformation.
  euler.rotation.pitch,                RigidBodyData.pRigidData[
  uRigidCnt].transformation.          euler.rotation.roll );
} /* for */
} /* for */
```

**Figure 6-1: Retrieving and Displaying Rigid Body Data**

## 6.2 Changing Rigid Body Settings

This section shows you how to add a rigid body and change the settings for the added body.

The routine `RigidBodyAddFromFile` adds a rigid body to the list of rigid bodies that the Optotrak System tracks while the routine `RigidBodyChangeSettings` is used to customize some parameters.

The routine `RigidBodyChangeSettings` can be used to customize some parameters:

- minimum number of markers
- cut off angle for marker inclusion
- maximum 3D marker error
- maximum raw sensor error
- maximum 3D RMS marker error
- maximum raw sensor RMS error

If you do not specify these parameters using `RigidBodyChangeSettings`, the default values will be used. `RigidBodyChangeSettings` can also be used to change these parameters for any rigid body that needs unique settings. A complete listing of all the locations where parameters are set for rigid bodies can be found [“Flags and Settings Associated with Rigid Bodies” on page 295](#).

## Sample Program That Changes Some Defaults In Rigid Body Settings

The sample code fragment in [Figure 6-2 on page 57](#) changes the settings for the transformation format, the algorithm used to determine the transformation and the parameters listed above

The routine `DataGetLatestTransforms` retrieves the latest frame of rigid body transformation data.

---

**Note** `DataGetLatestTransforms` always includes the 3D marker data. The related routine, `DataGetLatestTransforms2`, which optionally includes the 3D data, can be used with some minor changes to the program's buffer allocations. Refer to ["DataGetLatestTransforms" on page 163](#) and ["DataGetLatestTransforms2" on page 165](#) for a complete description of these routines.

---

The data display in [Figure 6-2 on page 57](#) is different from [Figure 6-1 on page 54](#). First, since the rigid body transformation data is in quaternion format, it is displayed as such. Second, this code fragment displays the 3D marker data used to determine the rigid body transformations after the rigid body transformations have been displayed. 3D marker data is available to the application program if required and in this sample they are simply printed to the screen. This code fragment is similar to Sample Program 10 on the API CD.

To add a rigid body, change default rigid body parameter settings and retrieve the data:

1. Add the rigid body to the Optotrak System's rigid body tracking list.
2. Change some of the default rigid body parameter settings.
3. Get the latest frame of transformation data.
4. Display the rigid body transformation data.
5. Display associated 3D marker data.

```

/* STEP 1
 * Add rigid body 1 for tracking to the OPTOTRAK System from a .RIG file*/

if( RigidBodyAddFromFile(
    RIGID_BODY_ID, /* ID associated with this rigid body.*/
    1,             /* First marker in the rigid body.*/
    "plate",      /* RIG file containing rigid body coordinates.*/
    0 ) )        /* Flags. */
{
    goto ERROR_EXIT;
}

```

```
    } /* if */

/* STEP 2
*Change the default settings for this rigid body 1
*/

if( RigidBodyChangeSettings(
    RIGID_BODY_1, /* ID associated with this rigid body. */
    4,           /* Minimum number of markers that must be seen before
    performing rigid body calculations.*/
    60,         /* Cut off angle for marker inclusion in calcs.*/
    (float)0.25, /* Maximum 3D marker error for this rigid body.*/
    (float)1.0, /* Maximum raw sensor error for this rigid body. */
    (float)1.0, /* Max 3D RMS marker error for this rigid body. */
    (float)1.0, /* Max raw sensor RMS error for this rigid body. */
                                OPTOTRAK_QUATERN_RIGID_FLAG
    | OPTOTRAK_RETURN_QUATERN_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
*Get and display ten frames of rigid body data
*/
fprintf( stdout, "Rigid Body Data Display\n" );
for( uFrameCnt = 0; uFrameCnt < 10; ++uFrameCnt )
{
    /*
    * STEP 3
    * Get a frame of data*/

    if( DataGetLatestTransforms( &uFrameNumber, &uElements,
        &uFlags, &RigidBodyData ) )
    {
        goto ERROR_EXIT;
    } /* if */

/* STEP 4
*Print out the rigid body transformation data
*/

    fprintf( stdout, "\n" );
    fprintf( stdout, "Rigid Body Transformation Data\n\n" );
    fprintf( stdout, "Frame Number: %8u\n", uFrameNumber );
    fprintf( stdout, "Transforms   : %8u\n", uElements );
    fprintf( stdout, "Flags       : 0x%04x\n", uFlags );
    for( uRigidCnt = 0; uRigidCnt < uElements; ++uRigidCnt )
    {
```



```

fprintf( stdout, "Rigid Body %u\n",
         RigidbodyData.pRigidData[ uRigidCnt].RigidId );
fprintf( stdout, "XT = %8.2f YT = %8.2f ZT = %8.2f\n",
         RigidbodyData.pRigidData[ uRigidCnt].transformation.
         quaternion.translation.x,          RigidbodyData.pRigidData[
         uRigidCnt].transformation.
         quaterion.translation.y,          RigidbodyData.pRigidData[
         uRigidCnt].transformation.
         quaterion.translation.z );
fprintf( stdout, "Q0 = %8.2f QX = %8.2f QY = %8.2f QZ =
         %8.2f\n",RigidbodyData.pRigidData[ uRigidCnt].transformation.
         quaterion.rotation.q0,          RigidbodyData.pRigidData[
         uRigidCnt].transformation.      quaterion.rotation.qx,
         RigidbodyData.pRigidData[ uRigidCnt].transformation.
         quaterion.rotation.qy,          RigidbodyData.pRigidData[
         uRigidCnt].transformation.      quaterion.rotation.qz
         );
} /* for */

/*STEP 5
*Print out the 3D data
*/

fprintf( stdout, "\nAssociated 3D Marker Data\n\n" );
for( uMarkerCnt = 0; uMarkerCnt < 6; ++uMarkerCnt )
{
    fprintf( stdout, "Marker %u X %f Y %f Z %f\n", uMarkerCnt + 1,
            RigidbodyData.p3dData[ uMarkerCnt].x,          RigidbodyData.p3dData[
            uMarkerCnt].y,          RigidbodyData.p3dData[ uMarkerCnt].z );
} /* for */
} /* for */

```

**Figure 6-2: Changing Default Rigid Body Parameter Settings**

## 6.3 Changing the Rigid Body Coordinate System

This section will show you how to change the coordinate system for a rigid body. This procedure uses some of the routines discussed in the previous section.

The rigid body transformations returned by the Optotrak System are in the coordinate system defined by the camera parameters that were first sent to the Optotrak System in the Program Initialization Code section of the application program (see [“Initializing, Retrieving System Status and Exiting from the Optotrak System” on page 24](#)). Some applications may need to have the transformations expressed relative to either a different coordinate system or to a constantly changing coordinate system. [Figure 6-3 on page 60](#) includes a sample code fragment that demonstrates the required code to accommodate a constantly changing coordinate system.

To illustrate how this is done, two rigid bodies are required. The first rigid body, given the ID “RIGID\_BODY\_1,” defines the coordinate system in which all other

transformations are measured. This rigid body is added to the Optotrak System's rigid body tracking list using the routine `RigidBodyAddFromFile`.

The second rigid body, given the ID "RIGID\_BODY\_2," is then added to the Optotrak System's tracking list using the routine `RigidBodyAdd`. Instead of using a rigid body file to define the rigid body, the application program specifies the required parameters to this routine, namely the number of markers defining the rigid body, and an array of 3D positions for these markers.

The routine `RigidBodyChangeFOR` rotates the transformations determined for RIGID\_BODY\_2 into the coordinate system defined by RIGID\_BODY\_1. The parameters to this routine specify which rigid body defines the coordinate system, and the mode for changing the coordinate system.

### Sample Program that Allows For a Constantly Changing Coordinate System

After all of the rigid bodies have been added and the coordinate system specified, the code fragment in [Figure 6-3 on page 60](#) requests and displays the rigid body transformation data, as well as the associated 3D marker data. With each frame of rigid body data, the Optotrak System also returns the 3D marker data used to determine the rigid body transformations. If the transformations are being expressed in a different coordinate system, the 3D data returned with the transformations are also transformed into this coordinate system.

---

**Note** This only applies to the 3D data retrieved using the routine `DataGetLatestTransforms`, and not to 3D data retrieved using the routine `DataGetLatest3D`. [Figure 6-3 on page 60](#) illustrates an example session on how to express transformations in a different coordinate system.

---

In this case, the mode for changing the coordinate system is `OPTOTRAK_CONSTANT_RIGID_FLAG`, which means that the coordinate system is determined each time a frame of rigid body transformation data is produced. This code fragment is similar to Sample Program 11 on the API CD.

To express transformations in a different coordinate system:

1. Add the first rigid body to the Optotrak System's rigid body tracking list.
2. Add the second rigid body to the Optotrak System's rigid body tracking list.
3. Change the coordinate system for rigid body transformation data.
4. Get the latest frame of transformation data.
5. Display the rigid body transformation data.
6. Display the transformed 3D marker data.

```

/*
 * STEP 1
 * Add rigid body 1 for tracking to the OPTOTRAK System from a .RIG file.
 */
if( RigidBodyAddFromFile(
    RIGID_BODY_1, /* ID associated with this rigid body.*/
    1, /* First marker in the rigid body.*/
    "plate", /* RIG file containing rigid body coordinates.*/
    OPTOTRAK_QUATERN_RIGID_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 2
 * Add rigid body 2 for tracking to the OPTOTRAK System from an array of *3D
 * points.
 */
if( RigidBodyAdd(
    RIGID_BODY_2, /* ID associated with this rigid body. */
    7, /* First marker in the rigid body. */
    6, /* Number of markers in the rigid body. */(float
    *)RigidBody2, /* 3D coords for each marker in the body. */
    NULL, /* no normals for this rigid body. */
    OPTOTRAK_QUATERN_RIGID_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 3
 * Change the default coordinate system to be defined by rigid body one.
 */

if( RigidBodyChangeFOR( RIGID_BODY_1, OPTOTRAK_CONSTANT_RIGID_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Get and display ten frames of rigid body data.
 */
fprintf( stdout, "Rigid Body Data Display\n" );
for( uFrameCnt = 0; uFrameCnt < 10; ++uFrameCnt )
{
    /*
     * STEP 4
     * Get a frame of data.
     */

```

```
if( DataGetLatestTransforms( &uFrameNumber, &uElements, &uFlags,
    &RigidBodyData ) )
{
    goto ERROR_EXIT;
} /* if */

/*
* STEP 5
* Print out the rigid body transformation data.
*/
fprintf( stdout, "\n" );
fprintf( stdout, "Rigid Body Transformation Data\n\n" );
fprintf( stdout, "Frame Number: %8u\n", uFrameNumber );
fprintf( stdout, "Transforms : %8u\n", uElements );
fprintf( stdout, "Flags      : 0x%04x\n", uFlags );
for( uRigidCnt = 0; uRigidCnt < uElements; ++uRigidCnt )
{
    fprintf( stdout, "Rigid Body %u\n",
        RigidBodyData.pRigidData[ uRigidCnt].RigidId );
    fprintf( stdout, "XT = %8.2f YT = %8.2f ZT = %8.2f\n",
        RigidBodyData.pRigidData[ uRigidCnt].transformation.
        euler.translation.x,
        RigidBodyData.pRigidData[
        uRigidCnt].transformation.
        euler.translation.y,
        RigidBodyData.pRigidData[ uRigidCnt].transformation.
        euler.translation.z );
    fprintf( stdout, "Y = %8.2f P = %8.2f R = %8.2f\n",
        RigidBodyData.pRigidData[ uRigidCnt].transformation.
        euler.rotation.yaw,
        RigidBodyData.pRigidData[
        uRigidCnt].transformation.
        euler.rotation.pitch,
        RigidBodyData.pRigidData[ uRigidCnt].transformation.
        euler.rotation.roll );
} /* for */

/*
* STEP 6
* Print out the 3D data.
*/
fprintf( stdout, "\nAssociated 3D Marker Data\n\n" );
for( uMarkerCnt = 0; uMarkerCnt < 12; ++uMarkerCnt )
{
    fprintf( stdout, "Marker %u X %f Y %f Z %f\n", uMarkerCnt + 1,
        RigidBodyData.p3dData[ uMarkerCnt].x,
        RigidBodyData.p3dData[ uMarkerCnt].y,
        RigidBodyData.p3dData[ uMarkerCnt].z );
} /* for */
} /* for */
```

**Figure 6-3: Expressing Transformations in a Different Coordinate System**

## 6.4 Transforming Previously Obtained Data

The sample code fragments discussed in this section will help you to understand how to transform previously obtained raw or 3D data into 6D data. The previous examples in this chapter obtain rigid body 6D data directly from the system, in real-time, by using the routine `DataGetLatestTransforms`.

The routine `OptotrakConvertRawTo3D` converts raw data to 3D position data on a frame-by-frame basis according to the current system camera parameters. `OptotrakConvertRawTo3D` returns the number of data elements (markers) comprising the frame as `uElements`, and the converted 3D positions in the array `pdtPositionData`. The raw data may have been collected by spooling directly to a memory buffer, or to a file that was later loaded to memory.

After the centroid data has been converted to 3D data, use the routine `OptotrakConvertTransforms` to transform the 3D position data to 6D data. After conversion, `OptotrakConvertTransforms` returns the number of elements (one rigid body in [Figure 6-4 on page 62](#)) and the transformed 6D data in the structure `dtRigidBodyData`.

### Sample Program to Convert and Transform Raw Data to 6D

[Figure 6-4 on page 62](#) is an example of how to convert previously collected full raw data to 3D data, and then transform the 3D data to 6D data on a frame by frame basis. This code fragment is similar to Sample Program 19 on the API CD. To convert and transform raw data to 6D data:

1. Convert each frame of raw data to its corresponding 3D position data.
2. Display the converted 3D position data.
3. Transform the converted 3D position data to 6D data.
4. Display the transformed 6D position data.

```
fprintf( stdout, "\n\n3D Position data:\n" );
pfRawDataCur  = pfRawData;
pdtPositionCur = pdtPositionData;
for( lnFrameCnt = 0; lnFrameCnt < lnFileFrames; lnFrameCnt++ )
{
    /*
    * STEP 1
    * Convert the raw data to 3D position data.
    */
    if( OptotrakConvertRawTo3D( &uElements,
                               pfRawDataCur,
                               pdtPositionCur ) )
    {
```

```
        goto ERROR_EXIT;
    } /* if */
    /*
    * STEP 2
    * Display the converted 3D position data.
    */
    for( uElementCnt = 0; uElementCnt < uElements; uElementCnt++ )
    {
    fprintf( stdout,
            "\n%5ld %3u %12.5f %12.5f %12.5f",
            lnFrameCnt,
            uElementCnt,
            pdtPositionCur[uElementCnt].x,
            pdtPositionCur[uElementCnt].y,
            pdtPositionCur[uElementCnt].z );
        } /* for */
    pfRawDataCur += uElements * NUM_SENSORS;
    pdtPositionCur += uElements;
} /* for */
fprintf( stdout, "\n\nRigid body data:\n" );
pdtPositionCur = pdtPositionData;
for( lnFrameCnt = 0; lnFrameCnt < lnFileFrames; lnFrameCnt++ )
{
    /*
    * STEP 3
    * Transform the 3D position data to 6D.
    */
    if( OptotrakConvertTransforms( &uElements,
                                   &dtRigidBodyData,
                                   pdtPositionCur ) )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
    * STEP 4
    * Display the transformed 6D data.
    */
    fprintf( stdout, "Frame %04u\n", lnFrameCnt );
    fprintf( stdout,
            "XT = %12.6f YT = %12.6f ZT =
            %12.6f\n", dtRigidBodyData.transformation.euler.translation.x, dtRigidBo
            dyData.transformation.euler.translation.y, dtRigidBodyData.transforma
            tion.euler.translation.z );
    fprintf( stdout,
            "Y = %12.6f P = %12.6f R =
            %12.6f\n", dtRigidBodyData.transformation.euler.rotation.yaw, dtRigidBo
            dyData.transformation.euler.rotation.pitch, dtRigidBodyData.transforma
            tion.euler.rotation.roll ); pdtPositionCur += uElementCnt;
} /* for */
```

**Figure 6-4: Transforming Previously Collected Raw Data to 6D Data**

## 6.5 Checking for Undetermined Transforms

This section discusses how to check that all transformations were successfully determined.

There are two flags located in the file specified by the element `pdatadest` in the routines `DataGetLatestTransforms`, `DataGetLatestTransforms2`, `DataReceiveLatestTransforms` and `DataReceiveLatestTransforms2`. These two flags, `OPTOTRAK_UNDETERMINED_FLAG` AND `OPTOTRAK_RIGID_ERR_MKR_SPREAD`, indicate if either the transformation could not be determined or the spread in the markers went out of bounds from the set values.

For more information on rigid body flags and error settings, see “[Flags and Settings Associated with Rigid Bodies](#)” on page 295, “[Description of Rigid Body Real-time Data Elements](#)” on page 243, and “[DataGetLatestTransforms](#)” on page 163.

---

**Note** The flag `OPTOTRAK_UNDETERMINED_FLAG` must be examined for every transformation.

---

### Sample Program to Check For Undetermined Transformations and Marker Spread Errors

[Figure 6-5 on page 64](#) shows how to check for two possible causes of an undetermined transformation: the `OPTOTRAK_UNDETERMINED_FLAG` from the `OptotrakRigidStruct` structure and the `OPTOTRAK_RIGID_ERR_MKR_SPREAD` that results from rigid body marker spread errors.

This sample fragment follows these steps:

1. Get the latest frame of transformation data.
2. Check the flags member for indication of an undetermined transform.
3. Further check the flags member for indication of marker spread errors.

```

/*
 * STEP 1
 * Get a frame of data.
 */

if( DataGetLatestTransforms( &uFrameNumber, &uElements, &uFlags,
                            &RigidBodyData ) )
{
    goto ERROR_EXIT;
} /* if */

```

```
/*
 * STEP 2
 * Check the returned flags member for improper transforms.
 */

if( RigidBodyData.pRigidData[0].flags & OPTOTRAK_UNDETERMINED_FLAG )
{
    fprintf( stdout, "\nUndetermined transform!" );

    /*
 * STEP 3
 * Further check the returned flags member for marker spread errors.
 */
    if( RigidBodyData.pRigidData[0].flags & OPTOTRAK_RIGID_ERR_MKR_SPREAD )
    {
        fprintf( stdout, " Marker spread error." );
    } /* if */
    continue;
} /* if */
```

**Figure 6-5: Checking for Undetermined Transformations**



---

## 7 Floating Point Programmer's Guide

This section will help you to understand:

- the file format used for all data files: the Northern Digital Floating Point (NDFP).
- how to convert raw and full raw data from an Optotrak System or raw data from an ODAU system to 3D data and voltages, respectively.
- how to read an existing data file, do calculations on the data and write the results to a new data file.

In an Optotrak System, both the Position Sensor and the ODAU produce raw data. This raw data can be converted with the NDI ToolBench software or with API routines. An application program can convert the raw data to a more usable format anytime after the data was originally collected.

---

**Note** To convert raw data to a more usable format, the camera file from the original raw data collection must be used. Otherwise, the conversion will produce incorrect values.

---

If either a centroid (raw) or full raw data file is specified, the output data format will be 3D marker data. If an ODAU raw data file is being converted, the output data format will be voltages. The raw data file input to this routine must be in NDFP file format.

## 7.1 The Northern Digital Floating Point Format

The Optotrak System is capable of producing several types of data files. All data files created by either the NDI ToolBench software or the Optotrak System application programs are written in the NDFP file format. This generalized file format allows all NDI data processing programs to read and manipulate data files.

This file format is used to describe frames of data consisting of any combination of character, integer and floating point elements. The naming convention for floating point files is:

<prefix>#<trialnumber>.<extension>

some standard prefixes are 'C' for 3D data, 'R' for raw data, and 'Ox' for data generated by any analog/digital device. The trial number and file extension are user-defined.

A floating point file consists of a header describing the file, followed by the data itself. The data is organized into frames, with each frame describing a single component of the data (such as the 3D positions of various markers at a particular point in time during a data collection). Each frame is then divided into a number of items, with each item describing a component in that frame (for example, the 3D positions of one marker). Finally, every item is subdivided into a number of subitems, with each subitem describing a component in that item, (for example, the value for the X coordinate in a 3D position). The data is stored after the header as a sequence of 4-byte, floating point numbers. Any data that is considered "missing", such as the 3D position of a marker that was not detected by the Position Sensor for a number of frames, is assigned a special value of constant BAD\_FLOAT, which equals -3.697314E28 (-3.69731x10<sup>28</sup>). These data points become placeholders and are subsequently ignored by analysis programs. Any custom programs should check for these special values by comparing against constant MAX\_NEGATIVE before doing any calculations.

### 7.1.1 File Header

The file header is fixed-size (256 bytes) and describes the organization of the file data section. The file data section, which contains the data in binary format, follows immediately after the header. There are six conventions that apply to the file header format:

- starting offset and length are measured in bytes
- integers are stored low byte first
- long integers are stored low word then high word
- strings are stored as characters, terminated by a null character (hex 0)

- floating point numbers follow IEEE 4 byte float conventions
- double precision floating point numbers follow IEEE 8 byte float conventions

Table 7-1 on page 67 describes the header.

**Table 7-1: Floating Point File Header Description**

| Field | Name         | Start | Size | Type         | Description                                               |
|-------|--------------|-------|------|--------------|-----------------------------------------------------------|
| 1     | filetype     | 0     | 1    | byte         | File type identifier (always 32 for floating point files) |
| 2     | items        | 1     | 2    | integer      | Number of items in each frame of data                     |
| 3     | subitems     | 3     | 2    | integer      | Number of floating point subitems in each frame           |
| 4     | numframes    | 5     | 4    | long integer | Number of frames of data                                  |
| 5     | frequency    | 9     | 4    | float        | Frequency at which data was collected                     |
| 6     | UserComments | 13    | 60   | string       | A user-supplied comment                                   |
| 7     | SysComments  | 73    | 60   | string       | May contain file information                              |
| 8     | DescripFile  | 133   | 30   | string       | Unused at present                                         |
| 9     | cutoff       | 163   | 2    | integer      | Cutoff frequency used if 3D data was filtered             |
| 10    | CollTime     | 165   | 10   | string       | The time at which the data was collected (hh:mm:ss)       |
| 11    | CollDate     | 175   | 10   | string       | The date on which the data was collected (mm/dd/yy)       |
| 12    | FrameStart   | 185   | 4    | long integer | Unused at present.                                        |

**Table 7-1: Floating Point File Header Description (Continued)**

| Field | Name            | Start | Size | Type    | Description                                                         |
|-------|-----------------|-------|------|---------|---------------------------------------------------------------------|
| 13    | Extended-Header | 189   | 2    | integer | Contains the value 12345 if the following four fields are also used |
| 14    | CharSubitems    | 191   | 2    | integer | Number of character subitems                                        |
| 15    | IntSubitems     | 193   | 2    | integer | Number of integer subitems                                          |
| 16    | DoubleSubitems  | 195   | 2    | integer | Number of double precision floating point subitems                  |
| 17    | ItemSize        | 197   | 2    | integer | Size of each item, including all its subitems (in bytes)            |
| 18    | padding         | 199   | 57   | string  | Unused                                                              |

### 7.1.2 Sample Formats

The following tables describe the floating point file format for the most common files generated by the Optotrak System.

**Table 7-2: Optotrak 3D Data Format**

| Items   | Marker 1       |                |                | Marker 2       |                |                | ... | Marker M       |                |                |
|---------|----------------|----------------|----------------|----------------|----------------|----------------|-----|----------------|----------------|----------------|
| frame 1 | X <sub>1</sub> | Y <sub>1</sub> | Z <sub>1</sub> | X <sub>1</sub> | Y <sub>1</sub> | Z <sub>1</sub> | ... | X <sub>1</sub> | Y <sub>1</sub> | Z <sub>1</sub> |
| frame 2 | X <sub>2</sub> | Y <sub>2</sub> | Z <sub>2</sub> | X <sub>2</sub> | Y <sub>2</sub> | Z <sub>2</sub> | ... | X <sub>2</sub> | Y <sub>2</sub> | Z <sub>2</sub> |
| ⋮       | ⋮              | ⋮              | ⋮              | ⋮              | ⋮              | ⋮              |     | ⋮              | ⋮              | ⋮              |
| frame N | X <sub>N</sub> | Y <sub>N</sub> | Z <sub>N</sub> | X <sub>N</sub> | Y <sub>N</sub> | Z <sub>N</sub> | ... | X <sub>N</sub> | Y <sub>N</sub> | Z <sub>N</sub> |

where:

Number of frames = N

Number of items/frame = M (number of markers)

Float subitems/item = 3

Item size = 12 bytes  
 X, Y and Z translations are in mm.

**Table 7-3: Optotrak Raw (Centroid) Data Format**

| Item    | Marker 1                                         | Marker 2                                         |     | Marker M                                         |
|---------|--------------------------------------------------|--------------------------------------------------|-----|--------------------------------------------------|
| frame 1 | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> | ... | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> |
| frame 2 | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> | ... | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> |
| ⋮       | ⋮ ⋮ ⋮                                            | ⋮ ⋮ ⋮                                            | ⋮   | ⋮ ⋮ ⋮                                            |
| frame N | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> | ... | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> |

where:

Number of frames = N  
 Number of items/frame = M (number of markers)  
 Float subitems/item = S (number of sensors)  
 Item Size = (S x 4) bytes

**Table 7-4: Optotrak Full Raw Data Format**

| Item    | Marker 1                                                                                                         |     | Marker M                                                                                                         |
|---------|------------------------------------------------------------------------------------------------------------------|-----|------------------------------------------------------------------------------------------------------------------|
| frame 1 | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> status <sub>1</sub> status <sub>2</sub> ... status <sub>s</sub> | ... | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> status <sub>1</sub> status <sub>2</sub> ... status <sub>s</sub> |
| frame 2 | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> status <sub>1</sub> status <sub>2</sub> ... status <sub>s</sub> | ... | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> status <sub>1</sub> status <sub>2</sub> ... status <sub>s</sub> |
| ⋮       | ⋮ ⋮ ⋮ ⋮ ⋮                                                                                                        | ⋮   | ⋮ ⋮ ⋮ ⋮ ⋮                                                                                                        |
| frame N | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> status <sub>1</sub> status <sub>2</sub> ... status <sub>s</sub> | ... | r <sub>1</sub> r <sub>2</sub> ... r <sub>s</sub> status <sub>1</sub> status <sub>2</sub> ... status <sub>s</sub> |

where:

Each status element contains the signal strength high, amplification, error code and signal strength low.  
 Number of frames = N  
 Number of items/frame = M (number of markers)  
 Float subitems/item = S (number of sensors)  
 Char subitems/item = 4 x S  
 Item size = (S x 8) bytes

**Table 7-5: ODAU Raw Data Format, Without Digital Input**

| Item    | Channel        |                |       |                |
|---------|----------------|----------------|-------|----------------|
| frame 1 | r <sub>1</sub> | r <sub>2</sub> | ..... | r <sub>C</sub> |
| frame 2 | r <sub>1</sub> | r <sub>2</sub> | ..... | r <sub>C</sub> |
| ⋮       | ⋮              | ⋮              |       | ⋮              |
| frame N | r <sub>1</sub> | r <sub>2</sub> | ..... | r <sub>C</sub> |

where:

Number of frames = N

Number of items/frame = 1

Integer subitems/item = C (number of channels)

Item size = (C x 2) bytes

**Table 7-6: ODAU Raw Data Format, With Digital Input**

| Item    | Channel        |                |       |                | Digital Input |
|---------|----------------|----------------|-------|----------------|---------------|
| frame 1 | r <sub>1</sub> | r <sub>2</sub> | ..... | r <sub>C</sub> | D             |
| frame 2 | r <sub>1</sub> | r <sub>2</sub> | ..... | r <sub>C</sub> | D             |
| ⋮       | ⋮              | ⋮              |       | ⋮              | ⋮             |
| frame N | r <sub>1</sub> | r <sub>2</sub> | ..... | r <sub>C</sub> | D             |

where:

Number of Frames = N

Number of items/frame = 1

Integer subitems/item = [C (number of channels) + 1]

Item Size = [(C x 2) + 2] bytes

**Table 7-7: ODAU Converted Data Format, Voltages,  
No Digital Input**

| Item    | Channel |       |       |       |
|---------|---------|-------|-------|-------|
| frame 1 | $v_1$   | $v_1$ | ..... | $v_C$ |
| frame 2 | $v_1$   | $v_2$ | ..... | $v_C$ |
| ⋮       | ⋮       | ⋮     |       | ⋮     |
| frame N | $v_1$   | $v_2$ | ..... | $v_C$ |

where:

Number of frames = N

Number of items/frame = 1

Float subitems/item = C (number of channels)

Item size = (Cx4) bytes

**Table 7-8: ODAU Converted Data Format, Voltages  
With Digital Input**

| Item    | Channel |       |       |       | Digital Input |
|---------|---------|-------|-------|-------|---------------|
| Frame 1 | $v_1$   | $v_2$ | ..... | $v_C$ | D             |
| Frame 2 | $v_1$   | $v_2$ | ..... | $v_C$ | D             |
| ⋮       | ⋮       | ⋮     |       | ⋮     | ⋮             |
| Frame N | $v_1$   | $v_2$ | ..... | $v_C$ | D             |

where:

Number of Frames = N

Number of items/frame = 1

Float subitems/item = C (number of channels)

Integer subitems/item = 1

Item size = ((Cx4) + 2) bytes

**Tables 7-9 to 7-11 refer to files that are generated by the Data Analysis Package.**

**Table 7-9: Data Analysis Package for One Rigid Body Euler Representation**

| Item    | Rotation Values (radians) |                |                | Translation Values (mm) |                |                | Error |
|---------|---------------------------|----------------|----------------|-------------------------|----------------|----------------|-------|
| frame 1 | R <sub>z</sub>            | R <sub>y</sub> | R <sub>x</sub> | T <sub>x</sub>          | T <sub>y</sub> | T <sub>z</sub> | error |
| frame 2 | R <sub>z</sub>            | R <sub>y</sub> | R <sub>x</sub> | T <sub>x</sub>          | T <sub>y</sub> | T <sub>z</sub> | error |
| ⋮       | ⋮                         | ⋮              | ⋮              | ⋮                       | ⋮              | ⋮              | ⋮     |
| frame N | R <sub>z</sub>            | R <sub>y</sub> | R <sub>x</sub> | T <sub>x</sub>          | T <sub>y</sub> | T <sub>z</sub> | error |

where:

- N = Number of frames
- Number of items/frame = number of rigid bodies
- Float subitems/item = 7
- Item size = 28 bytes

**Table 7-10: Data Analysis Package for One Rigid Body Quaternion Representation.**

| Item    | Rotation Quaternion Format |                |                |                | Translation Value (mm) |                |                | Error |
|---------|----------------------------|----------------|----------------|----------------|------------------------|----------------|----------------|-------|
| frame 1 | q <sub>0</sub>             | q <sub>x</sub> | q <sub>y</sub> | q <sub>z</sub> | T <sub>x</sub>         | T <sub>y</sub> | T <sub>z</sub> | error |
| frame 2 | q <sub>0</sub>             | q <sub>x</sub> | q <sub>y</sub> | q <sub>z</sub> | T <sub>x</sub>         | T <sub>y</sub> | T <sub>z</sub> | error |
| ⋮       | ⋮                          | ⋮              | ⋮              | ⋮              | ⋮                      | ⋮              | ⋮              | ⋮     |
| frame N | q <sub>0</sub>             | q <sub>x</sub> | q <sub>y</sub> | q <sub>z</sub> | T <sub>x</sub>         | T <sub>y</sub> | T <sub>z</sub> | error |

where:

- Number of frames = N
- Number of items/frame = number of rigid bodies
- Float subitems/item = 8
- Item size = 32 bytes



**Table 7-11: Data Analysis Package for One Rigid Body  
Rotation Matrix Representation**

| Item    | Rotation Matrix Values |                 |                 |                 |                 |                 |                 |                 |                 | Translation Values (mm) |                |                | Error |
|---------|------------------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-------------------------|----------------|----------------|-------|
| frame 1 | R <sub>00</sub>        | R <sub>01</sub> | R <sub>02</sub> | R <sub>10</sub> | R <sub>11</sub> | R <sub>12</sub> | R <sub>20</sub> | R <sub>21</sub> | R <sub>22</sub> | T <sub>x</sub>          | T <sub>y</sub> | T <sub>z</sub> | error |
| frame 2 | R <sub>00</sub>        | R <sub>01</sub> | R <sub>02</sub> | R <sub>10</sub> | R <sub>11</sub> | R <sub>12</sub> | R <sub>20</sub> | R <sub>21</sub> | R <sub>22</sub> | T <sub>x</sub>          | T <sub>y</sub> | T <sub>z</sub> | error |
| ⋮       | ⋮                      | ⋮               | ⋮               | ⋮               | ⋮               | ⋮               | ⋮               | ⋮               | ⋮               | ⋮                       | ⋮              | ⋮              | ⋮     |
| frame N | R <sub>00</sub>        | R <sub>01</sub> | R <sub>02</sub> | R <sub>10</sub> | R <sub>11</sub> | R <sub>12</sub> | R <sub>20</sub> | R <sub>21</sub> | R <sub>22</sub> | T <sub>x</sub>          | T <sub>y</sub> | T <sub>z</sub> | error |

where:

Number of frames = N

Number of items/frame = number of rigid bodies

Float subitems/item = 13

Item size = 52 bytes

## 7.2 Converting Optotrak and ODAU Raw Data Files

This section explains how to convert raw data files or full raw data files to 3D data files and ODAU raw data files to voltage files.

Use the routine FileConvert to convert raw data files once the Optotrak System and ODAU raw data files have been collected (see [“Retrieving Buffered Optotrak Data” on page 34](#) and [“Retrieving Buffered ODAU Data” on page 46](#)). The parameters for this routine are: the name of the input raw data file, the name of the output file for the converted data, and the type of raw data that is being converted.

To convert a raw or full raw data file with the routine FileConvert, the data is read from the input file, converted to 3D using previously loaded camera parameters, and written to the output file.

---

**Note** An application program converting raw data or full raw data must ensure that the correct camera parameters have been loaded.

---

Conversion of the ODAU raw data files is done on the host computer and requires no special parameters.

## Sample Program that Converts an Optotrak System Data File and an ODAU Raw Data File

Figure 7-1 on page 76 is a sample of how to convert an Optotrak System data file and an ODAU raw data file. This code fragment is similar to Sample Program 14 on the API CD.

To convert an Optotrak System data file and an ODAU raw data file, follow the steps listed below.

1. Collect the Optotrak System data file and the ODAU raw data file.
2. Convert the Optotrak System data file to a 3D format.
3. Convert the ODAU raw data file to a voltage format.

```
/*
 * STEP 1
 * Set up a collection for the ODAU.
 */
if( OdausetupCollection(
    ODAU1,          /* Id the ODAU parameters are for. */
    NUM_CHANNELS,  /* Number of analog channels to collect. */
    ODAU_GAIN,     /* Gain to use for the analog channels.*/
    ODAU_DIGITAL_INPB_INPA, /* Mode for the Digital I/O port.*/
    (float)100.0,  /* Frequency to collect data frames at. */
    (float)90000.0, /* Frequency to scan channels at. */
    1,            /* Stream mode for the data buffers. */
    (float)2.0,   /* Number of seconds of data to collect. */
    (float)0.0,   /* Number of seconds to pretrigger data. */
    0 ) )        /* Flags. */
{
    goto ERROR_EXIT;
} /* if */

/*
 * Set up a collection for the OPTOTRAK.
 */
if( OptotrakSetupCollection(
    NUM_MARKERS,   /* Number of markers in the collection. */
    (float)50.0,   /* Frequency to collect data frames at. */
    (float)2500.0, /* Marker frequency for marker maximum on-time. */
    30,           /* Dynamic or Static Threshold value to use. */
    160,         /* Minimum gain code amplification to use. */
    1,           /* Stream mode for the data buffers. */
    (float)0.20,  /* Marker Duty Cycle to use. */
    (float)6.0,   /* Voltage to use when turning on markers. */
    (float)2.0,   /* Number of seconds of data to collect. */
    (float)0.0,   /* Number of seconds to pre-trigger data by. */
```

```
        OPTOTRAK_BUFFER_RAW_FLAG ) )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
    * Activate the markers.
    */
    if( OptotrakActivateMarkers() )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
    * Initialize a file for spooling of the OPTOTRAK centroid data.
    */

    if( DataBufferInitializeFile( OPTOTRAK, "R#001.S14" ) )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
    * Initialize a file for spooling of the ODAU raw data.
    */
    if( DataBufferInitializeFile( ODAU1, "O1#001.S14" ) )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
    * Spool data to the previously initialized files.
    */
    fprintf( stdout, "Collecting Data Files\n" );
    if( DataBufferSpoolData( &uSpoolStatus ) )
    {
        goto ERROR_EXIT;
    } /* if */
    fprintf( stdout, "Spool Status: 0x%04x\n", uSpoolStatus );

    /*
    * De-activate the markers.
    */
    if( OptotrakDeActivateMarkers() )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
    * STEP 2
    * Convert the centroid data in the file 'R#001.S14' and write the
    * 3D data to the file 'C#001.S14'
    */
    fprintf( stdout, "Converting OPTOTRAK centroid data file\n" );
```

```

if( FileConvert( "R#001.S14", "C#001.S14", OPTOTRAK_RAW ) )
{
    goto ERROR_EXIT;
} /* if */
fprintf( stdout, "File conversion complete\n" );

/*
 * STEP 3
 * Convert the ODAU raw data file to a voltage data file.
 */

fprintf( stdout, "Converting ODAU Data File\n" );
if( FileConvert( "O1#001.S14", "V1#001.S14", ODAU_RAW ) )
{
    goto ERROR_EXIT;
} /* if */
fprintf( stdout, "Conversion Complete\n" );

```

**Figure 7-1: Converting Optotrak System Data Files and ODAU Raw Data Files**

## 7.3 Processing NDFP Format Files

Review this section for information on how to read data in an NDFP file, do calculations based on the data and write to an NDFP file.

A data file can be opened in one of two modes - read-only mode or read-write mode. When a file is opened, the application program must use a unique file ID so that it can identify the file in later operations. When a file is opened in read-only mode, the necessary values are read from the header (which describes the organization of the data in the file) into the parameters of the FileOpen routine for subsequent use by the application program. If a file is opened in read-write mode, the application program must correctly specify these parameters so that they can be written properly to the file header. For a detailed description of the file header parameters used in FileOpen, see “FileOpen” on page 222 and “FileOpenAll” on page 224.

If a file is opened in read-only mode, use the routine FileRead to read data, see “FileRead” on page 226. Data is read directly from the specified file into memory allocated by the application program. You can specify the number of frames to be read with a single invocation of FileRead.

If a file is opened in read-write mode, the application program can write data frames to it. The routine FileWrite (“FileWrite” on page 228) is used to write data to a specified file. Once the application program is finished with a data file, it must close the data file using the routine FileClose (“FileClose” on page 220).




---

**If the file is not closed properly before the program exits, data in the file could be lost.**

---

[Figure 7-2 on page 81](#) is a sample of how to read and write data to and from NDFP format files.

Collected files are processed in this order:

1. The collected file is opened in read-only mode.
2. The first frame of the input file is read and stored as the base frame for the calculation phase.
3. Each frame of data from the input file is read and has the values from the base frame subtracted from it.
4. The resultant data frame is written to the output file.
5. The input and output files are closed and the program exits.

### Sample Program to Read, Calculate and Write Data in NDFP Files

The first step in the sample program collects a 3D data file once marker 1 comes into view. A program similar to this was presented in [“Advanced Buffered Data Retrieval Without Blocking” on page 36](#) and is not discussed here.

Once the data file has been collected, the sample program shuts down the processor system in the Optotrak System, since communication with the Optotrak System is not required for the file processing routines presented here. The program then processes the data file. This sample code fragment is similar to Sample Program 15 on the API CD.

The code fragment in [Figure 7-2 on page 81](#) follows these steps:

1. Collect a 3D data file once marker 1 comes into view.
2. Open the 3D data file for input.
3. Open a second file for output.
4. Read the first frame of data from the input file.
5. Read the  $i$ th frame of data from the input file.
6. Do the calculations on the  $i$ th frame of data.
7. Write the  $i$ th frame of data to the output file.
8. Close the both input and output files.

```

/*
 * STEP 1
 * Set up a collection for the OPTOTRAK.
 */

```

```
if( OptotrakSetupCollection(
    NUM_MARKERS,          /*Number of markers in the collection. */
    FRAME_RATE,          /*Frequency to collect data frames at. */
    (float)2500.0,       /*Marker frequency for marker max on-time.*/
    30,                  /*Dynamic or Static Threshold value to use.*/
    160,                /*Minimum gain code amplification to use. */
    1,                  /*Stream mode for the data buffers. */
    (float)0.2,         /*Marker Duty Cycle to use. */
    (float)6.5,         /*Voltage to use when turning on markers.*/
    COLLECTION_TIME,    /*Number of seconds of data to collect */
    (float)0.0,         /*Number of seconds to pre-trigger data by.*/
    0 ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Activate the markers.
 */

if( OptotrakActivateMarkers() )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Loop until marker 1 comes into view.
 */

fprintf( stdout, "Waiting for Marker 1\n" );
do
{
    /*
     * Get a frame of 3D data.
     */

    if( DataGetLatest3D( &uFrameNumber, &uElements, &uFlags,
                        pBase3dData ) )
    {
        goto ERROR_EXIT;
    } /* if */
} /* do */
while( pBase3dData[ 0].x < MAX_NEGATIVE );

/*
 * Initialize a file for spooling of the OPTOTRAK 3D data.
 */

if( DataBufferInitializeFile( OPTOTRAK, "C#001.S15" ) )
```

```
{
    goto ERROR_EXIT;
} /* if */

/*
 * Spool data to the previously initialized file.
 */

fprintf( stdout, "Collecting Data File\n" );
if( DataBufferSpoolData( &uSpoolStatus ) )
{
    goto ERROR_EXIT;
} /* if */
fprintf( stdout, "Spool Status: 0x%04x\n", uSpoolStatus );

/*
 * De-activate the markers.
 */

if( OptotrakDeActivateMarkers() )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Shutdown the transputer message passing system.
 */

if( TransputerShutdownSystem() )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 2
 * Open the 3D data file we just collected as our input data file.
 */

fprintf( stdout, "Processing 3D data file\n" );
if( FileOpen( "C#001.S15",
              INPUT_FILE,
              OPEN_READ,
              &nFileItems,
              &nFileSubItems,
              &nFileFrames,
              &fFileFrequency,
              szFileComments,
              &pFileHeader ) )
{
    goto ERROR_EXIT;
} /* if */
```

```
/*
 * STEP 3
 * Open a new file as our output file.
 */

strcpy( szFileComments, "Normalized 3D data file" );
if( FileOpen( "NC#001.S15",
             OUTPUT_FILE,
             OPEN_WRITE,
             &nFileItems,
             &nFileSubItems,
             &lnFileFrames,
             &fFileFrequency,
             szFileComments,
             &pFileHeader ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 4
 * Read the first frame of 3D data from the input file and store it.
 */

if( FileRead( INPUT_FILE, 0L, 1, pBase3dData ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * Read each frame of the input file and subtract the base position.
 * Write the resultant frame to the output file.
 */

for( lnFrameCnt = 0L; lnFrameCnt < lnFileFrames; ++lnFrameCnt )
{

    /*
     * STEP 5
     * Read the current frame from the input file.
     */
    if( FileRead( INPUT_FILE, lnFrameCnt, 1, pInput3dData ) )
    {
        goto ERROR_EXIT;
    } /* if */
    /*
     * STEP 6
     * Perform the subtraction for each in item in the frame.
     */
    for( nItemCnt = 0; nItemCnt < nFileItems; ++nItemCnt )
    {
```



```

/*
 * If the input data for this marker is missing set the *output data
 * to missing as well.
 */
if( pInput3dData[ nItemCnt].x < MAX_NEGATIVE )
{
    pOutput3dData[ nItemCnt].x =
        pOutput3dData[ nItemCnt].y =
            pOutput3dData[ nItemCnt].z = BAD_FLOAT;
} /* if */
/*
 * Else perform the subtraction.
 */
else
{
    pOutput3dData[ nItemCnt].x = pInput3dData[ nItemCnt].x -
        pBase3dData[ nItemCnt].x
    pOutput3dData[ nItemCnt].y = pInput3dData[ nItemCnt].y -
        pBase3dData[ nItemCnt].y
    pOutput3dData[ nItemCnt].z = pInput3dData[ nItemCnt].z -
        pBase3dData[ nItemCnt].z

    } /* else */
} /* for */
/*
 * STEP 7
 * Write the calculated frame to the output file.
 */

if( FileWrite( OUTPUT_FILE, lnFrameCnt, 1, pOutput3dData ) )
{
    goto ERROR_EXIT;
} /* if */
} /* for */
/*
 * STEP 8
 * Close the input and output files.
 */
FileClose( INPUT_FILE );
FileClose( OUTPUT_FILE );
fprintf( stdout, "File processing complete\n" );

```

**Figure 7-2: Processing NDFP Format Files**



## 8 Retrieving Data With a Secondary Host Computer

This section will help you to understand how to run user-written application programs from either the primary or a secondary computer. It includes some basic information on using two computers with the Optotrak System and explains how to retrieve both real-time and buffered data from the secondary host.

You need to purchase the *Secondary PC Interface Kit*, or use an Ethernet connection if using an Optotrak Certus System, to allow you to connect two computers to your system. Generally, the primary computer runs the standard NDI ToolBench software, while the secondary computer runs an API program that retrieves and manipulates the data. It is also possible to run the NDI ToolBench software on both computers to test the system connections. Finally, you can write API programs to run on both computers.

The physical connection of the secondary host computer depends on the devices connected to the Optotrak System. Complete instructions on connecting the secondary host computer, are in the guide *Installing the Secondary PC Interface Kit*.



---

**The primary host program must load the Optotrak System and camera parameters and it must complete these operations before the secondary host sample application programs begin. If the primary host does not complete these operations before the secondary host program starts, it may terminate unsuccessfully.**

---

## 8.1 Retrieving Optotrak System Real-time Data on a Secondary Host Computer

The Program Initialization Code section for a secondary host application differs from a primary host application. The primary host is responsible for downloading the appropriate transputer programs and startup code and loading the camera parameters. The secondary host application simply initializes communications with the Optotrak System.

---

**Note** When the secondary host application invokes the `TransputerInitializeSystem` routine, the flag `OPTO_SECONDARY_HOST_FLAG` must be specified. This flag identifies the application program as a secondary host application to the Optotrak System

---

Once communication has been established, the secondary host application program requests and receives the system status with `OptotrakGetStatus`. This contains the configuration of the current Optotrak System collection set up by the primary host and allows the secondary host to allocate the required amount of memory for the real-time Optotrak System data. The size of memory required depends upon the type of data that is requested. There are several types of real-time data that can be requested from the Optotrak System; 3D data and raw data are the two main formats. Each data type has a unique frame size; calculate the frame size using:

$$\text{Size} = (\text{Number of elements}) \times (\text{size of elements}) \text{ bytes.}$$

To determine the size of the elements see [“Sample Formats” on page 68](#) and [“Real-time Data Types” on page 237](#). There are also methods that can be used to retrieve real-time data, a complete listing can be found in [“Real-time Data Retrieval Routines” on page 157](#).



**Warning!**

---

**The routines for data retrieval copy the data directly into memory allocated by the application program. You must ensure that enough memory is allocated for the frame of real-time data, to prevent the possibility of invalid memory regions being overwritten and possible system crashes.**

---

Use the routine `TransputerShutdownSystem` to remove the secondary host from the message-passing system. This should be the last operation executed by the secondary host.

### Sample Program that Retrieves and Displays Full Raw Data From a Secondary Host

[Figure 8-1 on page 88](#) is an example of a secondary host application program that retrieves and displays full raw data via a routine call. The amount of memory

required for full raw data is calculated using:

$$\text{Size} = (\text{nMarkers} \times \text{sizeof}[\text{FullRawDataType}]) \text{ bytes}$$

This sample code fragment is similar to Sample Program 16 on the API CD. Follow these steps:

1. Initialize communications with the system of processors.
2. Get the current Optotrak System status.
3. Allocate memory for storing the real-time full raw data.
4. Request/receive/display 10 frames of real-time full raw data.
5. Disconnect the application program from the system of processors.

```

/*****
Name:                main
Input Values:
    int
        argc          :Number of command line parameters.
    unsigned char
        *argv[]       :Pointer array to each parameter.

Output Values:
    None.
Return Value:
    None.

*****/
void main( int argc, unsigned char *argv[] )
{
    int
        nFlags,
        nNumSensors
        nNumOdaus,
        nRigidBodies,
        nMarkers,
        nThreshold,
        nMinimumGain,
        nStreamData,
        nSensorCode;
    float
        fVoltage,
        fDutyCycle,
        fCollectionTime,
        fFrameFrequency,
        fMarkerFrequency,
        fPreTriggerTime;

```

```
    unsigned int
        uFlags,
        uElements,
        uFrameCnt,
        uMarkerCnt,
        uSensorCnt,
        uFrameNumber;
    FullRawDataType
        *pFullRawData;
    char
        szNDErrorString[MAX_ERROR_STRING_LENGTH + 1];

/*
 * STEP 1
 * Initialize communications with the system of processors.
 */

if( TransputerInitializeSystem( OPTO_LOG_ERRORS_FLAG |
                               OPTO_SECONDARY_HOST_FLAG ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 2
 * Retrieve the OPTOTRAK System status.
 */
if( OptotrakGetStatus( &nNumSensors, &nNumOdaus, &nRigidBodies, &nMarkers,
                      &nThreshold, &nMinimumGain, &nStreamData, &fDutyCycle,
                      &fFrameFrequency, &fMarkerFrequency, &fVoltage, &fCollectionTime,
                      &fPreTriggerTime, &nFlags ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 3
 * Allocate memory for receiving the real-time OPTOTRAK raw data.
 */

pFullRawData = (FullRawDataType *)calloc(nMarkers, sizeof( FullRaw
DataTypes ) );
if( NULL == pFullRawData )
{
    fprintf( stdout, "Error: Unable to allocate required memory\n"
            );TransputerShutdownSystem();
    exit( 1 );
} /* if */

/*
```

```

* STEP 4
* Get and display ten frames of OPTOTRAK Raw data.
*/

for( uFrameCnt = 0; uFrameCnt < 10; ++uFrameCnt )
{
    /*
    * Get a frame of data.
    */

    if( DataGetLatestRaw( &uFrameNumber, &uElements, &uFlags,
                          pFullRawData ) )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
    * Print out the data.
    */

    fprintf( stdout, "Frame Number: %8u\n", uFrameNumber );
    fprintf( stdout, "Elements      : %8u\n", uElements );
    fprintf( stdout, "Flags          : 0x%04x\n", uFlags );
    for( uMarkerCnt = 0; uMarkerCnt < nMarkers; ++uMarkerCnt )
    {
        /*
        * Print out the current marker number.
        */
        fprintf( stdout, "Marker %u\t\tCentroid Peak  DRC
                    Code\n", uMarkerCnt + 1 );

        /*
        * Print out the information for each sensor.
        */

        for( uSensorCnt = 0; uSensorCnt<NUM_SENSORS; ++uSensorCnt )
        {
            /*
            * Print out the current sensor number.
            */

            fprintf( stdout, "\tSensor %u\t", uSensorCnt + 1 );

            /*
            * Print out the centroid. If it is bad print out *the string
            *'missing'.
            */

            if( pFullRawData[ uMarkerCnt].fCentroid[ uSensorCnt]
                <MAX_NEGATIVE )
            {

```

```
    fprintf( stdout, " missing " );
  } /* if */
else
  {
    fprintf( stdout, "%8.2f ",                pFullRawData[
      uMarkerCnt].fCentroid[ uSensorCnt] );
  } /* else */

/*
 * Print out the rest of this sensor's information.
 */

nSensorCode = pFullRawData[ uMarkerCnt].
  SensorData[ uSensorCnt].ucCode;
  fprintf( stdout, "%4d %4d %s\n",
    pFullRawData[ uMarkerCnt]
      .SensorData[ uSensorCnt].ucPeak,
    pFullRawData[ uMarkerCnt]
      .SensorData[ uSensorCnt].ucDRC,
    pSensorStatusString[ nSensorCode] );
  } /* for */
} /* for */

/*
 * STEP 5
 * Shutdown the message-passing system.
 */

if( !TransputerShutdownSystem() )
{
  goto ERROR_EXIT;
} /* if */

/*
 * Exit the program.
 */

exit( 0 );

ERROR_EXIT:
if( OptotrakGetErrorString( szNDErrorString,
  MAX_ERROR_STRING_LENGTH + 1 ) == 0 )
{
  fprintf( stdout, szNDErrorString );
} /* if */
TransputerShutdownSystem();
exit( 1 );
} /* main */
```

**Figure 8-1: Retrieving Optotrak Full Raw Data on a Secondary Host Computer**



## 8.2 Retrieving Buffered Data on a Secondary Host Computer

When you retrieve buffered data on a secondary host computer, the primary computer has already loaded the appropriate transputer programs and startup code, loaded the correct camera parameters, and configured a collection on the Optotrak System. The secondary host must determine the number of markers before it can retrieve data. After receiving the Optotrak System status, the application program uses this information to allocate the required memory for the real-time data and a file is initialized for spooling the data before the spooling variables are defined.

### Sample Program that Retrieves Buffered Data to a Secondary Host

Figure 8-2 on page 91 is a code fragment that initiates communication with the Optotrak System and retrieves the Optotrak System status. The application waits for marker 1 to come into view and then starts a data file collection. Data are collected until all the data has been spooled or marker 1 goes out of view. For a detailed explanation of the steps that relate to the markers coming into or out of view, see “Advanced Buffered Data Retrieval Without Blocking” on page 36 for a detailed explanation of steps 4 - 7. This code fragment is similar to Sample Program 17 on the API CD.

To retrieve buffered data to a secondary host:

1. Allocate memory for the real-time 3D data.
2. Initialize a data file for spooling the 3D buffered data.
3. Initialize the spooling variables to the required starting state.
4. Wait for marker 1 to come into view.
5. Start spooling buffered data to the secondary host computer.
6. Retrieve real-time 3D data. If marker 1 goes out of view, stop spooling data.
7. Write any incoming buffered data to the appropriate spool destination and check for spool complete.

```

/*
 * STEP 1
 * Allocate memory for receiving the real-time OPTOTRAK 3D data.
 */

p3dData = (Position3d *)calloc( nMarkers, sizeof( Position3d ) );if( NULL
    == p3dData )
{
    fprintf( stdout, "Error: Unable to allocate required memory\n" );

```

```
    TransputerShutdownSystem();
    exit( 1 );
} /* if */

/*
 * STEP 2
 * Initialize a file for spooling of the OPTOTRAK 3D data.
 */

if( DataBufferInitializeFile( OPTOTRAK, "C#001.S17" ) )
{
    goto ERROR_EXIT;
} /* if */

/*
 * STEP 3
 * Initialize the necessary spooling variables.
 */

uSpoolStatus      =
uSpoolComplete    =
uRealtimeDataReady = 0;

/*
 * STEP 4
 * Loop until marker 1 comes into view.
 */

fprintf( stdout, "Waiting for Marker 1\n" );
do
{
    /*
     * Get a frame of 3D data.
     */

    if( DataGetLatest3D( &uFrameNumber, &uElements, &uFlags, p3dData ) )
    {
        goto ERROR_EXIT;
    } /* if */
} /* do */
while( p3dData[ 0 ].x < MAX_NEGATIVE );

/*
 * STEP 5
 * Start the OPTOTRAK spooling data to us.
 */

if ( DataBufferStart() )
```

```

{
    goto ERROR_EXIT;
} /* if */
fprintf( stdout, "Collecting Data File\n" );

/*
 * Loop around spooling data to file until marker 1 goes out of view.
 */

do
{
    /*
     * STEP 6
     * Get a frame of 3D data.
     */

    if( DataGetLatest3D(&uFrameNumber, &uElements, &uFlags, p3dData ) )
    {
        goto ERROR_EXIT;
    } /* if */

    /*
     * Check to see if marker 1 is out of view and stop the OPTOTRAK
     * from spooling data if this is the case.
     */

    if( p3dData[ 0].x < MAX_NEGATIVE )
    {
        if( DataBufferStop() )
        {
            goto ERROR_EXIT;
        } /* if */
    } /* if */

    /*
     * STEP 7
     * Write data if there is any to write.
     */

    if( DataBufferWriteData( &uRealtimeDataReady, &uSpoolComplete,
                            &uSpoolStatus ) )

    {
        goto ERROR_EXIT;
    } /* if */
} /* do */
while( !uSpoolComplete );
fprintf( stdout, "Spool Status: 0x%04x\n", uSpoolStatus );

```

**Figure 8-2: Retrieving Optotrak Buffered Data on a Secondary Host Computer**



## 9 Optotrak API Routines

### 9.1 Overview

This section is your reference for all the routines in the Optotrak API. The function, inputs and outputs are described. The routines are grouped by functionality in this order:

- Optotrak Specific
- Optotrak Device Handles
- ODAU Specific
- Real-time Data Retrieval
- Buffered Data Retrieval
- Rigid Body Specific
- Rigid Body Related
- File Processing
- Registration & Alignment

Within each function, the routines are listed alphabetically. A summary of all the routines is given in [Table 9-1 on page 96](#).

Most of the routines in the Optotrak API return an integer value that indicates whether the routine has been completed successfully. The application program uses this value to decide if the program should continue or exit. If a routine returns a non-zero value, then an error has occurred in the execution of that routine. If a routine returns `NDI_NO_ERROR_CODE` (zero) then the routine has completed successfully.

If a routine returns unsuccessfully, then the API routine `OptotrakGetErrorString` can be used to determine the last error. To see all of the error messages:

1. specify the flag `OPTO_LOG_ERRORS_FLAG` in the routine `TransputerInitializeSystem`
2. look in the file `opto.err`

A complete list of the error messages is provided in [“Error Messages and Constants” on page 273](#). See [“Flags and Settings Associated with Rigid Bodies” on page 295](#) for details on rigid body error flags and settings.

## 9.2 Optotrak Certus Specific Routines

### 9.2.1 Device Handles

The Optotrak Certus System Control Unit is capable of automatically identifying tools and strobers that are connected to the system. Each tool or strober is assigned a unique ID, called a device handle, which persists until the tool is unplugged, the application releases it, or the session is ended. Device handles, which are defined in “[Optotrak Device Handle Routines](#)” on page 135, may assume any of the following states:

**FREE** A device handle in the FREE state is not associated with a device. Initially, all device handles are in the FREE state. A device handle can also enter the FREE state from the UNOCCUPIED state when instructed by the application program.

**OCCUPIED** A device handle enters the OCCUPIED state when the Optotrak Certus System Control Unit detects that a new device has been connected to the system. The Optotrak Certus System Control Unit is responsible for informing the application program that a device handle has entered the OCCUPIED state.

**INITIALIZED** When the application program recognizes an OCCUPIED device, the associated device handle changes to the INITIALIZED state. While a device handle is in the INITIALIZED state, the Optotrak Certus System Control Unit may receive commands from the host computer to change the configuration of the device handle. Yet, the device handle will not be able to activate markers, monitor switches, or set visible LED information.

**ENABLED** When instructed by the application program, a device handle enters the ENABLED state from the INITIALIZED state. Devices in the ENABLED state are able to activate markers, monitor switches, and set visible LED information.

**UNOCCUPIED** A device handle enters the UNOCCUPIED state when the Optotrak Certus System Control Unit detects that a device has been disconnected from the system. The Optotrak Certus System Control Unit is responsible for informing the application program that a device handle has entered the UNOCCUPIED state.

If the Optotrak Certus System Control Unit detects that an INITIALIZED or ENABLED device has been disconnected from the system, it will mark the device handle as UNOCCUPIED. The host application program will then be responsible for instructing the Optotrak Certus System Control Unit to change the status of the device handle from UNOCCUPIED to FREE.

## 9.2.2 Device Handle Properties

Information about tools and strobers is conveyed using device handle properties. Device handle properties may describe the physical attributes of a device or describe parameters set by the application.

Currently supported properties include:

**DH\_PROPERTY\_NAME** identifies the device name.

**DH\_PROPERTY\_MARKERSTOFIRE** identifies the number of markers to activate.

**DH\_PROPERTY\_MAXMARKERS** identifies the maximum number of markers on the device.

**DH\_PROPERTY\_STARTMARKERPERIOD** identifies the index of the first marker activated by the device within the system activation order.

**DH\_PROPERTY\_SWITCHES** identifies the number of switches available on the device.

**DH\_PROPERTY\_VLEDS** identifies the number of visible LEDs on the device.

**DH\_PROPERTY\_PORT** identifies the System Control Unit strober port to which the device is plugged into.

**DH\_PROPERTY\_ORDER** identifies the index within the port.

**DH\_PROPERTY\_SUBPORT** identifies the sub-port to which the device is plugged into.

**DH\_PROPERTY\_FIRINGSEQUENCE** identifies the marker activation sequence.

**DH\_PROPERTY\_HAS\_ROM** identifies whether or not the device contains a Read Only Memory (ROM).

**DH\_PROPERTY\_TOOLPORTS** identifies the number of ports on a strober.

**DH\_PROPERTY\_3020\_CAPABILITY** identifies whether or not the system has the ability to activate Optotrak 3020 markers.

**DH\_PROPERTY\_3020\_MARKERSTOFIRE** identifies the number of Optotrak 3020 System markers to activate.

**DH\_PROPERTY\_3020\_STARTMARKERPERIOD** identifies the index of the first marker activated by the device within the system activation order.

**DH\_PROPERTY\_STATUS** identifies the status of the device handle.

---

**Note** Additional properties may be added to the Optotrak Certus Firmware in a future revision.

---

## 9.3 Optotrak API Routines

**Table 9-1: Optotrak API Routines**

| Functionality          | Routine Name                     | Page     |
|------------------------|----------------------------------|----------|
| Optotrak               | TransputerDetermineSystemCfg     | page 101 |
|                        | TransputerInitializeSystem       | page 102 |
|                        | TransputerLoadSystem             | page 103 |
|                        | TransputerShutdownSystem         | page 105 |
|                        | OptotrakActivateMarkers          | page 105 |
|                        | OptotrakChangeCameraFOR          | page 106 |
|                        | OptotrakConvertRawTo3D           | page 108 |
|                        | OptotrakConvertTransforms        | page 109 |
|                        | OptotrakDeActivateMarkers        | page 110 |
|                        | OptotrakGetCameraParameterStatus | page 111 |
|                        | OptotrakGetErrorString           | page 112 |
|                        | OptotrakGetNodeInfo              | page 113 |
|                        | OptotrakGetStatus                | page 114 |
|                        | OptotrakLoadCameraParameters     | page 116 |
|                        | OptotrakSaveCollectionToFile     | page 117 |
|                        | OptotrakSetCameraParameters      | page 118 |
|                        | OptotrakSetProcessingFlags       | page 120 |
|                        | OptotrakSetStroberPortTable      | page 123 |
|                        | OptotrakSetupCollection          | page 124 |
|                        | OptotrakSetupCollectionFromFile  | page 128 |
| OptotrakStopCollection | page 133                         |          |



**Table 9-1: Optotrak API Routines**

| <b>Functionality</b>       | <b>Routine Name</b>                     | <b>Page</b>              |
|----------------------------|-----------------------------------------|--------------------------|
| Optotrak<br>Device Handles | OptotrakDeviceHandleEnable              | <a href="#">page 135</a> |
|                            | OptotrakDeviceHandleFree                | <a href="#">page 135</a> |
|                            | OptotrakDeviceHandleGetNumberProperties | <a href="#">page 136</a> |
|                            | OptotrakDeviceHandleGetProperties       | <a href="#">page 137</a> |
|                            | OptotrakDeviceHandleGetProperty         | <a href="#">page 138</a> |
|                            | OptotrakDeviceHandleSetBeeper           | <a href="#">page 139</a> |
|                            | OptotrakDeviceHandleSetProperties       | <a href="#">page 139</a> |
|                            | OptotrakDeviceHandleSetVisibleLED       | <a href="#">page 140</a> |
|                            | OptotrakGetDeviceHandles                | <a href="#">page 141</a> |
|                            | OptotrakGetNumberDeviceHandles          | <a href="#">page 142</a> |

**Table 9-1: Optotrak API Routines**

| <b>Functionality</b>     | <b>Routine Name</b>         | <b>Page</b> |
|--------------------------|-----------------------------|-------------|
| ODAU Related Routines    | OdauGetStatus               | page 143    |
|                          | OdauSaveCollectionToFile    | page 144    |
|                          | OdauSetAnalogOutputs        | page 145    |
|                          | OdauSetDigitalOutputs       | page 146    |
|                          | OdauSetTimer                | page 148    |
|                          | OdauSetupCollection         | page 149    |
|                          | OdauSetupCollectionFromFile | page 152    |
| Real-time Data Retrieval | DataGetLatestCentroid       | page 157    |
|                          | DataGetLatest3D             | page 158    |
|                          | DataGetLatestOdauRaw        | page 160    |
|                          | DataGetLatestRaw            | page 161    |
|                          | DataGetLatestTransforms     | page 163    |
|                          | DataGetLatestTransforms2    | page 165    |
|                          | DataGetNext3D               | page 167    |
|                          | DataGetNextCentroid         | page 168    |
|                          | DataGetNextOdauRaw          | page 170    |
|                          | DataGetNextRaw              | page 171    |
|                          | DataGetNextTransforms       | page 172    |
|                          | DataGetNextTransforms2      | page 174    |
|                          | DataIsReady                 | page 177    |
|                          | DataReceiveLatest3D         | page 178    |
|                          | DataReceiveLatestCentroid   | page 179    |
| DataReceiveLatestOdauRaw | page 180                    |             |

**Table 9-1: Optotrak API Routines**

| <b>Functionality</b>                 | <b>Routine Name</b>          | <b>Page</b> |
|--------------------------------------|------------------------------|-------------|
| Real-time Data Retrieval, continued. | DataReceiveLatestRaw         | page 182    |
|                                      | DataReceiveLatestTransforms  | page 183    |
|                                      | DataReceiveLatestTransforms2 | page 185    |
|                                      | ReceiveLatestData            | page 187    |
|                                      | RetrieveSwitchData           | page 187    |
|                                      | RequestLatest3D              | page 188    |
|                                      | RequestLatestCentroid        | page 189    |
|                                      | RequestLatestOdauRaw         | page 190    |
|                                      | RequestLatestRaw             | page 191    |
|                                      | RequestLatestTransforms      | page 192    |
|                                      | RequestNext3D                | page 193    |
|                                      | RequestNextCentroid          | page 194    |
|                                      | RequestNextOdauRaw           | page 195    |
|                                      | RequestNextRaw               | page 196    |
| RequestNextTransforms                | page 197                     |             |
| Buffered Data Retrieval              | DataBufferAbortSpooling      | page 198    |
|                                      | DataBufferInitializeFile     | page 199    |
|                                      | DataBufferInitializeMem      | page 200    |
|                                      | DataBufferSpoolData          | page 201    |
|                                      | DataBufferStart              | page 202    |
|                                      | DataBufferStop               | page 203    |
|                                      | DataBufferWriteData          | page 204    |

**Table 9-1: Optotrak API Routines**

| <b>Functionality</b>     | <b>Routine Name</b>          | <b>Page</b> |
|--------------------------|------------------------------|-------------|
| Rigid Body Specific      | RigidBodyAdd                 | page 206    |
|                          | RigidBodyAddFromDeviceHandle | page 208    |
|                          | RigidBodyAddFromFile         | page 208    |
|                          | RigidBodyChangeFOR           | page 210    |
|                          | RigidBodyChangeSettings      | page 211    |
|                          | RigidBodyDelete              | page 213    |
| Rigid Body Related       | CombineXfrms                 | page 214    |
|                          | CvtQuatToRotationMatrix      | page 215    |
|                          | CvtRotationMatrixToQuat      | page 216    |
|                          | DetermineEuler               | page 216    |
|                          | DetermineR                   | page 217    |
|                          | InverseXfrm                  | page 218    |
|                          | TransformPoint               | page 219    |
| File Processing          | FileClose                    | page 220    |
|                          | FileCloseAll                 | page 221    |
|                          | FileConvert                  | page 221    |
|                          | FileOpen                     | page 222    |
|                          | FileOpenAll                  | page 224    |
|                          | FileRead                     | page 226    |
|                          | FileReadAll                  | page 227    |
|                          | FileWrite                    | page 228    |
|                          | FileWriteAll                 | page 230    |
| Registration & Alignment | nOptotrakAlignSystem         | page 231    |
|                          | nOptotrakCalibrigSystem      | page 232    |
|                          | nOptotrakRegisterSystem      | page 234    |

## 9.4 Optotrak Specific Routines

See “Retrieving Real-time Optotrak Data” on page 29, “Retrieving Buffered Optotrak Data” on page 34, and “Initializing, Retrieving System Status and Exiting from the Optotrak System” on page 24 for general a discussions of Optotrak System routines.

### 9.4.1 TransputerDetermineSystemCfg

#### Function

Determines the current configuration of the networks of internal processors in the Optotrak System and writes the information to the default network information file, `system.nif`. This routine also retrieves the camera parameters from each Position Sensor and stores that information in the default camera parameter file, `standard.cam`.

#### Prototype

```
int TransputerDetermineSystemCfg( char *pszInputLogFile )
```

#### Parameters

`pszInputLogFile` points to a null terminated string containing the name of a log file that can be used to log any status information, messages and errors encountered while determining and writing the network information file. `pszInputLogFile` can be set to NULL if logging is not required.

#### Description

`TransputerDetermineSystemCfg` determines the current configuration of the Optotrak System network of processors, in its default operation mode, and writes the information to the default network information file, `system.nif`. It also stores default camera parameters in `standard.cam`.

`TransputerDetermineSystemCfg` is functionally equivalent to the preliminary initialization of the Optotrak System done with the command-line utility programs `optset32.exe` (Windows NT/2000/XP) or `buildnif` (Linux, SGI, and Sun). This routine should be called once at the beginning of the application program if the preliminary Optotrak System initialization has not been done with the command-line utility programs.

After an application program calls `TransputerDetermineSystemCfg`, it must call `TransputerLoadSystem` to load the network of processors. You only need to invoke `TransputerLoadSystem` once after power-up and determining the system configuration. After the network of processors has been loaded, the application

program can connect and disconnect from the processor system as desired by using the routines `TransputerInitializeSystem` and `TransputerShutdownSystem` respectively.



---

**You must re-determine the system configuration and re-initialize the Optotrak System if either the cabling connections to the communication ports at the back of the System Control Unit or the Position Sensor are changed, or if the order of the Position Sensors is changed. The cabling connection has changed once a cable is disconnected. It is a good practice to re-initialize your setup and re-determine the system configuration if you are not certain whether the cabling connections have been changed.**

---

In an alternative operation mode, `TransputerDetermineSystemCfg` can store the network information file information internally within the API instead of externally in the default system network information file, `system.nif`. To do this, set the flag `OPTO_USE_INTERNAL_NIF` with a call to `OptotrakSetProcessingFlags`, prior to calling `TransputerDetermineSystemCfg`. The subsequent call to `TransputerLoadSystem` will signal access to the internally stored network information file information by passing an empty string argument.

### See Also

`TransputerLoadSystem`, `TransputerInitializeSystem`, `TransputerShutdownSystem`, `OptotrakSetProcessingFlags`

## 9.4.2 `TransputerInitializeSystem`

### Function

Initializes the system so the application program can communicate with the Optotrak System.

Verifies that the code running on all Optotrak components is compatible with the current version of the API.

### Prototype

```
int TransputerInitializeSystem( unsigned int uFlags )
```

### Parameters

`uFlags` specifies the mode in which the application program is to connect to the Optotrak System. This parameter passes information using bit flags. To pass two or more flags to the routine, separate them with the logical OR operator “|”.

Values:

**OPTO\_LOG\_ERRORS\_FLAG** writes error information to the file `opto.err` in the current directory.

**OPTO\_LOG\_MESSAGES\_FLAG** writes a detailed summary of all messages sent and received from the Optotrak System to the file `opto.err` in the current directory.

**OPTO\_SECONDARY\_HOST\_FLAG** registers the computer as the secondary host computer. If this flag is not set the computer is registered as the primary host computer.

## Description

`TransputerInitializeSystem` initializes the system so the application program can communicate with the Optotrak System. It verifies that the code running on all Optotrak components is compatible with the API. This routine must be called if the application program has just downloaded system code with the routine `TransputerLoadSystem`. This routine must also be called to re-establish communication with the Optotrak System if the `TransputerShutdownSystem` routine was previously called.



**Warning!**

---

**You must re-determine the system configuration and re-initialize the Optotrak System if either the cabling connections to the communication ports at the back of the System Control Unit or the Position Sensor are changed, or if the order of the Position Sensors is changed. The cabling connection has changed once a cable is disconnected. It is a good practice to re-initialize your setup and re-determine the system configuration if you are not certain whether the cabling connections have been changed.**

---

See Sample Program 1 on the API CD for an example that uses this routine.

## See Also

`TransputerLoadSystem`, `TransputerShutdownSystem`

### 9.4.3 `TransputerLoadSystem`

#### Function

Loads the Optotrak System network of processors with the system code according to the specified Network Information File (`.nif`).

## Prototype

```
int TransputerLoadSystem( char *pszNifFile )
```

## Parameters

**pszNifFile** is a pointer to a null terminated string that specifies the network information file to use when loading the network of processors in the Optotrak System. Do not include the file extension when specifying the file name in this parameter. For example, use “system” and not “system.nif”.

If **pszNifFile** is set to NULL, the standard network information file, **system.nif**, will be used. If the API flag **OPTO\_USE\_INTERNAL\_NIF** is set, the internal network configuration will be used. The network configuration can be stored internally by calling **TransputerDetermineSystemCfg** with the API set to use internal network configuration information.

## Description

**TransputerLoadSystem** downloads the appropriate transputer program and startup code according to the system configuration specified in the network information file. Once this routine has completed, the application program must call **TransputerInitializeSystem** to establish a communications connection with the Optotrak System.

You only need to invoke **TransputerLoadSystem** once after power-up. After the network of processors has been loaded successfully, the application program can connect and disconnect from the processor system as desired by using the routines **TransputerInitializeSystem** and **TransputerShutdownSystem** respectively.

---

**Note** After a call to the routine **TransputerLoadSystem**, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may be generated.

---



---

**You must re-determine the system configuration and re-initialize the Optotrak System if either the cabling connections to the communication ports at the back of the System Control Unit or the Position Sensor are changed, or if the order of the Position Sensors is changed. The cabling connection has changed once a cable is disconnected. It is a good practice to re-initialize your setup and re-determine the system configuration if you are not certain whether the cabling connections have been changed.**

---

See Sample Program 1 on the API CD for an example of code that uses this routine.



## See Also

TransputerDetermineSystemCfg, TransputerInitializeSystem,  
TransputerShutdownSystem

### 9.4.4 TransputerShutdownSystem

#### Function

Informs the Optotrak System that the application program will no longer be communicating with it.

#### Prototype

```
int TransputerShutdownSystem( void )
```

#### Parameters

None.

#### Description

TransputerShutdownSystem disconnects the application program from the Optotrak System. This routine should be called either before an application program exits, or once it has finished using the Optotrak System. You must call TransputerInitializeSystem if TransputerShutdownSystem has been invoked and the program needs to resume communications with the Optotrak System.

See Sample Program 1 on the API CD for an example of code that uses this routine.

## See Also

TransputerLoadSystem, TransputerInitializeSystem

### 9.4.5 OptotrakActivateMarkers

#### Function

Activates the IRED markers.

#### Prototype

```
int OptotrakActivateMarkers( void )
```

## Parameters

None.

## Description

OptotrakActivateMarkers turns the markers on.

Markers are automatically activated after the routines OptotrakSetupCollection or OptotrakSetupCollectionFromFile are called without passing the flag OPTOTRAK\_NO\_FIRE\_MARKERS\_FLAG. However, if the routine OptotrakDeActivateMarkers has been called since the last system download, then the routines will not automatically activate the markers.

It is best not to rely on the automatic activation of markers. You should call OptotrakActivateMarkers before the beginning of a collection and use OptotrakDeActivateMarkers when the collection has finished.

If there is a significant amount of elapsed time between a call to OptotrakSetupCollection or OptotrakSetupCollectionFromFile and the time when the markers are viewed, then it is a good practice to use OptotrakDeActivateMarkers after calling the collection setup routines. This reduces the wear on the markers. The markers should then be activated only when needed using OptotrakActivateMarkers.

There may be a short propagation delay between when OptotrakActivateMarkers is called and when the markers become activated. Some initial frames may be reported as MISSING if data is requested immediately after a call to this routine. Use a sleep routine of two seconds after calling OptotrakActivateMarkers before requesting data. The sleep time may need to be extended to three to five seconds for markers on axon strobbers.

See Sample Program 2 on the API CD for an example of code that uses this routine.

## See Also

OptotrakDeActivateMarkers, OptotrakSetupCollectionFromFile, OptotrakSetupCollection

### 9.4.6 OptotrakChangeCameraFOR

#### Function

Changes the Optotrak System's measurement coordinate system.

## Prototype

```
int OptotrakChangeCameraFOR( char      *pszInputCamFile,
                             int       nNumMarkers,
                             Position3d *pdtAlignedPositions,
                             Position3d *pdtMeasuredPositions,
                             char      *pszAlignedCamFile,
                             Position3d *pdt3dErrors,
                             float     *pfRmsError )
```

## Parameters

**pszInputCamFile** points to a null terminated string specifying the name of the current camera parameter file. If it is NULL then the default file, standard.cam, is used.

**nNumMarkers** specifies the number of markers used.

**pdtAlignedPositions** points to an array containing the positions of the markers in the desired coordinate system.

**pdtMeasuredPositions** points to an array containing the measured positions of the markers in the current coordinate system.

**pszAlignedCamFile** points to a null terminated string that specifies the name of the new camera parameter file.

**pdt3dErrors** points to an array that contains the 3D errors of the alignment transformation.

**pfRmsError** points to the RMS distance error of the alignment transformation.

## Description

OptotrakChangeCameraFOR changes the Optotrak System's measurement coordinate system. This routine produces a new camera parameter file that puts the markers in aligned positions. The new file is produced by comparing a measured and aligned (or desired) view.

A typical use of this routine applies to rigid bodies, defined by markers, and placed in the field of view. The routine takes the measured positions of the markers and sets the desired positions from the rigid body file description of the object.

You must call OptotrakLoadCameraParameters after calling OptotrakChangeCameraFOR. Use the new camera parameter filename specified in pszAlignedCamFile with the new coordinate system.

You may also use the `nOptotrakAlignSystem` routine to change the Optotrak System's measurement coordinate system but it requires a data file to be collected using a rigid body. `nOptotrakAlignSystem` is intended for users who are familiar with the 'align' command-line utility.

See Sample Program 18 on the API CD for an example that uses this routine.

## See Also

`OptotrakLoadCameraParameters`, `nOptotrakAlignSystem`

## 9.4.7 OptotrakConvertRawTo3D

### Function

Converts a frame of raw data to its corresponding 3D positions, according to the current camera parameters.

### Prototype

```
int OptotrakConvertRawTo3D( unsigned int *puElements,  
                           void          *pSensorReadings,  
                           Position3d   *pdt3DPositions )
```

### Parameters

`puElements` is set to the number of data elements (markers) that were converted.

`pSensorReadings` points to the frame of raw data.

`pdt3DPositions` points to memory set aside by the application program that will receive the converted 3D positions.

### Description

`OptotrakConvertRawTo3D` can be used to convert previously obtained raw position data to their corresponding 3D positions. The camera parameters currently loaded in the system are used for the conversion. You can obtain real-time raw data frames at high frequencies, store them, and convert them to 3D positions at a later time with this routine.



**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 241. This routine copies the 3D positions directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

See Sample Program 19 on the API CD for an example that uses this routine.

## See Also

DataGetLatestRaw

## 9.4.8 OptotrakConvertTransforms

### Function

Transforms an input frame of 3D position data to the corresponding 6D data, according to the current rigid body definitions.

### Prototype

```
int OptotrakConvertTransforms( unsigned int          *puElements,
                               struct OptotrakRigidStruct *pDataDest6D,
                               Position3d             *pDataDest3D )
```

### Parameters

**puElements** is set to the number of data elements (rigid bodies) that were used in the transformation.

**pDataDest6D** points to the memory set aside by the application program that will contain the transformed 6D data.

**pDataDest3D** points to the array of 3D positions to be transformed.

### Description

OptotrakConvertTransforms converts previously obtained 3D positions to their corresponding 6D data. The rigid body definitions currently loaded in the system are used to determine the transformations. The routine also returns the number of rigid body transformations in the frame. This routine, in conjunction with OptotrakConvertRawTo3D, allows an application program to obtain real-time

centroid data frames at high frequencies, store them, and convert them to 3D positions and 6D data at a later time.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size (See “[Optotrak 3D Data](#)” on page 241 and “[Optotrak Rigid Body Transformation Data](#)” on page 242.) This routine copies the transformation data directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

See Sample Program 19 on the API CD for an example of code that uses this routine.

### See Also

DataGetLatestRaw, OptotrakConvertRawTo3D, OptotrakSetProcessingFlags

## 9.4.9 OptotrakDeActivateMarkers

### Function

De-activates the IRED markers.

### Prototype

```
int OptotrakDeActivateMarkers( void )
```

### Parameters

None.

### Description

OptotrakDeActivateMarkers turns the markers off.

De-activate the markers after using either `OptotrakSetupCollectionFromFile` or `OptotrakSetupCollection` if there is a delay between setting up the collection and viewing the markers. This reduces the wear on the markers. Activate the markers when needed using `OptotrakActivateMarkers`.

See Sample Program 2 on the API CD for an example of code that uses this routine.

### See Also

OptotrakActivateMarkers

## 9.4.10 OptotrakGetCameraParameterStatus

### Function

Lists all of the camera parameter sets available in a camera parameter file, as well as the current camera parameter set marker type, wavelength type and model type.

### Prototype

```
int OptotrakGetCameraParameterStatus( int *pnCurrentMarkerType,
                                     int *pnCurrentWaveLength,
                                     int *pnCurrentModelType,
                                     char *szStatus,
                                     int nStatusLength )
```

### Parameters

**pnCurrentMarkerType** is set to the current camera parameter set marker type.

**pnCurrentWaveLength** is set to the current camera parameter set wavelength type.

**pnCurrentModelType** is set to the current camera parameter set lens model type.

**szStatus** points to a character buffer that will contain the identifying information for all of the available camera parameter sets as a null terminated string.

---

**Note** Since the number of parameter sets contained within an extended camera parameter file has no defined upper limit, the buffer size should be large enough to handle large numbers of sets, otherwise the status string is truncated to fit the buffer.

---

**nStatusLength** is the length of the buffer **szStatus**.

### Description

**OptotrakGetCameraParameterStatus** is an API routine that determines all of the camera parameter sets available to an application program in an extended camera parameter file. The extended file has been loaded with the API routine **OptotrakLoadCameraParameters**. The parameter set information is returned as a null terminated string. A typical example is given below for one camera with 3 lens parameter sets:

```
C3-03269
  Model 0 NDI Metal Base Marker 950 nm
  Model 1 NDI Metal Base Marker 950 nm
  Model 2 NDI Metal Base Marker 950 nm
```

The routine can return the marker type, marker wavelength and the lens model type of the parameter set currently being used by the API. Application programs can use this information to select another camera parameter set by calling the API routine `OptotrakSetCameraParameters` with the desired marker type and wavelength and lens model type as arguments.

## See Also

`OptotrakLoadCameraParameters`, `OptotrakSetCameraParameters`

### 9.4.11 `OptotrakGetErrorString`

#### Function

Obtains the most recent API error message.

#### Prototype

```
int OptotrakGetErrorString( char* szErrorString, int nBufferSize )
```

#### Parameters

`szErrorString` points to a character buffer that contains the error message as a null terminated string.

`nBufferSize` is the size of the buffer. Set `nBufferSize` to `MAX_ERROR_STRING_LENGTH + 1` to ensure that the buffer will contain the largest message. If the buffer is too small, the message is truncated.

#### Description

`OptotrakGetErrorString` obtains error information from the API whenever an API routine returns unsuccessfully. If several routines fail in succession, only the last error message of the last failed routine will be available.

`OptotrakGetErrorString` is a replacement for the global variable `szNDErrorString`, which was used in early versions of the API to obtain error messages directly. A complete list of the error messages is provided in [“Error Messages and Constants” on page 273](#).

See Sample Program 1 on the API CD for an example of code that uses this routine.

## See Also

None.



## 9.4.12 OptotrakGetNodeInfo

### Function

Requests and returns the node information for a specified node. A node is any device connected in the Optotrak System.

### Prototype

```
int OptotrakGetNodeInfo( int nNodeId,
                        struct OptoNodeInfoStruct *pdtNodeInfo )
```

### Parameters

**nNodeId** is the node identification number for the node from which information is being requested. The valid range is from 0 to 6. You may also use the pre-defined node identifiers corresponding to the node number:

**Table 9-2: Predefined Node Identifiers**

| Node Name       | Node Number |
|-----------------|-------------|
| OPTOTRAK        | 0           |
| DATA_PROPRIETOR | 1           |
| ODAU1           | 2           |
| ODAU2           | 3           |
| ODAU3           | 4           |
| ODAU4           | 5           |
| SENSOR_PROP1    | 6           |

**pdtNodeInfo** points to a structure that will receive the requested information. The structure is defined as:

```
struct OptoNodeInfoStruct
{
    unsigned long ulHWType,
                  ulHWRev,
                  ulMemSize;
    char          szFreezeId[ 80],
                  szSwDesc[ 80],
                  szSerialNo[ 32];
};
```

The following hardware types are defined for the element `ulHWType`:

- OPTO\_HW\_TYPE\_SENSOR
- OPTO\_HW\_TYPE\_SU
- OPTO\_HW\_TYPE\_ODAU
- OPTO\_HW\_TYPE\_REALTIME
- OPTO\_HW\_TYPE\_CERTUS\_SENSOR
- OPTO\_HW\_TYPE\_CERTUS\_SU

---

**Note** A System Control Unit equipped with a Real-time option will be reported as `OPTO_HW_TYPE_REALTIME` instead of `OPTO_HW_TYPE_SU`.

---

### Description

OptotrakGetNodeInfo obtains node information from the API for a specified node, such as ODAU1.

### See Also

None.

## 9.4.13 OptotrakGetStatus

### Function

Retrieves the current status information of the Optotrak System.

### Prototype

```
int OptotrakGetStatus( int    *pnNumSensors,
                      int    *pnNumOdaus,
                      int    *pnNumRigidBodies,
                      int    *pnMarkers,
                      float  *pfFrameFrequency,
                      float  *pfMarkerFrequency,
                      int    *pnThreshold,
                      int    *pnMinimumGain,
                      int    *pnStreamData,
                      float  *pfDutyCycle,
                      float  *pfVoltage,
                      float  *pfCollectionTime,
                      float  *pfPreTriggerTime,
                      int    *pnFlags )
```

---

## Parameters

---

**Note** You can pass a NULL pointer in a parameter if you do not require the data from it.

---

**pnNumSensors** is set to the number of sensors that are currently in the Optotrak System.

**pnNumOdaus** is set to the number of ODAUs that are currently active in the Optotrak System.

**pnNumRigidBodies** is set to the number of rigid bodies for which the Optotrak System determines and returns transformations when an application program requests the latest set of rigid body transformations.

**pnMarkers** is set to the number of markers in the current collection running on the Optotrak System.

**pfFrameFrequency** is set to the frame frequency at which the Optotrak System is collecting data.

**pfMarkerFrequency** is set to the marker frequency the Optotrak System is using to strobe the individual markers.

**pnThreshold** is set to the threshold value the Optotrak System is using to process sensor pixel data.

**pnMinimumGain** is set to the minimum gain value for the current Optotrak System collection.

**pnStreamData** indicates whether the Optotrak System is to send buffered data back by request only, or to automatically send buffered data back once data spooling is initiated.

Values:

\***pnStreamData** = 0: Send buffered data back at request only.

\***pnStreamData** = 1: Send buffered data back automatically.

**pfDutyCycle** is set to the duty cycle value for the current Optotrak System collection.

**pfVoltage** is set to the voltage the Optotrak System is using when turning on the IRED markers.

**pfCollectionTime** is the duration, in seconds, for buffered data collections.

**pfPreTriggerTime** is not currently supported and will be 0.

**pnFlags** is set to the flag values used when the last Optotrak System collection was configured. This includes the version of the System Control Unit. If the System

Control Unit is an Optotrak Certus, `pnFlags` returns the value `OPTOTRAK_CERTUS_FLAG`. If the System Control Unit is an Optotrak 3020, `pnFlags` returns `OPTOTRAK_3020_FLAG`.

## Description

`OptotrakGetStatus` requests and receives the current status information of the Optotrak System. This routine returns information about the current collection running on the Optotrak System, as well as general system information (e.g. the number of sensors and ODAU devices in the system, the Optotrak System variant, etc.).

See Sample Program 1 on the API CD for an example of code that uses this routine.

## See Also

`OdauGetStatus`, `OptotrakSetupCollection`

### 9.4.14 OptotrakLoadCameraParameters

#### Function

Loads the parameter settings specified in the camera parameter file to the Optotrak System.

#### Prototype

```
int OptotrakLoadCameraParameters( char *pszCamFile )
```

#### Parameters

`pszCamFile` points to a null terminated string that specifies which camera parameter file to use. If `NULL` is passed into this parameter, the default camera parameter file “standard.cam” will be used. Do not include the file extension when specifying the camera parameter file to this routine. For example, use “standard” and not “standard.cam”.

#### Description

`OptotrakLoadCameraParameters` reads the contents of the specified camera parameter file, packs the camera parameters into the appropriate format, and sends the information to the Optotrak System.



---

**You must ensure that the Optotrak System is using the correct set of camera parameters. Without the camera parameters the Optotrak System is unable to properly convert the raw data into 3D position data.**

---

Camera parameter files for a single Optotrak System are created with `TransputerDetermineSystemCfg`.

Camera parameter files for single Optotrak Systems can also be created using the alignment routines, or by using a command line utility. For information on registration and alignment routines, please see [“Registration and Alignment Routines” on page 231](#). For instructions on how to create camera parameter files from the command line, please see [“Initializing the Optotrak System” on page 15](#).

To create camera parameter files for multiple systems, use the registration facility in the NDI ToolBench software (see “Using the Calibrate Interface Window” in the *“NDI ToolBench User Guide”*) or use the `nOptotrakRegisterSystem` routine. The command line utility `register.exe` or `calibrig.exe` may also be used to create camera parameter files for multiple systems.

The routine `OptotrakLoadCameraParameters` searches for the specified camera parameter file in a number of directories. The current directory is searched first. If the file is not in the current directory, the routine searches the “realtime” subdirectory under the path specified in the environment variables `ND_USER_DIR` and `ND_DIR`.

As an example, suppose `ND_USER_DIR` was set to `c:\ndigital` and `ND_DIR` was set to `d:\ndigital`. `OptotrakLoadCameraParameters` would first look for the file in the current directory, then in the directory `c:\ndigital\realtime`, and finally in the directory `d:\ndigital\realtime`.

See Sample Program 1 on the API CD for an example of code that uses this routine.

## See Also

None.

## 9.4.15 OptotrakSaveCollectionToFile

### Function

Saves the current collection parameters to a collection parameter file in ASCII format that can be read by `OptotrakSetupCollectionFromFile` and other programs.

## Prototype

```
int OptotrakSaveCollectionToFile( char *pszCollectFile )
```

## Parameters

**pszCollectFile** points to a null terminated string specifying the name of the collection parameter file.

## Description

Use this routine to save the current collection parameter values to an ASCII file.

## See Also

OptotrakSetupCollectionFromFile

### 9.4.16 OptotrakSetCameraParameters

## Function

Specifies the camera parameter set used by the API according to the given marker type, wavelength type, and model type.

## Prototype

```
int OptotrakSetCameraParameters( int nMarkerType,  
                                int nWaveLength,  
                                int nModelType )
```

## Parameters

**nMarkerType** is the type of marker for which the camera parameter set has been optimized. The marker types currently supported are:

**Table 9-3: Supported Types of Markers**

| <b>nMarkerType</b> | <b>Description</b>           |
|--------------------|------------------------------|
| 1                  | Metal base markers (default) |
| 2                  | Ceramic base markers         |
| 3                  | Unknown marker type          |

**nWaveLength** is the type of marker wavelength for which the camera parameter set has been optimized. The marker wavelength types currently supported are:

**Table 9-4: Supported Marker Wavelength**

| <b>nWaveLength</b> | <b>Wavelength (nm)</b> |
|--------------------|------------------------|
| 0                  | 950 (default)          |
| 1                  | 880                    |

**nModelType** is the type of lens model for which the camera parameter set has been optimized. The lens model types currently supported are:

**Table 9-5: Supported Types of Lens Model**

| <b>nModelType</b> | <b>Description</b>                             |
|-------------------|------------------------------------------------|
| 0                 | The original lens model type (default).        |
| 1                 | A new lens model.                              |
| 2                 | A new lens model optimized for larger volumes. |

## Description

`OptotrakSetCameraParameters` specifies the camera parameters to be used by the API. The camera parameter sets are identified according to their marker type, marker wavelength, and lens model type. Multiple camera parameter sets can be stored in the extended camera parameter files. When an extended camera parameter file is loaded by the API routine `OptotrakLoadCameraParameters`, the default camera parameter set is automatically selected. Call `OptotrakSetCameraParameters` afterwards to select a non-default camera parameter set. The routine will return an error if the requested camera parameter set is not part of the loaded camera parameter sets. Use `OptotrakGetCameraParameterStatus` to list all of the available camera parameter sets.

## See Also

`OptotrakLoadCameraParameters`, `OptotrakGetCameraParameterStatus`

## 9.4.17 OptotrakSetProcessingFlags

### Function

Sets or clears bit flags that enable various processing options for the API. These flags can also be set in the initialization file, optotrak.ini.

### Prototype

```
int OptotrakSetProcessingFlags( unsigned int uFlags )
```

### Parameters

**uFlags** specifies the flags used to control various processing options. This parameter passes information using bit flags. To pass two or more flags to the routine, separate them using the logical OR operator “|”.

Values:

**OPTO\_LIB\_POLL\_REAL\_DATA** controls blocking in real-time data retrieval routines.

**OPTO\_CONVERT\_ON\_HOST** causes 3D conversions to be done on the host computer instead of on the system.

**OPTO\_RIGID\_ON\_HOST** causes rigid body 6D transformations to be done on the host computer instead of on the system.

**OPTO\_USE\_INTERNAL\_NIF** causes the API to use network information stored internally rather than the network information loaded from a network information file (the default file system.nif or some user specified equivalent file).

### Description

OptotrakSetProcessingFlags controls the way various data processing is done. If this routine is not called, the **OPTO\_LIB\_POLL\_REAL\_DATA** bit is set and **OPTO\_CONVERT\_ON\_HOST**, **OPTO\_RIGID\_ON\_HOST**, and **OPTO\_USE\_INTERNAL\_NIF** are not set in the system by default. As a result, all real-time data retrieval routines will not block. All 3D conversions and rigid body transformations will be done on the system and the API will get its network information from an external file such as system.nif.

When the OptotrakSetProcessingFlags routine is called, the bit flags settings in the parameters will override these default settings.



## Alternative Methods of Setting the First Three Flags

The first three bit flags (OPTO\_LIB\_POLL\_REAL\_DATA, OPTO\_CONVERT\_ON\_HOST, and OPTO\_RIGID\_ON\_HOST) can either be specified in a collection parameter file called by OptotrakSetupCollectionFromFile or during system initialization. To set these flags during system initialization, you must first specify the corresponding parameters in the [OPTOTRAK System] section of the parameter file, optotrak.ini:

```
[OPTOTRAK System]
bConvertOnHost    = FALSE
bRigidOnHost      = FALSE
PollRealDataFlag = TRUE
```

After setting the flags, call TransputerInitializeSystem. If you use TransputerInitializeSystem to read the values from the file optotrak.ini, the flag settings are not obvious from the program source code and may be changed by another program or user. The values are, however, easily viewed in the .ini file. If you use OptotrakSetProcessingFlags, then it is clear from the program source code what flags are used, but you cannot change the flag values without recompiling your program and the flag settings are not easily viewed if you do not have access to the source code. Choose the method that is most appropriate to your project.

## The Effect of OPTO\_LIB\_POLL\_REAL\_DATA on DataIsReady

In the default situation, the routine DataIsReady returns FALSE if data is not yet ready for retrieval, and TRUE when the requested data is waiting to be retrieved. If OPTO\_LIB\_POLL\_REAL\_DATA is not set, DataIsReady blocks until the data is actually ready, and then returns TRUE. The computer processor will be less busy because the program will avoid making multiple calls to DataIsReady. For users with the SCSI option, this saves one SCSI cycle for each group of function calls of the type {RequestLatest*type*, DataIsReady, DataReceiveLatest*type*}, where *type* represents a data type such as 3D, Transforms, and OdauRaw.

## Choosing an Appropriate Value for OPTO\_LIB\_POLL\_REAL\_DATA

Most application programs for real-time data retrieval use one of the following three approaches:

**First Approach:** With OPTO\_LIB\_POLL\_REAL\_DATA not set and using DataGetLatest*type*, call DataIsReady to request for available data. DataIsReady will block until the data is ready. This approach is the easiest method to use, but may result in slower processing than other methods.

```
{
    DataGetLatesttype (blocks until data is ready)
    some data processing code
}
```

**Second Approach** With `OPTO_LIB_POLL_REAL_DATA` not set, and using `RequestLatesttype` and `DataReceiveLatesttype`, the waiting time is reduced because the request for data is being processed by the system while the user program is doing data processing steps. In the first approach, the request for the next frame of data is not sent until the processing of the current frame of data is completed.

```
RequestLatesttype
{
    DataIsReady (blocks until data is ready)
        ReceiveLatesttype
        RequestLatesttype
        some data processing code
}
DataIsReady
ReceiveLatesttype
```

**Third Approach** With `OPTO_LIB_POLL_REAL_DATA` set and using `RequestLatesttype` and `DataReceiveLatesttype`, more processing cycles (and SCSI cycles, if applicable) are used, but there is no blocking while waiting for data to be ready for retrieval and your program can do other processing tasks while waiting for data.

```
{
    RequestLatesttype
        while( !DataIsReady )
        {
            some code (not for processing the current frame of data)
        }
    ReceiveLatesttype
    some data processing code for the current frame of data
}
```

### **Setting `OPTO_CONVERT_ON_HOST` and `OPTO_RIGID_ON_HOST`**

The bit flags `OPTO_CONVERT_ON_HOST` and `OPTO_RIGID_ON_HOST` provide you with the option of doing 3D conversions and rigid body 6D transformations on the host computer, instead of on the system. For computers with Pentium processors or faster, doing the conversions on the host computer can be faster than doing them on the system. Also, because the host computer does not have the same memory constraints as the system, the maximum number of rigid bodies that can be handled when on-host transformations are enabled is much greater (85 versus 10 for on-system conversions).

If neither `OPTO_CONVERT_ON_HOST` nor `OPTO_RIGID_ON_HOST` are set (the default), then all conversions and rigid body transformations are done on the system, as in the past, but if both are set, then both are done on the host computer.

---

**Note** You can set `OPTO_RIGID_ON_HOST` and not set `OPTO_CONVERT_ON_HOST` to have only the rigid body transformations done on the host computer, but the combination of setting `OPTO_CONVERT_ON_HOST`, but not `OPTO_RIGID_ON_HOST`, is not supported.

---

### Setting `OPTO_USE_INTERNAL_NIF`

To use internally generated and stored network information for downloading the transputer programs and startup code, use the bit flag `OPTO_USE_INTERNAL_NIF`. The default method of obtaining the network information is the external file `system.nif`, or a user-specified equivalent file. If this flag is set, then the application program must subsequently call `TransputerDetermineSystemCfg` to determine the network information and store it locally.

See Sample Program 19 on the API CD for an example of code that uses this routine.

### See Also

`DataIsReady`, `DataGetLatesttype`, `DataReceiveLatesttype`, `RequestLatesttype`, `TransputerDetermineSystemCfg`, `TransputerInitializeSystem`, `OptotrakSetupCollectionFromFile`

## 9.4.18 `OptotrakSetStroberPortTable`

### Function

Configures the strober port firing table for the Optotrak System according to the specified parameters.

### Prototype

```
int OptotrakSetStroberPortTable( int nPort1,
                                int nPort2,
                                int nPort3,
                                int nPort4 )
```

### Parameters

`nPort1` specifies the number of markers connected to strober port 1.

`nPort2` specifies the number of markers connected to strober port 2.

`nPort3` specifies the number of markers connected to strober port 3.

**nPort4** specifies the number of markers connected to strober port 4 (Must be set to 0 for Optotrak Certus System).

---

**Note** The Optotrak Certus System only has three ports. Port4 must be set to 0 when using the Optotrak Certus System.

---

## Description

OptotrakSetStroberPortTable is used for programs running on an Optotrak 3020 System or an Optotrak 2020 System. This function should be called before calling OptotrakSetupCollection or OptotrakSetupCollectionFromFile. The strober port table can also be specified in the [Strober Table] section of the parameter file used by OptotrakSetupCollectionFromFile.

This routine can also be used with an Optotrak Certus System, but only if no device handle routines are called.

---

**The total number of markers specified for the four ports must equal the number of markers specified by the NumberOfMarkers parameter in the section [OPTOTRAK System] of the parameter file used by OptotrakSetupCollectionFromFile, or nMarkers specified in OptotrakSetupCollection. If the two numbers do not match, you will receive an error message.**

---

## See Also

OptotrakSetupCollection, OptotrakSetupCollectionFromFile

## 9.4.19 OptotrakSetupCollection

### Function

Configures the Optotrak System collection using the specified parameters.

### Prototype

```
int OptotrakSetupCollection( int    nMarkers,  
                           float  fFrameFrequency,  
                           float  fMarkerFrequency,  
                           int     nThreshold,  
                           int     nMinimumGain,  
                           int     nStreamData,  
                           float  fDutyCycle,
```

```
float fVoltage,
float fCollectionTime,
float fPreTriggerTime,
int nFlags )
```

## Parameters




---

**Do not set marker voltage, duty cycle or marker power such that the markers become hot and cause personal injury and/or property damage.**

---

**nMarkers** specifies the number of markers for which data is collected.

Optotrak 3020 System Bounds: [1, 255]

Optotrak Certus System Bounds: [1, 512]

**fFrameFrequency** specifies the rate at which an entire frame of data (a position for each marker) is generated.

Optotrak 3020 System Bounds:  $\left[1, \frac{3500}{N+1}\right]$

Optotrak Certus System Bounds:  $\left[1, \frac{4600}{N+2}\right]$  (not for every system configuration)

where N = number of markers.

**fMarkerFrequency** specifies the frequency at which individual markers are strobed.

Optotrak 3020 System Bounds: [1000, 3500]

Optotrak Certus System Bounds: [1500, 4600]

**nThreshold** specifies the noise threshold under which all sensor data is ignored when the centroid is being calculated. The value passed via this parameter defines either a static threshold value or a dynamic threshold value, depending on the mode in which the Optotrak System is functioning.

Dynamic Threshold Bounds: [0, 100]

Static Threshold Bounds: [0, 3000]

**nMinimumGain** is the maximum amplification that can be applied to the signal received by the sensors. Bounds: [0, 255]

**nStreamData** indicates whether the Optotrak System is to either send buffered data back by request only, or automatically send buffered data back once data spooling is initiated.

Values:

**nStreamData** = 0: Send buffered data back at request only.

**nStreamData** = 1: Send buffered data back automatically.

**fDutyCycle** is the fraction of time that a marker is actually turned on during the marker period. Bounds: [0.1, 0.85] (lower upper bounds for some system configurations)

**fVoltage** is the voltage level used for strobing the markers. Bounds: [7.0, 12.0]

**fCollectionTime** is the duration time, in seconds, for buffered data collections. Bounds: [0, 99999]

**fPreTriggerTime** is not supported and must be 0.

**nFlags** indicates settings of low-level parameters for the collections. Generally the flags parameter is set to zero. However, different types of collections can be set up by setting individual bits in the flags parameter using the given constants and the logical OR “|” operator.

Values:

**OPTOTRAK\_NO\_INTERPOLATION\_FLAG** suppresses time interpolation of raw data within a frame.

**OPTOTRAK\_FULL\_DATA\_FLAG** collects and buffers full raw data rather than just raw data if it is used in combination with **OPTOTRAK\_BUFFER\_RAW\_FLAG**.

**OPTOTRAK\_BUFFER\_RAW\_FLAG** collects and buffers raw data instead of 3D data.

**OPTOTRAK\_NO\_FIRE\_MARKERS\_FLAG** suppresses the automatic activation of the markers upon completion of the `OptotrakSetupCollection` routine. If this flag is set, `OptotrakActivateMarkers` must be used to activate the markers before starting a collection. See [“OptotrakActivateMarkers” on page 105](#) for further details.

**OPTOTRAK\_STATIC\_THRESHOLD\_FLAG** signals the Optotrak System to use the `nThreshold` parameter as a static noise threshold value. If this flag is not set, a percentage of the waveform peak is used as the threshold. In this case, the parameter `nThreshold` specifies the percentage of the peak to use. The recommended value is 30. Bounds is [0, 255].

**OPTOTRAK\_AUTO\_DUTY\_CYCLE\_FLAG** dynamically adjusts the exposure time for each marker so that all of the markers have roughly the same perceived power. This flag is obsolete since Optotrak 3020 Systems and Optotrak Certus Systems use electronic shutter control.

**OPTOTRAK\_EXTERNAL\_CLOCK\_FLAG** indicates to the system that a TTL clock signal will be supplied to either pin 2 of the 9-pin D-shell connector on the back of the 3020 System Control Unit, or pin 3 of the 9-pin D-shell connector on

the back of the Certus System Control Unit. Set the collection frequency in your program to be as close as possible to the frequency of the clock signal that is in use.

**OPTOTRAK\_EXTERNAL\_TRIGGER\_FLAG** requires that a combination of software and hardware events must be used to start a collection. First the program must issue the call to start a collection in the usual way (using `DataBufferStart` or `DataBufferSpoolData`). This call will not start the collection, but will set up the system to prepare for the collection. The data collection actually starts when a high to low TTL signal is provided to either pin 3 of the 9-pin D-shell connector on the back of the 3020 System Control Unit, or pin 7 of the 9-pin D-shell on the back of the Certus System Control Unit.

**OPTOTRAK\_GET\_NEXT\_FRAME\_FLAG** indicates to the system that the next available frame should be returned instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. Instead of using this flag, you can use the routine `DataGetNexttype` for the type of data you require.

**OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG** indicates that switch data should be returned by the System Control Unit to the application program. If this flag is not set, all switch data is ignored. This flag is only used with Optotrak Certus Systems.

## Description

`OptotrakSetupCollection` uses the specified parameters to configure a collection on the Optotrak System. This routine only configures the Optotrak System collection parameters. If an ODAU is connected to the Optotrak System, then `OdaSetupCollection` must also be called.



**The last collection setup routine called must be for the Optotrak System. For an example of this, see “Sample Program 11” on page 260.**

**Note** After a call to the routine `OptotrakSetupCollection`, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may be generated.

See Sample Program 2 on the API CD for an example of code that uses this routine.

## See Also

OptotrakSetupCollectionFromFile, OdaSetupCollection, OptotrakStopCollection

### 9.4.20 OptotrakSetupCollectionFromFile

#### Function

Configures the Optotrak System collection parameters based on the contents of an collection parameter file.

#### Prototype

```
int OptotrakSetupCollectionFromFile( char *pszCollectFile )
```

#### Parameters

**pszCollectFile** points to a null terminated string specifying the collection parameter file to be used.

#### Description

OptotrakSetupCollectionFromFile reads a specified collection parameter file and sets up a collection based on its parameters. Collection parameter files can be created either by the NDI ToolBench software or by any ASCII editor. A sample parameter file is found in [Figure 9-1 on page 129](#).

```
[OPTOTRAK System]
NumberOfMarkers      = 10
FrameFrequency       = 40
MarkerFrequency      = 2500
DutyCycle            = 0.50
StroberVoltage       = 7.0
DynamicThresholdPercent = 30
InterpolationFlag    = TRUE
StaticThresholdFlag  = FALSE
StaticThresholdValue = 200
MinimumGain          = 160
CollectionTime       = 10

[RigidBody 0]
StartMarker          = 1
RigidFile            = Sample
PerformQuaternionFlag = 1
PerformIterativeFlag = 0
ReturnFormat         = EULER
```



```
[Strober Table]
Port1           = 6
Port2           = 4
Port3           = 0
Port4           = 0
```

**Figure 9-1: A Sample of a Parameter File**

## Parameter File Details

There are usually three types of sections in a collection parameter file. If the parameter file is generated by NDI ToolBench or another NDI program, it may contain some additional sections. There may also be sections for ODAU devices.

The first section type is called `[OPTOTRAK System]`. Within this section, values may be specified for any collection setup or system setup parameter. The beginning of the setup parameters section is indicated by the text “OPTOTRAK System”.

The second section type is called `[RigidBody n]` where  $n$  is the rigid body identifier. If real-time rigid body calculations are available (see “[Real-time Rigid Body Programmer’s Guide](#)” on page 51), then the RigidBody sections of the collection parameter file can be used.

For rigid body conversions done on the system, up to 10 rigid bodies can be specified, with rigid body identifiers ranging from 0 to 9, inclusive. With rigid body conversions done on the host computer, up to 85 rigid bodies can be specified.

The third section type is called `[Strober Table]`. This supports multiple strober ports.

The parameters that are contained in each section type are described below:

**The parameters for the `[OPTOTRAK System]` section are:**

---

**Note** You must name the parameters exactly as shown in the list below and you must always include the `FrameFrequency` and `NumberOfMarkers` parameters. You do not need to specify the other parameters: defaults are used if they are not specified.

---

`NumberOfMarkers` is the number of markers to setup for the collection.

Optotrak 3020 System Bounds: [1, 255]

Optotrak Certus System Bounds: [1, 512]

`FrameFrequency` is the frame rate for the collection.

Optotrak 3020 System Bounds: [1,  $\frac{3500}{N+1}$ ]

Optotrak Certus System Bounds:  $[1, \frac{4600}{N+2}]$  (not for every system configuration)

where N = number of markers.

**MarkerFrequency** is the marker rate for the collection.

Optotrak 3020 System Bounds: [1500, 3500]

Optotrak Certus System Bounds: [1500, 4600]

**DutyCycle** is the fraction of the marker period that the marker is turned on.

Bounds: [0.1, 0.85] (lower upper bounds for some system configurations)

Default value is 0.5.

**StroberVoltage** is the voltage applied to the markers. Bounds: [7.0, 12.0]. Default value is 7.0 Volts.



---

**Do not set marker voltage, duty cycle or marker power such that the markers become hot and cause personal injury and/or property damage.**

---

**CollectionTime** is the duration of the file collections, in seconds. Bounds: [0, 99999]. Default value is 2 seconds.

**StreamData** is set to 1 if the data is to be streamed back to the application. If set to 0, then the application must request the data explicitly. Bounds: [0, 1]. Default value is 1.

**PreTriggerTime** is not currently supported. Set to 0 if specified. Default value is 0.

**MinimumGain** is rarely required in the parameter file. The value indicates the minimum gain to be used by the system. The gain values affect the amount of signal amplification on a sensor. Bounds: [0, 255]. Default value is 160.

**DynamicThresholdPercent** is rarely required in the parameter file. It indicates the percentage of the signal peak to be considered when determining the centroid. Bounds: [0, 100]. Default value is 30.

**StaticThresholdFlag** is rarely required in the parameter file. Set to 1 only if static thresholding is required rather than dynamic thresholding in the determination of the centroids. Bounds: [0, 1]. Default value is 0.

**StaticThresholdValue** is rarely required in the parameter file. If static thresholding is set, then this value is used as the threshold for determining centroids. Bounds: [0, 3000]. Default value is 200.

**InterpolationFlag** is rarely required in the parameter file. Set to 0 if no marker interpolation is desired; set to 1 if marker interpolation is required. Bounds: [0, 1]. Default value is 1 (interpolation on).

**Buffer3dFlag** is set to 1 if 3D data is to be buffered. Bounds: [0, 1]. Default value is 0 (off).

**BufferFullRawFlag** is set to 1 if full raw data is to be buffered during data collections and 0 if raw data is to be buffered. Bounds: [0, 1]. Default value is 0 (off).

**GetNextFrameFlag** is set to 1 to indicate to the system that the next available frame should be returned instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. Bounds: [0, 1]. Default is 0 (off).

**ExternalSyncFlag** is set to 1 to indicate to the system that a TTL clock signal will be supplied to either pin 2 of the 9-pin D-shell connector on the back of the 3020 System Control Unit, or pin 3 of the 9-pin D-shell connector on the back of the Certus System Control Unit. The collection frequency should be set in the application program as closely as possible to the frequency of the clock signal that is used. Bounds: [0, 1]. Default is 0 (off).

**ExternalTriggerFlag** is set to 1 to require that a combination of software and hardware events be used to start a data file collection. The application program must first start a collection in the usual manner (using `DataBufferStart` or `DataBufferSpoolData`), which will not actually start the collection, but will set up the system to prepare for the collection. The data collection will start when a high-to-low TTL signal is provided to either pin 3 of the 9-pin D-shell connector on the back of the 3020 System Control Unit, or pin 7 of the 9-pin D-shell connector on the back of the Certus System Control Unit. Bounds: [0, 1]. Default is 0 (off).

**PollRealDataFlag** is set to 1 to have the application program poll the system to determine if data is ready to be sent from the system. This flag can also be set by calling `OptotrakSetProcessingFlags` with the `OPTO_LIB_POLL_REAL_DATA` bit set. Bounds: [0, 1]. Default is 1 (on).

**bConvertOnHost** is set to 1 to have the conversion of centroid data to 3D data done on the host computer instead of on the system. This flag can also be set by calling `OptotrakSetProcessingFlags` with the `OPTO_CONVERT_ON_HOST` bit set. Bounds: [0, 1]. Default is 0 (off).

**bRigidOnHost** is set to 1 to have the transformation of 3D position data to 6D rigid body data done on the host computer instead of on the system. This flag can also be set by calling `OptotrakSetProcessingFlags` with the `OPTO_RIGID_ON_HOST` bit set. Bounds: [0, 1]. Default is 0 (off).

**bUseColParmOnlyFlag** specifies that only the following parameters will be used: `FrameFrequency`, `MarkerFrequency`, `DutyCycle`, `StroberVoltage`, and `CollectionTime`. This parameter is only used with Optotrak Certus Systems.

**bUseSwitchesFlag** indicates that switch data should be returned by the System Control Unit to the application. If this flag is not set, all switch data is ignored. This parameter is only used with Optotrak Certus Systems.

**ModelType** indicates the camera parameter set model type currently being used by the API to convert full raw data to 3D positions. Bounds: formally [0, 64564]; the actual upper bound is dependent on the number camera parameter model types currently supported. Default is ModelType 0.

**MarkerType** indicates the marker type that the API has been optimized for in the current camera parameter set. Bounds: formally [1, 64564]; the actual upper bound is dependent on the number of marker types currently supported. Default is MarkerType 1 (metal base).

**WaveLength** indicates the marker wavelength type that the API has been optimized for the current camera parameter. Bounds: formally [0, 64564]; the actual upper bound is dependent on the number of marker wavelengths currently supported. Default is Wavelength type 0 (950 nm).

**The parameters for the [RigidBody] section are:**

---

**Note** **The RigidFile and StartMarker for the rigid body must be specified to correctly load and use the rigid body.**

---

**RigidFile** is the name of the file defining the rigid body (without the file extension). The system searches for the file <RigidFile>.rig under the current directory and under the “rigid” subdirectory under the path specified in the ND\_DIR and ND\_USER\_DIR environment variables. No default is provided.

**StartMarker** is the number of the first marker for the rigid body. Bounds: [1, NumberOfMarkers]. No default is provided.

**PerformQuaternionFlag** is set to 1 if the rigid body determination is to be done using the quaternion algorithm. This is a high-speed rigid body algorithm that requires 3D position data for at least three markers in the rigid body. This algorithm can be used in conjunction with the iterative algorithm. If neither the quaternion nor the iterative algorithm are specified explicitly, the system will use the quaternion algorithm only by default. Bounds: [0, 1]. Default value is 0.

**PerformIterativeFlag** is set to 1 if the rigid body determination is to be done using the iterative algorithm. This is a rigid body algorithm that requires sensor values (centroids) for the markers in the rigid body rather than 3D values. This algorithm can be used in conjunction with the quaternion algorithm. If neither the quaternion nor iterative algorithm are specified explicitly, the system will use the quaternion algorithm only by default. Bounds: [0, 1]. Default value is 0.

**ReturnFormat** is the desired return format for real-time rigid bodies. There are three formats available; Euler, rotation matrix, and quaternion. Euler indicates that Euler parameters are returned when real-time rigid bodies are requested. The Euler parameters are roll ( $R_z$ ), pitch ( $R_y$ ), yaw ( $R_x$ ), and translation values  $T_x$ ,  $T_y$ ,  $T_z$ . Rotation matrix indicates that a  $3 \times 3$  rotation matrix is returned as well as the translation values ( $T_x$ ,  $T_y$ ,  $T_z$ ). Quaternion signals that quaternion format is returned by the system. A quaternion is a four dimensional vector ( $q_0$ ,  $q_x$ ,  $q_y$ , and  $q_z$ ) providing orientation in addition to the translation values ( $T_x$ ,  $T_y$ ,  $T_z$ ).

---

**Note** The return format is independent of the algorithm used to determine the rigid body. All three formats can be returned when using either the Quaternion algorithm or the iterative algorithm. Valid values: [EULER, MATRIX, QUATERNION]. Default value is EULER.

---

The parameters in the [Strober Table] section contain the number of markers connected to each strober port.

**Port1** is the number of markers connected to strober port 1.

**Port2** is the number of markers connected to strober port 2.

**Port3** is the number of markers connected to strober port 3.

**Port4** is the number of markers connected to strober port 4 (option not available with Optotrak Certus Systems).

---

**Note** The Optotrak Certus System only has three ports. Port4 must be set to 0 when using an Optotrak Certus System.

---

See the OptotrakSetupCollection routine for further details on the specified values.

---

**Note** The total number of markers specified for the four ports must equal the number of markers specified by the NumberOfMarkers parameter in the section [OPTOTRAK System]. If the two numbers do not match, you will receive an error message.

---

## See Also

OptotrakSetupCollection, OdaSetupCollection, OptotrakSaveCollectionToFile

## 9.4.21 OptotrakStopCollection

### Function

Stops the collection executing in the Optotrak System.

## Prototype

```
int OptotrakStopCollection( )
```

## Parameters

None

## Description

OptotrakStopCollection stops the collection executing in the Optotrak Certus System. This function should be called before changing a device handle's status.

When a device is either added to, or removed from an Optotrak Certus System, that device handle's status is changed within the Optotrak Certus System Control Unit. The uFlags parameter—found in numerous Optotrak API routines (e.g. OptotrakGetLatestRaw)—informs the system that a change has occurred. However, if a device is disconnected from the system during a collection, that collection will return MISSING for the disconnected device.

The application program should call OptotrakStopCollection and re-initialize the device handles if a device is either added to, or removed from a system. After all devices have been re-initialized you may begin a new collection.

## See also

OptotrakSetupCollection, OptotrakSetupCollectionFromFile

## 9.5 Optotrak Device Handle Routines

The following routines are commonly used for identifying and manipulating device handles. Device handles are used in Optotrak Certus Systems only, please see [“Optotrak Certus Specific Routines” on page 94](#) for a general discussion of Optotrak Certus device handles.

### 9.5.1 OptotrakDeviceHandleEnable

#### Function

Enables the device associated with a specified device handle.

#### Prototype

```
int OptotrakDeviceHandleEnable( int nDeviceHandleID )
```

#### Parameters

`nDeviceHandleID` identifies the device.

#### Description

A device must be enabled before the application can retrieve data from that device. Once a device is enabled it may activate markers, monitor switches, or set visible LED information.

### 9.5.2 OptotrakDeviceHandleFree

#### Function

Removes a device handle no longer used by the system.

#### Prototype

```
int OptotrakDeviceHandleFree ( int nDeviceHandleID )
```

#### Parameters

`nDeviceHandleID` identifies the device.

## Description

When a device is disconnected from the system the associated device handle is set to `DH_STATUS_UNOCCUPIED` by the Optotrak Certus System Control Unit. The application program should perform all required cleanup (e.g. free memory associated with the device) and then call `OptotrakDeviceHandleFree` to indicate that the device is no longer required by the system. When the Optotrak Certus System Control Unit receives this command the device handle state changes from `UNOCCUPIED` to `FREE`.

If the same device is later reconnected to the system it will receive a new device handle.

### 9.5.3 OptotrakDeviceHandleGetNumberProperties

#### Function

Determines the number of properties associated with a device handle. The application should call this routine before calling `OptotrakDeviceHandleGetProperties` to retrieve the device properties.

#### Prototype

```
int OptotrakDeviceHandleGetNumberProperties( int nDeviceHandleID,  
   int *pnProperties )
```

#### Parameters

`nDeviceHandleID` identifies the device.

`pnProperties` is set to the number of properties associated with a specified device.

#### Description

After a device is initialized the application may request the device properties, including the physical location, the number of switches, and the device capabilities.

Each of these properties are stored in a `DeviceHandleProperty` structure. The application should allocate enough memory for all the device handle properties before calling the `OptotrakDeviceHandleGetProperties` routine. Each property requires the amount of memory needed for a `DeviceHandleProperty` structure.

See Optotrak Certus Sample Program 3 on the API CD for sample code that uses this routine.



## 9.5.4 OptotrakDeviceHandleGetProperties

### Function

Retrieves the device handle properties for a device, including the physical location, the number of switches, and the device capabilities.

### Prototype

```
int OptotrakDeviceHandleGetProperties( int                nDeviceHandleID,
                                     DeviceHandleProperty *grdtProperties,
                                     int                nProperties )
```

### Parameters

**nDeviceHandleID** identifies the device.

**grdtProperties** points to the buffer that will store the device handle properties.

**nProperties** is the number of properties expected in the buffer.

### Description

OptotrakDeviceHandleGetProperties retrieves property information from the specified device. Each property is stored in a DeviceHandleProperty structure. The application must allocate enough memory for all the device handle properties. Each property requires the amount of memory needed for a DeviceHandleProperty structure. To determine the number of properties associated with a device, the application should first call OptotrakDeviceHandleGetNumProperties. For a list of properties, please see [“Device Handle Properties” on page 95](#).

Values for dtPropertyType:

```
DH_PROPERTY_TYPE_STRING,
DH_PROPERTY_TYPE_INT,
DH_PROPERTY_TYPE_FLOAT,
DH_PROPERTY_TYPE_DOUBLE,
DH_PROPERTY_TYPE_CHAR,
```

```
typedef struct
{
    unsigned int
        uPropertyID;
    enum DeviceHandlePropertyType
        dtPropertyType;
    union
    {
        int
```

```
        nData;  
    float  
        fData;  
    double  
        dData;  
    char  
        cData;  
    char  
        szData[MAX_PROPERTY_STRING];  
    } dtData;  
} DeviceHandleProperty;
```

## 9.5.5 OptotrakDeviceHandleGetProperty

### Function

Retrieves a specified device handle property from a device.

### Prototype

```
int OptotrakDeviceHandleGetProperty( int          nDeviceHandleID,  
                                     DeviceHandleProperty *pdtProperty,  
                                     unsigned int    uPropertyID )
```

### Parameters

**nDeviceHandleID** identifies the device.

**pdtProperty** points to the buffer where the device handle property information will be written.

**uPropertyID** specifies the property to be retrieved from the device.

### Description

OptotrakDeviceHandleGetProperty retrieves a single property from a device. The actual property information—specified by the uPropertyID parameter—is stored at the location pointed to by pdtProperty. For a list of properties, please see [“Device Handle Properties” on page 95](#).

See Optotrak Certus Sample Program 3 on the API CD for sample code that uses this routine.

## 9.5.6 OptotrakDeviceHandleSetBeeper

### Function

Starts the beeper on a device for a specified period of time.

### Prototype

```
int OptotrakDeviceHandleSetBeeper( int          nDeviceHandleID,  
                                   unsigned int  uTime_ms )
```

### Parameters

**nDeviceHandleID** identifies the device.

**uTime\_ms** is the length of time that the beeper will be on, in milliseconds.

### Description

OptotrakDeviceHandleSetBeeper starts the beeper on the strober for the specified period of time.

See Optotrak Certus Sample Program 12 on the API CD for sample code that uses this routine.

## 9.5.7 OptotrakDeviceHandleSetProperties

### Function

Sets one or more properties for a device.

### Prototype

```
int OptotrakDeviceHandleSetProperties( int          nDeviceHandleID,  
                                       DeviceHandleProperty *grdtProperties,  
                                       int          nProperties )
```

### Parameters

**nDeviceHandleID** identifies the device.

**grdtProperties** is an array of DeviceHandleProperty structures specifying the properties to be set for the device.

**nProperties** is the size of the properties array.

## Description

OptotrakDeviceHandleSetProperties sends a set of properties to a specified device. Properties are stored in the grdtProperties array, which must contain at least nProperties entries. For information on the structure of the properties array, please see “[OptotrakDeviceHandleGetProperties](#)” on page 137.

See Optotrak Certus Sample Program 15 on the API CD for sample code that uses this routine.

## 9.5.8 OptotrakDeviceHandleSetVisibleLED

### Function

Sets the state of a visible light emitting diode (LED) for a device.

### Prototype

```
int OptotrakDeviceHandleSetVisibleLED ( int          nDeviceHandleID,  
                                       int          nLED,  
                                       unsigned int  uState )
```

### Parameters

**nDeviceHandleID** identifies the device.

**nLED** identifies the visible LED.

**uState** identifies the state of the visible LED. The three possible states are VLEDST\_ON, VLEDST\_OFF, and VLEDST\_BLINK.

### Description

OptotrakDeviceHandleSetVisibleLED sets the state of a visible LED for a device.

---

**Note** The first LED on a strober is reserved. The status of this LED cannot be set. Attempting to change the status of this LED will return an error.

---

See Optotrak Certus Sample Program 11 on the API CD for sample code that uses this routine.

## 9.5.9 OptotrakGetDeviceHandles

### Function

Retrieves a list of device handles assigned by the Optotrak System, including the status of each device.

### Prototype

```
int OptotrakGetDeviceHandles( DeviceHandle *grdtDeviceHandles,
                             int          nDeviceHandles,
                             unsigned int *puFlags )
```

### Parameters

**grdtDeviceHandles** is an array of DeviceHandle structures that will be filled by the System Control Unit with a list of existing device handles.

**nDeviceHandles** is the size of the array of DeviceHandle structures.

**puFlags** is the flags returned by the System Control Unit indicating the status of the device configuration.

### Description

OptotrakGetDeviceHandles requests a list of known device handles from the Optotrak System. The application program should first determine the number of device handles by calling the OptotrakGetNumberDeviceHandles routine. The application program should allocate enough memory to accept the expected number of device handles.

If the number of device handles requested does not match the number of device handles known by the System Control Unit, the puFlags variable will contain the flag OPTO\_TOOL\_CONFIG\_CHANGED\_FLAG.

The DeviceHandle structure contains the device handle ID and the status of the device handle.

Device status:

```
DH_STATUS_UNOCCUPIED,
DH_STATUS_OCCUPIED,
DH_STATUS_INITIALIZED,
DH_STATUS_ENABLED,
DH_STATUS_FREE
```

```
typedef struct
{  int
```

```
        nID;                /* device ID */
    DeviceHandleStatus
        dtStatus;          /* device status */
} DeviceHandle;
```

You should check the status of each device handle once the list of device handles has been determined. If the device status is `DH_STATUS_UNOCCUPIED`, the application needs to free the device handle using the `OptotrakDeviceHandleFree` routine. If the device status is `DH_STATUS_INITIALIZED`, the application needs to enable the device using the `OptotrakDeviceHandleEnable` routine.

### 9.5.10 OptotrakGetNumberDeviceHandles

#### Function

Retrieves the number of device handles currently assigned by the System Control Unit.

#### Prototype

```
int OptotrakGetNumberDeviceHandles( int *pnDeviceHandles )
```

#### Parameters

`pnDeviceHandles` is set to the number of devices assigned by the System Control Unit.

#### Description

`OptotrakGetNumberDeviceHandles` queries the Optotrak System for the number of assigned devices. The application should allocate enough memory for the assigned device handles before calling the `OptotrakGetDeviceHandles` routine.

## 9.6 ODAU Specific Routines

See the “[ODAU Programmer’s Guide](#)” on page 41 for a general discussion of ODAU routines.

### 9.6.1 `OdauGetStatus`

#### Function

Retrieves the current status from the specified ODAU device. This applies to both ODAU I and II. All parameters are set in the routines “[OdauSetupCollection](#)” on page 149 or “[OdauSetupCollectionFromFile](#)” on page 152.

#### Prototype

```
int OdauGetStatus( int      nOdauId,
                  int      *pnChannels,
                  int      *pnGain,
                  int      *pnDigitalMode,
                  float     *pfFrameFrequency,
                  float     *pfScanFrequency,
                  int      *pnStreamData,
                  float     *pfCollectionTime,
                  float     *pfPreTriggerTime,
                  unsigned  *puCollFlags,
                  int      *pnFlags )
```

#### Parameters

`nOdauId` specifies which ODAU to retrieve status information from. The possible values are constants `ODAU1`, `ODAU2`, `ODAU3` and `ODAU4`.

---

**Note** You may set any of the following parameters to **NULL** if the data provided is not required.

---

`pnChannels` is set to the number of analog channels that are being sampled on the ODAU device.

`pnGain` specifies the gain to be applied to the analog channels.

`pnDigitalMode` indicates the operation mode of the digital port.

`pfFrameFrequency` is set to the frequency at which the frames of data are being sampled.

`pfScanFrequency` (**ODAU II only**) is set to the frequency at which the channels are being sampled within a frame.

**pnStreamData** indicates whether the ODAU is to send buffered data back by request only, or to automatically send buffered data back once data spooling is initiated.

Values:

**nStreamData** = 0: Send buffered data back at request only.

**nStreamData** = 1: Send buffered data back automatically.

**pfCollectionTime** specifies the duration time for buffered data collections.

**pfPreTriggerTime** is not supported and will be 0.

**puCollFlags** is 0 if the analog data is being collected in single-ended mode, and is ODAU\_DIFFERENTIAL\_FLAG if the analog data is being collected in differential mode.

**pnFlags** is set to the flag values used when the last ODAU collection was configured.

## **Description**

OdauGetStatus requests and receives the current status information from the specified ODAU device. Status information returned by this routine consists of the current collection configuration for the specified ODAU device and the flags parameters.

## **See Also**

OtotrakGetStatus, OdauSetupCollection

## **9.6.2 OdauSaveCollectionToFile**

### **Function**

Saves the ODAU collection parameters to an ASCII file to be read by OdauSetupCollectionFromFile, or by other programs. This applies to both ODAU I and II.

### **Prototype**

```
int OdauSaveCollectionToFile( char *pszCollectFile )
```



## Parameters

**pszCollectFile** points to a null terminated string specifying the ODAU collection parameter file to be used. Do not specify the file extension. For example, use “expt1” not “expt1.ini”.

## Description

Use this routine to save the ODAU collection parameter values to an ASCII file. If ODAU units are currently connected and set up, then the current collection parameters are saved.

## See Also

OdaSetupCollectionFromFile

### 9.6.3 OdaSetAnalogOutputs

#### Function

Allows an application program to set or retrieve the voltage values on the analog output channels on an ODAU II. This applies only to ODAU II units.

#### Prototype

```
int OdaSetAnalogOutputs( int      nOdaId,  
                        float     *pfVoltage1,  
                        float     *pfVoltage2,  
                        unsigned   uChangeMask )
```

#### Parameters

**nOdaId** specifies on which ODAU unit to set the outputs. The possible values are constants ODAU1, ODAU2, ODAU3 and ODAU4.

**pfVoltage1** is the desired voltage for channel 1, and must be in range -5V to +5V. The API sets it to the actual value of the output channel upon return from this function.

**pfVoltage2** is the desired voltage for channel 2, and must be in range -5V to +5V. The API sets it to the actual value of the output channel upon return from this function.

**uChangeMask** is a mask used to indicate which channel(s) are to be updated:

- 0 = no channels
- 1 = channel 1 only
- 2 = channel 2 only
- 3 = channels 1 and 2

Use the value of 0 to obtain the current voltages.

## Description

Use `OdauSetAnalogOutputs` to retrieve the current ODAU II voltages (`uChangeMask` set to 0), or to set new voltage values for the analog output channels (`uChangeMask` set to 1, 2 or 3).

## See Also

None.

### 9.6.4 `OdauSetDigitalOutputs`

#### Function

Allows an application program to set the states of the digital output channels on an ODAU I or II.

#### Prototype

```
int OdauSetDigitalOutputs (int      nOdauId,  
                           unsigned *puDigitalOut,  
                           unsigned  uUpdateMask )
```

#### Parameters

**nOdauId** specifies on which ODAU unit to set the outputs. The possible values are constants `ODAU1`, `ODAU2`, `ODAU3` and `ODAU4`.

**puDigitalOut** points to a digital output word (unsigned value) containing the desired digital states. For an ODAU II, the lower 8 bits are used, while for an ODAU I, all 16 bits are used. Upon return from this call, `puDigitalOut` will contain the current states of the output bits.

**uUpdateMask** is a mask of bits (that is, it is an unsigned value containing the number 1 in each of the bits corresponding to the channel that should be updated – lower 8 bits for an ODAU II, 16 bits for an ODAU I) used to alter the digital output channels. Using a mask comprised of all bits set to zero allows an application program to obtain the current digital output values without changing them. ODAU I examples are in [Table 9-6 on page 147](#) while ODAU II examples are in [Table 9-7 on page 147](#).

**Table 9-6: Examples of ODAU I Digital Output and Update Masks**

|                                  | Hex    | Binary           | Decimal |
|----------------------------------|--------|------------------|---------|
| All channels                     | 0xFFFF | 1111111111111111 | 65535   |
| All of Port B only               | 0xFF00 | 1111111100000000 | 65280   |
| All of Port A only               | 0x00FF | 0000000011111111 | 255     |
| Port B Channel 1<br>(i.e. bit 8) | 0x0100 | 0000000100000000 | 256     |
| Port A Channel 1<br>(i.e. bit 0) | 0x0001 | 0000000000000001 | 1       |

**Table 9-7: Examples of ODAU II Digital Outputs and Update Masks**

|                    | Hex  | Binary   | Decimal |
|--------------------|------|----------|---------|
| All channels       | 0xFF | 11111111 | 255     |
| All of Port B only | 0xF0 | 11110000 | 240     |
| All of Port A only | 0x0F | 00001111 | 15      |
| Port B channel 1   | 0x10 | 00010000 | 16      |
| Port A channel 1   | 0x01 | 00000001 | 1       |

## Description

Use `OdauSetDigitalOutputs` to retrieve the current ODAU digital output values (`uUpdateMask` set to 0), or to set the states of the digital output channels (`uUpdateMask` set to non-zero values).

---

**Note** When ODAU II data is returned either as real-time or spooled data, the 8 digital channels (if requested) will be in the upper 8 bits of the 16 bit data word.

---

## See Also

None

### 9.6.5 OdauSetTimer

#### Function

Allows an application program to set the user timer on an ODAU II. This routine applies only to ODAU II units.

#### Prototype

```
int OdauSetTimer( int          nOdauId,
                  unsigned int  uTimer,
                  unsigned int  uMode,
                  unsigned long ulVal )
```

#### Parameters

**nOdauId** is the ODAU unit whose user timer is being set. The possible values are constants `ODAU1`, `ODAU2`, `ODAU3` and `ODAU4`.

**uTimer** is the timer to be set. The user timer is `ODAU_TIMER2`.

**uMode** is the user timer operating mode. Allowed values are in [Table 9-8 on page 148](#).

**Table 9-8: User Timer Operating Mode Values**

| <b>uMode Value</b> | <b>User Timer Operating Mode</b> |
|--------------------|----------------------------------|
| 1                  | <code>ODAU_TIMER_ONE_SHOT</code> |
| 2                  | <code>ODAU_TIMER_RATE_GEN</code> |

**Table 9-8: User Timer Operating Mode Values**

| <b>uMode Value</b> | <b>User Timer Operating Mode</b> |
|--------------------|----------------------------------|
| 3                  | ODAU_TIMER_SQUARE_WAVE           |
| 4                  | ODAU_TIMER_RETRIGGER             |

**uVal** is the value to be loaded into the user timer. Bounds are [0, 65535].

### Description

Use `OdauSetTimer` to setup the ODAU II user timer in one of the supported modes.

### See Also

None.

## 9.6.6 OdauSetupCollection

### Function

Configures the collection parameters for the specified ODAU based on the specified parameters. This routine applies to both ODAU I and ODAU II devices.

### Prototype

```
int OdauSetupCollection( int    nOdauId,
                        int    nChannels,
                        int    nGain,
                        int    nDigitalMode,
                        float   fFrameFrequency,
                        float   fScanFrequency,
                        int    nStreamData,
                        float   fCollectionTime,
                        float   fPreTriggerTime,
                        unsigned uFlags )
```

### Parameters

**nOdauID** specifies the ODAU device to which the collection parameters are to be applied. The possible values are constants ODAU1, ODAU2, ODAU3 and ODAU4.

**nChannels** specifies the number of analog channels to collect data from. Bounds are [1,256].

---

**Note** Without the use of multiplexer boards, the maximum values used should be 16 for single-ended collections and 8 for differential. Multiplexer boards can only be used with an ODAU II, not an ODAU I. Differential mode is supported only on an ODAU II.

---

**nGain** specifies the gain (amplification) to be applied to the analog channels. As the gain (amplification) increases, the voltage range that can be sampled decreases.

The allowed values for an ODAU I are:

- nGain** = 1: Voltage range -10 to 10 volts.
- nGain** = 2: Voltage range -5 to 5 volts.
- nGain** = 4: Voltage range -2.5 to 2.5 volts.
- nGain** = 8: Voltage range -1.25 to 1.25 volts.

while the allowed values for an ODAU II are:

- nGain** = 1: Voltage range -10 to 10 volts.
- nGain** = 5: Voltage range -2 to 2 volts.
- nGain** = 10: Voltage range -1 to 1 volts.
- nGain** = 100: Voltage range -0.1 to 0.1 volts.

**nDigitalMode** indicates the operation mode of the digital port. The digital port provides 16 pins of digital input/output for an ODAU I, and 8 pins for an ODAU II. The pins are grouped into two blocks: RegisterA and RegisterB. The specified digital mode determines the configuration of each register (input, output or disabled). The allowed values are set out [Table 9-9 on page 150](#). The Multiplexer values are supported only for the ODAU II.

**Table 9-9: Configuration of Digital Port**

| <b>nDigitalMode</b>      | <b>Port A</b> | <b>Port B</b> |
|--------------------------|---------------|---------------|
| ODAU_DIGITAL_PORT_OFF    | Off           | Off           |
| ODAU_DIGITAL_INPB_INPA   | Input         | Input         |
| ODAU_DIGITAL_OUTPB_INPA  | Input         | Output        |
| ODAU_DIGITAL_OUTPB_OUTPA | Output        | Output        |
| ODAU_DIGITAL_OFFB_MUXA   | Multiplexer   | Off           |

**Table 9-9: Configuration of Digital Port**

| <b>nDigitalMode</b>     | <b>Port A</b> | <b>Port B</b> |
|-------------------------|---------------|---------------|
| ODAU_DIGITAL_INPB_MUXA  | Multiplexer   | Input         |
| ODAU_DIGITAL_OUTPB_MUXA | Multiplexer   | Output        |

---

**Note** To collect digital data at the same time as using **Ports A and B** in output mode, include the digital input flag by setting **nDigitalMode** to: **ODAU\_DIGITAL\_OUTPB\_OUTPA | ODAU\_DIGITAL\_INPB\_INPA**. Also, if one or both ports are set for output, then they should not have external input signals applied to them. Even if they are set to also read input, they should only read back the values set by the output function.

---

**fFrameFrequency** specifies the sampling frequency at which the ODAU is to collect frames of data. Bounds are [1,27000] for an ODAU I, and [1,100000] for an ODAU II.

**fScanFrequency** [ODAU II only] sets the rate at which the channels are sampled within the frame. Bounds are [1, 100000]. Maximum recommended scan rates vary with the gain selected. See [Table 9-12 on page 155](#).

**nStreamData** indicates whether the ODAU is to send buffered data back by request only, or to automatically send buffered data back once data spooling is initiated.

Values:

**nStreamData** = 0: Send buffered data back at request only.

**nStreamData** = 1: Send buffered data back automatically.

**fCollectionTime** specifies the duration time for buffered data collections. Bounds: [1,99999]

**fPreTriggerTime** is not currently supported and must be 0.

**uFlags** should be set to 0 to indicate the collection of analog data in single-ended mode, or it should be set to **ODAU\_DIFFERENTIAL\_FLAG** if differential mode is to be used. For an ODAU I, differential mode is not supported, and **uFlags** should be set to 0.

## Description

**OdauSetupCollection** uses the specified parameters to set up a collection on the specified ODAU device. This routine only configures the collection parameters for the specified ODAU.



---

**The correct setup collection routine must be invoked for each device connected to the Optotrak System. For example, if there were two ODAUs in the Optotrak System configuration, the application program must invoke `OdauSetupCollection` for each ODAU, followed by an invocation of `OptotrakSetupCollection`. The Optotrak System must be setup last.**

---

See Sample Program 12 on the API CD for an example of code that uses this routine.

## See Also

`OptotrakSetupCollectionFromFile`, `OdauSetupCollection`

## 9.6.7 `OdauSetupCollectionFromFile`

### Function

Allows you to setup ODAU collection parameters from information stored in an ASCII file. This routine applies to both ODAU I and ODAU II devices.

### Prototype

```
int OdauSetupCollectionFromFile( char *pszCollectFile )
```

### Parameters

`pszCollectFile` points to a null terminated string specifying the ODAU collection parameter file to use.

### Description

This routine allows you to set up ODAU collection parameters from information in an ASCII file. These settings may be in the same file as the settings used for `OptotrakSetupCollectionFromFile`, or in a different file. The parameters should be in this format:

```
[Odau 01]
FrameFrequency = 2000.0
NumberOfChannels = 3
Gain = 1
DigitalModeFlag = 0
bDifferential = FALSE
ScanFrequency = 90000.0
```



## Parameter File Details:

---

**Note** The ASCII file must include a separate section for each ODAU unit in the system. Each section must start with the heading [OdaU 01] (or [OdaU 02], [OdaU 03] or [OdaU 04]). All parameters are optional and defaults are used if parameters are omitted. The order of the parameters in the file is not important.

---

**FrameFrequency** is the frequency for ODAU collection, in Hz. Bounds are [1,27000] for ODAU I, and [1,100000] for ODAU II. Default is 100Hz.

**NumberOfChannels** is the number of channels of analog data to collect. Bounds are [1,256]. Without the use of multiplexer boards, the maximum values should be 16 for single-ended collections and 8 for differential collections. Multiplexer boards can only be used with an ODAU II, not an ODAU I. Default is 8.

**Gain** is the gain used for the analog channels. Range is [1,100], but valid values depend on the model of ODAU in use. For ODAU I, valid values are 1, 2, 4 and 8. For the ODAU II, valid values are 1, 5, 10 and 100. Default is 1.

**DigitalModeFlag** indicates whether there is digital input and/or output, and, if using an ODAU II, whether or not multiplexer boards are in use. Range is [0,255]. Default is 0. Valid values are in [Table 9-10 on page 153](#).

**Table 9-10: Values for Digital Mode Flag**

| Digital Mode Flag | Port A       | Port B       |
|-------------------|--------------|--------------|
| 0                 | off          | off          |
| 17                | input        | input        |
| 33                | input        | output       |
| 34                | output       | output       |
| 51                | input/output | input/output |
| 4                 | multiplexer  | off          |
| 20                | multiplexer  | input        |
| 36                | multiplexer  | output       |


---

**Note** Output readings for ports A and B can be obtained with a DigitalModeFlag setting of either 34 or 51, but a setting of 51 must be used to obtain collected digital data as well. Also, if one or both ports are set for output, then they should not have external input signals applied to them. Even if they are set to also read input, they should only read back the values that are set by the output function.

---

### Description of Ports A and B:

For an ODAU I, port A is bits 0 to 7 (pins 18, 5, 17, 4, 16, 3, 15, and 2, respectively, on a 25-pin D-shell), and port B is bits 8 to 15 (pins 25, 12, 24, 11, 23, 10, 22, and 9, respectively). The pin layout of the 25-pin D-shell is displayed in [Figure 9-2 on page 154](#) and the pin assignments are indicated in [Figure 9-11 on page 154](#).



**Figure 9-2: Numbering on the ODAU I 25-pin D-shell Plug**

For an ODAU II, port A consists of pins 25, 27, 29 and 31 of the 50-pin front panel connector, while port B consists of pins 26, 28, 30 and 32. Pin 1 of the connector is in the upper left corner as you look at the front of the ODAU II. Pin 2 is beneath pin 1. All odd numbered pins are in the top row. The signal on each pin is shown in [Figure 9-13 on page 156](#).

**bDifferential** [ODAU II only] If set to FALSE (0), then data collection will be in single-ended mode. If set to TRUE (1), then data collection will be in differential mode. Default is 0.

**ScanFrequency** [ODAU II only] This frequency sets the rate at which the channels are sampled within the frame. Valid values are [1, 100000]. Maximum recommended scan rates vary with the gain selected. Default is 90000. See [Table 9-12, “Recommended Scan Rates for ODAU II,” on page 155](#).

**Table 9-11: ODAU I Digital I/O Port Pin Configuration**

| Pin Number | Signal Name | Pin Number | Signal Name |
|------------|-------------|------------|-------------|
| 1          | Ground      | 14         | Vcc(+5V)    |
| 2          | Bit 7       | 15         | Bit 6       |
| 3          | Bit 5       | 16         | Bit 4       |
| 4          | Bit 3       | 17         | Bit 2       |
| 5          | Bit 1       | 18         | Bit 0       |
| 6          | Unused      | 19         | Unused      |

**Table 9-11: ODAU I Digital I/O Port Pin Configuration (Continued)**

| Pin Number | Signal Name                  | Pin Number | Signal Name                  |
|------------|------------------------------|------------|------------------------------|
| 7          | Unused                       | 20         | Unused                       |
| 8          | Trigger for Collection Frame | 21         | Trigger for Collection Start |
| 9          | Bit 15                       | 22         | Bit 14                       |
| 10         | Bit 13                       | 23         | Bit 12                       |
| 11         | Bit 11                       | 24         | Bit 10                       |
| 12         | Bit 9                        | 25         | Bit 8                        |
| 13         | Ground                       |            |                              |

**StreamData** indicates whether the ODAU is to send buffered data back by request only, or to automatically send buffered data back once data spooling is initiated.

Values:

**nStreamData** = 0: Send buffered data back at request only.

**nStreamData** = 1: Send buffered data back automatically.

Default is 1.

**CollectionTime** is the length of data collection, in seconds. It should be set to the same value for both the ODAU and the Optotrak System. Default is 2.

**PreTriggerTime** Not currently supported - it must be 0.

**Table 9-12: Recommended Scan Rates for ODAU II**

| Gain | Maximum Recommended Scan Frequency (Hz) |
|------|-----------------------------------------|
| 1    | 90,000                                  |
| 5    | 90,000                                  |
| 10   | 70,000                                  |
| 100  | 20,000                                  |

## See Also

OptotrakSetupCollectionFromFile

Table 9-13: ODAU II Front Panel Port Pin Configuration

| Pin Number | Signal Name | Pin Number | Signal Name  |
|------------|-------------|------------|--------------|
| 1          | AI GND      | 2          | AI GND       |
| 3          | ACH0        | 4          | ACH8         |
| 5          | ACH1        | 6          | ACH9         |
| 7          | ACH2        | 8          | ACH10        |
| 9          | ACH3        | 10         | ACH11        |
| 11         | ACH4        | 12         | ACH12        |
| 13         | ACH5        | 14         | ACH13        |
| 15         | ACH6        | 16         | ACH14        |
| 17         | ACH7        | 18         | ACH15        |
| 19         | AI SENSE    | 20         | DAC0 OUT     |
| 21         | DAC1 OUT    | 22         | EXT REF      |
| 23         | AO GND      | 24         | DIG GND      |
| 25         | ADIO0       | 26         | BDIO0        |
| 27         | ADIO1       | 28         | BDIO1        |
| 29         | ADIO2       | 30         | BDIO2        |
| 31         | ADIO3       | 32         | BDIO3        |
| 33         | DIG GND     | 34         | + 5 V        |
| 35         | + 5 V       | 36         | SCANCLK      |
| 37         | !EXT STROBE | 38         | !START TRIG  |
| 39         | !STOP TRIG  | 40         | !EXTCONV     |
| 41         | TIMER CLK 2 | 42         | TIMER GATE 2 |
| 43         | TIMER OUT 2 | 44         | N/C          |
| 45         | N/C         | 46         | !FRAME CLK   |
| 47         | N/C         | 48         | N/C          |
| 49         | N/C         | 50         | FOUT         |

## 9.7 Real-time Data Retrieval Routines

See the “Real-time Rigid Body Programmer’s Guide” on page 51 and “Real-time Data Types” on page 237 for more general information on data types and rigid body routines.

### 9.7.1 DataGetLatestCentroid

#### Function

Requests and receives the latest frame of centroid (raw) data.

#### Prototype

```
int DataGetLatestCentroid( unsigned int *puFrameNumber,
                          unsigned int *puElements,
                          unsigned int *puFlags,
                          void          *pDataDest )
```

#### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of markers for which there is data in the received frame.

**puFlags** indicates the current status of the Optotrak System.

The `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specific to the Optotrak Certus System Control Unit. If this flag is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. In this case, if `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

DataGetLatestCentroid retrieves the latest frame of centroid data from the Optotrak System. The routine also returns the frame number associated with the data, the number of markers for which there are data in the frame, and the status of the Optotrak System data buffer.

The centroid data can be converted to 3D data or 6D transformation data by using conversion routines such as OptotrakCovertRawTo3D or OptotrakConvertTransforms. To optimize real-time capabilities, use the DataGetLatest*type* routines.

This routine copies the centroid data frame from the Optotrak System directly into a memory block set aside by the application program.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Sample Formats” on page 68. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

## See Also

DataGetNextCentroid, RequestLatestCentroid, RequestNextCentroid

### 9.7.2 DataGetLatest3D

#### Function

Requests and receives the latest frame of 3D data.

#### Prototype

```
int DataGetLatest3D( unsigned int    *puFrameNumber,  
                   unsigned int    *puElements,  
                   unsigned int    *puFlags,  
                   void            *pDataDest )
```

#### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of 3D positions (number of markers) in the frame of data received.

**puFlags** indicates the current status of the Optotrak System.

The `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specific to the Optotrak Certus System Control Unit. If this flag is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

`DataGetLatest3D` retrieves the latest frame of 3D data from the Optotrak System. The routine also returns the frame number associated with the data, the number of 3D marker positions in the frame, and status of the Optotrak System data buffer.

The conversion of the raw data to 3D data can be processed in the Optotrak System or in the host computer depending on the status `OPTO_CONVERT_ON_HOST` flag set in the routine `OptotrakSetProcessingFlags`.



**Warning!**

---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 242 and “[Sample Formats](#)” on page 68. This routine copies the 3D positions directly into a memory block. If the memory block size is too small, the 3D data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

The 3D data can be converted to 6D transformation data by using the conversion routine `OptotrakConvertTransforms`. To optimize real-time capabilities, use the `DataGetLatesttype` routines.

See Sample Program 6 on the API CD for an example of code that uses this routine.

## See Also

`DataGetNext3D`, `RequestLatest3D`, `RequestNext3D`, `OptotrakSetProcessingFlags`

### 9.7.3 DataGetLatestOdaRaw

#### Function

Requests and receives the latest frame of raw data from the specified ODAU.

#### Prototype

```
int DataGetLatestOdaRaw( int          nOdaId,  
                        unsigned int *puFrameNumber,  
                        unsigned int *puElements,  
                        unsigned int *puFlags,  
                        void          *pDataDest )
```

#### Parameters

**nOdaId** specifies the ODAU device from which to retrieve the frame of raw data. Use the defined constants ODAU1, ODAU2, ODAU3, and ODAU4.

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of data elements in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is specified during a call to the OptotrakSetupCollection routine, the puFlags parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, puFlags will contain OPTO\_TOOL\_CONFIG\_CHANGED\_FLAG. If there is new switch data available, puFlags will contain OPTO\_SWITCH\_DATA\_CHANGED\_FLAG. To retrieve the switch data, call RetrieveSwitchData.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is not specified, puFlags will be set to indicate the current status of the Optotrak System data buffer. If puFlags is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received from the ODAU.

#### Description

DataGetLatestOdaRaw retrieves the latest frame of raw data from the specified ODAU device. This data frame consists of an integer for each analog channel being sampled, followed by an integer for the digital port data, if it is being sampled. The routine also returns



- the frame number associated with the data
- the number of data elements in the frame
- the status of the ODAU data buffer.

The number of data elements in the frame equals the number of channels if digital data is not sampled, or the number of channels plus 1, if digital data is sampled.




---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 247. This routine copies the raw data directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---



---

**Note** If digital data is requested from an ODAU II, it is stored in the upper 8 bits of the last unsigned integer subitem.

---

See Sample Program 12 on the API CD for an example of code that uses this routine.

## See Also

DataGetNextOdaRaw, RequestLatestOdaRaw, RequestNextOdaRaw

### 9.7.4 DataGetLatestRaw

#### Function

Requests and receives the latest frame of full raw data.

#### Prototype

```
int DataGetLatestRaw( unsigned int    *puFrameNumber,
                    unsigned int    *puElements,
                    unsigned int    *puFlags,
                    void             *pDataDest )
```

#### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of markers for which there is data in the received frame.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

### Description

`DataGetLatestRaw` retrieves the latest frame of full raw data (centroid data and status) from the Optotrak System. The routine also returns the frame number associated with the data, the number of markers for which there are data in the frame, and the status of the Optotrak System data buffer.



**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 241. This routine copies the full raw data directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

See Sample Program 3 on the API CD for an example of code that uses this routine.

### See Also

`DataGetNextRaw`, `RequestLatestRaw`, `RequestNextRaw`

## 9.7.5 DataGetLatestTransforms

### Function

Requests and receives the latest frame of 6D rigid body transformation data, as well as the associated 3D position data.

### Prototype

```
int DataGetLatestTransforms( unsigned int   *puFrameNumber,
                           unsigned int   *puElements,
                           unsigned int   *puFlags,
                           void           *pDataDest )
```

### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of rigid body transformations in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

### Description

`DataGetLatestTransforms` retrieves the:

- latest frame of rigid body transformation 6D data
- 3D data used in the calculation of the rigid body transformations
- frame number associated with the data

- number of rigid body transformations in the frame
- status of the Optotrak System data buffer

The number of valid 3D markers returned from `DataGetLatestTransforms` may be fewer than the actual number of markers visible to the Optotrak System. The algorithm may have omitted one or more markers from the best fit routine.

The 6D data is followed immediately by the 3D data, see “[Organization of Rigid Body Transformation Data](#)” on page 242. If you only need the rigid body transformation data, or wish to collect the 6D and 3D data in separate buffers, use the routine `DataGetLatestTransforms2` instead.

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines `RigidBodyAdd` or `RigidBodyAddFromFile`.

The conversion of the 3D data to 6D transformation can be processed in the Optotrak System or in the host computer depending on the status `OPTO_RIGID_ON_HOST` flag set in the routine `OptotrakSetProcessingFlags`.



**Warning!**

---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 246 and “[Sample Formats](#)” on page 68. This routine copies the rigid body positions directly into a memory block. If the memory block size is too small, the rigid body data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** The parameter `pdatadest` includes the bit flag `OPTOTRAK_UNDETERMINED_FLAG`. This flag will be set if the transformation could not be determined. This flag must be checked for each transformation.

---

See “[Flags and Settings Associated with Rigid Bodies](#)” on page 295 for information on flag settings and where to set error parameters.

See Sample Program 9 on the API CD for an example of code that uses this routine.

## See Also

DataGetLatestTransforms2, DataGetNextTransforms, RequestLatestTransforms, RequestNextTransforms, OptotrakSetProcessingFlags, RigidBodyAdd, RigidBodyAddFromFile

### 9.7.6 DataGetLatestTransforms2

#### Function

Requests and receives the latest frame of 6D rigid body transformation data, and the associated 3D position data, in separate buffers.

#### Prototype

```
int DataGetLatestTransforms2(unsigned int          *puFrameNumber,
                           unsigned int          *puElements,
                           unsigned int          *puFlags,
                           struct OptotrakRigidStruct *pDataDest6D,
                           Position3d           *pDataDest3D )
```

#### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of rigid body transformations in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is specified during a call to the OptotrakSetupCollection routine, the puFlags parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, puFlags will contain OPTO\_TOOL\_CONFIG\_CHANGED\_FLAG. If there is new switch data available, puFlags will contain OPTO\_SWITCH\_DATA\_CHANGED\_FLAG. To retrieve the switch data, call RetrieveSwitchData.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is not specified, puFlags will be set to indicate the current status of the Optotrak System data buffer. If puFlags is set to a non-zero value, then a buffering error has occurred.

**pDataDest6D** points to memory set aside by the application program to store the 6D rigid body transformation data returned by the Optotrak System.

`pDataDest3D` points to memory set aside by the application program to store the associated 3D position data returned by the Optotrak System. If `pDataDest3D` is set to NULL, then the 3D data is not returned.

## Description

Retrieves data from the Optotrak System. It is similar to `DataGetLatestTransforms`, except that it stores the retrieved data into two buffers. The latest frame of rigid body transformation data, and the associated 3D data used in the calculation of the rigid body transformations are in separate buffers for easier access. You can also choose not to have the 3D data returned by passing a NULL value for `pDataDest3D`. The routine also returns the frame number associated with the data, the number of rigid body transformations in the frame, and status of the Optotrak System data buffer.

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to TRUE in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines `RigidBodyAdd` or `RigidBodyAddFromFile`.

The conversion of the 3D data to 6D transformation can be processed in the Optotrak System or in the host computer depending on the status `OPTO_RIGID_ON_HOST` flag set in the routine `OptotrakSetProcessingFlags`.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See [“Size Calculation” on page 242](#), [“Size Calculation” on page 246](#) and [“Sample Formats” on page 68](#). This routine copies the 6D rigid body transformation data and 3D position data directly into a memory blocks. If the memory block sizes are too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** The parameter `pdatadest6D` includes the bit flag `OPTOTRAK_UNDETERMINED_FLAG`. This flag will be set if the transformation could not be determined. This flag must be checked for each transformation.

---

See [“Flags and Settings Associated with Rigid Bodies” on page 295](#) for information on flag settings and where to set error parameters.

## See Also

DataGetNextTranforms2, DataGetLatestTransforms, RequestLatestTransforms, RequestNextTransforms, OptotrakSetProcessingFlags, RigidBodyAdd, RigidBodyAddFromFile

### 9.7.7 DataGetNext3D

#### Function

Requests and receives the next frame of 3D data.

#### Prototype

```
int DataGetNext3D( unsigned int *puFrameNumber,
                  unsigned int *puElements,
                  unsigned int *puFlags,
                  void          *pDataDest )
```

#### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of 3D positions (number of markers) in the frame of data received.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

DataGetNext3D is similar to DataGetLatest3D, except that it retrieves the next available frame of 3D data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. The routine also returns the frame number associated with the data, the number of 3D marker positions in the frame, and the status of the Optotrak System data buffer.

DataGetNext3D is the same as DataGetLatest3D with the flag OPTOTRAK\_GET\_NEXT\_FRAME\_FLAG set in the collection.

The conversion of the raw data to 3D data can be processed in the Optotrak System or in the host computer depending on the status OPTO\_CONVERT\_ON\_HOST flag set in the routine OptotrakSetProcessingFlags.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 242. This routine copies the 3D positions directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

The 3D data can be converted to 6D transformation data by using the conversion routine OptotrakConvertTransforms. To optimize real-time capabilities, use the DataGetLatest*type* routines.

## See Also

DataGetLatest3D, RequestLatest3D, RequestNext3D, OptotrakSetProcessingFlags

### 9.7.8 DataGetNextCentroid

#### Function

Requests and receives the next frame of centroid data.

#### Prototype

```
int DataGetNextCentroid( unsigned int *puFrameNumber,  
                        unsigned int *puElements,  
                        unsigned int *puFlags,  
                        void          *pDataDest )
```



## Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of markers for which there is data in the received frame.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

`DataGetNextCentroid` is similar to `DataGetLatestCentroid`, except that it retrieves the next available frame of centroid data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. The routine also returns the frame number associated with the data, the number of markers for which there is data in the frame, and the status of the Optotrak System data buffer.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 241 and “Sample Formats” on page 68. This routine copies the centroid data frame directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

## See Also

`DataGetLatestCentroid`, `RequestLatestCentroid`, `RequestNextCentroid`,

## 9.7.9 DataGetNextOdaRaw

### Function

Requests and receives the next frame of raw data from the specified ODAU.

### Prototype

```
int DataGetNextOdaRaw( int          nOdaId,  
                      unsigned int *puFrameNumber,  
                      unsigned int *puElements,  
                      unsigned int *puFlags,  
                      void          *pDataDest )
```

### Parameters

**nOdaId** specifies which ODAU device to retrieve the frame of raw data. Use the defined constants ODAU1, ODAU2, ODAU3, and ODAU4.

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of data elements in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is specified during a call to the OptotrakSetupCollection routine, the puFlags parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, puFlags will contain OPTO\_TOOL\_CONFIG\_CHANGED\_FLAG. If there is new switch data available, puFlags will contain OPTO\_SWITCH\_DATA\_CHANGED\_FLAG. To retrieve the switch data, call RetrieveSwitchData.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is not specified, puFlags will be set to indicate the current status of the Optotrak System data buffer. If puFlags is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received from the ODAU.

### Description

DataGetNextOdaRaw is similar to DataGetLatestOdaRaw, except that it retrieves the next available frame of raw data from the specified ODAU device instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the ODAU frame

rate. The data frame consists of an integer for each analog channel being sampled, followed by an integer for the digital port data, if it is being sampled. The routine also returns the frame number associated with the data, the number of data elements in the frame, and status of the ODAU data buffer.

DataGetNextOdaRaw is the same as DataGetLatestOdaRaw with the flag OPTOTRAK\_GET\_NEXT\_FRAME\_FLAG set in the collection.

If digital data is requested from an ODAU II, it is stored in the upper 8 bits of the last unsigned integer subitem.




---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 247 and “Sample Formats” on page 68. This routine copies the ODAU raw data directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

## See Also

RequestLatestOdaRaw, RequestNextOdaRaw, DataGetLatest3D

### 9.7.10 DataGetNextRaw

#### Function

Requests and receives the next frame of full raw data.

#### Prototype

```
int DataGetNextRaw( unsigned int *puFrameNumber,
                  unsigned int *puElements,
                  unsigned int *puFlags,
                  void          *pDataDest )
```

#### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of markers for which there is data in the received frame.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

`pDataDest` points to memory set aside by the application program to store the data received by the Optotrak System.

### Description

`DataGetNextRaw` is similar to `DataGetLatestRaw`, except that it retrieves the next available frame of full raw data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. The routine also returns the frame number associated with the data, the number of markers for which there is data in the frame, and the status of the Optotrak System data buffer.

`DataGetNextRaw` is the same as `DataGetLatestRaw` with the flag `OPTOTRAK_GET_NEXT_FRAME_FLAG` set in the collection.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 241 and “[Sample Formats](#)” on page 68. This routine copies the raw data frame directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

### See Also

`DataGetLatestRaw`, `RequestLatestRaw`, `DataGetNextCentroid`,  
`DataGetLatestCentroid`

## 9.7.11 DataGetNextTransforms

### Function

Requests and receives the next frame of 6D rigid body transformation data, as well as its associated 3D position data.

## Prototype

```
int DataGetNextTransforms( unsigned int  *puFrameNumber,
                          unsigned int  *puElements,
                          unsigned int  *puFlags,
                          void          *pDataDest )
```

## Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of rigid body transformations in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

`DataGetNextTransforms` is similar to `DataGetLatestTransforms`, except that it retrieves both the next available frame of rigid body transformation 6D data and the 3D data used in the calculation of the rigid body transformations, instead of the latest frame. This data is received from the Optotrak System. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate.

The 6D data is followed immediately by the 3D data. If you only need the rigid body transformation data, or wish to collect the 6D and 3D data in separate buffers, use the routine `DataGetNextTransforms2` instead.

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak

3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

DataGetNextTransforms is the same as DataGetLatestTransforms with the flag OPTOTRAK\_GET\_NEXT\_FRAME\_FLAG set in the collection.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines RigidBodyAdd or RigidBodyAddFromFile.

The conversion of the 3D data to 6D transformation data can be processed in the Optotrak System or in the host computer depending on the status OPTO\_RIGID\_ON\_HOST flag set in the routine OptotrakSetProcessingFlags.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 246 and “Sample Formats” on page 68. This routine copies the rigid body transformation data directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** The parameter pdatadest includes the bit flag OPTOTRAK\_UNDETERMINED\_FLAG. This flag will be set if the transformation could not be determined. This flag must be checked for each transformation.

---

See “Flags and Settings Associated with Rigid Bodies” on page 295 for information on flag settings and where to set error parameters.

### See Also

DataGetNextTransforms2, DataGetLatestTransforms, OptotrakSetProcessingFlags, RigidBodyAdd, RigidBodyAddFromFile, DataReceiveLatestTransforms, RequestLatestTransforms.

## 9.7.12 DataGetNextTransforms2

### Function

Requests and receives the next frame of 6D rigid body transformation data, and its associated 3D position data, in separate buffers.

### Prototype

```
int DataGetNextTransforms2( unsigned int                               *puFrameNumber,
```

```

unsigned int          *puElements,
unsigned int          *puFlags,
struct OptotrakRigidStruct *pDataDest6D,
Position3d           *pDataDest3D )

```

## Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of rigid body transformations in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest6D** points to memory set aside by the application program to store the 6D rigid body transformation data returned by the Optotrak System.

**pDataDest3D** points to memory set aside by the application program to store the associated 3D position data returned by the Optotrak System. If `pDataDest3D` is set to `NULL`, then the 3D data is not returned.

## Description

`DataGetNextTransforms2` is similar to `DataGetLatestTransforms2`, except that it retrieves the next available frame of rigid body transformation 6D data and the 3D data used in the calculation of the rigid body transformations, instead of the latest frame. The data is retrieved from the Optotrak System. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate.

`DataGetNextTransforms2` is also similar to `DataGetNextTransforms`, except that the 6D and 3D data is retrieved in separate buffers for easier access. You can also choose not to have the 3D data returned by passing a `NULL` value for `pDataDest3D`. The routine also returns the frame number associated with the data, the number of rigid body transformations in the frame, and status of the Optotrak System data buffer.

The 6D data is followed immediately by the 3D data. If you only need the rigid body transformation data, or wish to collect the 6D and 3D data in separate buffers, use the routine `DataGetLatestTransforms2` instead.

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines `RigidBodyAdd` or `RigidBodyAddFromFile`.

The conversion of the 3D data to 6D transformation data can be processed in the Optotrak System or in the host computer depending on the status `OPTO_RIGID_ON_HOST` flag set in the routine `OptotrakSetProcessingFlags`.

`DataGetNextTransforms2` is the same as `DataGetLatestTransforms2` with the flag `OPTOTRAK_GET_NEXT_FRAME_FLAG` set in the collection.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 242, “[Size Calculation](#)” on page 246 and “[Sample Formats](#)” on page 68. This routine copies the 6D rigid body transformation data and 3D position data directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** The parameter `pdatadest6D` includes the bit flag `OPTOTRAK_UNDETERMINED_FLAG`. This flag will be set if the transformation could not be determined. This flag must be checked for each transformation.

---

See “[Flags Affecting Rigid Bodies](#)” on page 297 for information on flag settings and where to set error parameters.

### See Also

`DataGetLatestTransforms`, `RequestLatestTransforms`, `RigidBodyAdd`, `RigidBodyAddFromFile`, `RequestNextTransforms`



### 9.7.13 DataIsReady

#### Function

Determines if there is a frame of data to be received.

#### Prototype

```
int DataIsReady( void )
```

#### Parameters

None.

#### Return Value

**Non-zero:** if there is a data frame waiting to be received.

**Zero:** if there is no data waiting to be received.

#### Description

DataIsReady is used by an application program to determine if a frame of requested data is ready to be received. If an application program attempts to retrieve data using one of the RequestNext*type* or RequestLatest*type* routines, it must call DataIsReady to determine whether or not the data is ready to be received, and then it should invoke the appropriate DataReceiveLatest*type* routine to receive the waiting data.

DataIsReady is different from all other routines; its return value indicates the presence of data to be received. If data is ready to be received, then DataIsReady returns a non-zero value, otherwise it returns a value of zero.

If the bit flag OPTO\_LIB\_POLL\_REAL\_DATA has been set in the OptotrakSetProcessingFlags routine, DataIsReady blocks until the data is ready, rather than returning FALSE.

#### See Also

DataReceiveLatest3D, DataReceiveLatestRaw, DataReceiveLatestTransforms, DataReceiveLatestOdaRaw, OptotrakSetProcessingFlags

## 9.7.14 DataReceiveLatest3D

### Function

Receives the currently waiting frame of 3D data.

### Prototype

```
int DataReceiveLatest3D(unsigned int *puFrameNumber
                      unsigned int *puElements,
                      unsigned int *puFlags,
                      void          *pDataDest )
```

### Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of 3D positions (number of markers) in the frame of data received.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

### Description

`DataReceiveLatest3D` receives data that was requested previously using `RequestLatest3D`. Invoke this routine only after a call to the `DataIsReady` routine returns `TRUE`. The routine also returns the frame number associated with the data, the number of 3D marker positions in the frame, and status of the Optotrak System data buffer.



Warning!

---

You must ensure that the memory block reserved for the API to store the data is the correct size. “[Size Calculation](#)” on page 242 and “[Sample Formats](#)” on page 68. This routine copies the 3D positions directly into a memory block. If the memory block size is too small, the 3D data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.

---



---

**Note** Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns TRUE, the application program does not know which data type it is receiving.

---

See Sample Program 5 on the API CD for an example of code that uses this routine.

## See Also

`DataIsReady`, `RequestLatest3D`

### 9.7.15 `DataReceiveLatestCentroid`

#### Function

Receives the currently waiting frame of centroid (raw) data.

#### Prototype

```
int DataReceiveLatestCentroid( unsigned int *puFrameNumber,
                             unsigned int *puElements,
                             unsigned int *puFlags,
                             void          *pDataDest )
```

#### Parameters

`puFrameNumber` is set to the frame number associated with the frame of data received.

`puElements` is set to the number of markers for which there is data in the received frame.

`puFlags` indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value

indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, puFlags will contain OPTO\_TOOL\_CONFIG\_CHANGED\_FLAG. If there is new switch data available, puFlags will contain OPTO\_SWITCH\_DATA\_CHANGED\_FLAG. To retrieve the switch data, call RetrieveSwitchData.

If OPTOTRAK\_SWITCH\_AND\_CONFIG\_FLAG is not specified, puFlags will be set to indicate the current status of the Optotrak System data buffer. If puFlags is set to a non-zero value, then a buffering error has occurred.

pDataDest points to memory set aside by the application program to store the data received by the Optotrak System.

### Description

DataReceiveLatestCentroid is used by an application program to receive centroid data that was requested previously using RequestLatestCentroid. This routine should be invoked only after a call to the DataIsReady routine returns TRUE. The routine also returns the frame number associated with the data, the number of markers for which there are data in the frame, and the status of the Optotrak System data buffer.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 241 and “Sample Formats” on page 68. This routine copies the raw data frame directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine DataIsReady returns TRUE, the application program does not know which data type it is receiving.

---

### See Also

DataIsReady, RequestLatestCentroid, RequestNextCentroid

## 9.7.16 DataReceiveLatestOdaRaw

### Function

Receives the currently waiting frame of raw data from the specified ODAU.

## Prototype

```
int DataReceiveLatestOdaRaw( unsigned int *puFrameNumber,
                           unsigned int *puElements,
                           unsigned int *puFlags,
                           int          *pDataDest )
```

## Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of data elements in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received from the ODAU.

## Description

`DataReceiveLatestOdaRaw` receives data that was requested previously using `RequestLatestOdaRaw`. Invoke this routine only after a call to the `DataIsReady` routine returns `TRUE`. The routine also returns the frame number associated with the data, the number of data elements in the frame, and status of the Optotrak System data buffer.

---

**Note** If digital data is requested from an ODAU II, it is stored in the upper 8 bits of the last unsigned integer subitem.

---



Warning!

You must ensure that the memory block reserved for the API to store the data is the correct size. See [“Size Calculation” on page 247](#) and [“Sample Formats” on page 68](#). This routine copies the data from the ODAU directly into a memory block. If the memory block size is too

**small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DatalsReady` returns `TRUE`, the application program does not know which data type it is receiving.

---

## See Also

`DatalsReady`, `RequestLatestOdaRaw`

### 9.7.17 `DataReceiveLatestRaw`

#### Function

Receives the currently waiting frame of full raw data.

#### Prototype

```
int DataReceiveLatestRaw( unsigned int *puFrameNumber,  
                        unsigned int *puElements,  
                        unsigned int *puFlags,  
                        void          *pDataDest )
```

#### Parameters

`puFrameNumber` is set to the frame number associated with the frame of data received.

`puElements` is set to the number of markers for which there is data in the received frame.

`puFlags` indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

`pDataDest` points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

`DataReceiveLatestRaw` receives full raw data that was requested previously using `RequestLatestRaw`. Invoke this routine only after a call to the `DataIsReady` routine returns `TRUE`. The routine also returns the frame number associated with the data, the number of markers for which there are data in the frame, and the status of the Optotrak System data buffer.




---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “Size Calculation” on page 241 and “Sample Formats” on page 68. This routine copies the raw data frame directly into a memory block. If the memory block size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---



---

**Note** Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns `TRUE`, the application program does not know which data type it is receiving.

---

## See Also

`DataIsReady`, `RequestLatestRaw`

### 9.7.18 DataReceiveLatestTransforms

#### Function

Receives the currently waiting frame of 6D rigid body transformation data and its associated 3D position data.

#### Prototype

```
int DataReceiveLatestTransforms( unsigned int *puFrameNumber,
                                unsigned int *puElements,
```

```
unsigned int *puFlags,  
void        *pDataDest )
```

## Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of rigid body transformations in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest** points to memory set aside by the application program to store the data received by the Optotrak System.

## Description

`DataReceiveLatestTransforms` receives the latest rigid body transformation 6D data that was requested previously using `RequestLatestTransforms`, immediately followed by the 3D data used in the calculation of the rigid body transformations. This data is received from the Optotrak System. If you only need the rigid body transformation data, or wish to collect the 6D and 3D data in separate buffers, use the routine `DataReceiveLatestTransforms2` instead. Invoke this routine only after a call to the `DataIsReady` routine returns `TRUE`. The routine also returns the frame number associated with the data, the number of rigid body transformations in the frame, and status of the Optotrak System data buffer.



**Warning!**

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 246 and “[Sample Formats](#)” on page 68. This routine copies the rigid body transformation data directly into a memory block. If the memory block**



---

**size is too small, the data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

**Note** Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns `TRUE`, the application program does not know which data type it is receiving.

---

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routine `RigidBodyAdd` or `RigidBodyAddFromFile`

## See Also

`DataIsReady`, `RequestLatestTransforms`, `RigidBodyAdd`, `RigidBodyAddFromFile`, `DataReceiveLatestTransforms2`

## 9.7.19 DataReceiveLatestTransforms2

### Function

Receives the currently waiting frame of 6D rigid body transformation data, and its associated 3D position data, in separate buffers.

### Prototype

```
int DataReceiveLatestTransforms2 (
    unsigned int      *puFrameNumber,
    unsigned int      *puElements,
    unsigned int      *puFlags,
    struct OptotrakRigidStruct *pDataDest6D,
    Position3d        *pDataDest3D )
```

## Parameters

**puFrameNumber** is set to the frame number associated with the frame of data received.

**puElements** is set to the number of rigid body transformations in the received frame of data.

**puFlags** indicates the current status of the Optotrak System.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, the `puFlags` parameter will contain a value indicating that the device configuration has changed or that switch data is available. If the device configuration has changed, `puFlags` will contain `OPTO_TOOL_CONFIG_CHANGED_FLAG`. If there is new switch data available, `puFlags` will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. To retrieve the switch data, call `RetrieveSwitchData`.

If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is not specified, `puFlags` will be set to indicate the current status of the Optotrak System data buffer. If `puFlags` is set to a non-zero value, then a buffering error has occurred.

**pDataDest6D** points to memory set aside by the application program to store the 6D rigid body transformation data returned by the Optotrak System.

**pDataDest3D** points to memory set aside by the application program to store the associated 3D position data returned by the Optotrak System. If `pDataDest3D` is set to `NULL`, then the 3D data is not returned.

## Description

`DataReceiveLatestTransforms2` is similar to `DataReceiveLatestTransforms`, except that it receives the waiting rigid body transformation data, and the associated 3D data used in the calculation of the rigid body transformations, in separate buffers for easier access. The data is received from the Optotrak System. You can also choose not to have the 3D data returned by passing a `NULL` value for `pDataDest3D`. Invoke this routine only after a call to the `DataIsReady` routine returns `TRUE`. The routine also returns the frame number associated with the data, the number of rigid body transformations in the frame, and status of the Optotrak System data buffer.



---

**You must ensure that the memory block reserved for the API to store the data is the correct size. See “[Size Calculation](#)” on page 242, “[Size Calculation](#)” on page 246, and “[Sample Formats](#)” on page 68. This routine copies the 6D rigid body transformation data and 3D position data directly into a memory block. If the memory block size is too small, the 3D data is copied to an invalid memory area on the host computer, causing unpredictable behaviour. This may include a system crash.**

---

---

**Note** Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns `TRUE`, the application program does not know which data type it is receiving.

---

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines `RigidBodyAdd` or `RigidBodyAddFromFile`.

### See Also

`DataIsReady`, `RequestLatestTransforms`, `RigidBodyAdd`, `RigidBodyAddFromFile`, `DataReceiveLatestTransforms`

## 9.7.20 ReceiveLatestData

---

**Note** This function is now obsolete. Use `DataReceiveLatest3D`, `DataReceiveLatestRaw`, `DataReceiveLatestTransforms`, `DataReceiveLatestTransforms2`, `DataReceiveLatestCentroid` or `DataReceiveLatestOdaRaw` instead.

---

## 9.7.21 RetrieveSwitchData

### Function

Requests and retrieves the data generated from the switch on an Optotrak Certus strober that accompanied the latest real-time data received from the Optotrak System.

### Prototype

```
int RetrieveSwitchData( int      nNumSwitches,
                      boolean *pbSwitchData )
```

## Parameters

**nNumSwitches** indicates the number of switches.

**pbSwitchData** is the data buffer to which switch data will be written.

## Description

The Optotrak System returns switch data to the API with every frame of real-time data. If `OPTOTRAK_SWITCH_AND_CONFIG_FLAG` is specified during a call to the `OptotrakSetupCollection` routine, every API routine that retrieves real-time data will indicate whether or not new switch data is available. If there is new switch data available, the `puFlags` variable in the `DataGetLatest` routine will contain `OPTO_SWITCH_DATA_CHANGED_FLAG`. Use the `RetrieveSwitchData` routine to retrieve the new switch data if required.

---

**Note** If an application program ends with a switch pressed, the System Control Unit will continue to monitor the changes in switch status. If you start a new application program without re-initializing and re-determining the configurations (using `TransputerLoadSystem` and `TransputerDetermineSystemCfg`), the application program will assume a change in switch status.

---

See Optotrak Certus Sample Program 14 on the API CD for an example of code that uses this routine.

### 9.7.22 RequestLatest3D

#### Function

Requests the latest frame of 3D data from the Optotrak System.

#### Prototype

```
int RequestLatest3D( void )
```

#### Parameters

None.

## Description

RequestLatest3D requests that the Optotrak System determine the latest frame of 3D data and send it back to the application program. Unlike the related routine DataGetLatest3D, this routine does not block and wait for the data to be returned.

The routine DataIsReady determines if there is 3D data available to be received. When the data is ready to be received, DataIsReady returns a non-zero value (TRUE). The application program should then use the routine DataReceiveLatest3D to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine DataIsReady returns TRUE, the application program does not know which data type it is receiving.**

---

## See Also

RequestNext3D, DataIsReady, DataReceiveLatest3D, DataGetLatest3D

### 9.7.23 RequestLatestCentroid

#### Function

Requests the latest frame of centroid data.

#### Prototype

```
int RequestLatestCentroid ( void )
```

#### Parameters

None.

#### Description

RequestLatestCentroid requests the Optotrak System to return the latest frame of centroid data. Unlike the related routine DataGetLatestCentroid, this routine does not block and wait for the data to be returned. The application program must determine if the raw data is ready, and then receive it.

The routine `DataIsReady` determines if there is raw data available to be received. When the data is ready to be received, `DataIsReady` returns a non-zero value (TRUE). The application program should then use the routine `DataReceiveLatestCentroid` to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



---

**Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns TRUE, the application program does not know which data type it is receiving.**

---

## See Also

`RequestNextCentroid`, `DataIsReady`, `DataReceiveLatestCentroid`,  
`DataGetLatestCentroid`

## 9.7.24 RequestLatestOdaRaw

### Function

Requests the latest frame of raw data from the specified ODAU device.

### Prototype

```
int RequestLatestOdaRaw( int nOdaId )
```

### Parameters

`nOdaId` specifies the ODAU device from which to request the frame of raw data. Use one of the defined constants `ODAU1`, `ODAU2`, `ODAU3` and `ODAU4`.

### Description

`RequestLatestOdaRaw` requests that the specified ODAU return the latest frame of raw data to the application program. Unlike the related routine `DataGetLatestOdaRaw`, this routine does not block and wait for the data to be returned. The application program determines if the ODAU data is ready, and then receive it.

The routine `DataIsReady` determines if there is data that can be received. When data is ready to be received, `DataIsReady` returns a non-zero value (TRUE). Use the routine

---

DataReceiveLatestOdaRaw to receive the data from the ODAU. Never request data without a followup call to receive the data.



---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine DataIsReady returns TRUE, the application program does not know which data type it is receiving.**

---

## See Also

DataIsReady, DataReceiveLatestOdaRaw, DataGetLatestOdaRaw, RequestNextOdaRaw

### 9.7.25 RequestLatestRaw

#### Function

Request the latest frame of full raw data.

#### Prototype

```
int RequestLatestRaw( void )
```

#### Parameters

None.

#### Description

RequestLatestRaw requests that the Optotrak System return the latest frame of full raw data. Unlike the related routine DataGetLatestRaw, this routine does not block and wait for the data to be returned. The application program must determine if the raw data is ready, and then receive it.

Use DataIsReady to determine if there is full raw data that can be received. When the data is ready to be received, DataIsReady returns a non-zero value (TRUE). Use the routine DataReceiveLatestRaw to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns `TRUE`, the application program does not know which data type it is receiving.**

---

## See Also

`DataIsReady`, `DataReceiveLatestRaw`, `DataGetLatestRaw`, `RequestNextRaw`

## 9.7.26 RequestLatestTransforms

### Function

Request the latest frame of rigid body transformation data.

### Prototype

```
int RequestLatestTransforms( void )
```

### Parameters

None.

### Description

`RequestLatestTransforms` requests that the Optotrak System determine the latest frame of rigid body 6D transformation data and return it to the application program. Unlike the related routines `DataGetLatestTransforms` and `DataGetLatestTransforms2`, this routine does not block and wait for the data to be returned. The application program determines if the transformation data is ready and then receives it.

The routine `DataIsReady` determines if there is 6D data to be received. When data is ready to be received, `DataIsReady` returns a non-zero value (`TRUE`). The application program should then use the routines `DataReceiveLatestTransforms` or `DataReceiveLatestTransforms2` to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other,**



---

**and the routine `DataIsReady` returns `TRUE`, the application program does not know which data type it is receiving.**

---

You may use this routine if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines `RigidBodyAdd` or `RigidBodyAddFromFile`.

### See Also

`DataIsReady`, `DataReceiveLatestTransforms`, `DataReceiveLatestTransforms2`, `DataGetLatestTransforms`, `DataGetLatestTransforms2`, `RequestNextTransforms`

## 9.7.27 RequestNext3D

### Function

Requests the next frame of 3D data from the Optotrak System.

### Prototype

```
int RequestNext3D( void )
```

### Parameters

None.

### Description

`RequestNext3D` is similar to `RequestLatest3D`, except that it requests the next available frame of 3D data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. Unlike the related routine `DataGetNext3D`, this routine does not block and wait for the data to be returned. The application program determines if the 3D data is ready, and then receives it.

The routine `DataIsReady` determines if there is 3D data available to be received. When the data is ready to be received, `DataIsReady` returns a non-zero value (TRUE). The application program should then use the routine `DataReceiveLatest3D` to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



---

**Do not send out another request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns TRUE, the application program does not know which data type it is receiving.**

---

### See Also

`DataIsReady`, `DataReceiveLatest3D`, `RequestLatest3D`, `DataGetNext3D`

## 9.7.28 RequestNextCentroid

### Function

Request the next frame of centroid data.

### Prototype

```
int RequestNextCentroid( void )
```

### Parameters

None.

### Description

`RequestNextCentroid` is similar to `RequestLatestCentroid`, except that it requests the next available frame of centroid (raw) data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. Unlike the related routine `DataGetNextCentroid`, this routine does not block and wait for the data to be returned. The application program determines if the centroid data is ready, and then receives it.

The routine `DataIsReady` determines if there is centroid data available to be received. When the data is ready to be received, `DataIsReady` returns a non-zero value (TRUE). The application program should then use the routine `DataReceiveLatestRaw` to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns TRUE, the application program does not know which data type it is receiving.**

---

## See Also

`DataIsReady`, `DataReceiveLatestCentroid`, `RequestLatestCentroid`,  
`DataGetNextCentroid`

## 9.7.29 RequestNextOdaRaw

### Function

Requests the next frame of raw data from the specified ODAU device.

### Prototype

```
int RequestNextOdaRaw( int nOdaId )
```

### Parameters

`nOdaId` specifies the ODAU device from which to request the frame of raw data.

### Description

`RequestNextOdaRaw` is similar to `RequestLatestOdaRaw`, except that it requests that the specified ODAU return the next available frame of raw data to the application program instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the ODAU frame rate. Unlike the related routine `DataGetLatestOdaRaw`, this routine does not block and wait for the data to be returned. The application program determines if the raw data is ready, and then receives it.

The routine `DataIsReady` determines if there is data available to be received. When data is ready to be received, `DataIsReady` returns a non-zero value (TRUE). The application program should then use the routine `DataReceiveLatestOdaRaw` to receive the data from the ODAU. Never request data without a followup call to receive the data.



---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns TRUE, the application program does not know which data type it is receiving.**

---

## See Also

`DataIsReady`, `DataReceiveLatestOdaRaw`, `RequestLatestOdaRaw`,  
`DataGetLatestOdaRaw`

### 9.7.30 RequestNextRaw

#### Function

Request the next frame of full raw data.

#### Prototype

```
int RequestNextRaw( void )
```

#### Parameters

None.

#### Description

`RequestNextRaw` is similar to `RequestLatestRaw`, except that it requests the next available frame of full raw data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. Unlike the related routine `DataGetNextRaw`, this routine does not block and wait for the data to be returned. The application program first determines if the raw data is ready, then receives it.

The routine `DataIsReady` determines if there is full raw data available to be received. When the data is ready to be received, `DataIsReady` returns a non-zero value (TRUE).

The application program should then use the routine `DataReceiveLatestRaw` to receive the data from the Optotrak System. Never request data without a followup call to receive the data.




---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine `DataIsReady` returns `TRUE`, the application program does not know which data type it is receiving.**

---

## See Also

`DataIsReady`, `DataReceiveLatestRaw`, `RequestLatestRaw`, `DataGetNextRaw`

### 9.7.31 RequestNextTransforms

#### Function

Requests the next frame of rigid body transformation data.

#### Prototype

```
int RequestNextTransforms( void )
```

#### Parameters

None.

#### Description

`RequestNextTransforms` is similar to `RequestLatestTransforms`, except that it requests the next available frame of rigid body 6D transformation data from the Optotrak System instead of the latest frame. This prevents multiple instances of the same frame from being returned by the system if the data is being requested more frequently than the frame rate. Unlike the related routines `DataGetLatestTransforms` and `DataGetLatestTransforms2`, this routine does not block and wait for the data to be returned. The application program determines if the raw data is ready, then receives it.

The routine `DataIsReady` determines if there are 6D data available to be received. When data is ready to be received, `DataIsReady` returns a non-zero value (`TRUE`). The application program should then use either `DataReceiveLatestTransforms` or

DataReceiveLatestTransforms2 to receive the data from the Optotrak System. Never request data without a followup call to receive the data.



---

**Do not send out a new request for data until the data from an earlier request has been received. This is particularly important if the program is requesting different types of data. For example, if an application requests two types of data, one immediately after the other, and the routine DataIsReady returns TRUE, the application program does not know which data type it is receiving.**

---

You may use this routine if conversions are performed on-host by setting bConvertOnHost and bRigidOnHost to TRUE in the Optotrak System parameter file or with the routine OptotrakSetProcessingFlags. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use this routine.

Before any rigid body transformations can be determined and returned, the application program must load the required rigid body parameters using the routines RigidBodyAdd or RigidBodyAddFromFile.

### See Also

DataIsReady, DataReceiveLatestTransforms, DataReceiveLatestTransforms2, DataGetLatestTransforms, DataGetLatestTransforms2

## 9.8 Buffered Data Retrieval Routines

### 9.8.1 DataBufferAbortSpooling

#### Function

Nullifies all bindings between devices and spool destinations.

#### Prototype

```
int DataBufferAbortSpooling( void )
```

#### Parameters

None.

## Description

DataBufferAbortSpooling nullifies any relationship between a device and a spool destination made using the routines DataBufferInitializeFile or DataBufferInitializeMem.



Warning!

---

**Do not confuse this routine with the routine DataBufferStop, which terminates the current spool operation. DataBufferAbortSpooling only affects the relationship between devices and spool destinations.**

---

Invoke this routine to associate a different spool destination with a device's data spooling. Once DataBufferAbortSpooling is called, all data spooling relationships must be reset using either DataBufferInitializeFile or DataBufferInitializeMem. Do not invoke DataBufferAbortSpooling if the data spooling procedure has already started.

## See Also

DataBufferInitializeFile, DataBufferInitializeMem

## 9.8.2 DataBufferInitializeFile

### Function

Initializes a Northern Digital Floating Point format file as a destination for buffered data.

### Prototype

```
int DataBufferInitializeFile( unsigned int uDataId, char *pszFileName )
```

### Parameters

**uDataId** identifies the device from which buffered data is to be spooled.

Values:

**OPTOTRAK** specifies the Optotrak System Control Unit's data buffer.

**ODAU1** specifies the first ODAU's data buffer.

**ODAU2** specifies the second ODAU's data buffer.

**ODAU3** specifies the third ODAU's data buffer.

ODAU4 specifies the fourth ODAU's data buffer.

**pszFileName** points to a null terminated string that specifies the name of the file to which the data is spooled.

## Description

DataBufferInitializeFile creates a Northern Digital Floating Point format file to spool data to from the specified device's data buffer. If the file already exists, it is replaced by the new one, even if spooling is aborted. Once the spooling process has begun, data is retrieved from the specified device's data buffer and written to the initialized file. When the spooling process is complete, the file's header is updated with the appropriate information, and the file is closed.

See Sample Program 3 on the API CD for an example of code that uses this routine.

## See Also

DataBufferInitializeMem, DataBufferAbortSpooling

### 9.8.3 DataBufferInitializeMem

#### Function

Initializes a block of memory as a destination for buffered data.

#### Prototype

```
int DataBufferInitializeMem( unsigned int    uDataId,  
                           SpoolPtrType   pMemory )
```

#### Parameters

**uDataId** identifies the device from which buffered data is to be spooled.

Values:

**OPTOTRAK** specifies the Optotrak System Control Unit's data buffer.

**ODAU1** specifies the first ODAU's data buffer.

**ODAU2** specifies the second ODAU's data buffer.

**ODAU3** specifies the third ODAU's data buffer.

**ODAU4** specifies the fourth ODAU's data buffer.



**pMemory** specifies the memory area to which the buffered data is spooled.

## Description

`DataBufferInitializeMem` does the required initialization to spool data from the specified device's data buffer to an area of memory reserved by the application program. Once the spooling process has begun, data is retrieved from the specified device's data buffer and written to the appropriate memory location.



**Warning!**

---

**You must ensure the specified memory block is large enough to store the entire contents of the data collection. If the memory block reserved by the application program is too small, data is written to an invalid memory location. This can cause unpredictable results for the application program. This may include a system crash.**

---

To determine the required amount of memory the application program must know how much data will be spooled.

Memory = [collection time x frame frequency] x (size of 1 frame) bytes

To determine the size of one frame, see [“Real-time Data Types” on page 231](#).

See Sample Program 4 on the API CD for an example of code that uses this routine.

## See Also

`DataBufferInitializeFile`, `DataBufferAbortSpooling`

### 9.8.4 `DataBufferSpoolData`

#### Function

Spools an entire data buffer.

#### Prototype

```
int DataBufferSpoolData( unsigned int *puSpoolStatus )
```

#### Parameters

**puSpoolStatus** indicates if there was a data buffer error during the data spooling.

## **Description**

DataBufferSpoolData does everything required to spool buffered data to the initialized spool destinations. This routine starts all devices spooling buffered data at the current frame. It then receives and writes the spooled data to the appropriate spool destinations, either memory or file.

If a spool destination was not initialized for the device's data, then the data spooled from this device is lost. Once all the data has been spooled to the application program, this routine fills in the puSpoolStatus parameter and does the required clean-up (e.g. close data files, etc.).

The status returned via the puSpoolStatus parameter indicates if a data buffering error occurred on any of the devices spooling data. If the returned spool status is non-zero, then an error occurred during the spooling period. The lower byte of this parameter contains a code indicating the error, and the upper byte contains the ID of the device on which the error occurred, i.e. the Optotrak System or one of the ODAU devices

See Sample Program 3 on the API CD for an example of code that uses this routine.

## **See Also**

DataBufferInitializeFile, DataBufferInitializeMem

### **9.8.5 DataBufferStart**

#### **Function**

Starts all data buffers spooling data back to the host computer.

#### **Prototype**

```
int DataBufferStart( void )
```

#### **Parameters**

None.

#### **Description**

DataBufferStart starts spooling all initialized data buffers at the current frame. For every device that you want to spool data from, either a file or memory block must have been initialized. Use this routine with DataBufferWriteData. By using these routines together, an application program can spool data, without blocking for the

entire length of the collection time. Once `DataBufferStart` has been invoked, the routine `DataBufferWriteData` must be called repeatedly until it signals that data spooling has completed.

An application program can stop spooling buffered data before the collection time has elapsed by invoking the routine `DataBufferStop`. This routine stops all devices currently spooling buffered data at the current frame. However, it is still necessary to stay in the data writing loop until the `DataBufferWriteData` routine indicates that the data spooling has completed.

See Sample Program 5 on the API CD for an example of code that uses this routine.

## See Also

`DataBufferStop`, `DataBufferWriteData`, `DataBufferSpoolData`

## 9.8.6 `DataBufferStop`

### Function

Stops all devices spooling buffered data back to the host computer.

### Prototype

```
int DataBufferStop( void )
```

### Parameters

None.

### Description

`DataBufferStop` stops spooling data from the data buffers from all devices to the application program at the current frame. This routine allows an application program to terminate data buffer spooling before the entire collection time has elapsed. However, it is still necessary to stay in the data writing loop until the `DataBufferWriteData` routine indicates that the data spooling has completed.

See Sample Program 7 on the API CD for an example of code that uses this routine.

## See Also

`DataBufferStart`, `DataBufferWriteData`

## 9.8.7 DataBufferWriteData

### Function

Receives data buffer data and writes them to the appropriate destination.

### Prototype

```
int DataBufferWriteData( unsigned int   *puRealtimeData,  
                        unsigned int   *puSpoolComplete,  
                        unsigned int   *puSpoolStatus,  
                        unsigned long  *pulFramesBuffered )
```

### Parameters

**puRealtimeData** indicates if there is a frame of real-time data to be received.

**puSpoolComplete** indicates if the data spooling is complete.

**puSpoolStatus** indicates if there was a data buffer error during the data spooling.

**pulFramesBuffered** indicates the number of frames of data received since the latest call to `DataBufferStart`, including those received by the current call.

### Description

`DataBufferWriteData` receives and writes spooled buffer data to the appropriate, initialized, spool destination. The routine first checks if there is any real-time data to be received that was requested using one of the “Request” prefixed routines. If there is data to be received, the routine sets the `puRealtimeData` variable to a non-zero value (TRUE) and then returns immediately. At this point the application program must receive the real-time data frame using the corresponding “Receive” prefixed routine.

If there is no real-time data to receive, but there is spooled buffer data to receive, then the `DataBufferWriteData` routine receives the data and writes it to the appropriate spool destination, or discards the data if no spool destination was initialized for this device. Once all the data for the collection time has been received, this routine sets the `puSpoolComplete` variable to a non-zero value (TRUE), fills in the `puSpoolStatus` variable, and does the required clean up (e.g. closes data files, etc.).

The status returned via the `puSpoolStatus` variable indicates if a data buffering error occurred on any of the devices spooling data. If the returned spool status is non-zero, then an error occurred during the spooling period. The lower byte of this variable contains a code indicating the error and the upper byte contains the ID of the device on which the error occurred, e.g. the Optotrak System or ODAU1.

The `pulFramesBuffered` value can be used by the application program to determine the degree to which the data spooling has been completed after each call to the routine. If `pulFramesBuffered` is set to `NULL` when the function is called, then the variable is not updated with the number of frames buffered.

See Sample Program 5 on the API CD for an example of code that uses this routine.

### **See Also**

`DataBufferStart`, `DataBufferStop`

## 9.9 Rigid Body Specific Routines

You may use these routines if conversions are performed on-host by setting `bConvertOnHost` and `bRigidOnHost` to `TRUE` in the Optotrak System parameter file or with the routine `OptotrakSetProcessingFlags`. If you wish to perform the conversion internally, you must have an Optotrak Certus System or an Optotrak 3020 System along with the Optotrak Real-time Rigid Body Option to use these routines.

See [“Real-time Rigid Body Programmer’s Guide” on page 51](#) and [“Flags and Settings Associated with Rigid Bodies” on page 295](#) for more information on rigid bodies.

### 9.9.1 RigidBodyAdd

#### Function

Adds a rigid body to the list of rigid bodies that the Optotrak System tracks.

#### Prototype

```
int RigidBodyAdd( int    nRigidBodyId,  
                int    nStartMarker,  
                int    nNumMarkers,  
                float  *pRigidCoordinates,  
                float  *pNormalCoordinates,  
                int    nFlags )
```

#### Parameters

**nRigidBodyId** specifies a unique identifier to be associated with the rigid body.

**nStartMarker** specifies the marker number of the first IRED marker in the rigid body.

**nNumMarkers** specifies the number of IRED markers in the rigid body.

**pRigidCoordinates** specifies an array of 3D positions that define the locations of the IRED marker positions of the rigid body in its home position.

**pNormalCoordinates** specifies an array of 3D positions, each of which defines the normal vectors for an IRED marker of the rigid body in its home position. If `NULL` is specified for this parameter, then no normal coordinates are added for the specified rigid body.

**nFlags** contains the bit flags that configure the manner in which the Optotrak System processes the rigid body. For a complete description of these flags, see [“RigidBodyChangeSettings” on page 211](#).

## Description

RigidBodyAdd adds a rigid body to the Optotrak System rigid body tracking list according to the specified parameters. Once the rigid body has been added, the Optotrak System determines transformations for this rigid body whenever it is requested.

The positions of the rigid body markers are passed to this routine via an array of 3D (x, y, z) positions. The array must be organized as follows:

**Table 9-14: Position Array for Rigid Body Markers**

| Marker     | 3D Position |       |       |
|------------|-------------|-------|-------|
| Marker 1   | $X_1$       | $Y_1$ | $Z_1$ |
| Marker 2   | $X_2$       | $Y_2$ | $Z_2$ |
| ⋮          | ⋮           | ⋮     | ⋮     |
| Marker $M$ | $X_M$       | $Y_M$ | $Z_M$ |

The RigidBodyAdd routine assumes the X, Y, Z coordinates are defined as C-type floats. The X, Y, Z variables for the normal vectors are in the same form as the rigid body marker positions, but they specify a 3D vector instead of a 3D position.

Normal vectors are used to exclude a marker's data from the rigid body transformation calculations whenever the marker is viewed at very oblique angles.

---

**Note** After a call to the routine RigidBodyAdd, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may be generated.

---

See Sample 11 on the API CD for an example of code that uses this routine.

## See Also

RigidBodyAddFromFile, RigidBodyChangeSettings, DataGetLatestTransforms, DataGetLatestTransforms2

## 9.9.2 RigidBodyAddFromDeviceHandle

### Function

Adds the rigid body defined in a device to the list of rigid bodies tracked by the Optotrak System.

### Prototype

```
int RigidBodyAddFromDeviceHandle( int  nDeviceHandle,  
                                int  nFlags )
```

### Parameters

**nDeviceHandle** identifies the device.

**nFlags** contains the flags that configure the way in which the Optotrak System processes the rigid body. For a complete description of these flags, see [“RigidBodyChangeSettings” on page 211](#).

### Description

If the specified device contains rigid body information, `RigidBodyAddFromDeviceHandle` informs the System Control Unit that the device should be tracked as a rigid body.

## 9.9.3 RigidBodyAddFromFile

### Function

Uses the specified rigid body file to add a rigid body to the list of rigid bodies that the Optotrak System tracks.

### Prototype

```
int RigidBodyAddFromFile( int  nRigidBodyId,  
                        int  nStartMarker,  
                        char *pszRigFile,  
                        int  nFlags )
```

### Parameters

**nRigidBodyId** specifies a unique identifier to be associated with the rigid body.



**nStartMarker** specifies the marker number of the first IRED marker in the rigid body.

**pszRigFile** specifies the rigid body file that defines the rigid body being added. When specifying the rigid body file to this routine, the application program must not use the file extension. For example, the application program should specify “frame” instead of “frame.rig.”

**nFlags** contains the bit flags that configure the way in which the Optotrak System processes the rigid body. For a complete description of these flags, see “RigidBodyChangeSettings” on page 211.

## Description

RigidBodyAddFromFile adds a rigid body to the Optotrak System rigid body tracking list according to the specified parameters. This routine differs from RigidBodyAdd in that the rigid body coordinates are read from a rigid body file. Once the rigid body has been added, the Optotrak System determines the rigid body’s transformations whenever it is requested.

This routine searches for the specified rigid body file in a number of directories. The first directory searched is the current directory. If the file is not found in the current directory, the routine makes use of the ND\_USER\_DIR and ND\_DIR environment variables to search for the file. The default location for a rigid body file is in the rigid subdirectory under the directory where the Northern Digital software was installed, e.g. c:\ndigital.

As an example, suppose ND\_USER\_DIR is set to c:\ndigital and ND\_DIR is set to d:\ndigital. RigidBodyAddFromFile would first look for the file in the current directory, then in the directory c:\ndigital\rigid, and finally in the directory d:\ndigital\rigid.

---

**Note** After a call to the routine RigidBodyAddFromFile, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may be generated.

---

See Sample Program 9 on the API CD for an example of code that uses this routine.

## See Also

RigidBodyAdd, RigidBodyChangeSettings, DataGetLatestTransforms

## 9.9.4 RigidBodyChangeFOR

### Function

Changes the default coordinate system for the rigid body transformation calculations.

### Prototype

```
int RigidBodyChangeFOR( int nRigidBodyId, int nRotationMethod )
```

### Parameters

**nRigidBodyId** specifies the ID of the rigid body that determines the coordinate system. If -1 is specified as the rigid body ID, then the coordinate system is changed to the default coordinate system.

**nRotationMethod** specifies the method used to rotate the coordinate system. The options are:

**OPTOTRAK\_STATIC\_RIGID\_FLAG** indicates that the Optotrak System is to use the current position and orientation of the specified rigid body to define the static coordinate system.

**OPTOTRAK\_CONSTANT\_RIGID\_FLAG** indicates that the Optotrak System is to define the coordinate system based on the specified rigid body's transformation each time the rigid body transformations are determined. This is equivalent to measuring one moving rigid body with respect to another moving rigid body.

### Description

RigidBodyChangeFOR changes the coordinate system used to determine the rigid body transformations. The coordinate system can be changed so that the determined rigid body transformations are expressed in a coordinate system based on the position and orientation of one of the rigid bodies. A parameter to the routine specifies which rigid body to use. If the specified rigid body has an ID of -1, the default coordinate system is used.

There are two methods of determining the coordinate system in which the other rigid body transformations are expressed. If **OPTOTRAK\_STATIC\_RIGID\_FLAG** is specified, the coordinate system is based on the current position and orientation of the specified rigid body. If **OPTOTRAK\_CONSTANT\_RIGID\_FLAG** is specified, the coordinate system is computed each time rigid body transformations are requested, based on the position and orientation of the specified rigid body at that

time. If the coordinate system is being constantly computed, the transformation of the specified rigid body is expressed in the default coordinate system. This provides the application program with the transformation needed to determine the new coordinate system.

See Sample Program 11 on the API CD for an example of code that uses this routine.

## See Also

None.

## 9.9.5 RigidBodyChangeSettings

### Function

Changes parameters affecting the calculations of the specified rigid body.

### Prototype

```
int RigidBodyChangeSettings( int    nRigidBodyId,
                           int    nMinMarkers,
                           int    nMaxMarkerAngle,
                           float  fMax3dError,
                           float  fMaxSensorError,
                           float  fMax3dRmsError,
                           float  fMaxSensorRmsError,
                           int    nFlags )
```

### Parameters

**nRigidBodyId** specifies the ID of the rigid body that is being changed.

**nMinMarkers** specifies the minimum number of markers that must be viewed to determine a rigid body transformation.

**nMaxMarkerAngle** specifies the maximum angle, in degrees, that a marker can rotate from a Position Sensor, and still have that data included in the determination of the rigid body transformation. To use this parameter, normal vectors must also be loaded for the rigid body.

**fMax3dError** specifies the maximum calculated 3D error, in millimeters, allowed for a single marker. If a marker has an error greater than this value, then the rigid body transformation is recalculated with that marker's data ignored.

**fMaxSensorError** specifies the maximum calculated raw sensor error allowed for a each marker/sensor pair. If a marker/sensor pair has an error greater than this value,

then the rigid body transformation is recalculated with that marker/sensor pair's data ignored.

**fMax3dRmsError** specifies the maximum allowable 3D RMS error, in millimeters. If the calculated 3D RMS error is greater than **fMax3dRmsError**, then the transformation is deemed to be undetermined.

**fMaxSensorRmsError** specifies the maximum allowable raw sensor RMS error. If the calculated raw sensor RMS error is greater than the **fMaxSensorRmsError**, then the transformation is deemed to be undetermined.

**nFlags** contains the bit flags that configure the way the Optotrak System processes the rigid body. The following list details the flags and their significance:

**OPTOTRAK\_NO\_RIGID\_CALCS\_FLAG** to stop the Optotrak System from determining rigid body transformations, but keep the rigid body in the tracking list.

**OPTOTRAK\_DO\_RIGID\_CALCS\_FLAG** to restart the Optotrak System determining rigid body transformations if they are turned off for this rigid body.

**OPTOTRAK\_QUATERN\_RIGID\_FLAG** to force the Optotrak System to use the quaternion algorithm when determining the rigid body transformation.

The quaternion algorithm is a closed-form solution of absolute position and orientation of a rigid body based on proprietary enhancements made to the original algorithm developed by Berthold K.P. Horn.

**OPTOTRAK\_ITERATIVE\_RIGID\_FLAG** to force the Optotrak System to use the Iterative Euler Angle algorithm when determining the rigid body transformation.

The Iterative Euler Angle algorithm calculates an initial estimate using the quaternion algorithm and further refines the result using sensor information. In theory, this algorithm generates more accurate results than the quaternion algorithm. In practice, the improvement is negligible and requires longer time to complete.

**OPTOTRAK\_RETURN\_QUATERN\_FLAG** to force the Optotrak System to return the rigid body transformation in quaternion format.

**OPTOTRAK\_RETURN\_MATRIX\_FLAG** to force the Optotrak System to return the rigid body transformation in Rotation Matrix format.

**OPTOTRAK\_RETURN\_EULER\_FLAG** to force the Optotrak System to return the rigid body transformation in Euler Angle format in radians.

## Description

RigidBodyChangeSettings uses the specified parameters to control the way in which the Optotrak System processes information for the specified rigid body. This routine can be used to suppress rigid body calculations for a rigid body, change the algorithm used to determine the rigid body transformation, and select the format for the rigid body transformation.

Some of the settings (flags) can be combined while others may not. If both the OPTOTRAK\_NO\_RIGID\_CALCCS\_FLAG and the OPTOTRAK\_DO\_RIGID\_CALCCS\_FLAGS are specified together, the first flag is ignored and the calculations are performed.

When specifying the transformation determination algorithm, the flags may be used individually or together. If just the OPTOTRAK\_QUATERN\_RIGID\_FLAG is specified, then just the Quaternion algorithm is used. If just the OPTOTRAK\_ITERATIVE\_RIGID\_FLAG is specified, then just the Iterative Euler Angle algorithm is used. If both flags are specified, then the Quaternion algorithm is used, and its result is then used as an initial estimate for the Iterative Euler Angle algorithm. Finally, the flags for the transformation return format must be used exclusively.

The routines RigidBodyAdd and RigidBodyAddFromFile each include a flags parameter that uses the above flag definitions. If no flags are specified when adding a rigid body, a set of default values is used. The default for the minimum markers is 3, the default for the maximum 3D error is 0.25 mm, and the default combination for the flags is OPTOTRAK\_DO\_RIGID\_CALCCS\_FLAG, OPTOTRAK\_QUATERN\_RIGID\_FLAG, OPTOTRAK\_ITERATIVE\_RIGID\_FLAG, and OPTOTRAK\_RETURN\_EULER\_FLAG.

See Sample Program 10 on the API CD for an example of code that uses this routine.

## See Also

RigidBodyAdd, RigidBodyAddFromFile

## 9.9.6 RigidBodyDelete

### Function

Removes the specified rigid body from the Optotrak System rigid body tracking list.

### Prototype

```
int RigidBodyDelete( int nRigidBodyId )
```

## Parameters

**nRigidBodyId** specifies the ID of the rigid body that is to be removed.

## Description

RigidBodyDelete completely removes the rigid body associated with the specified ID from the Optotrak System rigid body tracking list. After a rigid body is removed, no more transformations are determined for that body.

---

**Note** Do not use this routine to remove the rigid body that is being used to define a constantly changing coordinate system; the rigid body must first be removed from the coordinate system calculation before it can be deleted.

---

## See Also

RigidBodyAdd, RigidBodyAddFromFile, RigidBodyChangeFOR

## 9.10 Rigid Body Related Routines

The following routines are commonly used for manipulating rigid body data structures and transformation data.

### 9.10.1 CombineXfrms

#### Function

Combines two Euler angle format transformations specified in roll-pitch-yaw geometry.

#### Prototype

```
void CombineXfrms( transformation *pdtXfrm1,  
                  transformation *pdtXfrm2,  
                  transformation *pdtNewXfrm )
```

#### Parameters

**pdtXfrm1** points to the first transformation.

**pdtXfrm2** points to the second transformation.

**pdtNewXfrm** points to the combined transformation.

## Description

CombineXfrms combines two Euler angle format transformations specified by pdtXfrm1 and pdtXfrm2 into a combined transformation, pdtNewXfrm. The routine is an implementation of the following relation:

$$T_{13} = T_{23} \cdot T_{12}$$

$T_{12}$  is the first transformation (pdtXfrm1) and is defined as the transformation from frame 1 to frame 2

$T_{23}$  is the second transformation, (pdtXfrm2) and is defined as the transformation from frame 2 to frame 3

$T_{13}$  is the transformation from frame 1 to frame 3.

## See Also

DetermineR

### 9.10.2 CvtQuatToRotationMatrix

#### Function

Determines the rotation matrix corresponding to the specified quaternion.

#### Prototype

```
void CvtQuatToRotationMatrix( QuatRotationStruct   *pdtQuatRot,
                             RotationMatrixType   pdtRotMatrix)
```

#### Parameters

**pdtQuatRot** points to the quaternion format rotation being converted.

**pdtRotMatrix** points to the  $3 \times 3$  rotation matrix.

#### Description

CvtQuatToRotationMatrix converts the input quaternion format rotation specified by pdtQuatRot to its corresponding rotation matrix, which is returned in pdtRotMatrix. This routine complements CvtRotationMatrixToQuat, which converts an input rotation matrix to its corresponding quaternion format.

## See Also

`CvtRotationMatrixToQuat`

### 9.10.3 `CvtRotationMatrixToQuat`

#### Function

Determines the quaternion format rotation corresponding to the specified rotation matrix.

#### Prototype

```
void CvtRotationMatrixToQuat( QuatRotationStruct *pdtQuatRot,  
                             RotationMatrixType  pdtRotMatrix)
```

#### Parameters

`pdtQuatRot` points to the quaternion defining the rotation.

`pdtRotMatrix` points to the  $3 \times 3$  rotation matrix being converted.

#### Description

`CvtRotationMatrixToQuat` converts the input rotation matrix specified by `pdtRotMatrix` to its corresponding quaternion format, which is returned in `pdtQuatRot`. This routine complements `CvtQuatToRotationMatrix`, which converts an input quaternion to its corresponding rotation matrix.

## See Also

`CvtQuatToRotationMatrix`

### 9.10.4 `DetermineEuler`

#### Function

Determines the Euler angles from the specified  $3 \times 3$  rotation matrix. Values are returned in radians.

#### Prototype

```
void DetermineEuler( RotationMatrixType  pdtRotMatrix,
```



```
rotation          *pdtEulerRot )
```

## Parameters

**pdtRotMatrix** points to the  $3 \times 3$  rotation matrix being converted.

**pdtEulerRot** points to the resultant Euler angle format rotation.

## Description

DetermineEuler converts the input rotation matrix specified by **pdtRotMatrix** to its corresponding Euler angle format, which is returned in **pdtEulerRot**. This routine complements DetermineR, which converts input Euler angles to their corresponding rotation matrix.

## See Also

DetermineR

### 9.10.5 DetermineR

## Function

Determines the  $3 \times 3$  rotation matrix from the specified Euler angle format rotation.

## Prototype

```
void DetermineR( rotation          *pdtEulerRot,
                 RotationMatrixType pdtRotMatrix )
```

## Parameters

**pdtEulerRot** points to the Euler angle format rotation being converted.

**pdtRotMatrix** points to the resultant  $3 \times 3$  rotation matrix.

## Description

DetermineR converts the input Euler angle format rotation specified by **pdtEulerRot** to its corresponding rotation matrix, which is returned in **pdtRotMatrix**. The routine is an implementation of the following relation:

$$R(r, p, y) = R(r) \cdot R(p) \cdot R(y)$$

**R** is a  $3 \times 3$  rotation matrix

*r* is the roll, or the rotation about the *z*-axis

*p* is the pitch, or the rotation about the *y*-axis

*y* is the yaw, or the rotation about the *x*-axis

This routine complements DetermineEuler, which converts the input rotation matrix to the corresponding Euler angles.

### **See Also**

DetermineEuler

## **9.10.6 InverseXfrm**

### **Function**

Determines the inverse to the specified Euler angle format transformation.

### **Prototype**

```
void InverseXfrm( transformation *pdtXfrm,  
                transformation *pdtInverseXfrm )
```

### **Parameters**

**pdtXfrm** points to the input transformation.

**pdtInverseXfrm** points to the resultant inverse transformation.

### **Description**

InverseXfrm calculates the inverse of the Euler angle format transformation specified by pdtXfrm and outputs it as pdtInverseXfrm. This routine implements the relation:

$$T_{\text{inverse}} = T^{-1}$$

where **T** = the Euler angle format transformation.

### **See Also**

None.

## 9.10.7 TransformPoint

### Function

This routine calculates the 3-D position of a given point after it has been transformed by the specified rotation and translation.

### Prototype

```
void TransformPoint( RotationMatrixType   pdtRotationMatrix,
                    Position3d          *pdtTranslation,
                    Position3d          *pdtOriginalPositionPtr,
                    Position3d          *pdtTransformedPositionPtr )
```

### Parameters

**pdtRotationMatrix** points to the 3 x 3 rotation matrix defining the rotational part of the transformation.

**pdtTranslation** points to the structure defining the translational part of the transformation.

**pdtOriginalPositionPtr** points to the structure defining the original position of the point being transformed.

**pdtTransformedPositionPtr** points to the structure containing the transformed position of the point. The original contents of this structure will be overwritten in a safe manner so that it can be the same as pdtOriginalPositionPtr if the original position does not need to be retained.

### Description

TransformPoint determines the 3-D position of a given point after being transformed by the specified rotation and translation according to the following relation:

$$p_1 = R p_0 + T$$

where

$p_1$ =the position after transformation

$p_0$ =the position before transformation

$R$ =the matrix rotation component of the transformation

$T$ =the translation component of the transformation

This routine does not require separate memory to be allocated for the output position if the input position need not be retained – calling TransformPoint with the

values for `pdtOriginalPositionPtr` and `pdtTransformedPositionPtr` will overwrite the original position with the transformed position.

### See Also

None.

## 9.11 File Processing Routines

### 9.11.1 FileClose

#### Function

Closes the specified file.

#### Prototype

```
int FileClose( unsigned int uFileId )
```

#### Parameters

`uFileId` identifies the file that is to be closed.

#### Description

`FileClose` closes the file corresponding to a specified file ID that was opened with `FileOpen`. This allows the application program to open another file with that ID. Once the file is closed, the associated file ID remains invalid until a new file is opened with that ID.

Files that were opened with `FileOpenAll` should be closed with `FileCloseAll`, not `FileClose`.

See Sample 15 on the API CD for an example of code that uses this routine.

### See Also

`FileOpen`, `FileOpenAll`, `FileCloseAll`

## 9.11.2 FileCloseAll

### Function

Closes a specified file that was opened with FileOpenAll.

### Prototype

```
int FileCloseAll( unsigned int uFileId )
```

### Parameters

**uFileId** identifies the file that is to be closed.

### Description

Use FileCloseAll to close files opened with FileOpenAll. This routine correctly handles files having all subitem types. Closing the file allows the application program to open another file with that ID. Once the file is closed, the associated file ID remains invalid until a new file is opened with that ID.

Files that were opened with FileOpen should not be closed with FileCloseAll, but should be closed with FileClose instead.

### See Also

FileOpen, FileClose, FileCloseAll

## 9.11.3 FileConvert

### Function

Converts raw data from the specified input file to the appropriate type and writes the converted data to the specified output file. Within this routine, raw data includes full raw data, centroid data and ODAU raw data. The camera file that has been loaded using the OptotrakLoadCameraParameters will be used for the data conversion.

### Prototype

```
int FileConvert( char          *pszInputFilename,  
                char          *pszOutputFilename,  
                unsigned int  uFileType )
```

## Parameters

**pszInputFilename** is a null terminated string specifying the name of the input file containing the raw data to be converted.

**pszOutputFilename** is a null terminated string specifying the name of the output file to which the converted data is to be written.

**uFileType** specifies the type of raw data that is being converted. This parameter may have one of the following values:

Values:

**OPTOTRAK\_RAW** indicates that the input file contains raw data, and should be converted to 3D data.

**ANALOG\_RAW** indicates that the input file contains ODAU raw data, and should be converted to voltages.

## Description

FileConvert converts the raw data stored in the specified input file to the appropriate type, and writes the converted data to the specified output file. The two types of raw data are raw data, and analog raw data. These data files are converted to 3D data and voltage data, respectively. An analog raw data file is created by spooling data from an ODAU's data buffer.



---

**Do not invoke this routine if the Optotrak System is currently spooling data. Doing so will corrupt the input file.**

---

See Sample 8 on the API CD for an example of code that uses this routine.

## See Also

TransputerInitializeSystem, OptotrakLoadCameraParameters

### 9.11.4 FileOpen

#### Function

Opens a Northern Digital Floating Point format file in the specified mode. This routine only allows access to float subitems.

## Prototype

```
int FileOpen( char          *pszFilename,
              unsigned int  uFileId,
              unsigned int  uFileMode,
              int           *pnItems,
              int           *pnSubItems,
              long int      *pLnFrames,
              float         *pfFrequency,
              char          *pszComments,
              void          **pFileHeader )
```

## Parameters

**pszFilename** is a null terminated string specifying the name of the Northern Digital Floating Point format file that is to be opened.

**uFileId** specifies a unique ID that is to be used to identify the newly opened file.

**uFileMode** specifies the mode in which the file is to be opened. Two modes are possible:

**OPEN\_READ** indicates that the file is to be opened in read-only mode.

**OPEN\_WRITE** indicates that the file is to be opened in read-write mode.

**pnItems** is the number of items in the file.

**pnSubItems** is the number of floating point type subitems in the file.

**pLnFrames** is the number of frames in the file.

**pfFrequency** is the frequency at which the data in the file was collected.

**pszComments** is a null terminated string, of up to 59 characters in length, containing user comments for the file.

**pFileHeader** is set to a pointer to the Northern Digital file header for the specified file.

## Description

`FileOpen` opens the specified Northern Digital Floating Point format file, allowing an application program to read or write frames of data from or to the file. The application program can open up to 16 files, using the file IDs 0 - 15 to identify each of the files. If the file is opened in read mode, then data can only be read from it, but if the file is opened in read-write mode, then data can be written to it as well.

When a file is opened, the contents of the routine parameters, `pnItems`, `pnSubItems`, `plnFrames`, `pfFrequency`, and `pszComments` are either read from the file's header or written to it.

The parameter `pFileHeader` is set to point to the file header for the file. While it is not recommended that an application program change the contents of the file header itself, this option is available if necessary.

`FileOpen` only allows access to float type subitems. Use the related routine, `FileOpenAll` to access all subitem types. Any files opened in an application program should be closed before the program exits. Use `FileClose`, not `FileCloseAll`, to close files opened with `FileOpen`.

See Sample 15 on the API CD for an example of code that uses this routine.

## See Also

`FileClose`, `FileRead`, `FileWrite`, `FileOpenAll`

## 9.11.5 FileOpenAll

### Function

Opens a Northern Digital Floating Point format file in the specified mode, allowing access to all subitem types.

### Prototype

```
int FileOpenAll( char          *pszFilename,
                unsigned int   uFileId,
                unsigned int   uFileMode,
                int             *pnItems,
                int             *pnSubItems,
                int             *pnCharSubItems,
                int             *pnIntSubItems,
                int             *pnDoubleSubItems,
                long int        *plnFrames,
                float           *pfFrequency,
                char            *pszComments,
                void            **pFileHeader )
```

### Parameters

**pszFilename** is a null terminated string that specifies the name of the Northern Digital Floating Point format file to be opened.

**uFileId** specifies a unique ID that identifies the newly opened file.



**uFileMode** specifies the mode in which the file is to be opened. Two modes are possible:

**OPEN\_READ** indicates that the file is to be opened in read-only mode.

**OPEN\_WRITE** indicates that the file is to be opened in read-write mode.

**pnItems** is the number of items in the file.

**pnSubItems** is the number of floating point type subitems in the file.

**pnCharSubItems** is the number of character type subitems in the file.

**pnIntSubItems** is the number of integer type subitems in the file.

**pnDoubleSubItems** is the number of double type subitems in the file.

**plnFrames** is the number of frames in the file.

**pfFrequency** is the frequency at which the data in the file was collected.

**pszComments** is a null terminated string, of up to 59 characters in length containing user comments for the file.

**pFileHeader** is set to a pointer to the Northern Digital file header for the specified file.

## Description

FileOpenAll is similar to FileOpen, but allows access to all of the subitem types (float, char, int and double). This routine opens the specified Northern Digital Floating Point format file, allowing an application program to read or write frames of data from or to the file. Up to 16 files can be opened, using the file IDs 0 - 15 to identify each of the files. If the file is opened in the read mode, then data can only be read from it, but if the file is opened in the read-write mode, then data can be written to it as well.

When a file is opened, the contents of the routine parameters, pnItems, pnSubItems, pnCharSubItems, pnIntSubItems, pnCharSubItems, plnFrames, pfFrequency, and pszComments are either read from the file's header or written to it.

The parameter pFileHeader is set to point to the file header for the file. While it is not recommended that an application program change the contents of the file header itself, this option is available. Any files opened in an application program should be closed before the program exits. Use FileCloseAll, not FileClose, to close files opened with FileOpenAll.

## See Also

FileOpen, FileCloseAll, FileRead, FileWrite

## 9.11.6 FileRead

### Function

Reads the specified frames for the specified file. This routine will only access float subitems.

### Prototype

```
int FileRead( unsigned int   uFileId,
              long int      lnStartFrame,
              unsigned int   uNumberOfFrames,
              void           *pDataDest )
```

### Parameters

**uFileId** specifies the file ID that identifies the file from which the data frames are to be read.

**lnStartFrame** specifies the starting frame of the data to be read. Frame numbers begin indexing at 0, so for files having  $N$  frames, the frames are identified as frames 0 to  $N-1$ .

**uNumberOfFrames** specifies the number of frames of data to be read.

**pDataDest** is the block of memory that the data frames are to be read into.

### Description

FileRead reads data frames from a file previously opened with FileOpen. The application program may specify the frame to start reading data, and the number of data frames to be read. The data is read from the file directly into the block of memory specified by the application program.

FileRead only allows access to float type subitems. To access all subitem types, use the related routine, FileReadAll.



**You must specify a memory block with pDataDest that is large enough to store all the data for the number of frames being read. If the memory block reserved by an application program is too small, data is written to an invalid memory location. This can cause unpredictable results for the application program, including a system crash.**

---

To determine the minimum amount of memory you must reserve, the application program must know the size of the data frames in the file being read.

Size of One Frame = [(Number of items in file) x (Number of subitems) x 4] bytes  
 where 4 is the size of float type subitems.

Amount of Memory = (Number of frames) x (Size of One Frame) bytes

See Sample 15 on the API CD for an example of code that uses this routine.

## See Also

FileOpen, FileReadAll, FileWrite

### 9.11.7 FileReadAll

#### Function

Reads the specified frames for the specified file, allowing access to all subitem types.

#### Prototype

```
int FileReadAll( unsigned int  uFileId,
                long int      lnStartFrame,
                unsigned int  uNumberOfFrames,
                void          *pDataDestFloat,
                void          *pDataDestChar,
                void          *pDataDestInt,
                void          *pDataDestDouble )
```

#### Parameters

**uFileId** specifies the file ID identifying the file from which the data frames are to be read.

**lnStartFrame** specifies the starting frame of the data to be read. Frame numbers begin indexing at 0, so for files having  $N$  frames, the frames are identified as frames 0 to  $N-1$ .

**uNumberOfFrames** specifies the number of frames of data to be read.

**pDataDestFloat** is the block of memory that the float type subitem data is to be read into.

**pDataDestChar** is the block of memory that the character type subitem data is to be read into.

**pDataDestInt** is the block of memory that the integer type subitem data is to be read into.

`pDataDestDouble` is the block of memory that the double type subitem data is to be read into.

## Description

`FileReadAll` is similar to `FileRead`, but recognizes all subitem types (float, char, int and double) rather than just float type subitems. `FileReadAll` reads data frames from a file previously opened with `FileOpenAll`. The application program may specify at which frame to start reading data, and the number of data frames to be read. The data is read from the file directly into the block of memory specified by the application program.

`FileReadAll` currently requires that the last four parameters (the destination data buffers) be non-null values. Even if there are no subitems of a particular type, you must still provide valid pointers to all four destination buffers.



---

**You must specify memory blocks, using the destination buffers, that are large enough to store all the data for the number of frames being read. If the memory blocks reserved by an application program are too small, data is written to invalid memory locations. This can cause unpredictable results for the application program, including a system crash.**

---

To determine the required amount of memory, the application program must know the size of the data frames for each type in the file being read.

Size of One Frame = [(Number of items in file) x (Number of subitems) x 8] bytes

Amount of Memory = (Number of frames) x (Size of One Frame) bytes

where 8 is the size of the a double subitem.

The size of a double subitem is used as it is the largest of all the subitem types; any file size calculated using this value will be large enough for all subitem types.

## See Also

`FileOpenAll`, `FileRead`, `FileWriteAll`

### 9.11.8 FileWrite

#### Function

Writes the specified frames of float type subitem data to the specified file.

## Prototype

```
int FileWrite( unsigned int  uFileId,
              long int      lnStartFrame,
              unsigned int  uNumberOfFrames,
              void          *pDataSrc )
```

## Parameters

**uFileId** specifies the file to which the data frames are to be written.

**lnStartFrame** specifies the starting frame of the data to be written. Frame numbers begin indexing at 0, so  $N$  frames are identified as frames 0 to  $N-1$ .

**uNumberOfFrames** specifies the number of frames of data to be written.

**pDataSrc** is the block of memory containing the data to be written.

## Description

FileWrite writes data frames to a file previously opened using the routine FileOpen. The routine can specify at which frame to start writing data and the number of data frames to be written. The data is written from the block of memory specified by the application program directly to the file.



**Warning!**

---

**When the application program first opened the file in read-write mode, it specified the number of frames that were to be written to the file. The application program must write exactly that number of frames to the file, since that is the number that is written to the file header, which governs subsequent file access. If too few frames are written, then another program reading from the file either fails or reads erroneous data. If extra frames are written to the file, then other programs will not be unable to read them, since only the number of frames specified in the file header are read.**

---

FileWrite writes float type subitems only. Use FileWriteAll to write all of the subitem types.

See Sample 15 on the API CD for an example of code that uses this routine.

## See Also

FileOpen, FileRead, FileWriteAll

## 9.11.9 FileWriteAll

### Function

Writes the specified frames of all subitem type data to the specified file.

### Prototype

```
int FileWriteAll( unsigned int  uFileId,
                 long int      lnStartFrame,
                 unsigned int  uNumberOfFrames,
                 void          *pDataSrcFloat
                 void          *pDataSrcChar
                 void          *pDataSrcInt
                 void          *pDataSrcDouble )
```

### Parameters

**uFileId** specifies the file to which the data frames are to be written.

**lnStartFrame** specifies the starting frame of the data to be written. Frame numbers begin indexing at 0, so  $N$  frames are identified as frames 0 to  $N - 1$ .

**uNumberOfFrames** specifies the number of frames of data to be written.

**pDataSrcFloat** is the block of memory containing the float type subitem data to be written.

**pDataSrcChar** is the block of memory containing the character type subitem data to be written.

**pDataSrcInt** is the block of memory containing the integer type subitem data to be written.

**pDataSrcDouble** is the block of memory containing the double type subitem data to be written.

### Description

FileWriteAll is similar to FileWrite, but can be used for all subitem types (float, char, int and double) instead of just float type subitems. FileWriteAll writes data frames to a file previously opened using the routine FileOpenAll. An application program can specify which frame to start writing data at, as well as the number of data frames to be written. The data is written from the memory blocks specified by the application program directly to the file.




---

**When the application program first opened the file in read-write mode, it specified the number of frames that were to be written to the file. The application program must write exactly that number of frames to the file, since that is the number that is written to the file header, which governs subsequent file access. If too few frames are written, then another program reading from the file either fails or reads erroneous data. If extra frames are written to the file, then other programs will not be unable to read them, since only the number of frames specified in the file header are read.**

---

## See Also

FileOpenAll, FileWrite, FileReadAll

## 9.12 Registration and Alignment Routines

### 9.12.1 nOptotrakAlignSystem

#### Function

Creates a camera parameter file that allows Position Sensors to report positions in a user determined coordinate system.

#### Prototype

```
int nOptotrakAlignSystem( AlignParms dtAlignParms,
                        float*      pfRMSError )
```

#### Parameters

**pfRMSError** points to a float structure where the RMS error for the alignment procedure will be stored.

**dtAlignParms** an AlignParms structure containing the registration parameters. The structure is defined as:

```
typedef struct AlignParametersStruct
{
    char
        szDataFile[_MAX_FNAME],
        szRigidBodyFile[_MAX_FNAME],
        szInputCamFile[_MAX_FNAME],
        szOutputCamFile[_MAX_FNAME],
        szLogFileName[_MAX_FNAME];
    boolean
```

```
        bVerbose;  
    } AlignParms;
```

The AlignParametersStruct contains the following parameters:

**szRawDataFile** Optotrak System data file collected using unaligned camera parameters.

**szRigidBodyFile** file containing the rigid body description for the collected data.

**szInputCamFile** camera parameter file containing unaligned camera parameters.

**szOutputCamFile** camera parameter file that will contain the newly aligned camera parameters.

**szLogFileName** optional filename for the log file. The log file will contain details about the alignment operation.

**bVerbose** indicates if the routine should output processing details to stdout.

### Description

A user may define a coordinate system by using a rigid body's local coordinate system. The nOptotrakAlignSystem routine takes an Optotrak System data file for a specific rigid body obtained using the unaligned Position Sensor(s), along with the appropriate camera parameter file, and determines the camera parameters required to transform data into the rigid body coordinate system.

In order to retrieve data from the Optotrak System in the aligned coordinate system, call the OptotrakLoadCameraParameters routine before any data retrieval or data buffering routine.

---

**Note** The data file must be collected using the specified rigid body.

---

### 9.12.2 nOptotrakCalibrigSystem

#### Function

Generates a set of camera parameters that allow multiple Position Sensors to report positions in a common coordinate system. This routine uses a static set of data to determine the new camera parameters.

#### Prototype

```
int nOptotrakCalibrigSystem( CalibrigParms dtCalibrigParms,
```



```
float*          pfRMSError )
```

## Parameters

**pfRMSError** points to a float structure where the RMS error for the registration procedure will be stored.

**dtCalibrigParms** is a CalibrigParms structure containing the registration parameters. The structure is defined as:

```
typedef struct CalibrigParametersStruct
{
    char
        *pszRawDataFile,
        *pszRigidBodyFile,
        *pszInputCamFile,
        *pszOutputCamFile,
        *pszLogFileName;
    boolean
        bVerbose;
} CalibrigParms;
```

The CalibrigParametersStruct contains the following parameters:

**szRawDataFile** is a raw data file collected using the unregistered camera parameters.

**szRigidBodyFile** is a file containing the rigid body description for the collected data.

**szInputCamFile** is a camera parameter file containing unregistered camera parameters.

**szOutputCamFile** is a camera parameter file that will contain the newly registered camera parameters.

**szLogFileName** is an optional filename for the log file. The log file will contain details about the registration operation.

**bVerbose** indicates if the routine should output processing details to stdout.

## Description

When using more than one Position Sensor simultaneously, the Position Sensors must report the marker positions in a common coordinate system. The nOptotrakCalibrigSystem routine takes a static raw data file for a specific rigid body

obtained from the unregistered Position Sensors, along with the appropriate camera parameter file, and determines a common coordinate system.

In order to retrieve data from the Optotrak System in the common coordinate system, call the `OptotrakLoadCameraParameters` routine before any data retrieval or data buffering routine.

To generate the input data file, collect a raw or full raw data file using the specified rigid body. Place the rigid body at a location where all Position Sensors can see at least 3 markers. Ensure that the rigid body is not moved during data collection.

---

**Note** The data file must be collected using the specified rigid body.

---

### 9.12.3 `nOptotrakRegisterSystem`

#### Function

Generates a set of camera parameters that allows multiple Position Sensors to report marker data in a common coordinate system. This routine uses a dynamic set of data to determine the new camera parameters.

#### Prototype

```
int nOptotrakRegisterSystem( RegisterParms dtRegisterParms,
                             float*       pfRMSError )
```

#### Parameters

`pfRMSError` points to a float structure where the RMS error for the registration procedure will be stored.

`dtRegisterParms` a `RegisterParms` structure containing the registration parameters. The structure is defined as:

```
typedef struct RegisterParamsStruct
{
    char
        zRawDataFile[_MAX_FNAME],
        szRigidBodyFile[_MAX_FNAME],
        szInputCamFile[_MAX_FNAME],
        szOutputCamFile[_MAX_FNAME],
        szLogFileName[_MAX_FNAME];
    float
        fXfrmMaxError,
        fXfrm3dRmsError,
        fSpread1,
```

```

        fSpread2,
        fSpread3;
    int
        nMinNumberOfXfrms,
        nLogFileLevel;
    boolean
        bCheckCalibration,
        bVerbose;
} RegisterParams;

```

The RegisterParamsStruct contains the following parameters:

**szRawDataFile** is a raw data file collected using the unregistered camera parameters.

**szRigidBodyFile** is a file containing the rigid body description for the collected data.

**szInputCamFile** is a camera parameter file containing unregistered camera parameters.

**szOutputCamFile** is a camera parameter file that will contain the newly registered camera parameters.

**szLogFileName** is an optional filename for the log file. The log file will contain details about the registration operation.

**fXfrmMaxError** is the maximum allowable error in the rigid body transformation. When an Optotrak System performs a collection, it locates the visible markers on a rigid body and stores that measurement information in a frame of data. The Optotrak System then performs a best-fit for each visible marker, by comparing the collected information to the stored information in a rigid body file. For each best-fit, there is some residual error. Using the residual error values, the Optotrak System performs an RMS error calculation. If the RMS error calculation returns a value less than that specified by fXfrmMaxError, the collected frame of data is used. If the RMS error calculation exceeds the value specified by fXfrmMaxError, the collected frame of data is rejected. Although, appropriate values for fXfrmMaxError depend on the type of rigid body you are measuring, values ranging from 0.2 - 0.5 mm may be considered typical.

**fXfrm3dRmsError** is the maximum allowable 3D RMS error. When multiple Optotrak Systems are used to perform a collection, each system creates points in space which represent a rigid body. The measurement information for these points in space are stored in a frame of data. The Optotrak Systems then perform a best-fit for each point in space, by comparing each other's collected

information. For each best-fit, there is some residual error. Using the residual error values, the Optotrak Systems perform a 3D RMS error calculation. If the 3D RMS error calculation returns a value less than that specified by `fXfrm3dRmsError`, each system's collected frame of data is used. If the 3D RMS error calculation exceeds the value specified by `fXfrm3dRmsError`, each system's collected frame of data is rejected. Although, appropriate values for `fXfrm3dRmsError` depend on the type of rigid body you are measuring, values ranging from 0.5 - 0.75 mm may be considered typical.

**fSpread1**, **fSpread2**, **fSpread3** is the minimum length for each side of the registration volume. The raw data used for registration must cover a volume that has, at minimum, the size specified by the spread parameters. Set the spread parameters to 0.0 if there are no minimum dimensions required.

**fMinNumberOfXfrms** is the minimum number of valid transformations required for a valid registration.

**nLogFileLevel** indicates the level of detail contained in the log file. The log file level must be one of the following: `REG_NO_LOG_FILE`, `REG_SUMMARY_LOG_FILE`, or `REG_DETAILED_LOG_FILE`.

**bCheckCalibration** indicates if the registration function should only check an existing calibration or generate a new camera parameter file.

**bVerbose** indicates if the routine should output processing details to stdout.

## Description

When using more than one Position Sensor simultaneously, the Position Sensors must report the marker positions in a common coordinate system. The `nOptotrakRegisterSystem` routine takes a dynamic raw data file for a specific rigid body obtained from the unregistered Position Sensors, along with the appropriate camera parameter file, and determines a common coordinate system.

To retrieve data from the Optotrak System in the common coordinate system, call the `OptotrakLoadCameraParameters` routine before any data retrieval or data buffering routine.

To generate the input data file, collect a raw or full-raw data file using the specified rigid body. Slowly move the rigid body in the desired volume. Ensure that the volume covered meets the minimum dimensions specified by the spread parameters.

---

**Note** The data file must be collected using the specified rigid body.

---

## 10 Real-time Data Types

This chapter will help you understand the types of real-time data returned by the devices in the Optotrak System. These data types are returned in response to requests for real-time data from an application program. The types of data your Optotrak System is capable of returning depend on the devices present in the system configuration, and the options purchased with the system.

An Optotrak System primarily returns raw (centroid) and 3D data. Full raw data may also be returned; it contains raw data in addition to information on the amplification, error codes, and signal strength.

An Optotrak System equipped to do rigid body calculations (see the [“Real-time Rigid Body Programmer’s Guide” on page 51](#)) can return data in three formats: Euler Representation, Quaternion Representation, and Rotation Matrix Representation.

An Optotrak System equipped with an ODAU can also return two additional types of data; raw analog and raw analog plus digital. (see [“ODAU Programmer’s Guide” on page 41](#)).

---

**Note** Routines that have *raw* in the routine name (for example, `DataGetNextRaw`) are accessing *full raw* data. Routines that have *centroid* in the routine name (for example, `DataReceiveLatestCentroid`) are accessing *raw (centroid only)* data.

---

### 10.1 “Missing” Marker Constants

There are two constants that are commonly used by NDI routines to indicate that a value is “missing”: `BAD_FLOAT` and `MAX_NEGATIVE`.

If an element of data is “missing” (i.e., it could not be determined), it is set to `BAD_FLOAT`. For example, if a sensor could not produce a valid centroid for a certain marker, then the sensor would assign `BAD_FLOAT` to that marker’s centroid value.

An application program can determine if an element of data is valid by comparing it to the `MAX_NEGATIVE` constant. Since the `BAD_FLOAT` value is less than the `MAX_NEGATIVE` value, the data element will be less than `MAX_NEGATIVE` if it is “missing.”

These constants are defined in the include file, `ndtypes.h` as follows:

```
#define BAD_FLOAT(float)      -3.697314E28
#define MAX_NEGATIVE(float) -3.0E28
```

## 10.2 Optotrak Raw and Full Raw Data

An Optotrak System primarily returns raw (centroid) and 3D data. Full raw data may also be returned: it contains raw data in addition to information on the amplification, error codes and signal strength.

### 10.2.1 Organization of Optotrak Raw and Full Raw Data

The Optotrak System organizes the full raw or raw (centroid) data by grouping according to the marker. For example, marker one's data is followed by marker two's data, and so on. The raw data for each marker is divided into two parts. First, there is a centroid for each sensor in the Optotrak System. Second, there is sensor status information for each sensor. This includes the signal strength high (peak), signal amplification (DRC), any error codes (marker status) and the signal strength low (peak nibble).

Table 10-1 on page 238 represents the data organization of one frame of data for Position Sensors reporting full raw data for M markers. The organization of one frame of centroid data does not include the sensor status. This information is presented for a series of frames in "Sample Formats" on page 68.

**Table 10-1: Organization of One Frame of Full Raw Sensor Data**

| <b>Marker Number</b> | <b>Centroid Data for Each Sensor</b> |       |     |       | <b>Sensor Status for Each Sensor</b> |                     |     |                     |
|----------------------|--------------------------------------|-------|-----|-------|--------------------------------------|---------------------|-----|---------------------|
| Marker 1             | $r_1$                                | $r_2$ | ... | $r_S$ | status <sub>1</sub>                  | status <sub>2</sub> | ... | status <sub>S</sub> |
| Marker 2             | $r_1$                                | $r_2$ | ... | $r_S$ | status <sub>1</sub>                  | status <sub>2</sub> | ... | status <sub>S</sub> |
| ⋮                    | ⋮                                    | ⋮     |     | ⋮     | ⋮                                    | ⋮                   |     | ⋮                   |
| Marker M             | $r_1$                                | $r_2$ | ... | $r_S$ | status <sub>1</sub>                  | status <sub>2</sub> | ... | status <sub>S</sub> |

There are M markers producing raw data for each sensor,  $r_S$ , and sensor status information, status<sub>S</sub>. The sensor status information describes the signal strength high (peak), signal amplification (DRC), any error codes (marker status) and the signal strength low (peak nibble).

## 10.2.2 Description of Optotrak Raw and Full Raw Data Elements

The information that is reported by a Position Sensor for each marker always contains the centroid element:

**Centroid:** Centroid is the positional component of the raw (centroid) data for the sensor. If the sensor could not determine the marker centroid for any reason, the centroid value is set to the constant value BAD\_FLOAT

The Position Sensor reports the following elements when full raw data is requested:

**Signal Strength High (Peak):** Signal Strength High is written as peak within the routines. This value is a number between 0 and 255 and indicates the marker signal strength measured by the sensor.

**Amplification (DRC):** The amount of amplification applied to the marker signal is called DRC within the routines. This value is scaled between 0 and 255, where 0 is the lowest possible amplification of the marker signal.

**Error Code (Marker Status):** The error code is also referred to as the marker status. Possible values are:

**CENTROID\_OK:** the sensor successfully calculated the marker centroid.

**CENTROID\_WAVEFORM\_TOO\_WIDE:** the marker wave form was too wide to determine a valid centroid.

**CENTROID\_PEAK\_TOO\_SMALL:** the marker signal was too weak to determine a valid centroid.

**CENTROID\_PEAK\_TOO\_LARGE:** the marker signal was too strong to determine a valid centroid.

**CENTROID\_WAVEFORM\_OFF\_DEVICE:** the marker wave form was partially out of the field of view of the sensor.

**CENTROID\_FELL\_BEHIND:** the sensor failed to process the current marker's signal before required to process the next marker's signal.

**CENTROID\_LAST\_CENTROID\_BAD:** the sensor was unable to determine the last marker's centroid and therefore cannot perform interpolation for the current marker's centroid. This would only occur if the currently executing collection was configured to perform interpolation on the raw sensor values.

**CENTROID\_BUFFER\_OVERFLOW:** the internal centroid (raw) data buffer on the sensor was full, so it was unable to store this current marker's data.

**CENTROID\_MISSED\_CCD:** the sensor has determined that it has missed a time control signal from the Optotrak System Control Unit, probably because of a communication error.

**CENTROID\_BAD\_CRC:** a communication error occurred in transmission of the marker's data.

**Signal Strength Low (PeakNibble):** The signal strength low is called peak nibble within the routines and is an extra four bits of resolution for the marker signal strength. The Certus, 3020 and 2020 Position Sensor has 12-bit resolution for a digitized marker signal. The peak element of the full raw data stores the most significant 8 bits of the digitized marker signal and the PeakNibble stores the least significant 4 bits. If Certus, 3020 or 2020 Position Sensors are being used, then PeakNibble is a value between 0 and 15.

### 10.2.3 C - Type Definition of Optotrak Full Raw Data

Use the C-type definition of sample type to access the full raw data returned by the Optotrak System. There are two C-types defined for full raw data: `SensorDataType` and `FullRawDataType`. `SensorDataType` accesses the sensor status information, while `FullRawDataType` provides access to both the centroid and the sensor status information.

---

**Note** You must have a fixed number of sensors in the Optotrak System to use the definition `FullRawDataType`.

---

```
typedef struct SensorDataStruct
{
    unsigned char ucPeak;
    unsigned char ucDRC;
    unsigned char ucCode;
    unsigned char ucPeakNibble;
} SensorDataType;

typedef struct FullRawDataStruct
{
    float          fCentroid[ NUM_SENSORS];
    SensorDataType SensorData[ NUM_SENSORS];
} FullRawDataType;
```



## 10.2.4 Size Calculation

The full raw data size is the product of the numbers of markers, number of sensors and the size of the FullRawDataType structure (8 bytes):

$$\text{Frame Size} = \text{Number of Markers} \times \text{Number of Sensors} \times 8.$$

## 10.2.5 C - Type Definition of Optotrak Centroid Data

Use the C-type definition of sample type to access the centroid data returned by the Optotrak System.

```
float fcentroid[ NUM_SENSORS]
```

## 10.2.6 Size Calculation

The centroid data size is the product of the numbers of markers and the size of the centroid data (4 bytes):

$$\text{Frame Size} = \text{Number of Markers} \times \text{Number of Sensors} \times 4$$

## 10.3 Optotrak 3D Data

A 3D data frame contains the X, Y, and Z coordinates for each marker in the current collection.

### 10.3.1 Organization of Optotrak 3D Data

The information is grouped by marker; the X, Y, and Z positions for marker one are followed by the X, Y, and Z positions for marker two and so on. If the Optotrak System is unable to determine a 3D position for a certain marker, then the X, Y, and Z coordinates for that marker are set to BAD\_FLOAT. [Table 10-2 on page 242](#) is the format of a 3D data frame for M markers. The structure for multiple frames in a Northern Digital Floating Point File is described in the [“Floating Point Programmer’s Guide” on page 65](#).

**Table 10-2: An Optotrak 3D Data Frame for M Markers**

| <b>Marker</b> | <b>Coordinate</b> |                |                |
|---------------|-------------------|----------------|----------------|
| Marker 1      | X <sub>1</sub>    | Y <sub>1</sub> | Z <sub>1</sub> |
| Marker 2      | X <sub>2</sub>    | Y <sub>2</sub> | Z <sub>2</sub> |
| ⋮             | ⋮                 | ⋮              | ⋮              |
| Marker M      | X <sub>M</sub>    | Y <sub>M</sub> | Z <sub>M</sub> |

### 10.3.2 C - Type Definition of Optotrak 3D Data Structure

The sample C-type definition, Position3d, is defined in the header file ndtypes.h. An application program can use this definition to access the elements of the 3D positions.

```
typedef struct Position3dStruct
{
    float x;
    float y;
    float z;
} Position3d;
```

### 10.3.3 Size Calculation

The 3D data size is the product of the numbers of markers and the size of the Position3d structure (12 bytes):

$$\text{Frame Size} = \text{Number of Markers} \times 12.$$

## 10.4 Optotrak Rigid Body Transformation Data

### 10.4.1 Organization of Rigid Body Transformation Data

The following is an example of a table format for a rigid body transformation data frame for B rigid bodies and M markers in Euler angle representation. The structure returned by the Optotrak System is OptotrakRigidStruct.

---

**Note** The rotation components are different for each of the three transformations: Euler, quaternion and rotation matrix.

---

The header portion of the returned data includes the elements Rigid Id, Flags, Quaternion Error and Iterative Error. The C-type definition for the 3D marker positions contained in this data frame was introduced in “C - Type Definition of Optotrak 3D Data Structure” on page 242.

**Table 10-3: An Optotrak Rigid Body Transformation Data Frame for B Rigid Bodies and M Markers (Euler Angle Format)**

|              | Header   | Rotational Component |       |       | Translational Component |       |       |
|--------------|----------|----------------------|-------|-------|-------------------------|-------|-------|
| Rigid Body 1 | Header 1 | $R_z$                | $R_y$ | $R_x$ | $T_x$                   | $T_y$ | $T_z$ |
| Rigid Body 2 | Header 2 | $R_z$                | $R_y$ | $R_x$ | $T_x$                   | $T_y$ | $T_z$ |
| ⋮            |          | ⋮                    | ⋮     | ⋮     | ⋮                       | ⋮     | ⋮     |
| Rigid Body B | Header B | $R_z$                | $R_y$ | $R_x$ | $T_x$                   | $T_y$ | $T_z$ |
| Marker 1     |          | $X_1$                |       | $Y_1$ | $Z_1$                   |       |       |
| Marker 2     |          | $X_2$                |       | $Y_2$ | $Z_2$                   |       |       |
| ⋮            |          | ⋮                    |       | ⋮     | ⋮                       |       |       |
| Marker M     |          | $X_M$                |       | $Y_M$ | $Z_M$                   |       |       |

## 10.4.2 Description of Rigid Body Real-time Data Elements

### Header

**RigidId:** The RigidId is the identifier assigned to the rigid body when added using the RigidBodyAdd or RigidBodyAddFromFile routine.

**flags:** Flags are the bit status flags for the rigid body transformation. If the bit flag `OPTOTRAK_UNDETERMINED_FLAG` is set, then the Optotrak System was unable to determine a transformation for this rigid body. For a complete description of all the flags that can be set in the element, refer to the documentation for the routine “RigidBodyChangeSettings” on page 211.

**QuaternionError:** The QuaternionError is the RMS quaternion error for the determined rigid body transformation, if the quaternion algorithm was used.

**IterativeError:** The IterativeError is the RMS error for the determined rigid body transformation, if the iterative Euler angle algorithm was used.

### Rotation

The rotation element varies with the transformation type. For a description of the rotation values see [“Sample Formats” on page 68](#).

### Translation

The translation element is reported for all three representations and describes the motion of the origin of the rigid body.

## 10.4.3 C - Type Definition of Rigid Body Transformation Data Structure

There are several pre-defined C-structures for accessing the transformation data found in the include files, ndopto.h and ndtypes.h. This include file is available on the API CD. The OptotrakRigidStruct structure allows an application program to access the transformation.

```
struct OptotrakRigidStruct
{
    long int          RigidId;
    long int          flags;
    float             QuaternionError,
    float             IterativeError;
    union TransformationUnion transformation;
};
```

The transformation returned by the Optotrak System is stored and accessed using a C-union. A C-union is used because it can be expressed in one of three formats: Euler Angle, Quaternion, and Rotation Matrix. Each of these formats consists of a rotational component and a translational component. The rotational component is expressed differently among the different transformation formats. The following tables illustrate each transformation format and how to access their elements.

**Table 10-4: Transformation Formats**

| <b>Format Accessed</b> | <b>Structure Member</b>        |
|------------------------|--------------------------------|
| Euler Angle            | Xfrm.transformation.euler      |
| Quaternion             | Xfrm.transformation.quaternion |
| Rotation Matrix        | Xfrm.transformation.rotation   |

**Table 10-5: Accessing Euler Angle Elements**

| <b>Euler Angle Element</b>   | <b>Structure Member</b> |
|------------------------------|-------------------------|
| rotation about the x-axis    | euler.rotation.yaw*     |
| rotation about the y-axis    | euler.rotation.pitch*   |
| rotation about the z-axis    | euler.rotation.roll*    |
| translation along the x-axis | euler.translation.x     |
| translation along the y-axis | euler.translation.y     |
| translation along the z-axis | euler.translation.z     |

\* angles are in radians

**Table 10-6: Accessing Quaternion Elements**

| <b>Quaternion Element</b>                    | <b>Structure Member</b>  |
|----------------------------------------------|--------------------------|
| first element in the quaternion orientation  | quaternion.rotation.q0   |
| second element in the quaternion orientation | quaternion.rotation.qx   |
| third element in the quaternion orientation  | quaternion.rotation.qy   |
| fourth element in the quaternion orientation | quaternion.rotation.qz   |
| translation along the x-axis                 | quaternion.translation.x |
| translation along the y-axis                 | quaternion.translation.y |
| translation along the z-axis                 | quaternion.translation.z |

**Table 10-7: Accessing Rotation Matrix Elements**

| Rotation Matrix Elements                             | Structure Member       |
|------------------------------------------------------|------------------------|
| a 3 x 3 floating point array;<br>the rotation matrix | rotation.matrix[i][i]  |
| translation along the x-axis                         | rotation.translation.x |
| translation along the y-axis                         | rotation.translation.y |
| translation along the z-axis                         | rotation.translation.z |

#### 10.4.4 Size Calculation

To determine the size of the data frame, multiply the number of rigid bodies in your collection by the size of the rigid body transformation data (64 + 16 bytes header data), and add to that the number of markers in the collection multiplied by the size of a 3D position.

$$\begin{aligned} \text{Frame Size} &= (\text{Number of Rigid Bodies} \times \text{Rigid Body Size}) + \\ &(\text{Number of Markers} \times 12) \\ &= (B \times (24 + 16)) + (M \times 12) \quad \text{Euler Transformation} \\ &= (B \times (28 + 16)) + (M \times 12) \quad \text{Quaternion Transformation} \\ &= (B \times (48 + 16)) + (M \times 12) \quad \text{Matrix Transformation} \end{aligned}$$

### 10.5 ODAU Raw Data

An ODAU raw data frame contains raw data for each analog channel being sampled, followed by digital data, if the Digital I/O port is configured for input. Because both the analog and the digital data are of the pre-defined C-type integer, it is not necessary for an application program to define a special type to access the data cleanly. See [Table 10-8 on page 247](#) for the organization of the elements in a general ODAU raw data frame.

**Table 10-8: The Organization of the Elements in a General ODAU Raw Data Frame**

| Channel              | Data           |
|----------------------|----------------|
| Channel 1            | a <sub>1</sub> |
| Channel 2            | a <sub>2</sub> |
| ⋮                    | ⋮              |
| Channel <sub>C</sub> | a <sub>C</sub> |
| Digital Data Input   | D              |

where:

- a = raw analog data
- D = digital data
- C = number of channels

If only raw data is sampled, the digital information is not returned.

A simple method for accessing the ODAU raw data is to cast it as an array of integers. The array would need one element for each analog channel being sampled, plus one element for the digital data if it were being sampled.

### 10.5.1 C - Type Definition for ODAU Raw Data

```
int OdauData[ ANALOG_CHANNELS + DIGITAL];
```

### 10.5.2 Size Calculation

Determine the size of the data frame multiplying the number of analog channels being sampled by the size of a 16-bit integer (2 bytes), and add another 2 bytes if digital data are being collected:

$$\text{Frame Size} = \text{Number of Analog Channels} \times 2 [+ 2 \text{ if digital input}]$$

---

**Note** Use the 2 byte integer size for both 32-bit and 16-bit applications, even though the size is specific to 16-bit applications. The data is formatted to this size by the system.

---





---

## Appendix A Libraries and Sample Application Programs

This appendix is divided into three sections. The first section describes the structure of the API installation CD. The remaining two sections describe the sample programs included with the API, as well as describe the operations these sample programs perform. This information is intended to explain the programmer's interface capabilities.

### A.1 API Installation CD

The API installation CD contains all the files required to build the sample application programs, as well as your own application programs. One CD is provided for both the PC version and Workstation version of the API.

To install the Optotrak API files on a Windows platform, insert the CD into your CD-ROM drive and click the Install button of the Optotrak Application Programmer's Interface dialog.

To install the Optotrak API files on another platform, consult the appropriate readme.txt file on the API installation CD.

Upon installing the Optotrak API files, three new directories will be created:

- **NDIoapi** this folder contains the library files and the Optotrak API sample programs.
- **ndlib** this folder contains the Optotrak API library files for the compiler(s) selected during installation. This folder also contains the library include files.
- **samples** this folder contains sample programs that show how the Optotrak API routines are used to communicate with an Optotrak System.

### A.2 Sample Programs for All Optotrak Systems

The following sample programs can be run on any Optotrak System. They are separated into five functional groups: Optotrak, Rigid Body, ODAU, File Processing, and Secondary Host.

- Optotrak specific samples demonstrate the capabilities of the API in conjunction with the standard Optotrak System.
- Rigid Body specific samples demonstrate how to use the API routines and data types if the Optotrak System includes either the Optotrak Real-time Rigid Body Option, the Accelerated Processing/Extended Buffering Option, or when on-host conversions are enabled.

- ODAU samples are examples of possible application programs when the configuration of the Optotrak System includes an ODAU device.
- File processing specific samples show how to use the API routines that manipulate Northern Digital Floating Point files.
- Secondary Host specific samples show how to use the API to interact with the Optotrak System when two computers are connected to the system.

---

**An ODAU device must not be connected to the Optotrak System to run these sample programs, except for the ODAU specific programs (Sample Programs 12 and 13)**

---

Table A-1: Sample Program Descriptions

| Functional Group                  | Sample Program                                                                  | Description                                                                                                                                                                     | Location |
|-----------------------------------|---------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| Optotrak Specific Sample Programs | 1                                                                               | The simplest application program provided, illustrating the steps required to initialize the Optotrak System, and then retrieve and display the Optotrak System status.         | page 252 |
|                                   | 2                                                                               | Shows the steps required to retrieve and display real-time 3D data from an Optotrak System.                                                                                     | page 253 |
|                                   | 3                                                                               | Illustrates the retrieval/display of real-time data, and the spooling of Optotrak System buffered data to a Northern Digital Floating Point File.                               | page 253 |
|                                   | 4                                                                               | Shows you how to spool buffered raw sensor data to a block of memory allocated by the application program.                                                                      | page 254 |
|                                   | 5                                                                               | Shows how to use the non-blocking method of retrieving real-time data, while at the same time spooling buffered data to file.                                                   | page 254 |
|                                   | 6                                                                               | Demonstrates the non-blocking method of spooling buffered data to a memory block, while at the same time retrieving and displaying real-time 3D data using the blocking method. | page 255 |
|                                   | 7                                                                               | Demonstrates a practical application of the non-blocking method for spooling buffered data.                                                                                     | page 255 |
|                                   | 8                                                                               | Shows how to collect a file of raw sensor data and convert it to a 3D data file.                                                                                                | page 256 |
|                                   | 18                                                                              | Shows how to change the camera's measurement coordinate system to allow the Optotrak System to compute 3D marker positions in a new frame.                                      | page 256 |
|                                   | 20                                                                              | Demonstrates TransputerDetermineSystemCfg.                                                                                                                                      | page 257 |
| 21                                | Demonstrates OptotrakGetCameraParameter-Status and OptotrakSetCameraParameters. | page 257                                                                                                                                                                        |          |

Table A-1: Sample Program Descriptions

| Functional Group                         | Sample Program | Description                                                                                                                              | Location                 |
|------------------------------------------|----------------|------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|
| Rigid Body Specific Sample Programs      | 9              | Shows how to load a rigid body, then retrieve and display real-time rigid body transformation data.                                      | <a href="#">page 258</a> |
|                                          | 10             | Shows how to load a rigid body, then retrieve and display real-time rigid body transformation data.                                      | <a href="#">page 259</a> |
|                                          | 11             | Shows how to load several rigid bodies, and use one of them to specify a new coordinate system for the rigid body calculations.          | <a href="#">page 260</a> |
|                                          | 19             | Shows how to convert previously collected raw data to their corresponding 3D positions, and then transforms the 3D positions to 6D data. | <a href="#">page 260</a> |
| ODAU Specific Sample Programs            | 12             | Shows how to configure Optotrak System collections when an ODAU device is connected.                                                     | <a href="#">page 261</a> |
|                                          | 13             | Shows how to configure Optotrak System collections when an ODAU device is connected.                                                     | <a href="#">page 261</a> |
| File Processing Specific Sample Programs | 14             | Shows how to collect an Optotrak System raw data file or an ODAU raw data file, and convert the file to the appropriate format.          | <a href="#">page 262</a> |
|                                          | 15             | Shows how to use the file processing capabilities of the Optotrak API.                                                                   | <a href="#">page 263</a> |
| Secondary Host Specific Sample Programs  | 16             | Shows the simplest form of a secondary host application program.                                                                         | <a href="#">page 263</a> |
|                                          | 17             | Presents a secondary host application that spools 3D data to file.                                                                       | <a href="#">page 264</a> |

## A.2.1 Optotrak Specific Sample Programs

### Sample Program 1

Sample Program 1 is the simplest sample application program provided. This sample shows you the steps required to initialize the Optotrak System, and then retrieve and display the Optotrak System status. The sample program follows these steps:

1. Loads the system of processors with the appropriate transputer programs and startup code.
2. Initiates communication with the system of processors.
3. Loads the appropriate camera parameters.
4. Request/receives/displays the current Optotrak System status. Passes NULL for those status variables that are not requested.
5. Disconnects the application program from the system of processors.

## Sample Program 2

Sample Program 2 shows you the steps required to retrieve and display real-time 3D data from the Optotrak System. This sample uses the blocking method of data retrieval by following these steps:

1. Loads the system of processors with the appropriate transputer programs and startup code.
2. Initiates communications with the system of processors.
3. Sets the optional processing flags to do the 3D conversions on the host computer.
4. Loads the appropriate camera parameters.
5. Sets up an Optotrak System collection.
6. Activates the IRED markers.
7. Requests/receives/displays 10 frames of real-time 3D data.
8. De-activates the markers.
9. Disconnects the application program from the system of processors.

## Sample Program 3

Sample Program 3 illustrates two main functions of the Optotrak System: the retrieval/display of real-time data, and the spooling of buffered data to a Northern Digital Floating Point File. The real-time data is raw sensor data. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.

5. Activates the IRED markers.
6. Requests/receives/displays 10 frames of real-time full raw data.
7. Initializes a file for spooling of data.
8. Collects and spools 3D data to disk.
9. De-activates the markers.
10. Disconnects the application program from the system of processors.

### **Sample Program 4**

Sample Program 4 shows you how to spool buffered raw sensor data to a block of memory allocated by the application program. The raw data is spooled to the memory block, and then displayed on the monitor. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Initializes a memory block for spooling of data.
7. Spools raw data to the memory block.
8. De-activates the markers.
9. Prints the contents of the memory block to the screen.
10. Disconnects the application program from the system of processors.

### **Sample Program 5**

Sample Program 5 shows you how to use the non-blocking method of retrieving real-time data, while at the same time spooling buffered data to file. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Initializes a data file for spooling data.
7. Starts spooling raw sensor data.

8. Spools raw data to file while at the same time requesting and displaying real-time 3D data.
9. De-activates the markers.
10. Disconnects the application program from the system of processors.

### **Sample Program 6**

Sample Program 6 demonstrates the non-blocking method of spooling buffered data to a memory block, while at the same time retrieving and displaying real-time 3D data using the blocking method. The following list details the steps performed by this sample program:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Initializes an area of memory for spooling data.
7. Starts spooling 3D data.
8. Spools the 3D data to memory while, at the same time, requesting and displaying real-time 3D data.
9. De-activates the markers.
10. Disconnects the application program from the system of processors.

### **Sample Program 7**

Sample Program 7 is a practical application of the non-blocking method for spooling buffered data. The application prepares the system to spool buffered data to disk, then waits for Marker 1 to come into view of the Optotrak System. The application program then starts spooling the buffered data to file. Buffered data is spooled until the entire trial has been spooled or the application program detects that Marker 1 has gone out of view. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Initializes a file for spooling 3D data.

7. Starts spooling when Marker 1 is in view.
8. Spools 3D data to file while, at the same time, requesting and examining 3D data.
9. Stops the pool once Marker 1 goes out of view or after 100 seconds of data are spooled.
10. De-activates the markers.
11. Disconnects the application program from the system of processors.

### **Sample Program 8**

Sample Program 8 show you how to collect a file of raw sensor data and convert it to a 3D data file. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Initializes a file for spooling raw data.
7. Spools the raw data to file.
8. De-activates the IRED markers.
9. Converts the raw data file to a 3D data file.
10. Disconnects the application program from the system of processors.

### **Sample Program 18**

This program shows you how to change the camera's measurement coordinate system to allow the Optotrak System to compute 3D marker positions in a new frame. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Requests and averages 50 frames of real-time 3D data.
7. De-activates the markers.
8. Creates a new camera parameter file using the averaged real-time 3D data.



9. Prints out the alignment transformation fit errors.
10. Loads the new camera parameters into the Optotrak System.
11. Activates the IRED markers.
12. Requests/receives/displays 10 frames of 3D data in the new coordinate system.
13. De-activates the IRED markers.
14. Disconnects the application program from the system of processors.

## Sample Program 20

This program demonstrates the routine `TransputerDetermineSystemCfg`, which can be used to determine the system configuration from within an application program. This is an alternative to determining the system configuration beforehand with one of the command-line utility programs (`optset32.exe` or `buildnif`). The program follows these steps:

1. Determines the system configuration and generates the external network information configuration file, `system.nif`, which is stored in the default location in the `ndigital\realtime` subdirectory.
2. Loads the system of processors with the appropriate transputer programs and startup code based on the newly generated external file, `system.nif`.
3. Initiates communications with the system of processors.
4. Disconnects the application program from the system of processors.
5. Sets the API to use the internally generated and stored system network information file configuration.
6. Determines the system configuration and store the network information internally.
7. Loads the system of processors with the appropriate transputer programs and startup code based on the internal network information file configuration.
8. Initiates communications again with the system of processors.
9. Disconnects the application program from the system of processors again.

## Sample Program 21

This program demonstrates the routines `OptotrakGetCameraParameterStatus` and

OptotrakSetCameraParameters. These can be used both to obtain status information about the camera parameter sets stored in the extended camera parameter files and to select a new camera parameter set to use for raw sensor data to 3D conversions. The program follows these steps:

1. Sets the processing flags to do the 3D conversions on the host computer. Unlike Sample Program 2, this is not optional. This must be done to use different camera parameter sets.
2. Loads the system of processors with the appropriate transputer programs and startup code.
3. Initiates communications with the system of processors.
4. Loads the appropriate camera parameters.
5. Obtains the status information for the parameter sets stored in the camera parameter files, as well as status information for the current parameter set.
6. Selects a new camera parameter set based on the given marker type, marker wavelength, and camera model type.
7. Sets up an Optotrak System collection.
8. Activates the IRED markers.
9. Requests/receives/displays 10 frames of real-time 3D data.
10. De-activates the markers.
11. Disconnects the application program from the system of processors.

### A.2.2 Rigid Body Specific Sample Programs

#### Sample Program 9

Sample Program 9 shows you how to use the rigid body specific routines to load a rigid body, then retrieve and display real-time rigid body transformation data. This sample assumes:

- there is a six marker rigid body connected to the Optotrak System.
- the rigid body file plate.rig, which defines the rigid body in its home position, exists in the appropriate directory.

The program follows these steps:

1. Loads the system of processors with the appropriate transputer programs and startup code.
2. Initiates communications with the system of processors.
3. Sets the optional processing flags to do the 3D conversions on the host computer.
4. Loads the appropriate camera parameters.
5. Sets up an Optotrak System collection.
6. Adds a rigid body to the Optotrak System tracking list using a rigid body file.
7. Activates the IRED markers.
8. Requests/receives/displays 10 frames of rigid body transformations.
9. De-activates the IRED markers.
10. Disconnects the application program from the system of processors.

### Sample Program 10

Sample 10 shows you how to use the rigid body specific routines to load a rigid body, then retrieve and display real-time rigid body transformation data.

The routine `RigidBodyChangeSettings` is also used to alter a number of the parameters which affect the determination of the rigid body transformations. For example, the return format for the rigid body is changed from the default Euler angle format to quaternion format.

This sample assumes:

- there is a six marker rigid body connected to the Optotrak System.
- the rigid body file `plate.rig`, which defines the rigid body in its home position, exists in the appropriate directory.

The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Adds a rigid body to the Optotrak System tracking list using a rigid body file.
7. Changes the default settings for the rigid body just added.

8. Requests/receives/displays 10 frames of rigid body transforms in Quaternion format; also display the attached 3D marker position data.
9. De-activates the IRED markers.
10. Disconnects the application program from the system of processors.

### Sample Program 11

Sample 11 shows you how to use several rigid bodies. After the rigid bodies have been loaded, one of them is used to specify a new coordinate system for the rigid body calculations.

This sample assumes that there are two rigid bodies connected to the Optotrak System, as indicated by the RigidBodyAdd routines. The program's source code must be modified to load the definitions of the connected rigid bodies before transformations are determined correctly. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Adds a rigid body to the Optotrak System's tracking list using an array of 3D points.
7. Adds a rigid body to the Optotrak System tracking list using a rigid body file.
8. Changes the coordinate system for the rigid body calculations.
9. Requests/receives/displays 10 frames of rigid body transformations; also displays attached 3D marker position data.
10. De-activates the IRED markers.
11. Disconnects the application program from the system of processors.

### Sample Program 19

This program shows you how to convert previously collected raw data to their corresponding 3D positions, and then transforms the 3D positions to 6D data. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Adds a rigid body to the Optotrak System tracking list using an array of 3D points.

6. Activates the IRED markers.
7. Initializes a file for spooling raw data.
8. Spools the raw data to the file.
9. De-activates the markers.
10. Opens the raw data file as an input file (for reading).
11. Reads the raw data from the input file.
12. Converts the raw data to 3D positions and prints them to the screen.
13. Transforms the 3D positions to 6D data and prints them to the screen.
14. Disconnects the application program from the system of processors.

### **A.2.3 ODAU Specific Sample Programs**

#### **Sample Program 12**

Sample 12 shows you how to configure the collections when an ODAU device is connected to the Optotrak System. Once the collections have been configured, this sample requests real-time ODAU raw data and displays it as voltages.

This sample assumes that the ODAU device is properly connected, and that the required network information file has been generated. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an ODAU collection.
5. Sets up an Optotrak System collection.
6. Activates the IRED markers.
7. Requests/receives/displays 10 frames of ODAU data.
8. De-activates the markers.
9. Disconnects the application program from the system of processors.

#### **Sample Program 13**

This sample shows you how to configure Optotrak System collections when an ODAU device is connected. After the collections are configured, buffered 3D data

and buffered ODAU raw data is spooled to the files c#001.s13 and o1#001.s13, respectively.

This sample assumes that the ODAU device is properly connected and the required network information file has been generated. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an ODAU collection.
5. Sets up an Optotrak System collection.
6. Activates the IRED markers.
7. Initializes files for spooling the Optotrak System and ODAU data.
8. Spools the Optotrak System and ODAU data to the disk files.
9. De-activates the markers.
10. Disconnects the program from the system of processors.

### A.2.4 File Processing Specific Sample Programs

#### Sample Program 14

This sample shows you how to collect both an Optotrak System raw data file and an ODAU raw data file, and convert the files to the appropriate formats. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an ODAU collection.
5. Sets up an Optotrak System collection.
6. Activates the IRED markers.
7. Initializes files for spooling the Optotrak System and ODAU data.
8. Spools the Optotrak System and ODAU data to the disk files.
9. De-activates the markers.
10. Converts the Optotrak System raw data file to a 3D format file.
11. Converts the ODAU raw data file to a voltage format file.
12. Disconnects the application program from the system of processors.

## Sample Program 15

This sample shows you how to use the file processing capabilities of the Optotrak API. A 3D data file is collected when marker 1 comes into view. After the file has been collected, it is opened as input and a new file is opened for output. The first frame of data in the input file is subtracted, on a marker by marker basis, from all frames in the input file. The resultant frames are written to the output file. The program follows these steps:

- 1-3. See Sample Program 1.
4. Sets up an Optotrak System collection.
5. Activates the IRED markers.
6. Waits until marker 1 comes into view.
7. Initializes a file for spooling 3D data.
8. Spools the 3D data to the file.
9. De-activates the markers.
10. Disconnects the application program from the system of processors.
11. Opens the new 3D file as an input file (for reading).
12. Opens another file as an output file (for writing).
13. Writes the data from the input file to the output file after performing the calculation.
14. Closes both the input and the output files.

### A.2.5 Secondary Host Specific Sample Programs

#### Sample Program 16

Sample Program 16 shows you the simplest form of a secondary host application program. Its main function is to retrieve and display real-time raw sensor data once a collection has been configured by the primary host computer. The program follows these steps:

1. Initializes the communication with the system of processors.
2. Gets the current Optotrak System status.
3. Allocates memory for storing the real-time raw data.

4. Requests/receives/displays 10 frames of real-time raw data.
5. Disconnects the application program from the system of processors.

### Sample Program 17

Sample Program 17 presents a secondary host application program that spools 3D data to file.

This secondary host sample application assumes that the primary host has configured the collection such that the Optotrak System buffers 3D data and not raw sensor data. The program follows these steps:

1. Initializes the communication with the system of processors.
2. Gets the current Optotrak System status.
3. Allocates memory for storing the real-time 3D data.
4. Initializes a file for spooling the 3D data.
5. Starts the spooling when Marker 1 is in view.
6. Spools 3D data to file while at the same time requesting and examining 3D data.
7. Stops the spool if Marker 1 goes out of view, or after all the buffered data has been spooled.
8. Disconnects the PC application program from the system of processors.

## A.3 Sample Programs for Optotrak Certus Systems

The following sample programs, which can only be run on an Optotrak Certus System, are separated into two functional groups: Optotrak and Rigid Body.

- Optotrak specific samples show the capabilities of the API in conjunction with the standard Optotrak Certus System.
- Rigid Body specific samples show how to use the API routines and data types if the Optotrak Certus System includes either the Optotrak Real-time Rigid Body Option, the Accelerated Processing/Extended Buffering Option, or when on-host conversions are enabled.



Table 10-9: Sample Program Descriptions

| Functional Group                    | Optotrak Certus Sample Program | Description                                                                                                                                          | Location |
|-------------------------------------|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| Optotrak Specific Sample Programs   | 1                              | Is the simplest sample application program provided. This sample uses basic calls to communicate with an Optotrak Certus System.                     | page 266 |
|                                     | 2                              | Demonstrates how to free an allocated device handle.                                                                                                 | page 266 |
|                                     | 3                              | Demonstrates how to get the number of properties associated with a handle, as well as retrieve one of the handle's properties.                       | page 267 |
|                                     | 4                              | Demonstrates how to get a device handle's properties, set one of those properties to a new value, and verify that the handle received the new value. | page 267 |
|                                     | 11                             | Demonstrates how to make a visible LED— associated with a specified handle—blink.                                                                    | page 268 |
|                                     | 12                             | Demonstrates how to control a beeper device.                                                                                                         | page 268 |
|                                     | 13                             | Demonstrates how to retrieve switch data.                                                                                                            | page 269 |
|                                     | 15                             | Demonstrates how to retrieve and edit a device handle properly.                                                                                      | page 269 |
|                                     | 16                             | Demonstrates how to use the blocking and non-blocking methods for retrieving real-time data.                                                         | page 270 |
|                                     | 17                             | Demonstrates how to retrieve and display real-time 3D data.                                                                                          | page 270 |
|                                     | 18                             | Demonstrates how to loop the real-time 3D data retrieval process.                                                                                    | page 271 |
| Rigid Body Specific Sample Programs | 14                             | Demonstrates how to collect rigid body data.                                                                                                         | page 271 |
| ODAU Specific Sample Programs       | 19                             | Demonstrates how to collect ODAU raw data.                                                                                                           | page 272 |

## A.3.6 Optotrak Specific Sample Programs

### Optotrak Certus Sample Program 1

Optotrak Certus Sample Program 1 is the simplest sample application program provided. This sample program uses basic calls to communicate with an Optotrak Certus System. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Sets up the collection.
7. Activates the markers.
8. Retrieves 3D data.
9. De-activates the markers.
10. Shuts down the system of processors.

### Optotrak Certus Sample Program 2

Optotrak Certus Sample Program 2 illustrates how to free an allocated device handle. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Frees an allocated device handle.
7. Attempts to access that device handle.

8. Shuts down the system of processors.

### **Optotrak Certus Sample Program 3**

Optotrak Certus Sample Program 3 illustrates how to get the number of properties associated with a handle, as well as retrieve one of the handle's properties. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Gets the number of properties associated with a device.
7. Gets a specific property value.
8. Shuts down the system of processors.

### **Optotrak Certus Sample Program 4**

Optotrak Certus Sample Program 4 illustrates how to get a device handle's properties, set one of those properties to a new value, and verify that the handle received the new value. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Gets number of properties associated with a device.
7. Gets all of the properties associated with that device.
8. Sets a specific property value.

9. Gets all of the properties associated with the device to ensure the new property value was received.
10. Shuts down the system of processors.

### **Optotrak Certus Sample Program 11**

Optotrak Certus Sample Program 11 illustrates how to make a visible LED—associated with a specified handle—blink. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Sets the visible LED status to VLEDST\_BLINK.
7. Shuts down the system of processors.

### **Optotrak Certus Sample Program 12**

Optotrak Certus Sample Program 12 illustrates how to control a beeper device. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Sets the beeper device to 'on'.
7. Sets the beeper device to 'off'.
8. Shuts down the system of processors.

---

## Optotrak Certus Sample Program 13

Optotrak Certus Sample Program 13 illustrates how to retrieve switch data. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Sets up the collection.
7. Activates the markers.
8. Retrieves 3D data.
9. Retrieves switch data.
10. De-activates the markers.
11. Shuts down the system of processors.

## Optotrak Certus Sample Program 15

Optotrak Certus Sample Program 15 illustrates how to retrieve and edit a device handle properly. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Initializes the devices.
5. Retrieves information for all devices connected.
6. For each device, retrieves and sets the name property value for the device handle.
7. Shuts down the system of processors.

## **Optotrak Certus Sample Program 16**

Optotrak Certus Sample Program 16 illustrates both the blocking and non-blocking methods of real-time data retrieval. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Initializes the devices.
4. Loads the Camera Parameters.
5. Sets system to perform data conversion on the host computer.
6. Sets up the collection.
7. Activates the markers.
8. Retrieves and displays a frame of raw data using the non-blocking method.
9. Retrieves and displays a frame of 3D data using the blocking method.
10. De-activates the markers.
11. Shuts down the system of processors.

## **Optotrak Certus Sample Program 17**

Optotrak Certus Sample Program 17 illustrates how to retrieve and display real-time 3D data. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Initializes the devices.
4. Sets system to perform data conversion on the host computer.
5. Loads the Camera Parameters.
6. Sets up the collection.
7. Activates the markers.
8. Retrieves and displays a frame of 3D data using the blocking method.
9. De-activates the markers.
10. Shuts down the system of processors.

---

## Optotrak Certus Sample Program 18

Optotrak Certus Sample Program 18 illustrates how to loop the real-time 3D data retrieval process. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Initializes the devices.
4. Sets system to perform data conversion on the host computer.
5. Loads the Camera Parameters.
6. Sets up the collection.
7. De-activates the markers.
8. Loops the retrieval and display of 3D data using the blocking method.
9. Activate the markers.
10. Stops loop if the strober configuration or switch data changes.
11. De-activates the markers.
12. Shuts down the system of processors.

### A.3.7 Rigid Body Specific Sample Programs

## Optotrak Certus Sample Program 14

Optotrak Certus Sample Program 14 illustrates how to collect rigid body data. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Gets the current Optotrak System status to retrieve the system version. If the system version is that of an Optotrak Certus, the program continues.
4. Loads the Camera Parameters.
5. Initializes the devices.
6. Sets up the collection.
7. Activates the markers.

8. Retrieves 3D data.
9. Adds rigid body information from a device handle.
10. Gets the rigid body data to ensure new data accepted.
11. De-activates the markers.
12. Shuts down the system of processors.

### **A.3.8 ODAU Specific Sample Programs**

#### **Optotrak Certus Sample Program 19**

Optotrak Certus Sample Program 19 illustrates how to retrieve and display ODAU raw data. The sample program follows these steps:

1. Determines the system configuration.
2. Initializes the communication with the system of processors.
3. Initializes the devices.
4. Sets system to perform data conversion on the host computer.
5. Loads the Camera Parameters.
6. Sets up the ODAU collection.
7. Sets up the Optotrak collection.
8. Activates the markers.
9. Loops the retrieval and display of 3D data and ODAU raw data.
10. Stops loop if the strober configuration or switch data changes.
11. De-activates the markers.
12. Shuts down the system of processors.



---

## Appendix B Error Messages and Constants

This section of the guide explains the various error messages and constants used throughout the Optotrak Application Programmer's Interface.

Each section lists the error message/constant and a brief description of the cause of the failure. The error messages are grouped according to the area of operation they are related to, and listed alphabetically within each group.

In general, remember that an error message may be generated much later in the program than the original coding error. Also, communication problems may be caused by hardware; check that all connections are secure and the power is on.

### B.1 Error Constants

All routines in the Optotrak API use a return value to indicate whether the routine has completed successfully. If a routine is successful, a zero value is returned, otherwise, one of the following codes is returned.

#### OPTO\_SYSTEM\_ERROR\_CODE

This value indicates that the routine failed because of a system related problem. For example, the application program tried to send a message to an ODAU device that did not respond.

#### OPTO\_USER\_ERROR\_CODE

This value indicates that the routine failed because of an invalid parameter. For example, the application program attempted to open a new file with an ID that was already in use.

### B.2 Error Messages

All error messages are recorded in the file `opto.err` if you have included the `OPTO_LOG_ERRORS_FLAG` in the routine `TransputerInitializeSystem` (see sample code from [“Program Sample Showing How to Initialize and Retrieve System Status” on page 25](#)).

Alternatively, if you are only interested in the most recent error, use the routine `OptotrakGetErrorString`.

For example, if the `OptotrakGetStatus` routine was unsuccessful because of a communications problem, looking in the `opto.err` file might reveal the information shown below, depending on the exact failure.

**ERROR: MessageSystem: Failed sending message preamble**

**ERROR: Optotrak: Failed requesting status**

If the routine `OptotrakGetErrorString` was used (method 2,) the null terminated string “Optotrak: Failed requesting status” would be returned.

---

**Note** After a call to the routine `TransputerLoadSystem`, `OptotrakSetupCollection`, `RigidBodyAdd` and `RigidBodyAddFromFile`, it is advisable to include a sleep routine to allow enough time for the routine to finish. The length of time required will depend on the speed of the host computer — the sample programs use a one second delay. If the sleep time is too short, the routine will fail and error messages may be generated.

---

### B.3 Message System Related Error Messages

**MessageSystem: Destination XX is not present: message ID YYYY.**

An invalid destination was specified when sending a message. XX is the destination device ID for the message. YYYY is the ID of the message the application program was attempting to send.

**MessageSystem: Failed receiving data. Last Message XXXX.**

Failed to receive the data portion of a message meant for the application program. XXXX is the value of the last message ID sent by the application program. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed receiving message header.**

The host computer has failed to successfully receive the header portion of a message. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed receiving message preamble.**

The host computer has failed to receive the preamble portion of a message. The message preamble is the first four bytes of a message that identifies it as a standard network message to the Optotrak System. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed receiving unsuccessful data: last ID XXXX.**

Failed receiving the data portion of an unsuccessful message. XXXX is the value of the last message ID sent by the application program. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed receiving YYYY bytes of data: last ID XXXX.**

Failed to receive the data portion of a message sent to the application program. XXXX is the value of the last message ID sent by the application program. YYYY is the size in bytes of the data portion of the message the host computer failed to receive. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed sending message data.**

Failed sending the data portion of the message. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed sending message header.**

Failed sending the header portion of the message. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Failed sending message preamble.**

Failed sending the message preamble. The message preamble is the first four bytes of a message that identifies it as a standard network message to the Optotrak System. There is a communication problem between the host computer and the Optotrak System.

**MessageSystem: Invalid link address XXX.**

The currently specified link adapter address is invalid. XXX is the decimal value of the currently specified address. This error may be caused either by a REALTIME environment variable address that does not match the actual address on the Northern Digital Interface adapter card, or there is another card in the PC that is conflicting with the Northern Digital Interface adapter card.

**MessageSystem: Received invalid message preamble XXXX.**

The host computer has received a message preamble that it is either invalid or not meant for the host computer. XXXX is the value the host computer

received as a message preamble. There is a communication problem between the host computer and the Optotrak System.

### **MessageSystem: Receiving invalid data at link address XXX.**

The host computer is currently receiving invalid data at the currently specified link adapter address. XXX is the decimal value of the currently specified link adapter address. This may be caused either by a REALTIME environment variable address which does not match the actual address on the Northern Digital Interface adapter card, or there is another card in the PC that is conflicting with the Northern Digital Interface adapter card.

### **MessageSystem: Too many ODAUs defined.**

Too many ODAU devices are connected in the Optotrak System. The maximum number of ODAUs that can be used in the Optotrak System is four.

### **MessageSystem: Too many sensors defined.**

Too many sensors are connected in the Optotrak System. The maximum number of sensors that can be used in the Optotrak System is 24.

### **MessageSystem: Unable to initialize communications with the Optotrak.**

The host computer could not open a communications link with the Optotrak System while attempting to set up the message system. There is a communication problem between the host computer and the Optotrak System.

### **MessageSystem: Unable to open OPTO.ERR file.**

Failed to open the file stream, opto.err, in the current directory. All error messages are logged to this file, if requested by the application program.

### **MessageSystem: Unknown message ID XXXX when getting configuration.**

Received an unknown message when retrieving configuration information from the Optotrak System. XXXX is the value of the message ID received. This is an internal message that should never appear. This may indicate an incompatibility between the API code and the TLD files in use. Contact NDI for assistance in resolving this error.

### **MessageSystem: Unsuccessful received: last ID XXXX.**

An unsuccessful message was received by the application program, indicating that the last request made by the host computer could not be carried

out. XXXX is the value of the last message ID sent by the application program.

## B.4 Transputer Related Error Messages

### **Transputer: Error downloading transputer code.**

The program failed to successfully download the transputer code in the specified network information file. The specified network information file may be incorrect, or there may be a hardware communications problem.

### **Transputer: Error shutting down collections.**

The program failed sending shutdown collection messages to the Optotrak System and ODAU devices. This error condition can occur as the `TransputerShutdownSystem` routine attempts to shutdown the Optotrak System and ODAU collections before discontinuing communication with the Optotrak System.

### **Transputer: REALTIME environment variable not defined.**

The REALTIME environment variable is not set to the address of the Northern Digital Interface adapter card. Refer to the Optotrak System Guide for details.

## B.5 Optotrak Related Error Messages

### **Optotrak: bRigidOnHost must be enabled if bConvertOnHost is enabled; bRigidOnHost set TRUE.**

An attempt was made to set the flag `bRigidOnHost` FALSE while the flag `bConvertOnHost` is set TRUE. This is an invalid state, and the API corrects it by setting `bRigidOnHost` TRUE. `bRigidOnHost` can only be set FALSE if `bConvertOnHost` is also FALSE.

### **Optotrak: Error shutting down collections.**

Failed sending shutdown collection messages to the Optotrak System and ODAU devices. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed finding desired camera parameter files for all cameras.**

You have requested a parameter set that is not part of the loaded camera parameter set.

**Optotrak: Failed loading camera parameter file.**

The specified camera parameter file was found but could not be loaded.

**Optotrak: Failed requesting node information.**

Failed sending a node information request to the Optotrak Administrator. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed requesting rigid status.**

Failed sending a rigid body status request to the Optotrak Administrator. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed requesting status.**

Failed sending a status request to the Optotrak Administrator. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed sending activate markers.**

Failed attempting to send the activate markers message to the Optotrak Administrator. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed sending camera parameters.**

Failed attempting to send the camera parameters to the Optotrak Administrator. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed sending collection parameters.**

Failed sending the standard set of collection parameters. There is a communication problem between the host computer and the Optotrak System, or one of the collection parameters was specified with an invalid value.

**Optotrak: Failed sending de-activate markers.**

Failed attempting to send the de-activate markers message to the Optotrak Administrator. There is a communication problem between the host computer and the Optotrak System.

**Optotrak: Failed sending strober port firing table.**

Failed to send the strober table message to the Optotrak. There is either a communication problem between the host computer and the Optotrak System, or an invalid table entry (i.e. an incorrect number of markers).

**Optotrak: Failed sending system parameters.**

Failed to send the system the set of collection parameters. There is either a communication problem between the host computer and the Optotrak System, or one of the system parameters was specified with an invalid value.

**Optotrak: Failed to determine FOR transform.**

The Optotrak System could not determine the required transformation while attempting to change the camera coordinate system.

**Optotrak: Invalid node number specified.**

An attempt was made to obtain node information on a device with an invalid node number.

**Optotrak: Memory block too small for camera parameters.**

The camera parameters in the specified camera parameter file exceed the size of the temporarily allocated memory buffer. This is an internal error message and should never appear. Contact NDI for assistance in resolving this error.

**Optotrak: Out of memory for camera parameters.**

Failed to allocate a temporary memory buffer for the camera parameters. The host computer has insufficient available memory for this operation.

**Optotrak: Strober port table not fully defined.**

Some, but not all of strober ports in the port table were correctly specified in the Strober Table section of the parameter file.

**Optotrak: Unable to find camera parameter file.**

The specified camera parameter file could not be found in any of the Northern Digital search directories.

**Optotrak: Unable to load the collection parameter file.**

The system was unsuccessful in attempting to load a collection parameter file. This error may have occurred because the file could not be found or it could not be opened for reading.

**Optotrak: Unable to load the collection parameter section XXXX.**

The section specified by XXXX could not be found in the collection parameter file. The section might be either missing, or incorrectly labelled.

**Optotrak: Unable to load the strober table from INI file.**

The Strober Table section in the collection parameter file could not be found. The section might be either missing, or incorrectly labelled.

**Optotrak: Unable to save the collection parameter file.**

The system was unsuccessful in its attempt to save a collection parameter file. This may be because the file could not be opened for writing.

**Optotrak: Unknown rigid return format in the collection parameter file.**

The ReturnFormat variable in a RigidBody section of the collection parameter file was not recognized. The three valid formats are Euler, quaternion, and matrix.

## **B.6 ODAU Related Error Messages**

**Odau: Error shutting down collections.**

Failed to successfully send shutdown collection messages to the Optotrak System and ODAU devices. There is a communication problem between the host computer and the Optotrak System.

**Odau: Failed loading ODAU Control FiFo.**

Failed while sending the FiFo control table to the specified ODAU device. There is a communication problem between the host computer and the



Optotrak System. An ODAU I also returns this message since it does not support the control FiFo option.

**Odau: Failed requesting status.**

Failed to send status request to the specified ODAU device. There is a communication problem between the host computer and the Optotrak System.

**Odau: Failed sending ODAU collection parameters.**

Failed sending collection parameters to the specified ODAU device. There is either a communication problem between the host computer and the Optotrak System, or one of the collection parameters was specified with an invalid value.

**Odau: Failed setting analog outputs.**

The ODAU analog outputs could not be set to the requested values. There is a communication problem between the PC and the Optotrak System. An ODAU I also returns this message since it does not support analog outputs.

**Odau: Failed setting digital outputs.**

The ODAU digital outputs could not be set to the requested values. There is a communication problem between the host computer and the Optotrak System.

**Odau: Failed setting ODAU Timer.**

The ODAU timer could not be set with the requested values. There is either a communication problem between the host computer and the Optotrak System, or one of the parameters (such as the mode or timer) was specified with an invalid value. An ODAU I also returns this message since it does not support timer setting.

**Odau: Invalid Gain.**

An invalid gain value was specified to the ODAU. Valid values for the ODAU II are 1, 5, 10, and 100. Valid values for the ODAU I are 1, 2, 4, and 8.

**Odau: Length of FiFo table does not match collection.**

The number of FiFo table entries specified does not match the number specified in the collection.

**Odau: Unable to load the collection parameter section XXXX, using defaults.**

The section specified by XXXX, containing settings for the ODAU, could not be found in the collection parameter file. The section might either be missing or mislabelled.

**Odau: Unable to save the collection parameter section XXXX.**

The section specified by XXXX could not be saved to the collection parameter file. There is a file writing problem on the host computer.

## **B.7 Real-time Related Error Messages**

**DatalsReady: Failed receiving message header.**

Failed to receive the header portion of a message sent to the application program. There is a communication problem between the host computer and the Optotrak System.

**DataReceiveLatestOdauRaw: Wrong type of data received.**

The Optotrak System expected to receive ODAU data, but some other type of data was received instead. Ensure that a request for ODAU data is sent before trying to receive it. Programs should not have multiple data requests outstanding. Each type of request should be followed by the appropriate type of receive.

**Receive: Data buffer was not allocated.**

A required data buffer has not been allocated. This is an internal error that is probably caused either by calling certain API routines in an improper sequence, or there is insufficient memory on the host computer.

**Receive: Failed allocating required buffer for latest transformed data.**

Failed to allocate the buffer space needed for transforming the data. The host computer has insufficient memory available for this operation.

**Receive: Failed receiving analog settings.**

An attempt to obtain the current ODAU analog output channel settings has failed. For ODAU II units, there is a communication problem between the host computer and the Optotrak System. A request made to an ODAU I also returns this message, since modifying output voltages is not supported on the ODAU I.

**Receive: Failed receiving digital settings.**

An attempt to obtain the current ODAU digital output channel states has failed. There is a communication problem between the host computer and the Optotrak System.

**Receive: Failed receiving latest data header.**

Failed to receive the data header portion of a real-time data frame. There is a communications problem between the host computer and the Optotrak System.

**Receive: Failed receiving latest data.**

Failed to receive the data portion of a real-time data frame. There is a communication problem between the host computer and the Optotrak System.

**Receive: Failed receiving latest transformed data.**

Failed to receive the data portion of a 6D data frame. There is a communication problem between the host computer and the Optotrak System.

**Receive: Failed resizing temporary buffer.**

Failed to reallocate needed buffer space. The host computer has insufficient available memory for this operation.

**Receive: Failed while discarding data.**

Failed to receive real-time buffer overflow data. There is a communication problem between the host computer and the Optotrak System.

**Receive: Unknown rigid body type.**

The rigid body type received with the transformation data was invalid. There is probably a communication problem between the host computer and the Optotrak System.

**Receive: Wrong type of data received.**

The Optotrak System was expected to send host computer a certain type of data, but a different type of data was received instead. This error is caused by an application having previously requested, but not received, a different data type.

**Request: Failed requesting latest 3D data.**

A failure occurred while sending a request for the latest frame of real-time 3D data. There is a communication problem between the host computer and the Optotrak System.

**Request: Failed requesting latest ODAU raw data.**

A failure occurred while sending a request for the latest frame of real-time ODAU raw data. There is a communication problem between the host computer and the Optotrak System.

**Request: Failed requesting latest raw data.**

A failure occurred while sending a request for the latest frame of real-time raw data. There is a communication problem between the host computer and the Optotrak System.

**Request: Failed requesting latest transformation data.**

A failure occurred while sending a request for the latest frame of real-time rigid body transformation data. There is either a communication problem between the host computer and the Optotrak System, or the system does not support transformation data. See the [“Real-time Rigid Body Programmer’s Guide”](#) on page 51.

**Request: Failed requesting next 3D data.**

A failure occurred while sending a request for the next frame of real-time 3D data. There is a communication problem between the host computer and the Optotrak System.

**Request: Failed requesting next ODAU raw data.**

A failure occurred while sending a request for the next frame of real-time ODAU raw data. There is a communication problem between the host computer and the Optotrak System.

**Request: Failed requesting next raw data.**

A failure occurred while sending a request for the next frame of real-time raw data. There is a communication problem between the host computer and the Optotrak System.

**Request: Failed requesting next transformation data.**

There was a failure while sending a request for the next frame of real-time rigid body transformation data. There is either a communication problem between the host computer and the Optotrak System, or the system does not support transformation data. See the [“Real-time Rigid Body Programmer’s Guide”](#) on page 51.

## B.8 Data Buffer Spooling Related Error Messages

**DataBuffer: Failed receiving message header.**

Failed to receive the header portion of an incoming message. There is a communication problem between the host computer and the Optotrak System.

**DataBuffer: Failed sending data buffer request.**

Failed while sending a request that buffered data be spooled to the application program. There is a communication problem between the host computer and the Optotrak System.

**DataBuffer: Failed sending start data buffer.**

Failed while sending a message that starts all data buffers spooling data back to the application program. There is a communication problem between the host computer and the Optotrak System.

**DataBuffer: Failed sending stop data buffer.**

Failed while sending the message that stops data buffers spooling data back to the application program. There is a communication problem between the host computer and the Optotrak System.

**DataBuffer: Out of memory for spool buffer.**

Failed to allocate memory for temporary storage that receives buffered data sent to the application program. The host computer has insufficient available memory for this operation.

**DataBuffer: Spool destination already defined for device.**

A spool destination has already been defined for the specified device's buffered data. To change spool destination assignments, first remove the current associations with `DataBufferAbortSpooling`.

**DataBuffer: Unable to open spool file.**

Failed to open the specified file for the buffered data. The file may already exist as a read-only file, or the application may not have write permission for the file destination.

**ReceiveDataBufferBlock: Failed receiving data block.**

Failed to receive the data portion of the spooled buffer data block. There is a communication problem between the host computer and the Optotrak System.

**ReceiveDataBufferBlock: Failed receiving data header.**

Failed to receive the data header portion of the spooled buffer data block. There is a communication problem between the host computer and the Optotrak System.

**ReceiveDataBufferBlock: Failed sending data buffer request.**

Failed while sending a request for the next block of buffered data. There is a communication problem between the host computer and the Optotrak System.

**SpoolDataToFile: Failed receiving data block.**

Failed while receiving buffered data to the temporary application memory buffer. There is a communication problem between the host computer and the Optotrak System.

**SpoolDataToFile: Failed writing data block.**

Failed to write data from the host computer memory buffer to the destination spool file.

**WriteFileHeader: Failed writing file header.**

Failed to write the Northern Digital file header data to the file at the completion of the spooling operation.

## B.9 Rigid Body Related Error Messages

### **RigidBody: Failed sending add normals.**

Failed while sending the message that adds normal vectors to a rigid body that is already in the rigid body tracking list. There is either a communication problem between the host computer and the Optotrak System, or the system does not support transformation data. See the “[Real-time Rigid Body Programmer’s Guide](#)” on page 51.

### **RigidBody: Failed sending add rigid.**

Failed while sending the message that adds a rigid body to the rigid body tracking list. There is either a communication problem between the host computer and the Optotrak System, or the system does not support transformation data. See “[Real-time Rigid Body Programmer’s Guide](#)” on page 51.

### **RigidBody: Failed sending delete rigid body.**

Failed while sending the message that removes a specified rigid body from the rigid body tracking list. There is probably a communication problem between the host computer and the Optotrak System.

### **RigidBody: Failed sending rotate FOR.**

Failed while sending the message that instructs the Optotrak System to express rigid body transformations in the specified coordinate system. There is probably a communication problem between the host computer and the Optotrak System.

### **RigidBody: Failed sending stop rotating FOR.**

Failed while sending the message that stops the Optotrak System from expressing rigid body transformations in a coordinate system other than the default one. There is probably a communication problem between the host computer and the Optotrak System.

### **RigidBody: Failed setting rigid body status.**

Failed while sending the message that changes settings for the specified rigid body. There is probably a communication problem between the host computer and the Optotrak System.

**RigidBody: Out of memory for rigid data.**

Failed to allocate memory to temporarily store the rigid body information. The host computer has insufficient memory available for this operation.

**RigidBody: Rigid body ID out of range.**

The specified rigid body ID exceeds the maximum limit.

**RigidBody: Unable to allocate required memory for transformation.**

Failed to allocate memory required for rigid body transformations. The host computer has insufficient memory available for this operation.

**RigidBody: Unable to find rigid body file.**

Failed to locate the specified rigid body file in any of the Northern Digital search directories.

**RigidBody: Unable to load rigid body file.**

Failed to successfully load the information in the specified rigid body file.

## **B.10 File Processing Related Error Messages**

**FileClose: Unable to write file header.**

Failed to successfully write the Northern Digital file header data to the file prior to closing the file.

**FileConvert: System currently spooling data.**

Tried to convert raw data while the system was spooling buffered data.

**FileConvert: Unable to allocate required memory for conversion.**

Failed to allocate a temporary memory buffer for storage and conversion of raw data. The host computer has insufficient memory available for this operation.

**FileConvert: Unable to allocate required memory for unpacking file data.**

Failed to allocate temporary memory buffers needed for unpacking the file data. The host computer has insufficient memory available for this operation.



**FileConvert: Unable to open input file.**

Failed to open the input raw data file that is to be converted. Ensure that the file exists and that the application has read permission.

**FileConvert: Unable to open output file.**

Failed to open a new file for the converted data. Ensure that the application has write permission. This may also indicate that the file already exists and is read-only.

**FileConvert: Unable to read input file header.**

Failed to read the Northern Digital file header from the input raw data file. Check that the file is a valid Northern Digital Floating Point (NDFP) format file.

**FileConvert: Unable to read raw data.**

Failed to read data from raw data file that is currently being converted.

**FileConvert: Unable to write converted data.**

Failed to write the converted data to the output data file.

**FileConvert: Unable to write data header.**

Failed to write the Northern Digital file header to the output data file.

**FileOpen: File ID already in use.**

Tried to open a file with a file ID that is currently assigned to an opened file.

**FileOpen: Unable to allocate memory for frame buffer.**

Failed to allocate a temporary memory buffer for storage of a frame buffer. The host computer has insufficient memory available for this operation.

**FileOpen: Unable to open file.**

Failed to open a file stream for the specified file. Ensure that the application has the appropriate permissions.

**FileOpen: Unable to read file header.**

Failed to read the Northern Digital file header from the specified file. Ensure that the file is a valid NDFP format file.

**FileOpen: Unable to write file header.**

Failed to write the Northern Digital file header to the specified file.

**FileOpenAll: File ID already in use.**

Tried to open a file with a file ID that is currently assigned to an opened file.

**FileOpenAll: Unable to allocate memory for frame buffer.**

Failed to allocate a temporary memory buffer for storage of a frame buffer. The host computer has insufficient memory available for this operation.

**FileOpenAll: Unable to open file.**

Failed to open a file stream for the specified file. Ensure that the application has the appropriate permissions

**FileOpenAll: Unable to read file header.**

Failed to read the Northern Digital file header from the specified file. Check that the file is a valid NDFP format file.

**FileOpenAll: Unable to write file header.**

Failed to write the Northern Digital file header to the specified file.

**FileRead: File not opened.**

Attempted to read from a file using a file ID that has no open file assigned to it. Before calling FileRead, open the file with FileOpen.

**FileRead: Specified destination pointer is null.**

The destination buffer pointer argument for FileRead has a null value. This error is caused by either incorrectly allocating, or failing to allocate, the required destination buffer in the application program before calling FileRead.

**FileRead: Specified frames are out of range.**

Attempted to read frames from the file that do not exist in the file.

**FileRead: Unable to read specified frames.**

Failed while attempting to read the specified frames from the specified file.

**FileRead: Unable to seek required position.**

Failed to move the file pointer to the start of the specified frame read range.

**FileReadAll: File not opened.**

Attempted to read from a file using a file ID that is not assigned to the file. Before calling FileReadAll, call FileOpenAll.

**FileReadAll: Specified destination char pointer is null.**

The pointer to the buffer for the character subitems read from the file is null. Before calling FileReadAll, the application must allocate a sufficiently large buffer for the character subitems. That application must provide a non-null pointer even if the file contains no character subitems.

**FileReadAll: Specified destination double pointer is null.**

The pointer to the buffer for the double subitems read from the file is null. Before calling FileReadAll, the application must allocate a sufficiently large buffer for the double subitems. The application must provide a non-null pointer even if the file contains no double subitems.

**FileReadAll: Specified destination float pointer is null.**

The pointer to the buffer for the float subitems read from the file is null. Before calling FileReadAll, the application must allocate a sufficiently large buffer for the float subitems. The application must provide a non-null pointer even if the file contains no float subitems.

**FileReadAll: Specified destination int pointer is null.**

The pointer to the buffer for the integer subitems read from the file is null. Before calling FileReadAll, the application must allocate a sufficiently large buffer for the integer subitems. That application must provide a non-null pointer even if the file contains no integer subitems.

**FileReadAll: Specified frames are out of range.**

Attempted to read frames from the file that do not exist in the file.

**FileReadAll: Unable to read specified frames.**

Failed while attempting to read the specified frames from the specified file.

**FileReadAll: Unable to seek required position.**

Failed to move the file pointer to the start of the specified frame read range.

**FileWrite: File not opened.**

Attempted to write to a file using a file ID that has no open file assigned to it. Use FileOpen to open the file before calling FileWrite.

**FileWrite: File not opened for write.**

Attempted to write to a file opened in read-only mode.

**FileWrite: Specified frames are out of range.**

Attempted to write non-existent frames to the file.

**FileWrite: Specified source pointer is null.**

The source buffer pointer argument for FileWrite has a null value. The application program must provide a valid pointer to a buffer containing the data to be written to the file.

**FileWrite: Unable to seek required position.**

Failed to move the file pointer to the start of the specified frame write range.

**FileWrite: Unable to write specified frames.**

Failed while attempting to write the specified frames to the specified file.

**FileWriteAll: File not opened.**

Attempted to write to a file using a file ID that has no open file assigned to it. Use FileOpenAll to open the file before calling FileWriteAll.

**FileWriteAll: File not opened for write.**

Attempted to write to a file opened in read-only mode.

**FileWriteAll: Specified frames are out of range.**

Attempted to write non-existent frames to the file.

**FileWriteAll: Specified source char pointer is null.**

The source buffer pointer argument for FileWriteAll (character type) has a null value. The application must provide a non-null pointer to a data buffer even if the output file contains no character subitems.

**FileWriteAll: Specified source double pointer is null.**

The source buffer pointer argument for FileWriteAll (double type) has a null value. The application must provide a valid pointer to a data buffer even if the output file contains no double subitems.

**FileWriteAll: Specified source float pointer is null.**

The source buffer pointer argument for FileWriteAll (float type) has a null value. The application must provide a valid pointer to a data buffer even if the output file contains no float subitems.

**FileWriteAll: Specified source int pointer is null.**

The source buffer pointer argument for FileWriteAll (integer type) has a null value. The application must provide a valid pointer to a data buffer even if the output file contains no integer subitems.

**FileWriteAll: Unable to seek required position.**

Failed while moving the file pointer to the start of the specified frame write range.

**FileWriteAll: Unable to write specified frames.**

Failed while attempting to write the specified frames to the specified file.



## Appendix C    **Flags and Settings Associated with Rigid Bodies**

This section is a brief description of rigid bodies as they are defined and used with the Optotrak System. Look in the *Optotrak Certus Rigid Body and Tool Design Guide* for more information on constructing rigid bodies and defining the rigid body within the Optotrak System. This section also summarizes information on the flags and error settings used within the API that apply to rigid body calculations.

### C.1    **Rigid Body Concepts and Terms**

A rigid body is an object that has at least three markers that are fixed in place so that they cannot move relative to each other. Rigid bodies require three fixed markers, but accuracy may be improved if more markers are used. Refer to the *Optotrak Certus Rigid Body and Tool Design Guide* and *Technical Bulletin TB0021 - Design and Manufacturing Tools Incorporating IRED Markers* for further information on rigid body design.

Markers can be real or imaginary. Imaginary markers are used to identify points where you do not want to place a marker. For example, it may be necessary to know the location of a tool tip on a robot. Using an optional digitizing probe, the tip of the tool can be defined relative to the remainder of the rigid body.

After a rigid body has been constructed, use the Optotrak System to accurately measure the positions of the markers and record that information in 3D data files. To do this, either several static positions or a slowly rotating dynamic position of the rigid body can be measured. NDI 6D Architect uses the 3D files to average the results to a single frame or view. This data is then used to generate a rigid body file (.rig) and align the rigid body in a user-defined coordinate system.

### C.2    **Accessing the Rigid Body with the API**

The API can be used to

- add rigid bodies to the rigid body tracking list (RigidBodyAdd or RigidBodyAddFromFile)
- delete rigid bodies from the rigid body tracking list (RigidBodyDelete)
- set some of the parameters that affect the accuracy of rigid body measurements (RigidBodyChangeSettings).
- set the type of transformation that will be calculated (RigidBodyChangeSettings)
- change the coordinate system (RigidBodyChangeFOR)

- obtain 6D transformation data (DataGetLatestTransforms, DataGetLatestTransforms2, DataGetNextTransforms, DataGetNextTransforms2, DataReceiveLatestTransforms, DataReceiveLatestTransforms2, RequestLatestTransforms and RequestNextTransforms)
- convert the data from one type of transformation to another (CombineXfrms, CvtQuatToRotationMatrix, CvtRotationMatrixToQuat, DetermineEuler, DetermineR, InverseXfrm, and TransformPoint)

Each of these routines is described in detail in “[Rigid Body Specific Routines](#)” on page 206 and “[Rigid Body Related Routines](#)” on page 214. Examples of rigid body code are located in the “[Real-time Rigid Body Programmer’s Guide](#)” on page 51.

The real-time data retrieval routines *type*Transforms have a prototype of the form:

```
int DataGetLatestTransforms( unsigned int  *puFrameNumber
                           unsigned int  *puElements
                           unsigned int  *puFlags
                           void          *pDataDest )
```

The real-time data retrieval routines *type*Transforms2 have a prototype of the form:

```
int DataGetLatestTransforms2 (unsigned int          *puFrameNumber
                              unsigned int          *puElements
                              unsigned int          *puFlags
                              struct OptotrakRigidStruct *pDataDest6D
                              Position 3D          *pDataDest3D )
```

The structure of pDataDest is shown in [Table 10-3 on page 243](#).

### OptotrakRigidStruct

Rigid body data consists of an array of the structure OptotrakRigidStruct. If the routine is of the *type*Transforms, there will also be an array of position3D data frames. Transformation data, error flags and error settings are all found in the structure OptotrakRigidStruct. This structure, shown in detail in “[Organization of Rigid Body Transformation Data](#)” on page 242, is comprised of a header, the rotation data and the translation data. The 3D marker data that the transformation data was calculated from may follow the 6D data depending on the API routines used to obtain the 6D data.

The c-type structure, available on the API CD in ndopto.h, is:

```
struct OptotrakRigidStruct
{
    long int          RigidId;
    long int          flags;
```



```

float          QuaternionError,
float          IterativeError;
union TransformationUnion transformation;
};

```

### C.3 Flags Affecting Rigid Bodies

Within the API, there are three locations for flags affecting rigid bodies.

Buffering errors are reported with the status flag puFlags in many of the real-time data retrieval routines. The element puFlags reports the current status of the Optotrak System buffer. If a buffering error has occurred, puFlags will be set to a non-zero number.

nFlags are instruction flags that set the way in which transformations are calculated in routines and can be set in either the routine RigidBodyChangeSettings or RigidBodyAdd/RigidBodyAddFromFile.

Lastly, flags that warn of errors in the rigid body transformations are reported within the flags element of OptotrakRigidStruct. These two flags are OPTOTRAK\_UNDETERMINED\_FLAG and OPTOTRAK\_RIGID\_ERR\_MKR\_SPREAD.

#### OPTOTRAK\_UNDETERMINED\_FLAG

The OPTOTRAK\_UNDETERMINED\_FLAG is set in the flags field of the OptotrakRigidStruct of each rigid body for which the Optotrak System was unable to determine the transformation. This is the most important flag to check to determine the quality of rigid body transformations.

---

**Note** You must check the flags element each time transformation data is received to see if the bit flags have been set.

---

#### OPTOTRAK\_RIGID\_ERR\_MKR\_SPREAD

OPTOTRAK\_RIGID\_ERR\_MKR\_SPREAD is set if the spread in the markers is smaller than the values set in the rigid body file. To make the Optotrak System look for errors in minimum marker spread, load the rigid body files with the routine RigidBodyAddFromFile. Three minimum distances (in mm) can be specified: the first is the magnitude of the largest vector found connecting any two visible markers, the second is the largest component of the vectors orthogonal to the first vector, and the third is the longest vector found orthogonal to the first two vectors. If any of the rigid body vectors are smaller than their corresponding spread parameter, then the transformation is considered undetermined, and the OPTOTRAK\_UNDETERMINED\_FLAG and

OPTOTRAK\_RIGID\_ERR\_MKR\_SPREAD bit flags are set. Specifying a value of zero for any of the minimum spread parameters will indicate to the Optotrak System to ignore that parameter in the check. For example, adding the following lines to the input rigid body file:

```
MinSpread1
100

MinSpread2
10

MinSpread3
0
```

will result in the Optotrak System flagging any transformation as undetermined if the longest vector found connecting any two markers is less than 100 mm, or if the largest component of any vector orthogonal to the first is less than 10 mm; the third spread parameter will be ignored.

## C.4 Error Settings for Rigid Bodies

These errors are used to determine the quality of the transformation.

**QuaternionError:** The QuaternionError is the RMS quaternion error for the determined rigid body transformation, if the quaternion algorithm was used.

**Iterative Error:** The IterativeError is the RMS error for the determined rigid body transformation, if the iterative Euler angle algorithm was used.

### C.4.1 Settings Changed in the Routine RigidBodyChangeSettings

There are six elements that define the rejection criteria within this routine. These are:

- nMinMarkers
- nMaxMarkerAngle
- nMax3dError
- fMaxSensorError
- fMax3dRMSError
- fMaxSensorRmsError

See “[RigidBodyChangeSettings](#)” on page 211 for a description of these elements. If the values set in these parameters are exceeded, the transformation will be set to undetermined.

# Index

## Numerics

6D Architect software, 7

## A

Additional Optotrak Manuals, 5  
 amplification, 29, 69, 238, 239  
 API *see* Application Programmer's Interface  
 Application Programmer's Interface, 8

## B

BAD\_FLOAT Constant, 237  
 buffered data  
   collection time, 115, 126, 144, 151  
   error messages, 285  
   from ODAU, 41, 47  
   from Optotrak  
     blocking, 34, 254  
     non-blocking, 36, 254, 255  
   initializing files for, 199  
   initializing memory blocks for, 200  
   on a secondary computer, 83, 89  
   set streaming, 115, 125, 144, 151, 155  
   spooling, 202, 203, 253  
 Buffered Data Retrieval Routines, 198

## C

camera  
   elements, 118, 132  
   errors, 119  
   extended files, 18  
   files, 7, 15, 16, 18, 111, 119  
   marker type, 111  
   marker wavelength, 111  
   parameter file, creating, 117  
   parameter files, 18  
   parameter set, 111, 112, 116  
   parameter sets, 119  
   parameters affecting transformations, 108, 111, 117  
   pointer, 107, 111, 116  
   routines, 18, 106, 107, 111, 118  
   two computers, 19

Centroid, 239  
 centroid data, 29  
   c-type definition, 241  
   routines, 157, 168, 179, 189, 194  
   size calculation, 241  
   variable definition, 2  
 Changing Rigid Body Settings, 54  
 Changing the Rigid Body Coordinate System, 57  
 characterized measurement volume, 11  
 Checking for Undetermined Transforms, 63  
 Client-Server Model, 14  
 COLLECT  
   data files, 66  
   for raw data conversion, 65  
   to build collection parameter files, 128  
 Collecting Buffered Optotrak Data, 34  
 Collection  
   Advanced Method, 36  
   Analog Data, 46  
   ODAU Data, 46  
   Optotrak Data, 34  
   Real Time Data, 29, 43, 51, 84  
   Secondary Host, 89  
   collection parameter  
     fCollectionTime, 126  
     fDutyCycle, 126  
     fFrameFrequency, 125  
     fVoltage, 126  
     nFlags, 126  
     nMarkerFrequency, 125  
     nMarkers, 125, 126  
     nMinimumGain, 125  
     nStreamData, 125  
     nThreshold, 125  
 Collection Parameter File, 128  
 CombineXfrms, 214  
 Connecting Two Host Computers to the Optotrak, 19  
 constants  
   BAD\_FLOAT, 237, 239  
   hexadecimal, 43  
   MAX\_NEGATIVE, 37, 237  
   OPTO\_SYSTEM\_ERROR\_CODE, 273  
   OPTO\_USER\_ERROR\_CODE, 273  
 converting  
   ODAU Raw Data Files, 73  
   Optotrak Raw Data Files, 73  
 coordinate system, 57, 106  
   changing, 210  
   Optotrak, 107

- c-type definition
    - ODAU, 247
    - Optotrak
      - 3D, 240
      - centroid, 241
      - full raw, 242
      - Rigid Body Transformation Data, 240
      - sensor, 240
  - c-type structure
    - OptotrakRigidStruct, 296
  - CvtQuatToRotationMatrix, 215
  - CvtRotationMatrixToQuat, 216
- D**
- data
    - converting raw files, 73
  - Data Buffer Spooling Related Error Messages, 285
  - Data Conversions and Transformations on the Host Computer, 17
  - data format
    - ODAU raw, 246
    - Optotrak 3D, 241
    - Optotrak raw sensor, 238
    - Optotrak rigid body transformation, 242
  - Data Proprietor, 14
  - data, collecting
    - Optotrak, 34
    - real time, 29
    - real time 3D, 29
      - on Secondary Host, 89
    - real time analog, 43
    - real time ODAU, 43
    - real time Optotrak, 29
    - real time rigid body transformations, 51
  - DataBufferAbortSpooling, 198
  - DataBufferGetLatest3D, 30
  - DataBufferInitializeFile, 34, 46, 199
  - DataBufferInitializeMem, 200
  - DataBufferSpoolData, 34, 46
  - DataBufferStart, 36, 202
  - DataBufferStop, 36, 203
  - DataBufferWriteData, 36, 204
  - DataGetLatest3D, 30, 158, 167
  - DataGetLatestOdaRaw, 43, 160
  - DataGetLatestRaw, 157, 161, 168, 171
  - DataGetLatestTransforms, 52, 55, 61, 63, 163
  - DataGetLatestTransforms2, 55, 63, 165
  - DataGetNextOdaRaw, 170
  - DataGetNextTransforms, 172
  - DataGetNextTransforms2, 174
  - DataIsReady, 32, 177
  - DataReceiveLatest3D, 178
  - DataReceiveLatestOdaRaw, 180
  - DataReceiveLatestRaw, 182
  - DataReceiveLatestTransforms, 183
  - DataReceiveLatestTransforms2, 63, 185
  - DataView, 8
  - DetermineEuler, 216, 219
  - DetermineR, 217
  - digital data, 151, 161, 246
    - format, 44
  - digital output, 146
  - digital port, 44, 171
    - configuration, 150
    - operation mode, 143
  - Digitize software, 9
  - duty cycle, 115
- E**
- error code, 29, 69, 239
  - error constants, 273
  - error messages, 273
    - data buffer spooling, 285
    - file processing, 288
    - messagesystem, 274
    - ODAU, 280
    - Optotrak, 277
    - real-time, 282
    - rigid body related, 287
    - transputer, 277
  - Error Settings Within the API, 297
  - euler transformation, 133
    - angle elements, 245
    - combining, 214
    - from rotation matrix, 216
    - inverse transformation, 218
    - rigid body format, 52
    - rotation values, 2
    - to rotation matrix, 217
- F**
- file format
    - Northern Digital Floating Point, 65, 66
    - samples, 68
  - File Processing Related Error Messages, 288
  - File Processing Routines, 220
  - FileClose, 76, 220
  - FileCloseAll, 221
  - FileConvert, 46, 73, 221

- FileOpen, 76, 222
  - FileOpenAll, 224
  - FileRead, 226
  - FileReadAll, 227
  - files
    - camera, 18, 117, 258
    - closing, 221
    - converting, 73, 222
    - converting raw data, 73
    - default, 15, 16, 24
      - camera, 18
    - floating point, 65, 66
    - manipulating, 250
    - naming, 29, 66
    - opening, 220, 223, 225
    - parameter, 128
    - processing, 76
    - processing NDFP data, 76
    - read-write, 77, 226, 227
    - size, calculating, 228, 241, 242, 246, 247
    - type, 29
  - FileWrite, 76, 228
  - FileWriteAll, 230
  - flags
    - low level, for collections, 126
    - marker spread, 63
    - rigid body, 52, 63, 212, 243, 295
    - rigid body error, 63
    - to set file conversion on host, 17
    - to set mode to connect to Optotrak, 102
    - to set processing, 258
    - to set processing options, 120, 253, 259
  - frame size
    - calculating
      - general, 84
      - ODAU, 247
      - Optotrak, 241, 242
      - rigid body, 246
- H**
- hardware types, 113
- I**
- initializing, 15
    - during program runtime, 17
    - from the command line, 17
    - within the API, 16, 17, 24
  - Initializing, Retrieving System Status and Exiting from the Optotrak, 24
  - initiating collection
    - ODAU, 44
    - Optotrak, 30
  - InverseXfrm, 218
  - iterative algorithm, 212
    - setting, 132
    - used with quaternion algorithm, 132
- L**
- LatestTransforms, 183
  - LatestTransforms2, 185
- M**
- marker
    - 3D error, 56, 211
    - 3D RMS error, 54, 211
    - cut off angle, 54, 211
    - data, 3D, 58
    - displaying positions, 30
    - duty cycle, 30, 126, 130
    - error codes, 239
    - frequency, 30, 115, 125, 130
    - missing, 237
    - sensor error, 211
    - spread errors, 63, 84
    - spread parameters, 297
    - type, 111, 118, 132
    - waiting for, 255
    - waiting to come into view, 37
    - wavelength, 112, 119, 132
  - matrix
    - rotation values, 2
  - MAX\_NEGATIVE Constant, 37, 237
  - measurement volume, 11
  - Message System Related Error Messages, 274
  - Missing Marker Constants, 237
  - mixed system capability, 13
- N**
- Naming Conventions Within Routines, 3
  - NDFP, *see* Northern Digital Floating Point
  - NDI PCI interface card, 7
  - NDI ToolBench software, 7
  - network information file, 15, 24, 101, 103, 257, 277
  - node identifiers, 113

nOptotrakAlignSystem, 231  
 nOptotrakCalibrigSystem, 232  
 Northern Digital Floating Point Format, 66

## O

### ODAU

- about, 3
- analog data, 43
- buffered data, 262
  - retrieving, 46
- collection parameter file, 145, 149, 152, 261
- configuration of digital port, 154, 156
- constants identifying, 42
- converted data
  - format, 71
- converting raw data to voltages, 65, 73
- digital data, 43
- digital output, 146
- digital port, 47
  - see also* digital port
- error message, 276, 277, 280–282
- gain, 47
- order of initializing, 42
- raw data
  - data frame, 246
  - format, 70
  - request/receive, 160, 170, 190, 196, 261
- retrieving data, 43, 46
- routines, 143, 143–156
  - listing, 98
- sample files, 250, 261
- setting up collection, 41, 42, 44
- status, 143
- types of, 41
- update masks, 147
- user timer, 148
- voltages, 43
  - see also* voltages

ODAU *see* Optotrak Data Acquisition Unit

OdauGetStatus, 143

OdauSaveCollectionToFile, 144

OdauSetAnalogOutputs, 145

OdauSetDigitalOutputs, 146

OdauSetTimer, 148

OdauSetupCollection, 42, 149

OdauSetupCollectionFromFile, 152

on-host

- conversion
  - enabling, 110

- number of rigid bodies, 52
- enabling rigid body, 120, 129
- OPTO\_CONVERT\_ON\_HOST, 121
- OPTO\_LIB\_POLL\_REAL\_DATA, 121
- OPTO\_RIGID\_ON\_HOST, 121
- OPTO\_SYSTEM\_ERROR\_CODE Constant, 273
- OPTO\_USER\_ERROR\_CODE Constant, 273
- Optoscope Administrator, 14
- Optotrak
  - 3D data, 188, 193, 241, 242
  - 6D data, 197
  - additional manuals, 5
  - buffered data, 34, 36
  - camera parameter file, 116
  - collection time, 155
  - configuring, 29
  - coordinate system, 106
  - data conversion, 17, 65, 66, 73
  - data format, 68
  - data to secondary host, 84
  - data types, 237
    - general size calculation, 241
  - disconnecting API from, 105
  - error messages, 277
  - flags
    - see also* flags
    - to set rigid body calculations, 212
  - initializing, 15, 24, 102, 103
  - low level parameters, 126
  - node identifier, 113
  - number of rigid bodies, 51
  - processors network, 103
  - raw data, 238, 240, 241
  - retrieving data, 29
  - rigid body data, 242, 244
  - routines, 101–133
  - sample programs, 252–257
  - sensor data
    - c-type definition, 240
  - sensor status, 238
  - set up, 124, 128
  - status, 114
  - strober table, 123
  - system configuration, 7
  - tracking list, 209
- Optotrak Data Acquisition Unit, 8
- OptotrakActivateMarkers, 105
- OptotrakChangeCameraFOR, 106
- OptotrakConvertRawTo3D, 61, 108
- OptotrakConvertTransforms, 109
- OptotrakDeActivateMarkers, 110
- OptotrakDeviceHandleEnable, 135
- OptotrakDeviceHandleFree, 135

- OptotrakDeviceHandleGetNumberProperties, 136
  - OptotrakDeviceHandleGetProperties, 137
  - OptotrakDeviceHandleGetProperty, 138
  - OptotrakDeviceHandleSetBeeper, 139
  - OptotrakDeviceHandleSetProperties, 139
  - OptotrakDeviceHandleSetVisibleLED, 140
  - OptotrakGetCameraParameterStatus, 19, 257
  - OptotrakGetDeviceHandles, 141
  - OptotrakGetErrorString, 111, 112, 118, 273
  - OptotrakGetNodeInfo, 113
  - OptotrakGetNumberDeviceHandles, 142
  - OptotrakGetStatus, 84, 114, 273
  - OptotrakLoadCameraParameters, 18, 24, 116
  - OptotrakRigidStruct, 63, 242, 244, 296, 297
  - OptotrakSaveCollectionToFile, 117
  - OptotrakSetCameraParameters, 19, 258
  - OptotrakSetProcessingFlags, 17, 120
  - OptotrakSetStroberPortTable, 123
  - OptotrakSetupCollection, 124
  - OptotrakSetupCollectionFromFile, 128
  - OptotrakStopCollection, 133
- P**
- PCI interface card, 7
  - port pin configuration, 154–156
  - Processing NDFP Data Files, 76
- Q**
- quaternion algorithm, 212
  - quaternion transformation
    - corresponding rotation matrix, 215
    - elements, 245
    - format, 52
    - specifying, 132
    - used with iterative algorithm, 132
    - values, 2
  - quick guide, 20
- R**
- raw data
    - about, 237
    - converting
      - to 3D, 73
    - converting ODAU raw to voltages, 43, 73
    - c-type definition, 240, 247
    - elements, 239
    - error codes, 239
    - error messages, 284
    - file format, 69, 70, 238, 246
    - file name convention, 66
    - file size, 241
    - from ODAU, 43
      - buffered, 46
    - from Optotrak
      - blocking, 29
    - on secondary computer, 84
    - routines, 108, 160, 161, 170, 171, 180, 182, 189, 190, 191, 195, 196, 221
    - specifying
      - buffered, 131
      - file type, 222
      - full raw, 126
  - read-write, 76, 223, 225, 229, 231
  - real-time data
    - ODAU, 43
    - Optotrak, 29
    - rigid body transformations, 51
    - routine based interface, 43
    - secondary host, 84
  - Real-time Data Retrieval Routines, 157
  - Real-time Related Error Messages, 282
  - ReceiveLatestData, 187
  - RequestLatest3D, 32, 188
  - RequestLatestOdaRaw, 190
  - RequestLatestRaw, 189, 191
  - RequestLatestTransforms, 192
  - RequestNext3D, 193
  - RequestNextOdaRaw, 195
  - RequestNextRaw, 196
  - RequestNextTransforms, 197
  - RetrieveSwitchData, 187
  - retrieving
    - camera parameters, 111
    - data
      - buffered ODAU, 46
      - buffered Optotrak, 34
      - on secondary host, 84, 89
      - real-time ODAU, 43
      - real-time Optotrak, 29
      - rigid body real-time, 51
    - Optotrak System Status, 24
  - rigid body
    - 3D to 6D, 109
    - adding, 206, 208
    - adding a rigid body, 51
    - changing settings, 54, 211
    - coordinate system, 57
    - c-type definition, 244
    - data frame size calculation, 246

- data organization, 242
- default settings, 213
- deleting, 213
- error messages, 287
- excluding marker data, 207
- file location, 209
- identification, 52
- iterative algorithm, 132
- marker array, 207
- number of, 115
- programmer's guide, 51
- retrieving data
  - real time, 51
  - transformations, 52, 163, 165, 172, 174, 183, 185, 192, 197
- return format, 133
- routines, 206
  - data retrieval, 163, 165, 172, 174, 183, 185, 192, 197
  - rigid body related, 214–220
  - rigid body specific, 206–214
- sample programs, 258
- settings, 54
- start marker, 132
- transformations, 61, 63
- transformations on-host, 18
- undetermined transforms, 63
- RigidBodyAdd, 58, 206, 243, 260
- RigidBodyAddFromDeviceHandle, 208
- RigidBodyAddFromFile, 51, 58, 208, 243, 297
- RigidBodyChangeFOR, 58, 210
- RigidBodyChangeSettings, 54, 211, 259, 297
- RigidBodyDelete, 213
- rotation matrix, 52, 133
  - convert from quaternion, 215
  - convert to quaternion, 216
  - corresponding to euler, 216, 217
  - corresponding to quaternion, 215
  - c-union, 244
  - file format, 73
  - pointer, 215, 217, 219
  - returning from Optotrak, 212
  - structure member, 244, 246
- Routines
  - CombineXfrms, 214
  - CvtQuatToRotationMatrix, 215
  - CvtRotationMatrixToQuat, 216
  - DataBufferAbortSpooling, 198
  - DataBufferGetLatest3D, 30
  - DataBufferInitializeFile, 199
  - DataBufferInitializeMem, 200
  - DataBufferSpoolData, 201
  - DataBufferStart, 202
  - DataBufferStop, 203
  - DataBufferWriteData, 204
  - DataGetLatest3D, 158, 167
  - DataGetLatestOdauRaw, 160
  - DataGetLatestRaw, 157, 161, 168, 171
  - DataGetLatestTransforms, 163
  - DataGetLatestTransforms2, 165
  - DataGetNextOdauRaw, 170
  - DataGetNextTransforms, 172
  - DataGetNextTransforms2, 174
  - DataIsReady, 177
  - DataReceiveLatest3D, 178
  - DataReceiveLatestOdauRaw, 180
  - DataReceiveLatestRaw, 179, 182
  - DataReceiveLatestTransforms, 183
  - DataReceiveLatestTransforms2, 185
  - DetermineEuler, 216, 219
  - DetermineR, 217
  - FileClose, 220
  - FileCloseAll, 221
  - FileConvert, 221
  - FileOpen, 222
  - FileOpenAll, 224
  - FileRead, 226
  - FileReadAll, 227
  - FileWrite, 228
  - FileWriteAll, 230
  - InverseXfrm, 218
  - OdauGetStatus, 143
  - OdauSaveCollectionToFile, 144
  - OdauSetAnalogOutputs, 145
  - OdauSetDigitalOutputs, 146
  - OdauSetTimer, 148
  - OdauSetupCollection, 149
  - OdauSetupCollectionFromFile, 152
  - OptotrakActivateMarkers, 105
  - OptotrakChangeCameraFOR, 106
  - OptotrakConvertRawTo3D, 108
  - OptotrakConvertTransforms, 109
  - OptotrakDeActivateMarkers, 110
  - OptotrakDeviceHandleEnable, 135
  - OptotrakDeviceHandleFree, 135
  - OptotrakDeviceHandleGetNumberProperties, 136
  - OptotrakDeviceHandleGetProperties, 137
  - OptotrakDeviceHandleGetProperty, 138
  - OptotrakDeviceHandleSetBeeper, 139
  - OptotrakDeviceHandleSetProperties, 139
  - OptotrakDeviceHandleSetVisibleLED, 140
  - OptotrakGetDeviceHandles, 141
  - OptotrakGetErrorString, 111, 112, 118
  - OptotrakGetNodeInfo, 113
  - OptotrakGetNumberDeviceHandles, 142
  - OptotrakGetStatus, 114
  - OptotrakLoadCameraParameters, 116



OptotrakSaveCollectionToFile, 117  
 OptotrakSetProcessingFlags, 120  
 OptotrakSetStrobePortTable, 123  
 OptotrakSetupCollection, 124  
 OptotrakSetupCollectionFromFile, 128  
 OptotrakStopCollection, 133  
 ReceiveLatestData, 187  
 RequestLatest3D, 188  
 RequestLatestOdaRaw, 190  
 RequestLatestRaw, 189, 191  
 RequestLatestTransforms, 192  
 RequestNext3D, 193  
 RequestNextOdaRaw, 195  
 RequestNextRaw, 196  
 RequestNextTransforms, 197  
 RigidBodyAdd, 206  
 RigidBodyAddFromDeviceHandle, 208  
 RigidBodyAddFromFile, 208  
 RigidBodyChangeFOR, 210  
 RigidBodyChangeSettings, 211  
 RigidBodyDelete, 213  
 TransputerInitializeSystem, 102  
 TransputerLoadSystem, 101, 103  
 TransputerShutdownSystem, 105

## S

### Sample Programs, 252–272

Certus Sample 1, 266  
 Certus Sample 11, 268  
 Certus Sample 12, 268  
 Certus Sample 13, 269  
 Certus Sample 14, 271  
 Certus Sample 15, 269  
 Certus Sample 16, 270  
 Certus Sample 17, 270  
 Certus Sample 18, 271  
 Certus Sample 19, 272  
 Certus Sample 2, 266  
 Certus Sample 3, 267  
 Certus Sample 4, 267  
 File Processing Specific, 262  
 ODAU Specific, 261  
 Optotrak, 252  
 Rigid Body Specific, 258  
 Sample 1, 252  
 Sample 10, 259  
 Sample 11, 260  
 Sample 12, 261  
 Sample 13, 261  
 Sample 14, 262  
 Sample 15, 263

Sample 16, 263  
 Sample 17, 264  
 Sample 18, 256  
 Sample 19, 260  
 Sample 2, 253  
 Sample 3, 253  
 Sample 4, 254  
 Sample 5, 254  
 Sample 6, 255  
 Sample 7, 255  
 Sample 8, 256  
 Sample 9, 258  
     Secondary Host Specific, 263  
 scan rates, 151, 154  
 SCSI, 8  
 SCU *see* System Control Unit  
 secondary host  
     collecting real time data, 84  
     connecting, 19  
     order of operation, 83  
     raw data, 84  
     removing from system, 84  
     retrieving data from, 83, 84, 89  
     sample programs, 263  
     specifying a computer as, 103  
 Setting Up Data Collection from the ODAU, 42  
 signal strength, 29, 69, 238, 239, 240  
 sleep routine, 24, 104, 127, 207, 209, 274  
 software  
     6D Architect, 7  
     Digitize, 9  
     NDI ToolBench, 7  
     Optotrak Application Programmer's Interface, 8  
 spool mapping  
     ODAU, 46  
     Optotrak, 34  
 status  
     Optotrak System, 24  
 strobers  
     error message, 279, 280  
     port-firing table, 123, 129, 279, 280  
 symbols, definitions of, 1  
 system control unit, 114, 126, 131  
     time control signal, 240  
 system status, 24, 89, 252, 264

## T

TLD files, 276  
 ToolBench software, 7

- transformations, 51, 259
    - changing coordinate system, 57, 58, 210, 260
    - checking, 63
    - combining, 214
    - converting, 109
    - error messages, 287, 288
    - on host computer, 17, 18, 120
    - quaternion, 55
    - raw to 3D, 17
    - retrieving, 163, 165, 173, 174, 184, 186, 192, 197, 259
  - Transforming Previously Obtained Data, 61
  - Transforms, 172, 192
  - Transforms2, 165, 174
  - transputer, 12, 17
    - loading programs, 25, 84, 104
  - Transputer Related Error Messages, 277
  - TransputerDetermineSystemCfg, 16, 257
  - TransputerInitializeSystem, 24, 84, 102, 273
  - TransputerLoadSystem, 16, 24, 101, 103
  - TransputerShutdownSystem, 25, 84, 105, 277
  - two computers
    - sample programs, 250
    - where to locate programs, 19, 83
- V**
- variables, definitions of, 2
  - voltages
    - calculating, 43
    - converting raw data to, 65, 73
    - data format, 71
    - retrieving, 146