

# Panda Programming Guide

Ahmad AlAttar

March 12, 2019

## Table of Content

<b>1</b>	<b>Introduction</b>	<b>3</b>
<b>2</b>	<b>Installation Guide</b>	<b>3</b>
2.1	Workspace Setup . . . . .	3
2.2	Installation of libfranka . . . . .	4
2.3	Installation of franka_ros . . . . .	4
2.4	Installation of panda.moveIt_config . . . . .	5
<b>3</b>	<b>Turn On/Off Panda</b>	<b>5</b>
3.1	Turn On Panda . . . . .	5
3.2	Turn Off Panda . . . . .	5
<b>4</b>	<b>Panda Gripper Control</b>	<b>6</b>
4.1	Launch Gripper Node . . . . .	6
4.2	Gripper Action Servers . . . . .	6
4.2.1	Gripper Homing Action . . . . .	6
4.2.2	Gripper Move Action . . . . .	7
4.2.3	Gripper Stop Action . . . . .	7
4.3	Gripper Control using MoveIt . . . . .	8
4.3.1	Launch MoveIt . . . . .	8
4.3.2	Launch RViz . . . . .	8
4.3.3	Setup . . . . .	8
4.3.4	Gripper Control . . . . .	9
<b>5</b>	<b>Franka Control</b>	<b>10</b>
5.1	Launch Nodes . . . . .	10
5.1.1	Launch Control Node . . . . .	10
5.1.2	Launch MoveIt . . . . .	10
5.1.3	Launch RViz . . . . .	10
5.2	Panda Control using MoveIt . . . . .	11
5.2.1	Setup . . . . .	11

5.2.2	Basic Information Print . . . . .	12
5.2.3	Simple Move . . . . .	12
5.2.4	Planning to Pose Goal . . . . .	12

---

## 1 Introduction

Franka Emika's Panda robot is a 7 DOF manipulator arm, shown below. The Panda Programming Guide will help you install the required packages, turn On/Off the robot, and how to program it using Python.



Figure 1: Franka Emika's Panda robot

---

## 2 Installation Guide

If a workspace exists with `libfranka`, `franka_ros`, and `panda_moveit_config` installed, **then skip the Installation Guide.**

---

### 2.1 Workspace Setup

In order to create ROS packages later on, we first need a catkin workspace. Follow these steps:

1. Go to your home directory:

terminal: `cd`

2. Create a directory called `Franka_ws`. This will be your workspace:

terminal: `mkdir Franka_ws`

3. Go inside the workspace:

terminal: `cd Franka_ws`

4. Create a folder called `src`:

terminal: `mkdir src`

5. Build your workspace:

terminal: `catkin_make`

---

## 2.2 Installation of libfranka

The library libfranka is the C++ implementation of the client side of the Franka Control Interface (FCI). Follow these steps to install libfranka:

1. Go into the src folder of your workspace:

terminal: `cd Franka_ws/src`

2. Execute this command to install required dependencies:

terminal: `sudo apt install build-essential cmake git libpoco-dev  
libeigen3-dev`

3. Clone libfranka with this command:

terminal: `git clone --recursive https://github.com/frankaemika/libfranka`

4. Go back to the root of your workspace:

terminal: `cd ..`

5. Build your workspace:

terminal: `catkin_make`

---

## 2.3 Installation of franka\_ros

The franka\_ros metapackage integrates libfranka into ROS. Follow these steps to install libfranka:

1. Go into the src folder of your workspace:

terminal: `cd Franka_ws/src`

2. Clone libfranka with this command:

terminal: `git clone --recursive  
https://github.com/frankaemika/franka_ros`

3. Go back to the root of your workspace:

terminal: `cd ..`

4. Build your workspace:

terminal: `catkin_make`

---

## 2.4 Installation of panda\_moveit\_config

The panda\_moveit\_config enabled the use of MoveIt to do motion planning for the Panda robot. Follow these steps to install panda\_moveit\_config:

1. Go into the src folder of your workspace:

terminal: `cd Franka_ws/src`

2. Clone panda\_moveit\_config with this command:

terminal: `git clone -b kinetic-devel  
https://github.com/ros-planning/panda_moveit_config.git`

3. Go back to the root of your workspace:

terminal: `cd ..`

4. Build your workspace:

terminal: `catkin_make`

---

## 3 Turn On/Off Panda

This section details how to turn On/Off the Panda robot.

---

### 3.1 Turn On Panda

Before programming, it is important to have the Panda robot in the Ready mode. Follow these steps to do so:

1. Connect all the wires needed.
2. Turn on the controller (you will see a blinking yellow light).
3. Open a browser.
4. Go to 192.168.0.88 using the browser.
5. Unlock the joint brakes (when done, you will see white or purple light).
6. If not already, set the E-Stop to Ready mode (you will see blue light).

---

### 3.2 Turn Off Panda

After you are done, make sure to switch off the Panda following these steps. **DO NOT DIRECTLY SWITCH OFF THE CONTROLLER.**

1. In the web interface click Shut down.
2. Wait for shutdown to complete (when done, panda lights will switch off).
3. Switch off the controller.

---

## 4 Panda Gripper Control

This section shows how to control the Panda robot gripper using action servers and MoveIt.

---

### 4.1 Launch Gripper Node

To be able to control the gripper, we launch the `franka_gripper` node (**NOT NEEDED IF `franka_control` IS ALREADY LAUNCHED**). Follow these steps:

1. Go to the workspace:

```
terminal 1: cd Franka_ws
```

2. Launch the `franka_gripper` node in a terminal:

```
terminal 1: roslaunch franka_gripper franka_gripper.launch
           robot_ip:=192.168.0.88
```

3. In another terminal, you can check that the node was launched:

```
terminal 2: rosnodetool list
```

---

### 4.2 Gripper Action Servers

After launching the `franka_gripper` node, these action servers will be available to you:

1. `MoveAction(width, speed)`
2. `GraspAction(width, epsilon_inner, epsilon_outer, speed, force)`
3. `HomingAction()`
4. `StopAction()`

All of the action servers can be used in Python. To send commands to the action servers, we need action clients.

#### 4.2.1 Gripper Homing Action

The following lines of code show how to create a simple action client for the `HomingAction` action server:

```

import rospy
import actionlib
from franka_gripper.msg import HomingAction, HomingGoal

if __name__ == '__main__':
    rospy.init_node('Franka_gripper_homing_action')
    client = actionlib.SimpleActionClient('/franka_gripper/homing', HomingAction)
    client.wait_for_server()
    goal = HomingGoal()
    client.send_goal(goal)
    client.wait_for_result(rospy.Duration.from_sec(5.0))

```

This action server homes the gripper and updates the maximum width given the mounted fingers.

#### 4.2.2 Gripper Move Action

The following lines of code show how to create a simple action client for the MoveAction action server:

```

import rospy
import actionlib
from franka_gripper.msg import MoveGoal, MoveAction

if __name__ == '__main__':
    rospy.init_node('Franka_gripper_move_action')
    client = actionlib.SimpleActionClient('/franka_gripper/move', MoveAction)
    client.wait_for_server()
    goal = MoveGoal(width = 0.08, speed = 0.08)
    client.send_goal(goal)
    client.wait_for_result(rospy.Duration.from_sec(5.0))

```

This action server moves to a target width with the defined speed.

#### 4.2.3 Gripper Stop Action

The following lines of code show how to create a simple action client for the StopAction action server:

```

import rospy
import actionlib
from franka_gripper.msg import StopAction, StopGoal

if __name__ == '__main__':
    rospy.init_node('Franka_gripper_stop_action')
    client = actionlib.SimpleActionClient('/franka_gripper/stop', StopAction)
    client.wait_for_server()

```

```
action = StopGoal()
client.send_goal(action)
client.wait_for_result(rospy.Duration.from_sec(5.0))
```

This action server aborts a running action. This can be used to stop applying forces after grasping.

---

### 4.3 Gripper Control using MoveIt

**MAKE SURE EITHER `franka_control` NODE OR `franka_gripper` NODE IS LAUNCHED.** Next we will need to launch MoveIt and RViz (RViz is optional).

#### 4.3.1 Launch MoveIt

Follow these steps to launch MoveIt:

1. Go to the workspace:

```
terminal 2: cd Franka_ws
```

2. Launch the MoveIt in a terminal:

```
terminal 2: roslaunch panda_moveit_config panda_moveit.launch
            controller:=position
```

#### 4.3.2 Launch RViz

RViz is used for visualization. Follow these steps to launch RViz (optional):

1. Go to the workspace:

```
terminal 3: cd Franka_ws
```

2. Launch the MoveIt in a terminal:

```
terminal 3: roslaunch panda_moveit_config moveit_rviz.launch
```

#### 4.3.3 Setup

Follow these setup steps (in a new Python file):

1. Import the following libraries:

```
import sys
import rospy
import moveit_commander
import moveit_msgs.msg
import geometry_msgs.msg
```



2. Initialize `moveit_commander` and `rospy`:

```
moveit_commander.roscpp_initialize(sys.argv)
rospy.init_node('move_group_python_interface', anonymous=True)
```

3. Instantiate a `RobotCommander` object which is the outer-level interface to the robot:

```
robot = moveit_commander.RobotCommander()
```

4. Instantiate a `PlanningSceneInterface` object which is an interface to the world surrounding the robot:

```
scene = moveit_commander.PlanningSceneInterface()
```

5. Instantiate a `MoveGroupCommander` object which is an interface to one group of joints:

```
group_name = "hand"
group = moveit_commander.MoveGroupCommander(group_name)
```

6. Create a `DisplayTrajectory` publisher which is used to publish trajectories for RViz to visualize:

```
display_trajectory_publisher =
rospy.Publisher('/move_group/display_planned_path',
moveit_msgs.msg.DisplayTrajectory, queue_size=20)
```

#### 4.3.4 Gripper Control

The following lines of code uses `MoveIt` to open the gripper:

```
joint_goal = group.get_current_joint_values()
joint_goal[0] = 0.03
joint_goal[1] = 0.03
```

```
group.go(joint_goal, wait=True)
group.stop()
```

The following lines of code uses `MoveIt` to close the gripper:

```
joint_goal = group.get_current_joint_values()
joint_goal[0] = 0.00
joint_goal[1] = 0.00
```

```
group.go(joint_goal, wait=True)
group.stop()
```

---

## 5 Franka Control

In this section, we will control the Panda arm using MoveIt in Python. It is important to first launch the required files: `franka_control`, `MoveIt`, and `RViz` (`RViz` is optional).

---

### 5.1 Launch Nodes

#### 5.1.1 Launch Control Node

The `franka_control` node exposes ROS services for controlling the Panda robot. Follow these steps to launch the `franka_control` node:

1. Go to the workspace:

```
terminal 1: cd Franka_ws
```

2. Launch the `franka_control` node in a terminal:

```
terminal 1: roslaunch franka_control franka_control.launch  
robot_ip:=192.168.0.88 load_gripper:=true
```

3. In another terminal, you can check that the node was launched:

```
terminal 2: roslaunch list
```

#### 5.1.2 Launch MoveIt

Follow these steps to launch `MoveIt`:

1. Go to the workspace:

```
terminal 2: cd Franka_ws
```

2. Launch the `MoveIt` in a terminal:

```
terminal 2: roslaunch panda_moveit_config panda_moveit.launch  
controller:=position
```

#### 5.1.3 Launch RViz

`RViz` is used for visualization. Follow these steps to launch `RViz` (optional):

1. Go to the workspace:

```
terminal 3: cd Franka_ws
```

2. Launch the `MoveIt` in a terminal:

```
terminal 3: roslaunch panda_moveit_config moveit_rviz.launch
```

---

## 5.2 Panda Control using MoveIt

This is where the robot actually move! **MAKE SURE YOU LAUNCHED franka\_control, MoveIt, AND RViz (optional)**. Run your code in terminal 4.

### 5.2.1 Setup

Follow these setup steps (in a new Python file):

1. Import the following libraries:

```
import sys
import rospy
import moveit_commander
import moveit_msgs.msg
import geometry_msgs.msg
from math import pi
```

2. Initialize moveit\_commander and rospy:

```
moveit_commander.roscpp_initialize(sys.argv)
rospy.init_node('move_group_python_interface', anonymous=True)
```

3. Instantiate a RobotCommander object which is the outer-level interface to the robot:

```
robot = moveit_commander.RobotCommander()
```

4. Instantiate a PlanningSceneInterface object which is an interface to the world surrounding the robot:

```
scene = moveit_commander.PlanningSceneInterface()
```

5. Instantiate a MoveGroupCommander object which is an interface to one group of joints:

```
group_name = "panda_arm"
group = moveit_commander.MoveGroupCommander(group_name)
```

6. Create a DisplayTrajectory publisher which is used to publish trajectories for RViz to visualize:

```
display_trajectory_publisher =
rospy.Publisher('/move_group/display_planned_path',
moveit_msgs.msg.DisplayTrajectory, queue_size=20)
```

### 5.2.2 Basic Information Print

To get some basic information, you can run the following:

```
# To print the reference frame of robot:
planning_frame = group.get_planning_frame()
print planning_frame

# To print the name of the end-effector link for this group:
eef_link = group.get_end_effector_link()
print eef_link

# To print a list of all the groups in the robot:
group_names = robot.get_group_names()
print group_names

# To print the entire state of the robot:
robot_state = robot.get_current_state()
print robot_state
```

### 5.2.3 Simple Move

To move the robot, run the following piece of code:

```
joint_goal = group.get_current_joint_values()
joint_goal[0] = 0
joint_goal[1] = -pi/4
joint_goal[2] = 0
joint_goal[3] = -pi/2
joint_goal[4] = 0
joint_goal[5] = pi/3
joint_goal[6] = 0

group.go(joint_goal, wait=True)
group.stop() # makes sure no residual movements left
```

### 5.2.4 Planning to Pose Goal

To plan a simple pose and move the robot, run the following piece of code:

```
pose_goal = geometry_msgs.msg.Pose()
pose_goal.orientation.w = 1.0
pose_goal.position.x = 0.4
pose_goal.position.y = 0.1
pose_goal.position.z = 0.4
group.set_pose_target(pose_goal)

plan = group.go(wait=True)
```

```
group.stop()  
group.clear_pose_targets() # clear targets after planning
```