

An Open-Source Plate reader

SOFTWARE USER GUIDE

Contents

Getting Started	S1
Requirements	S1
Installation.....	S1
Windows/Mac.....	S1
Raspberry-Pi.....	S3
Repository Structure Guide	S4
Hardware	S4
Software.....	S4
Before Running a Protocol.....	S5
Initial Function Testing.....	S5
Testing Light Sources	S5
Testing Linear Actuators	S5
Calibration Guide.....	S6
Loading a Plate	S6
Software Structure.....	S7
Structure Overview	S7
Structure Guide.....	S8
System.py.....	S8
GUI.py	S9
PlateConfiguration.py	S10
Protocol.py.....	S10
Absorbance.py	S11
Fluorescence.py	S11
Auxiliary.py.....	S12
Shake.py	S12
Kinetic.py.....	S12
Machine.py	S13
Arduino.py.....	S13
Spec.py.....	S14
System Level Control.....	S14
GUI Level Control.....	S15
Running the GUI.....	S15

The Intro Screen.....	S15
The Settings Screen.....	S16
The Plate/Well Selection.....	S17
The Protocol Selection	S18
Kinetic Toggling.....	S19
Absorbance Settings	S20
Fluorescence Settings	S20
Auxiliary Settings.....	S21
Review Protocol Sequence	S22

Getting Started

Requirements

The OSP software is written in Python 2.7 and utilizes the PyQt4 package to create a functional user interface. Below is a list of all the python packages required to successfully use the OSP device.

- PyQt 4
- Seabreeze
- Numpy
- Pickle
- Platform
- Xlsxwriter
- Matplotlib
- Time
- Serial
- Warnings

In addition, the Arduino IDE will have to be installed in order to load the OSP Arduino script (provided in this repository) on the device's Arduino UNO board.

The installation section below will guide you through installing all python packages and any programs required.

Installation

Windows/Mac

Installation of OSP software requirements on any Windows/Mac device is accomplished utilizing Anaconda, which should come prepackaged with all requirements except PyQt4, Seabreeze, Pyserial and Git. The following instructions require you to write in command line. If you have never used command line visit the following [LINK \(Windows\)](#) or [LINK \(Mac\)](#) for a quick intro.

Step 1. Install Anaconda – Visit the following [LINK](#) to download the Anaconda software. Be sure to download the version for Python 2.X and your specific OS environment (Windows/Mac). When prompted whether to **add Anaconda to your path** select yes.

Step 2. Install Git– Visit the following [LINK](#) to download the Git software. When prompted whether to **add Git to your path** select yes.

Step 3. Clone the OSP repository – Open the Terminal on your computer as an administrator. To do this, find the Terminal application, *right-click* it and select *Run as administrator*. To clone the repository onto your system, type the following lines of code (Make sure to use the ones specific to your OS):

Windows

```
cd %USERPROFILE%  
git clone https://github.com/brianchowlab/OSP.git
```

MAC

```
cd ~  
git clone https://github.com/brianchowlab/OSP.git
```

Step 4. Run the OSP installation script – In the open Terminal window type the following lines of code and leave the window open:

Windows

```
.\OSP\Software\install_pc_reqs.sh
```

MAC

```
./OSP/Software/install_pc_reqs.sh
```

Step 5. Install Arduino IDE – Visit the following [LINK](#) to download the latest Arduino IDE.

Step 6. Upload OSP script to Arduino UNO – Using the Arduino IDE, open up the OSP Arduino script located in the cloned repository folder *OSP\Software\Arduino Files\OSP_Serial_Communication*. Make sure the Arduino is plugged into the USB hub in the OSP device and that the OSP device is connected to your computer via USB. Click the **UPLOAD** button in the Arduino IDE to compile and upload the code to the Arduino UNO.

Step 7. (ONLY WINDOWS USERS) Setting up Seabreeze Drivers – Visit the following [LINK](#) and click the download button. Extract the .zip file to a known location. Open *Device Manager* on your Windows machine. Connect your computer to the OSP device via USB. In the *Device Manager's* list of devices make sure that there is now a tab named *Ocean Optics USB Devices*. In that tab, there should be a device named *Ocean Optics STS (WinUSB)*. If this device is present you can now move on to step 9.

If for some reason, you cannot find *Ocean Optics USB Devices* in the *Device Manager* you will have to manually set-up the drivers. To do this, in the *Device Manager* find the tab labeled *Other Devices*. In that tab there should be a device named *STS*. Right-click on this device name and select *Update drivers* from the drop-down menu. Next, a window will pop up asking you where the system should look for the drivers. Select the *Browse my computer for the driver software* option. On the next screen click the

browse button and direct the system to the *windows-driver-files* folder which you downloaded & extracted at the beginning of this step.

Step 8. Opening the user interface – In order to start utilizing the OSP device, make sure the device is **turned on** and plugged into the computer via USB. In your open Terminal type the following lines of code:

Windows

```
activate osp  
cd "OSP\Software\Python Files"  
python GUI.py
```

MAC

```
source activate osp  
cd "OSP/Software/Python Files"  
python GUI.py
```

Raspberry-Pi

The simplest way to install the OSP software on a Raspberry-Pi is downloading and installing the OSP.img file provided. In order to do this, you will need to have a completely empty 16GB Micro SD Card. If the card is not empty, all of its contents will be deleted in the installation of the OSP.img file. In addition, the first 3 installation steps need to be performed on a computer/laptop with an SD card slot.

Step 1. Download & Unzip the OSP.img – Visit the following [LINK](#) and click the **DISK IMAGE FILE** link and download button. Once downloaded, unzip the file to a known location.

Step 2. Download & Install the Etcher program – Visit the following [LINK](#) and download the *Etcher* program for your specific operating system.

Step 3. Install the OSP.img onto the SD card – Insert your empty SD card into the computer. Start up the *Etcher* software and follow the on-screen instructions in order to install the OSP.img onto the card. **This process can take up to 10 minutes.**

Step 4. Update the OSP repository – Once the OSP.img file is installed, insert it into your Raspberry-Pi (which should be connected to the **turned on** OSP device, touch-screen/monitor, keyboard & mouse). Turn the Raspberry-Pi on and when it is booted open the Terminal. Type the following and close the window when it is finished:

```
cd OSP  
git pull origin master
```

Step 5. Running the User Interface – To activate the user interface, make sure the OSP device is turned on and connected to the Raspberry-Pi via USB. Open up the Terminal and type the following:

```
cd "OSP/Software/Python Files"  
sudo python GUI.py
```

Repository Structure Guide

Below you will find a detailed explanation of how the repository is structured, including folder and file descriptions.

Hardware

All the files required for assembly of the OSP device are in this folder. The full *Assembly Guide* and *Parts List* can be found here as a PDFs.

CAD Files

All files necessary to 3D-print custom parts for the OSP device are found here. Each part is labeled according to its part number in the *Parts List*

Laser-Cut Files

All files necessary to laser-cut custom parts for the OSP device are found here. Each part is labeled according to its part number in the *Parts List*

PCB Files

All files necessary to order custom-made PCB board for the OSP device are found here. Each part is labeled according to its part number in the *Parts List*.

Software

All the files for installation and execution of the OSP software are in this folder. The full *Software Guide* can be found here as a PDF.

Arduino Files

The main Arduino script is found in this folder.

Graphical Designer Files

All raw GUI Designer files are found in this folder. GUI graphics are derived from these files.

Python Files

Data

This directory is where all saved data from experiments are saved in the form of .xlsx files.

Images

This directory contains all the image files referenced in the GUI.

[Protocols](#)

This directory is where all protocol sequences are saved from the GUI

[Screens](#)

This directory contains the Python equivalent of the graphical files found in *Software/Graphical Files*

[System Settings](#)

This directory contains the setting files such as calibrated well positions, excitation LED wavelengths, and the number of scans to average during a measurement.

[All Other Files](#)

These rest of the files are Python scripts containing the class objects from which the GUI is accessed. An explanation of each of these classes is provided in the next section.

[Before Running a Protocol](#)

[Initial Function Testing](#)

Once the OSP device is fully built and connected to a computer or Raspberry-Pi, it is useful to test the functionality of the linear actuators and light sources. To do so, open up the Arduino IDE, click on the Tool tab in the menu bar, and select the Serial Monitor from the drop-down menu. Follow the instructions below and enter the commands to test the functionality of individual components:

[Testing Light Sources](#)

To test LED 1 type: **L1;**

This should result in the LED 1 turning on. To turn it off, type the same command in again.

To test LED 2 type: **L2;**

This should result in the LED 2 turning on. To turn it off, type the same command in again.

To test LED 3 type: **L3;**

This should result in the LED 3 turning on. To turn it off, type the same command in again.

To test LED 4 type: **L4;**

This should result in the LED 4 turning on. To turn it off, type the same command in again.

To test LED 5 type: **L5;**

This should result in LED 5 & LED 6 turning on. To turn them off, type the same command in again.

[Testing Linear Actuators](#)

To extend Linear Actuators type: **M1500.001500.00;**

This should result in both linear actuators extending to about half their limit.

To contract the linear actuators type: **M1000.001000.00;**

This should result in both linear actuators fully contracting.

Calibration Guide

The software provided in this repository comes with pre-loaded calibrated positions for a 24-well plate and a 96-well plate. However, it is recommended that an initial calibration is performed after the OSP device is built to account for any small difference

Before performing any sort of calibration, you will need to prepare the calibration lid for the specific plate-type that is to be calibrated. Refer to the *Calibration Lid Assembly* at the end of the *General Assembly Guide* for instructions on how to do this.

The following steps should be followed every time a calibration is performed:

Step 1. Attach fiber to photodiode – Detach the fiber from the STS spectrophotometer and screw it into the photodiode mount.


Step 2. Adjust the iris size – As part of the *Top Optics Assembly* there is an iris. Using the markings on the rim of the iris, set the size of it to 2.

Step 3. Insert plate into OSP device – Place the calibration lid onto an empty micro-well plate and insert it into the device (refer to the next section about proper orientation).

Step 4. Initialize calibration – On the GUI press the *Calibrate* button and when prompted input the size of the micro-well plate (24 or 96). The machine will begin the calibration process and alert you when it is finished. **This process takes around 3 minutes.**

Loading a Plate

The orientation in which you load the plate is important not only for the calibration process but also for any sort of measurement performed on the OSP device. Please refer to the below picture to ensure you are loading the plate correctly. The plate should be held firmly in place with the attached springs.




Software Structure

The OSP software is object oriented and composed of a total of twelve class objects. The below sections will explain the structure of the software and how information is passed between the objects.

Structure Overview

The diagram below depicts the general flow of information between the 6 main software objects. The entire structure can be divided into three separate categories based on functionality.



Software workflow . (a) Simplified information flow between the device user and individual device components through the Raspberry Pi microcomputer. (b) Python-based API structure based on the information flow. The user can program the instrument through the GUI or a direct script.

System.py and *GUI.py* are the primary classes that allow the user to interface with the plate reader. They track and store user input. *System.py* is the object that contains all the vital functionality of translating the user input into commands to change device settings. This object is responsible for loading-in user settings (such as calibrated well coordinates and excitation led wavelengths) and for preparing the system for execution of input protocols. The user can call this object to set-up and execute protocol sequences from a Python script (more information about this and examples in *System Level Control*). Alternatively, the user can control the device via the provided GUI (*GUI.py*) object. The *GUI.py* object (the visual user interface) has access to all the functions and variables of the *System.py* object.

User input (into either the *GUI.py* or *System.py* constructs) results in the creation of a *PlateConfiguration.py* object containing a list of *Protocol.py* objects. This combination of classes organizes the desired sequence of protocols to be executed. The *PlateConfiguration.py* object contains a list of the selected wells and a separate list of the protocols to be run. As mentioned before, the list of protocols is a list of protocol objects (e.g. *Absorbance.py* is an object for an absorbance protocol and it inherits its attributes from the *Protocol.py* object). Each protocol object contains all the necessary parameters (spectrophotometer settings, excitation LED, acquisition type, etc.) for it to be executed.

The user is not limited to one *PlateConfiguration.py* structure. Performing separate protocols on different sets of wells would simply require multiple *PlateConfiguration.py* objects to be created. For example, to measure fluorescence from the set of wells [A1, A2, A3] using a 470nm LED and then measure fluorescence from the set of wells [D1, D2, D3] using a 610nm LED, the system would create two *PlateConfiguration.py* objects (one for each set of selected wells); the first *PlateConfiguration.py* object would represent the set of wells [A1, A2, A3] and contains a single *Fluorescence.py* object (for the scan with 470nm excitation), and the second *PlateConfiguration.py* object would represent the set of wells [D1, D2, D3] and also contains a single *Fluorescence.py* object (for the scan with 610nm excitation).

Once the *PlateConfiguration.py* object(s) is/are created, a list of these objects is sent to the *Machine.py* construct. This object is the control center of the entire software. *Machine.py* is responsible for parsing through the list of *PlateConfiguration.py* objects and their respective *Protocol.py* objects. As it parses these protocols, it sends the respective commands to the *Arduino.py* and *Spec.py* structures to perform the necessary operations to successfully execute the inputs. In other words, *Machine.py* is the object responsible for sending commands to the Arduino to tell it to move the linear actuators or to turn on/off the LEDs or external devices (except for the detector). In addition, it extracts spectra from the spectrophotometer and outputs the extracted data to a properly formatted Excel sheet. As a result, there is constant bi-directional communication between the *Machine.py*, *Arduino.py* & *Spec.py* objects.

A more detailed description of the functions and attributes of the involved class objects are below.

[Structure Guide](#)

This section describes each Python object and the functions which are available to call from it. For syntax details and examples of usage please refer to the System Level Control section as well as the commented source code.

[System.py](#)

[Description:](#)

This class object's main purpose is to track user inputs and adjust the system wide settings accordingly. It contains all the necessary functions to add plate objects to the protocol sequence, load calibration data, add protocols, and retrieve information about existing protocols.

[Functions:](#)

- load_calibration_data()** – Used to load calibrated well-position data into system.
- load_settings()** – Used to load the system settings from the settings file. This includes the labeled wavelengths of the excitation LEDs.
- save_settings()** – Used to save system settings to the settings file.
- load_well_labels()** – Used to make a list of well labels when the object is made.
- set_plate_type()** – Creates a new plate object and sets it as either 24 or 96 well.
- set_current_plate()** – Used to select a plate object so that it can be edited.
- get_current_plate_type()** – Returns whether the current plate is 24 or 96 well.
- get_plate_count()** – Returns the total number of plate objects.
- remove_plate_object()** – Removes the plate object.
- get_selected_wells()** – Returns the labels of the selected wells for a specific plate object.
- get_settings()** – Returns a list containing the system settings.
- select_wells()** – Select a specified set of wells for a specific plate object.
- clear_wells()** – Unselect any selected wells from a specific plate object.
- select_all_wells()** – Select all wells for a specific plate object.
- add_absorbance_protocol()** – Add an absorbance protocol to a specific plate object.
- add_fluorescence_protocol()** – Add a fluorescence protocol to a specific plate object.
- add_auxiliary_protocol()** – Add an auxiliary protocol to a specific plate object.
- turn_kinetic_tracking_on()** – Used to mark the start of a kinetic cycle.

add_kinetic_absorbance_protocol() – Used to add an absorbance protocol to a kinetic cycle.
add_kinetic_fluorescence_protocol() – Used to add a fluorescence protocol to a kinetic cycle.
add_kinetic_auxiliary_protocol() – Used to add an auxiliary protocol to a kinetic cycle.
get_all_protocol_brief_descriptions() – Returns brief descriptions (Labels & Index) of all protocols in a plate object.
get_all_protocol_full_description() – Returns detailed descriptions (Labels & Settings) of all protocol in a plate object.
initialize_machine() – Opens an instance of the Machine object.
initialize_datasheet() - Creates the data output file.
start_program() – Executes the protocol sequence.
close_machine() – Closes the instance of the Machine object.

GUI.py

Description:

This object contains all the code that comes together to create the User Interface. Running this object will bring up an instance of the GUI. For details about how to use the interface see the [GUI Level Control](#) section. This class object has access to all the functions that are in System.py.

Functions:

move_plate_out() – Moves plate holder into position by the front door panel.
initialize_calibration() – Runs the calibration sequence
change_settings() – Used to change system settings such as excitation LED wavelengths.
update_settings() – Used to update the system settings file.
start_plate_selection() – Opens the plate selection screen of the GUI
start_well_selection() – Opens the well selection screen of the GUI
start_protocol_selection() – Opens the protocol selection menu of the GUI
load_protocol_sequence() – Loads in a saved protocol sequence from a file.
save_protocol_sequence() – Used to save the currently selected protocol sequence to a file.
add_new_plate_configuration() – Adds a new plate object and brings up well selection menu.
back_function() – Used to return to previous screen.
open_protocol() – Opens the menu for a selected protocol.
add_protocol() – Adds the protocol and its settings to the protocol sequence.
review_protocols() – Brings up the review screen of the GUI
initialize_program() – Starts the execution of the protocol sequence
reset_system() – Resets system. Clears the protocol sequence and returns to main menu.
add_items() – Function to add object to the QTreeWidgetItem (PyQt)
get_current_item() - Function to get selected object from the QTreeWidgetItem (PyQt)
remove_plate() – Deletes a plate object and all its protocols from the protocol sequence.
remove_protocol_selection() – Deletes a protocol object from the protocol sequence.
edit_protocol_selection() – Brings up the menu of the chosen protocol and allows editing.
update_protocol() – Sets the edits made to a protocol and updates the protocol sequence
edit_knetic_protocol_selection() – Used to edit a protocol part of a kinetic cycle
update_kinetic_protocol() - Sets into place the edits made to a protocol that is a part of a kinetic cycle and updates the protocol sequence.
add_parent() – Function to add object to the QTreeWidgetItem (PyQt)

add_child() – Function to add object to the QTreeWidgetItem (PyQt)
add_subchild() – Function to add object to the QTreeWidgetItem (PyQt)

PlateConfiguration.py

Description:

This object contains all the information about a specific plate configuration input by the user. A plate configuration is defined as a set of wells to run specific protocols on. Each plate object contains the following information:

- (1) The plate type (24 or 96)
- (2) List of protocol of objects (the protocol sequence to perform on the selected wells)
- (3) List of selected wells.
- (4) Overall ordering if there is more than one plate configuration input by the user.

Functions:

add_protocol() – Used to add a protocol object to the list of protocols.
get_protocol_count() – Returns the current number of protocols in the protocol list
initialize_wells() – Creates a list of well labels when the plate object is created.
get_protocol_brief_description() – Get a short description of the specified protocol.
get_protocol_full_description() – Get a detailed description of the specified protocol.
remove_protocol() – Removes a protocol from the protocol list.
get_selected_wells() – Returns the labels of the wells selected for the plate object.
get_well_status() – Returns list describing whether a well is selected or not.
clear_wells() – Used to unselect all wells.
selected_all_wells() – Used to select all wells.
select_wells() – Used to select a specified set of wells.

Protocol.py

Description:

This object is a general parent class for all possible protocols. It contains the following general information about the selected protocol:

- (1) The protocol label (Absorbance, Fluorescence, Auxiliary, or Kinetic)
- (2) The plate index to which the protocol is assigned to
- (3) The order number of the protocol in the protocol list of the assigned plate object.

Functions:

get_brief_description() – Returns a short description (Label & Index in protocol list) of the protocol

Absorbance.py

Description:

This object contains all the system settings and necessary information to perform an absorbance measurement with settings defined by the user:

- (1) The exposure time.
- (2) The list of steps that need be run by the machine to perform the measurement.
- (3) The excel sheet layout of the measurement results.

Functions:

get_full_description() – Returns detailed description (Label & Setting Information) of protocol
get_exposure_time() – Returns exposure time of protocol
set_exposure_time() – Sets the exposure time of protocol
initialize_excel_section() – Creates the header of the protocol in the data output file.
get_dark() – Includes list of Machine.py commands to execute to perform a dark measurement.
start() – Includes list of Machine.py commands to execute to perform an Absorbance measurement.

Fluorescence.py

Description:

This object contains all the system settings and necessary information to perform a fluorescence measurement with settings defined by the user:

- (1) The exposure time.
- (2) The list of excitation LEDs to trigger.
- (3) The list of steps that need be run by the machine to perform the measurement.
- (4) The excel sheet layout of the measurement results.

Functions:

get_full_description() - Returns detailed description (Label & Setting Information) of protocol
get_exposure_time() - Returns exposure time of protocol
get_led_index() – Returns a string describing which excitation LED(s) were chosen
get_wavelength() – Returns the wavelength(s) of the selected LED(s)
set_exposure_time() – Sets the exposure time of the protocol
set_excitation_source() – Sets which excitation LED(s) will be used.
initialize_excel_section() - Creates the header of the protocol in the data output file.
get_dark() – Includes list of Machine.py commands to execute to perform a dark measurement.
start() – Includes list of Machine.py commands to execute to perform a Fluorescence measurement.

Auxiliary.py

Description:

This object contains all the system settings and necessary information to perform a measurement using the auxiliary port with settings defined by the user:

- (1) The exposure time.
- (2) The list of auxiliary ports to trigger.
- (3) The list of steps that need be run by the machine to perform the measurement.
- (4) The excel sheet layout of the measurement results.

Functions:

get_full_description() - Returns detailed description (Label & Setting Information) of protocol
get_duration() – Returns the duration in milliseconds
get_aux_index() – Returns which auxiliary port is being used.
set_duration() – Used to set the duration in milliseconds.
start() - Includes list of Machine.py commands to execute Auxiliary functionality.

Shake.py

Description:

This object contains all the system settings and necessary information to initialize and perform shaking of the plate:

- (1) The list of steps that need be run by the machine to perform the measurement.

Functions:

get_full_description() - Returns detailed description (Label & Setting Information) of protocol
start() - Includes list of Machine.py commands to execute Auxiliary functionality.

Kinetic.py

Description:

This object is for time-course measurements. This object represents a set of protocols which need to be performed on each individual well. The following information is contained within this object:

- (1) The order number of the object in the list of protocols of the assigned plate configuration
- (2) The type of kinetic method: Interval (regular interval over specified time) or Repeat (repeated set number of time)
- (3) List of protocol of objects (the protocol sequence to perform in the Kinetic cycle)
- (4) The layout of the excel sheet of the results of the kinetic cycle

Functions:

get_full_description() - Returns detailed description (Label & Setting Information) of protocol
add_protocol() - Add protocol to the Kinetic cycle.
remove_protocol() – Remove protocol from kinetic cycle.
get_protocol_count() - Returns the current number of protocols in the kinetic cycle.
get_protocol_brief_description() - Returns a short description (Label & Index in kinetic cycle) of the selected protocol.

get_protocol_full_description() - Returns a detailed description (Label & Setting information) of the selected protocol.
get_protocol() – Returns the selected protocol object.
initialize_excel_section_start() - Creates the header of the protocol in the data output file that represents the start of the kinetic cycle.
initialize_excel_section_stop() - Creates the header of the protocol in the data output file that represents the stop of the kinetic cycle.

Machine.py

Description:

This object contains all the functions to interface with the Arduino, Ocean Optics spectrophotometer, and Excel sheet to which the resulting data is being written to. Each protocol object contains a list of commands that the machine object needs to perform in order to execute the specific protocol.

Functions:

Initialize_arduino() – Creates an Arduino object and begins Serial communication.
initialize_spec() – Creates a Spec object.
move_plate() – Used to move plate to specific position.
turn_led_on() – Used to turn LED on
turn_led_off() – Used to turn LED off
trigger_aux() – Used to trigger an auxiliary port.
set_exposure_time() – Used to set the exposure time of the spec.
set_scans_to_avg() – Used to set the number of scans to average on the spec.
set_boxcar_width() – Used to set the boxcar width of the spec.
calibrate() – Used to run the calibration sequence
start_program() – Used to start protocol sequence execution.

Arduino.py

Description:

This object is used to interface with the Arduino. It allows the user to write and read to and from the Arduino's Serial.

Functions

write() – Used to write to the Serial port of the Arduino.
read() – Used to read from the Serial port of the Arduino.
close() – Disconnects system from Arduino.

Spec.py

Description:

This object is used to interface with the Ocean Optics Spectrophotometer through their API. It allows the user to set exposure time, boxcar width, scans to average, and to capture spectra.

Functions

For functions, see the linked [Seabreeze](#) documentation.

System Level Control

System level control of the OSP device refers to setting up and executing protocols via a manual Python script. An example script which shows the basics of how to write a Python script to execute OSP protocols is provided in the *Python Files* directory of the repository and can be opened with any text editor. Refer to the comments in the script for details about functions and their use. The script can also be run by typing the following lines of code in the terminal:

Windows

```
activate osp
cd %HOMEPROFILE%
cd "OSP\Software\Python Files"
python system_ex1.py
```

MAC

```
source activate osp
cd ~
cd "OSP/Software/Python Files"
python system_ex1.py
```


GUI Level Control

Running the GUI

If working from a Raspberry-Pi, the best way to initialize OSP is to use the Terminal, direct the Pi to the 'Python Files' directory, and run *GUI.py* from there. If not running on a Raspberry-Pi, simply run *GUI.py* either from the Terminal (refer to the *Running the User Interface* steps in the *Installation* section).

!!! Important Note !!!: The size of the GUI has been optimized for the Raspberry-Pi Touch Screen. If you have a high-resolution display (e.g. 4K), you may need to adjust resolution to display the GUI properly (recommended 1280x720). Users with standard resolution displays (e.g. 1920x1080) should be fine with their default settings and should not have to change the resolution.

The Intro Screen



1. The **START** button leads to the plate selection screen, where the user can input whether they will be using a 24-well plate or a 96-well plate for the protocol sequence.
2. The **CALIBRATE** button allows the user to re-calibrate the well positions in the system for either a 24-well plate or a 96-well plate. When pressed, a pop-up message shows up asking the user to input the plate-type (24 or 96) that they wish to calibrate the positions for.
3. The **PLATE OUT** button allows the user to have the device move the plate holder to the edge of the front door for plate loading and removal.
4. The **SETTINGS** button brings the user to the Settings Window, where they are allowed to change specific settings of the system.

The Settings Screen




Here, the user has access to change the system settings, specifically, the LED wavelength labels and the number of scans to average.

The screenshot shows a settings interface with a blue background. At the top left is an orange 'Back' button. The title 'SETTINGS' is centered at the top. Below the title is the instruction 'Set the wavelength of individual LEDs for excitement of wells.' To the left of the settings is a diagram of a well plate with three green LEDs labeled LED 1, LED 2, and LED 3. To the right of the diagram are three wavelength input fields: LED 1 (490 nm), LED 2 (490 nm), and LED 3 (490 nm). Further right is a section titled 'Measurement Options:' containing a '# Scans to Average:' label and a numeric input field. An orange 'SET' button is at the bottom center. Two arrows point to the interface: arrow '1' points to the '# Scans to Average:' input field, and arrow '2' points to the 'SET' button.

1. Here, the user can enter an integer to change the number of scans that the spectrophotometer averages when performing a measurement.
2. The three LEDs correspond to the three excitation light sources. Here the user can change the wavelengths to describe the LEDs connected to a specific LED socket position (or digital switch). These wavelengths are then used as reference for data output.
3. Once the user clicks the **SET** button the displayed settings are saved to '*.../System Settings/savedSettings.txt*'.
4. Hitting the **BACK** button will take the user to the home screen or the protocol selection screen (depending on which screen was up/active before the settings were accessed).

The Plate/Well Selection


The user has the option to select which type of plate they will perform the protocol sequence on. A protocol sequence can only be run on one type of plate. Plate selection leads the user to the well selection window for the chosen plate type.



*Click & Drag to select wells. Pressing the '**BACK**' button will return the user to the plate selection window.*

[The Protocol Selection](#)

From this screen, the user can choose to add specific protocols to their protocol sequence. Settings may also be changed from this screen.



1. The **ABSORBANCE** button takes the user to the absorbance protocol settings screen. This allows the user to set the exposure time.
2. The **AUXILIARY CONTROL** button takes the user to the auxiliary control settings screen. Here, the user can choose which auxiliary port to turn on and for how long.
3. The **FLUORESCENCE** button takes the user to the fluorescence protocol settings screen. Here the user chooses which excitation LED to use and sets the exposure time of the measurement.
4. The **KINETIC TOGGLE** switch allows the user to turn on/off a kinetic cycle (time-course).
5. The **SHAKE** button adds a protocol to shake the plate for 5 seconds.
6. The **NEW PLATE CONFIGURATION** button allows the user to create a new selection of wells for a different protocol sequence. Going down this path creates a new *PlateConfiguration.py* object.
7. The **LOAD PROTOCOL SEQUENCE** button allows the user to load a pre-existing (custom made) protocol sequence that has been saved to the repository. The loaded sequence will replace the current protocol sequence.
8. The **REVIEW** button takes the user to the review screen where they can see the protocol sequence resulting from their inputs, and edit any of the protocols individually or remove them.
9. The **RESET** button will take the user back to the main screen and clear the system of any protocols and well selections.

[Kinetic Toggling](#)

“Kinetics” allows the user to iteratively perform a specific protocol for time-course measurements without having to manually repeat the inputs. OSP allows two methods of kinetics: (1) Interval and (2) Repeat. Note that the user can only select one method, described below:

Back

KINETIC SETTINGS

☐ Interval

Duration : 1 ms

Time Interval : 1 ms

☐ Repeat

Number of Repeats : 1


SET

[Interval Kinetics:](#)


This method performs a protocol or a sequence of protocols in intervals (based on the interval specified by the user) for a specific duration (e.g. absorbance of a sample every 10 seconds for 2 minutes). ***The user must ensure that the time it takes to perform the specific protocol or sequence of protocols does not exceed the interval time.***

[Repeat Kinetics:](#)

This method simply repeats a protocol or a sequence of protocols the specified number of times. Once the SET button is pressed on the ‘Kinetic Settings’ screen, the user is taken back to the protocol selection menu, which will then have a background that is bright red. The kinetic state is set to ON as long as the background is red. To toggle the kinetic cycle off, press the kinetic button again and the screen will return to its original (blue) color. ***The user will not be able to review the protocol or add a new plate configuration until the kinetic cycle is turned off.***




Kinetics OFF



Kinetics ON

Absorbance Settings

The only tunable parameter here is exposure time. The light source for absorbance measurements by default is the top LED light source (by default, a 430nm LED + White LED).

The image shows a software interface for 'Absorbance' settings. It has a blue background. At the top left is an orange 'Back' button. Below it is a horizontal line, followed by the title 'ABSORBANCE' in black capital letters. Another horizontal line follows. Below this line, the text 'Exposure Time :' is on the left, and a white input field with '1 ms' and a vertical scrollbar is on the right. A third horizontal line is below the input field. At the bottom center is an orange 'SET' button.

Back


ABSORBANCE

Exposure Time : 1 ms

SET

Fluorescence Settings

The user must specify which excitation LED to use. The wavelength of the LED is referenced from the saved settings. The system will display a warning message if the user attempts to press the set button without checking one LED box. Multiple LEDs can be selected.

The image shows a software interface for 'Fluorescence' settings. It has a blue background. At the top left is an orange 'Back' button. Below it is a horizontal line, followed by the title 'FLOURESCENCE' in black capital letters. Another horizontal line follows. Below this line, the text 'Excitation LED :' is on the left, followed by four checkboxes: 'LED 1 (470nm)', 'LED 2 (525nm)', 'LED 3 (610nm)', and 'FIBER'. Below the checkboxes, the text 'Exposure Time :' is on the left, and a white input field with '1 ms' and a vertical scrollbar is on the right. A third horizontal line is below the input field. At the bottom center is an orange 'SET' button.

Back

FLOURESCENCE


Excitation LED : ☐ LED 1 (470nm) ☐ LED 2 (525nm) ☐ LED 3 (610nm) ☐ FIBER

Exposure Time : 1 ms

SET


[Auxiliary Settings](#)

Here, the user can select which auxiliary port (see “**Top Optical Train**” Assembly) they want to turn on by digital logic, and for what duration of time. Note that only one port can be turned on at one time. One port is specifically for unpowered LEDs. The other three communication ports are more customizable; to access their trigger ports, see the DB9 port pin-out diagram on the next page. By default, triggering a port also moves the motors to position the selected well under the top auxiliary optics attachment.




1. The AUX LED PORT (LED #4) controls an extra unpowered LED of choice, and the default limiting resistor sets the current to ~20mA (assuming a ~3V LED bias voltage). Setting the ‘duration’ sets the amount of time this auxiliary LED is on once the well is positioned under it.
2. Any of the three generic ‘ports’ allows direct control of external electrical connections that can be accessed through the DB9 header (see below diagram for pin-out). Each port contains two data pins (1 IN and 1 OUT) for digital control and inputs to/from the Arduino. The user can customize an Arduino script for automating external devices. For custom series of events where movement under the auxiliary piece is not required, this default can be eliminated by commenting out the movement function (`machine.move_plate()`) for the corresponding port number of the Auxiliary.py script.

When accessing the auxiliary port +12V power (pin 3), the total current draw should be less than 3 amps to ensure sufficient power for all internal components. When accessing the auxiliary port +5V power (pin 3), the total current draw should be less than 1.5 amps to prevent damage to the internal switch.



- | | |
|--------------------|--------------------|
| 1. +5V | 6. Port 1: Trigger |
| 2. Port 2: Trigger | 7. Port 2: Data In |
| 3. +12V | 8. Port 3: Trigger |
| 4. GND | 9. Port 3: Data In |
| 5. Port 1: Data In | |

[Review Protocol Sequence](#)



1. Window displaying the selected plates & protocols in their desired order. Each selection in the window is selectable. By highlighting one of the protocols, the user can edit/remove it.
2. The **SAVE LIST** button allows the user to save the current set of protocols in their order as a 'protocol' object. **NOTE:** When saving a sequence of protocols, the system saves the protocols and all their settings in shown order, and the plate type. ***The loaded sequence plate type cannot be changed unless the system is reset.***
3. The **REMOVE** button allows the user to remove one of the protocols OR an entire plate configuration from the protocol sequence. If a plate configuration (e.g. *Plate #1*) is selected, that entire plate object and all of its assigned protocols will be deleted. **NOTE: *The user can only delete one protocol or plate at a time.***
4. The **EDIT** button allows the user to edit a single protocol or a single selection of wells. In order to edit the selected wells, the user needs to highlight the selected wells tab.