

cellSens functions

|                                 |   | Dimension                              | Standard                  | Entry           |
|---------------------------------|---|--|---------------------------|-----------------|
| Layout                          | User experience customization   | ✓                                      | ✓                         | ✓               |
|                                 | Overlay multiple images   | ✓                                      | ✓                         | ✓               |
|                                 | Document groups for side-by-side image comparison                         | ✓                                      | ✓                         | ✓               |
|                                 | Movie playback  | ✓                                      | ✓                         | ✓               |
| View                            | Tile view (multiple images in a single data set shown side by side)       | ✓                                      | ✓                         | ✓               |
|                                 | Slice view for orthogonal plane viewing of 3D or time-lapse data sets     | ✓                                      | ✓                         | ✓               |
|                                 | Voxel view for isosurface and volumetric rendering of 3D and 4D data sets | ✓                                      | ✓                         | ✓               |
|                                 | Snapshot/multi-acquisition  | ✓                                      | ✓                         | ✓               |
| Image Acquisition               | Time-lapse at specified interval  | ✓                                      | ✓                         | ✓               |
|                                 | Automated multi-wavelength  | ✓                                      | Multi-channel Acquisition | ✓               |
|                                 | Z-Stack   | ✓                                      | ✓                         | ✓               |
|                                 | Multi-dimensional (xyz and wavelength)                                    | ✓                                      | ✓                         | ✓               |
|                                 | Graphical Experiment Manager  | ✓                                      | ✓                         | ✓               |
|                                 | Manual assisted panoramic imaging (manual MIA)                            | ✓                                      | Manual Process            | ✓               |
|                                 | Multiposition acquisition and stage navigator                             | Multi-position                         | ✓                         | ✓               |
|                                 | Automated panoramic imaging (auto MIA, requires motorized stage)          | Multi-position                         | ✓                         | ✓               |
|                                 | Instant EPI image (manual or motorized Z)                                 | ✓                                      | Manual Process            | ✓               |
|                                 | Simultaneous multi-color imaging (Image splitter needs)                   | Ratio/High End Device                  | ✓                         | ✓               |
|                                 | Live deblurring   | ✓                                      | ✓                         | ✓               |
|                                 | High Dynamic Range Imaging (HDR)  | ✓                                      | ✓                         | ✓               |
| Image Processing                | Multi-well Plate Acquisition  | Well Plate Navigator and Multiposition | ✓                         | ✓               |
|                                 | Geometry/combine/filter processing  | ✓                                      | ✓                         | ✓               |
|                                 | Fluorescence unmixing   | ✓                                      | ✓                         | ✓               |
|                                 | Brightfield unmixing  | ✓                                      | ✓                         | ✓               |
| Image Analysis                  | Deblurring (No/Nearest Neighbor, Wiener Filter)                           | ✓                                      | ✓                         | ✓               |
|                                 | Kymograph   | ✓                                      | ✓                         | ✓               |
|                                 | 2D deconvolution (constrained iterative deconvolution)                    | ✓                                      | ✓                         | ✓               |
|                                 | 3D deconvolution (constrained iterative deconvolution)                    | CI Deconvolution                       | ✓                         | ✓               |
|                                 | Region and line measurements  | ✓                                      | ✓                         | ✓               |
|                                 | Phase analysis  | ✓                                      | ✓                         | ✓               |
|                                 | Object analysis and classification  | Count & Measure                        | ✓                         | ✓               |
|                                 | Interactive measurement   | ✓                                      | ✓                         | ✓               |
|                                 | Intensity plot over time/z  | ✓                                      | ✓                         | ✓               |
|                                 | Colocalization  | ✓                                      | ✓                         | ✓               |
| Documentation and Collaboration | Object Counting (Manual)  | ✓                                      | ✓                         | ✓               |
|                                 | Online Ratio and Kinetics   | Ratio                                  | ✓                         | ✓               |
|                                 | Ratio analysis (off-line)   | ✓                                      | ✓                         | ✓               |
|                                 | Automatically compose Word reports  | ✓                                      | ✓                         | ✓               |
| Remoting                        | Database image and data management solution for microscopy                | Database Core                          | Database Core             | Database Client |
|                                 | Save and load images/documents from Database                              | Database Client                        | Database Client           | Database Client |
|                                 | Remote Live Image Viewing   | NetCam                                 | NetCam                    | NetCam          |

\* Three points angle, four points angle, arbitrary line, closed polygon, polyline and perpendicular line only.

Products with confirmed functionality

|                                |  | Dimension  | Standard              | Entry |
|--------------------------------|--|--|-----------------------|-------|
| Olympus                        | Camera                                     | DP20 <sup>1)</sup> , DP21, DP22, DP25 <sup>2)</sup> , DP26, DP27, DP70 <sup>1)</sup> , DP71 <sup>3)</sup> , DP72 <sup>2)</sup> , DP73 <sup>3)</sup> , DP80 <sup>3)</sup> | ✓                     | ✓     |
|                                | Microscope                                 | BX43, BX53, BX63, BX61, BX61WI, IX83, IX73, IX81, SZX16A   | ✓                     | ✓     |
|                                | Peripherals                                | IX81-ZDC, IX81-ZDC2, IX3-ZDC   | ✓                     | ✓     |
|                                | Motorized XY stage                         | BX-DSU, IX3-DSU, IX2-DSU, U-CBF  | Multi-position        | ✓     |
| Olympus Soft Imaging Solutions | Camera                                     | CC12, F-View II, Colorview I, Colorview II, Colorview III, Colorview IIIu, XM10, XC10, XC30, XC50, UC30, UC50, SC20, SC30, SC100   | ✓                     | ✓     |
|                                | Peripherals                                | cell^TIRF (multi-line, single line), MT20, USB-ODB converter, Real Time Controller (U-RTC and U-RTCE), U-FCB   | ✓                     | ✓     |
| Hamamatsu                      | Camera                                     | Orca R2 (C10600-10B), Orca O3 (C8484-03G), Orca O5 (C8484-05G), Orca ER (C4742-95-12ER), Orca Flash 2.8  | ✓                     | ✓     |
|                                |  | ImagEM C9100-13, ImagEMX2(C9100-23B), ORCA-Flash 4.0 V2(C11440-22CU), ORCA-Flash 4.0 LT  | High-End Camera       | ✓     |
| Q-Imaging                      | Camera                                     | MicroPublisher 3.3 RTV, MicroPublisher 5 RTV   | ✓                     | ✓     |
|                                |  | Monochrome: Xi Blue/Aqua, RETIGA (Exi, SRV, 2000R, 2000RV, 4000R, 4000RV, 6000, mono) QClick plus RGB slider   | ✓                     | ✓     |
| Photometrics                   | Camera                                     | Color : Xi Aqua  | High-End Camera       | ✓     |
|                                |  | OptiMOS, Rolera Thunder  | ✓                     | ✓     |
| Andor                          | Image Splitter                             | CoolSNAP HQ2   | High-End Camera       | ✓     |
|                                | Camera                                     | Evolve 512 Delta   | Ratio/High End Device | ✓     |
| Jenoptik                       | Camera                                     | Dual View DV2 /QuadView QV2  | High-End Camera       | ✓     |
|                                | Shutter                                    | iXon X3 897, iXon Ultra 897, Zyla4.2 (Camera-link), Zyla5.5(USB3.0)  | High-End Camera       | ✓     |
| Vincent Associates             | Shutter                                    | ProgRes C3, ProgRes C5   | ✓                     | ✓     |
|                                | Light Source                               | Uniblitz shutter (VCM-D1, VMM-D1, VMM-D3)  | ✓                     | ✓     |
| CoolLED                        | Light Source                               | precisExcite (pE-1, pE-2)  | ✓                     | ✓     |
|                                | Light Source                               | X-Cite 120 PC, X-Cite exacte, X-Cite XLED1   | ✓                     | ✓     |
| Sutter                         | Light Source                               | Lambda DG4   | ✓                     | ✓     |
|                                | Shutter, FW                                | Lambda 10-3/10-B   | ✓                     | ✓     |
| Prior                          | Motorized XY stage                         | Proscan (I, II, III), Optiscan   | Multi-position        | ✓     |
|                                | Shutter, FW, Z-drive                       | Proscan (I, II, III), Optiscan II  | ✓                     | ✓     |
| Ludl                           | Piezo Z (control via Real Time Controller) | NanoScanZ NZ100  | High-End Device       | ✓     |
|                                | Motorized XY stage                         | Mac 6000   | Multi-position        | ✓     |
| Objective Imaging              | Shutter, FW, Z-drive                       | Mac 6000   | ✓                     | ✓     |
|                                | Motorized XY stage controller              | Oasis 4i   | Multi-position        | ✓     |
| Märzhäuser                     | Z-drive controller                         | Oasis 4i   | ✓                     | ✓     |
|                                | Motorized XY stage                         | Tango  | Multi-position        | ✓     |
| Physik Instrumente             | Z-drive controller                         | Tango  | ✓                     | ✓     |
|                                | Piezo Z (control via Real Time Controller) | PIFOC P-721  | High-End Device       | ✓     |
| Yokogawa                       | CSU  | CSU-X1   | High-End Device       | ✓     |
|                                |  | CSU-X1   | High-End Device       | ✓     |

<sup>1)</sup> DP20/70 does not support Windows 7 64bit, Windows 8/8.1 32bit/64bit. <sup>2)</sup> DP25/DP71/DP72 does not support Windows 8/8.1 32bit/64bit. <sup>3)</sup> DP73/80 support only Windows 7/8/8.1 64bit.

Compatible image formats

|                |  |
|----------------|--|
| Read and write | JPEG, JPEG2000, TIFF, BMP, AVI, PNG, VSI (Virtual slide image),  |
| Read only      | GIF, PSD (Adobe Photoshop), TIFF (DP-BSW, FSX100, MetaMorph), OIF/OIB (Fluoview format), Cell, STK (MetaMorph), MRC (Medical Research Council) |

Recommended system requirements

|              |  |
|--------------|--|
| OS           | Microsoft Windows 8.1 Pro (32-bit/64-bit), Microsoft Windows 8 (32-bit/64-bit) Pro, Microsoft Windows 7 (32-bit/64-bit) Ultimate with SP1, Microsoft Windows 7 (32-bit/64-bit) Professional with SP1   |
| OS Language  | English, Simplified Chinese, Japanese, German, Russian (only for Entry and Standard) and all others with English like alphabet   |
| CPU          | Intel Core i5, Intel Core i7, Intel Xeon Recommended for high speed image acquisition: QuadCore  |
| RAM          | 4 GB Recommended for high speed image acquisition: 8GB or more only on Windows 7 64-bit operating system   |
| Graphic card | 1280x1024 (min. 1024 x768) monitor resolution with 32-bit-video card with separate graphics memory (no integrated graphics processor with shared memory)   |
| Port         | USB 2.0 port to connect devices to the system Fire Wire A to connect devices to the system (BX61, IX81, SZX2-MDCU, IX3-DSU etc...) Serial (RS232) to connect devices to the system (BX61, IX81, SZX2-MDCU etc...) Additional PCI/PCIe slots as necessary to connect third party peripherals (principally third party cameras) with proprietary interface cards |
| HDD          | 1 GB for installation Performance of hard disk is a limiting factor for image acquisition speed Recommended for high speed image acquisition: Solid State Drive (SSD)  |
| Drive        | DVD drive (Read: DVD-R DL)   |
| Web Browser  | Recommended for Windows 7: Microsoft Internet Explorer 8.0, 9.0, Recommended for Windows 8: Microsoft Internet Explorer 10, Recommended for Windows 8.1: Microsoft Internet Explorer 11  |

Image data courtesy of:

Hiroo Ueno, Ph.D.  
Department of Stem Cell Pathology, Kansai Medical University  
(cover page)

- OLYMPUS CORPORATION is ISO14001 certified.
- OLYMPUS CORPORATION is FM553994/ISO9001 certified.
- All company and product names are registered trademarks and/or trademarks of their respective owners. • Images on the PC monitors are simulated.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.



**OLYMPUS CORPORATION**  
Shinjuku Monolith, 2-3-1, Nishi Shinjuku, Shinjuku-ku, Tokyo, 163-0914, Japan  
**OLYMPUS EUROPA SE & CO. KG**  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
**OLYMPUS SCIENTIFIC SOLUTIONS AMERICAS CORP.**  
48 Woerd Avenue Waltham, MA, 02453, USA.  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373  
**OLYMPUS AUSTRALIA PTY. LTD.**  
3 Acacia Place, Notting Hill VIC 3168, Australia

**OLYMPUS LATIN AMERICA, INC.**  
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.  
**OLYMPUS (CHINA) CO., LTD.**  
A8F, Ping An International Financial Center, No. 1-3, Xinyuan South Road, Chaoyang District, Beijing, 100027 P.R.C.  
**OLYMPUS KOREA CO., LTD.**  
8F Olympus Tower, 446 Bongseunsa-ro, Gangnam-gu, Seoul, 135-509 Korea




Your Vision, Our Future

Imaging Software

cellSens

Seamless Workflow. Intuitive Operation.



# ADD SIMPLICITY TO EXPERIMENT DESIGN... LEAVE MORE TIME FOR RESEARCH

Olympus cellSens gives you a simpler way to work.


Enjoy full control over the user interface, with functions that are where you want them, when you need them.

Seamless operation, from image capture to report creation means more results with less effort.

Spend less time with your software. Have more time for research.




## Imaging



### Image Capture

Capture multi-color, time lapse, and z-stack images with ease. Just select the appropriate capture button, add relevant parameters, and click "Start". The Process Manager or Experiment Manager make it easy to capture multidimensional image.


## Processing



### Viewing and Processing

Automatically view your data in the colors and layout you choose. Take advantage of an array of advanced image processing functions, such as stitching, extended focus, deconvolution, and unmixing.


## Analyzing



### Measurement and Analysis

Make measurements using an intuitive interface. cellSens offers region of interest, phase analysis, and cell count capability. Export raw measurement data to MS Excel or a cellSens workbook with a single click.

## Reporting



### Collaboration and Communication

Actively collaborate with colleagues and coworkers with special tools including Database and Reporting functions. These functions make it simple to manage, share, and distribute your own image and data reports.

### Microscopy Research With a Personal Touch

With microscope optics pushing the boundaries of resolution and size at all magnifications and microscope design enabling new techniques, it is important to be able to efficiently capture and process the images produced. In addition, an increasing number of researchers are imaging using a microscope and it is therefore essential that imaging and analysis are both flexible and user-centric. The Olympus cellSens software family fulfils all these requirements with its unique personalisation concept.

# REDUCE CLUTTER AND CONFUSION BY DISPLAYING ONLY THE TOOLS AND WINDOWS YOU NEED

## It's Time to Get Personal

Olympus has been at the forefront of microscopy for over 90 years and has developed microscopes and systems for a broad spectrum of applications. As a result, we know that each researcher has individual requirements that can't all be met by fixed solutions. The cellSens software family consists of three packages, all featuring a peerless user-definable interface. As a result, each user can define what they want cellSens to show them within the defined work areas.

## Dynamic Interface


Creating an efficient workflow requires careful definition of the tasks and tools at each stage. With the cellSens platform's dynamic GUI, the same is true—the tools you need for each stage are clearly available, without clutter or the need to search. Olympus has created a number of interface layouts, which are developed with capabilities appropriate to the users needs.

- **Acquisition Layout**—for selecting between different acquisition processes and adjusting the camera settings
- **Processing Layout**—for post-acquisition functions such as image processing, execution of measurements, collection of data, presentation of resulting statistics
- **Count & Measure Layout**—for manual and automated measurement and object counting
- **Reporting Layout**—for generating reports to document and share results.
- **Create Layout**—a user can define his or her own layout in various arrangements




## Camera Control Panel

The most important microscope component that requires software control when imaging is the digital camera. Modern cameras feature a number of functions that can be changed to enhance or perfect an image; for example, exposure time and pixel binning. The cellSens Entry and Standard packages control such features on all Olympus digital microscopes and cameras. The Dimension package, in addition, controls such features on high-end research cameras as well. As a result, scientists can maximize the quality of their images.



## Dark Application Skin


The Dark Application Skin reduces computer monitor-generated ambient light and allows cellSens users to adapt to darkened environments; icon contrast remains high for easy recognition and quick selection.



Dark skin

## Arrange Windows as You Like

Organize the tools and windows for the job at hand to create a functional layout that works best for you.



Full screen

Floating panels

Docked panels

**Need Help? Online Help is Just a Click Away**

**Display Only Those Functions You Need on the Toolbar**

**Common Functions can be Grouped in a Single Tab**  
All necessary functions are placed where you want, when you need them. Layout tabs allow easy selection of functions according to your workflow. For instance, display camera control features in your Acquisition layout, and then remove them from view when you switch to the Processing layout.

**Display or Hide Windows as You Require, or Use Auto-hide for Clean Operation**

**Graphical Experiment Manager (GEM)**  
GEM enables the design of complex experiments by simply dragging and dropping icons onto the canvas.

**Create Flexible Workflow Toolbars for Repetitive Operations**  
cellSens lets you create custom toolbars for your most frequently used functions and then save them to the My Functions window. Custom buttons are also easy to use, with convenient tab access that further enhances workflow efficiency. Furthermore, appearance of each toolbar in the My Functions window can be customized by choosing an icon and/or text from Button Appearance window.

**Functional Panels are Grouped in Tabs for Easier Access**

# EMPOWERED TO DO WHAT YOU WANT

## What Researchers Wanted

## Our Solutions


Complex experimental procedure with flexible design

**Graphical Experiment Manager (GEM)**  
This function allows experiments to be designed with even greater versatility. Furthermore, image acquisition is available for up to 6 dimensions (XYZTλ multipoint).




Flexible well plate image capture

**Well Plate Navigator**  
Capture well plate samples automatically by using the well plate navigator in combination with the motorized stage. There is also enhanced flexibility to allow multiple experiments to be executed within a single well plate.




Intensity analysis

**Intensity Analysis**  
Visualize changes in intensity over time, and save this information for later analysis. Ratio Analysis function allows calibration, display and analysis of live/stored data reflecting changes in the intensity ratio between two acquisition channels.



Improved image detail


**Deconvolution**  
Choose between included 2D blind deconvolution and optional 3D blind deconvolution. This proprietary and highly efficient post-processing tool for both CCD and Confocal imaging enhances the ability to differentiate between imaged objects.



Kei Ito, Ph. D.  
Institute of Molecular and Cellular Biosciences, University of Tokyo

Unified task order management

**Macro Manager**  
Perform tasks, from imaging to processing and analysis, as a single composite procedure. Batch processing is also available, enabling multiple images to be subjected to preferred processes as a continuous series for a significant improvement in workflow efficiency.




## What Medical Researchers Wanted

## Our Solutions


Retention of intact observed images

**Live/Snapshot Function with White Balance Adjustment**  
Simply align the focus and select the appropriate white balance to capture images with true-to-life quality.




Observe large sample at once

**Panoramic Imaging**  
Create clear and seamless wide area images by automatic correction of mismatching between each images, even when using the manual stage. A fully functional wide-area focus map enables improved clarity in panoramic imaging.




Simultaneously monitoring of multiple images

**Image Comparison (simultaneous image windows)**  
Display images side by side for accurate comparison, with simultaneous zooming and movement.




Cell counting by hands

**Object Counting**  
Perform manual counts with self-set classes. Counts and proportions can then be undertaken for each class through simple mouse operation.



Nuclei counting with variable thresholding

**Particle Analysis**  
Set threshold levels for nuclei counts, or calculate parameters such as tissue slice total area and area ratios.




# AN ARRAY OF EASY-TO-USE FUNCTIONS TO TURN RESEARCH FINDING INTO COMPELLING PRESENTATIONS

## Image Capture

### Graphical Experiment Manager (GEM)

Dimension


Achieve a high flexibility in the design of experiments, with capabilities such as changing imaging conditions. Furthermore, using the High-end Device Solution provides compatibility with image splitting and piezo devices helps simultaneous two-color imaging and high-speed z-stack image acquisition.



### Well Plate Navigator

Dimension + Multiposition + Well Plate Navigator

The Well Plate Navigator Solution allows you to automatically scan and acquire images from different plate formats, either standard or customized. All acquired images can be saved into a structured database for easy access, together with their well position and user comments. Settings for imaging conditions can also be varied for individual wells, by column, by row or arbitrarily.



### Extended Focus Imaging


Dimension  
or  
Standard + Manual Process

By recording image data while the user gradually focuses through its sample the EFI function automatically creates a single all-in-focus image. The EFI process can be fully automated when cellSens Dimension is integrated into a motorized microscope. Such EFI composites can also be created by combining collections of previously captured images.

### Capture Multidimensional Images

Dimension + Multiposition


In combination with a motorized microscope, the Process Manager makes it easy to capture multi-color and multidimensional images. With the optional Multiposition Solution you can automatically capture multi-point and large area images.



### Panoramic Imaging

Dimension + Multiposition  
or  
Standard + Manual Process

The manual multiple image alignment function creates a single montage image as you scan the specimen. Multiple saved images with adjoining edges can also be combined into a single montaged image. Wide area imaging can be completely automated when cellSens Dimension and its optional Multiposition Solution are combined with a motorized microscope. This function can also be used in combination with a motorized z-focus to enable the capture of images auto-corrected for sample distortion and tilting. With the release of cellSens v.1.11, a multi-point focus map is now available to enable automated focusing across wide image areas.




## Viewing and Processing

### Unmixing

Dimension


With the linear unmixing algorithm in cellSens Dimension, you can readily separate fluorochromes which overlap in their emission spectra—such as GFP and YFP—to produce crosstalk-free fluorescent images. This linear unmixing tool can also separate autofluorescence-related background. Brightfield image unmixing is also available as part of cellSens Dimension.



### Best Focus Extraction

Dimension


Extract the best focus from images, including z-stack, time-lapse images. This function is effective in creating T-series images with the best focus possible, even when working with defocused time-lapse images.



### Deconvolution

Dimension + CI Deconvolution

cellSens Dimension includes a Live 2D deblurring algorithm for live preview and acquisition, to enable better focusing on thick specimens. Additional deconvolution techniques are available in cellSens to reassign out-of-focus light. The optional CI Deconvolution Solution employs the latest in Constrained Iterative Deconvolution technology to produce improved resolution, contrast and dynamic range with industry-leading speed.



### High Dynamic Range Imaging (HDRI)

Dimension

By automatically capturing many images at different exposures the HDRI function creates a final image with a much greater dynamic range, where low intensity signals are clearly visible without overexposing the bright areas of the sample.

# AN ARRAY OF EASY-TO-USE FUNCTIONS TO TURN RESEARCH FINDING INTO COMPELLING PRESENTATIONS

## Measurement and Analysis

### Manual Measurement

Dimension


or

Standard

or

Entry


Depending on the cellSens package different measurements are easily accessible, including distance between points, areas, intensity measurements and morphological parameters. Measurement data is saved as an image layer that can be exported to MS Excel and cellSens workbook formats, or viewed using OlyVia the free image viewer software.



### Automatic Object Measurement and Classification

Dimension + Count & Measure


cellSens Dimension has an extensive set of manual measurements that can be further expanded with the Count & Measure Solution. Easily perform automatic object measurement and classification in an interactive interface where recognized objects are always linked with their measurements.



### Intensity Analysis

Dimension

Graphically depict intensity and channel ratios, and export values to Excel or Workbook by simply setting the region of interest (ROI) on multi-color images captured via FRET or Ca2<sup>+</sup> imaging. Finer details of cell structures can also be brought into clear view through the use of ratio display, thanks to the intensity modulated display (IMD) that displays ratios and intensity in terms of hues and brightness. Furthermore, the ROI can be moved to capture measurements in line with cell movements, and online analysis is made possible through selection of the ratio option.




### Manual Count

Dimension

Standard

Perform manual counts with self-set classes. Counts and proportions can then be undertaken for each class through simple mouse operation.




## Collaboration and Communication

### Reporting

Dimension

A convenient Reporting tool combines images with image property data, measurement data and your own customized fields into a report template with easy drag-and-drop operation. These Microsoft Word\* reports will let you quickly and easily collaborate with colleagues and communicate your results.

\*Requires Microsoft Word version 2003 or later




### Database

Dimension + Database Core or Database Client

Standard + Database Core or Database Client

Entry + Database Client

The Database Core Solution allows the creation of user-defined databases, with full access control, which can be shared across a network. The database not only collects images but also all associated image properties, user comments and any kind of related file, like spreadsheets for other documents. An interactive query tool makes it easy to find the desired data, with automatic preview of the found images. With the Database Client Solution you can then conveniently deploy the capability to read and write to the shared database across many different stations.




### Remote Live Image

Dimension + NetCam

Standard + NetCam

The cellSens NetCam Solution lets any authorized network user see your live image in real time via a web browser.



### Solution

Each cellSens Package can be expanded towards a specific application by using optional "Solutions"

Dimension available solutions:  
 CI Deconvolution    Multiposition    Well Plate Navigator  
 Count & Measure    Ratio    Database Core  
 Database Client    NetCam    Photo Manipulation

Standard available solutions:  
 Multichannel Acquisition    Manual Process    Database Core  
 Database Client    NetCam

Entry available solution:  
 Database Client