BEDIENUNGSANLEITUNG  
INSTRUCTION MANUAL  
MODE D'EMPLOI  
HASZNALATI UTASÍTÁS  
LIBRETTO ISTRUZIONI  
MANUAL DE INSTRUCCIONES  
MANUAL DE INSTRUÇÕES  
INSTRUKCJA OBSŁUGI  
GEBRUIKSAANWIJZING  
NAVOD K POUŽITÍ  
NAVOD NA POUŽITIE  
KULLANIMI KÄÄVÄT  
INSTRUCȚIUNI DE UTILIZARE  
ИНСТРУКЦИИ ЗА ЕКСПЛОАТАЦИЯ  
ИНСТРУКЦИЯ ПО ЭКСПЛУАТАЦИИ

DIGITALE OBERARM-BLUTDRUCKMESSGERÄT  
DIGITAL DEVICE FOR READING BLOOD PRESSURE ON THE UPPER ARM  
APPAREIL DIGITAL POUR MESURE DE LA TENSION AU BRAS  
PÉLIKARI VEINYOMÁS MÉRÉSÉRE ÁLKALMAS KÉSZÜLÉK  
DISPOSITIVO DIGITALE PER LA MISURAZIONE DELLA PRESSIONE DEL BRACCIO  
APARATO DIGITAL DE MEDICIÓN DE LA PRESION SANGUINEA PARA EL BICEPS  
APARELHO DIGITAL PARA A MEDIÇÃO DA PRESSÃO SANGÜINEA NA PARTE SUPERIOR DO BRAÇO  
CYFROWE URZĄDZENIE DO MIERZENIA CIŚNIENIA KRWI NA RAMIENIU  
DIGITALE BOVENARM-BLOEDDRUKMETER  
DIGITALNI TONOMETR NA PAZI  
DIGITÁLNE ZARIADENIE BM122 NA MERAŇIE TLAKU KRVÍ V RAMENE  
DIGITAL KÖL ÖSTÖ TANSZYON ÖLÇME ALETİ  
DISPOZITIV DIGITAL PENTRU MĂSURAREA TENSIUNII ARTERIALE PE BRĂȚ  
Е ДИГИТАЛЕН АПАРАТ ЗА ИЗМЕРВАНЕ НА КРЪВНОТО НАТЯГАНЕ НА ГОРНАТА ЧАСТ НА РЪКАТА  
ЦИФРОВОЕ УСТРОЙСТВО BM122 ПРЕДНАЗНАЧЕНО ДЛЯ ИЗМЕРЕНИЯ КРОВЯНОГО ДАВЛЕНИЯ В ПЛЕЧЕВОЙ ЧАСТИ

# D

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Manufacturer:  
HEALTH o LIFE Co., Ltd.  
9F, No. 186, Jian Yi Road, Chung Ho City 235, Taipei, Taiwan

Authorized Representative in the European Community:  
Innovative business promotion GmbH Botzstr. 6 07743 Jena, Germany

# GB

## INSTRUCTION MANUAL BM122 DIGITAL DEVICE FOR READING BLOOD PRESSURE ON THE UPPER ARM

Dear Customer,

Please read this manual carefully before use in order to avoid damages resulting from misuse. Please pay particular attention to the safety information. If you pass the device on to a third party, these Operating Instructions must also be handed over.

### BLOOD PRESSURE MEASUREMENT

This unit allows you to measure the systolic and diastolic blood pressure as well as the heartbeat frequency oscillometrically. All readings can be read on the LCD display. This blood pressure measurement unit is meant only for measuring on the upper arm of adults. It is for home use and only for non-invasive blood pressure measuring. If the blood pressure is measured by a doctor in the hospital, then the stress or other effects this may have on the patient can result in higher blood pressure. Blood pressure is subject to constant change depending on circumstances. Therefore, one should not draw conclusions from a single measurement. It is better to take measurements at home and schedule regular readings. This will give you a more accurate knowledge of your actual blood pressure and help your doctor to make important decisions for diagnosis or treatment. If you are suffering from heart or liver problems, arterial deposition, diabetes, or similar problems, please seek advice from your doctor before carrying out the measurements.

### WHAT IS BLOOD PRESSURE?

Blood pressure is the pressure which the blood exerts to the arteries. The blood pressure inside the arteries is always fluctuating during your heart cycle. The highest pressure during this cycle is called systole and the lowest pressure is called diastole. With these two blood pressure values, the systole and the diastole, the doctor can determine the blood pressure of the patient. Your blood pressure can be influenced by many different factors, like physical effort, fear, the time of day, etc. The blood pressure is usually low in the morning and increases in the afternoon until the evening. It is generally lower in summer and higher in winter.

# GB

## WHAT IS THE PURPOSE OF MEASURING YOUR BLOOD PRESSURE AT HOME?

Measuring blood pressure during a check-up by a doctor in a hospital or clinic can make a person nervous, which results in a higher blood pressure. Various factors can influence a blood pressure measurement, therefore one single measurement does not necessarily give an accurate reading. The measured blood pressure of a relaxed patient, after getting up in the morning and before eating breakfast, is called fundamental blood pressure. Therefore measuring the blood pressure at home is important, because there the blood pressure is closest to this status.

### A. Classification of blood pressure according to WHO

The standards to measure high and low blood pressure, regardless of age, were determined by the World Health Organization (WHO), as shown in picture 1.


Fig. 1

### B. Fluctuations of blood pressure

Blood pressure fluctuates depending on the time of day and season. The fluctuations are more marked in people with high blood pressure. Normally blood pressure is higher when one is working and then lowers during the night. The fluctuations of blood pressure during a day are illustrated in picture 2. These measurements were taken every 5 minutes. The thick line represents the period of sleep. The rise of the blood pressure at 4 o'clock in the morning (shown in the illustration Fig. 2 as A) and at noon (shown in the illustration as B) represent a pain attack.


Fig. 2

# GB

## POINTS TO NOTE BEFORE USING THE MEASURING INSTRUMENT UNIT

**CAUTION!** Read the instructions carefully before using the device.

Please follow instructions carefully.

Do not use device near mobile phones and microwaves to avoid electromagnetic disturbances.


- Blood pressure should be taken by a doctor or an expert who knows your medical history. By making regular use of the blood pressure measuring unit and recording your readings, you can keep your doctor informed of your blood pressure tendencies.
- Connect the tube with the appliance's tube connector.
- Tie the sleeve securely around your upper arm. The sleeve should be at the same height as your heart. If necessary, support the arm on a table by placing an object (e.g. a cushion) underneath.
- The measuring instrument should not be moved or jostled during the measurement, because this can lead to incorrect measurement of the blood pressure.
- Do the measurement in a quiet place and in a relaxed position without talking.
- Do not twist the seal around the sleeve of a jacket or a shirt, because it is not possible to get an accurate measurement this way.
- Remember that the blood pressure fluctuates naturally during the daytime and other factors such as smoking, drinking alcohol, bathing, taking medicine, and physical effort influence blood pressure levels.
- The values of the blood pressure measurement available in this measuring appliance correspond with the values of a trained specialist, who measures with a sleeve while listening with a stethoscope. The measurement values are between the limit values which were determined by the national standardization for electronic and automated sphygmomanometers of the United States.


# GB

## DESCRIPTION OF THE DEVICE

### A. Part names

### B. LCD display


# GB

## CORRECT USE OF THE SLEEVE

1. Find the brachial artery (Arteria Brachialis) of your left arm, which is located 25 mm above the elbow. Test where your pulse is strongest.


2. Now, place the sleeve around your left upper arm. The hook and loop adhesive fabric must face outwards and the tube should be directed towards the hand. The whole arm should be uncovered. There must not be any pressure on veins from clothing, jewellery, etc., as this can alter the measurement. Also, do not place the sleeve over clothing.

3. Place the sleeve so that it is positioned about 2 - 3 cm above the crease of the elbow. The tube must be in line with the hand, as if the tube were dividing the hand down the middle. Close the sleeve so that it surrounds the upper arm firmly. Check that the hook and loop adhesive fabric is securely fastened.

4. The person whose blood pressure is being taken should now be sitting on a chair with the left elbow supported on a table.

**The tube must not be bent or twisted!**

# GB

## CORRECT BODY POSTURE DURING MEASUREMENT

To obtain as accurate a reading as possible, always remember to hold the sleeve at the same height as the heart. If the sleeve is higher or lower than the heart this causes deviations from your actual blood pressure.

### Measuring blood pressure while sitting:



- Place your elbow on the table or other surface.
- The sleeve should be at the same height as your heart. If necessary, place a cushion or a book underneath your lower arm. Keep your arm and the rest of your body in a relaxed position.
- Any exertion on the body, such as drinking, eating, talking, making or receiving telephone calls, etc., leads to incorrect measurements and should therefore be avoided.

## DESCRIPTION OF THE SYMBOLS IN THE DISPLAY

Brachial Artery

Symbol in the display	Condition/ Cause	Corrective measures
	This symbol appears when the pulse value is detected (blinking) together with a sound.	
	This icon appears when the batteries are almost empty or have been inserted incorrectly.	Replace all old batteries with new ones. Always observe the correct polarity (+ and -) of the batteries while inserting them.
<b>ERROR</b>	This symbol appears when an accurate measurement is impossible.	Repeat the measurement. Make sure that the sleeve is attached properly and your body posture is accurate.


# GB

## INSTRUCTIONS FOR MEASUREMENT

- **Inserting the batteries:**

Remove the cover of the battery compartment. Observe the correct polarity while inserting the batteries. At first all symbols will appear on the LCD display for 3 seconds, then the date 1 (month) – 1 (day) and the time 12:00.

- **Set the date and the time**

Press the button **MODE** and the month indicator will flash.

Set the month with the button **SET**.

Press the button **MODE** again and the day indicator will flash.

Set the day with the button **SET**.

Press the button **MODE** again and the hour indicator will flash.

Set the hour (in 12-hour format, for the afternoon and night time a small "P" will appear under the hour indicator) with the button **↓**.

Press the button **MODE** again and the minute indicator will flash.

Set the minute with the button **SET**.

Press the button **MODE** once again to return to the time and date indicator.

The device is now ready to operate.

- **Steps to measure the blood pressure:**

a) Secure the sleeve around your left upper arm (see chapter *Correct Use of the Sleeve*).

b) Sit on a chair and assume an upright body position (see chapter *Correct Body Posture During Measurement*).

c) Select with the button **SET** the desired number where you would like to save the measurement. Press the button **START/STOP** to start the appliance. The pumping process starts automatically. The pressure is raised automatically to 160 mmHg. If there is a tendency of high blood pressure, the air pressure will automatically be raised to the required level. **The inflation process can be stopped at any time by pressing the START/STOP button.** There should not be more than 2 measurements within a short period of time as this could alter the measured data.

d) After inflation, the air will escape slowly during the measurement process. The heart symbol for the pulse will appear on the display.

## GB

- e) When the blood pressure measurement has been taken, the systolic, diastolic and pulse values will be shown for 1 minute in the display and the data will automatically disappear after this time. During the display time, the appliance can be switched off by pressing the **START/STOP** button.
- f) If no button is pressed for 1 minute and no function is used during this time, the energy saving function will switch the appliance off automatically.
- **To save, check and delete the measured data:**
  - a. **Save data:**

After each blood pressure measurement, the systolic, diastolic, and pulse data of the respective measurement will be saved automatically.  
The available database can store up to 99 measurements.  
If the number of measurements exceeds 99, the oldest measurement will be deleted automatically.
  - b. **Check data:**

By pressing the **MEMORY** button, the database can be selected and, simultaneously, the number of full memory spaces will be shown in the upper edge of the LCD display.  
If no memory space is taken, NONE is indicated.  
Press **MEMORY** button to check the most recent measurement.  
Systolic and diastolic blood pressure, pulse value, time, and date are shown.  
Press **MEMORY** button again to check the other measurement values, sorted by save order.  
To return to normal display, press **START/STOP** button.
  - c. **Delete data:**

**Delete all data:**  
Remove batteries from the battery compartment and reinsert them. All data will be deleted.  
Press **MEMORY** button, to check if the data has been deleted. If done correctly, no data will be shown in each database.

## GB

### • Changing the batteries:

If  appears in the display or the display is dark, the batteries are dead. Replace the batteries. Until they are replaced, no further blood pressure measurements are possible.

Change all batteries with battery type R6/LR6/AA.

Only use top quality identical batteries. Do not use new and expired batteries simultaneously. If the batteries have leaked, remove them from the battery compartment using a cloth and dispose them according to regulations. Wipe battery compartment and other compartments carefully with a cloth. Be sure not to come in direct contact with the battery acid..

**Batteries do not belong in your household waste.**

### MAINTENANCE

- To clean the casing, use a dry cloth moistened with water or a mild cleaning agent, and wipe the unit using a dry cloth. If the sleeve is soiled and needs to be cleaned, use a dry or slightly damp cloth only.
- Do not use any spray cleansers or strong cleaning agents.
- If the unit is not going to be used for an extended period of time, remove the batteries (expired batteries can damage the unit).

### STORAGE

- Always store the appliance in its tray after use.
- Protect the measuring unit and its accessories from direct sunlight, high temperatures, humidity, and dust.
- Do not subject the unit to extremely low (below -20°C) or high (over 50 °C) temperatures.
- Do not fold the sleeve.

### PREVENTION OF MALFUNCTIONING

- Do not twist the sleeve.
- Do not pump air into the sleeve unless it is wound around the arm.
- Do not try to disassemble the unit or to change the components of the appliance or sleeve.
- Do not drop the measuring unit or allow it to suffer any heavy impact.

# GB

## TAKE NOTE

- **Caution:** The sleeve size only covers a size range of 24 to 32 cm in diameter and is only to be used by adults.
- Do not use the measuring unit if your arm is injured.
- If the measuring appliance should continue inflating, open the sleeve immediately or press the **START/STOP** button.
- By pressing **START/STOP** during inflation, the inflation process is stopped and the air is immediately released from the sleeve.
- Should by any means pain occur during the measurement, stop immediately by pushing the **START/STOP** button or by opening the sleeve.

## TECHNICAL DATA

Measuring procedure	Oscillometrical
Measurement range	Pressure 0 ~ 300 mmHg, Pulse 40 ~ 199 beats per minute
Measurement accuracy	Pressure +/- 3mmHg, Pulse max. +/- 5%
Inflation	Electric rolling pump
Deflation	Magnetic valve
Display	Systole, Diastole, pulse value, year/date/time, memory database
Memory Database	with 99 memory spaces
Automatic switch off	1 minute after last pressing of any button
Weight	Approx. 265g with batteries
Sleeve size	Arm diameter of about 24 to 32 cm
Operating temperature range	+10°C to +40°C, <85% relative humidity
Storage temperature range	-20°C to +70°C, <85% relative humidity
Power supply	6V DC 4 x 1,5 V R6/LR6/AA alkaline batteries
Battery life	300 measurements (2 per day)
Measures	80 x 120 x 40 mm

# GB

## INSTRUCTIONS ON ENVIRONMENT PROTECTION

Do not dispose of this product in the usual household garbage at the end of its life cycle; hand it over at a collection point for the recycling of electrical and electronic appliances. The symbol on the product, the instructions for use or the packing will inform about the methods for disposal.

The materials are recyclable as mentioned in its marking. By recycling, material recycling or other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

## GUARANTEE AND CUSTOMER SERVICE

Before delivery our devices are subjected to rigorous quality control. If, despite all care, damage has occurred during production or transportation, please return the device to your dealer. In addition to statutory legal rights, the purchaser has an option to claim under the terms of the following guarantee:

For the purchased device we provide 2 years guarantee, commencing from the day of sale. During this period we will remedy all defects free of charge, which can be demonstrably attributed to material or manufacturing defects, by repair or exchange.

Defects which arise due to improper handling of the device and malfunctions due to interventions and repairs by third parties or the fitting of non-original parts, are not covered by this guarantee.

BM122/0708

ELTA GmbH  
Carl-Zeiss-Str. 8  
D- 63322 Rödermark/  
Ober-Roden, Germany

 0197

# GB

## NOTE

This blood pressure measuring device conforms to European regulations and carries the CE-Symbol "CE 0197".

The quality of this device has been proven and conforms to the requirements of the European Community regulation 93/42/EEC (EN IEC 60601-1 – general safety requirements, EN IEC 60601-1-2:2001 – electromagnetic compatibility; requirements and screening processes) from July 14, 1993 according to medical appliance's and service standards as follows:

- EN 1060-1: Non-invasive blood pressure device – general requirements
- EN 1060-3: Non-invasive blood pressure device – additional requirements for electromagnetic blood pressure appliances.
- EN 1060-4: Non-invasive blood pressure device – screening process to determine the measurement accuracy of automatic non-invasive blood pressure appliances.

### Classifications:

- Device with internal power supply
- BF-type device
- IPX0
- Device is not suitable for use near inflammable anaesthetics mixtures with air or oxygen or nitrogen.
- Continuous operation with short term charging

Manufacturer:  
HEALTH & LIFE Co., Ltd.  
9F, No. 186, Jian Yi Road, Chung Ho City 235, Taipei, Taiwan

Authorized Representative in the European Community:  
Innovative business promotion GmbH Botzstr. 6 07743 Jena, Germany

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