GYENNO WATCH

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This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

Warning: Changes or modifications to this unit not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

FCC Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

CE 2200  RoHS  Bluetooth® FCC ID:2ACGF-GH210001

GYENNO WATCH

Product Manual

| Product Description

The GYENNO Watch, a new generation of customized smart watch that is developed based on Android, continues the development philosophy of GYENNO healthcare and strengthens sensory interaction, family care, health traction, family interaction based on the standard smart watch functions (as a health and information center).

- 1 Watch Screen
Touchscreen display can display information and receive interactive input
- 2 Watch Body
- 3 Strap Connection
The same as the connection of a traditional watch
- 4 Power Button
- 5 Strap
- 6 Heart Rate Measurement Sensor
Heart rate measurement
- 7 Charging Contact
Adopts the GYENNO standard five point rechargeable magnetic interface and contact charging
- 8 GYENNO SafePin Cable
- 9 Reset Button

Package List

GYENNO Watch × 1, Watch Strap × 1, GYENNO SafePin Cable × 1, Product Manual × 1, Brochure × 1.


Main Function

* All pictures are representative schematics, for the actual display effect please refer to the real product.


Exterior

A variety of dial styles can be selected according to your own preferences, meeting your different style and matching demands.


Life Information

In addition to displaying the traditional time and calendar, the GYENNO Watch also displays weather and moon phase information, so the changes of nature are always in your grasp. When the battery is low, the scale point on the dial will show a red warning.


Notification Reminder

The GYENNO Watch can simultaneously display important information received from a mobile phone, including SMS, phone calls, e-mail, QQ and WeChat. These messages will be sorted according to the latest time and displayed on the screen.


Interactive Reminders

If you find that someone else is in an irregular state, such as having forgotten to eat, having a heart rate too high or too low, not having had water for a long time (requires the GYENNO Cup), or having sat down for too long, you can send a tactile reminder to that person.


Time Tools

A timer, alarm clock and stopwatch allow you to freely manage and control your own time.


Heart Rate Monitor

The GYENNO Watch can detect the blood flow under the lower skin via an Optical CMOS Sensor to sense changes in heart rate, thus comprehensively monitoring your health and exercise effectiveness.


Sleep Monitor

The GYENNO Watch can help you completely record your deep/light sleep information at night, including time to fall asleep and sleep distribution, so as to improve your sleeping habits based on the data.


Motion Detection

The GYENNO Watch's built-in motion detection function enables you to monitor your motion in real-time, automatically distinguishing between running and walking up and down stairs, as well as recording and calculating the calories consumed, and urges you to maintain an active status.


Heartbeat Transmission

With the GYENNO Watch, you can transmit your own heartbeat to the people you miss in real-time and convey your heart.

Health Index

With the function of a high and low ranking health index, it can show the health status and trends of contacts in the Watch, so you know their health status. For different health statuses of your relatives and friends, you can send reminders and care via the GYENNO Watch.

Exercise Mode

In exercise mode, the watch can detect the amount of movement, enabling you to burn calories more efficiently while exercising.

Material and design

Appearance

- (1) Case: Moonlight Silver / Space Gray aerospace aluminum alloy
- (2) Strap: Light Gray / Dark Brown calfskin, Space Grey / Silver stainless steel
- (3) External Display: Gorilla Glass 3
Strengthens the protective effect of the display, and the adoption of unique compact technology allows a lower degree of depression than common tempered glass when it is damaged.
- (4) Internal Display: Multi-touch TFT display
An excellent display effect that still displays clearly and thoroughly even when used in sunlight.

Sensors

- (1) Tactile Engine
A linear motor with a unique structure provides you with a real touch sense.
- (2) Heart Rate Sensor
Detects blood flow under the skin through the Optical CMOS Sensor and senses changes in heart rate.
- (3) Other Sensors: 1. Accelerometer 2. Magnetometer 3. Gyroscope
Automatically distinguishes between running and walking up and down stairs, and accurately records your exercise and sleep.

GYENNO SafePin Cable (Magnetic Charging)

The magnetic charging connector of the GYENNO SafePin Cable is designed to quickly dock with the GYENNO Watch. To keep it safe, the connector is held in place magnetically so that if someone trips over the cord, it will pull out of the socket without damaging the connector or the power socket, or pulling the watch off the surface on which it is placed. In addition, the magnet is also designed to prevent the connector from fraying or weakening over time.

Instructions

*For more information on how to use, please visit the GYENNO website and check the *GYENNO Watch User Guide*.

1. Charging

Charging takes about 1.5 hours.

Please charge this product before first use or when the battery is insufficiently charged.

- 1.1 Connect the GYENNO SafePin Cable (Magnetic Charging) ⑧ to the charger or the power transmitting USB interface of an electronic device;
- 1.2 Connect the magnet end of the GYENNO SafePin Cable (Magnetic Charging) ⑧ to the charging contact ⑦;
- 1.3 When charging is in progress, you can see the charging indicator on the watch's display ①;
- 1.4 When it is fully charged, the charging indicator on the watch display ① will show that charging is completed.

2. Power

Long press the power button ④ to switch on. The dial will show the QR Code of this device if it has not been connected with the GYENNO Health APP.

3. Connect to GYENNO Health

* Please make sure your phone system is Android 4.3 or iOS 7.0 or above, and supports Bluetooth 4.0

- 3.1 Download the GYENNO Health APP: If you are an iPhone user, go to the Apple APP STORE and search for GYENNO Health, if you are an Android user, go to the Google Play store, SnapPea or other mainstream application market and search for GYENNO Health and download; or go to www.gyenno.com and open the application download page to download;
- 3.2 Complete the user registration process according to the APP user registration wizard;
- 3.3 Enter the APP settings page and choose New device, and complete the connection between the GYENNO WATCH and GYENNO Health according to the setup wizard.

4.Adjust and Fit

Wear the GYENNO Watch to get started. In order to obtain the best experience, please adjust the tightness of the strap to fit your wrist.

5.Touch-screen Operation

 <p>Tap</p>	 <p>Swipe from left to right</p>	 <p>Swipe up and down</p>
<p>Confirm/Enter</p>	<p>Return to the previous screen</p>	<p>Move up and down to view the screen interface</p>

* For more information on how to use, please visit our GYENNO website and check the *GYENNO Watch User Guide*.

6.Switch Screen Interface, Restart, Shutdown and Restore Factory Settings

Short press the power button  to switch on/off the watch screen; long press the power button  to enter the system interface choose restart, shutdown, or reset factory settings.

7.Watch Upgrade

The GYENNO Watch firmware can be upgraded. In the GYENNO Health APP, enter "Settings - My device - My Watch - Watch Upgrade" to upgrade. Before upgrading, please make sure the watch and phone is fully charged.

Maintenance and Repair

1. The GYENNO Watch has a splash-proof feature but cannot be immersed in water. You can use it when you are exercising, in the rain for a short duration and amid domestic water spluttering, but it cannot be used when you are swimming, taking a shower or soaking. For example, you can wear and use the GYENNO Watch while you exercising, washing your hands or in the rain, but you cannot immerse it in water. In addition, the leather strap cannot come into contact with water.
2. If water splashes on the watch, please clean and dry it with a soft, lint-free cloth. If you need to clean the dial, please clean it when it has been shut down with a cleaning cloth that has been dampened with water. You cannot use any liquid other than water to clean it. Do not clean when the GYENNO Watch is charging.
3. Use GYENNO brand straps or straps authorized by GYENNO.
4. To clean the strap, please detach it from the GYENNO Watch. For the leather section of the strap, please use a soft and lint-free cloth that has been moistened with a small amount of water (if needed) and then clean. After cleaning, let the strap dry completely. When storing the leather strap, keep it out of direct sunlight, high temperatures or humidity. Please do not immerse the leather strap in water; the leather strap is not water resistant.

This product benefits from a nationwide quality assurance and warranty service. Its warranty period is one year.
If you need maintenance or repairs, please call GYENNO Customer Services:
4008-033-037

Notes

The GYENNO Watch is not a medical device. The battery life and charging times may vary due to specific use and settings. The device color may vary and change over time. In order to give users a better experience, GYENNO will continue to optimize and improve the product. During the process of use, actual results may be inconsistent with the display effect in the pictures and presentation. Please refer to the real product.

For more information about the GYENNO Watch, please visit:
<http://www.gyenno.com/watch.html>

Contact Information

If you need services or information, or have any questions, please do not hesitate to contact us:

WEB: www.gyenno.com

Service Email: service@gyenno.com

Service Tel: 4008-033-037

Manufacturer: GYENNO Technologies CO., LTD. Authorized Products

Address: Room B-201, Shenzhen Institute of Industry-University-Research of Nanjing University, 3rd YueXing Road, Science Park, Nanshan District, ShenZhen, P.R. China

Product Certification

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

Warning: Changes or modifications to this unit not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

FCC Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

CE 2200  RoHS  Bluetooth® FCC ID:2ACGF-GH210001