
Setting Up the MedSense Beacon

Step 1: Registering Devices

Log into MedSense HQ and follow prompts to register each component of the MedSense system. You will need the serial number found on each device in order to complete the registration process.

Step 2: Installing and Calibrating a Beacon

- Locate the Beacon on the wall behind the patient's bed along the centerline of the bed. Keep a clear line of sight between the Beacon and the Badges.

- Install brackets at heights between 1.5 and 2.5 meters, using either the screws or adhesive provided.
- Supply power to the Beacon using only the provided power supplies, either by plugging into the wall or using batteries.
- Certain patient rooms will require more than one Beacon and placement will vary depending on room layout.

Using the MedSense Beacon

MedSense Beacon

POWER

When the Beacon is fully charged and working properly, you will see a blue blinking light. If it is battery powered, batteries will last approximately six months. A blinking red light means low power and maintenance is needed promptly.

To change the batteries, remove the Beacon from the wall and unscrew the battery door. Replace the batteries with six AA batteries, either Alkaline or rechargeable.

CARE INSTRUCTIONS

To clean the Beacon, wipe down in place with a cloth.

Legal Statements

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

To assure continued compliance, any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

Radiation Exposure Statement

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment and it also complies with Part 15 of FCC RF Rules. Operation is subject to the following two conditions:

- This device may not cause interference, and
 - This device must accept any interference, including interference that may cause undesired operation of the device.
-

Caution!

The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user authority to operate the equipment.

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