
# **Wi-Mic** HB-1303 User Manual





# Identifying parts



- 1. Power button
- 2. Mute button
- 3. Charging/Power LED
- 4. External Mic port
- 5. Internal Mic
- 6. Hook for neck loop
- 7. Reset button
- 8. Charging port for Base Unit
- 9. Micro-USB charging port
- 10. Wireless pair button

# Charging the Wi-Mic

- 1. Turn off the Wi-Mic.
- Plug the micro-USB charger into micro-USB charging port of the Wi-Mic. RED charging/power LED lights up during charging.
  (Optional: The Wi-Mic can be charged when placed it into the charging port of the Base Unit.)
- 3. The charging/power LED is off when charging is completed.

#### **Basic operation**

#### Power On

Press and hold the power button for 2 seconds until GREEN power LED lights up.

## Power off

Press and hold the power button for 2 seconds until power LED is off.

#### Mute On

To mute the internal/external Mic, press the mute button once. The GREEN power LED is flashing when Wi-Mic is in mute mode.

## Mute off

To cancel the mute status of Wi-Mic, press mute button once. The GREEN power LED lights up and stops flashing.


- 1. Turn on both Wi-Mic and PAE-300.
- 2. First, press the TV function button on all of the PAE-300 needed for wireless pairing to set them into TV mode. Then, press and hold the Wireless Pair Button on those PAE-300 until the GREEN power LED on those PAE-300 is flashing.
- 3. Second, press and hold the Wireless Pair Button of the Wi-Mic until the GREEN power LEDs on both PAE-300 and Wi-Mic are fast flashing for pairing.
- 4. The GREEN power LEDs on both Wi-Mic and PAE-300 are lit solid when the pairing succeeds.

#### Note:

please follow the above procedures; otherwise, it will result in pairing failure.

#### **How Wi-Mic Works**

- Spokesman can speak to the internal mic on Wi-Mic. The voice of the spokesman will go through the Wi-Mic, and then is transmitted to all PAE-300 which are paired with the Wi-Mic. The PAE-300 users are able to listen to the spokesman's voice through the PAE-300 headset.
- External mic can be used by inserting it into the external mic port on the Wi-Mic. Spokesman can speak to the external mic instead of the internal mic without holding the Wi-Mic all the time.



## **Reset Wi-Mic**

When the Wi-Mic cannot be operated in any situation, press the reset button by a small pin object to reset the Wi-Mic.



www.vitasound.com info@vitasound.com 1-877-967-0742