



Contents

1	<i>SAFETY</i>	2
1.1	Warnings	2
1.2	Symbols for Safety	3
1.3	Function Symbols	4
1.4	Remarks for Safe Use	6
2	<i>Overview</i>	12
2.1	Description of MiroCam Capsule Endoscope System	12
2.2	System Overview	13
2.3	Feature of MiroCam Capsule Endoscope System	14
2.4	Function Explanation	15
2.5	Intened Puposue and Side Effect	18
3	<i>Product</i>	19
3.1	Pictures	19
3.2	Supply Items	21
4	<i>Installation</i>	22
4.1	Installation	22
5	<i>Technical Data</i>	24
5.1	Overview	24
5.2	Hardware	24
5.3	Software	25
5.4	COMPLIANCE/APPROVALS	25
6	<i>Environment condition for storage and transport.</i>	27
6.1	Safety environment condition for storage	27
6.2	Environment condition for transport	27
7	<i>Troubleshooting</i>	28
7.1	Introduction	28
7.2	Who should perform repairs	28
7.3	Obtain replacement parts	28
7.4	Troubleshooting guide	28

1 SAFETY

1.1 Warnings

The MiroCam® has been manufactured to conform to the International Standard for Medical Electrical Equipment: General Requirements for Safety IEC 60601-1, together with the Collateral Standard for Electromagnetic Compatibility Requirement and Tests IEC 60601-1-2. The MiroCam® has been manufactured to conform Medical Equipment with respect to electric shock fire and mechanical hazard only CAN/CSA C22.2 NO.601.1.

To guarantee safety and reliability of the equipment, special attention should be paid to all the precautions directed throughout this manual, such as " **WARNING**", " **CAUTION**", and " **NOTE**". Also, to prevent incorrect operations and a failure in maintenance, it is essential that you should be familiar with the functions, operation and maintenance instructions by thoroughly reading this manual before operating the equipment.

WARNING

"**WARNING**" is used to alert the user of a potential serious hazard that may lead to death or serious injury of patients and personnel.

CAUTION

"**CAUTION**" is used to indicate information the user should know to avoid a hazard that may lead to device malfunction, device failure, damage to the device, or damage to other property.

NOTE

"**NOTE**" is used to indicate important information on installation, operation, or maintenance which you should to know for the safe and effective use of the equipment.

1.2 Symbols for Safety

This section describes a set of symbols that IEC (The International Electrotechnical Commission) has established for medical electronic equipment to classify a connection or warning of any potential hazards.



IEC 348: Notice for the user to pay special attention to the following details



IEC 878-02-03: Indicates that this is classified into Type BF equipment



EN 980: Denotes Date of Manufacture



EN 980: Denotes Date of Manufacture Address



IEC60601, ANNEX D: Denotes "ON" status of main power switch



IEC60601, ANNEX D: Denotes "OFF" status of main power switch

SN

EN 980: Denotes serial number



IEC 417-5031: Denotes DC (Direct Current)



IEC 417-5032: Denotes AC (Alternating Current)

A

Denotes Ampere, the unit of current

V

Denotes Volt, the unit of Voltage

Hz

Denote Herz, the unit of Frequency



IEC 417-5021: Denotes potential equalization terminal



























Single Use Only















The Available Period




1.3 Function Symbols

1.3.1 Application Function Symbols

Symbol	Description	Symbol	Description
	Functions & Windows for Receiver Unit		Functions & Windows for Diagnosis Report
	Function & Windows for Patient Database Listing		Function & Windows for Export Patient Image Data
	Function & Windows for Image Display and Diagnosis		Function & Windows for Backup Patient Image Data
	Capturing current patient image to the thumbnail		Landmark current time of patient image to the thumbnail
	Single previous image		Single next image
	Previous Selected image		Next Selected image
	Display image continuously		Display image continuously reverse direction
	Pause image display		Enlarge image display
	Single image display mode		Dual image display mode
	Quad image display mode		Display full image
	Display selected image		Display red image
	Grouping image display in Thumbnail		Single image display in Thumbnail
	Minimize application software		Close application software
	Mark circle in current image		Mark arrow in current image
	Erase circle and arrow in current image		Select marker color

Symbol	Description	Symbol	Description
	Undo mark in current image		Redo mark in current image
	Doctor login		Cancel login
	Patient's diagnosis history		Delete current patient information
	Progress reporting current patient		Complete reporting current patient
	Complete Export current patient		Complete backup current patient
	Doctor logout		Display full patient list

1.3.2 Receiver Function Symbols

	Symbol	Description
SIG		Status of capsule signal receiving Green : Good Yellow : Bad
INI		Status of initialization of Receiver Unit Green : Good Yellow : Bad
BAT		Status of Battery Green : Full charged Yellow : Not charged


1.4 Remarks for Safe Use


- Follow these safety instructions and the clinically adopted precautions.
- The manufacturer is not liable for harm or damage caused by improper, unauthorized, unprofessional or inexperienced use of the device
- IntroMedic.Co.Ltd., is NOT responsible for physical harm or equipment problems caused by the user's careless operation or mismanagement.
- The user MUST have read and understood the user's manual. ONLY trained and qualified persons may operate the system.
- The user's manual must ALWAYS be with the equipment. It is the USER'S RESPONSIBILITY to ensure this!
- CAUTION is advised that the foreign substances, ironware or fluids, such as water or disinfecting fluid, not enter the equipment.
- ONLY authorized personnel may perform any type of repair. Never attempt to open the covers, panels or casings.
- DO NOT crease, bend, fold or twist the probe cables and take care to guard them against mechanical stress (e.g. wheels or heels)!
- The probes must not be exposed to mechanical shock (e.g. by dropping). Any damage caused that way invalidates the warranty.
- CAUTION: Damage/injuries to the probe or cable may lead to safety hazard. It is HIGHLY advised to have them repaired IMMEDIATELY.
- DO NOT handle fluids in the vicinity of the system.
- When using a cart purchased elsewhere, implement the latch guard method to prevent the wheels from rolling.
- DO NOT USE in moist or damp places.
- DO NOT operate the equipment with wet hands.
- Avoid using the equipment in extreme temperatures or humid environments.
- DO NOT keep the equipment or carry out the procedure in places such as areas exposed to direct sunlight, vicinity of heaters, vicinity of chemical materials or gases, areas moist/damp or dusty, or poorly ventilated areas.
- DO NOT disassemble or open the equipment without permission. The free repair warranty will be automatically invalidated.
- DO NOT carry out the procedure in areas with high vibrations or in environments where high electro-magnetic waves are generated.
- DO NOT pull out the power cord by grabbing the cable. When disconnecting the power cord, grab the plug and then pull out. This prevents short-circuits, disconnection, or cord damage.


- CAUTION: Verify that the power voltage supplied from the power receptacle, matches with the voltage the system requires. Check Voltage and Frequency on the AC/DC adaptor.
- CAUTION: Verify that all connection terminals are securely connected to the system.
- CAUTION: Turn off the power switch before connecting the probe.
- DO NOT discard probes, cables and connectors with general waste. Separately discard as industrial waste or medical waste.
- CAUTION: Discard the battery according to the regulations of industrial waste. DO NOT discard with general waste.
- DO NOT carry out the procedure simultaneously with other procedures using medical products or equipment.
- DO NOT carry out the procedure on patients with pacemakers and/or defibrillators.
- DO NOT use for purposes other than medical treatment.
- DO NOT connect the USB cable to the receiver while the receiver's patch probes are still connected
- DO NOT charge the receiver while the receiver's patch probes are still connected.
- Connect USB cable to receiver only after mounting it on charger.
- DO NOT install any other programs onto the workstation utilized for review and diagnosis of patient image data.
- DO NOT carry out the procedure on patients with known or suspected gastrointestinal tract obstruction, perforation, stricture or fistular.
- DO NOT carry out the procedure on patients who have difficulty swallowing food or pills(dysphagia).
- DO NOT carry out the procedure on patient who have difficulty communicating.
- DO NOT carry out the procedure on patient with indigestion or delay in gastric emptying.
- DO NOT carry out the procedure on patient who may be affected by electromagnetic radiation, such as pregnant women, infants, and patients with heart disease or epileptics.
- DO NOT carry out the procedure on patients with diverticulosis in the Small Bowel.
- DO NOT carry out the procedure on patients who are recommended against having the procedure by a Physician.
- The capsule is disposable and should not be reused.


1.4.1 Environmental Conditions for Operation

- Temperature : -10 - +70 (Operating : +10 - +40)
- Relative humidity : 10% - 80% (Operating : 45% - 75%)
- Atmospheric pressure : 700hPa to 1060hPa


 **WARNING** DO NOT operate the equipment in the vicinity of generators, power stations, X-ray devices, and broadcasting stations where high levels of electro-magnetic waves are generated. The electro-magnetic waves can cause equipment malfunctions.

 **CAUTION** If the equipment has been brought in from a cold environment (stock room, airfreight) into a warm room, initial activation should take place after a few hours, to allow for temperature adjustment and balance and evaporation of condensed humidity.


 **WARNING** DO NOT operate the equipment in the vicinity of heat sources, strong electric or magnetic fields (close to a transformer), or near instruments generating high-frequency signals.


 **WARNING**


- Do not use this device when the output contains d.c. components.
- Do not use the unit with h.f. surgical equipment. It may result in burns at the site of the electrodes and possible damage to the capsule and receiver.
- Do not use the unit in close radius (within 1 m) of short wave or microwave therapy equipment. It may produce instability in the captured image.


 **WARNING** This device is a Class B device according to EN60601-1-2 standards. This equipment can cause radio interference in residential areas. In this case, the owner (or operator) can be held responsible to take appropriate measures or take proper measures for compensation.

1.4.2 Safety Precautions

-  **CAUTION** - Make sure the environment is without interference from electromagnetic fields.
- Make sure the environment is without noise and vibration.
 - DO NOT carry out the procedure while using other equipments, devices or products.
 - The instruction for use of the probes MUST be observed.
 - DO NOT use on patients with pacemakers or defibrillators.

 **NOTE** This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of U.S. FCC Rules.

 **CAUTION** EQUIPOTENTIAL BONDING:
In hospital, doctors and patients are subjected to dangerous, uncontrollable compensating currents. These currents are due to the potential differences between connected equipment and touchable conducting parts as found in medical rooms. The safest solution to the problem is consistent equipotential bonding. Medical equipment is connected with connecting leads made up with angle sockets to the equipotential bonding network in medical rooms.

 **CAUTION** DO NOT use the capsule if the package is unsealed. If capsule package is unsealed, it can be infected.

- DO NOT reuse a used capsule.
- To prevent unexpected accidents like fire or explosion, do not use it near or in the presence of inflammable or ignitable substances.
- DO NOT disassemble the equipment case nor open the cover. In case of problems in equipment or need for other service, please contact our head office or a customer support staff immediately.
- Only the accessories authorized and designed by our company should be used with this equipment. Faults resulting from the usage of unapproved or inappropriate accessories are not guaranteed against.
- This equipment may have an effect on other products or be effected by them.
- Follow your doctor's instructions and abide by the guidelines in the user manual.
- DO NOT try to upload the data while the patch probes are still connected to the receiver.
- DO NOT charge the rechargeable battery in the receiver while the patch probes are still connected to the receiver.
- Stay away from high frequency radiation sites (such as high voltage, radar, installation power plants, MRI, CT or electric blankets etc.) during your capsule endoscope procedure. (It may result in serious side effects requiring an emergency operation.)

- In case of any symptoms of abdominal pain, vomiting, fever, heart trouble, dizziness or seizure during or after the capsule endoscope procedure, please notify your doctor.
- Always check the connection between the receiver and the patch probes.
- Always check that the battery in the receiver is fully charged before use.
- DO NOT use the capsule if the package is unsealed.
- After ingesting the capsule, always check whether the capsule has been excreted.
- Prior to undergoing the capsule endoscopy procedure, patients with diabetes must control their insulin doses.
- For more accurate data and better analysis, patients can have solid food for lunch, but must have liquid food as dinner on the day before the procedure and need to fast for 12 hours before the procedure.
- 12 hours before the procedure, take a laxative such as Sodium phosphate or PEG solution.
- DO NOT smoke for 12 hours prior to the capsule endoscope procedure.
- DO NOT apply body lotion to your body before the procedure.
- DO NOT bite the capsule before swallowing.
- Avoid excessive physical activity during the capsule endoscope procedure.
- 2 hours after swallowing the capsule, patients may drink water, and after 4 hours, patients may have liquid food.
- When undergoing the capsule endoscope procedure, DO NOT make physical contact with another person undergoing this same procedure.
- During operation of the receiver, DO NOT touch it or get it wet.
- Use the provided batteries only and never take it out of the receiver body during the procedure.
- During upload of the data recorded in the receiver to the PC, avoid disconnecting the USB (It may damage the patient's data).
- Always confirm that the USB is connected by checking on the screen of the MiroCam application program.

1.4.3 Cleaning and maintenance

1.4.3.1 System and accessories

- Other units and accessories should be cleanly maintained. For cleaning, rub them lightly with a soft cloth wet with warm water at least once a week. Do not use organic solvents such as lacquer, thinner, ethylene and oxide because they can damage the equipment. Be careful that foreign substances do not enter the main system when cleaning. Please maintain the probe according to the disinfection instructions.
- The vaginal manometry probe must not be cleaned with a wet towel before usage. Use the vaginal probe ONLY after putting a disposable latex condom on it, and then discard only the condom after use.
- The patch probe for abdominal attachment must not be rubbed or wiped.
- ALWAYS operate the equipment under sanitary environmental conditions. DO NOT use heat or gas for the disinfection of the probes and cables.
- CLEAN vaginal probes, anal probes, patch probes and vaginal manometry probe should ALWAYS be used.

1.4.3.2 Service documents

If required, or upon request, the supplier may provide block diagrams, list of spare parts, descriptions, adjustment instructions or other related information which may help qualified technical personnel in repairing specified parts of the equipment which have been defined repairable by the manufacturer.

1.4.3.3 Moving the Equipment

- CAUTION when moving equipment.
- WARNING: Excessive impact/shock causes internal damage.
- If wiring is connected/disconnected when moving, check the exact wiring status after moving.
- If damage to the equipment is found out after moving, immediately contact the company or dealer.

**WARNING**

Before moving the system, always make sure to disconnect the monitor from the main system, and then safely move the main system and monitor separately. Connect the main system and monitor only after moving is finished and the systems structure is settled and stabilized.

2 Overview

2.1 Description of MiroCam Capsule Endoscope System

MiroCam is Capsule Endoscope Systems and accessories to be used to diagnosis of patients. It consists of the capsule and receiver.

2.1.1 Photos

2.1.1.1 Capsule



2.1.1.2 Receiver




2.1.1.3 Accessory



2.2 System Overview

MiroCam Capsule Endoscope Systems consist of Capsule Receiver and Image Processing Workstation (Software).

MiroCam System Block Diagram



2.2.1 Capsule

MiroCam Capsule Endoscope System' Capsule Unit consists of optical dome, white LED, image sensor, lens and gold courted cage.

2.2.2 Receiver

MiroCam Capsule Endoscope System's Receiver unit consist of 8 channel connector, differential amplifier, filter, data demodulator, data storage memory and USB communication channel.

2.2.3 IPW

MiroCam Capsule Endoscope System's Software is running on the Microsoft Windows XP using the IBM compatible PC platform and provide easy-to-use icon menu. It consist of Image reconstruction part and Image display part.

2.3 Feature of MiroCam Capsule Endoscope System

2.3.1 Main Feature


- Capsule
 - High Resolution Image (320 X 320)
 - Widest Field of view (150 degrees)
 - Smallest Capsule size (11 X 24mm)
 - Longest operation time
- Receiver
 - Built-in circuit protection functions
 - Small, light and easy to carry around
- IPW
 - User Friendly software feature
 - Red Color Detection
 - Color of images are virtually natural
 - Dual and quad image views
- System
 - Over Intel P4 3.0GHz
 - Over 1G Byte Memory
 - Over 100GByte Hard Disk

2.4 Function Explanation


2.4.1 Display flow





Image Review







Report View




2.4.2 Function of Control and Display of Signal

- All operations start after completing or selecting the patient registration on the default screen.
- If select Patient in Patient List for data review and press  button, Review window is displayed.
- Press  button for playing image data.


- Press  button for stop playing image data.
- Press  button for save current image to thumbnail.
- Press  button for recording diagnosis data and printing paper.
- Press  button for operation for receiver unit initialization and image data upload.

2.4.3 Sequence of Operation


- Initialize Receiver Unit

Connect USB Cable between Receiver unit with Workstation → Click  button
→ Click **Format** button Click **Input Patient Info** button → Input
Patient information → Click **Save** button


- Upload Image Data

Connect USB Cable between Receiver unit with Workstation → Click  button
→ Click **Upload data** button

- Display Image Data and Diagnosis

Select Patient in Patient List for data review → Click  button → Click  button

- Print Diagnosis Data

Select Patient in Patient List for data review → Click  button → Input
diagnosis result → Insert diagnosis point in patient image data

2.4.4 Connection and Disconnection

Part	Where	When	Method
Patch Probe	Connector of upper side of receiver unit	Before process capsule endoscope	<p>Connection : Hold the connector on the end of the receiver's patch probe and put it into the socket of the top of the patch.</p> <p>Disconnection : Hold the connector on the end of the receiver's patch probe and pull it out by hands..</p>
USB Cable	USB Connector of left side of receiver unit	After process capsule endoscope	<p>Connection : Hold the USB cable connector and push it into the USB connector on the side of the receiver.</p> <p>Disconnection : Hold the connector on the tip of the USB cable and pull it out by hands.</p>
Cradle	Lower side of receiver unit	After process capsule endoscope	<p>Connection : Place down the receiver into the Cradle as its slot. .</p> <p>Disconnection : Pull out the receiver from the Cradle vertically.</p>
Adaptor	Rear side of cradle	System Installiation	<p>Connection : Put the DC cable of the charger into the DC jack of the Cradle..</p> <p>Disconnection : Pull out the connector of the DC cable of the charger by hands.</p>

2.5 Intened Pupose and Side Effect

2.5.1 Intended Purpose

- **MiroCam** is intended for the detection of pathologies as an adjunctive tool in the diagnosis of small bowel gastrointestinal disorders and diseases.

3 Product

3.1 Pictures

3.1.1 Main Device

- Capsule



- Receiver



3.1.2 Accessories

■ Patch Cable



- Model : 21-1-19-1, 10P
- Manufacturer : Mednis

■ USB Cable



- Model : AV22201-06
- Manufacturer : BELKIN

■ Cradle



- Model : MR1000-C
- Manufacturer : JEC Korea

■ Adaptor



- Model : JMW128KAXXXNXX(Compliance with requirements of IEC601-1)
- Manufacturer : JEC Korea

3.2 Supply Items

3.2.1 Basic Supply Item

No.	Model	MiroCam Product Code	Quantity	Description	Replace Period
1	Capsule	MC1000-C	1 ea		-
2	Receiver	MR1000-R	1 ea		-
3	Cradle	MR1000-C	1 ea		-
4	Adaptor	MR1000-A	1 ea		-
5	Power Code	MR1000-P	1 ea		-
6	Patch Probe	MR1000-S	1 pk		Single use
7	Ground Pad	MR1000-E	1 ea		Single use
8	Patch Cable	MR1000-D	1 ea		6 month
9	Software	MW1000-S	1 ea		-
10	Sack	MR1000-K	1 ea		-
11	PC	MW1000-P	1 ea		-
12	Monitor	MW1000-L	1 ea		-
13	Printer	MW1000-I	1 ea		-
14	Keyboard	MW1000-K	1 ea		-
15	Mouse	MW1000-M	1 ea		-
16	Back	MR1000-G	1 ea		-
17	Cart	MW1000-C	1 ea		-
18	User Manual	MW1000-U	1 ea		-

3.2.2 Optional Supply Item


No.	Model	MiroCam Product Code	Quantity	Description	Replace Period
1	Battery	MR1000-B	1 ea		6 month
2					
3					
4					

4 Installation


4.1 Installation

4.1.1 Specification of packing

■ MiroCam Capsule Endoscope System Box



■ MiroCam Capsule Endoscope Receiver Set Box



■ MiroCam Capsule Endoscope Battery Charger Box

Rev#	DATE	REVISION/NOTE	CHECK
△	..		
△	..		

No.	DESCRIPTION	MATL	QTY.	REMARK
1	Charger Box	Pulp	1	---
SCALE	NIS	MODEL	Box	
DATE	2017.02.22	ITEM	Charger Box	
DWG. Sgn	SH/KM	DWG. NO.	DWG. NO. 1/1	
CHE. Sgn				
ENG. APP. Sgn				
USR. APP. Sgn				

■ MiroCam Capsule Endoscope Accessory Box

Rev#	DATE	REVISION/NOTE	CHECK
△	..		
△	..		

No.	DESCRIPTION	MATL	QTY.	REMARK
1	Accessory Box	Pulp	1	---
SCALE	NIS	MODEL	Box	
DATE	2017.02.22	ITEM	Accessory Box	
DWG. Sgn	SH/KM	DWG. NO.	DWG. NO. 1/1	
CHE. Sgn				
ENG. APP. Sgn				
USR. APP. Sgn				

5 Technical Data

5.1 Overview

Product Name : MiroCam Capsule Endoscope System

Model Name : MiroCam

Serial Number : Refer to the label which is attaching in the reverse side of the product.

5.2 Hardware

5.2.1 Capsule

- Size : 10.8 X 24mm
- Weight : 3.25g
- Material : Human Compliance Plastic
- Light : 6 white LED
- View Angle : 150°(in Image)
- View Depth : 3 cm
- Enlargement Ratio : 1:8
- Detectable Range : under 0.1mm
- Sampling Ratio : 2.9 fps
- Working time : Over 11 Hours
- Mechanical Safety : Compatible ISO60601-1-1
- Chemical Safety : Safe in pH=2 ~ pH=8
- Battery Type : Silver Oxide Cell(3Vdc, 70mA)
- Operation Temperature : 20 ~ 40
- Storage Temperature : 0 ~ 50

5.2.2 Receiver

- Operation System : Firmware
- Recording Time : 11 Hours
- Weight : 350g, include battery
- Operation Voltage : 3.7Vdc, 0.45A
- Battery Type: Lithium Ion Battery (3.7Vdc, 8.8A)
- Battery Weight : 215g
- Operation Temperature : 0 ~ 40
- Storage Temperature : 0 ~ 55
- Category : Type BF

5.2.3 Image Workstation

- Operating System : Windows XP Professional

- CPU : Core 2 Duo E6300(1.86GHz/2M)
- Memory : DDR II 1GByte(667MHz)
- Display Adaptor : Geforce 7600GT 256MB
- Hard Disk : SATA II 160GB, SATA II 320GB
- ODD : DVD-RW
- Monitor Resolution : 1280 X 1024
- Monitor Contrast : 700 : 1
- Printer Resolution : 4800 X 1200 dpi
- Printer Paper : A4

5.2.4 Cradle

- Input Voltage : 110~220VAC
- Input Current : 3A
- Output Voltage ; 4.2VDC
- Output Current : 4A
- Operation Display : LED Display

5.3 Software

5.3.1 Version

- MiroCam Application Software Version 1.00

5.3.2 Specification

- Language : English
- Data Export : JPEG Image, AVI Video Clip, PDF Data Report
- Data Display : Single or Multi Image, Time Bar, Color Bar, Diagnosis Data
- Event Marker : Small Image with Explanation
- Display Ratio : 5 ~ 30 fps
- Display Mode : Single View, Dual View, Quad View
- Running Mode : Normal Mode, Fast Mode

5.3.3 Feature

- Color status Display
- User Friendly software feature
- Automated detection of GI tract bleeding
- Color of images are virtually natural

5.4 COMPLIANCE/APPROVALS

The MiroCam capsule endoscope system and accessories complies with the Medical Device Directive 93/42/EEC (CE₀₈₄₃).

In addition, the product complies with

- IEC 60601-1:1988 + A1:1991 + A2:1995 (EN 60601-1:1990 + A1:1992 + A2:1995)
- EN 60601-1-2
- EN 10993-1
- EN 10993-10

The possibility of hazards arising from software errors was minimized in compliance with EN1441 and EN60601-1-4.

5.4.1 FCC STATEMENT

- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
- CAUTION : Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
- NOTE : This equipment has been tested and found to comply with the limit for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help

6 Environment condition for storage and transport.

6.1 Safety environment condition for storage

6.1.1 Receiver unit

- Environment condition for storage
 - Temperature : -10 - +70
 - Relative humidity : 10% - 80%
 - Atmospheric pressure : 700hPa to 1060hPa
- Keep in the place which is not contacted with water.
- Keep in the place where the direct light dose not shines.
- Keep in the place which is not contacted with something that is a hazardous to the human body.
- Keep in the place which the kid's hand dose not closes.
- Don't keep it in the (storing) place for chemicals or near gas.
- When you clean up the outside of the product, only use gauze with water. Other liquid like alcohol should not be allowed.
- When you charge the battery, please use the cradle.
- If you don't use this product for so long, please separate the battery from the body of it.

6.1.2 Image Workstation

- Treat this based upon the way how to handle a general PC.

6.2 Environment condition for transport

- Temperature : -10 - +70
- Relative humidity : 10% - 80%
- Atmospheric pressure : 700hPa to 1060hPa

7 Troubleshooting

7.1 Introduction

This section explains how to troubleshoot the MiroCam capsule endoscope system if problems arise. Tables are supplied that list of possible difficulties, probable cause, and recommended actions to correct the difficulties.

7.2 Who should perform repairs

Only qualified service personnel should open the HMT2100 housing, remove and replace components, or make adjustments.

7.3 Obtain replacement parts

IntroMedic Technical Services provides technical assistance information and replacement parts. To obtain replacement parts, contact IntroMedic Co.,Ltd. or your local representative. Refer to part names and part numbers listed on Section 3.2, *Supply Items*.

7.4 Troubleshooting guide

Condition	Recommended Action
1. The MiroCam capsule endoscope system workstation fails to power-up when the I/O power switch is on.	<ol style="list-style-type: none"> 1. Ensure that the MiroCam capsule endoscope workstation main system is plugged into an operational AC outlet in accordance with the input specification rated on the side panel of the workstation main system. 2. If the condition persists, contact IntroMedic Co.,Ltd..
2. The MiroCam capsule endoscope system receiver unit fails to power-up when the I/O power switch is on.	<ol style="list-style-type: none"> 1. Recharging battery of receiver unit. 2. Replacement battery of receiver unit. 3. If the condition persists, contact IntroMedic Co.,Ltd..
3. It's not blinking when it's out of the case.	<ol style="list-style-type: none"> 1. DO NOT use this capsule and retry another capsule. 2. Contact IntroMedic Co.,Ltd. for replacement capsule.
4. INI LED, on receiver unit, is lighted as yellow.	<ol style="list-style-type: none"> 1. Initialization receiver unit in MiroCam application software. 2. If the condition persists, contact IntroMedic Co.,Ltd..
5. BAT LED, on receiver unit, is lighted as yellow.	<ol style="list-style-type: none"> 1. Recharging battery of receiver unit. 2. Replacement battery of receiver unit. 3. If the condition persists, contact IntroMedic Co.,Ltd..