



Fig. 131

Cleaning bottle
Assembly

Order No.: 14 0514 57248



Fig. 132

Suction cups
1 set, consisting of 2 pcs.

Order No.: 14 3000 00403



Fig. 133

Storage rails for the unload drawer
1 set, consisting of 3 pcs.

Order No.: 14 0514 55967



Fig. 134

Cut-resistant gloves

1 pair, size M

Order No.:**14 0514 55967**

Fig. 135

Tool kit HistoCore SPECTRA CV

Consisting of:

- Screwdriver, 5.5x150
- Leica brush
- T16A fuse

Order No.:**14 0514 54189****Consumables**

Designation	Order No.
mounting medium	
HistoCore SPECTRA X1 (1 packaging unit, 4 bottles with 125 ml each)	380 1733
Coverglasses	
HistoCore SPECTRA CV premium coverglass 1x24x50 (2x300 pcs.)	380 0152

10. Warranty and Service

Warranty

Leica Biosystems Nussloch GmbH guarantees that the contractual product delivered has been subjected to a comprehensive quality control procedure based on the Leica in-house testing standards, and that the product is faultless and complies with all warranted technical specifications and/or agreed-upon characteristics.

The scope of the warranty is based on the contents of the contract concluded. The warranty terms of your Leica sales organization or the organization from which you have purchased the contractual product shall apply exclusively.

Service information

If you are in need of technical customer support or spare parts, please contact your Leica representative or the Leica dealer where you purchased the instrument.

The following information about the instrument is required:

- Model name and serial number of the instrument.
- Location of the instrument and name of a contact person.
- Reason for the service call.
- The delivery date of the instrument.

11 Decommissioning and Disposal

11. Decommissioning and Disposal



Warning

The instrument or parts of the instrument must be disposed of according to existing applicable, local regulations. All objects contaminated with spilled reagents must be disinfected immediately with a suitable disinfectant to prevent spreading to other areas of the lab or to lab personnel.

Please refer to the Chapter ([→ P. 109 – 7. Cleaning and Maintenance](#)) and the Chapter Decontamination Certificate ([→ P. 163 – 12. Decontamination Certificate](#)) at the end of these Instructions for Use for notes on cleaning the HistoCore SPECTRA CV coverslipper.

The instrument can be contaminated when using biohazardous specimens. Thorough disinfecting is required before re-commissioning or disposal (e.g. multiple cleaning steps, disinfection or sterilization). Dispose of the instrument in compliance with the applicable laboratory specifications.

Contact your Leica representative for more information.

12. Decontamination Certificate

Any product that is to be returned to Leica Biosystems or serviced on site must be properly cleaned and decontaminated. The associated decontamination certificate template can be found on our website www.LeicaBiosystems.com in the product menu. This template must be used to enter all required data.

If a product is returned, a copy of the completed and signed decontamination certificate must either be enclosed or handed over to a service technician. The user shall be responsible for products that are returned without a completed decontamination certificate or with a missing decontamination certificate. Return shipments that are classified by the company as a potential hazard source will be returned to the sender at his/her own expense and risk.

www.LeicaBiosystems.com



Leica Biosystems Nussloch GmbH
Heidelberger Str. 17 - 19
69226 Nussloch
Germany

Phone: +49 6224 - 143 0
Fax: +49 6224 - 143 268
Web: www.LeicaBiosystems.com