Implantable Insulin Pump System
Sistema de bomba de insulina implantable
Le système pompe à insuline implantable
Implantierbare Insulinpumpe
Pompa impiantabile per insulina
Implanteerbaar insulinepompstelsysteem van
Implanterbara insulinpumpssystem

options to be considered for diagnosing and restoring a failed or impaired MIP System.

Physicians must be completely familiar with the function of the MIP System prior to its implantation and use.

Patients must be fully trained on how to use the MIP System, and must demonstrate proficiency in programming their PPC and responding to alarms prior to using the MIP System on their own.

Individual components of the MIP System should NOT be used if damaged prior to or during the implantation procedure. The Pump should not be placed in direct contact with other metal implants.

When refilling the Pump, never push on the syringe plunger. When the Refill Needle is properly seated in the Pump fill port, vacuum in the Pump reservoir will automatically draw in insulin from the syringe. When refilling the Pump, never release the syringe-locking ring while there is a vacuum. The handle of the syringe may snap back forcibly and damage the syringe.

When measuring the Pump stroke volume, always check for leaks at the Pump - Pipette connection. If a leak occurs, the measurement is inaccurate and should be repeated.

ADVERSE EVENTS
Although the safety of the MIP System has been well demonstrated during clinical evaluations spanning more than 10 years, the following adverse events have been reported:

Table listing adverse events such as Abdominal Pain, Inflammation at Refill Site, Infection, Necrosis, Retinal Disorder, Abnormal Liver Function, Ileus, Hyperglycemia, and Hypoglycemia, along with their descriptions.

* A small number of patients may develop elevated titers of anti-insulin antibodies in their bloodstream. Clinical evaluations have shown that elevated anti-insulin antibodies alone do not interfere with the maintenance of glycemic control provided with this treatment.

STERILE DEVICE
The Pump and Infusion Catheter are sterilized by ethylene oxide and provided in tamper evident packaging. Do not use if the package has been opened or damaged.

WIRELESS COMMUNICATION NOTICE
This device complies with part 15 of U.S. FCC regulations and European standards for the wireless communication of medical devices.

Table showing regulatory information for the MIP System, including FCC and CE markings with corresponding numbers.

English

** NOTE **Please refer to the physician and patient manuals for detailed instructions for use and information regarding the Medtronic MiniMed Implantable Insulin Pump System

Table mapping MIP models (MMT-2007D, MMT-3150, etc.) to their respective components like Pump, Catheter, and Refill Kit.

DESCRIPTION
The Medtronic MiniMed Implantable Insulin Pump System ("MIP System") is a programmable medical device, designed to administer insulin for the treatment of Type 1 or insulin-requiring Type 2 diabetes mellitus.

INDICATIONS AND USAGE
The MIP System is indicated for the long-term, continuous, intraperitoneal administration of exogenous insulin, to maintain glycemic control in patients with Type 1 or insulin-requiring Type 2 diabetes mellitus.

- CONTRAINDICATIONS
* Patients unwilling or unable to:
- Test their blood glucose (BG) four or more times per day
- Use their BG results to update Pump programming
- Switch to subcutaneous insulin injections, if their Pump fails
- Comply with the guidance of their primary care physician
- Return to their physician every 2-3 months for insulin refills
* Patients who reside or travel above 8,000 ft. (2,400 m.) elevation
* Patients with a history of hypersensitivity reactions to titanium alloy or silicone materials that are used in the MIP System

WARNINGS
The Pump can only be used with the MMT-3150 PPC and either the MMT-4024 or MMT-4027 Infusion Catheter. The Pump can only be refilled or rinsed using the MMT-4102 Refill Needle and MMT-4105 Refill Kit.

PRECAUTIONS
Patients must always have supplies available for taking subcutaneous insulin injections, in the event of a MIP System malfunction. If a malfunction is suspected, patients should immediately discontinue the use of their MIP System, and contact their primary care physician for instructions and MIP System diagnosis.

Physicians should review the Physician's Manual for a full description of pump rinsing, catheter flushing and software reset

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