



Beyond The Basics: Performing Radial STEMI Procedures

Rajkumar K. Sugumaran, MD
Cardiac Solutions
Phoenix, AZ



Terumo Medical Corporation Disclaimer






-The opinions expressed by Dr. Sugumaran during this presentation are his own and do not necessarily reflect those of Terumo Medical Corporation. Unattributed data, device selection and procedural guidance is a matter of physician preference and opinion derived from Dr. Sugumaran's own observations and experiences and should be treated accordingly.

-All visual representations and other imagery contained in this presentation are used with permission from Dr. Sugumaran.



TIS-255-04082016

Objective


-  Why Transradial PCI for STEMI?
-  Perfusion Test and Patient Prep
-  Radial Access
-  Catheter Selection and Manipulation
-  Hemostasis



TIS-255-04082016


Vascular Access Complications

Hospital Resources Consumed in Treating Complications associated with Percutaneous Coronary Interventions



TERUMO TIS-255-04082016
 Kogutson, A, Cohen, D, Brown, MD, Silver, A, Becker, E, and Culler, S. American Journal of Cardiology 2010; 107: 202-207

Reduction in Bleeding Complications



*The Association of Arterial Access Site At Angioplasty With Transfusion And Mortality: The MORTAL Study. (Mortality benefit of Reduced Transfusion After PCI for the first 90 days). JACC: Cardio Interv. 2014; 8: 1165-1174. Authors: @Bharathini, W, Freer, S, Silver, A, Cohen, D, Culler, S
 †Effects of Radial Versus Femoral Artery Access in Patients With Acute Coronary Syndromes With or Without ST-Segment Elevation. Shante R, Mehta, MD, MD, JACC: Cardio Interv. 2014; 8: 1175-1184.
 ‡Radial Versus Femoral Randomized Investigation in ST-Segment Elevation Acute Coronary Syndrome: The RIFLE-STEACS (Radial Versus Femoral Randomized Investigation in ST-Elevation Acute Coronary Syndrome) Study. Ferrero Hernandez, MD, PhD, et al. European Society of Cardiology, 2014

TERUMO TIS-255-04082016

Radial Continues to be Supported by Evidence

Minimizing Adverse Haemorrhagic Events by Transradial Access Site and Systemic Implementation of Angiox (the MATRIX Trial)¹

Study Methods

- Randomised, Superiority trial
- 8404 Patients, 74 Centers
- Co-primary End Points
 - MACE
 - NACE
- Secondary End Points
 - Individual components of Composite outcomes
 - All cause Mortality
 - Stroke
 - MI
 - Bleeding

Results

- Radial garners superior outcomes
 - Radial as compared with femoral access reduces NACE through a reduction of bleeding and all-cause mortality
- Co-Primary End Points
 - 15% relative reduction in MACE
 - 17% relative reduction in NACE
- Secondary End Points
 - 28% Reduction in all-cause mortality
 - 33% reduction in Bleeding

TERUMO TIS-255-04082016
 Q100 2014.0n.17

Transradial PCI for STEMI

The Prevalence and Outcomes of Transradial PCI for STEMI

◆ Key Points:

- ◆ The authors concluded that the wider usage of TRI for STEMI may significantly improve patient outcomes.
- ◆ TRI patients were also significantly less likely to have vascular complications than the femoral PCI population.
- ◆ TRI was associated with a lower risk of bleeding and in-hospital mortality while there was no difference in procedural success.




TIS-255-04082016

Bellotti D, Galassi G, et al. J Am Coll Cardiol. 2013;61:2010-2016

Transradial and Bivalrudin

This figure demonstrates the rate of percutaneous coronary intervention-associated bleeding in 501,017 patients grouped by vascular access and anticoagulation.






TIS-255-04082016

Bellotti D, et al. Clin Cardiol. 2014; 36(3):157-162

Common Hurdles with Transradial PCI for STEMI





TIS-255-04082016

Transradial Advantages

- ◆ Benefits of radial access for all patients, especially STEMI
- ◆ When to start performing radial STEMI?
- ◆ Why to perform radial STEMI?
- ◆ Discussing with your Cath Lab Staff in advance of radial STEMI case



TIS-255-04082016

Room Set-up and Prep



TIS-255-04082016

Room Set up – Right Radial



Some physicians will access radial artery with arm at 90 degrees. Once sheath is inserted and secured, bring right arm in near right groin site.

Other physicians work in same plane as wrist with equipment easily accessible.



TIS-255-04082016

Room Set up – Left Radial

Begin with the arm adducted 90 degrees

To work from the patient's right, use pillows, straps or blankets to elevate arm

Assess position for patient comfort and compliance prior to prepping site



TIS-255-04082016

Room Set up – Left Radial



Access the left radial with arm positioned at 90 degrees to patient

Position the arm near patient's left groin if working from the patient's right



TIS-255-04082016

Room Set up – Considerations

Shave and prep medial to lateral; from the mid forearm to the mid palm

Shave the upper arm if RHC is to be performed

Consider shave and or Prep femoral access site.

Femoral access can prove useful for emergent access of the femoral artery or vein



TIS-255-04082016

Administering Heparin to Avoid RAO

Comparison of the Effect of Intra-Arterial Versus Intravenous Heparin on Radial Artery Occlusion After Transradial Catheterization

Key Points:

- ◆ No statistical difference between intra-arterial and systemic heparin administration in regards to RAO.
- ◆ 500 consecutive patients
- ◆ Early RAO 5.6% (ia) vs. 6% (systemic)
- ◆ Late RAO 4% (ia) vs 3.2% (systemic)




TIS-255-04082016


Small Procedure: Volume 136, Issue 8, Pages 1883-1891 (15 October 2008)

Radial Artery Spasm

Reasons



What Happens




Radial Artery Spasm occurs in 2%-6% of patients



TIS-255-04082016

Pain's Role of Transradial Intervention, The Basics and Beyond © 2012 by Tejas Patel

Radial Artery Spasm Prevention



TIS-255-04082016

Pain's Role of Transradial Intervention, The Basics and Beyond © 2012 by Tejas Patel

Antispasmodic Cocktail

Verapamil is a very acidic drug. To reduce burning effect, Physician should consider diluting cocktail with the patient's blood in a 10-20mL syringe. Other substitutions would be Nicardipene or Cardene




TIS-255-04082016

OPTITORQUE™

◆ Coronary diagnostic catheter

- ◆ Available in 5 Fr and 6 Fr catheter sizes
- ◆ Radial shapes are designed to eliminate catheter swaps
- ◆ Shaft with 2-ply stainless steel braid designed for 1-to-1 torque and accurate placement
- ◆ Large lumen for high contrast flow
- ◆ Atraumatic soft tip
- ◆ Designed to provide greater visibility around the ostium and lower contrast pressure from the end hole




Optitorque™
Coronary Diagnostic Catheters



TIS-255-04082016

HEARTRAIL™ Guide Catheters




- Provides multiple points of contact against the contralateral wall
- Ikari left offers versatility – for use in LCA and RCA
- Optimal STEMI guide cath to reduce DTB




TIS-255-04082016

CASE 1

◆ 57 y/o morbidly obese WF with a h/o CAD s/p Anterior MI LAD stents 2014, HTN and HLD presented to hospital for worsening chest pain. Troponin 0.28




TIS-255-04082016



TIG 4.0
CATHETER




TIS-255-04082016



IKARI 3.5 L
GUIDE
RUFATHROUGH®
NE Coronary
GuideWire




TIS-255-04082016




ASPIRATION THROMBECTOMY



TIS-255-04082016




TIS-255-04082016



EXTRA SUPPORT WITH/UP GUIDELNER



TIS-255-04082016



TIS-255-04082016

POLLING QUESTION

- ◆ Which is the strongest independent predictor for radial artery spasm?
- A. Female gender
- B. Diabetes
- C. Hypertension
- D. Small radial artery diameter
- E. Unsuccessful access at first attempt



TIS-255-04082016

Answer


- D. Small radial artery diameter



TIS-255-04082016

CASE 2

◆ 48 y/o Hispanic female with a h/o Diabetes mellitus, HTN, and HLD developed shortness of breath and crushing chest pain after dinner. Troponin 0.20




TIS-255-04082016




KARI 3.5 L
GUIDE
PUSH THROUGH[®]
NS Coronary
GuideWire




TIS-255-04082016



TIS-255-04082016




TIS-255-04082016



DES 3.0X16




TIS-255-04082016



IKARI 3.5L GUIDE



TIS-255-04082016



KARI 3.5 L
GLIDE
RATHROUGH®
NS Coronary
Guidance
PRIMARY
STENTING
DES 2.5 X12



TIS-255-04082016



3.0 X10 NC POST
DILATION



TIS-255-04082016



TIS-255-04082016

POLLING QUESTION

- ◆ Which of the following is a predictor of RAO (Radial arterial occlusion)?
- ◆ A. Female gender
- ◆ B. Length of procedure
- ◆ C. Lack of blood flow during compression D. Insufficient anticoagulation
- ◆ E. Ratio of artery diameter/ diameter of sheath < 1
- ◆ F. All of the above



TIS-255-04082016

Answer

- ◆ F. All of the above




TIS-255-04082016

Radial Artery Occlusion (RAO)

Remaining Challenges and Opportunities for Improvement in Percutaneous Transradial Coronary Procedures

Key Points:

- The authors grouped strategies aimed at minimizing the risk of arterial occlusion into three categories: Proven to Reduce Risk, May Reduce Risk, or Not Shown to Reduce Risk.
- The Four "Proven to Reduce Risk" strategies were:






TIS-255-04082016

Rao, S.V., Bernal, J., Barendse, O.F., Eur Heart J. 2012; Oct. 33(20):2501-6

Catheter Selection & Manipulation

 TIS-255-04082016

Radial Hemostasis

 TIS-255-04082016

Radial Hemostasis

Prevention of Radial Artery Occlusion – Patent Hemostasis Evaluation Trial (PROPHET Study)
A Randomized Comparison of Traditional Versus Patency Documented Hemostasis after Transradial Catheterization

Conclusion

Patent hemostasis is successful in significantly lowering the incidence of radial artery occlusion after TRA, without compromising hemostatic efficacy





Fig. 1. Patent hemostasis leads to a significant decrease in the incidence of radial artery occlusion at 24-hr and 30-day follow-up.

Radial Hemostasis

Impact of Two Different Hemostatic Devices on Radial Artery Outcomes after Transradial Catheterization

Conclusion

- TR BAND® Radial Compression Device provides equivalent hemostatic efficacy and a lower incidence of radial artery occlusion after transradial catheterization compared to the HemoBand.
- A device with a lower incidence of this complication is desirable over other available choices.



TIS-255-04082016

Furukawa T, J Invasive Cardiol 2009 Mar; 21(3): 101-4

Radial Hemostasis

At the conclusion of the procedure, a hemostasis device is recommended for access site management



TIS-255-04082016

Radial Hemostasis

- Important to understand the "Patent" hemostasis technique.
- Close attention to placement of the TR BAND® Radial Compression Device
- Monitor duration of compression



© 2014 TERUMO MEDICAL CORPORATION



TIS-255-04082016

Radial Hemostasis

Once hemostasis is achieved transport the patient with the inflation syringe attached either to the patients chart or taped to the patient.



TERUMO

TIS-255-04082016

Thank you! Questions?



TERUMO

TIS-255-04082016
