


Tour the Scope
















Overview

The scope is designed to view high-frequency oscillating data or single high-speed events that the graph display may not be able to capture.

Maximum Sampling Rate: The maximum supported sampling rate is a function of the computer-sensor interface that is connected to the computer (for example: 850 Universal Interface → 10 MHz, ScienceWorkshop 750 → 250 kHz).

Features	 See
Continuous Mode	"What is Continuous Mode?": 165
Fast Monitor Mode	"What is Fast Monitor Mode?": 166
Single shot	"What does Single Shot do?": 175
Multiple triggers	"How do I select the scope trigger?": 138
Pre-trigger	"How do I pre-trigger?": 141
Image freeze	"How do I freeze and record the Scope display?": 142

Toolbar

Icon	Description
	Adjust y-axis scale to fit data
	Activate and control scope trigger
	Stop collection after one trace
	Automatically adjust sample rate based on time-axis scale
	Activate to view multiple runs; Select visible run(s)
	Increase trace offset
	Set trace offset value to zero
	Decrease trace offset
	Show data coordinates and access Delta Tool
	Creates a data set from active traces
	Increase number of data points in trace
	Decrease number of data points in trace
	Add new y-axis to scope display
	Remove active element or axis
	Allow rearrangement of axes