




Medtronic

PEAK PlasmaBlade™ for Pacemaker/ICD Implants and Revisions



Conclusion

Use of the PEAK PlasmaBlade for pacemaker/ICD implants and revisions may reduce the likelihood of transvenous lead damage compared to traditional electrosurgery. Although the overall risk of transvenous lead injury is relatively low, the patient and financial consequences are serious and the use of additional measures to reduce potential risk, such as the PEAK PlasmaBlade, should be considered.



PEAK PlasmaBlade: Setting 5



Traditional ESU: 30 Watts

Copolymer lead, perpendicular device orientation, cut mode

Selected References*

1. Loh SA, Carlson GA, Chang EI et al. Comparative healing of surgical incisions created by the PEAK PlasmaBlade, conventional electrosurgery and a scalpel. *Plas Reconstr Surg* 2009; 124(6):1849-1859. Chronic wound healing study conducted in living porcine model.
2. Data on file. ETR-00023.
3. Ruidiaz ME, Messmer D, Huang EJ et al. Comparative healing of human cutaneous surgical incisions created by the PEAK PlasmaBlade, conventional electrosurgery, and a standard scalpel. *Plas Reconstr Surg* 2011; 128(1): 104-111. Chronic wound healing study conducted in subjects undergoing abdominoplasty.
4. Chang E, Carlson GA, Vose JG, et al. Comparative healing of rat fascia following incision with three surgical instruments. *J Surg Res* 2011; 167:e47-e54.
5. Data on file. VR-0055.
6. Weisberg IL, Desai S, Davison P et al. Effects of pulsed RF energy compared to standard electrosurgery on transvenous lead insulation materials. Heart Rhythm Society Annual Meeting 2010. Denver, Colorado.
7. Lim K, Reddy S, Desai S et al. Effects of electrocautery on transvenous lead insulation materials. *J Cardiovasc Electrophysiol* 2009; 20(4):429-435.
8. Reynolds MR, Cohen DJ, Kugelmass AD et al. The frequency and incremental cost of major complications among medicare beneficiaries receiving implantable cardioverter-defibrillators. *J Am Coll Cardiol* 2006; 47:2493-2497.

* Performance has not been specifically established in all surgical procedures.

Rx only. For a listing of indications, contraindications, precautions, and warnings, please refer to the Instructions For Use (IFU) that accompany PEAK PlasmaBlade disposable devices and/or the PEAK Surgery System User Guide.

For further information, please call 866-777-9400 or 603-742-1515.

You may also consult our website:

www.medtronicadvancedenergy.com

Medtronic Advanced Energy LLC

180 International Drive
Portsmouth, NH 03801
USA
www.medtronicadvancedenergy.com
Tel: (866) 777-9400
Fax: (866) 222-0900

International Telephone Numbers

Adriatic Regional Office 385-1-488-1120
Australia 1800-668-670
Baltic Regional Office 37-1-67560226
Belgium 32-2456-09-09
Canada 1800-217-1617
China 86-21-50800998
Czech Republic 420-2-9657-9580
France 33-470-679-800
Germany 49-2159-8149-209

Greece 30-210-67-79-099
Hong Kong 852-2919-1312
Hungary 36-30-5052987
India 91-22-26836733
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-6-4795-1506
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09

Netherlands 31-45-566-8800
Poland 48-22-465-6942
Russian Federation 7-495-580-73-77
Singapore 65-6776-6255
South Africa 27-11-466-1820
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
UK 44-1923-205-166
USA 1-603-742-1515

Significant Reduction in Transvenous Lead Damage When Compared to Traditional Electrosurgical Devices

The PEAK PlasmaBlade is a surgical device that uses very brief (40µs range), high frequency pulses of RF energy to induce electrical plasma along the edge of a thin (12.5µm), 99.5% insulated electrode. Due to the low duty cycle from RF pulsing and proprietary TPS insulating technology, the PEAK PlasmaBlade uses less total energy and operates at significantly lower temperatures than traditional electrosurgical technology (40 – 170°C vs. 200 – 350°C).^{1,2} Comparatively, PEAK PlasmaBlade incisions demonstrate 74% less thermal injury depth than traditional electrosurgical devices. (p<0.05).³

Pre-clinical and clinical studies with the PEAK PlasmaBlade have demonstrated that the improved thermal injury profile of the device, compared to traditional electrosurgery, results in the following benefits:

- Improves wound strength, wound healing and the cosmetic appearance of scars^{1,3,4}
- Decreases inflammatory cell counts³ and serous drainage⁵
- Enables patients to return to a normal diet earlier⁵
- Reduces patient intra-operative narcotic consumption by 22% (p=0.07)⁵
- Reduces patient post-operative narcotic consumption by 28% (p=0.59)⁵


Materials and Methods⁶

A series of ten polyurethane, silicone, and silicone-urethane copolymer transvenous leads were superficially tunneled into chicken breasts maintained at 37°C. These leads were then subjected to simulated surgical extraction using traditional electrosurgery or the PEAK PlasmaBlade. Extraction was performed with either parallel or perpendicular-to-lead technique using purely Cut or Coag mode at 3 second power outputs of either 20W or 30W. Lead damage was numerically characterized (0 to 3 scale, by severity) in a blinded fashion by visual and microscopic inspection.


Results⁶

Using traditional electrosurgery, significant lead damage was noted in all polyurethane leads, with more damage occurring with 30W vs. 20W, Cut vs. Coag mode, and perpendicular vs. parallel orientation. Considering alternative materials, silicone leads demonstrated less damage than polyurethane, and co-polymer leads demonstrated the greatest amount of damage of all three lead types with traditional technology. **Comparatively, the PEAK PlasmaBlade did not damage the silicone or polyurethane lead in Coag mode with either parallel or perpendicular technique. Using Cut mode, only minimal damage was demonstrated with perpendicular technique in the polyurethane and co-polymer leads.** Of the three insulation materials, silicone lead insulation demonstrated the highest tolerance to electrosurgery, regardless of technique or energy mode.


Thermal Injury Profile¹



Scalpel




PEAK PlasmaBlade Cut




Electrosurgery Cut

Device Effect on Transvenous Lead Insulation


PEAK PlasmaBlade




Traditional ESU



The PEAK PlasmaBlade: Pacemaker/ICD Implants and Revisions



PEAK PlasmaBlade: 40 - 170°C¹



Traditional ESU: 200 - 350°C¹

Discussion


Silicone, polyurethane, and copolymers are widely used in transvenous pacing leads due to their favorable flexibility, insulative, and tunneling characteristics. While they are readily available and cost-effective, polyurethane and copolymers in particular are also susceptible to damage from high-temperature electro-surgical instruments.^{6,7} Surgeons are advised by lead manufacturers to use low power settings and maintain operative vigilance during dissection to prevent damage. However, this

is not always possible. The consequences of lead damage during generator or battery replacement carry significant morbidity and mortality risk, including increased length of stay and death, and serious financial implications – averaging between \$5,000 and \$20,000 per incident.⁸ The PEAK PlasmaBlade's lower operating temperature and thermal spread resulted in reduced damage regardless of orientation or mode when compared to traditional electro-surgery.




Temperature Profile²

PEAK PlasmaBlade Cutting Temperature vs. Traditional ESU



■ The PEAK PlasmaBlade demonstrated an average 64% reduction in blade temperature compared to traditional electro-surgery for similar Cut settings

PEAK PlasmaBlade Coagulation Temperature vs. Traditional ESU



■ The PEAK PlasmaBlade demonstrated an average 40% reduction in blade temperature compared to traditional electro-surgery for similar Coag settings