**Manufacturer:**

Phonak AG  
Laubisrütistrasse 28  
CH-8712 Stäfa  
Switzerland

**EU Representative:**

Phonak GmbH  
D-70736 Fellbach-Oeffingen  
Germany

[www.phonak.com/phonak-roger](http://www.phonak.com/phonak-roger)

029-0265-02/NZ 00/2013-08/visu1 © Phonak AG All rights reserved

# Roger MyLink

## User Guide

roger

**PHONAK**  
life is on

---

## Contents

<b>1. Welcome</b>	<b>5</b>	<b>4. Using headphones</b>	<b>21</b>
<b>2. Getting to know your Roger MyLink</b>	<b>6</b>	<b>Step 1. Plug the headphones into Roger MyLink's headphone socket</b>	<b>21</b>
2.1 Compatibility	7	<b>Step 2. Turn Roger MyLink on</b>	<b>21</b>
2.2 Device description	7	<b>Step 3. Put the headphones on</b>	<b>22</b>
2.3 Indicator light	9	<b>5. More detailed information</b>	<b>23</b>
<b>3. Getting started</b>	<b>14</b>	<b>6. Troubleshooting</b>	<b>24</b>
<b>Step 1. Charge your Roger MyLink</b>	<b>14</b>	<b>7. Important information</b>	<b>26</b>
<b>Step 2. Detach the neckloop</b>	<b>16</b>	7.1 Hazard warnings	26
<b>Step 3. Hang Roger MyLink around the neck and reattach the loop</b>	<b>16</b>	7.2 Product safety information	28
<b>Step 4. Switch Roger MyLink on</b>	<b>17</b>	7.3 Other important information	31
<b>Step 5. Choose how to wear Roger MyLink</b>	<b>17</b>	7.4 International use	32
<b>Step 6. Turn on your Roger microphone</b>	<b>18</b>	<b>8. Service and warranty</b>	<b>33</b>
<b>Step 7. Switch your hearing aid(s) to the T-Coil program</b>	<b>18</b>	<b>9. Information and explanation of symbols</b>	<b>35</b>
<b>Step 8. Connect Roger MyLink to your Roger microphone</b>	<b>18</b>		
<b>Step 9. Adjust Roger MyLink's volume</b>	<b>20</b>		

---

## 1. Welcome

Congratulations on choosing Roger MyLink from Phonak. Roger MyLink is an easy-to-use and affordable Roger receiver suitable for use with any hearing aid featuring a T-Coil (from Phonak or any other manufacturer). Roger MyLink is compatible with most Roger microphones from Phonak.

Your Roger MyLink is a Swiss premium-quality product developed by Phonak, the world's leading hearing healthcare company.

Please read this User Guide thoroughly to benefit from all the possibilities your Roger MyLink offers.

If you have any questions, please consult your hearing care professional or local Phonak support center.

Phonak – life is on

[www.phonak.com](http://www.phonak.com)


---

## 2. Getting to know your Roger MyLink

Your hearing system consists of three main elements: Roger MyLink, your hearing aids and a Roger microphone, such as Roger Pen.

The Roger microphone transmits the speaker's voice wirelessly to your hearing aids via the Roger MyLink. You only need to activate the T-Coil in your hearing aid(s) to hear the speaker's voice loud and clear in your ear(s).

---


### 2.1 Compatibility

Your hearing aids must feature a T-Coil to receive speech signals from Roger MyLink. Please ask your hearing care professional if you are unsure whether or not your hearing aids have a T-Coil.

A compatible Roger microphone is also required to pick up the voice of the speaker. In this User Guide, Roger Pen is used as an example.

### 2.2 Device description

- ① Fixed side of neckloop
- ② Detachable side of neckloop
- ③ Green /orange status indicator

## 2. Getting to know your Roger MyLink

- ① Volume control
- ② On/off

- ① Charging socket
- ② Headphone socket (2.5 mm)

## 2.3 Indicator light

Green light

Orange light

- ① Battery almost exhausted
- ② Battery completely exhausted
- ③ Charging
- ④ Charging complete
- ⑤ Volume control: minimum/maximum position
- ⑥ Volume control: middle position

## 2. Getting to know your Roger MyLink

### General

#### Switching Roger MyLink on

The green light is visible for 3 seconds

### Meaning

Roger MyLink has been switched on and is ready for use.

### User action

Switch on the Roger microphone and set the hearing aids to the T-Coil program.

### Battery status

a)  60 seconds 

3 orange flashes (and 3 beeps), repeated every 60 seconds

Battery is almost empty.

Charge Roger MyLink as soon as possible (within an hour).\*

b) 

10 orange flashes (and 10 beeps)

Battery is completely exhausted. Roger MyLink can no longer be used without first being charged. Roger MyLink will automatically switch off after ten flashes.

Charge Roger MyLink immediately (within an hour). \*

\* We recommend you charge Roger MyLink and your microphone at the same time.

## 2. Getting to know your Roger MyLink

### Charging

c) 

Green double-flashes every 5 seconds

d) 

The green light is constantly visible

### Volume

e)  One orange flash (and a double beep)

f)  One green flash (and a single beep)

### Meaning

Roger MyLink is connected to the charger and is being charged.

Roger MyLink's battery is fully charged.

The volume control is set to either maximum or minimum position.

The volume control is in the middle position.

### User action

Wait until battery is charged. Charging a completely exhausted battery takes around two hours. The battery will be 80% charged within one hour.

Unplug the charger. The green light will disappear, meaning Roger MyLink is ready for use.

### 3. Getting started


#### Step 1. Charge your Roger MyLink

Fully charge Roger MyLink and your Roger microphone before using the system for the first time.

An empty Roger MyLink battery is:

- 80% charged after one hour
- 100% charged after two hours

Plug the charging unit into an easily accessible power socket.

Voltage: 7.5V  
Current: 250 Ma

⚠ Only use original Phonak charging units.

Plug the connector from the charger into the Roger MyLink charging socket.

- Once fully charged, Roger MyLink can be used for 10 or more hours.
- Roger MyLink cannot be used when the charger is still attached.


---

### 3. Getting started

#### Step 2. Detach the neckloop

#### Step 3. Hang Roger MyLink around the neck and reattach the loop


---

#### Step 4. Switch Roger MyLink on

#### Step 5. Choose how to wear Roger MyLink

---

### 3. Getting started

**Step 6.** Turn on your Roger microphone

**Step 7.** Switch your hearing aids to the T-Coil program

① It may be necessary to contact your hearing care professional to enable the T-Coil program in your hearing aids.

**Step 8.** Connect Roger MyLink to your Roger microphone

Hold Roger MyLink within 10 cm / 4 inches of your Roger microphone and press Connect on the Roger microphone.

---

① Step 8 is only required the first time you use Roger MyLink. Roger MyLink will then remain connected with your microphone even after Roger MyLink has been rebooted.

**Ready!**

### 3. Getting started

#### Step 9. Adjust Roger MyLink's volume

Volume up

Volume down

**Short button press:** The volume will increase/decrease only slightly.

**Long button press:** The volume will continue to increase/decrease until you stop pressing the button.

Once you have chosen your preferred volume, Roger MyLink will automatically store this setting for next time.

- ① If the wirelessly transmitted signal is not strong enough, increase the volume setting on Roger MyLink.

### 4. Using headphones

People who do not wear hearing aids can use Roger MyLink with headphones. To benefit from maximum hearing comfort, we recommend using official Roger MyLink headphones.

#### Step 1. Plug the headphones into Roger MyLink's headphone socket


#### Step 2. Turn Roger MyLink on

---

#### 4. Using headphones

##### Step 3. Put the headphones on

- The neckloop is automatically switched off when the headphones are attached (to save battery power).
- Use Roger MyLink's volume control to increase or decrease the volume of the transmitted signal (see page 20).

---

#### 5. More detailed information

Operating range: Up to 20 meters / 66 feet

Battery: Rechargeable Lithium-Polymer accumulator

Power supply: 7.5 VDC / 250 mA

Wireless: 2.4 GHz

## 6. Troubleshooting

The following is a summary of frequently asked questions about Roger MyLink and our answers to them.

<i>Most probable cause</i>	<i>Solution</i>
<b>The green light is not visible when Roger MyLink is switched on</b>	
■ Battery is completely exhausted	■ Charge Roger MyLink
<b>I cannot hear sound from the Roger microphone</b>	
■ One or all devices are switched off	■ Make sure all devices (Roger MyLink, Roger microphone, hearing aids) are turned on and fully operational
■ The T-Coil is not activated	■ Switch your hearing aids to the T-Coil program
■ Roger MyLink is not connected to your Roger microphone	■ Connect Roger MyLink with your Roger microphone by pressing 'Connect' on your microphone. Note that Roger MyLink needs to be within 10cm/4 inches of your Roger microphone
■ The Roger microphone is muted	■ Unmute the Roger microphone (refer to its User Guide)
■ Roger MyLink is out of range of the Roger microphone	■ Move closer to the Roger microphone to move within its operating range

<i>Most probable cause</i>	<i>Solution</i>
<b>The signal is interrupted during audio streaming</b>	
■ Roger MyLink is too far away from the Roger microphone	■ Move closer to the Roger microphone to move within its operating range
■ Roger MyLink is shielded by an obstacle (e.g. human body)	■ Adjust the device orientation and make sure you are in line of sight to the Roger microphone
<b>The sound of Roger MyLink is bad</b>	
■ Electromagnetic fields from other equipment may create an audible hum through the T-Coil of your hearing aids	■ Move away from electronic equipment until the hum stops
■ Roger microphone is not worn correctly	■ Make sure the Roger microphone is worn as indicated in its User Guide
<b>The volume is too loud or too soft</b>	
■ Volume is not correctly set	■ Adjust the volume on Roger MyLink until it suits your preference

---

## 7. Important information

- 📖 Please read the information on the following pages before using your Phonak Roger MyLink.

### 7.1 Hazard warnings

- ⚠️ Persons with pacemakers or other medical devices should ALWAYS consult their physician, the manufacturer of their device BEFORE using the Roger MyLink. The use of Roger MyLink with a pacemaker or other medical devices should ALWAYS be in accordance with the safety recommendations of the physician responsible for your pacemaker or the pacemaker manufacturer.
- ⚠️ Keep this device out of reach of children under 3 years of age.
- ⚠️ Switch Roger MyLink off when you enter a room with an activated loop system.
- ⚠️ Do not use water or fluids for cleaning. It is recommended to use Phonak CleanLine products. Your hearing care professional will be glad to advise you.

- ⚠️ Opening the Roger MyLink might damage it. If problems occur that cannot be resolved by following the guidelines in the Troubleshooting section of this User Guide, consult your hearing care professional.
- ⚠️ Your Roger MyLink should only be repaired by an authorized service center. Changes or modifications to this device that were not explicitly approved by Phonak are not permitted.
- ⚠️ Please dispose of electrical components in accordance with your local regulations.
- ⚠️ Do not use your Roger MyLink in areas where electronic equipment is prohibited. If in doubt, please ask the person in charge.
- ⚠️ Roger MyLink should not be used on aircrafts unless specifically permitted by flight personnel.
- ⚠️ External devices should only be connected if they have been tested in accordance with corresponding IECXXXXX standards.

---

## 7. Important information

- ⚠ Only use accessories that are approved by Phonak AG.
- ⚠ When operating machinery, ensure that no parts of your Roger MyLink can get caught in the machine.
- ⚠ Do not charge Roger MyLink while it is on your body.

### 7.2 Product safety information

- ⓘ Protect the Roger MyLink from excessive moisture (bathing, swimming), heat (radiator, car dashboard), and direct contact with skin when sweating (workout, fitness, sport).
- ⓘ X-ray radiation, CT or MRI scans may destroy or adversely affect the correct functioning of the Roger MyLink.
- ⓘ Protect the Roger MyLink from excessive shock and vibration.
- ⓘ Do not expose your Roger MyLink to temperatures and humidity outside of the recommended conditions stated in this User Guide.

- ⓘ Do not use excessive force when connecting your Roger MyLink to its different cables.
- ⓘ Protect all openings (microphones, audio input and charger) from dirt and debris.
- ⓘ If the Roger MyLink has been dropped or damaged, if it overheats during charging, has a damaged cord or plug, or has been dropped into liquid, stop using your Roger MyLink and contact your hearing care professional.
- ⓘ Never use a microwave or other heating devices to dry the Roger MyLink.
- ⓘ Clean the Roger MyLink using a damp cloth. Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the Roger MyLink.
- ⓘ Only recharge your Roger MyLink using chargers supplied by Phonak or by stabilized chargers providing 7.5VDC and  $\geq 250$  mA.

---

## 7. Important information

- ① When the Roger MyLink is not in use, turn it off and store it safely.
- ① The Roger MyLink has a built-in rechargeable battery. This battery can be charged and discharged hundreds of times but may eventually wear out. The battery may then only be replaced by authorized personnel.

---

## 7.3 Other important information

- High-powered electronic equipment, larger electronic installations and metallic structures may impair and significantly reduce the Roger MyLink's operating range.
- Roger MyLink collects and stores internal technical data. This data may be read by a hearing care professional in order to check the device, as well as to help you use your device correctly.
- The digital signal sent from the Roger microphone to a connected receiver cannot be overheard by other devices which are not in the transmitter's network.
- Be aware that wireless transmission from transmitter to receiver may be interrupted due to interference or moving out of range, and that words which are then spoken by the transmitter wearer may not be heard by the receiver wearer.

---

## 7. Important information

- Avoid exposing the receiver to hair spray or other cosmetic chemicals.
- The sensitivity of the hearing aid microphone may be reduced by the use of your Roger MyLink and therefore environmental awareness may be reduced.

### 7.4 International use

- The Roger MyLink is operating in the 2.4GHz ISM Band, which is unlicensed and worldwide available.

---

## 8. Service and warranty

### 8.1 Local warranty

Please ask the hearing care professional, where you purchased your Roger MyLink, about the terms of the local warranty.

### 8.2 International warranty

Phonak offers a one year limited international warranty, valid as of the date of purchase. This limited warranty covers manufacturing and material defects. The warranty is valid only if proof of purchase is shown.

The international warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer goods.

---

## 8. Service and warranty

### 8.3 Warranty limitation

This warranty does not cover damage from improper handling or care, exposure to chemicals, immersion in water or undue stress. Damage caused by third parties or non-authorized service centers renders the warranty null and void. This warranty does not cover any services performed by a hearing care professional in his/her office.

Serial number:

\_\_\_\_\_

Authorized hearing  
care professional  
(stamp/signature):

Date of purchase:

\_\_\_\_\_

\_\_\_\_\_

---

## 9. Information and explanation of symbols

With the CE symbol, Phonak AG confirms that this Phonak product meets the requirements of the Medical Devices Directive 93/42/EEC as well as the R&TTE Directive 1999/5/EC on radio and telecommunications equipment.  
**CE mark applied: 2013**

This symbol indicates that Roger MyLink products meet the requirements of a Type B component according to EN 60601-1.

This symbol indicates that it is important for the user to read and take into account the relevant information in this user guide.

---

## 9. Information and explanation of symbols

This symbol indicates that it is important for the user to pay attention to the relevant warning notices in these user guides.

Important information for handling and product safety.

Australian EMC and Radiocommunications compliance label.

---

### Operating conditions

This device is designed such that it functions without problems or restrictions if used as intended, unless otherwise noted in these user guides.

Temperature: -20° to +60° Celsius  
(-4° to +140° Fahrenheit)

Humidity transportation: Up to 90%  
(non condensing)  
Humidity storage: 0% to 70%, if not  
in use.

Atmospheric pressure: 500hPA to  
1100 hPA.

---

## 9. Information and explanation of symbols

The symbol with the crossed-out garbage bin is to make you aware that this device may not be thrown away as normal household waste. Please dispose of old or unused device, at waste disposal sites intended for electronic waste, or give your device to your hearing care professional for disposal. Proper disposal protects the environment and health.

This product has an embedded, non-replaceable battery. Do not attempt to open the product or remove the battery as this may cause injury and damage the product.

Please contact the local recycling facility for removal of the battery.

Important Safety Information  
Declaration of Conformity  
Service and Warranty

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**Section 3.0** Roger + FM microphone/transmitter, microphone/émetteur

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**Section 5.0** Roger SoundField loudspeaker, haut-parleurs

**Section 6.0** Service and Warranty, Service et Garantie

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Roger DynaMic

Roger AudioHub

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Roger MyLink, Roger 10, Roger 11, Roger 13, Roger 15, Roger 16, Roger X, Roger receiver Guide-U

### Roger SoundField loudspeaker, haut-parleurs

Roger DigiMaster 5000

Roger DigiMaster 7000

Roger DigiMaster X

The information contained herein is subject to change without notice.

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## Section 1.0 All devices

### 1.1 Important Safety Information

Please read the information on the following pages before using your Phonak device.

This symbol on the product or its packaging means that it should not be disposed of with your other household waste. It is your responsibility to dispose of your waste equipment separately from the municipal waste stream. The correct disposal of your old appliance will help prevent potential negative consequences for the environment and human health.

#### 1.1.1 Hazard warnings

-  Keep this device out of reach of children aged under 3 years.
-  Opening the device might damage it. If problems occur that cannot be resolved by following the guidelines in the Troubleshooting section of this User Guide, consult your hearing care professional.
-  Your device should only be repaired by an authorized service center. Changes or modifications to this device that were not explicitly approved by Phonak are not permitted.
-  Please dispose of electrical components in accordance with your local regulations.
-  Do not use your device in areas where electronic equipment is prohibited. If in doubt, please ask the person in charge.
-  device should not be used on aircrafts unless specifically permitted by flight personnel.
-  External devices should only be connected if they have been tested in accordance with corresponding IECXXXXX standards.
-  Only use accessories that are approved by Phonak AG.

### 1.1.2 Product safety information

- ① Protect the device from excessive shock and vibration.
- ① Do not expose your device to temperatures and humidity outside of the recommended conditions stated in this User Guide.
- ① Never use a microwave or other heating devices to dry the device.
- ① Clean the device using a damp cloth. Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the device.
- ① Only recharge your device using chargers supplied by Phonak or by stabilized chargers providing 5VDC and  $\geq 500$  mA.
- ① Do not connect an USB or audio cable to your device that exceeds 3 meters (9ft) in length.
- ① When the device is not in use, turn it OFF and store it safely.

### 1.1.3 Other important information

- Persons with pacemakers or other medical devices should ALWAYS consult their physician, the manufacturer of their device BEFORE using the device. The use of the device with a pacemaker or other medical devices should ALWAYS be in accordance with the safety recommendations of the physician responsible for your pacemaker or the pacemaker manufacturer.

## 1.2 Declaration of Conformity

For the devices used in body worn operation, the devices have been tested and found in compliance with the requirements of EN 62209-2, FCC rule paragraph 2.1093, the ANSI/IEEE C95.1:1999, the NCRP Report Number 86 for uncontrolled environment, according to the Health Canada's Safety Code 6 and the Industry Canada Radio Standards Specification RSS-102 for General Population/Uncontrolled exposure.

For the devices, which are not intended to be used in body worn operation, the antenna(s) used for these devices must be installed to provide a separation distance

of at least 20cm from all persons and must not be collocated or operating in conjunction with any other antenna or transmitter to comply with FCC and Industry Canada RF exposure limits for general population/uncontrolled exposure.

## Section 2.0 Roger microphone/transmitter

### 2.1 Hazard warnings

- ⚠ Using your device's cables in any way contradictory to their intended purpose (e.g., wearing the USB cable around the neck) can cause injury.
- ⚠ When operating machinery, ensure that no parts of your device can get caught in the machine.
- ⚠ Do not charge device while it is on your body.
- ⚠ Caution: electric shock. Do not insert plug alone into electrical outlets.

#### 2.1.1 Product safety information

- ⓘ Protect the device from excessive moisture (bathing, swimming), heat (radiator, car dashboard), and direct contact with skin when sweating (workout, fitness, sport).
- ⓘ X-ray radiation, CT or MRI scans may destroy or adversely affect the correct functioning of the device.
- ⓘ Do not use excessive force when connecting your device to its different cables.
- ⓘ Protect all openings (microphones, audio input and charger) from dirt and debris.
- ⓘ If the device has been dropped or damaged, if it overheats during charging, has a damaged cord or plug, or has been dropped into liquid, stop using your device and contact your hearing care professional.
- ⓘ The device has a built-in rechargeable battery. This battery can be charged and discharged hundreds of times but may eventually wear out. The battery may then only be replaced by authorized personnel.

### 2.1.2 Other important information

- Device collects and stores internal technical data. This data may be read by a hearing care professional in order to check the device, as well as to help you use your device correctly.
- The digital signal sent from the device to a connected receiver cannot be overheard by other devices which are not in the microphone's network.

### 2.2 Declaration of Conformity

Hereby, Phonak Communications AG, Länggasse 17, CH-3280 Murten, declares that these Phonak products are in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. To obtain a copy of the Declaration of conformity, please contact the manufacturer or the local Phonak representative whose address can be taken from the list on [www.phonak.com](http://www.phonak.com) (worldwide locations).

This device complies with Part 15 of the FCC Rules and with RSS210 of Industry Canada. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. The device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. The device generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.

- 
- 
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help.

**FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.**

---

## Section 3.0 Roger + FM microphone/transmitter

### 3.1. Hazard warnings

- 
-  Using your device's cables in any way contradictory to their intended purpose (e.g., wearing the USB cable around the neck) can cause injury.
  -  When operating machinery, ensure that no parts of your device can get caught in the machine.
  -  Do not charge device while it is on your body.
  -  Caution: electric shock. Do not insert plug alone into electrical outlets.
- 

#### 3.1.1 Product safety information

-  Protect the device from excessive moisture (bathing, swimming), heat (radiator, car dashboard), and direct contact with skin when sweating (workout, fitness, sport).
  -  X-ray radiation, CT or MRI scans may destroy or adversely affect the correct functioning of the device.
  -  Do not use excessive force when connecting your device to its different cables.
  -  Protect all openings (microphones, audio input and charger) from dirt and debris.
  -  If the device has been dropped or damaged, if it overheats during charging, has a damaged cord or plug, or has been dropped into liquid, stop using your device and contact your hearing care professional.
-

- ① The device has a built-in rechargeable battery. This battery can be charged and discharged hundreds of times but may eventually wear out. The battery may then only be replaced by authorized personnel.

### 3.1.2 Other important information

- Device collects and stores internal technical data. This data may be read by a hearing care professional in order to check the device, as well as to help you use your device correctly.
- Your device is approved to operate at a radio frequency that is specific to your country or region and might not be approved for use outside your country or region. Be aware that operation during international travel may cause interference with other electronic devices, or other electronic devices may cause interference with your device.
- Be aware that your device's FM radio signals might also be picked up and overheard by other FM receivers.

### 3.2 Declaration of Conformity

Hereby, Phonak Communications AG, Länggasse 17, CH-3280 Murten, declares that these Phonak products are in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. To obtain a copy of the Declaration of conformity, please contact the manufacturer or the local Phonak representative whose address can be taken from the list on [www.phonak.com](http://www.phonak.com) (worldwide locations).

This device complies with Part 15 and part 95 of the FCC Rules and with RSS210 of Industry Canada. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. The device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. The device generates uses and can radiate radio frequency energy and, if

not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.**

This transmitter is authorized by rule under the Low Power Radio Service (47 CFR part 95) and must not cause harmful interference to TV reception or United States Navy SPASUR installations. You do not need an FCC license to operate this transmitter. This transmitter may only be used to provide: auditory assistance to persons with disabilities, persons who require language translation, or persons in educational settings; health care services to the ill; law enforcement tracking services under agreement with a law enforcement agency; or automated maritime telecommunications system (AMTS) network control communications.

Two way voice communications and all other types of uses not mentioned above are expressly prohibited. This device may not interfere with TV reception or Federal Government radar. For Canada: if TV channel 13 is used in the area, the installer shall reduce or adjust the radio frequency radiated power so that nearby TV channel 13 receivers do not receive radio interference from the system installed. Suggestions: A test with a TV receiver equipped with a "rabbit-ear antenna" and tuned to channel 13 should be conducted at the perimeter of the user's intended coverage

area and should not overlap other users' areas without the latter's consent. If this does not solve the problem, a channel near the 217 MHz edge and not near 216 MHz should be tried. This Class B digital apparatus complies with Canadian ICES-003.

---

## Section 4.0 Roger receivers

### 4.1 Product safety information

- ① Protect the device from excessive moisture (bathing, swimming), heat (radiator, car dashboard), and direct contact with skin when sweating (workout, fitness, sport).
- ① X-ray radiation, CT or MRI scans may destroy or adversely affect the correct functioning of the device.

#### 4.1.1 Other important information

- The digital signal sent from the device to a connected receiver cannot be overheard by other devices which are not in the microphone's network.
- Be aware that wireless transmission from microphone to receiver may be interrupted due to interference or moving out of range, and that words which are then spoken by the microphone wearer may not be heard by the receiver wearer.
- Avoid exposing the receiver to hair spray or other cosmetic chemicals.
- The sensitivity of the hearing aid microphone may be reduced by the use of your device and therefore environmental awareness may be reduced.

Universal receivers may sound a loud beep at start-up, which might be annoying for cochlear implant users. Please contact your hearing care professional to discuss deactivating this beep function.

### 4.2 Declaration of Conformity

Roger receivers are products of Phonak AG, which have been tested and confirmed to be in compliance with the essential requirements and other relevant provisions of the Medical Devices Directive 93/42/EEC and Directive 1999/5/EC. To obtain a copy

of the Declaration of conformity, please contact the manufacturer or the local Phonak representative whose address can be taken from the list on [www.phonak.com](http://www.phonak.com) (worldwide locations).

This device complies with Part 15 of the FCC Rules and with RSS210 of Industry Canada. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. The device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. The device generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.**

---

## Section 5.0 Roger SoundField loudspeaker

### 5.1 Other important information

- The digital signal sent from the device to a connected receiver cannot be overheard by other devices which are not in the microphone's network.

### 5.2 Declaration of Conformity

Hereby, Phonak Communications AG, Länggasse 17, CH-3280 Murten, declares that these Phonak products are in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. To obtain a copy of the Declaration of conformity, please contact the manufacturer or the local Phonak representative whose address can be taken from the list on [www.phonak.com](http://www.phonak.com) (worldwide locations).

This device complies with Part 15 of the FCC Rules and with RSS210 of Industry Canada. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. The device generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**FCC Caution:** Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.

---

## 6.0 Service and Warranty

### Local warranty

Please ask the hearing care professional, where you purchased your device, about the terms of the local warranty.

### International warranty

Phonak offers a one year limited international warranty, valid as of the date of purchase. This limited warranty covers manufacturing and material defects. The warranty is valid only if proof of purchase is shown.

### Warranty limitation

This warranty does not cover damage from improper handling or care, exposure to chemicals, immersion in water or undue stress. Damage caused by third parties or non-authorized service centers renders the warranty null and void. This warranty does not cover any services performed by a hearing care professional in his/her office.

Serial number: \_\_\_\_\_

Purchase date: \_\_\_\_\_

Authorized hearing  
care professional  
(stamp/signature): \_\_\_\_\_

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029-3149/V1.00/2013-04/visu/ Printed in Vietnam © Phonak AG All rights reserved