


Rondish Company Limited
Unit G & H, 4/F, Block 1, Kwai Tak Ind. Ctr,
15-33 Kwai Tak St., Kwai Chung, N.T.,
Hong Kong

Early Alert Bed Sensor Pad (WPAS-10) User Instructions

Sensor Pad and Monitor Instructions

Select Pad Setting:

- 1 Open the pocket to access the 3-way switch on the transmitter module. The default position is OFF
- 2 "TX-only" will activate Early-alert function only on the transmitter side of the sensor pad. Alarm goes off when a patient is sitting on the "TX" side edge of the sensor pad and about to leave their bed. The other side remains in a normal setting. Alarm goes off when the patient leaves their bed.
"Both-sides" will enable Early-Alert function on both sides of the sensor pad.




Pad with Monitor:

- 1 Press the "Program" button on the wireless monitor unit – you will hear a beeping from the monitor
(NOTE: for CordFree monitor press Program + Reset)
- 2 Press and release the centre zone of sensor pad to add it into memory. You will see a visual and audible feedback from the monitor to confirm the signal was received
- 3 Test that the sensor pad is working by pressing and releasing



Secure the Pad:

- 1 Position the sensor pad across the bed beneath the hips of the patient. The neck of the sensor pad should exit out the side of the bed.
- 2 Optional clip sets are available to prevent pad movement. They loop through the eyelets as shown below
- 3 They can be clipped to the bed from to hold the sensor pad in place.



Care & Maintenance Instructions:

- Use disinfectant wipes or antibacterial cleansers to clean the pad
- **DO NOT** immerse the sensor pad in any liquid
- **DO NOT** expose the sensor pad to high temperature sanitizers

Storage:

- Always keep the sensor pad flat when in storage
- Operating temperature: 10°C ~ 40°C (50°F ~ 104°F)
- Storage temperature: 0°C ~ 50°C (32°F ~ 122°F)

Rondish Company Limited
Unit G & H, 4/F, Block 1, Kwai Tak Ind. Ctr,
15-33 Kwai Tak St., Kwai Chung, N.T.,
Hong Kong

FCC & CE warning statements

Any changes or modifications not expressly approved by the party responsible for compliance could void your authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

RF exposure statements

1. This Transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. This equipment complies with FCC radiation exposure requirement set forth for an uncontrolled environment. End user must follow the specific operating instructions for satisfying RF exposure compliance.

Hereby, Rondish Company Limited declares that the radio equipment type, WPAS-10 is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: www.rondish.com

RONDISH COMPANY LIMITED

Unit G/H, 4/F, Block 1, Kwai Tak Industrial Centre, 15-33 Kwai Tak St., Kwai Chung, N.T., Hong Kong.

Tel.: (852) 2541 6685 Fax: (852) 2541 7411 Email: rcf@rondish.com Website: www.rondish.com