ReliOn®

HEMOGLOBIN A1c COLLECTION INSTRUCTIONS

Experts in Lab Results by Mail

PLEASE READ THOROUGHLY

Remember to complete all your personal information on the Hemoglobin A1c Test Authorization Form. Please read and provide your signature and date under the Patient Consent in order to approve your sample for testing. Heritage Labs will not complete any testing without your signature and date.

Toll-free ReliOn customer support:
1-888-764-2384

A1c Importance

Daily plasma glucose monitoring is important; though additionally a more long-term picture is necessary in order to assess the effectiveness of your diabetes treatment plan. The A1c test shows your average plasma glucose level over a two to three month period. By testing these levels regularly and using the results, along with your daily plasma glucose records, you can determine how well your diabetes is being controlled.

Count on A1c to help you stay in control

The American Diabetes Association (ADA) recommends that people with diabetes monitor glycemic control with an A1c test every two to three months. Also recommended are an annual test of microalbumin levels to verify proper kidney function and a cholesterol test, including LDL and HDL levels, to help determine whether you are at risk for heart disease or stroke.

Use only as indicated

The ReliOn® collection card is intended to collect blood for analysis of A1c. The Appraise collection kit is mailed to Heritage Labs for analysis and the results are released to you in the manner and address as directed by you on the Hemoglobin A1c Test Authorization Form.


1. Check The Contents Of Your Kit

Your Appraise kit contains:

1. Collection Instructions
2. Hemoglobin A1c Test Authorization & Collection Form
3. Lancet
4. Alcohol Pad
5. Gauze Pad
6. Adhesive Bandage
7. Postage Paid Return Envelope


2. Fill Out The Hemoglobin A1c Test Authorization Form

Fill out the Hemoglobin A1c Test Authorization Form and complete the required information.

3. Rinse Hands

Rinse hands in warm tap water. Clean the selected puncture site with an alcohol pad and dry it with a gauze pad.


4. Stimulate Blood Flow

Stimulate blood flow to the selected finger by letting your hand hang at your side for 15-20 seconds. Shake your hand back and forth several times.


5. Use The Lancet Provided

Using the lancet provided, remove the cap, place palm up and position lancet on finger. Press down on lancet to puncture site. Wipe off the first blood drop with the alcohol pad.


6. Blood Should Begin To Flow

Blood should begin to flow freely. Place a **large** free-flowing drop of blood on each circle of the Collection Form as indicated. Do not place one drop of blood on top of the other. If blood flow stops, wipe with alcohol pad again to assist blood flow.

7. Let The Blood Spot Air-Dry

Let the blood spots air-dry for at least 30 minutes. Make sure specimen is dry prior to folding Hemoglobin A1c Test Authorization Form.

8. Tear Off Top Portion Of The Test Authorization Form

Tear off top portion of the Hemoglobin A1c Test Authorization Form. Retain for your records. Your specimen number will be necessary to retrieve results electronically.

9. Mailing Instructions

1. Place the Hemoglobin A1c Test Authorization Form into the postage paid return envelope.
2. Mail the sample within **three days** of completion of the collection.
3. Results will be released to you in the manner as directed by you on the Hemoglobin A1c Test Authorization Form.

Good daily plasma glucose control makes for good A1c results

People with frequent high plasma glucose readings may have a high A1c test result. To maintain an A1c level of less than 7%, your daily plasma glucose should rarely go above 170 mg/dL in any given two to three month period. Use the chart below to understand how your daily plasma glucose results relate to A1c levels.

Comparing Plasma Glucose with A1c

Average Glucose, mg/dL	% A1c
100	5
135	6
170	7
205	8
240	9
275	10
310	11
345	12

Use A1c to help reduce the risk of complications

A1c is a type of hemoglobin that interacts with plasma glucose. The A1c test gives you a picture of your average plasma glucose level over the previous two to three months and can be used to predict your risk for complications. The higher the A1c result, the higher your risk.


A1c measurements should be a regular part of your diabetes care.

ReliOn® delivers accurate results without the hassle

With ReliOn®, you can collect a sample anytime and anywhere. What's more, the accuracy of the ReliOn® A1c test averages 98% or more when compared with standardized laboratory methods using whole blood collected in a tube. As with any lab test, improper specimen collection may affect results.

Please review your test results with your doctor or health care professional.

Dried Blood Spot vs. Whole Blood Hemoglobin A1c

Store and handle your ReliOn® kit properly

The ReliOn® collection kit may be stored at room temperature (65°-78° F). When properly stored, the kit is stable until the expiration date printed on your ReliOn® kit. Do not use the ReliOn® kit past the expiration date.

Limitations

As with any laboratory test, there are factors that can affect test results. Hemoglobin A1c (A1c) values should be used in conjunction with information from other clinical and medical evaluations. A1c is not a diagnosis for detecting diabetes. Improper collection of specimens may affect laboratory results. If you have any questions regarding the performance of A1c testing, please consult your doctor.

A word of warning

Persons with hemophilia or anti-coagulation therapy should consult their doctor or health care professional prior to using this kit.

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