- PATIENT PROGRAMMER
USER MANUAL
- PATIENT PROGRAMMER
GEBRUIKERSHANDLEIDING
- PATIENTEN-PROGRAMMIERGERÄT
GEBRAUCHSANLEITUNG
- PROGRAMMATEUR PATIENT
MANUEL D'UTILISATION
- PROGRAMADOR DEL PACIENTE
MANUAL DEL USUARIO

PATIENT PROGRAMMER
USER MANUAL

LP0563-01 Rev A, 2014-02

Patient Programmer

Model MN20600-02

User Manual

Copyright © 2014 Spinal Modulation, Inc. All rights reserved. Unauthorized duplication or use is prohibited. AXIUM, SPINAL MODULATION and the Man in Circle Design are trademarks and/or registered trademarks of Spinal Modulation, Inc.


Authorization to affix the CE mark
granted in 2011

Contents

Explanation of Symbols on Product or Package Labeling	1
Glossary	2
Introduction	3
Indications for Use.....	3
Description.....	3
Warnings.....	3
Precautions – For Your Programmer and Your Stimulator.....	4
Precautions – For Your Therapy.....	5
Programmer Overview.....	7
Programmer Features	7
Charging the Battery.....	8
Programmer Power Up	8
Main Menu.....	9
Stimulator Binding	10
Connect	10
Magnet.....	10
Connecting with Your Stimulator	10
Pain Control Screen	11
ID Heading.....	11
Select Group	11
Back to Main Menu	11
Programmer Status Bar.....	12
Adjusting Your Stimulator Settings.....	12
Turn OFF All Stimulation.....	12
Pain Control Tab	13
Turn Stimulation On or Off for a Body Region.....	13
Adjust Stimulation to a Body Region	13
Device, Physician and Clinic Information.....	14
My Info Tab	14
Device Tab and Physician Tab.....	14
Stimulator Identification Information	14
Physician Information	14
Clinic Tab	15

To Use with Your TNS Device	15
To Use with Your INS Device.....	16
Cleaning your Programmer	16
Disposal of your Programmer	16
Environmental Conditions for Storage and Operation.....	16
Appendix I: Troubleshooting.....	20

Explanation of Symbols on Product or Package Labeling

REF	Model Number
	Serial Number
	Read the Manual
	Consult the Manual
	Contents of Package are Non-Sterile
	Manufacturing Date
	Manufacturer
	Protected against Electric Shock
	Not waterproof. Applies to the Programmer when it is not in its carrying case.
	Limited waterproof. Applies to the TNS. Applies to the Programmer in its carrying case.
	Turns the Programmer ON and OFF. Turns stimulation OFF on the TNS.
	Keep Dry
	Store between -10°C and 50°C (14°F and 122°F)
	Store between 0 and 93% humidity
	The device is a radio transmitter
	Magnet. Shows the location of the Programmer magnet.
	Authorized European Representative

Glossary

Lead – Surgical wire: takes electrical signals from the neurostimulator to the stimulation area

Stimulation – Small electrical pulses: produces a tingling sensation and replaces pain signals

Stimulator – Device that makes electrical pulses that will stimulate the nerves in your spine: can refer to either the Trial Neurostimulator or Implantable Neurostimulator

Trial Neurostimulator (TNS) – External Stimulator that clips onto your belt: attaches to the connector cable, which is connected to the leads that are implanted in the area near your spine

Implantable Neurostimulator (INS) – Stimulator implanted in your back or abdomen: attaches to leads implanted in the area near your spine

Connector Cable – Cable that connects the leads to your Trial Neurostimulator

Programmer – Portable, hand-held device: allows you to adjust your stimulation settings

Clinical Programmer – Portable, hand-held device: allows the physician to program your Stimulator.

Computer Tomography (CT) Imaging – Computerized X-ray imaging: produces electronic images of tissues and organs

Diathermy – High energy heat: used to cut or cauterize during surgery or a type of therapy

Electromagnetic Interference (EMI) – Electrical signals that interfere with the device function

Magnetic Resonance Imaging (MRI) – Medical imaging: produces electronic images of tissues and organs

Paresthesia – Tingling sensation felt during therapy delivery: produced by spinal cord stimulation

Precaution – Situation that could cause uncomfortable stimulation and possible damage to the Stimulator or Programmer

Program – Instructions or changes to stimulation settings that are put into the Programmer and transmitted to the Stimulator

Stimulation Level – Amount of stimulation: can be increased or decreased within a range set up by your doctor

Warning – Potentially serious hazard that could cause injury or death

Introduction

Your Axium® Patient Programmer is used to program your Trial or Implantable Stimulator. This User Manual gives detailed instructions on how to safely use your Programmer and your Stimulator. See your doctor if you have any questions.

Indications for Use

The Spinal Modulation Neurostimulator System is indicated for the management of chronic intractable pain.

Description

The Axium Patient Programmer is a portable, hand-held device. It is used to communicate with your Stimulator. Your Programmer contains information about your Stimulator. You will use it to adjust your stimulation settings or turn stimulation off. It can run on its rechargeable battery. It can also be plugged into a power outlet.

A Clinical Programmer is used by your doctor to initially set up your Stimulator and make adjustments later if needed. Only your doctor and/or Spinal Modulation clinical personnel use the Clinical Programmer.

Warnings

- Do not use your Programmer or the Stimulator until your doctor has trained you.
- Do not use your Programmer until your doctor has set up your Stimulator.
- Do not have any elective magnetic resonance imaging (MRI). If MRI is necessary, your doctor should remove any leads and disconnect the TNS or INS device. Using MRI in the area of the leads may move the leads. It can also damage the TNS or INS device. Using MRI may cause uncomfortable ("jolting" or "shocking") levels of stimulation. MRI may cause heating at the tip of the lead. Unintended stimulation could result in tissue damage.
- Before having a CT scan, tell your doctor that you have an implanted device. All stimulation for your device should be turned OFF before the procedure. After the scan, your doctor should turn it back on and make sure the system is working properly.
- Do not have any diathermy (high energy heat) procedures. Diathermy may cause injury to you or damage your Stimulator.
- Do not remove your leads or Connector Cable by yourself. This may cause serious injury and could cause an infection.
- Do not open or modify the Programmer or TNS. Keep them closed to protect them. Modifications to the device may cause improper operation.
- Do not transport the Programmer outside of its carrying case. Operate it only in a moisture-free environment. The Programmer may malfunction if it becomes wet.

- Magnetic fields may interfere with the Stimulator. This can affect the stimulation level and cause discomfort. Avoid theft detection devices at store and library exits. Avoid security screeners at airports. Inform security personnel of your devices. They will assist you. Do not stand near the screening equipment.
- Avoid areas with warning notices preventing entry by patients wearing a pacemaker. Ask your doctor for more information.
- Power generators, arc welders and large magnetized speakers may cause interference. Do not stand near these or similar devices.
- Be aware of where you place your Charger. Pets, children or you can become entangled in the cord, which could cause a fall or strangulation.
- If contact with the Stimulator System causes a rash, report this to your doctor. If your throat or tongue starts to swell, get emergency aid immediately.

Precautions – For Your Programmer and Your Stimulator

Follow these precautions to maintain proper function of your Programmer and Stimulator.

- Do not drop or mishandle your Programmer or Stimulator. Physical damage to the devices may keep them from working properly.
- Do not wash the Programmer or TNS device. Excess water may keep them from working properly. Use a soft damp cloth to gently wipe the devices if needed.
- Do not use abrasive or caustic cleaning products on your Programmer or TNS device.
- Avoid contact with body fluids for the TNS and Programmer. Contamination may cause damage to the devices.
- Do not shower or bathe with the TNS device. You may take a sponge bath. You must take care not to get the TNS device wet.
- Do not use any equipment or accessories that are not supplied with your Programmer.
- Do not place your Programmer close to cards with magnetic strips. The Programmer has a magnet that could demagnetize your cards. Keep the Programmer away from computer hard drives and magnetic storage devices.
- Do not operate the Programmer or Stimulator outside the temperature range of 5°C to 40°C (23°F to 113°F). Rapid temperature changes may affect device operation.
- Do not store the Programmer outside the temperature range of -10°C to 50°C (14°F to 122°F).
- Do not leave the Programmer in a car or other places where temperatures can reach 50°C (122°F).
- An unlikely failure of your Stimulator System is possible due to random component failure. If any part of your Stimulator System stops working or changes how it works, turn stimulation OFF. Contact your doctor during normal business hours.
- Return your Programmer and your TNS to your doctor at the end of the trial period, or when no longer being used. Do not discard or burn the TNS or Programmer. Fire may cause the internal batteries to explode.
- Do not attempt to dispose of the TNS or Programmer yourself.
- Do not replace the TNS or Programmer battery by yourself, even if it does not seem to be working. Only Spinal Modulation personnel should replace the TNS or Programmer batteries.

- Do not use any other company's device to program your Stimulator. Use only the Programmer provided by Spinal Modulation.
- Do not allow unauthorized use of your Programmer. This may cause unwanted changes in the programming.
- Do not use the Programmer or Stimulator near explosive or flammable gases. This may cause serious injury.
- Do not use the Charger if the power cord is damaged. This may cause injury or damage your Stimulator.
- To remove power from the Charger when not in use, unplug it from the power outlet.
- Avoid unnecessary programming of your device. Frequent usage will wear the battery down faster.
- Do not open or shred the battery.
- Do not short-circuit the battery.
- Do not subject the battery to mechanical shock.
- In the event of a battery leak, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- In case of a non-responding unit, Spinal Modulation personnel will observe the plus (+) and minus (-) marks on the battery and equipment and ensure correct use. They will also wipe the battery terminals with a clean dry cloth if they become dirty.
- Keep the battery out of the reach of children.
- Seek medical advice immediately if a battery has been swallowed.
- Keep the Programmer clean and dry to prevent affecting the battery.

Precautions – For Your Therapy

Follow these precautions to maintain appropriate therapy:

- Follow the proper wound care techniques given to you by your doctor.
- Do not rub or press on the implant site. This may cause the leads to move or your skin to erode. It may also cause the INS to move.
- Wear the TNS on the outside of your clothing or on a belt.
- Changes in your body position can increase pain or cause uncomfortable stimulation. Use the Programmer to adjust stimulation levels. Turn OFF stimulation, if needed.
- Avoid excessive bending, twisting and stretching. Do not lift objects over ten pounds. These activities may cause the leads to move. You may experience higher or lower levels of stimulation.
- Avoid driving a car or operating other potentially dangerous machinery while stimulation is ON. You could be distracted from driving or device operation if sudden changes in stimulation occur.
- Your Stimulator may affect the operation of other implantable devices such as pacemakers or implantable cardiac defibrillators. Tell your doctor about any other implantable devices you have or are scheduled to get.

- Tell your regular doctors or healthcare providers that you have a Stimulator. Do not have any elective medical procedures without first discussing them with your doctor. Some medical devices or therapies, such as those listed below, may interfere with your Stimulator:
 - Electrocautery – Uses an electric probe to cauterize blood vessels and stop bleeding during surgery.
 - Lithotripsy – Uses high-output shock waves to break up gallstones and kidney stones.
 - Therapeutic Radiation – Uses ionizing radiation to destroy cancer cells.
 - High-output ultrasound – Uses high frequency sound waves to treat bone and muscle injuries, or to stimulate muscle or improve blood flow.
 - RF Ablation – Uses radio frequency energy to cause controlled tissue damage.
 - Microwave Ablation – Uses high speed alternating electric field to cause controlled tissue damage.
 - Dental procedures, electrolysis, static field therapeutic magnets and diagnostic X-ray.
- Appoint a family member or friend to tell emergency medical personnel that you have a Stimulator, in case you need emergency care. You will be given a Medical Alert Card to carry with you. This card will inform emergency medical personnel that you have a Stimulator.

If you have any concerns about your stimulator device, contact your doctor during normal business hours.

RF OPERATING FREQUENCIES

Nearby equipment emitting strong magnetic fields can interfere with RF communication, even if the other equipment complies with CISPR emission requirements. The operating characteristics are as follows:

MICS band: 402-405 MHz

The effective radiated power is below the limits as specified in

Europe: EN ETSI 301 839-2

USA FCC 47 CFR Part 95; 95.601-95.673 Subpart E, 95.1201-95.1219

FCC ID: Y8L-MN20600-02

This device may not interfere with stations operating in the 400.150–406.000 MHz band in the Meteorological Aids, Meteorological Satellite, and Earth Exploration Satellite Services and must accept any interference received, including interference that may cause undesired operation.

Programmer Overview

Your Axium Patient Programmer is a portable, hand-held device. It is powered by an internal rechargeable battery. The Programmer can also be plugged into a power outlet for use or for recharging. Your Programmer communicates with your Stimulator to control your stimulation. Your doctor will explain how to use the Programmer to adjust stimulation for your best pain relief.

Keep your Programmer near you at all times. This will allow stimulation adjustment, if needed. Carry your Programmer in its carrying case. The case provides protection from water.

Your Programmer System consists of:

- Programmer (and internal magnet)
- Auxiliary Magnet
- Carrying Case
- Patient Information Sheet
- Stylus
- Charger
- Patient Medical Alert Card
- Patient Programmer User Manual (this document)

Programmer

Stylus

Programmer Charger

Programmer Features

With your Axium Patient Programmer, you can:

- Turn stimulation ON or OFF for each body region
- Adjust the stimulation level for each body region
- Change the Group for stimulation - See "Select Group" under the "Pain Control Screen" section
- Turn OFF all stimulation if needed
- View your Stimulator ID information
- View your ID number
- View your lead implant date
- View your doctor's name, clinic name and contact information

Charging the Battery

The Programmer System comes with a Charger. The Programmer must be fully charged before using it for the first time. It takes approximately two to four hours to fully charge the battery. The Programmer Status Bar at the bottom of the screen shows the battery charge level.

1. Plug the Charger into a power outlet.
2. Connect the Charger to your Programmer.

When the battery is charging, the battery icon on the screen shows "AC". When charging is complete, the amber light becomes green.

Your Programmer will operate when it is connected to a power outlet. It does not use battery power when connected to an outlet. Connect your Programmer to the Charger and plug it in to an outlet regularly to keep it charged. Do not leave Charger plugged into a power outlet overnight.

NOTE: Use only the Model MN23400-U Charger. It was included with your Programmer. Use of any other charger could damage your Programmer.

The battery can be expected to last at least 500 discharge cycles with normal use. Once it fails to hold a charge, it should be replaced by Spinal Modulation personnel only. Do not attempt to change the battery. Your Programmer and Charger can be expected to last up to two years with normal use.

Programmer Power Up

Press the  button to turn ON your Programmer screen. The Main Menu will display.

NOTE: If your Programmer does not turn ON, charge the battery and try again.

The device is charging when the amber light is on.

Press here to turn on Programmer screen.

Main Menu

The Programmer Main Menu displays two main functions:

- **Connect:** Allows you to connect to your Stimulator; also allows you to adjust stimulation settings.
- **Programmer Setup:** Allows you to set your Programmer date and time, and to view information about your Stimulator.

The Main Menu also shows your physician, and the clinic phone number.

The Programmer status bar is at the bottom of the Main Menu. It displays your Programmer–Stimulator connection status, the battery charge level and the time. See the “Programmer Status Bar” section in this User Manual for more detail.

You can change the time and date and access the Programmer Info screen from the Setup screen. You can also view your Stimulator serial number and your patient ID.

The Programmer Info screen displays your programmer serial number, programmer software version, firmware version and manufacturing date.

Stimulator Binding


Your doctor will bind and unbind your Stimulator to your Programmer. You cannot edit this information.

Connect

Use the “Connect” button on the Main Menu to connect with your Stimulator and to adjust your Stimulator settings. See the “Connecting with Your Stimulator” section below.

Magnet

A magnet is located under the magnet symbol on the back side of the Programmer. Place the magnet over the Stimulator to check for connection between the Programmer and Stimulator. See the “Connecting with Your Stimulator” section below for more detail.

Connecting with Your Stimulator

Use your Programmer to communicate with your Stimulator.

- Turn ON your Programmer. The Main Menu will display. Your Programmer must be within three feet of your Stimulator.
- Press the “Connect” button. The Programmer will begin searching for the Stimulator. An icon shows on the screen to show that it is busy.
- Hold the magnet on the Programmer over your Stimulator and move it around in a circular motion to start communication.
- The Programmer chimes when it is connected to your Stimulator. The Pain Control screen will then display. “Connected” shows in the status bar at the bottom left of the screen. If your Programmer cannot connect to your Stimulator, an error message displays. “Disconnected” shows in the status bar.
- If your Programmer cannot connect to your Stimulator, go back to the Main Menu. Press “Connect” again. Move the magnet in a circular fashion over your Stimulator. Repeat until the Programmer connects to the Stimulator.


NOTE: *If after two minutes your Programmer cannot connect to the Stimulator, place the Programmer over the Stimulator again. An error message may display. The message asks you to confirm that your Programmer is near enough to the Stimulator. After confirming that your Programmer is within three feet of the Stimulator, press “OK”. Press “Connect” again.*

Pain Control Screen

The ID Heading is at the top of the Pain Control Screen. The “Turn OFF All Stimulation” button is just below the ID Heading.

Two tabs are below the “Turn OFF All Stimulation” button: the “Pain Control” tab and the “My Info” (Information) tab. See the “Adjusting Your Stimulator Settings” section in this User Manual for more detail.

The “Exit” button at the bottom right side of the screen returns you to the Main Menu.

ID Heading

ID: J052122 **Stimulator CB0708**

Located at the top of the screen, the ID Heading displays the following information:

- **ID:** Displays your identification number (ID).
- **Stimulator Serial Number:** Displays your Stimulator's serial number.


Select Group

The “Group” button is in the center of the Pain Control screen. Press “Group” to display a drop down menu. The drop down menu has up to four groups defined by your doctor. When you select a group name, the Stimulator switches settings to the new group.

Back to Main Menu

The “Exit” button closes the Pain Control window. The session ends, and returns to the Main Menu.

NOTE: When programming is complete, select the “Exit” button. Turn off the programmer to conserve power.

Programmer Status Bar

The Programmer Status Bar is located at the bottom of the Programmer screen. The status bar shows the Programmer-Stimulator connection status, the battery charge level and the time.

- **Programmer-Stimulator Connection Status:**
 - Shows "Connecting" when the Programmer is trying to connect to the Stimulator.
 - Shows "Connected" when the Programmer is connected to the Stimulator.
 - Shows "Disconnected" when the Programmer is disconnected from the Stimulator.
- **Battery Level:** Shows the Programmer battery charge level.
- **Clock:** Shows the time. See the "Main Menu" section in this User Manual for more detail.


Adjusting Your Stimulator Settings

You can adjust your Stimulator settings from the Pain Control screen. Stimulation can be turned ON or OFF for up to four regions of your body. You can also adjust the stimulation level for any of those regions.

Turn OFF All Stimulation

- Press the "Turn OFF All Stimulation" button to stop all stimulation therapy. A window appears, asking you to confirm that you want to turn OFF all stimulation.

NOTE: After turning OFF all stimulation, you can restore stimulation therapy for each of the body regions individually. See "Turn Stimulation On or Off for a Body Region" section below.

Pain Control Tab

Select the "Pain Control" tab on the "Pain Control" screen. From the "Pain Control" tab, you can turn stimulation ON or OFF for each body region. You can also adjust the stimulation level for each body region.

Turn Stimulation On or Off for a Body Region

Your Programmer shows the names of one to four designated body regions that your leads affect. To turn stimulation ON or OFF for a body region:

- Select the body region by pressing the desired tab.
- Press the "OFF" button to stop stimulation to that region. When stimulation is OFF, the "OFF" button is black.
- Press the "ON" button to start stimulation to that region. When stimulation is ON, the "ON" button is green.

Adjust Stimulation to a Body Region

Select the correct body region tab on the Pain Control screen. Be sure that "Enable" is "ON".


- Press the "-" button to decrease the stimulation level.
- Press the "+" button to increase the stimulation level.

Stimulation Level Indicator:

The stimulation level indicator is between the "-" and "+" buttons. The indicator moves up or down as you adjust the stimulation level for the selected body region. The indicator shows the current stimulation level as compared to the maximum set by your doctor.

NOTE: The indicator bar is completely green when you have reached the maximum stimulation level.

NOTE: Be sure that "Enable" is ON when adjusting the Stimulation Level. If "Enable" is OFF, the "+" button will not respond.



Device, Physician and Clinic Information

My Info Tab

The "My Info" (Information) tab contains three tabs, the "Device" tab, the "Physician" tab, and the "Clinic" tab.

Device Tab and Physician Tab

The "Device" tab and "Physician" tab display the following information:

Stimulator Identification Information

- Stimulator Firmware Version
- Stimulator Voltage Information

NOTE: The battery information is for an INS device . It is not for an external TNS device.

- Implant Date: The date the INS device was implanted. This information is not for the external TNS device.


Physician Information

- Name: Your doctor's name
- Phone: Your doctor's contact phone number
- Email: Your doctor's email contact

Clinic Tab

"Clinic" tab displays the following information:

- **Clinic Name:** Your clinic's name
- **Address:** Your clinic's address
- **After Hours Contact:** A phone number to call in case of an emergency

To Use with Your TNS Device

To connect your Programmer to your TNS, push the "Connect" button on the Programmer. Move the Programmer magnet over the TNS in a circular motion. The Programmer will chime when connection is made. You may then use the Programmer to adjust your stimulation settings.

To quickly turn OFF stimulation, press the red button on the TNS for more than two seconds or push the "Turn OFF all Stimulation" button on the Programmer screen. To enable stimulation after pressing either button, you must connect with your Programmer and turn the stimulation back on.

To Use with Your INS Device

To connect your Programmer to your INS, push the "Connect" button on the Programmer. Move the Programmer magnet over the implant location in a circular motion. The Programmer will chime when connection is made. You may then use the Programmer to adjust your stimulation settings.

Cleaning your Programmer

To clean your Programmer, wipe carefully with a damp – but not wet – soft cloth. Do not get the Programmer wet. Always follow the Precautions listed above. They instruct you not to wash the Programmer or TNS. Excess water may keep them from working properly.

Disposal of your Programmer

Return your Programmer and your TNS to your doctor at the end of the trial period or when no longer being used. Do not discard or burn the TNS or Programmer. Fire may cause the internal batteries to explode. Do not attempt to dispose of the TNS or Programmer yourself.

Environmental Conditions for Storage and Operation

Storage Temperature/Humidity: -10°C to 50°C (14°F to 122°F) at relative humidity up to 93%

Storage Altitude: 700 hPa to 1060 hPa

Operating Temperature/Humidity: 5°C to 40°C at relative humidity from 15% to 93%

Operating Altitude: 700 hPa to 1060 hPa

Guidance and Manufacturer's Declarations

GUIDANCE AND MANUFACTURER'S DECLARATION <i>Electromagnetic Emissions</i>		
The Spinal Modulation Neurostimulator System is intended for use in the electromagnetic environment specified below. The customer or the user of the Spinal Modulation Neurostimulator System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic Environment – Guidance
RF Emissions 1	Group 2	The Spinal Modulation Neurostimulator System must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	
		The Spinal Modulation Neurostimulator System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
CISPR 14-1	Complies	The Patient Programmer is not intended to be connected to other equipment except the Programmer Charger.

GUIDANCE AND MANUFACTURER'S DECLARATION
Electromagnetic Emissions

The Spinal Modulation Neurostimulator System is intended for use in the electromagnetic environment specified below. The customer or the user of the Spinal Modulation Neurostimulator System should assure that it is used in such an environment.

Immunity	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD)	IEC 61000-4-2	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Pass	Mains power quality should be that of a typical commercial or home environment
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth		Mains power quality should be that of a typical commercial or home environment
Voltage dips, short interruptions and voltage variations on power supply	input lines IEC 61000-4-11 <5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s NOTE UT is the a.c. mains voltage prior to application of the test level.		Mains power quality should be that of a typical commercial or home environment
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m		Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial, hospital, or home environment.

GUIDANCE AND MANUFACTURER'S DECLARATION
Electromagnetic Immunity

The Spinal Modulation Neurostimulation System is intended for use in the electromagnetic environment specified below. The customer or the user of the Spinal Modulation Neurostimulation System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of Spinal Modulation Neurostimulation System, than 0.2 meter, based on transmitters of 80 MHz to 2.5 GHz. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Recommended separation distances between portable and mobile RF communications equipment and the Spinal Modulation Neurostimulation System

The Spinal Modulation Neurostimulation System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Spinal Modulation Neurostimulation System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the System.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12m	0.12m	0.23m
0.1	0.37m	0.37m	0.74m
1	1.17m	1.17m	2.33m
10	3.70m	3.70m	7.37m
100	11.70m	11.70m	23.30m

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Appendix I: Troubleshooting

Pop Up Message	Possible Solution
<p>Connection with your Stimulator was lost. Please reconnect.</p>	<ul style="list-style-type: none"> • Move your Programmer over your Stimulator to establish connection and reconnect to the device. • If you still have difficulty connecting, bring your Programmer closer to your Stimulator and continue to move the magnet over your Stimulator. • Move to another location, as there may be interference in your current location.
<p>Unable to connect to your Stimulator. Please contact your physician during normal business hours.</p>	<ul style="list-style-type: none"> • Hit "OK" and try to reconnect to your Stimulator. • Move your Programmer over your Stimulator to establish connection and reconnect to the device. • If you still have difficulty connecting, bring your Programmer closer to your Stimulator and continue to move the magnet over your Stimulator. • Move to another location as there may be interference in your current location. • Contact your doctor during normal business hours if the problem does not go away.
<p>Unable to connect to your Stimulator. Please try again.</p>	<ul style="list-style-type: none"> • Move your Programmer over your Stimulator to establish connection and reconnect to the device. • If you still have difficulty connecting, bring your Programmer closer to your Stimulator and continue to move the magnet over your Stimulator. • Move to another location, as there may be interference in your current location.
<p>Your Stimulator battery is low. It will need to be replaced soon. Please contact your physician during normal business hours.</p>	<ul style="list-style-type: none"> • Contact your doctor during normal business hours to set up an appointment. Your battery has reached the Elective Replacement Interval (ERI).

Pop Up Message	Possible Solution
<p>Your Stimulator battery needs to be replaced. Stimulation has been turned OFF permanently. Please contact your physician during normal business hours.</p>	<ul style="list-style-type: none"> • Contact your doctor during normal business hours to set up an appointment. Your battery has reached End of Service (EOS) and will not stimulate until it has been replaced.
<p>Stimulation for one or more leads has been turned OFF. Please contact your physician during normal business hours.</p>	<ul style="list-style-type: none"> • Move your Programmer over your Stimulator to establish connection and reconnect to the device. • Go to the Pain Control screen and turn on the lead that has been turned off. • If it does not allow you to turn it back on, contact your doctor during normal business hours.
<p>All stimulation has been turned OFF. Please contact your physician during normal business hours.</p>	<ul style="list-style-type: none"> • Contact your doctor during normal business hours.
<p>Stimulation has been turned OFF due to a magnet. Please use your Programmer to restore stimulation.</p>	<ul style="list-style-type: none"> • Move your Programmer over your Stimulator to establish connection and reconnect to the device. • Go to the Pain Control screen and turn on each lead one at a time.
<p>You have turned OFF all Stimulation. Please use your Programmer to restore stimulation.</p>	<p>You have turned off the device by pressing the switch on the TNS or by pushing the "Turn OFF All Stimulation" button on your Programmer.</p> <ul style="list-style-type: none"> • Move your Programmer over your Stimulator to establish connection and reconnect to the device. • Go to the Pain Control screen and turn on the leads that have been turned off.
<p>Programmer battery is low. Please recharge.</p>	<ul style="list-style-type: none"> • The battery has reached 30% on your Programmer and needs to be recharged.

Spinal Modulation hereby declares that Axium Neurostimulators and Programmers are in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. The Declaration of Conformity may be obtained upon request from the Manufacturer. Manufacturer contact information may be found on the back cover of this manual.

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Copyright © 2014 Spinal Modulation, Inc. Alle rechten voorbehouden. Onbevoegde duplicatie of gebruik is verboden. AXIUM, SPINAL MODULATION en het Man-in-Circle-ontwerp zijn handelsmerken en/of gedeponeerde handelsmerken van Spinal Modulation, Inc.
Copyright © 2014 Spinal Modulation, Inc. Alle Rechte vorbehalten. Die unbefugte Vervielfältigung oder Verwendung ist verboten. AXIUM, SPINAL MODULATION und das Mensch-in-Kreis-Design sind Marken und/oder eingetragene Marken von Spinal Modulation, Inc.
Copyright © 2014 Spinal Modulation, Inc. Tous droits réservés. Toute reproduction ou utilisation non autorisée est interdite. AXIUM, SPINAL MODULATION et le logo Man in Circle sont des marques et/ou des marques déposées de Spinal Modulation, Inc.
Copyright © 2014 Spinal Modulation, Inc. Todos los derechos reservados. Queda prohibido el uso o duplicado no autorizado. AXIUM, SPINAL MODULATION y el diseño del hombre en el círculo son marcas comerciales o marcas comerciales registradas de Spinal Modulation, Inc.
Consultants B.V.**
Spinal Modulation Pty Ltd