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Base de Communication E₉₀

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E₉₀ Communication Base User manual

Copyright © 2006 Télémédic Inc. All rights reserved.

This equipment is protected by the Canadian copyright laws. It cannot be reproduced, distributed or modified in any ways by other entities than Télémédic, except under written consent by Télémédic.

All the requisite efforts have been done to ensure that the information contained in this document was complete and precise at the moment of printing. However, the information can be modified without notice.

Declaration of conformity

This device is registered as CSA International certified. To see the certificates of official declaration, consult the online documentation on Télémédic's website: <http://www.telemedic.ca>

This device complies with Industry Canada's low power radio devices norms RSS-210. Certification number: 6418A-E90001

This device complies with the FCC (Federal Communications Commission) intentional radio frequencies transmission.

Identification number: FCC ID : T3W-E90-001

Warranty

Consult the warranty booklet included with this product to learn about the limited warranty given by Télémédic on that product.

Trademark

Télémédic, Télémédic's logo and VieGie Med are registered trademarks. Other brands and products names mentioned in this document are properties of the respective holders.

Ordering information

You can order updates of this document directly on Télémédic's website: <http://www.telemedic.ca>.

Comments

Télémédic appreciates your comments relating to this manual. Please send your comments directly to info@telemedic.ca.

Support

For all technical questions or defectiveness, please contact Télémédic's user support service at (450) 250-1414 or at info@telemedic.ca

Meaning of symbols used in this manual:

	Warning: Failure to these instructions may lead to unexpected operation or defects of the system.
	Important: Important element to consider in the use of this device.

Meaning of device symbols:

	Warning: Consult accompanying documents. This manual and the associated sections on Télémédic's website (see manual for references) should be consulted before using this device.
	This equipment complies with general requirements for safety of medical electrical equipment CAN/CSA C22.2 No.601.1 and CAN/CSA C22.2 No.60601-1-1, UL 2601-1. This equipment also complies with electromagnetic compatibility requirements for medical electrical equipment CAN/CSA C22.2 No 60601-1-2
	This device complies with FCC (Federal Communications Commission) rules part 15 (Radio frequency devices) for the United States.
	This equipment includes radio frequency transmitter which generate non-ionizing radiation.

The E₉₀ communication base

The E₉₀ communication base is a device that acts as a relay between medical peripherals (M₉₃ portable unit or device equipped with the I_{9X} interface card) and the main server located at the hospital or health professional's office, depending on the case.

This device, located at the patient's home, receives data by wireless radio frequency, then ensures their periodic transmission. This transmission is made by telephone connection or by network connection, such as the Internet.

It is highly recommended to read this manual before using the E₉₀ communication base.

In this section, we describe:

- the main components;
- the device's functions and terms of use;
- the precautionary measures;
- the maintenance;
- the functional preventive check-ups;
- the optimal operating performance.

Main components

E90 communication base Indicators

Indicator No.

- 1 On/off switch
- 2 Forced server communication switch
- 3 Power input connector
- 4 RS-232 connector (reserved for technical staff use)
- 5 RJ-45 network socket (LAN)
- 6 RJ-11 telephone socket
- 7 Power indicator LED (green)

Indicator No.

- 8 Radio frequency communication indicator LED (yellow)
- 9 Main server communication indicator LED (yellow)
- 10 Green indicator LED showing that internal modem is operational
- 11 Green indicator LED showing that internal network card is operational
- 12 Yellow indicator LED reserved for technical staff use
- 13 Wireless telecommunication antenna

Illustration 2 E90 communication base components

The device functions and-terms of use

Turning on power

Connect the power supply provided by TéléMédic (**TéléMédic Part no 31-004 / Manufacturier CUI Inc. # DPR090080-P5P-SZ**) to the E₉₀ communication base power input connector (**illustr. 1 – no. 3**). The power supply should then be connected to a standard electrical outlet (120V, 60Hz).

Activate the power switch (**illustr. 1 – no. 1**). The indicator lights 7, 8, 9 and 12 of **illustration 1** should flash and turn off some time later, indicating that the device is now ready to receive and send data. The indicator 7 of **illustration 1** stays on to indicate that power is on. Notice that this indicator flashes when supplying is made from the internal batteries of the E₉₀ communication base. In this case, unless an electrical failure occurs, the connection of the power connector should be checked. Notice that it is recommended to plug all cords before turning power on.

Vital signs receipt from M₉₃ portable unit or from the I_{9x} interface card

Sending vital signs towards the E₉₀ communication base is done automatically at predetermined hours so you don't have to intervene.

The yellow indicator flashes (**illustr. 1 - indicator 8**) during receipt of information by radio frequency.

Transmission

To use the transmission by telephone connection, connect one end of a standard telephone cord in the telephone inlet of the E₉₀ communication base (**illustr. 1 – no. 6**) and the other end in the telephone wall socket. Note that you can use a telephone socket multiplier if the socket is already used by another phone. During a telephone call, indicator 10 goes on.

In the case of a network transmission, such as by Internet, you have to connect the end of a standard network cable in the network socket (**illustr. 1 - no 5**) on the E₉₀ communication base and the other end in an available socket of your networking equipment (router, modem, concentrator, etc.). Notice that the network has to dynamically assign a network address (DHCP) and has to allow the communication on the TCP/IP port selected; this information should be provided by your supplier in order to keep the device working properly through this communication path. During a network communication, indicator 11 goes on.

You can simultaneously connect both the network socket and the telephone socket. In this case, the E₉₀ communication base will first try to make the transmission by network and will automatically use the telephone line as a second resort; for example, in the case of a network breakdown.

The E₉₀ communication base is configured to transmit data to the server in accordance with predetermined hours or frequencies. If necessary, you can press the switch located at the back of the communication base (**illustr. 1 – no. 2**) to force a transmission to the server.

When the E₉₀ communication base tries to contact the server, the yellow indicator light (**illustr. 1 - indicator 9**) goes on. The same indicator flashes when data is transferred to the server.

The yellow indicator (**illustr. 1 - indicator 12**) is reserved to the use of TéléMédic's technical personnel. If this indicator flashes, please contact TéléMédic's technical support (refer to the **Support** section at the end of the present document).

Using the telephone

If you are using your phone and the E₉₀ communication base tries to use your telephone line, you will hear a dialling tone in your handset. In this case, it is not necessary to end your call: the E₉₀ communication base will try to do its transfer later, without losing information.

If the E₉₀ communication base uses the telephone line and you pick up the handset, you won't hear no tone and it will be impossible for you to make a phone call. You just have to wait until the E₉₀ communication base terminates its transfer, indicator 9 will be off.

If you have high speed Internet ADSL connection (on a telephone line), make sure to connect your E₉₀ communication base through a telephone line filter provided by your Internet service provider.

TéléMédic recommends that you use a telephone socket multiplier (T adapter), telephone cable (modular line cord) and networking cable CSA approved.

Precautionary measures

The E₉₀ communication base is a precise scientific device that requires some basic precautions during use. Following the advice listed below will ensure your communication base's longevity, durability and optimal operating performance.

Never:

- open or dismantle the device;
- hit or throw;
- subject the device to sudden temperature changes;
- put the device in an area exposed to high temperature, high moisture or dusty contents;
- clean the device with a vaporizing liquid, solvent or detergent;
- immerse the device.

Maintenance

- Use a soft lightly damp cloth.

Functional preventive check-ups

To ensure your E₉₀ communication base's optimal operating performance, do the following check-ups:

- The power switch (**illustr. 1** no. 1) has to be up;
- Indicator 7 of **illustration 1** has to be on;
- Make sure the power supply is connected to a standard electrical outlet (120V, 60Hz) and that the plug is connected to the E₉₀ communication power supply input connector (**illustr. 1** – no. 3).

Optimal operating performance

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. This medical equipment needs special precautions regarding electromagnetic compatibility (EMC). Installation and services shall take care of following precautions:

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

The use of telephone cord or network cable other than those provided or specified by Télémédic may affect the electromagnetic compatibility.

For installation adjacent, stacked or very close to electrical equipment, it is important to verify operation of either product. For E₉₀ Communication base refer to functional preventive check-ups above. For the other equipment refer to the manufacturer operating manual.

Also, portable and mobile RF communication equipment ca affects the E₉₀ communication base.

Technical support

If you detect or suspect any malfunction or failure of the E₉₀ communication base, or if you have questions about its installation or use, please contact Télémédic's technical support at **(450) 250-1414**.

Technical specifications

- Radio communication: Transmission power = 750µW (-1.25 dBm)
Transmission range = 100 metres/328 ft.
Frequency = 916MHz (frequency band ISM)
TéléMédic's proprietary protocol
- Standby generation: 3 rechargeable internal batteries AA, 1.2V, Ni-MH
Autonomy = 24 hours
 Batteries are not intended to be changed by the operator.
- Power consumption: 200mA
- Memory capacity: 512Ko up to 40000 vital signs values (device for 3 vital signs)
- Data retention: 40 years
- Inputs/Outputs: Serial = connector RS232 D-SUB DB9
Telephonic = connector RJ-11
Network = connector RJ-45
- Dimensions: 12.5 cm x 18 cm x 3.6 cm / 5 in x 7 in x 1.4 in
- Weight: 362 grams/12.77 oz.
- Temperature: Operation = 10-35 °C (50-95 °F)
Storage = 5-40 °C (41-104 °F)
- Humidity: Operation = 20-80% / Storage = 10-90%
- Security: CAN/CSA C22.2 No.601.1, CAN/CSA C22.2 No.60601-1-1 et UL 2601-1.
- Configurable frequency and period of data transmission
- Automatic detection of communication port

Conformity

This device complies with Part 15 of the FCC rules (US) and RSS-210 (Canada). Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by TéléMédic could void the user's authority to operate the equipment.

6700 Choquette Avenue
St-Hyacinthe (Quebec) J2S 8L1
Tel: (450) 250-1414
Fax: (450) 250-1416
www.telemedic.ca

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Copyright © 2006 TéléMédic Inc. Tous droits réservés.