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Carte d'interface I_{9x}

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I_{9x} Interface card User Manual

Copyright © 2006 Télémédic Inc. All rights reserved.

This equipment is protected by the Canadian copyright laws. It cannot be reproduced, distributed or modified in any ways by others entities than Télémédic, except if a written consent is made by Télémédic.

All the requisite efforts have been done to ensure that the information contained in this document was complete and precise at the moment of printing. However, the information can be modified without notice.

Declaration of conformity

This device is registered as CSA International certified. To see the certificates of official declaration, consult the online documentation on Télémédic's website:

<http://www.telemedic.ca>

This device complies with Industry Canada's low power radio devices norms RSS-210. Certification number: 6418A-I9X001

This device complies with the FCC (Federal Communications Commission) transmission of intentional radio frequencies.

Identification number FCC ID : T3W-I9X-001

Warranty

Consult the warranty booklet included with this product to learn about the limited warranty given by Télémédic on that product.

Trademark

Télémédic, Télémédic's logo and VieGie Med are registered trademark. Other brands and products names mentioned in this document are properties of the respective holders.

Ordering information

You can order updates of this document directly on Télémédic's website: <http://www.telemedic.ca>.

Comments

Télémédic appreciates your comments relating to this manual. Please send your comments directly to info@telemedic.ca.

Support

For all technical questions or defectiveness, please contact Télémédic's user support service at (450) 250-1414 or at info@telemedic.ca

Meaning of symbols used in this manual:

	Warning: Failure to these instructions may lead to unexpected operation or defects of the system.
	Important: Important element to consider in the use of this device.

Meaning of device symbols:

	Warning: Consult accompanying documents. This manual and the associated sections on Télémédic's website (see manual for references) should be consulted before using this device.
	This equipment complies with general requirements for safety of medical electrical equipment CAN/CSA C22.2 No.601.1 and CAN/CSA C22.2 No.60601-1-1, UL 2601-1. This equipment also complies with electromagnetic compatibility requirements for medical electrical equipment CAN/CSA C22.2 No 60601-1-2
	This device complies with FCC (Federal Communications Commission) rules part 15 (Radio frequency devices) for the United States.
	This equipment includes radio frequency transmitter which generate non-ionizing radiation.

The I_{9x} interface card

Some patients need several vital signs detection devices. The I_{9x} interface card allows the incorporation of medical peripherals with a data communication interface such as a blood glucose monitor, blood pressure monitor, scale, etc... to Télémédic's system.

The I_{9x} interface card is therefore used in conjunction with the device you have to use. It is used to capture data and to register data collected by your device to ensure radio frequency (RF) data transmission to the E₉₀ communication terminal.

It is highly recommended to read this manual before using the I_{9x} interface card.

In this section, we describe:

- the main components;
- the device's functions and terms of use;
- the precautionary measures;
- the maintenance;
- the functional preventive check-ups;
- the optimal operating performance.

Main components

I9X interface card indicators

Indicator No.

- 1 Patient identification key connector (type IDX)
- 2 Radio frequency communication indicator LED
- 3 Three-coloured indicator LED displaying device's status

Indicator No.

- 4 RS-232 connector for medical devices
- 5 Connector for Bayer glucometer
- 6 Power input connector


Illustration 2 - I9x Interface card components

The device functions and terms of use

The I_{9x} interface card captures vital signs collected by your medical devices (blood glucose monitor, scales, blood pressure monitor, etc.) that integrate a data communication output. As a patient, you just have to use your medical devices as usual and the I_{9x} interface card will automatically proceed with the data collection and transmission.

The device's status

The three-coloured indicator light (**illustr. 1- no. 3**) displays the I_{9x} interface card status. The following table explains the different states of this indicator.

Interface card status	Indicator light	Description
Low battery	Yellow - flashing every 6 seconds	Battery is low. The I _{9x} card is not functioning and is awaiting connection to recharge.
Battery recharging	Orange - flashing every 6 seconds	Battery is recharging; everything is well connected and the I _{9x} interface card is still functional.
Acquisition	Continuous green followed by green flashing (when data present)	The I _{9x} card is in the process of reading vital signs from a connected medical device.
Undetected material or absence of data	Continuous green followed by absence of flashing or orange flashing	The I _{9x} card was not able to read vital signs from the connected medical device.
Non functional I _{9x} device	Continuous yellow flashing	The I _{9x} card has to be replaced. Contact TéléMédic's technical support.

Table 3 Indicator lights description

Battery recharge

Battery must be recharged for approximately 4 hours before first use.

Connect the power supply plug (**illustr. 1 - no. 6**), provided by TéléMédic (**TéléMédic Part number 50-15001 / Manufacturer CUI Inc. DPD060050-P12P-SZ**). Then connect the power supply to a standard electrical outlet (120V, 60Hz).

The subsequent rechargeable time is approximately 4 hours. Battery recharge should be done once a week or according to use.

Notice that the I_{9x} interface card still functions normally even if connected to the power supply.

Medical devices connection

To use the I_{9x} interface card, first select the appropriate cable to connect your medical device and insert the proper I_{9x} interface card connector. See table 2 below for the appropriate cabling and connectors to use in accordance with the medical device type. Note that the two I_{9x} interface card medical devices connectors can be used simultaneously.

Device Type	Model	Manufacturer	Appropriate TéléMédic cable	Connector to use on the I _{9x} interface card (reference illustration 1)
Scale	UC-321P	Life Source	31-001	Connector number 4
Blood pressure monitor	UA-767PC	AND	31-003	Connector number 4
Blood glucose monitor	6118 (Breeze)	Bayer	31-002	Connector number 5

Table 4 List of supported medical equipment

Patient key connector

The patient key connector (**illustr. 1** - no. 1) is used to identify patients when the I_{9x} interface card and the medical device connected to it are shared by several patients. In that situation, just insert the patient key (IDX type) in the connector. Note that the I_{9x} interface card emits a sound signal during insertion and removal of the key. This indicates that the key is properly functioning. Notice that the insertion and removal of a patient key forces the reading of vital signs stored in the medical devices. Any vital signs values read at these moments will be automatically associated to the patient assigned to the key.

Data acquisition and transmission

The acquisition of vital signs values stored in the medical device connected to the I_{9x} interface card is automatically processed. However, it is important to note that the medical devices should be connected for proper data reading. As soon as the I_{9x} interface card detects new vital signs, it emits a sound signal and forwards them by radio frequency to the E₉₀ communication terminal. The wireless communication indicator light (**illustr. 1** - no. 2) flashes several times during the transmission.

Precautionary measures

The I_{9x} interface card is a precise scientific device that requires some basic precautions during use. Following the advice listed below will ensure your I_{9x} interface card longevity, durability and optimal operating performance.

Never:

- open or dismantle the device;
- hit or throw;
- subject the device to sudden temperature changes;
- put the device in a area exposed to high temperature, high moisture or dusty contents;
- clean the device with a vaporizing liquid, solvent, or detergent;
- immerse the device.

Maintenance

- Use a soft lightly damp cloth.

Functional preventive check-ups

To ensure your I_{9x} interface card optimal operating performance, do the following check-ups:

- Connect the power supply and verify that the three-coloured indicator light (**illustr. 1** - no. 3) flashes orange.
- Connect the I_{9x} interface card to your medical device and collect vital signs as recommended by the device manufacturer. The I_{9x} interface card detects the device and the three-coloured indicator light (**illustr. 1** - no. 3) becomes green and flashes when data is available.
- Introduce a patient key in the patient key connector (**illustr. 1** - no. 1). See paragraph above for patient key info.

Optimal operating performance

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. This medical equipment needs special precautions regarding electromagnetic compatibility (EMC). Installation and services shall take care of following precautions:

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

The use of cables other than those provided or specified by Télémédic may affect the electromagnetic compatibility.

For installation adjacent, stacked or very close to electrical equipment, it is important to verify operation of either product. For I_{9x} interface card refer to functional preventive check-ups above. For the other equipment refer to the manufacturer operating manual.

Also, portable and mobile RF communication equipment ca affects the I_{9x} interface card.

Technical support

If you detect or suspect any malfunction or failure of the I_{9x} interface card, or if you have any questions concerning its installation or use, please contact Télémédic's technical support at **(450) 250-1414**.

Technical specifications

- Data acquisition: Acquisition frequency: 2 minutes to 24 hours
Acquisition time: 30 minutes increment
- Radio communication: Transmission power = 750µW (-1.25 dBm)
Transmission range = 100 metres /328 ft.
Frequency = 916MHz (frequency band ISM)
Télémédic's proprietary protocol
- Alimentation: Internal rechargeable battery 3.7V Lithium Ion
Normal mode autonomy = 504 hours (21 days)
Continuous mode autonomy = 96 hours (4 days)
 Batteries is not intended to be changed by the operator
- Power consumption: Typical 30mA, 500mA while recharging battery
- Memory capacity: 128 Ko, up to 11000 vital signs values
- Data retention: 40 years
- Inputs/Outputs: **Standard serial communication (illustr. 1- no. 4)**
3.5mm /0.138 in. Stereo jack connector.
Serial communication only for Bayer devices
(illustr. 1 - no. 5) 2.5mm /0.098 in. Stereo jack connector
Patient key (illustr. 1 - no. 1)
- Dimensions: 6 cm x 4.5 cm x 2.6 cm /2.36 in. X 1.77 in. X 1.02 in.
- Weight: 45 grams /1.59 oz.
- Temperature: Operation = 10-35 °C (50 °F - 95 °F) /
Storage = 5-40 °C(41 °F - 104 °F)
- Humidity: Operation = 20-80% / Storage = 10-90%
- Security: CAN/CSA C22.2 No.601.1, CAN/CSA C22.2 No.60601-1-1
et UL 2601-1.
- Configurable frequency, period and duration of vital signs acquisition
- Configurable frequency and period of vital signs transmission
- Remote configuration

Conformity

This device complies with Part 15 of the FCC rules (US) and RSS-210 (Canada). Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by TéléMédic could void the user's authority to operate the equipment

6700, Choquette Avenue
St-Hyacinthe (Quebec) J2S 8L1
Tel: (450) 250-1414
Fax: (450) 250-1416
www.telemedic.ca

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Copyright © 2006 Télémédic Inc. Tous droits réservés.