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Unité Portable M₉₃

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M₉₃ Portable Unit User Manual

Copyright © 2006 Télémédic Inc. All rights reserved.

This equipment is protected by the Canadian copyright laws. It cannot be reproduced, distributed or modified in any ways by other entities than Télémédic, except under written consent by Télémédic.

All the requisite efforts have been done to ensure that the information contained in this document was complete and precise at the moment of printing. However, the information can be modified without notice.

Declaration of conformity

This device is registered as CSA International certified. To see the certificates of official declaration, consult the online documentation on Télémédic's website:

<http://www.telemedic.ca>

This device complies with Industry Canada's low power radio devices norms RSS-210. Certification number: 6418A-M93

This device complies with the FCC (Federal Communications Commission) transmission of intentional radio frequencies.

Identification number FCC ID : T3W-M93-XXX

Warranty

Consult the warranty booklet included with this product to learn about the limited warranty given by Télémédic on that product.

Trademark

Télémédic, Télémédic's logo and VieGie Med are registered trademarks. Other brands and products names mentioned in this document are properties of the respective holders.

Ordering information

You can order updates of this document directly on Télémédic's website: <http://www.telemedic.ca>.

Comments

Télémédic appreciates your comments relating to this manual. Please send your comments directly to info@telemedic.ca.

Support

For all technical questions or defectiveness, please contact Télémédic's user support service at (450) 250-1414 or at info@telemedic.ca

Meaning of symbols used in this manual:

	Warning: Failure to these instructions may lead to unexpected operation or defects of the system.
	Important: Important element to consider in the use of this device.

Meaning of device symbols:

	Warning Consult accompanying documents. This manual and the associated sections on Télémédic's website (see manual for references) should be consulted before using this device.
IPX4	This symbol indicates that this device is splash proof.
	This symbol indicates that this device provides a certain level of protection against electrical shock.
	This equipment complies with general requirements for safety of medical electrical equipment CAN/CSA C22.2 No.601.1 and CAN/CSA C22.2 No.60601-1-1, UL 2601-1. This equipment also complies with electromagnetic compatibility requirements for medical electrical equipment CAN/CSA C22.2 No 60601-1-2
	This device complies with FCC (Federal Communications Commission) rules part 15 (Radio frequency devices) for the United States.
	This equipment includes radio frequency transmitter which generate non-ionizing radiation.

M₉₃ portable unit

The M₉₃ portable unit performs periodic vital signs acquisition and transmission.

By wearing this device, you benefit from a distant and precise scientific monitoring of your medical status without compromising your mobility or autonomy. The vital signs acquired by the device are:

- Heart rate;
- Blood oxygen saturation;
- Skin surface body temperature.

It is highly recommended to read the present manual before using the M₉₃ portable unit.

In this section, we describe:

- the main components;
- the device's functions and terms of use;
- the precautionary measures;
- the maintenance;
- the functional preventive check-ups;
- the optimal operating performance.

Main components

M93 Portable unit indicators

Indicator No.

- 1 Mainframe
- 2 Date, time and messages display screen (LCD)
- 3 Vital signs reading indicator – ♥
- 4 Indicator that shows vital sign acquisition is affected by too much patient's movement – ㄥㄥ
- 5 Radio frequency communication indicator – ㄩㄩㄩ
- 6 Continuous vital signs reading mode indicator – ♥
- 7 Equipement type identification – *M93*

Indicator No.

- 8 Battery recharge connector
- 9 Vital signs sensor
- 10 Display activation button (time, battery power, message selection)
- 11 Integrated battery
- 12 Wrist band
- 13 Continuous acquisition mode start and stop button
- 14 Radio frequency vital signs transmission integrated antenna

Illustration 12 M93 portable unit components

M₉₃ functions and terms of use

Installation

M₉₃ portable unit must be worn around the wrist, such as a regular watch, with your thumb inside the ring where sensors are located. (**illustr. 1** - no. 9).

The device's mainframe should be facing up on the right or left arm, depending on the model. The position and close fit can be adjusted by using the strap.

Illustration 13 M93 portable unit installation

Turning on power

M₉₃ portable unit is permanently on and set on reduced power consumption mode. To activate the device, push the display activation switch (**illustr.1** - no. 10) or the continuous mode on /off switch (**illustr. 1** – no. 13).

Switch activation

The display activation switch (**illustr. 1** – no.10) enables you to display the following messages (on a rotational basis, by pressing on the switch):

Date and time:

- First line: MM/DD
- Second line: HH : MM

Battery power:

- First line: Battery
- Second line: Percentage indicating the energy level available (65% in this example).

Illustration 14 Time and battery power display

When going to daylight saving time or back to standard time, it is possible that a slight time lag can be observed between the actual time and the time displayed on the M₉₃ portable unit. The actual time will automatically be adjusted during the next complete communication with the communication base and server.

Vital signs display (According to configuration)

French		English	
			
Fréquence Cardiaque = 80 bpm prise à 12 :00		Heart Rate = 80 ppm taken at 12 :00	
French		English	
			
Saturation en Oxygène = 99% prise à 12 :00		Oxygen Saturation = 99% taken at 12 :00	
French		English	
			
Température = 30 °C prise à 12 :00		Temperature = 80 °F taken at 12 :00	

Display closing

In order to maximise battery use, the M₉₃ portable unit crystal liquid display automatically shuts off after a 10 second idle period. Just press the switch located on the right side (**illustr. 1** – no. 10) to activate display.

Data acquisition

Vital signs acquisition is done automatically by the M₉₃ portable unit without you having to intervene. When the M₉₃ portable unit is set on acquisition mode, the acquisition indicator is displayed on the screen (**illustr. 1** - no. 3).

At that moment, if your hand is in movement, the movement indicator will be displayed on the screen (**illustr. 1** – no. 4). This happens when your hand is moving too much for the device to be able to read your vital signs. If this situation lasts more than 30 seconds and your M₉₃ unit is not set on continuous mode, you will hear a sound signal indicating to stop your activities (hand movement) for the duration of the acquisition time; The default setting being two minutes. Note that the sound signal is turned off by default between 8pm and 8am to avoid being disturbed during the night. These parameters and default settings can be modified by your health professional.

Continuous mode vital signs acquisition

The continuous mode acquisition is a means of manually reading vital signs. This mode is activated by pressing the continuous mode activation switch, located on the left side (**illustr. 1** – no.13). The portable unit will then immediately read vital signs and the continuous reading indicator will be displayed on the M₉₃ portable unit screen (**illustr. 1** – no. 6).

The continuous mode acquisition will stop when pressing the continuous reading switch a second time. It is possible to notice a slight delay before the continuous mode actually cuts itself off, this only means that the last data acquisition has to be completed. The data collected will then be automatically transferred to the E₉₀ communication terminal. Moreover, if you forget to turn the continuous mode off, it will automatically shut itself off after a 20 minute data acquisition period (by default).

Transmission

Programmed in advance, the collected data transmission is automatically carried out by the M₉₃ portable unit at 6:30am, 11am and 8:30pm every day, by default. No intervention is required.

When the M₉₃ portable unit is set on radio frequency mode, the vital signs transmission indicator is displayed on the screen (**illustr. 1** – no. 5). The transmission is carried out to the E₉₀ communication terminal.

If transmission fails after a few seconds, an "X" will appear on the screen and disappear after a few seconds. Do not worry, the M₉₃ portable unit has huge storing capacities and can retain several days worth of acquired data. This means no data loss while your are away from the E₉₀ communication terminal during the scheduled transmission period. The transmission will proceed during the next scheduled period, when you will be in proximity to the E₉₀ communication terminal (less than 100 meters/330 feet).

Battery recharge

Battery must be recharged for approximately 4 hours before first use.

Connect the power supply plug provided by Télémédic (**Télémédic Part no 50-15001 / Manufacturer CUI Inc. DPD060050-P12P-SZ**) to the battery recharge connector (**illustr. 1** – no. 8), then connect the power supply to a standard electrical outlet (120V, 60Hz). Display will show "Charging" when properly connected.

French	English
	
Illustration 18 Battery recharge message	

When the portable unit is in "charging" mode, the acquisition and transmission functionality are deactivated. The portable unit awaits deconnection to resume its normal operations.

French	English
	

Illustration 19 Charging completed message

When battery is low (15% of its capacity), the "battery low" message will be displayed on the M₉₃ portable unit.

French	English
	

Illustration 20 Battery low message

When battery is too low for the portable unit to operate, the "battery empty" message appears and will not disappear until battery is recharged. If this situation lasts for a few days, the portable unit will not have enough power to display the message and it will disappear.

French	English
	

Illustration 21 Battery empty message

 Subsequent charging time is of approximately 4 hours during which the patient must take off the portable unit. Battery recharge should be done once a week or according to use.

M₉₃ portable unit will initialize when charging empty battery. During this 5 second period, the following message will be displayed on the screen.

Once initialization is completed, the software version will be displayed for a few seconds.
Ex.: 106-615-100-R010

Illustration 22 Initializing and software messages

Precautionary measures

M₉₃ portable unit is a precise scientific device that requires some basic precautions. Following the advices listed below will ensure your portable unit's longevity, durability and optimal operating performance.

Never:

- Open or dismantle the device;
- hit or throw;
- subject the device to sudden temperature changes;
- put the device in an area exposed to high temperature, high moisture or dusty contents;
- clean the device with a vaporizing liquid, a solvent, or a detergent;
- immerse the device.

The wrist band (**illustr. 1** – no. 12) as well as the sensors (**illustr 1** – no. 9) can be detached. If you detect or suspect malfunction or part failure or defect, **do not try** dismantling, please contact technical support.

Maintenance

- Use a soft lightly damp cloth to clean the portable unit parts such as the wrist band, mainframe and ring.

Functional preventive check-ups

To ensure your M₉₃ portable unit's optimal operating performance, do the following check-ups:

- Press the message selection button (**illustr. 1** – no. 10). A message should reappear on the display screen.
- If no messages appear, battery could be completely empty: initiate the battery recharge process.

Optimal operating performance

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. This medical equipment needs special precautions regarding electromagnetic compatibility (EMC). Installation and services shall take care of following precautions:

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

For installation adjacent, stacked or very close to electrical equipment, it is important to verify operation of either product. For M₉₃ portable unit refer to functional preventive check-ups above. For the other equipment refer to the manufacturer operating manual.

Also, portable and mobile RF communication equipment can affects the M₉₃ portable unit.

Technical support

If you detect or suspect any malfunction or failure of the M₉₃ portable unit, or if you have any questions concerning its installation or use, please contact Télémédic's technical support at **(450) 250-1414**

Technical specifications

- Vital signs acquisition:
 - SPO₂ (Oxygen saturation)
(0% to 100% ± 3 digits)
 - Heart rate (0 to 511 PPM)
 - Body temperature (28 °C to 37 °C ± 0,10 °C)
(82.4 °F to 98.6 °F±0,10°F)
- Programmable acquisition:
 - Acquisition frequency: 2 minutes to 24 hours
 - Acquisition time: 30 minutes increment
 - Acquisition duration: 1 minute to 24 hours
 - Sampling frequency: 1 second
 - Continuous mode average: 1 to 86400 samples
 - Standard mode average: 1 to 86400 samples
- Radio communication:
 - Transmission power = 750µW (-1.25 dBm)
 - Transmission range = 100 meters/328 ft.
 - Frequency = 916MHz (frequency band ISM)
 - TéléMédic's proprietary protocol
- Electrical power:
 - Rechargeable internal battery: 3.7V Lithium Ion
 - Standard mode autonomy = 504 hours (21 days)
 - Continuous mode autonomy = 42 hours
 -  Battery is not intended to be changed by the operator.
- Power consumption:
 - Typical 20mA, 500mA while charging battery.
- Storing capacity:
 - 128 Ko, up to 11000 vital signs values
- Data retention:
 - 40 years
- Dimensions:
 - 5.2 cm x 6.8 cm x 2.7 cm/2.05 in. x 2.67 in. x 1.06 in.
- Weight:
 - 136 grams/4.79 oz.
- Temperature:
 - Operation = 10-35 °C (50-95 °F)
 - Storage = 5-40 °C (41-104 °F)
- Humidity:
 - Operation = 20-80% / Storage = 10-90%
- Security:
 - IPX4, Type BF patient applied part
 - CAN/CSA C22.2 No.601.1, CAN/CSA C22.2 No.60601-1-1 et UL 2601-1.
- Configurable frequency, period and duration of vital signs acquisition
- Configurable frequency and period of vital signs transmission
- Remote configuration

Conformity

This device complies with Part 15 of the FCC rules (US) and RSS-210 (Canada). Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by TéléMédic could void the user's authority to operate the equipment.

6700 Choquette Avenue
St-Hyacinthe (Quebec) J2S 8L1
Tel: (450) 250-1414
Fax: (450) 250-1416
www.telemedic.ca

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Copyright © 2006 TéléMédic Inc. Tous droits réservés.