TRANSTEK

User Manual

Blood Pressure Monitor TMB-1591-BS



FCC ID:OU9TMB1591B2



- Thank you very much for selecting TRANSTEK Blood Pressure Monitor TMB-1591-BS.
- To use the monitor correctly and safely, please read the manual thoroughly.
- Please keep this manual well in order to reference in future.

Table of Contents

INTRODUCTION
BEFORE YOU START
MEASUREMENT
INFORMATION FOR USER
ABOUT BLOOD PRESSURE
TROUBLESHOOTING. 18 SPECIFICATIONS. 19 AUTHORIZED COMPONENT 20 CONTACT INFORMATION. 20 FCC STATEMENT. 21 EMC GUIDANCE. 22

INTRODUCTION INTRODUCTION

General Description

Thank you for selecting TRANSTEK arm type blood pressure monitor (TMB-1591-BS). The monitor features blood pressure measurement, pulse rate measurement and the result storage. The design provides you with two years of reliable service.

Readings taken by the TMB-1591-BS are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

This manual contains important safety and care information, and provides step by step instructions for using the product. Read the manual thoroughly before using the product.

Features:

- · 68mm*90mm Digital LCD display with White backlight
- · Maximum 99 records
- Measuring during inflation technology

♥ Indications for Use

The Transtek Blood Pressure Monitor is digital monitors intended for use in measuring blood pressure and heartbeat rate with arm circumference ranging from 22cm to 42cm(about 83/4"-161/2"). It is intended for adult indoor use only.

♥ Contraindications

- 1. The device is not suitable for use on may be pregnant women or pregnant women.
- 2. The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.

Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

♥ Safety Information

The signs below might be in the user manual, labeling or other component. They are the requirement of standard and using.

$\overline{}$,
③	Refer to instruction manual/booklet To signify that the instruction manual/booklet must be read.
†	Symbol for "Type BF applied part"
À	Caution Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
X	The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
<u>6</u> 8	Symbol for "Recycle"
SN	Symbol for "Serial Number"
==	Symbol for "Direct Current"
	Symbol for "Manufacturer"
쎄	Symbol for "Date and Country of manufacture"
	For indoor use only
	Symbol for "Class II Equipment"
_	

INTRODUCTION

- ACAUTION

- * This device is intended for indoor, home use.
- * This device is not intended for public use.
- * This device is portable, but it is not intended for use during patient transport.
- * This device is not suitable for continuous monitoring during medical emergencies or operations.
- *This device is intended for no-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm, or for any purpose other than obtaining a blood pressure measurement
- * This device is for adults. Do not use this device on neonates or infants. Do not use it on children unless otherwise instructed by a medical professional.
- * Do not use on the women in pregnant, including pre-eclamptic, patients.
- * The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.
- * The effectiveness of this device has not been established for use:
- -on users with common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation
- -on users with peripheral arterial disease,
- -on users undergoing intravascular therapy, or with arteriovenous (AV) shunt.
- Consult a medical professional before use.
- * Do not use this device for diagnosis or treatment of any health problem or disease. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or health care professional.
- * If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.
- * This device may be used only for the intended use described in this manual, the manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse
- * Report any unexpected operation or events to the manufacturer.
- * Do not apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.
- * Warning: Do not kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result injury.
- * Warning: Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- * Warning: Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.
- * Warning: Do not place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt, which could disrupt blood circulation and cause injuries.
- * Do not place the cuff on the arm on the same side of a mastectomy (especially when lymph
- nodes have been removed). it is recommended to take measurements on the unaffected side.
- * Do not wrap the cuff on the same arm to which another monitoring device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.
- * Please check that the operation of the device do not result in prolonged impairment of patient blood circulation.
- *Warning: On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.
- * Warning: Do not use this device with high-frequency (HF) surgical equipment at the same time.


CAUTION

- * Warning: This device is not AP/APG equipment. Do not use the device where flammable anesthetic are present, or in environments mixture with air of with oxygen or nitrous oxide.
- * The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast translentburst signal.
- *Wireless communication equipment, such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies may cause interference that may affect the accuracy of measurements. A minimum distance of 1 foot (30 cm) should be kept from such devices during a measurement.
- * You can use this device to take your own measurement, no third-party operator is required.
- * Please use the device under the environment which is provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
- * The device may require up to 30 minutes to warm up / cool down from the minimum / maximum storage temperature before it is ready for use.
- * Warning: Excessive cuff tube lengths could cause strangulation if you don't manage them properly.
- * Warning: Do not touch output of the batteries/adapter and the user simultaneously.
- Adapter is specified as a part of ME EQUIPMENT.
- * Warning: The power cord is considered the disconnect device for isolating this equipment from supply mains. Do not position the equipment so that it is difficult to reach or disconnect.
- * The blood pressure monitor, its adapter, and the cuff are suitable for use within the patient environment.
- * Warning: Do not use this device if you are allergic to polyester, nylon, or plastic.
- * Warning: Only use accessories approved by manufacturer. Using unapproved accessories might cause damage to the unit and injure users.
- * Warning: If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Power button immediately to release the air from the cuff.
- No calibration is required within two years of reliable service.
- * Do not attempt to repair the unit yourself if it malfunctions. Only have repairs carried out by authorized service centers.
- * At the request of authorized service personnel, circuit diagrams, component part lists, descriptions, and calibration procedures will be made available by the manufacturer or distributor.
- * It is recommended that the performance should be checked after repair, maintenance, and every two years of use, by retesting the requirements in limits of the error of the cuff pressure
 - indication and air leakage (testing at least at 50 mmHg and 200 mmHg).

 Warning: Do not use the device while under maintenance, or being serviced.
- * Store your device, cuff and adapter in a clean and dry place, protect it against extreme moisture,
- heat, lint, dust and direct sunlight. Never place any heavy objects on it.
- * Make sure the rubber tube of the cuff is not squeezed, stretched, or kinked during storage.
- * Warning: Keep the device, cuff, and batteries away from children as they may pose a risk of choking or strangulation if used improperly.
- * Clean both device and cuff with a soft, dry cloth. If necessary use a dampened cloth and natural detergent. Do not use alcohol, benzene, or other harsh chemicals.
- * Do not wash the cuff in a washing machine or dishwasher!
- * The service life of the cuff may vary by the frequency of washing, skin condition, and storage state. The typical service life is 10000 times.
- Dispose of accessories, detachable parts, and the device according to the local guidelines.

INTRODUCTION

▼ LCD Display Signal




SYMBOL	DESCRIPTION	EXPLANATION		
SYS	Systolic blood pressure	High pressure result		
DIA	Diastolic blood pressure	Low pressure result		
PULSE	Pulse	Pulse/minute		
mmHg	mmHg	Measurement Unit of the blood pressure		
((<u>\(\)</u>))	Irregular heartbeat	Irregular heartbeat detection		
	Battery Indicator	Indicate the current battery		
	Grade	The grade of the blood pressure		
A E	Shocking reminder	Shocking will result in inaccurate		
*!	Data transmission error	Data transmission error		
\bigcirc	Heartbeat	Heartbeat dectetion during measurement		

♥ Monitor Components



Component list of pressure measuring system

- 1 Cuff
- 2 Air pipe 3 PCBA
- 4 Pump
- 5 Valve



♥ List

1. Blood Pressure Monitor 2. Cuff (about 22cm~42cm) (TMB-1591-BS)







3. 4*AA alkaline batteries

4. User manual




5. AC adaptor (UE08WCP-060100SPA)

▼ The Choice of Power Supply

- 1.Battery powered mode: 6VDC 4*AA alkaline batteries
- 2.AC adaptor powered mode: 6V === 1A

(Can be supplied by AC adaptor model UE08WCP-060100SPA!)

Right picture is the hole in for power adaptor.



ACAUTION

In order to get the best effect and protect your monitor, please use the right battery and special power adaptor. The power adapter is a part of the device After using, please pull out the adaptor plug insulates from the main supply. Do not position the device in a position where it is difficult to disconnect from the supply mains.

Installing and Replacing the Batteries

- 1. Slide off the battery cover.
- 2. Install the batteries by matching the correct polarity, as shown.
- 3. Replace the cover.



Replace the batteries whenever the below happen

- The to + m shows
- •The display dims
- The display does not light up


ACAUTION .

- Remove batteries if the device is not likely to be used for some time.
- The old batteries are harmful to the environment, do not dispose with other daily trash.
- Remove the old batteries from the device and follow your local recycling quidelines.

♥ Pairing

Turn on Bluetooth and APP. Make sure both are ON when pairing is proceeding.

When the blood pressure monitor is off, press On/Standby button for 2 full seconds to enter bluetooth pairing mode. The symbol flashes, indicating the pairing is proceeding.



Then please select the user ID you want to connect with your smartphone on the app to continute the pair-up.

If succeed, symbol will turn off.



If fail, symbol turns off.

will flash all the time until the blood pressure

Bluetooth Module No.: LS51802

RF Frequency Range: 2402 MHz to 2480 MHz


Output Power Range: -3.21 dbm Supply Voltage: 2V-3.6 V

Transmitting Distance: 10 meters

▼ Tie the Cuff

- Tie the cuff on your upper arm, then position the tube off-center toward the inner side of arm in line with the little finger.
- The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.
- **3**.Sit comfortably with your arm resting on a flat surface.
- 4. Correct position:
 - Bare your arm or wear tights only when starting measurement.
 - Sit comfortably with legs uncrossed, feet flat on the floor, back and arm supported. The center of the cuff should be at the same level as the right atrium of the heart.
- Rest for 5 minutes before measuring.
- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same upper arm, or as directed by a physician.








♥ Start Measurement

When the blood pressure monitor is off, press On/Standby button to turn it on, it will finish the whole measurement.

LCD display








Inflating and Measuring

Display and save the result.





The blood pressure monitor will proceed to data transmission after measurement. The bluetooth symbol flashes on the LCD indicates data is transmitting.



After successful transfer, the device powers off the Bluetooth radio and the icon (and rectangle) are removed.

If the user presses and releases the On/Standby button, another reading is initiated.

If the user presses and holds the On/Standby button for 2 seconds, the device powers down.

Or if there is no operation, after a 10 seconds of inactivity, the device powers down.


Notes: During inflation, the heart icon in the upper left blinks in accordance with the user's pulse rate.

Additionally, during inflation the progress metre to the right of the digits builds vertically up as the pressure increases according to the table below.

Segments	Pressure
1	>= 0
2	>= 40
3	>= 80
4	>= 120

Segments	Pressure
5	>= 140
6	>= 160
7	>= 180
8	>= 200

During the measurement, if you press the On/Standby button to stop the measurement, the numerics are cleared. It will display as below:



If you press and release the On/Standby button, another reading is initiated

If you press and hold the On/Standby button for 2 seconds, the device powers down.

Or if there is no operation, after a 10 seconds of inactivity, the device powers down.

-/i\CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the oldest record (99) is dropped from the list.

▼ Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances.



wait at least 1 hour after dinner or drinking





Wait at least 20 minutes after taking a bath



In a very cold environment



When talking or moving your fingers



When you want to discharge urine

▼ Maintenance

In order to get the best performance, please follow the instructions below.



Put in a dry place and avoid the sunshine



Avoid intense shaking and collisions



Using wet cloths to remove dirt



Avoid touching water, clean it with a dry cloth in case.



temperature environment



Avoid washing the cuff

What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.





■ What is the standard blood pressure classification?

The chart on the right is the standard blood pressure classification published by American Heart Association (AHA).

Blood Pressure Category	Systolic mmHg (upper#)		Diastolic mmHg (lower#)	
Normal	less than 120	and	less than 80	
Elevated	120-129	and	less than 80	
High Blood Pressure (Hypertension) Stage 1	130-139	or	80-89	
High Blood Pressure (Hypertension) Stage 2	140 or higher	or	90 or higher	
Hypertensive Crisis insult your doctor immediately	Higher than 180	and/or	Higher than 120	



Please consult a physician if your measuring result falls outside the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

♥ Irregular Heartbeat Detector


An irregular heartheat is detected when a heartheat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, the monitor records all the pulse intervals and calculate the average; if there are two or more pulse intervals, the difference between each interval and the average is more than the average value of ±25%, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of ±15%, the irregular heartbeat symbol appears on the display when the measurement results are appear.



The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why does my blood pressure fluctuate throughout the day?

 Individual blood pressure varies multiple times everyday. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.



3. Wait at least 3 minutes for another measurement.

Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc, Also, there is the "white coat" effect, which means blood pressure usually increases in clinical settings.

Is the result the same if measuring on the right arm?

It is ok for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.



What you need to pay attention to when you measure your blood pressure at home:

If the cuff is tied properly.

If the cuff is too tight or too loose.

If the cuff is tied on the upper

If you feel anxious.

Taking 2-3 deep breaths before beginning will be better for measuring.

Advice: Relax yourself for 4-5 minutes until you calm down.



16 arm every time.

TROUBLESHOOTING SPECIFICATIONS

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the products not operating as you think it should, check here before arranging for servicing.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY		
	Display	Batteries are exhausted.	Replace with new batteries		
No power	will not light up.	Batteries are inserted incorrectly.	Insert the batteries correctly		
		AC adaptor is inserted incorrectly.	Insert the AC adaptor tightly		
Low batteries	The display indicates the "BAT LO" message, pauses for 3 seconds. The battery icon shows empty (does not flash.)	Batteries are low.	Replace with new batteries		
Error message	*! shows	Unsuccessful pairing.	Check if both the APP and Bluetooth are on, operate and send the data again.		
	E 1 shows	The cuff is not secure.	Readjust the cuff and relax for a moment and then measure again.		
Warning message	"out " shows	Out of measurement range	Relax for a moment. Refasten the cuff and ther measure again. If the problem persists, contact your physician.		

	Battery powered mode:			
	6VDC 4*AA alkaline batteries			
Power supply	AC adaptor powered mode:6V===1A			
	(Can be supplied by AC adaptor model			
	UE08WCP-060100SPA!)			
Display mode	Digital LCD V.A.68mm*90mm			
Measurement mode	Oscillographic testing mode			
Measurement range	Rated cuff pressure:			
weasurement range	0mmHg~299mmHg(0kPa ~ 39.9kPa)			
	Measurement pressure:			
	SYS: 60mmHg~230mmHg (8.0kPa~30.7kPa)			
	DIA: 40mmHg~130mmHg (5.3kPa~17.3kPa)			
	Pulse value: (40-199)beat/minute			
Accuracy	Pressure:			
	5℃-40℃ within±3mmHg			
	pulse value:±5%			
Normal working condition	A temperature range of :+5°C to +40°C			
rtormar working condition	A relative humidity range of 15% to 90%,			
	non-condensing, but not requiring a water			
	vapour partial pressure greater than 50 hPa			
	An atmospheric pressure range of :			
	700 hPa to 1060 hPa			
Storage & transportation	Temperature:-20°C to +60°C			
condition	A relative humidity range of ≤ 93%, non-condensing			
condition	at a water vapour pressure up to 50hPa			
Measurement perimeter				
of the upper arm	About 22cm~42cm			
Net Weight	Approx.300g(Excluding the dry cells)			
External dimensions	Approx.92mm*140mm*46mm			
Attachment	4*AA alkaline batteries,user manual,AC adapter			
Mode of operation	Continuous operation			
Degree of protection	Type BF applied part			
Protection against ingress of water	IP22			
Software Version	V01			

FCC STATEMENT

▼ Authorized Component

please use the TRANSTEK authorized adapter.

Adapter

Type: UE08WCP-060100SPA Input: 100~240V, 50~60Hz,400mA

Output: 6V == 1A

(Conforms to UL certificate)

♥ Contact Information

For more information about our products, please visit www.transtekcorp.com.you can get customer service, usual problems and customer download, transtek will serve you anytime.

Manufactured by: Guangdong Transtek Medical Electronics Co., Ltd.
Company: Guangdong Transtek Medical Electronics Co., Ltd.
Address: Zone A, No.105, Dongli Road, Torch Development District,
528437 Zhongshan, Guangdong, China

▼ FCC Statement

FCC ID:OU9TMB1591B2

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution. The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- -- Reorient or relocate the receiving antenna.
- -- Increase the separation between the equipment and receiver.
- -- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- -- Consult the dealer or an experienced radio/TV technician for help. FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

EMC GUIDANCE EMC GUIDANCE

▼ EMC Guidance

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment TMB-1591-BS including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

- All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
- 2. Guidance and manufacturer's declaration-electromagnetic emissions and Immunity.

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions				
Emissions test	Compliance			
RF emissions CISPR 11	Group 1			
RF emissions CISPR 11	Class [B]			
Harmonic emissions IEC 61000-3-2	Class A			
Voltage fluctuations / flicker emissions IEC 61000-3-3	Comply			

Table 2

Guidance and manufacturer's declaration – electromagnetic Immunity					
IEC 60601-1-2 Test level	Compliance level				
±8 kV contact	±8 kV contact				
±2 kV, ±4 kV, ±8 kV, ±15 kV air	±2 kV, ±4 kV, ±8 kV, ±15 kV air				
±2 kV for power supply lines	±2 kV for power supply lines				
±1 kV signal input/output	Not Applicable				
100 kHz repetition frequency	100 kHz repetition frequency				
±0.5 kV, ±1 kV differential mode	±0.5 kV, ±1 kV differential mode				
±0.5 kV, ±1 kV, ±2 kV common mode	Not Applicable				
0% UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle	0% Uτ; 0,5 cycle. At 0°, 45°, 90°, 135° 180°, 225°, 270° and 315°. 0% Uτ; 1 cycle and 70% Uτ; 2530 cycles; Single phase: at 0°. 0% Uτ; 250 / 300 cycle				
30 A/m	30 A/m				
50 Hz / 60 Hz	50 Hz / 60 Hz				
3 V	3 V				
0,15 MHz – 80 MHz	0,15 MHz – 80 MHz				
6 V in ISM and amateur radio bands	6 V in ISM and amateur radio bands				
between 0,15 MHz and 80 MHz	between 0,15 MHz and 80 MHz				
80% AM at 1 kHz	80% AM at 1 kHz				
10 V/m	10 V/m				
80 MHz – 2,7 GHz	80 MHz – 2,7 GHz				
80% AM at 1 kHz	80% AM at 1 kHz				
	IEC 60601-1-2 Test level ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air ±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency ±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode 0% UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0% UT; 1 cycle and 70% UT; 25/30 cycles; Single phase: at 0°. 0% UT; 250 / 300 cycle 30 A/m 50 Hz / 60 Hz 3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz 10 V/m 80 MHz – 2,7 GHz				

EMC GUIDANCE EMC GUIDANCE

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity								
Radiated RF IEC61000-4-3 (Test specifications	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
for ENCLOSURE PORT	385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27	27
IMMUNITY to RF wireless communicati-	450	430-470	GMRS 460, FRS 460	FM ± 5k Hz deviation 1 kHz sine	2	0.3	28	28
ons equipment)	710	704-787	LTE Band	Pulse	0.2	0.3	9	9
,	745		13, 17	modulation 217 Hz				
	780			Pulse	-			28
	810	800-960	800/900,	modulation 18 Hz	2	0.3	28	26
	870							
	930							
	1720	1700- 1990	CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS	217 Hz	2	0.3	28	28
	1845							
	1970							
	2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28	28
	5240	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9
	5500							
	5785							
	3100					<u> </u>		