

Keypad Lights

The CapTel keypad lights up for easy dialing in dark environments. You can set the keypad lights to go off when your CapTel is not in use. When set to “off” the keypad lights will go off after approximately 30 seconds of non-use. The default setting is Keypad Lights On.

 **Settings** →  **General Settings** →  **Keypad Lights On**

Language

You can set CapTel menus and captions to be in Spanish or English. The default setting is English.

 **Settings** →  **General Settings** →  **English** or  **Spanish**

Setting the Ringer Volume

You can turn the Ringer Volume up to a setting you hear best. The default setting is five (5).

 **Settings** →  **Sounds** →  **Ringer Volume**

Adjusting the Ringer Pitch

You can adjust the frequency of the CapTel ringer to a High, Medium, or Low range that you hear best. The default setting is Med.

 **Settings** →  **Sounds** →  **Ringer Pitch**

Tone Control

You can adjust the frequency of the sound coming through the CapTel handset to a frequency range that you hear best. Select from Low, Medium, High, or Customizable Frequency Ranges. To customize the frequency range, see page 29 for instructions.

 **Settings** →  **Sounds** →  **Tone**

Time & Date

You can adjust the time and date displayed on your CapTel phone screen. Touch the white area around the date/number that you want to change.

 **Settings** →  **Time and Date**

About Your CapTel Phone

Provides details about how your CapTel phone is set up to help with troubleshooting issues.

 **Settings** →  **About**

C A P T E L 2 4 0 0 i I N D E X

Answering a call.....	26	Captions.....	4, 27, 31-41, 46,71
with captions	26	button.....	6,9
on an extension phone	26	change font size.....	38
Answering Machine.....	49	change font color.....	39
erase messages	52, 53	erase.....	37
play messages	51	how captions scroll.....	41
play messages aloud/silent	58	In Spanish.....	13, 72
record personal greeting.....	54	review during a call.....	34
see captions.....	51, 52, 58	review after a call.....	35
set number of rings	57	review captions of previous call.....	46
turn on/off.....	50	save captions after hanging up.....	36
turn on/off.....	50	turn on/off.....	32
Audio Jacks	11	Contacts.....	59
Buttons.....	8,9	add new	60
captions	6,9	add picture to contact....	60, 62, 64, 65
CUST SERV.....	6,9	remove picture from contact.....	??
dialing pad.....	6,8	editing.....	62
flash.....	6,8	erase a contact	63
home.....	6,9	dial from.....	61
mute.....	6,9	Conversation Memory.....	IFC, 34
redial.....	6,8	Corrections	33
tone	6,8	change color	39
volume	6,8	Customer Service button	6, 9
Calls	23	calling for help.....	22, 43, BC
make/answer calls.....	24, 26	Spanish customer service.....	??
use Call History.....	47	Dialing	23
Call History.....	44	by picture	65
dial from.....	47	dialing 9 in an office setting	47
review calls	45	from Call History	47
review captions from call	46	from Contact list	61
erase.....	48	from Favorites.....	68
Call waiting	8	from number pad.....	24
		with Speed Dial.....	65
		Display Screen.....	6, 41
		adjust brightness.....	42, 70
		screen timeout	42, 70
		night mode.....	70
		Flash button	8

Flasher	71	Reset (problem solving)	xx
set the ring flasher	71	Ringer	72
Headset	11	ring flasher	71
use with CapTel	11,28	adjust volume	72
Indicator Lights	8, 9, 27	turn ringer off	72
on keypad	72	adjust pitch	72
mute	9, 27	Ring flasher	71
volume	8, 27	Scroll, captions across display	41
ring flasher	71	smooth scrolling	41
Internet	12, 15-18	word-by-word scrolling	41
connecting to the Internet	15	Settings	69-73
requirements	12	Setup	12
using a router	15,16	help with set up	22
using WiFi	17	requirements	12
Keypad	6	using a router	15,16
lights	72	using wireless/WiFi	17-18
dialing from	24	Spanish	13, 72
Language	13, 72	captions	13, 72
Spanish	13, 72	customer service in Spanish	??
English	13, 72	Speed dial	65
Messages, Answering Machine	51	adding contact to	65
play messages	51	adding picture to	64
erase messages	52	dialing from	65
play aloud/silent	58	erasing contact from	63
Mute button	6,9	Time & Date	19, 73
Neckloop	11	Tone Control	28, 73
use with CapTel	11,28	adjust during a call	28
Network	15	customize Tone setting	29
WiFi	17	Volume	8
change settings	17	adjust during a call	27
available networks	17	button	6, 8
view saved networks	18	save conversation volume	71
Reading captions	31	Ringer volume	72
after hang up	35	WiFi/Wireless network	17, 18
during a call	34	view saved networks	18
from previous call	46	enter network manually	17
Redial button	6,8	remove saved network	??

FCC Information

Compliance Statement (Part 15.19)

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

Warning (Part 15.21)

Changes or modification not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC Interference Statement (Part 15.105b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not

occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF Exposure (OET Bulletin 65)

To comply with FCC/IC exposure requirements for mobile transmitting devices, this transmitter should only be used or installed at location where there is at least 20 cm separation distance between the antenna and all persons.

To comply with FCC/IC RF exposure limits for general population/uncontrolled exposure, the antenna(s) used for this transmitter must be installed to provide a separation distance of at least 20 cm from all persons and must not be co-located or operating in conjunction with any other antenna or transmitter.

CapTel is the latest innovation from  **Ultratec**.

Ultratec, Inc.

450 Science Drive

Madison, WI 53711

(888) 269-7477 V/TTY

(866) 670-9134 (Spanish CapTel Customer Service)

Email: CapTel@CapTel.com • Fax: (608) 204-6167 • Web: www.CapTel.com

©2013 Ultratec, Inc. All rights reserved.

Ultratec and CapTel are registered trademarks of Ultratec, Inc.

Getting Help

If you have any questions as you use your new CapTel 2400i, we're here to help. Use any of these ways to contact us:

Pick up the handset and press  (**CUST SERV**) to automatically speed dial directly to our helpful Customer Service team. If you are in an office setting, you may need to dial a number such as "9" first to get an outside line first and then press  (**CUST SERV**). *(Available 24 hours a day, 7 days a week. Customer Service is closed on major holidays.)*

Help materials are also available online at our website.



Phone: 1-888-269-7477 *(Available 24 hours a day, 7 days a week)*

FAX: (608) 204-6167



Mail: CapTel Customer Service
450 Science Drive
Madison, WI 53711



Email: CapTel@CapTel.com

Online: www.CapTel.com
(Live Chat help available at our website)