

8.4 Subject Details

Shows all available subject details.

ERT - MasterScope CT

Current Subject
 SUBJECT INITIALS AB
 SUBJECT NUMBER 99333002
 SEX Female
 DATE OF BIRTH 08FEB1958

PROTOCOL CENTER # 99333

Visit V101 Visit Summary Visit Calendar Subject Details

Search/New Subject

V101

Started

Select Next Action

Pre ECG

FeNO

Pre PFT

Ipratropium Bromide

Reschedule Cancel

Subject Details

Subject Initials AB Patient initials must be 1-3 characters in uppercase

Subject Number 99333 002 Subject Number must be between 001 to 999

Date of Birth 08 Feb 1958 Age must be between 40 and 120

Age 56 Age will be calculated automatically

Sex Female Male Female or Male

Height 165 cm inches Height must be between 90cm and 230cm

Ethnicity

Asian Japanese

Black Mixed Ethnicity

Caucasian Native American

Chinese Pacific Islander

Hispanic or Latino Other

Indian (Indian subcontinent)

Save Undo

Center Number: 99333 | Systemid: 028904 | Study Version: 800286-800286-1.4.0-4.0.21 (20JUN2014 11:17:42) | Battery Power: 100% Remote: Disabled SYS#: 800286-028904-000004 | 23JUL2014 14:26:15

Update the values as described in Chapter “**System Overview**” - “**New Subject**”.



Click on <Save> to store the data as entered.

9. Spirometry Measurement

Spirometry consists of the Slow Spirometry and /or Forced Spirometry measurement.

9.1 Basic Conditions prior to starting a Measurement



In order to perform high-quality measurements, the following conditions should be fulfilled:

- A waiting time of 10 to 15 minutes prior to the first examination ensures normal ventilation as a prerequisite for a high reproducibility of the examination results.
- To improve motivation, cooperation and coordination of the subject, the measurement procedure should be explained or, if necessary, demonstrated prior to starting the measurement.
- The subject is always measured in a straight sitting position.
- He/she should keep his/her head straight or slightly extended.
- The nose-clip is put on the lower half of the nose.
- The subject should slightly bite on the mouthpiece.
- Make sure that the subject carefully seals his/her lips around the mouthpiece.
- The subject's tongue should be below the mouthpiece.
- Dentures do not influence the measurement results.
- Tongue piercings should be removed prior to the measurement.


Activities that should be avoided prior to starting a measurement:

- Smoking within at least 1 h of testing
- Consuming alcohol within 4 h of testing
- Performing vigorous exercise within 30 min of testing
- Wearing clothing that substantially restricts full chest and abdominal expansion
- Eating a large meal within 2 h of testing.


9.2 Slow Spirometry (VC and IC)

Definition of VC and IVC

The VC is the volume change at the mouth between the position of full inspiration and complete expiration, expressed in litres at BTPS. The slow VC can be derived in two ways. The expiratory vital capacity (EVC) is the maximal volume of air from the point of maximal inhalation. The IVC or VCin is the maximal volume of air inhaled from the point of maximal exhalation, achieved by an unforced exhalation from end-tidal inspiration. These maneuvers are unforced, except at the point of reaching RV or TLC, respectively, where extra effort is required.




Tracing of tidal breathing followed by an inspiratory maneuver to total lung capacity (TLC) followed by a full but not forceful expiration to residual volume (RV) to record EVC or VCex; FRC: functional residual capacity.



Tracing of tidal breathing followed by an expiratory maneuver to residual volume (RV), followed by a full inspiration to total lung capacity (TLC) to record IVC or VCin. FRC: functional residual capacity; ERV: expiratory reserve volume.

Definition of IC

Inspiratory capacity (IC) is volume change recorded at the mouth when taking a full inspiration with no hesitation, from a position of passive end-tidal expiration, i.e. FRC, to a position of maximum inspiration, expressed in liters at BTPS. IC is an indirect estimate of the degree of lung hyperinflation at rest, and is useful to assess changes in FRC with pharmacological interventions and physical exercise.



Tracing of tidal breathing followed by an inspiratory maneuver to total lung capacity (TLC) to record inspiratory capacity (IC).

9.3 Perform a Slow Spirometry Measurement

Please observe the instructions for hygiene of your system.

To ensure optimal subject safety, we recommend only the use of accessories which are distributed with the MasterScope.



Make the proper preparations.
(see "**Basic Conditions prior to starting a Measurement**")

ERT - MasterScope CT

Current Patient PATIENT NUMBER 999926001 SEX Male DATE OF BIRTH 02JAN1965

PROTOCOL SIMULATION STUDY 01 CENTER # 999926

ERT
 Getting It Done. Right.

Logged in as DrHouse

Log off

Study Help

Calibration

Send Data

Tools

New DCR

Shutdown

Visit V1

Visit Summary

Visit Calendar

Patient Details

Search/New Patient

V1

COPD Medication Washout
 Questionnaire Reminder
Select Next Action
 Pre PFT
 Dosing Salbutamol
 Post PFT
 Inclusion

Reschedule Cancel

Select Next Action

Date and time	Operator	Action Type	Action details 1	Action details 2	Action details 3	Action details 4
06FEB2015 09:01:06	DrHouse	COPD Medication Washout	COPD Medication Washout Check: Yes			
06FEB2015 09:01:19	DrHouse	Questionnaire Reminder	Questionnaires completed by patient?: Yes			

Print

Next Action

Pre PFT

Continue With 'Pre PFT'

Center # 999926 | SystemId: 028933 | Study Version: 800154-SIMULATION STUDY 01-0 2-4 0 53 (04FEB2015 08:59:20) | Battery Power: 100% Remote: Disabled SYS# 800154-028933-000001 | 06FEB2015 09:01:25



Start Spirometry by pressing the <**Continue With 'Pre PFT'**> icon.

Please follow the notes on the screen.

ERT - MasterScope CT

Current Patient
PATIENT NUMBER 999926001
SEX Male
DATE OF BIRTH 02JAN1965

999926002

PROTOCOL SIMULATION STUDY 01
CENTER # 999926

ERT
Getting It Done. Right.

Visit V1 Visit Summary Visit Calendar Patient Details

Logged in as DrHouse

Log off

Study Help

Calibration

Send Data

Tools

New DCR

Shutdown

Search/New Patient

V1

COPD Medication Washout

Questionnaire Reminder

Pre PFT

Dosing Salbutamol

Post PFT

Inclusion

Reschedule Cancel

PFT Preparation - Please confirm

Patient Details: 999926001, 50 yrs, Male, 171 cm

Patient Position: Seated

Patient Reminder: Please use nose clip

Ambient Conditions: Room Temperature (Celsius) 22
Relative Humidity (Percent) 33
Absolute Pressure (Hectopascal) 999

Current User: Jim House (DrHouse)

Confirm your user identity to continue or click [Change user]

Use fingerprint to continue with measurement

OPTIONAL ACTIONS Cancel Change user Change Ambient Conditions

Center #: 999926 | Systemid: 028933 | Study Version: 800154-SIMULATION STUDY 01-0.2-4.0.53 (04FEB2015 08:59:20) | Battery Power: 100% Remote: Disabled SYS#: 800154-028933-000001 | 06FEB2015 10:59:49

Change Ambient Conditions

Please check and, if necessary, correct the ambient conditions shown on the screen by clicking <**Change Ambient Conditions**>.

Put your finger on the fingerprint sensor to continue.

As an alternative, click the <**Use finger print to continue with measurement**> icon to continue.

Use finger print to continue with measurement

Prior to starting the Spirometry measurement, the authorized investigator has to verify again.

Study SIMULATION STUDY 01 - SIMULATION

Fingerprint Temporary Password Remote Support



User Confirmation

Please confirm your identity as IM by using your fingerprint again!


Please put your finger on the fingerprint sensor again or click [Other User] to change user.

Other User

Any measurement starts with a zero adjustment of the pneumotach.



Continue with <OK>.



Do not breathe through the pneumotach and do not move the pneumotach during zero adjustment.

If the zeroing fails, a corresponding message will be displayed.

When the zero adjustment is complete, instruct the subject to close his/her nose with the nose-clip and to approach the mouthpiece. Make sure his/her lips seal around the mouthpiece.

Instruct the subject to hold the sensor in one hand.




The results of the Spirometry measurement depend on the subject's cooperation. Explain the breathing exercise to your subject thoroughly.

Spirometry V7.0.2.8 [SP-01]

T1111001

F1

To start a Slow Spirometry measurement, click this icon once or press F1.



F7

To see the results for this measurement, click this icon once or press F7 on the keyboard.

F8

To re-zero the Spirometry Sensor, click this icon once or press F8 on the keyboard.

F10

Upon completing the recommended number of measurements, click this icon once or press F10.

13:46

Copyright © eResearchTechnology GmbH 2010

45% | 00:17




Start with **<F1>**.
The subject should breathe normally until stable tidal breathing is shown.

Spirometry V6.11.7 [SP-02]

999260004

Spirogram



	Pred	Best	% Pred
IC [L]	2.745		

* Breathe normally for at least 7 breaths
* When F1 becomes visible again, press F1 again to proceed with IC

F10

15:52


Activate ERV/IC measurement

00:30




After some breaths, the <F1> icon will re-appear.


Wait until the baseline is stable (indicated by the dFRC and dVT bars).
Stability is achieved when both bars are green.




Wait longer



Cont.with <F1>
(2nd best)



Continue with <F1>

Note:

The volume signal must run horizontally after the measurement was started (subject not yet breathing into device). If you observe curves going up or down (see example) during a measurement, the pneumotach has not been properly zeroed.



In order to fix this, press <F7> "**Result**" to enter the Result phase and then press <F8> "**Zero**" to start zero adjustment again.



Make certain that the pneumotach is not moved and the subject is not breathing through the mouthpiece.



Hint: Close the opening of the pneumotach during zero adjustment with the palm of your hand.



Two vertical lines will automatically be set in the spirogram window indicating the space to be filled by the breathing maneuver.


9.3.1 Recording of Different SVC Maneuvers on the MasterScope



Depending on the study protocol one of the following maneuver needs to be performed to record the parameters needed – please also see the protocol-specific Short Guide.



IC - Inspiratory Capacity

The measurement starts with tidal breathing (establishing a stable baseline may require more than 5 tidal breaths). From tidal breathing the subject takes in a maximum inspiration to total lung capacity (TLC). When TLC is achieved, the subject's exhalation should be relaxed and the subject should come back to tidal breathing (no maximal expiration to ERV).



VCEx or VEX - Expiratory Vital Capacity

The measurement starts with tidal breathing (establishing a stable baseline may require more than 5 tidal breaths). From tidal breathing the subject takes in a maximum inspiration to total lung capacity (TLC). When TLC is achieved, a full but relaxed exhalation to residual volume (RV) should occur, followed by an inspiration back to tidal breathing to record EVC or VCex.





The **Time Ex (expiration time) of the maneuver** should last at least 6 seconds. For visualization, a bar appears which is slowly filled up in red color during expiration. After 3 or 6 seconds of expiration (children or adult), the filling color turns green and this test quality criterion is fulfilled.

The **EOT (End of Test)** bar shows the change in volume in the last second. As the bar shrinks and changes from red to green, the end of test quality criterion is fulfilled.

VCin or IVC - Inspiratory Vital Capacity

The measurement starts with tidal breathing (establishing a stable baseline may require more than 5 tidal breaths). From tidal breathing the subject exhales fully to residual volume (RV), followed by a full inspiration to total lung capacity (TLC) and back to tidal breathing to record IVC or VCin.



The Time Ex (expiration time) of the maneuver should last at least 6 seconds. For visualization, a bar appears which is slowly filled up in red color during expiration. After 3 or 6 seconds of expiration (children or adult), the filling color turns green and this test quality criterion is fulfilled.

The **EOT (End of Test)** bar shows the change in volume in the last second. As the bar shrinks and changes from red to green, the end of test quality criterion is fulfilled.



Press <F7> to end the test and see the results.



The quality of the slow spirometry measurement depends on the subject's cooperation. In order to assess reproducibility (repeatability) and thus the quality of cooperation, it is recommended to perform at least 3 trials.



Press <F1> to record another trial.



This traffic light indicates the repeatability of the trial:

Red: Less than 3 efforts were performed.

Yellow: 3 or more accepted efforts were performed and the repeatability criteria were not met.


Green: 3 or more accepted efforts were performed and the repeatability criteria were met.



A green traffic light does not imply that all ATS criteria for the test are met. The green light only signifies that repeatability has been established.



The best breathing maneuver (in our example the red one) will automatically be calculated and displayed on the screen.



Explanation of icons from left to right:

- Show previous trial
- Show all trials
- Show next trial



Up to 8 efforts can be performed.




You may now continue with the Forced Spirometry measurement by pressing <F2> (depends on study setup. In screenshot above <F2> is not shown but might be available in your trial.).



It is not possible to add Slow Spirometry efforts once you have started a Forced Spirometry effort.

9.4 Forced Spirometry Measurement

Parameter Definition



FEV1	Forced expiratory volume after 1 second
FVC	Forced expiratory vital capacity
PEF	Peak expiratory flow (peak flow)
FEF25-75	Mean maximal expiratory flow between 25 and 75% of FVC
FEV3	Forced expiratory volume after 3 seconds
FEV6	Forced expiratory volume after 6 seconds
FEV1/FVC	FEV1 in % of forced expiratory vital capacity
FEV1/FEV6	FEV1 in % FEV6

9.5 Perform a Forced Spirometry Measurement

Please observe the instructions for hygiene of your system.

To ensure optimal subject safety, we recommend only the use of accessories which are distributed with the MasterScope.



Make the proper preparations.
(see “**Basic Conditions prior to starting a Measurement**”)

ERT - MasterScope CT

Current Patient PATIENT NUMBER 999926001 SEX Male DATE OF BIRTH 02JAN1965

PROTOCOL SIMULATION STUDY 01 CENTER # 999926

ERT® Getting It Done. Right.

Visit V1 Visit Summary Visit Calendar Patient Details

Logged in as DrHouse

Log off

Study Help

Calibration

Send Data

Tools

New DCR

Shutdown

Search/New Patient

V1

COPD Medication Washout

Questionnaire Reminder

Select Next Action

Pre PFT

Dosing Salbutamol

Post PFT

Inclusion

Reschedule Cancel

Select Next Action

Date and time	Operator	Action Type	Action details 1	Action details 2	Action details 3	Action details 4
06FEB2015 09:01:06	DrHouse	COPD Medication Washout	COPD Medication Washout Check: Yes			
06FEB2015 09:01:19	DrHouse	Questionnaire Reminder	Questionnaires completed by patient?: Yes			

Print

Next Action

Pre PFT

Continue With 'Pre PFT'

Center # 999926 | SystemId: 028933 | Study Version: 800154-SIMULATION STUDY 01-0.2-4.0.53 (04FEB2015 08:59:20) | Battery Power: 100% Remote: Disabled SYS#: 800154-028933-000001 | 06FEB2015 09:01:25

Continue With 'Pre PFT'

Start Spirometry by pressing the <Continue With 'Pre PFT'> icon.

Please follow the notes on the screen.

ERT - MasterScope CT

Current Patient
PATIENT NUMBER 999926001
SEX Male
DATE OF BIRTH 02JAN1965

999926002

PROTOCOL SIMULATION STUDY 01
CENTER # 999926

ERT
Getting It Done. Right.

Visit V1 Visit Summary Visit Calendar Patient Details

Search/New Patient

V1

COPD Medication Washout

Questionnaire Reminder

Pre PFT

Dosing Salbutamol

Post PFT

Inclusion

PFT Preparation - Please confirm

Patient Details: 999926001, 50 yrs, Male, 171 cm

Patient Position: Seated

Patient Reminder: Please use nose clip

Ambient Conditions: Room Temperature (Celsius) 22
Relative Humidity (Percent) 33
Absolute Pressure (Hectopascal) 999

Current User: Jim House (DrHouse)


Confirm your user identity to continue or click [Change user]

Use fingerprint to continue with measurement

OPTIONAL ACTIONS

Cancel Change user Change Ambient Conditions

Center # 999926 | SystemId: 028933 | Study Version: 800154-SIMULATION STUDY 01-0-2-4.0.53 (04FEB2015 08:59:20) | Battery Power: 100% Remote: Disabled SYS#: 800154-028933-000001 | 06FEB2015 10:59:49



Please check and, if necessary, correct the ambient conditions shown on the screen.

Put your finger on the fingerprint sensor to continue.

As an alternative, click the **<Use finger print to continue with measurement>** icon to continue.



Prior to starting the Spirometry measurement, the authorized investigator has to verify again.

Study SIMULATION STUDY 01 - SIMULATION

Fingerprint Temporary Password Remote Support



User Confirmation

Please confirm your identity as IM by using your fingerprint again!


Please put your finger on the fingerprint sensor again or click [Other User] to change user.

Other User

Any measurement starts with a zero adjustment of the pneumotach.



Continue with <OK>.



Do not breathe through the pneumotach and do not move the pneumotach during zero adjustment.
If the zeroing fails, a corresponding message will be displayed.

When the zero adjustment is complete, instruct the subject to close his/her nose with the nose-clip and to approach the mouthpiece. Make sure his/her lips seal around the mouthpiece.

Instruct the subject to hold the pneumotach in one hand.



The results of the Spirometry measurement depend on the subject's cooperation. Explain the breathing exercise to your subject thoroughly.


Spirometry V7.0.2.8 [SP-01]

T1111001

 To start a Forced Spirometry measurement, click this icon once or press F2 on the keyboard.



 To see the results for this measurement, click this icon once or press F7 on the keyboard.

 To re-zero the Spirometry Sensor, click this icon once or press F8 on the keyboard.

 Upon completing the recommended number of measurements, click this icon once or press F10.

13:47

Copyright © eResearchTechnology GmbH 2010 45% 00:10



Start with <F2>.

The subject should breathe normally until stable tidal breathing is shown.




From tidal breathing, the subject is instructed to inhale as deeply as possible. Without interruption, the subject should now immediately exhale as fast and as much (FEV1) and as long (FVC) as possible. The maneuver is usually completed by a maximal inhalation (VC IN) and continue to normal breathing.



Please note: During the whole examination, the subject must not leave the mouthpiece.

The upper window shows the spirogram in volume-time display.
The lower window shows the flow-volume curve.



The **Time Ex (expiration time)** of the forced expiratory maneuver should last at least:

- 3 seconds for patients <10 years
- 6 seconds for patients ≥ 10 years.


For visualization, a bar appears:

After 3 seconds (first bar for patients < 10 years) or 6 seconds (second bar for patients ≥ 10 years) of expiration, the filling color turns green and this test quality criterion is fulfilled.

The **EOT (End of Test)** bar shows the change in volume in the last second. As the bar shrinks and changes from red to green, the end of test quality criterion is fulfilled.



Screen display after pressing <F7> "Result".



Explanation of icons from left to right:

- Show expiration only
- Show inspiration only
- Show F/V loop
- Switch to Tiffeneau
- Show quality errors




The quality of the flow-volume curve depends on the subject's cooperation. In order to assess reproducibility (repeatability) and thus the quality of cooperation, it is recommended to perform at least 3 trials. The result of the best and the second best trial for FEV1 and FVC may differ by ≤ 150 mL or $< 5\%$. For $FVC \leq 1$ L a difference of ≤ 100 mL is valid*.

* Literature:

- MR Miller et al. Series "ATS/ERS Task Force: Standardisation of Lung Function Testing", Standardisation of Spirometry, Eur Respir J 2005; 319-338 Copyright © ERS Journals Ltd. 2005

Screen display after three reproducible trials:



This traffic light indicates the repeatability of the trial:



Red: Less than 3 efforts were performed.

Yellow: 3 or more accepted efforts were performed and the repeatability criteria were not met.

Green: 3 or more accepted efforts were performed and the repeatability criteria were met.




A green traffic light does not imply that all ATS criteria for the test are met. The green light only signifies that repeatability has been established.



Up to 8 efforts can be performed.



Clicking the Quality Errors Button opens the quality errors window shown below. Quality errors of each effort are highlighted and described in this window.



	Pred	Best	% Pred	1st	2nd	3rd
Deselected by						
FEV1 [L]	3.847	3.478	90.40	3.478	3.363	
FVC [L]	4.760		Measured at: 19JUL2019 08:41:48			
FEV1/FVC	0.81		Deselected by System			
PEF [L/min]		617		617	613	642
Error Code		AC			C	IEFD

With mouse over on the effort number (here 3rd effort) the date/time of the measurement will be displayed and if the effort was deselected by the system or by the user (if applicable).

Colors above the effort number matching the color of the flow/volume loop for easy identifying the loop vs. the parameters. Dashed color indicates the effort was deselected while a solid colored line indicates the effort is valid.

	Pred	Best	% Pred	1st	2nd
Deselected by					
FEV1 [L]	3.847	3.478	90.40	3.478	3.363
FVC [L]	4.760	4.951	104.02	4.951	4.770
FEV1/FVC	0.81	0.70	86.64	0.70	0.71
PEF [L/min]		617		617	613
Error Code		AC			C
A - Less than 3 acceptable measurements C - FVC is not repeatable					

With mouse over of the error codes (here in best column error code A and C) the description of the error code will be displayed.

9.6 ATS Criteria (ATS/ERS 2005)

Forced Spirometry: Acceptability of individual efforts

- **A plateau was detected at the end of the expiration**
The volume change in the last second of expiration is less than 25 mL.
- **Expiration time is long enough (6 sec)**
The duration of expiration is greater than 6 seconds.
- **Back-extrapolated volume is not excessive**
The back extrapolated volume is less than 150 mL or less than 5% of FVC.
- **No late peak flow was detected**
The peak flow is reached earlier than 120 ms after extrapolated start of the maneuver. Indicates a poor effort.
- **No coughing was detected in the first part of the expiration**
Coughing during the first second of the maneuver affects the measured FEV1 value and other parameters.

Forced Spirometry: Repeatability of test set

- **No repeatability: Less than 3 acceptable measurements recorded**
- **FEV1 repeatability is unacceptable**
Difference in FEV1 from best to second best is greater than 150 mL.
- **FVC repeatability is unacceptable**
Difference in FVC from best to second best is greater than 150 mL.
- **PEF repeatability is unacceptable**
Difference in PEF from best to second best is greater than 0.67 L/sec.

Slow Spirometry: Repeatability of test set

- **IC repeatability is unacceptable**
Coefficient of Variation (CV) is greater than 6%.
- **VC repeatability is unacceptable**
Difference in VC from highest to second highest is greater than 150 mL.

Messages for repeatability and/or acceptability are evaluated. This is also recalculated after every new measurement performed.

In the MasterScope, acceptability and repeatability errors are implemented as error codes. These error codes are shown on the screen and on the PFT reports:

Forced Spirometry:

E code	Description
A	No repeatability: Less than 3 acceptable forced measurements
B	FEV1 repeatability is unacceptable
C	FVC repeatability is unacceptable
D	Expiration time was too short
E	No plateau was detected at the end of the expiration
F	Back extrapolated volume was too large
G	PEF repeatability is unacceptable
H	Late peak flow detected
I	Coughing was detected in the first part of the expiration

Slow Spirometry:


E code	Description
L	IC repeatability is unacceptable
M	VC repeatability is unacceptable
N	No repeatability: Less than 3 acceptable slow measurements
O	Expiration time was too short
P	No plateau was detected at the end of the expiration
Q	End-expiratory level of the tidal breathing was too variable

9.7 Unacceptable Efforts

Immediately after each effort, messages are displayed, which indicate whether the last effort should be disregarded (not accepted) by the operator and how to prevent such errors in the future. MasterScope automatically discards efforts which show the following issues:


- **Back-extrapolated volume too large**
- **Coughing during first second of test**
- **Delayed PEF (possible poor effort)**

In such a scenario, a message like the following appears displaying all of the criteria that are not met:




Depending on the study setup the system either allows the user to decide if the efforts should be kept or it will not allow to deselect the effort for the following errors:

- **Expiration time shorter than 6 sec.**
- **No plateau at end of expiration detected**



The automatic removal of unacceptable efforts is based on ATS/ERS guidelines. The investigator can accept an effort again if indicated by clicking on the deselected dot with the left mouse button. To disregard an effort, click on the colored dot with the left mouse button.



Discarded efforts are not counted for the minimum three efforts.

9.8 Selection of Best Parameters

The best Forced Spirometry parameters will be selected for reporting and displayed as follows:

- All parameters, except for FEV1 and FVC, are taken from the measurement with the highest sum of FEV1+FVC (best test).
- The best FEV1 is the maximum value of all accepted maneuvers.
- The best FVC is the maximum value of all accepted maneuvers.



The best FEV1 and FVC may come from different test curves.

9.9 Measurement Summary Result and Comments




After acceptable efforts have been obtained and no additional efforts are indicated, proceed by pressing <F10>.

The Summary screen shows the best and predicted values of the parameters FEV1 and FVC.

Best = Best values from all trials

Pred = Predicted value

%Pred = Best value in % of predicted value



Operator comment

Enter an operator comment in the respective entry field or select a preset comment and move it to the entry field by clicking the arrow key.

It is recommended to enter comments in the text field when applicable.

A comment is expected if the subject, for instance, was not able to perform acceptable and/or repeatable efforts.

Comments should be meaningful and entered in English. This will help the Overread department to analyze the data.

If further flow volume measurements are required to achieve acceptable results, press <Cancel>. This returns you to the measurement screen.

When the measurement results satisfy the requirements, press <Ok> to save the test and exit.

Cancel

Ok

9.10 Review Results and Edit Comments

At any time, measurement details can be reviewed by clicking on the corresponding tab and Info icon in the Visit Summary screen.

ERT - MasterScope CT

Current Patient: PATIENT NUMBER: 999926001, SEX: Male, DATE OF BIRTH: 02JAN1965

PROTOCOL: SIMULATION STUDY 01, CENTER #: 999926

Visit V1 | Visit Summary | Visit Calendar | Patient Details

Select Next Action

Date and time	Operator	Action Type	Action details 1	Action details 2	Action details 3	Action details 4
06FEB2015 09:01:06	DrHouse	COPD Medication Washout	COPD Medication Washout Check: Yes			
06FEB2015 09:01:19	DrHouse	Questionnaire Reminder	Questionnaires completed by patient?: Yes			
06FEB2015 11:00:58	DrHouse	Pre PFT	Best FEV1 at 11:01:38	FEV1: 3.78 L FVC: 4.96 L FEV1/FVC: 0.8	%Pred FEV1: 105 % %Pred FVC: 109 %	

Next Action

Dosing Salbutamol

OPTIONAL ACTIONS


Center #: 999926 | SystemId: 028933 | Study Version: 800154-SIMULATION STUDY 01-0.2-4.0.53 (04FEB2015 08:59:20) | Battery Power: 100% Remote Disabled SYS#: 800154-028933-000001 | 06FEB2015 12:02:20

Typical information provided with the measurement results are:


1. Parameter values from measurements
2. Flow Volume curves
3. ATS error codes for single measurements and best test.
4. Comments for this test, which can be amended.

9.11 Amending Spirometry Measurements

Once a spirometry test has been saved, it can be reopened and amended if no other workflow steps (inclusion, randomization, data transfer) have been completed.



In order to amend the last test, click the icon <**Amend PFT**>. The spirometry module will allow performing additional measurements.



Amending PFT is only applicable until the measurement results have been transferred via Data Transfer to the ERT Backend Center.

10. Resting ECG



MasterScope ECG allows to record a 12-lead resting ECG measurement.

ERT - MasterScope CT

Current Subject
 SUBJECT INITIALS: AB
 SUBJECT NUMBER: 99333002
 SEX: Female
 DATE OF BIRTH: 08FEB1958

PROTOCOL CENTER # 99333

ERT
 Getting It Done, Right.

Visit V101 | Visit Summary | Visit Calendar | Subject Details

Search/New Subject

V101

Started

Select Next Action

Pre ECG

FeNO

Pre PFT

Ipratropium Bromide

Reschedule | Cancel

OPTIONAL ACTIONS

Optional FeNO

Amend Medication | Washout

Shut down

Logged in as IM

Log off

Study Help

Calibration

Send Data

Tools

New DCR

Print

Select Next Action

Date and time	Operator	Action Type	Action details 1	Action details 2	Action details 3	Action details 4
- No Visit Action performed -						

Next Action

Pre ECG

Help text for Pre ECG (VT)

Continue With 'Pre ECG'

Center Number: 99333 | SystemId: 028904 | Study Version: 800286-800286-1.4.0-4.0.21 (20JUN2014 11:17:42) | Battery Power: 100% Remote: Disabled SYS#: 800286-028904-000004 | 23JUL2014 14:25:01




Click the <Continue With 'Pre ECG'> icon.

Following screen appears:

Please follow the notes on the screen.


Put your finger on the fingerprint sensor to continue.

As an alternative, click the **<Use finger print to continue with measurement>** icon to continue.

 **Use finger print to continue with measurement**

Prior to starting the ECG measurement, the authorized investigator has to verify again.

After the user confirmation, please follow the instructions given in chapter **"Performing an ECG Recording"**.




10.1 Information on ECG Recording

An electrocardiogram (ECG) is a graphic recording of the changes occurring in the electrical potentials (millivolt changes) at defined sites on the skin. The continuously changing electrical fields are the result of depolarization and polarization of the heart and are distributed in the body without any delay.

The electrical fields are caused by the cardiac cells, which are electrically polarized. The ECG is a graphic recording of cardiac electrical activity but is not a measure for cardiac pumping capacity (muscle strength).

10.1.1 The Waveform

Willem Einthoven (1860-1927), Professor of Physiology and Winner of the 1924 Nobel Prize, developed the ECG Standard Leads I, II and III, which are named after their inventor.



Einthoven named the prominent waves alphabetically P, Q, R, S, T and U.

The flat amplitudes P, T and U are called waves, Q, R and S are called peaks.

The P-wave represents the wave of depolarization that spreads from the atrium.

The Q, R and S peaks, also referred to as QRS-complex, represent the wave of depolarization from the ventricle.

The T-wave represents the repolarizations of the ventricle.

The U-wave is undefined.

10.1.2 ECG Leads

To minimize artifacts, the skin of the defined lead positions has to be prepared thoroughly.

Preparing the subject's skin:


1. Identify the (10) electrode sites on the torso by referring to the picture and description below.
2. Remove any hair from the electrode site using a razor.
3. Wipe oils from the electrode sites with an alcohol prep pad.
4. Remove any dead skin from the electrode sites with an abrasive cleaner.
Two to three moderate rubs at each site should be sufficient.



TIP: Electrodes should be stored in an air-tight container. Electrodes will dry out if not stored properly which will cause loss of adhesion and conductivity. Please note the storage conditions indicated on the electrode packaging.


Correct electrode placement is important for acquiring a successful ECG recording.

10.1.3 Chest ECG



- V1 Fourth intercostal space at the right sternal border
- V2 Fourth intercostal space at the left sternal border
- V3 Midway between V2 and V4
- V4 Fifth intercostal space at the left of the midclavicular line
- V5 Anterior axillary line at same horizontal level as V4
- V6 Mid-axillary line on same horizontal level as V4 and V5
- LA Left deltoid
- RA Right deltoid
- RL Right thigh
- LL Left thigh

10.1.4 Limb ECG




When connecting the electrodes to MasterScope ECG, the tiny and fast potential differences originating from the heart can be detected on the surface of the body between either two individual electrodes or between one individual electrode and a group of combined electrodes and recorded by MasterScope ECG.


The different measurement setups are commonly referred to as leads.
For a standard 12-lead ECG, four electrodes are placed at the limbs and six at the chest.

The 12 leads are:


Three bipolar limb leads: I, II and III (according to **Einthoven**)




Three unipolar limb leads: aVR, aVL and aVF (according to **Goldberger**)



Six unipolar chest leads: V1, V2, V3, V4, V5, V6 (according to **Wilson**)



In contrast to the limb leads, the chest leads have to be positioned precisely. The lead positions are internationally standardized.



- V1: 4th intercostal space, right sternal border
- V2: 4th intercostal space, left sternal border
- V3: midway between V2 and V4
- V4: 5th intercostal space, left mid-clavicular line
- V5: between V4 and V6, left anterior axillary line
- V6: level with V4, left mid-axillary line

10.1.5 Basic Conditions for ECG Recording

For high-quality ECG recording, certain criteria have to be met:



- Mentally prepare the subject for the examination in order to eliminate pain and consequently tachycardia and muscle tremor.
- Ambient temperature should be at least 23 °C to avoid shivering; make sure that the subject is lying comfortably on a suitable couch or bed and eliminate all sources of noise.
- Check the condition of your equipment to ensure proper signal sampling.
- Make sure that the chest electrodes are positioned according to international standards and pay attention to polarity of the limb and chest electrode cables.
- The subject should try to avoid movement during the measurement because this can lead to motion artifacts.



Only use original electrode cables delivered by ERT. If the wrong cable is used the defibrillation energy delivered to the patient can decrease, the device can be damaged, or electric shock to the operator or other persons occur.

10.1.6 Preparing for the Measurement

1. Attach electrodes

Procedure:

- Clean the subject's skin with a skin-sensitive agent to remove probable fatty residues.
Make sure that the skin is dry before applying the electrodes.
Only the supplied disposable electrodes are to be used:
Remove protecting foil and attach the electrode to the skin.

2. Connect the ECG amplifier electrode cables with the electrodes.


10.2 Performing an ECG Recording

A screen displaying the correct placement of the electrodes for the ECG recording will appear.

10.2.1 Electrode Test



Before an ECG recording is started, an automatic electrode contact measurement is performed. This measurement is started by pressing <F1>.



The electrode test program automatically checks whether the electrodes are in good contact with the skin. If a contact is poor, the respective electrode will immediately be indicated on the screen and directly at the electrode.


- good contact (only on the screen)
- bad contact (flashing)



If an electrode contact is indicated as poor, please check the respective contact.

The lead placement check will start automatically as soon as all electrodes are in good contact with the skin.


To enable an easy check of the electrode placement, the system will illuminate the electrodes in a fixed sequence.



Depending on the used ECG, the electrodes are visually checked ("running lights"), starting with "V1" and ending with "LL":




Chest ECG electrode sequence:	"e" shape
Limb ECG electrode sequence:	mirrored "S" shape



Pressing <OK> will start the recording of the ECG signals.

10.2.2 ECG Recording



The recorded ECG waveforms are displayed on the screen.



Scaling and Filter




Press <F4> to change between the scalings 5 mm/mV and 10 mm/mV.




Press <F5> to turn on/off the notch filter. Turn this on if there is noise seen on the ECG signal.

Before using the notch filter ensure there is no other interference causing the noise, e.g. electric interference next to the ECG amplifier.

10.2.3 ECG Evaluation



Inspect the recording quality of the signals. If the signals are free of drift and noise, press <F2> to analyse the last 10 seconds of the recording and halt the acquisition. HES® Interpretation will be performed and the results will be displayed on the screen.



Interpretation according to HES®

MasterScope ECG provides the Hannover ECG System HES®, which has been developed together with leading cardiologists all over the world.

Today, the HES® algorithm is considered as quasi standard for ECG recording and interpretation.



WARNING

Automatic interpretation of the ECG is not possible for pediatric subjects with an age below 16 years and for pacemaker subjects.



WARNING

A qualified physician has to reassess all MasterScope ECG measurements. An interpretation by MasterScope ECG is only significant when considered together with other clinical findings. ECG interpretation statements made by MasterScope ECG represent partial qualitative and quantitative information on the subject's cardiovascular condition and no therapy or drugs can be administered solely on the interpretation statements.



It is extremely important that the ECGs submitted have clean baselines and are free of any drifts. See "Recording Quality ECGs" on the next page.



If the quality of the recorded data is not sufficient, a new ECG recording should be started prior to pressing <F10> (= save and exit ECG recording). For this, press <F1> to delete the recorded data and restart ECG recording.




Press <F10> to save the current ECG recording and exit the program.

Leave Measurement Program

Comment

Save Cancel



Enter a comment and click <Save>.

Following screen appears:

ERT - MasterScope CT | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING MODE | TRAINING

Training Patient SUBJECT INITIALS AAA SUBJECT NUMBER T00011001 SEX Female DATE OF BIRTH 1900-01-01

PROTOCOL CENTER # 0001

ERT *Getting It Done. Right.*

Visit V101 Visit Summary Visit Calendar Subject Details

Search/New Subject

V101

Started

Select Next Action

Pre PFT

Ipratropium Bromide

Salbutamol - Albuterol

Post PFT

Reschedule

Select Next Action

Date and time	Operator	Action Type	Action details 1	Action details 2	Action details 3	Action details 4
08MAY2014 12:30:48	IM	Medication Washout				
08MAY2014 12:32:46	IM	Pre ECG	HR: 60 1/min RR: 1000 msec QT: 418 msec QTcF: 418 msec dP: 100 msec dPD: 200 msec dQRS: 88 msec			

Print

Next Action

FeNo

Continue With 'FeNo'

OPTIONAL ACTIONS + Optional FeNo + Optional ECG + Optional PFT

Center Number: 00011 | SystemId: 999999 | Study Version: 800286-800286-0.22.1-4.0.13 (28APR2014 11:17:47) | Remote: Disabled | SYS#: 800286-999999-000001 | Battery Power: 100% | 08MAY2014 12:33:06



Click this icon to display the detailed ECG report.