

DMS Pro Quick Set Up

The Eversense DMS Pro software is intended for clinic use. This sheet shows a short cut for setting up a single DMSPro account for your entire clinic. This short cut is not recommended if you have multiple users in the clinic. For complete instructions, see the DMS Pro User Guide (eversensedisabetes.com).

1 Sign Up for DMS Pro

Go to <https://uspro.eversensedms.com>.
Click on **Sign Up**.

Click on **Sign-Up**

2 Request an Eversense DMS Pro Account

Complete the form by entering in your clinic's information.

- Under "First Name" and "Last Name" it is recommended to either put your clinic name, or a primary health care provider in the clinic. (This is what your patients will see when they share data.)
- Select a "Title" from the dropdown, and complete the remaining required fields.

Notes:

- If you want to add individual providers to your Eversense DMS Pro clinic, see the Eversense DMS Pro User Guide.
- You cannot use the same email address for both an Eversense DMS Pro account and a personal Eversense DMS account.

Click **Submit** when form is completed.

3 You Will Receive Two Emails

- Account registration confirmation.
- Account approval email (may take up to 24 hours). Click on the link* in the email to complete account activation.
- Record important account information in the box to the right.

***The link expires after 72 hours, but you can request that a new email be sent.**

Clinic Name:

Email:


Password:

Clinic ID#: (see Step 4)


4 Sign in to the clinic DMS Pro account

The assigned clinic ID is at the top.

Note: Patients will need the Clinic Number when they send a request to share their CGM data with the clinic.



The sign-in page for Eversense DMS Pro. It features the Eversense logo at the top, followed by a language selector (English) and a link to sign up. Below this is a form with fields for Email and Password, both marked as mandatory. A 'Sign In' button is highlighted with a red box. At the bottom, there are links for 'Forgot your password' and 'Problem logging in? Contact Us'.




The clinic selection page. At the top, there is a blue header with a menu icon. Below the header, the text 'XYZ Clinic - 9' is displayed and highlighted with a red box. Underneath is a 'List of Patients' section with a search bar and a 'Search Patient' button. A refresh icon is also visible.

5 Inviting Patients to the Clinic


- On the navigation bar on the left, click **Patient** then **Invite Patient**.
- Enter the patient's Eversense email (the one used when setting up their mobile app and creating their Eversense account).
- Click **Search Eversense**.

If the patient's email is registered, the field to the right will be displayed. Complete all mandatory information, then hit **Send**. This generates an invitation email to the patient.


Note: If the entered email is not registered with Eversense DMS, you will receive an error message. Confirm with the patient what email address was used when setting up their account.



The navigation bar on the left side of the Eversense DMS Pro interface. It includes links for Home, Staff, Patients, Invite Patient (highlighted with a red box and a mouse cursor), Pending Requests, and Help.




The search form for Eversense. It has an 'Email' field with 'sample@email.com' entered. A 'Search Eversense' button is highlighted with a red box. Below the search field are 'Send' and 'Reset' buttons.



The patient details form. It shows the email 'sensedemo2@gmail.com' and a message 'User already exist in Eversense DMS'. Fields for MRN, Last Name (Brown), Date Of Birth (Dec 31, 1989), and Gender (Male) are visible. A 'Send' button with a checkmark is highlighted with a red box.

Once the patient accepts the invitation, the **✗** under Data Share will turn to **✓**. Click on the patient's last name to view their reports.



The patient list table in the Eversense DMS Pro interface. It shows a list of patients with columns for Patient last name, Patient first name, MRN, Gender, DOB, Provider Name, Status, and Data Share. The 'Data Share' column for the patient 'Kasodimos' is highlighted with a red box and a checkmark.

Patient last name	Patient first name	MRN	Gender	DOB	Provider Name	Status	Data Share
Kasodimos	SENS DEMO-Tommy	12345	Male	Jun 24, 1993	Barbara Montgomery, MD	Active	✓

Note: Patients may request to share their data. Pending requests from patients appear in the left navigation bar. Click on **Patients > Pending Requests**.

For additional information or support on Eversense DMS Pro, refer to your DMS Pro User Guide, contact your Eversense Representative, or contact Eversense Customer Care at:

1-844-SENS4U (736-7348)
Support@eversensediababetes.com



Continuous Glucose Monitoring System

MANUFACTURED BY:
SENSEONICS, INC.
20451 Seneca Meadows Parkway | Germantown
MD 20876-7005 | USA
301.515.7260 | global.eversensediababetes.com



DISTRIBUTED BY:
ASCENSIA DIABETES CARE US, INC.
5 Wood Hollow Road | Parsippany
NJ 07054 | USA
+1.862.225.2902 | ascensia.com