

User's Guide



Shop online at

omega.com®

Ω OMEGA®

omega.com

e-mail: info@omega.com

For latest product manuals:

omegamanual.info

ISO 9001
CERTIFIED
CORPORATE QUALITY

STAMFORD, CT

ISO 9001
CERTIFIED
CORPORATE QUALITY

MANCHESTER, UK

TCL Series Portable Sample Cooling Device



OMEGAnet® Online Service
omega.com

Internet e-mail
info@omega.com

Servicing North America:

U.S.A.:
ISO 9001 Certified

Omega Engineering, Inc., One Omega Drive, P.O. Box 4047
Stamford, CT 06907-0047 USA
Toll Free: 1-800-826-6342
FAX: (203) 359-7700

TEL: (203) 359-1660
e-mail: info@omega.com

Canada:

976 Bergar
Laval (Quebec), Canada H7L 5A1
Toll-Free: 1-800-826-6342
FAX: (514) 856-6886

TEL: (514) 856-6928
e-mail: info@omega.ca

For immediate technical or application assistance:

U.S.A. and Canada: Sales Service: 1-800-826-6342/1-800-TC-OMEGA®
Customer Service: 1-800-622-2378/1-800-622-BEST®
Engineering Service: 1-800-872-9436/1-800-USA-WHEN®

Mexico:

En Español: 001 (203) 359-7803
info@omega.com.mx

FAX: (001) 203-359-7807
e-mail: espanol@omega.com

Servicing Europe:

Benelux:

Managed by the United Kingdom Office
Toll-Free: 0800 099 3344
FAX: +31 20 643 46 43

TEL: +31 20 347 21 21
e-mail: sales@omega.nl

Czech Republic:

Frystatska 184
733 01 Karviná, Czech Republic
Toll-Free: 0800-1-66342
FAX: +420-59-6311114

TEL: +420-59-6311899
e-mail: info@omegashop.cz

France:

Managed by the United Kingdom Office
Toll-Free: 0800 466 342
FAX: +33 (0) 130 57 54 27

TEL: +33 (0) 161 37 29 00
e-mail: sales@omega.fr

Germany/Austria:

Daimlerstrasse 26
D-75392 Deckenpfronn, Germany
Toll-Free: 0 800 6397678
FAX: +49 (0) 7056 9398-29

TEL: +49 (0) 7059 9398-0
e-mail: info@omega.de

United Kingdom:
ISO 9001 Certified

OMEGA Engineering Ltd.
One Omega Drive, River Bend Technology Centre, Northbank
Irlam, Manchester M44 5BD England
Toll-Free: 0800-488-488
FAX: +44 (0)161 777-6622

TEL: +44 (0)161 777-6611
e-mail: sales@omega.co.uk

It is the policy of OMEGA Engineering, Inc. to comply with all worldwide safety and EMC/EMI regulations that apply. OMEGA is constantly pursuing certification of its products to the European New Approach Directives. OMEGA will add the CE mark to every appropriate device upon certification.

The information contained in this document is believed to be correct, but OMEGA accepts no liability for any errors it contains, and reserves the right to alter specifications without notice.

WARNING: These products are not designed for use in, and should not be used for, human applications.

BEFORE USE:

Please read the following instructions:



Carefully read the operating manual first before operating the device.



For indoor use only



Ambient temperature range +5°C to +40°C



Use in a well-ventilated area.



Relative humidity not exceeding 80%



Power supply fluctuation not exceeding 10%

Overview

The **TCL-11 and TCL-12** are ideal for incubating any number of small laboratory sample vessels or micro tube sizes at sub-ambient temperatures for extended periods of time. They are an excellent alternative cooling device that replaces the need for those inconvenient and messy ice buckets that can often leave sample tubes floating in melting wet ice. This can lead to valuable sample loss and potential contamination. Each TCL-11 and TCL-12 are supplied with a pack of special ceramic coated metal beads designed to help maintain the cooling temperature while supporting the sample vessels during use. Beads are chemically inert and may be autoclaved (in a separate autoclavable container) if needed. Beads may also be pre-chilled prior to use for quicker cooling times. The TCL-11 and TCL-12 are manufactured to fit standard aluminum block tube inserts and includes a protective lid which helps to maintain the cool temperature.

Two models to choose from:

TCL-11 is the standard analog cooling device. Set to maintain temperature at 0°C.

TCL-12 is the digital (display & setting) cooling device. Temperature range can be easily adjusted by the user from 0°C up to ambient 25°C.

Simply assemble the electrical cord and adapter, plug in and turn on your TCL-12. Adjust the unit to its appropriate cooling temperature. Place sample vessels into the TCL-12.

No mess, no sample loss, no potential contamination.

For additional information on other quality labware accessories that might be used with your new TCL-11 and TCL-12 see below:

TCL-10-B1	Aluminum Block Insert 0.2ml tubes
TCL-10-B2	Aluminum Block Insert 0.5ml tubes
TCL-10-B3	Aluminum Block Insert 1.5ml tubes
TCL-10-B4	Aluminum Block Insert 2.0ml tubes
TCL-10-B5	Aluminum Block Insert Combination tubes
TCL-10-CB	Ceramic Beads




Figure 1: Overview of the TCL-11 (1) case; (2) lid; (3) front panel




Figure 2: Overview of TCL-12 (Digital). (1) case; (2) lid; (3) front panel with digital display

Operation

Before Use. Please make sure you have carefully read this operating manual. If there are any questions related to the proper use of TCL-11 and TCL-12 please contact Omega Engineering.

Contact with low (cold) temperatures and vessels may cause cold burns. Please use care and wear protective gloves to protect hands. Do not place liquids directly into the unit. The TCL-11 and TCL-12 are rated to operate at 90-230V, output 12V, 5A, 60 W.

TCL-11: simply assemble the electrical cord to its adapter. Plug into the Base Unit and power outlet. Turn on the TCL-11. It will cool to a final set temperature of 0°C.




Figure 3: The rear view of TCL-11 and TCL-12. (4) switch (5) power inlet

TCL-12. Once the power is on, the TCL-12 shows the actual temperature of the sample chamber. Press and hold “mode” key, it shows the set temperature, use the “up” or “down” key to increase or decrease the setting temperature. After programming desired temperature, release the “mode” key, and the actual temperature will be shown. The TCL-12 will begin to ramp to set temperature.




Figure 4: Overview of the front panel of TCL-12. (1) digital display; (2) mode (set key); (3) down key; (4) up key

Post Use. The ceramic coated beads and any block insert may remain cold for some time. Allow to warm up to ambient temperature before removing from the TCL-11 and TCL-12 and make sure to wear protective gloves. Avoid handling cold objects with wet hands.

Maintenance. Before cleaning always be sure to disconnect the TCL-11 and TCL-12 from the power source. The outer casing of the Units may be wiped clean with a damp cloth, mild detergent or ethanol. Do not immerse the TCL-11 or TCL-12 in water or solvents. Do not use aggressive solvents or abrasive cleaners. Please check chemical compatibility with wetted parts before using any other decontamination method or lab solvent. Always dispose of any contaminated cleaning articles properly & safely.

WARRANTY/DISCLAIMER

OMEGA ENGINEERING, INC. warrants this unit to be free of defects in materials and workmanship for a period of **36 months** from date of purchase. OMEGA's WARRANTY adds an additional one (1) month grace period to normal **three (3) year product warranty** to cover handling and shipping time. This ensures that OMEGA's customers receive maximum coverage on each product.

If the unit malfunctions, it must be returned to the factory for evaluation. OMEGA's Customer Service Department will issue an Authorized Return (AR) number immediately upon phone or written request. Upon examination by OMEGA, if the unit is found to be defective, it will be repaired or replaced at no charge. OMEGA's WARRANTY does not apply to defects resulting from any action of the purchaser, including but not limited to mishandling, improper interfacing, operation outside of design limits, improper repair, or unauthorized modification. This WARRANTY is VOID if the unit shows evidence of having been tampered with or shows evidence of having been damaged as a result of excessive corrosion; or current, heat, moisture or vibration; improper specification; misapplication; misuse or other operating conditions outside of OMEGA's control. Components in which wear is not warranted, include but are not limited to contact points, fuses, and triacs.

OMEGA is pleased to offer suggestions on the use of its various products. However, OMEGA neither assumes responsibility for any omissions or errors nor assumes liability for any damages that result from the use of its products in accordance with information provided by OMEGA, either verbal or written. OMEGA warrants only that the parts manufactured by the company will be as specified and free of defects. OMEGA MAKES NO OTHER WARRANTIES OR REPRESENTATIONS OF ANY KIND WHATSOEVER, EXPRESSED OR IMPLIED, EXCEPT THAT OF TITLE, AND ALL IMPLIED WARRANTIES INCLUDING ANY WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY DISCLAIMED. LIMITATION OF LIABILITY: The remedies of purchaser set forth herein are exclusive, and the total liability of OMEGA with respect to this order, whether based on contract, warranty, negligence, indemnification, strict liability or otherwise, shall not exceed the purchase price of the component upon which liability is based. In no event shall OMEGA be liable for consequential, incidental or special damages.

CONDITIONS: Equipment sold by OMEGA is not intended to be used, nor shall it be used: (1) as a "Basic Component" under 10 CFR 21 (NRC), used in or with any nuclear installation or activity; or (2) in medical applications or used on humans. Should any Product(s) be used in or with any nuclear installation or activity, medical application, used on humans, or misused in any way, OMEGA assumes no responsibility as set forth in our basic WARRANTY/DISCLAIMER language, and, additionally, purchaser will indemnify OMEGA and hold OMEGA harmless from any liability or damage whatsoever arising out of the use of the Product(s) in such a manner.

RETURN REQUESTS/INQUIRIES

Direct all warranty and repair requests/inquiries to the OMEGA Customer Service Department. BEFORE RETURNING ANY PRODUCT(S) TO OMEGA, PURCHASER MUST OBTAIN AN AUTHORIZED RETURN (AR) NUMBER FROM OMEGA'S CUSTOMER SERVICE DEPARTMENT (IN ORDER TO AVOID PROCESSING DELAYS). The assigned AR number should then be marked on the outside of the return package and on any correspondence.

The purchaser is responsible for shipping charges, freight, insurance and proper packaging to prevent breakage in transit.

FOR **WARRANTY** RETURNS, please have the following information available BEFORE contacting OMEGA:

1. Purchase Order number under which the product was PURCHASED,
2. Model and serial number of the product under warranty, and
3. Repair instructions and/or specific problems relative to the product.

FOR **NON-WARRANTY** REPAIRS, consult OMEGA for current repair charges. Have the following information available BEFORE contacting OMEGA:

1. Purchase Order number to cover the COST of the repair,
2. Model and serial number of the product, and
3. Repair instructions and/or specific problems relative to the product.

OMEGA's policy is to make running changes, not model changes, whenever an improvement is possible. This affords our customers the latest in technology and engineering.

OMEGA is a registered trademark of OMEGA ENGINEERING, INC.

© Copyright 2009 OMEGA ENGINEERING, INC. All rights reserved. This document may not be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form, in whole or in part, without the prior written consent of OMEGA ENGINEERING, INC.

Where Do I Find Everything I Need for Process Measurement and Control? **OMEGA...Of Course!**

Shop online at omega.comSM

TEMPERATURE

- ☑ Thermocouple, RTD & Thermistor Probes, Connectors, Panels & Assemblies
- ☑ Wire: Thermocouple, RTD & Thermistor
- ☑ Calibrators & Ice Point References
- ☑ Recorders, Controllers & Process Monitors
- ☑ Infrared Pyrometers

PRESSURE, STRAIN AND FORCE

- ☑ Transducers & Strain Gages
- ☑ Load Cells & Pressure Gages
- ☑ Displacement Transducers
- ☑ Instrumentation & Accessories

FLOW/LEVEL

- ☑ Rotameters, Gas Mass Flowmeters & Flow Computers
- ☑ Air Velocity Indicators
- ☑ Turbine/Paddlewheel Systems
- ☑ Totalizers & Batch Controllers

pH/CONDUCTIVITY

- ☑ pH Electrodes, Testers & Accessories
- ☑ Benchtop/Laboratory Meters
- ☑ Controllers, Calibrators, Simulators & Pumps
- ☑ Industrial pH & Conductivity Equipment

DATA ACQUISITION

- ☑ Data Acquisition & Engineering Software
- ☑ Communications-Based Acquisition Systems
- ☑ Plug-in Cards for Apple, IBM & Compatibles
- ☑ Datalogging Systems
- ☑ Recorders, Printers & Plotters

HEATERS

- ☑ Heating Cable
- ☑ Cartridge & Strip Heaters
- ☑ Immersion & Band Heaters
- ☑ Flexible Heaters
- ☑ Laboratory Heaters

ENVIRONMENTAL MONITORING AND CONTROL

- ☑ Metering & Control Instrumentation
- ☑ Refractometers
- ☑ Pumps & Tubing
- ☑ Air, Soil & Water Monitors
- ☑ Industrial Water & Wastewater Treatment
- ☑ pH, Conductivity & Dissolved Oxygen Instruments