GenePix[®] Professional 4200A

MICROARRAY SCANNER

User's Guide

Part Number 2500-0153 Rev C August 2005 Printed in USA

Copyright 2005 Axon Instruments / Molecular Devices Corp.

No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, microfilming, recording, or otherwise, without written permission from Molecular Devices Corp.

QUESTIONS? See Axon's Knowledge Base: **<http://support.axon.com>**

PLEASE READ!!!! SAFETY

There are important safety issues that you must take into account when using this instrument. Please carefully read the safety warnings starting on page 37 before you use this instrument.

VERIFICATION

This instrument is extensively tested and thoroughly calibrated before leaving the factory. Nevertheless, researchers should independently verify the basic accuracy of the instrument using suitable test samples.

WARNING

If this equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

DISCLAIMER

This equipment is not intended to be used, and should not be used, in human experimentation or applied to humans in any way.

Table of Contents

Chapter 1 Introduction	1
Chapter 2 Installation	3
About This User's Guide.....	3
Computer Requirements.....	3
Setting Up and Moving the GenePix scanner	4
Preparing to Install the GenePix Scanner.....	4
Quick Overview of Hardware and Software Installation	5
Installing the SCSI Card.....	5
Installing the GenePix Pro Software	6
Connecting the GenePix Scanner to the Computer	7
Windows 2000 and Windows XP Installation	7
Connecting the Security Key.....	8
Troubleshooting Installation	8
SCSI Card Installation Problems.....	8
Chapter 3 Operation	9
Inserting a Slide.....	9
Functional Checkout	10
Installing Additional Filters	10
Chapter 4 Imaging System Hardware	15
Discussion of Hardware Components	15
Optics	15
Mechanical	20
Electronic	21
General Items of Interest.....	22
Spatial Resolution	22

Dynamic Range and Detection Limit 25

Specifications 25

Technical Assistance..... 27

Warranty and Repair Service 29

Licensing Notice..... 33

Declaration of Conformity..... 35

Important Safety Information..... 37

Index 57

Chapter 1

Introduction

The GenePix[®] Professional 4200A is the latest microarray scanner from Axon Instruments. It has several enhanced features for high performance quantitative imaging set by the benchmark GenePix 4000 series scanners. Key features of the 4200A are:

- Up to four internally installed laser excitation sources;
- Gradient neutral density wheel permitting 5–100% laser transmission;
- Data Scan resolution from 5 μm to 100 μm ;
- Fast Preview Scan at 40 μm ;
- Uses 16-bit ultra-low noise digitization;
- Dynamic range of 10^4 ;
- Easy-access 16-position emission filter wheel with standard filters for each laser installed at the factory;
- Full integration with GenePix Pro and Acuity software.

Axon Instruments has consulted widely among the scientific pioneers of microarray techniques to produce an instrument that fulfills all their image acquisition and primary analysis needs.

Chapter 2

Installation

About This User's Guide

This User's Guide focuses on the technical and theoretical aspects of the GenePix Professional 4200A hardware. We recommend that you use this Guide in conjunction with the GenePix Pro Software User's Guide, and the on-line Help within GenePix Pro (press <F1>). For information on setting up your GenePix Professional 4200A scanner and software, please refer to Chapter 1. For a guided tutorial on scanning your first slide, see Chapter 4 of the GenePix Pro Software User's Guide.

Computer Requirements

Please refer to the Axon Instruments website for up-to-date information.

Should you choose to purchase your own computer, Axon recommends that due consideration be made for the selection of a quality supplier and/or components.

Setting Up and Moving the GenePix scanner

Placement

Place the GenePix Professional 4200A scanner onto a secure laboratory bench top or reinforced table. The scanner weighs approximately 45 kg (99 lb) and requires two persons to lift. *Carefully* place the scanner into its final position; do not slide it into position as this may result in internal damage or misalignment of the optics.

Moving the Scanner

Always use the Park Scanner utility prior to moving the scanner, whether making a short move to another location or when shipping the scanning long distance. The Park Scanner utility is accessed by opening the GenePix Utilities program, which can be found in the Start Menu / GenePix Pro menu or by directly running the GP_Uutilities.exe executable file located in the GenePix Pro installation directory.

Instructions:

1. Close GenePix Pro, if open and start GenePix Utilities program.
2. Press Park Scanner button.
3. Wait for status message indicating that scanner has been successfully parked.
4. Immediately turn off scanner power; do not reopen GenePix Pro as this will un-park the scanner.
5. Place one or two pieces of adhesive tape to secure the door closed.

Failure to park the scanner may result in misalignment or damage of the optical system.

Preparing to Install the GenePix Scanner

If you have purchased a GenePix scanner *without the optional host computer* please follow the instructions below. If you purchased a GenePix scanner with a bundled host computer, the necessary hardware and software are already installed. If so, skip to “Connecting the GenePix Professional 4200A to the Computer”.

Quick Overview of Hardware and Software Installation

Here is a brief outline of how to get started with your GenePix scanner. Detailed instructions follow immediately below.

If you have purchased a GenePix scanner without the optional host computer:

1. Install the SCSI card (see below.)
2. Install GenePix Pro and attach the security key to parallel or USB port.

For all users:

3. Attach the GenePix scanner to the computer with the SCSI cable.
4. Attach the power cord to the GenePix scanner.
5. Turn on the GenePix scanner.
6. Turn on the computer.

Installing the SCSI Card

Included with your GenePix scanner is a SCSI adapter card for your host computer. This card is used to transfer data and instructions between the host computer and the GenePix scanner.

The SCSI card that is included with the GenePix scanner is a PCI-based Plug'n'Play card from AdvanSys that will install easily on most computers. We recommend that you use the included SCSI card.

To install the AdvanSys SCSI card:

1. The computer must have an available PCI slot.
2. Turn off computer and unplug the power cord. Remove the computer cover according to the computer manufacturer's instructions.

3. Remove the AdvanSys SCSI adapter from bag; do not touch the metal traces or contacts, and carefully insert into the open PCI slot. Do not force into place, it should go in with a gentle push. Screw the endplate to the chassis.
4. Replace computer cover, plug in power cord and turn on the computer. Do not connect GenePix scanner yet.
5. Plug and Play operating systems such as Window 2000 and Windows XP will automatically recognize the Advansys SCSI adapter and install the correct driver.
6. Install GenePix Pro software (see below.)
7. Close Windows and turn off the computer. Attach the SCSI cable to the SCSI adapter and the GenePix scanner (the scanner should also be powered off; see below.)
8. Turn on GenePix scanner first and wait for approximately 15 seconds before turning on computer.

If you already have a SCSI card installed on your host computer and do not want to replace it with the included SCSI card, you must verify that it is a SCSI-2 device that is Win-ASPI compliant. You must also ensure that you have the proper cabling to connect the SCSI card to the GenePix scanner. While it is possible that an existing SCSI card may work with the GenePix scanner, technical support is provided only for the SCSI card that is bundled with the GenePix. If you are using the GenePix scanner on a computer with multiple SCSI devices please ensure that there are no SCSI ID conflicts. At the factory, the GenePix scanner SCSI ID is set to 4.

Installing the GenePix Pro Software

To install the GenePix Pro software, please consult the instructions in the GenePix Pro manual.

Connecting the GenePix Scanner to the Computer

Caution: Failure to follow these instructions in the correct order can lead to damage to your instrument.

Before connecting the GenePix scanner to your SCSI connector, please ensure that your computer has been turned off, and that the GenePix scanner power is turned off. Once the computer has been turned off, connect the 50-pin data cable between the SCSI adapter of your computer and the GenePix scanner hardware. The SCSI connector and power connector are located at the back of the GenePix scanner.

Once the SCSI connection has been made, connect the power cord to the GenePix scanner. Turn on the power supply, and then turn on your computer.

Windows 2000 and Windows XP Installation

If the connections have been made correctly *and* GenePix Pro software has already been installed, when you restart your computer Windows will show a message that indicates that “New Hardware” has been detected. Windows 2000 will automatically install the necessary driver for the GenePix scanner.

Windows XP requires that you manually install the driver using the *Add New Hardware* Wizard. This wizard starts any time a new Plug and Play device is detected. Working through the wizard, select the option to “Install the software automatically (Recommended)”. Press *Next* and then press *Finish* when prompted. Repeat the procedure seven additional times.

If you are connecting the GenePix scanner to a computer that was purchased with the GenePix scanner, the driver is already installed.

Connecting the Security Key

The GenePix Pro software requires a security key to operate. If the key is absent you will be able to open images in GenePix Pro, but you will not be able to save your settings files or save your analysis results. In addition, you can analyze only 1000 features without a key.

To install the GenePix Pro key, please consult the instructions in the GenePix Pro manual.

Troubleshooting Installation

SCSI Card Installation Problems

If you followed the above instructions and still require further assistance on installation of the included SCSI card, contact Technical Support. If you are using a SCSI card that was not provided with your GenePix scanner, you should contact the supplier of that SCSI card for technical support.

GenePix Pro automatically detects the GenePix Professional 4200A scanner. If this does not occur please verify the following items:

- The data cable between the computer and the scanner is securely connected.
- The GenePix scanner has power (the green power light at the front of the scanner will be glowing).
- The SCSI card is operational (see above).
- The GenePix SCSI ID 4 does not conflict with any other SCSI devices that you may have already installed on your computer (*e.g.*, a Jaz drive).
- The scanner power was turned on *before* booting the computer.

Chapter 3

Operation

Inserting a Slide

To scan a slide, first lift the door. The slide holder will move to its forward-most position for loading. Gently lift the tab on the front of slide holder to raise the top half. With one finger, hold open the flat-spring positioner on the inside left edge of the slide holder. Drop the slide *face down* into the slide holder frame. The slide should be supported by the three sapphire balls on the frame. A slight flick of the flat-spring positioner will ensure proper positioning of the slide. Do not push the slide into position. Gently close the slide holder lid, and close the door. You are now ready to scan your sample. See “Image Acquisition” in the Tutorial section of the *GenePix Pro User’s Guide* for scanning instructions.

Note: The slide holder is a precision component designed to ensure proper focusing and field uniformity. Improper handling may damage the slide holder and affect imaging performance. Never touch the slide holder while it is moving. Never force the slide holder closed, or apply significant pressure to it for any reason. Do not rub or scratch the sapphire balls.

Functional Checkout

The Functional Checkout procedure is to help you confirm that your GenePix Professional 4200A microarray scanner is performing according to the specifications set by Axon Instruments, Inc. A similar, but much more detailed procedure was performed at the factory prior to shipping this unit to you. We recommend that you follow the procedure once upon receiving the instrument and then any time after that if you wish to check the performance of the instrument.

The GenePix 4200A microarray scanner has been designed to provide robust, long-lived microarray scanning for years of normal use. However, if the images you acquire on the GenePix scanner do not appear correct and if you believe there is a problem with the scanner, it is important to follow these instructions before contacting Technical Support:

1. Start the GenePix Pro software and click the Report tab.
2. Click the Functional Checkout link.
3. Follow the Functional Checkout procedure.

GenePix Pro automatically logs hardware performance into the Hardware Diagnostics Report every two hours during operation.

Installing Additional Filters

Your GenePix 4200A Professional scanner comes with one standard filter for each installed laser. The standard filters are not accessible by the user. You can install additional filters in any of the twelve remaining filter positions.

The installation process involves two steps:

1. Install the filters in the scanner;
2. Record the information for the new filters in GenePix Pro.

1. Install the filter(s) in the scanner

Turn off power to the instrument before opening the filter wheel compartment.

1. Use a 5/32" hex wrench (provided with the optional filter kit) to access the filter wheel compartment on the front of the scanner (Figure 1). Turn the hex wrench counterclockwise to release the filter door latch. Pull the handle to remove the filter cartridge (Figure 2).

Figure 1: Opening the filter compartment.

Figure 2: Removing the filter cartridge.

2. Turn the wheel to align any open position with the hole in the top of the filter cartridge. The number of that filter position is visible in the smaller adjacent hole. Note the number of the filter position you will use: this number must be recorded in GenePix Pro (see below).
3. **Handle filters by the edges only.** The rim of the filter is marked with an arrow indicating the direction of light propagation through the filter. Drop the filter into the hole with the **arrow pointing down**.
4. Using the cylindrical wrench, screw a lock-ring into the hole to gently secure the filter. Do not over-tighten.

Figure 3: Tightening the lock-ring.

5. Repeat for any additional filters.
6. When you have installed all additional filters, insert the filter cartridge into the instrument. Close the door and turn the hex key clockwise to secure the latch.

2. Record the information for the new filters in GenePix Pro

1. Turn on the scanner and computer. If GenePix Pro is running, shut it down.
2. From the Windows Start menu, go to *Programs / Axon Laboratory / GenePix Pro / Filter Tool*.

Installed	Formal name	Alternate name
☒ Filter 1	670DF40	Cy5
☒ Filter 2	530DF30	Alexa 488
☐ Filter 3		
☐ Filter 4		
☐ Filter 5		
☐ Filter 6		
☐ Filter 7		
☐ Filter 8		
☒ Filter 9	575DF35	Cy3
☐ Filter 10		
☐ Filter 11		
☐ Filter 12		
☐ Filter 13		
☐ Filter 14		
☐ Filter 15		
☐ Filter 16		

Buttons: Save, Close, Help

3. In the *Filter Tool* program, click the checkbox for the filter position you are adding.
4. Enter the midpoint of the wavelength range that is transmitted by the filter. For example, a 530DF30 filter transmits a 30 nm band of light centered at 530 nm (*i.e.*, it transmits 515 – 545 nm), therefore the midpoint of this filter is entered as “530”. For a long-pass filter, enter the peak emission of the dye you will use with that filter.
5. In GenePix Pro you can display either the filter’s *Formal name* (*e.g.*, “530DF30”) or an *Alternate name* (*e.g.*, “Alexa 488” or “Blue Filter”). Enter both names here so that each user can choose their preferred nomenclature in GenePix Pro.
6. After you have entered this information for each new filter, click *Save*.
7. Close the *Filter Tool* program.

When you open GenePix Pro, the new filters are listed in the Hardware Settings dialog box.

Chapter 4

Imaging System Hardware

The GenePix Professional 4200A scanner is one of the most sensitive scanners currently available. This chapter introduces some of the fundamental hardware concepts needed to understand the operation of the GenePix scanner.

Discussion of Hardware Components

When discussing the hardware components of the GenePix scanner it is convenient to reduce the system to its optical, mechanical and electrical components.

Optics

Fundamentally, the GenePix scanner uses a laser-excitation-based fluorescence scanning and imaging system. The optical system can be considered in terms of excitation light, mirrors, lenses, filters and photo-detection.

The optical lens is located beneath the slideholder, therefore the specimen must be placed feature-side-down into the slide holder. GenePix Pro software automatically displays the image as if you were viewing the specimen from the top view; *i.e.*, the upper-left corner of the slide is also the upper-left corner in the acquired image.

Excitation

In the GenePix scanner standard configuration, laser excitation is provided by individual 635 nm and 532 nm lasers. These wavelengths correspond to the ideal wavelengths used to excite the fluorophores Cy3 and Cy5 (Amersham Biosciences), which are currently the most widely used fluorescent dyes for microarray imaging. The individual lasers have been selected because of their superior optical performance and reliability. Since the performance of such lasers is often sensitive to external temperature fluctuations, the GenePix scanner uses an active temperature stabilization design to minimize temperature-based laser fluctuations. Furthermore, each laser has stabilized feedback output to ensure that photon output is constant.

The GenePix Professional 4200A scanner can be configured with two additional lasers for a total of four installed laser light sources. A blue laser with 488 nm excitation can be used with fluorescein and similar dyes; and a yellow laser with 594 nm excitation can be used with Texas Red and similar dyes. The Alexa Fluor™ dyes from Molecular Probes offers a number of dyes compatible with all four lasers excitation sources available in the GenePix 4200A scanner.

Laser excitation sources are often used in scanning systems instead of conventional wide-band white light sources such as a xenon or a mercury light source. The reasons for this are quite simple: conventional wide-band light sources do not deliver high enough photon density at the appropriate wavelength for optimum signal-to-noise ratio. In addition, it is often difficult to use optical filters to obtain spectrally pure, narrow-band excitation light.

Optical Path

Excitation laser light is directed onto the slide after passing through a series of filters and mirrors. While laser light is by definition very narrow band, the GenePix scanner employs additional optical filters to ensure that no

spurious excitation light is directed onto the slide. If the laser light impinges on an appropriate fluorophore bound on the slide, emission light of a longer wavelength is emitted (the difference between the emission and excitation wavelengths is called the Stokes shift). These emission photons are directed back through the optical system where they pass through another set of bandpass filters (see below) before reaching the photodetector.

Emission Filters

Light from the slide is filtered before reaching the photodetectors, to ensure that as far as possible, the detectors measure only emission photons from the appropriate fluorescent dye. To achieve this, the filter's 'window' is matched as closely as possible to the peak emission spectrum of the selected dye, to maximize acceptance of the desired emission light and at the same time exclude stray laser excitation light, and light emitted by non-target fluorescence (autofluorescence, reagent background etc.).

The GenePix scanner comes with standard emission filters installed, one for each laser installed at the time of purchase. Each standard emission filter is optimized for the emission spectrum of the respective dye. The factory-installed Cy5 ('red') filter is a bandpass filter with cut-on frequency at a wavelength of 655 nm, and cut-off frequency at 695 nm. This is ideal for selecting the emission photons of the red dye Cy5, which has an emission spectrum centered on 667 nm, and excluding light from the 635 nm laser and the Cy3 dye. Likewise, the Cy3 ('green') filter transmits 550–600 nm light; the standard 'yellow' emission filter transmits 607–695 nm light; and the standard 'blue' filter transmits 508–560 nm. The spectra for the standard emission filters in the GenePix scanner have very sharp cut-on and cut-off frequencies, so the transmission coefficient is above 0.8 within the entire specified band, and close to zero outside it. The standard emission filters are custom-designed to reject laser light at optical density (OD) 8 for each respective laser excitation, and reject broad-spectrum light at OD 5.

Emission filters are located in a sixteen-position filter wheel that is readily accessible in the GenePix Professional 4200A. Four positions in the filter wheel are reserved for the standard emission filters (red, yellow, green, blue) and are installed at the factory and are not accessible. The remaining twelve positions are available to the user for custom filter installation (see below). Additional filters can be obtained directly from Axon, or a custom filter can be ordered directly through an optical filter company such as Chroma Technology Corp. (Brattleboro, VT); Omega Optical Inc. (Brattleboro, VT); or Semrock (Rochester, NY.) We recommend that you use latex gloves or finger cots when handling filters to avoid getting fingerprints on the optical surfaces.

If you decide to use different dye and/or filter combinations, it is important to match the optical characteristics of filter, dye, and laser. The transmission band of the filter should be centered as close as possible to the emission peak of the chosen dye, and must be narrow enough to exclude light from the laser, other dyes, etc., but not so narrow that too many emission photons are excluded, otherwise sensitivity will be lost. We also strongly recommend that a reject filter of OD 8 be employed to eliminate laser excitation.

Installing Filters

Before installing filters, turn the scanner power off. To gain access to the filter wheel, use the hex key to open the front access panel and gently pull the filter wheel cassette out by the handle.

The filter wheel cassette is numbered and there is an access port on the top of the cassette in which to install filters. There is an opening at the rear of the cassette that allows you to rotate the filter. Choose a filter position by moving the filter wheel so that the designed position is lined up with the top access port. There is an arrow on the mounting ring of the filter that indicates the propagation direction. The arrow on all filters obtained from Axon Instruments indicate the propagation path *to the PMT detector*; thus the arrow should be facing down (into the table) when installing a filter in

the filter wheel. Different manufacturers have different conventions, so if you obtain a filter from a third-party, consult with the manufacturer's documentation in order to determine the correct orientation. Once you determine the orientation of the filter, carefully drop it into position; a small hex key or similar object can be used to gently maneuver the filter into the seated position—do not touch the optical surface as this will irreparably damage the filter, but it is safe to touch the black mounting ring. Next, secure the filter in place with the threaded metallic ring using the filter ring tool. It is advisable to record the filter information and position number on a piece of paper. Install additional filters if desired, and replace filter wheel cassette back into the scanner. The access panel must be securely closed before scanning can be resumed. To uninstall a filter, perform this procedure in reverse order.

The filter information is stored in the PC board of the scanner (not the computer). After each filter installation or change is made, turn on the scanner and open the *Filter Tool* program, located in the Axon Laboratory / GenePix Pro 5 folder of the start menu. Activate the filter position by pressing on the check box, and enter a description that will help you identify the filter. There are two fields available for each position in which you can enter filter information. We recommend that you use the filter supplier's designation as the "Official name" and an easily recognizable description for "Alternative name" such as the dye name (*e.g.*, FITC, Texas Red.)

Photo-Detection

The GenePix scanner uses a high-sensitivity, low-noise photo-multiplier tube (PMT) to detect the emitted fluorescent light. Loosely defined, PMTs are optical components that convert incident photons into electrons via the photoelectric effect. When an incident photon impinges on the active surface of the PMT (the photocathode), an electron is generated. This electron flows through a series of electron multipliers (dynodes) to the

anode. The amount of current that flows from the anode is directly proportional to the amount of incident light at the photocathode.

The amount of amplification that a PMT can produce depends on the number of dynodes in the PMT, and the voltage that is applied to the PMT. With PMTs it is possible to achieve a gain of 10^7 . When you increase the PMT voltage setting in the GenePix 4200A software, you are increasing the gain of the PMTs.

Note: It is important to realize that when you increase the gain of the system, you are also increasing the sensitivity to non-specific fluorescence, thus increasing any electronic noise in the system. This increase in noise can often be overcome by line averaging but the cost is decreased acquisition speed. In general, the signal-to-noise ratio is not improved by increasing the PMT voltage.

The PMT in the GenePix scanner meets two essential criteria for high-quality imaging. First, since different photocathodes have varying sensitivities to specific wavelengths of light (quantum efficiency), the PMT has been carefully chosen for the typical wavelengths of fluorescence emission light produced by microarrays. Second, the PMT has been further selected for its reliability and optimal signal-to-noise performance.

Mechanical

The GenePix scanner is a mechanical device. Its laser-based scanning system uses a patent-pending mechanical design. The GenePix scanner operates by slowly scanning the slide in the *y*-direction, and rapidly scanning in the *x*-direction. The design of the GenePix scanner allows scanning up to four wavelengths of a 25 mm × 75 mm slide at 5- μ m resolution in about 30 minutes.

When using such a mechanical device, one must be able to monitor the variability in the performance of its mechanical components. The GenePix scanner has been designed with a number of optimizations to ensure that variability in the mechanical performance is monitored in real time during acquisition, and corrected. The

utmost attention has been given to selecting and implementing components to create a highly reliable, high-performance instrument.

Placing Slides into Slide Holder

Standard microscope slides (1" × 3"; 25 × 75 mm) can be used in GenePix scanners. Open the door of the scanner and lift up the top of the slide holder by lift on the clasp. Hold the slide in one hand with the thumb and index finger, positioned near the bottom of the slide. Carefully place the slide, feature-side-down into the slide holder and gently click the side spring 2 to 3 times to help position the slide into place, if necessary. Slowly close the top of the slide holder, ensuring that the clasp locks into place and close the scanner door.

Perform a Preview or Data Scan using the controls in the GenePix Pro software. Refer to the GenePix Pro user's manual for more information.

Proper Use of Slide Holder

The slide holder retracts into the scanner when the scanner door is closed. Opening the door or pressing the Eject button in GenePix Pro will return the slide holder to the loading position. It takes a few moments to fully eject the slide holder; **do not open the slide holder before it has come to a complete stop** as this can cause misalignment. This type of misuse is evident when scrapes or even paint from the scanner enclosure are found on the slide holder. Misalignment of the slide holder caused by prematurely opening the slide holder while the stage is still in motion can result in deterioration of scanner performance and may require servicing.

Electronic

While the optical components of the GenePix scanner are indispensable for generating the electrical current from the PMT, the electronic components are just as important. They govern the operation of the PMT and the conversion of their analog current signals into a digital signal that can be processed by the host computer.

Axon Instruments is the world leader in the design and construction of low-noise amplifiers. Our amplifiers are used around the world by researchers in the neuroscience community. We have 20 years of expertise creating low-noise amplifiers and high-performance computer hardware to digitize these signals. Our expertise in this area has been applied to the GenePix scanner. Integrated with the GenePix scanner is not only our high-performance analog-digital converter, but also the experience and knowledge to create low-noise, high-performance amplifiers. We have employed a unique design to ensure that the PMTs in GenePix scanner are linear over a wide dynamic range, and that the signal-to-noise level exceeds all other currently available microarray scanners.

General Items of Interest

Spatial Resolution

The GenePix Professional 4200A is a scanning instrument. This means that a single beam of light is rapidly scanned across the microarray, and a composite pixilated image is created from the digitized signals from the PMT. The spatial resolution of the system refers to the size of the pixel. It determines the minimum distance that can be distinguished between two points of light. In the GenePix scanner, the highest spatial resolution is 5 microns.


It is possible to change the resolution of Data Scans by changing the pixel size setting in the Hardware Settings dialog of the GenePix Pro software. It is important to note that this does not change the size of the scanning beam. The GenePix scanner implements changes in resolution by increasing the *y*-direction step size and averaging the *x*-direction samples over multiples of 5- μm spots. For example, to change the pixel size from 5 μm to 10 μm the instrument only scans every second *x*-direction line (by increasing each *y*-direction step from 5 to 10 μm —consequently a 10- μm scan is twice as fast as a scan at 5- μm) and outputs the average of each adjacent pair of 5- μm spots in the *x*-direction as the value of a single 10- μm ‘pixel’. Similarly, a scan at 20- μm pixel size is implemented by scanning every fourth *x*-direction line (and so is four times faster than at 5 μm) and


using the average of four contiguous 5- μm spots on each x-direction line as the value of a 20- μm ‘pixel’.

Data Scan pixel size can be increased beyond 20 μm in increments of 20 μm , but the scanning pattern (and therefore scan duration) remains the same. Every fourth x-direction line is scanned and the spots within the larger pixel size are averaged. For example, a 60- μm data scan will average the values of thirty-six 5- μm spots (three x-direction lines contribute twelve 5- μm spots each) to calculate the value for each 60- μm ‘pixel’. This is illustrated schematically in the figure below.


Changing the resolution of a Data Scan may often be useful if you have large features and you do not need to scan at maximum resolution. As you decrease the resolution, you will decrease the final image size (in megabytes). However, keep in mind that you need to acquire enough pixels for each feature in order to make accurate measurements.

For a Preview Scan, the same 5- μm spots are scanned. However the y-direction step size is increased to 40 μm , and each 40- μm Preview Scan ‘pixel’ is an average of eight contiguous 5- μm spots in the x-direction. This speeds image acquisition for balancing the PMT gains while minimizing photobleaching.

D) 60- μ m Data Scan

E) Preview Scan (40- μ m)

Representation of how the GenePix scanner samples 5- μ m spots and averages them to calculate the values of the image pixels for the various scan modes and pixel sizes. On the left side of each image pair is the raw data scan from the GenePix Scanner. Each square represents a 5 μ m $\times$ 5 μ m spot as scanned by the 5- μ m laser beam. On the right are the resultant image pixels at various resolutions. The shading indicates the spots that are scanned by the GenePix scanner as well as the corresponding pixels in the final image.

Dynamic Range and Detection Limit

The dynamic range of an instrument is the signal width that can be resolved above the noise the instrument produces. The dynamic range depends on all the components in a given system. For example, if you use a PMT with a dynamic range of 10^6 but process the signals through a sub-standard 8-bit analog-to-digital converter, you will end up with an instrument with a poor dynamic range.

Furthermore, while an instrument may have a large dynamic range, one may not be able to use the whole range reliably if the response is non-linear. Therefore when rating an instrument one must consider the range over which its response remains linear. For the GenePix scanner, the dynamic range is 10^4 .

The dynamic range of an instrument is often considered in conjunction with its detection limit. The detection limit is defined as the dye concentration for which the signal-to-noise ratio = 3. Signal-to-noise ratio is calculated as

$$(\text{Signal} - \text{Background}) / \text{Standard Deviation of Background}$$

Although you may be able to see features below this level, the quantitative accuracy diminishes significantly.

Specifications

Sample type:	Standard microscope slides (1" × 3", 25 × 75 mm, 26 × 76 mm; 0.9–1.2 mm thick)
Maximum scan area:	22 × 73 mm
Maximum resolution:	5 μm
Scan mode and order:	Sequential; fully user-configurable scan order
Detector:	Photomultiplier tube (PMT)
Digital resolution:	16-bit
Image format:	Single-image or multi-image TIFF
Computer interface:	SCSI (adapter included)
Scan time:	Four minutes per channel for full size scan area at 10 μm resolution

Dynamic range:	Four orders of magnitude
Rotational Uniformity:	Less than 4% non-uniformity
Dimensions:	43 cm (16.9") × 65 cm (25.6") × 34 cm (13.4") (W×D×H)
Weight:	45 kg (99 lb)
Power supply:	Line voltage: 85 to 264 VAC (110 to 340 VDC) universal voltage input Line frequency: 50–60 Hz Power: 150 W / Fuse: 2.0 A slow. 5 × 20 mm

Technical Assistance

If you need help to resolve a problem, there are several ways to contact Axon Instruments / Molecular Devices:

World Wide Web

www.axon.com

Phone

(800) 635-5577

Fax

+1 (510) 675-6300

E-mail

axontech@axon.com

Questions?

See Axon's Knowledge Base: <http://support.axon.com>

Warranty and Repair Service

Standard Warranty

Axon Instruments / Molecular Devices, Corp. warrants its non-consumable hardware products to be free from defects in materials and workmanship for 12 months from date of invoice or date of purchased installation visit, which ever is later. The warranty covers the cost of parts and labor to repair the product. Please keep the shipping container for future use. If you should need an additional container, one can be provided for an additional charge. Products returned to our factory for repair should be properly packaged with transportation charges prepaid. Axon Instruments / Molecular Devices will pay for the return shipping of the product to the customer. If the shipment is to a location outside the United States, the customer will be responsible for paying all duties, taxes and freight clearance charges if applicable.

The warranty is valid when the product is used for its intended purpose and does not cover products which have been modified without approval from Axon Instruments / Molecular Devices, or which have been damaged by abuse, accident or connection to incompatible equipment.

This warranty is in lieu of all other warranties, expressed or implied.

Repair Service

The company reserves the right to cease providing repair maintenance, parts and technical support for its non-consumable hardware products five years after a product is discontinued. Technical support for old versions of software products will cease 12 months after they are upgraded or discontinued.

If you purchased your instrument from a Distributor or OEM Supplier, contact them for repair service.

If you purchased your instrument from Axon Instruments/ Molecular Devices, contact our Technical Support Department. If it is determined your instrument must return to the factory for repair, the Technical Support Representative will issue a Service Request (SR) number. Our Logistic Coordinator will contact you with specific instructions.

Out-of-Warranty Repair Service

Out-of-warranty repair service is available. Contact the supplier for current rates or see our web site at: www.axon.com.

Optional Service Agreement

Purchasing an optional Service Agreement extends the coverage of the Standard Warranty. Contact the supplier for current rates.

Shipping

The GenePix 4200A is a solidly built instrument designed to survive shipping around the world. However, in order to avoid damage during shipping, the GenePix 4200A must be properly packaged.

In general, the best way to package the GenePix 4200A is in the original factory carton. If this is no longer available, we recommend that you carefully wrap the

GenePix 4200A in at least three inches (75 mm) of foam or "bubble-pack" sheeting. The wrapped instrument should then be placed in a sturdy cardboard carton. Mark the outside of the box with the word FRAGILE and an arrow showing which way is up.

We do NOT recommend using loose foam pellets to protect the GenePix 4200A. During shipping, there is a good chance that the instrument will shift within the loose pellet packing and be damaged.

If you need to ship the GenePix 4200A to another location, or back to the factory, and you do not have a means to adequately package it, Axon Instruments can ship the proper packaging material to you for a small fee. This may seem an expense you would like to avoid, but it is inexpensive compared to the cost of repairing an instrument that has sustained shipping damage.

It is your responsibility to package the instrument properly before shipping. If the packaging is inadequate, and the instrument is damaged during shipping, the shipper will not honor your claim for compensation.

Licensing Notice

Axon Instruments / Molecular Devices Corp. is not licensed under any patents owned by Oxford Gene Technology Limited ("OGT"), covering oligonucleotide arrays and methods of using them to analyze polynucleotides. The purchase of Axon Instruments / Molecular Devices products does not convey any license under any of OGT's patent rights, including any right to make or use oligonucleotide arrays under OGT's patents.

Customers may use Axon Instrument / Molecular Devices products to analyze oligonucleotide arrays according to OGT's patented methods if those arrays have either been purchased from OGT's licensed suppliers, or have been made by the customer under a license from OGT.

Please contact OGT to enquire about a license under OGT's patents at licensing@ogt.co.uk <<mailto:licensing@ogt.co.uk>>.

USE OF THIS INSTRUMENT WITH MICROARRAYS MAY REQUIRE A LICENSE FROM ONE OR MORE THIRD PARTIES THAT HAVE PATENTS RELEVANT TO SUCH USE. AXON INSTRUMENTS / MOLECULAR DEVICES DOES NOT SUGGEST OR PROMOTE THE USE OF THIS INSTRUMENT IN A MANNER THAT INFRINGES ON THE PATENT RIGHTS OF A THIRD PARTY. YOU ARE ENCOURAGED TO EVALUATE WHETHER A LICENSE IS REQUIRED FOR YOUR SPECIFIC APPLICATION OF THIS INSTRUMENT. COMPANIES THAT HAVE INTELLECTUAL PROPERTY RIGHTS IN THE POTENTIAL FIELD OF APPLICATION OF THIS INSTRUMENT INCLUDE WITHOUT LIMITATION, AFFYMETRIX, INC.

(“AFFYMETRIX”), AGILENT, AND OXFORD GENE TECHNOLOGY. THIS INSTRUMENT HAS NOT BEEN LICENSED OR APPROVED FOR DIAGNOSTIC APPLICATIONS.

THE USE OF THIS INSTRUMENT IN CONNECTION WITH MICROARRAYS MAY BE WITHIN THE SCOPE OF PATENTS HELD BY AFFYMETRIX. TO THE EXTENT THAT AFFYMETRIX PATENT RIGHTS ENCOMPASS THIS INSTRUMENT OR ITS USE, AFFYMETRIX HAS GRANTED A LIMITED PATENT LICENSE FOR RESEARCH USE ONLY AND NOT FOR USE IN DIAGNOSTIC PROCEDURES. SUCH LICENSE, IF APPLICABLE, IS LIMITED TO USE OF THIS INSTRUMENT WITH SPOTTED MICROARRAYS SEPARATELY LICENSED BY AFFYMETRIX. NO LICENSE IS CONVEYED, BY IMPLICATION, ESTOPPEL OR OTHERWISE, TO USE THIS INSTRUMENT WITH MICROARRAYS MADE USING *IN SITU* OR PHOTOLITHOGRAPHIC SYNTHESIS. NO OTHER LICENSE IS CONVEYED, BY IMPLICATION, ESTOPPEL OR OTHERWISE, UNDER ANY AFFYMETRIX PATENT OR OTHER INTELLECTUAL PROPERTY RIGHT.

This instrument is licensed by Affymetrix under the following patents: U.S. Patent Nos. 5,578,832; 5,631,734; 5,834,758; 5,936,324; 5,981,956; 6,025,601; 6,141,096; 6,171,793; 6,185,030; 6,201,639; 6,207,960; 6,218,803; 6,225,625; 6,252,236; 6,262,838; 6,335,824; 6,403,320; 6,403,957; 6,407,858; 6,472,671; 6,490,533; 6,545,264; 6,597,000; 6,643,015; and 6,650,411.

Declaration of Conformity

Manufacturer: Molecular Devices Corp. / Axon Instruments
3280 Whipple Road
Union City, CA 94587 USA

Type of Equipment: Microarray Scanner

Model Number: GenePix Professional 4200A

Serial Number: _____

Year of Manufacture: 2002

Application of Council Directives:
EC EMC Directive 89/336/EEC as amended
EC Low Voltage Directive 73/23/EEC as amended

Harmonized Standards to which Conformity is Declared:

EMC:

EN 61326 04: 1997 (Emissions)
EN 61326 04: 1997 (Immunity)

Low Voltage:

EN 61010-1:1993 +A1 +A2:1995
EN 60825-1: 1994 +A11:1996 (Laser Safety)

I, the undersigned, hereby declare that the equipment specified above conforms to the above Directives and Standards.

Authorized Signature and Date: (Signature on file)

Important Safety Information

The operator of the GenePix Professional 4200A Microarray Scanner instrument is assumed to be trained in the correct operation of the instrument and the safety issues. Throughout the *GenePix Professional 4200A Microarray Scanner User's Guide*, the word "you" refers to this trained operator.

Safety Text Used In This Guide

Make sure you follow the precautionary statements presented in this guide.

Laser Safety Statement: The Center for Devices and Radiological Health (CDRH) of the U.S. Food and Drug Administration has implemented regulations for laser products manufactured since August 1, 1976. Compliance is mandatory for products marketed in the United States. The scanner is certified as a "Class 1" laser product under the U.S. Department of Health and Human Services (DHHS) Radiation Performance Standard according to the Radiation Control for Health and Safety Act of 1968. Since radiation emitted inside the scanner is completely confined within protective housings and external covers, the laser beam cannot escape during any phase of normal user operation.

Warning: Using controls, making adjustments, or performing procedures other than those specified in the GenePix 4200A User's Guide may result in hazardous exposure to laser radiation, high voltage, or moving parts. Exposure to these hazards can cause severe or fatal injury.

Caution: Indicates that damage to the instrument, loss of data, or invalid data could occur if the user fails to comply with the advice given.

Important: Highlights information that is critical for optimal performance of the system.

Note: Identifies items of general interest.

The Protective Housing

The protective housing and instrument interlocks are designed to protect you from exposure to the laser radiation, high voltage, or moving parts. You do not need to access the interior of the instrument. There are no serviceable parts inside.

Warning: Do not defeat the interlock, open the protective housing, or try to gain access to the interior of the instrument through any other openings. Opening the housing can damage the instrument components and result in hazardous exposure to laser radiation, high voltage, or moving parts.

The GenePix 4200A is a Class I laser device. There are high power lasers embedded inside it, but the user cannot access these. The moving Mirror and Lens Assembly located under the Slide Holder directs and focuses these Class 3B laser beams directly up, onto the surface of the sample slide. Slides can be placed safely in the laser's scanning area because the slide-loading door has redundant interlocks that block the beams. There are two panels on the device that are not interlocked, but these are fastened with screws and require a tool to remove them. They are shown in the Non-Interlocked Panels section below. Do not remove any screws to open the housing!

38 • Important Safety Information

- The lasers inside can cause eye damage before you have a chance to blink!
- There are high voltages inside.
- There are no user serviceable parts inside.
- The unit will be damaged if you attempt to open the protective housing.

Interlock Failure Symptoms

If you experience any of the following symptoms, you may have an interlock failure. It is unsafe to continue to use the GenePix 4200A with only one interlock functioning. Please contact Technical Support immediately at either of the numbers indicated on the rear panel of the GenePix 4200A.

1. Mirror and lens assembly keeps moving below the slide after the door is opened.
2. The blue “Scanning” LED at the front of the GenePix 4200A doesn’t extinguish when the door is open.
3. You cannot scan any slides.
4. You can no longer hear the distinctive metal-on-metal sound that the mechanical interlock makes when it falls into place. This normally happens when the slide loading door has opened about ¼”.

Lifting Hazard

Warning: The GenePix Professional 4200 scanner weighs approximately 99 lbs. (45 kg) and requires two persons to lift. There are four molded handholds on the underside of the scanner (see figure). Use appropriate lifting technique to avoid injury.

Lasers

Cautions

- Do not stick anything into the openings or vents of the GenePix 4200A, other than loading the slide. You may cause an interlock to fail.
- Do not attempt to move any parts below the slideholder. You may cause an interlock to fail.

Embedded Laser Classification and Power

Wavelength	Power	Divergence	Duration	Embedded Laser's Class
488 nm	20 mW	<1.2 mrad	Continuous	Class IIIb
532 nm	14±4 mW	<1.2 mrad	Continuous	Class IIIb
594 nm	1-2 mW	0.92 mrad	Continuous	Class IIIb
635 nm	15 mW	0.7 mrad	Continuous	Class IIIb

Warning: Exposure to laser radiation can cause injury. For example, viewing the laser radiation directly can cause blindness.

High-Voltage Hazard

High-voltage electronics can be found in the GenePix 4200A Microarray Scanner instrument. Under normal operating conditions, you are protected from high voltage within the instrument. In the unlikely event the scanner receives a surge from the AC power line greater than 1 kV, the scanner may reset. In order to continue normal operation simply cycle power to both the computer and scanner.

Warning: Do not try to gain access to the interior of the instrument. Exposure to high-voltage electronics can cause severe or fatal injury.

High-voltage electronics can also be found in the scanner's controlling computer. See the computer manufacturer documentation on the high-voltage hazard warnings. Make sure you follow the instructions on the safe operation of the computer.

Moving Parts

The GenePix 4200A Microarray Scanner instrument contains moving parts that can cause injury. Under normal operating conditions, the instrument is designed to protect you from the moving parts. The interlock and protective housing are designed so that you cannot access the moving parts during a scan.

Warning: Do not try to gain access to the interior of the instrument. The moving parts inside the instrument can cause injury.

Fuses

There are two 2A/250V fuses located adjacent to the power switch at the rear of the unit.

Hazardous Material Precautions

Use standard laboratory procedures and cautions when working with chemicals.

Caution: Always follow the manufacturer's precautions when working with chemicals. Axon Instruments is not responsible or liable for any damages caused by or as a consequence of the use of any hazardous material.

Non-Interlocked Panels

View of the top and rear of GenePix Professional 4200A scanner showing the twenty screws that attach the protective housing.

DO NOT REMOVE THESE SCREWS!

DANGER
LASER RADIATION WHEN OPEN
AVOID DIRECT EXPOSURE TO BEAM

Don't remove these screws.

Maintenance and Service

No user maintenance is required. No user service is allowed. If there is a problem, please contact Technical Support.

List of Controls

This Guide constitutes a list of controls.

Caution: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

Safety Labels

If a label becomes illegible or is missing for any reason, please contact Microarray System Technical Support for a free replacement label. While waiting for the replacement label, copy the label shown in this section and attach the copy of the label to the instrument.

Symbol Explanations

		
Device standby. Ready to scan slide. Stage waiting at the start position.	Scanning.	Eject. Slideholder ready to be unloaded.

GenePix™ PROFESSIONAL 4200A
Axon Instruments

PATENTS PENDING

CLASS 1 LASER PRODUCT

This device complies with 21 CFR 1040.10 & 1040.11

LIFTING HAZARD, HEAVY OBJECT.

FUSE  **250V, T2A**

POWER 100 - 240V~	FREQUENCY 50 - 60Hz
------------------------------------	--------------------------------------

MAX OPERATING CURRENT: 1.25 AMP

CE 

MADE IN U.S.A.
Axon Instruments, Inc.
3280 Whipple Road • Union City
California 94587 • U.S.A.
Tel. +1 (510) 675-6200
Fax +1 (510) 675-6300

CAUTION
Visible laser radiation when case is open and interlocks defeated.
CLASS 1 LASER PRODUCT

DANGER
HAZARDOUS VOLTAGE.
To prevent electric shock and damage to the instrument, do not attempt to open unit.

WARNING
NO CUSTOMER SERVICEABLE COMPONENTS ARE CONTAINED WITHIN THIS UNIT. CONSULT THE MANUFACTURER FOR REPAIR/RETURN INSTRUCTIONS.

THE USER SHALL BE MADE AWARE THAT, IF THE EQUIPMENT IS USED IN ANY MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

Site Requirements

The GenePix Professional 4200A is designed to operate optimally if the site requirements are met. As with any precision optical instrument, care should be taken to maintain a low-dust, low-vibration environment. Temperature and humidity extremes may compromise performance.

Environmental Temperature: 50–86° F (10–30° C).

Environmental Humidity: 5–95% non-condensing.

Power Requirements: Universal Input 50-60 Hz, 100–240V~AC, 2A max.

Space Requirements: 25.6 in deep × 16.9 in wide × 13.4 in high (65 cm × 43 cm × 34 cm); plus additional 6 in (15.2 cm) clearance.

Weight Requirements: Sufficient to support 99 lbs (45 kg) with minimal vibration.

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**LIFTING HAZARD,
HEAVY OBJECT.**

**HAZARDOUS VOLTAGE.**
To prevent electric shock and damage to the instrument, do not attempt to open unit.

**WARNING**
NO CUSTOMER SERVICEABLE COMPONENTS ARE CONTAINED WITHIN THIS UNIT. CONSULT THE MANUFACTURER FOR REPAIR/RETURN INSTRUCTIONS.
THE USER SHALL BE MADE AWARE THAT, IF THE EQUIPMENT IS USED IN ANY MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

**CAUTION**
Visible laser radiation when case is open and interlocks defeated.
CLASS 1 LASER PRODUCT

**DANGER**
HAZARDOUS VOLTAGE.
To prevent electric shock and damage to the instrument, do not attempt to open unit.

**WARNING**
NO CUSTOMER SERVICEABLE COMPONENTS ARE CONTAINED WITHIN THIS UNIT. CONSULT THE MANUFACTURER FOR REPAIR/RETURN INSTRUCTIONS.
THE USER SHALL BE MADE AWARE THAT, IF THE EQUIPMENT IS USED IN ANY MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

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PATENTS PENDING
CLASS 1 LASER PRODUCT

This device complies with 21 CFR 1040.10 & 1040.11

LIFTING HAZARD, HEAVY OBJECT.

FUSE 250V, T2A
POWER 100 - 240V~
FREQUENCY 50 - 60Hz
MAX OPERATING CURRENT: 1.25 AMP

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WARNING
NO CUSTOMER SERVICEABLE COMPONENTS ARE CONTAINED WITHIN THIS UNIT. CONSULT THE MANUFACTURER FOR REPAIR/RETURN INSTRUCTIONS.
THE USER SHALL BE MADE AWARE THAT, IF THE EQUIPMENT IS USED IN ANY MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

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Index

- Alexa Fluor Dyes, 16
- Amplifiers, 22
- Custom Filters, 18
- Cy3, Cy5, 16
- Digitization, 1
- Emission Filter, 17, 18
- Filters
 - Installing, 10
- Fluorophore, 17. See Cy3, Cy5
- GenePix
 - Installation, 3
- Hardware, 15
 - Amplifiers, 22
 - Laser Excitation, 16
 - Quantum Efficiency, 20
- Help
 - On-Line, 3
- Hardware
 - Spatial Resolution, 22
- Lasers
 - Temperature Fluctuations, 16
- On-line Help, 3
- Photo-Multiplier Tube
 - Gain, 20
- Pixel Size, 22, 23, 24
- Preview Scan, 23, 24
- Spatial Resolution, 22, 23
- Stokes Shift, 17
- Temperature Fluctuations, 16