
Criterion™ Precast Gels

Instruction Manual and Application Guide



Bio-Rad Technical Support

For help and technical advice, please contact the Bio-Rad Technical Support department. In the United States, the Technical Support department is open Monday–Friday, 5:00 AM–5:00 PM, Pacific Time.

Phone: 1-800-424-6723

Fax: 1-510-741-5802

Email: LSG_TechServ_US@bio-rad.com (for U.S. and international customers)

Online technical support and worldwide contact information are available at www.consult.bio-rad.com.

Legal Notices

No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording, or any information storage or retrieval system, without permission in writing from Bio-Rad Laboratories.

Bio-Rad reserves the right to modify its products and services at any time. This user guide is subject to change without notice.

Although prepared to ensure accuracy, Bio-Rad assumes no liability for errors, or for any damages resulting from the application or use of this information.

Brij and Tween are trademarks of ICI Americas, Inc. Coomassie is a trademark of BASF Aktiengesellschaft. Ficoll is a trademark of Amersham Pharmacia Biotech. StrepTactin is a trademark of Institut für Bioanalytik GmbH. StrepTactin is covered by German patent application P 19641876.3. Bio-Rad Laboratories, Inc. is licensed by Institut für Bioanalytik GmbH to sell these products for research use only. SYBR is a trademark of Invitrogen Corporation. SYPRO is a trademark of Molecular Probes, Inc. Bio-Rad is licensed to sell SYPRO products for research use only, under U.S. Patent 5,616,502. Triton is a trademark of Union Carbide.

Criterion™ TGX Stain-Free™ and Criterion Stain Free™ precast gels are covered by U.S. Patent No. 7,569,130.

Purchase of Criterion™ XT Bis-Tris gels, XT MOPS running buffer, XT MES running buffer, XT MOPS buffer kit, and XT MES buffer kit is accompanied by a limited license under U.S. patents 6,143,154; 6,096,182; 6,059,948; 5,578,180; 5,922,185; 6,162,338; and 6,783,651 and corresponding foreign patents.

Contents

| | |
|---|----------|
| Chapter 1: Criterion™ Precast Gels | 1 |
| 1.1 Introduction | 1 |
| 1.2 Gel Formulations | 2 |
| 1.3 Comb Configurations | 2 |
| 1.4 Specifications | 2 |
| 1.5 Storage Conditions | 3 |
| 1.6 Important Notes | 3 |
| | |
| Chapter 2: Setup and Basic Operation | 4 |
| 2.1 Workflow Overview | 4 |
| 2.2 Required Materials | 5 |
| 2.3 Setting Up and Running Criterion Gels in the Criterion Cell | 5 |
| 2.4 Removing the Gel | 5 |
| | |
| Chapter 3: SDS-PAGE | 7 |
| 3.1 Introduction | 7 |
| 3.2 Criterion Gel Selection Guide for SDS-PAGE | 7 |
| 3.2.1 Criterion™ TGX™ and Criterion™ TGX Stain-Free™ Gels | 8 |
| 3.2.2 Criterion Tris-HCl and Criterion Stain Free™ Gels | 9 |
| 3.2.3 Criterion™ XT Bis-Tris Gels | 10 |
| 3.2.4 Criterion XT Tris-Acetate Gels | 11 |
| 3.3 SDS-PAGE Buffers | 12 |
| 3.3.1 Running Buffers | 12 |
| 3.3.2 Sample Buffers | 12 |
| 3.4 Sample Preparation | 12 |
| 3.5 Running Conditions | 13 |

| | |
|--|-----------|
| Chapter 4: Native PAGE | 14 |
| 4.1 Introduction | 14 |
| 4.2 Criterion Gel Selection Guide for Native PAGE | 14 |
| 4.2.1 Criterion TGX and Criterion TGX Stain-Free Gels | 14 |
| 4.2.2 Criterion Tris-HCl and Criterion Stain Free Gels | 15 |
| 4.2.3 Criterion XT Tris-Acetate Gels | 15 |
| 4.3 Native PAGE Buffers | 15 |
| 4.4 Sample Preparation | 16 |
| 4.5 Running Conditions | 16 |
| | |
| Chapter 5: Criterion Stain Free System | 17 |
| 5.1 Introduction | 17 |
| 5.2 Criterion Stain Free Workflow | 18 |
| 5.3 Electrophoresis with Criterion TGX Stain-Free Gels | 18 |
| 5.4 Using the Gel Doc™ EZ Imager | 18 |
| | |
| Chapter 6: Peptide Analysis | 19 |
| 6.1 Introduction | 19 |
| 6.2 Criterion Tris-Tricine/Peptide Gels | 19 |
| 6.2.1 Gel Composition | 19 |
| 6.2.2 Gel Selection Guide | 19 |
| 6.3 Tris-Tricine/Peptide Buffers | 19 |
| 6.4 Sample Preparation | 20 |
| 6.5 Running Conditions | 20 |
| | |
| Chapter 7: Isoelectric Focusing (IEF) | 21 |
| 7.1 Introduction | 21 |
| 7.2 Criterion IEF Gels | 21 |
| 7.2.1 Gel Composition | 21 |
| 7.2.2 Gel Selection Guide | 21 |
| 7.3 IEF Buffers | 21 |
| 7.4 Sample Preparation | 22 |
| 7.5 Running Conditions | 22 |

| | |
|--|-----------|
| Chapter 8: Protease Analysis by Zymogram PAGE | 23 |
| 8.1 Introduction. | 23 |
| 8.2 Criterion Zymogram Gels. | 23 |
| 8.2.1 Gel Composition | 23 |
| 8.2.2 Gel Selection Guide. | 23 |
| 8.3 Zymogram Buffers. | 23 |
| 8.4 Sample Preparation. | 24 |
| 8.5 Running Conditions. | 24 |
| | |
| Chapter 9: Nondenaturing Nucleic Acid PAGE | 25 |
| 9.1 Introduction. | 25 |
| 9.2 Criterion TBE Gels. | 25 |
| 9.2.1 Gel Composition | 25 |
| 9.2.2 Gel Selection Guide. | 25 |
| 9.3 Nondenaturing Nucleic Acid PAGE Buffers | 25 |
| 9.4 Sample Preparation. | 26 |
| 9.5 Running Conditions. | 26 |
| | |
| Chapter 10: Denaturing Nucleic Acid PAGE | 27 |
| 10.1 Introduction. | 27 |
| 10.2 Criterion TBE-Urea Gels | 27 |
| 10.2.1 Gel Composition | 27 |
| 10.2.2 Gel Selection Guide. | 27 |
| 10.3 TBE-Urea Buffers | 27 |
| 10.4 Sample Preparation. | 27 |
| 10.5 Running Conditions. | 28 |
| | |
| Chapter 11: 2-D Electrophoresis | 29 |
| 11.1 Introduction. | 29 |
| 11.2 Equilibration. | 29 |
| 11.3 Agarose Overlay | 29 |
| 11.4 Second-Dimension Electrophoresis. | 29 |

| | |
|--|-----------|
| Chapter 12: Detection | 30 |
| 12.1 SDS-PAGE and Native PAGE Detection | 30 |
| 12.2 Peptide Gel Staining | 31 |
| 12.3 IEF Gel Staining | 31 |
| 12.4 Zymogram Gel Staining | 32 |
| 12.5 TBE Gel Staining | 32 |
| 12.6 TBE-Urea Gel Staining | 32 |
| | |
| Chapter 13: Blotting | 33 |
| 13.1 Introduction | 33 |
| 13.2 Transfer | 33 |
| 13.2.1 Transfer Buffers | 33 |
| 13.2.2 Wet Transfer Using the Criterion Blotter | 33 |
| 13.2.3 Transfer Using the Trans-Blot® Turbo™ System | 34 |
| 13.2.4 Semi-Dry Transfer Using the Trans-Blot® SD Cell | 36 |
| 13.3 Total Protein Blot Stains | 36 |
| 13.4 Immunodetection | 37 |
| | |
| Chapter 14: Troubleshooting | 38 |
| | |
| Appendix A: Quick Start Guide | 40 |
| Appendix B: Buffers | 42 |
| Appendix C: Related Literature | 45 |
| Appendix D: Ordering Information | 46 |

1

Criterion™ Precast Gels

1.1 Introduction

Criterion precast gels are an effective system for performing polyacrylamide gel electrophoresis (PAGE). These 13.3 x 8.7 cm gels are wider and longer than traditional mini format gels, and their innovative, easy-to-use design produces excellent resolution while accommodating more samples per gel. Designed for use with the Criterion family of vertical electrophoresis cells, which includes the Criterion (2-gel capacity) and Criterion™ Dodeca™ (12-gel capacity) cells, Criterion precast gels allow separation of more samples than mini format gels and provide significant cost and time savings. Some of the unique features provided are:

- Integrated buffer chamber that eliminates buffer leaks
- Capacity for up to 26 samples per gel
- Compatibility with multichannel pipets (12+2 and 26-well)
- Outlined and numbered wells that simplify sample loading and identification
- Patented¹ J-foot design that eliminates post-run gel processing steps and improves gel drying and blotting results
- Criterion™ TGX Stain-Free™ formulations for rapid gel imaging without staining



¹ U.S. patent 6,093,301.

1.2 Gel Formulations

Criterion precast gels are available in a range of formulations for virtually every electrophoresis application (Table 1.1). All Criterion gels are composed of polyacrylamide with a bisacrylamide crosslinker, and they are available in a selection of single percentages and gradients.

Table 1.1. Criterion precast gel formulations.

| Application | Gel Formulation | Sample Buffer | Running Buffer |
|----------------------------|---------------------------|---------------|----------------------------|
| SDS-PAGE | Criterion Tris-HCl | Laemmli | Tris/glycine/SDS |
| | Criterion Stain Free | | |
| | Criterion TGX™ | | |
| | Criterion TGX Stain-Free | | |
| | Criterion™ XT Bis-Tris | | |
| Native PAGE | Criterion XT Tris-acetate | XT | XT MOPS or XT MES |
| | Criterion Tris-HCl | XT | XT Tricine |
| | Criterion Stain Free | Native | Tris/glycine |
| | Criterion TGX | | |
| | Criterion TGX Stain-Free | | |
| Criterion XT Tris-acetate | | | |
| Peptide analysis | Criterion Tris-Tricine | Tricine | Tris/Tricine/SDS |
| Isoelectric focusing (IEF) | Criterion IEF | IEF | Anode and cathode buffers |
| Protease detection | Criterion zymogram | Zymogram | Tris/glycine/SDS |
| dsDNA separation | Criterion TBE | Nucleic acid | Tris/boric acid/EDTA (TBE) |
| ssDNA and RNA separation | Criterion TBE-urea | TBE-urea | TBE |

1.3 Comb Configurations

| Comb Type | Well Volume |
|------------------------|--|
| 12+2 well ¹ | 45 µl with two 15 µl reference wells |
| 18-well | 30 µl |
| 26-well ¹ | 15 µl |
| Prep+2 well | 800 µl with two 15 µl reference wells |
| IPG+1 well | 11 cm ReadyStrip™ IPG strip (450 µl) with one 15 µl reference well |

1.4 Specifications

| | |
|-----------------------------|-------------------|
| Gel material | Polyacrylamide |
| Gel dimensions (W x L) | 13.3 x 8.7 cm |
| Gel thickness | 1.0 mm |
| Resolving gel height | 6.5 cm |
| Cassette dimensions (W x L) | 15.0 x 10.6 cm |
| Cassette material | Styrene copolymer |
| Comb material | Polycarbonate |
| Running buffer | 460 ml (per gel) |

¹ Multichannel pipet compatible.

1.5 Storage Conditions

Table 1.2. Storage conditions for Criterion precast gels. Store gels flat. Shelf life is from date of manufacture; expiration dates are printed on the cassettes.

| Storage Temperature | Gel Formulation | Shelf Life |
|---------------------|---------------------------|------------|
| Ambient | Criterion XT Bis-Tris | 12 months |
| 2–8°C | Criterion TGX | 12 months |
| | Criterion TGX Stain-Free | 12 months |
| | Criterion Tris-HCl | 12 weeks |
| | Criterion Stain-Free | 12 weeks |
| | Criterion XT Tris-acetate | 8 months |
| | Criterion Tris-Tricine | 12 weeks |
| | Criterion IEF | 24 weeks |
| | Criterion zymogram | 8 weeks |
| | Criterion TBE | 12 weeks |
| | Criterion TBE-urea | 8 weeks |

1.6 Important Notes

Use each Criterion precast gel as soon as possible after removing it from the storage pouch.

Improper storage of Criterion precast gels can produce numerous artifacts.

- Store gels flat
- Avoid prolonged storage at temperatures above those recommended
- Do not freeze gels
- If you suspect your gels have been stored improperly, discard them


Do not run more than one gel type in the same apparatus at the same time. Different gel percentages and formulations have different conductivities and different run times.

Use unstained standards with Criterion TGX Stain-Free and Criterion Stain Free gels, as some prestained standards are not detected by the Gel Doc™ EZ imager. To monitor electrophoresis, use 10 µl of a 1:1 mixture of Precision Plus Protein™ unstained and prestained standards.

2

Setup and Basic Operation

2.1 Workflow Overview



2.2 Required Materials

- Criterion™ precast gels
- Criterion or Criterion™ Dodeca™ cell
- PowerPac™ Basic or PowerPac HC power supply (or equivalent)
- Sample buffer
- Running buffer (460 ml per gel)

2.3 Setting Up and Running Criterion Gels in the Criterion Cell

1. Each Criterion gel is packaged in a plastic storage tray. Remove the cover of the tray by lifting the corner tab and pulling it diagonally across the package. Remove the gel from the package.
2. Remove the comb and rinse the wells with deionized water (diH_2O) or running buffer.
3. Remove the tape from the bottom of the cassette by pulling the tab across the gel.
4. Insert the cassette into one of the slots in the Criterion cell tank so that the integrated upper buffer chamber faces the center of the cell (**A**).


5. Fill each integrated upper buffer chamber with 60 ml running buffer. **A**

6. Fill each half of the lower buffer tank with 400 ml running buffer to the marked fill line.

7. Load samples using a syringe or a pipet with gel-loading tips.


Optional: place a sample loading guide on the outer edge of the cassette to help align pipet tips with the wells (this is particularly useful when using multichannel pipets).

8. Place the lid on the tank, aligning the color-coded banana plugs with corresponding jacks on the lid. See Chapters 3–10 for power supply settings.




2.4 Removing the Gel

1. After electrophoresis is complete, turn off the power supply and disconnect the electrical leads. **B**
2. Remove the lid from the tank and remove the gel(s) from the cell. Pour off and discard the upper running buffer.
3. Invert the cassette and place the integral buffer chamber over the cassette-opening tool built into the Criterion cell lid (**B**).



■ Criterion Precast Gels

4. Press down firmly to break the seals on both sides of the cassette. The cassette splits open approximately $\frac{1}{3}$ of the way. Alternatively, open the gel cassette by sliding the tapered back of the comb into the slits on either side of the cassette.
5. Pull the two halves of the cassette apart from the top to completely expose the gel (C).
6. Remove the gel by either floating the gel into a fixing or staining solution or by carefully lifting the bottom edge of the gel from the cassette.



3

SDS-PAGE

3.1 Introduction

Criterion™ precast gels provide versatile systems for the separation of proteins by either molecular weight (SDS-PAGE) or mass-to-charge ratio (native PAGE). (See Chapter 4 for native PAGE applications and protocols.) This versatility is possible because Criterion gels are made without SDS, allowing the sample buffer and running buffer to determine the separation mechanism.

SDS-PAGE relies on a discontinuous buffer system. Two ions differing in electrophoretic mobility form a moving boundary when voltage is applied. Proteins have an intermediate mobility that causes them to concentrate, or stack, into a narrow zone at the beginning of electrophoresis. As that zone moves through the gel, the sieving effect of the polyacrylamide gel matrix causes proteins of different molecular weights to move at different rates. This stacking effect is responsible for the high resolving power of SDS-PAGE: the sample is loaded in a relatively broad zone, and the moving boundary concentrates the proteins into sharp bands prior to separation.

Protein samples for SDS-PAGE are prepared using SDS and usually a thiol reducing agent such as β -mercaptoethanol or dithiothreitol (DTT). SDS forms complexes with proteins, giving them a rodlike shape and similar mass-to-charge ratio. The reducing agent disrupts disulfide bonds between and within proteins, allowing complete denaturation and dissociation. Heat treatment in the presence of SDS and reducing agent effectively eliminates the effects of native charge and higher order structure on electrophoretic mobility, so the migration distance depends primarily on molecular weight.

Molecular weight is estimated by plotting the logarithm of protein molecular weight vs. the relative mobility (R_f) of the protein (R_f = distance migrated by the protein/distance migrated by the dye front) or by using the point-to-point semilog interpolation method in Quantity One® or Image Lab™ software. Refer to bulletins 3133 and 3144 for more information.

3.2 Criterion Gel Selection Guide for SDS-PAGE

A number of Criterion gel types are available for SDS-PAGE (Table 3.1) in both single and gradient polyacrylamide percentages. Use the protein migration charts and tables to select the gel type that offers optimum resolution of your sample:

- Use single-percentage gels to separate bands of similar molecular weight. Optimum separation occurs in the lower half of the gel, so use a percentage in which the protein migrates to the lower half of the gel
- Use gradient gels to separate samples containing a broad range of molecular weights. Gradient gels allow resolution of both high and low molecular weight bands on the same gel. Larger pore sizes at the top of the gel permit resolution of larger molecules, and smaller pore sizes toward the bottom of the gel restrict excessive separation of small molecules

Table 3.1. Criterion precast gels for SDS-PAGE.

| Gel Formulation | Gels | Description |
|--------------------|----------------------------|--|
| TGX (Laemmli-like) | Criterion™ TGX™ | Laemmli-like, extended shelf life gels |
| | Criterion™ TGX Stain-Free™ | Laemmli-like, extended shelf life gels with trihalo compounds for rapid fluorescence detection |
| Tris-HCl (Laemmli) | Criterion Tris-HCl | Tris-HCl Laemmli gels |
| | Criterion Stain Free™ | Tris-HCl Laemmli gels with trihalo compounds for rapid fluorescence detection |
| Bis-Tris | Criterion™ XT Bis-Tris | Based on a Bis-Tris HCl buffer system (pH 6.4); use these gels with Criterion XT MES buffer for optimum resolution of small proteins |
| Bis-Tris | Criterion XT Bis-Tris | Based on a Bis-Tris HCl buffer system (pH 6.4); use these gels with Criterion XT MOPS buffer for optimum resolution of midsized proteins |
| Tris-acetate | Criterion XT Tris-acetate | Based on a Tris-acetate buffer system (pH 7.0) |

3.2.1 Criterion TGX and Criterion TGX Stain-Free Gels

Criterion TGX (Tris-Glycine eXtended shelf life) gels are Laemmli-like gels with a proprietary modification that extends their shelf life to 12 months and enhances separation characteristics relative to conventional gel types. The TGX formulation yields Laemmli-like separation patterns with short run times and exceptionally straight lanes and sharp bands. TGX gels offer excellent staining quality, greater transfer efficiency, and molecular weight estimation without the need for special, expensive buffers.

These gels are run using standard Laemmli sample buffer and Tris/glycine/SDS running buffer. Two types of TGX formulations are available:

- Criterion TGX — Laemmli-like, extended shelf life gels
- Criterion TGX Stain-Free — Laemmli-like, extended shelf life gels with trihalo compounds that allow rapid fluorescent detection of proteins with the Criterion Stain Free system, eliminating staining and destaining steps for faster results (see Chapter 5 for more details)

Gel Composition

| | |
|--------------|--|
| Crosslinker | 2.6% C |
| Stacking gel | 4% T, 2.6% C |
| Shelf life | ~12 months at 2–8°C; expiration date is printed on each cassette |

| Gel Percentage | Optimum Separation Range |
|----------------------|--------------------------|
| 7.5% | 40–200 kD |
| 10% | 30–150 kD |
| 12% | 20–120 kD |
| 18% | 10–50 kD |
| 4–15% | 20–250 kD |
| 4–20% | 10–200 kD |
| 8–16% | 20–120 kD |
| 10–20% | 10–100 kD |
| Any kD™ ¹ | 10–200 kD |

¹Any kD is a unique single-percentage formulation that provides a broad separation range and short running time.

3.2.2 Criterion Tris-HCl and Criterion Stain Free Gels


These Tris-HCl, Laemmli gels use discontinuous glycinate and chloride ion fronts to form moving boundaries to stack and then separate denatured proteins by size. They are run using standard Laemmli sample buffer and Tris/glycine/SDS running buffer.

- Criterion Tris-HCl — Tris-HCl formulation that offer a shelf life of 12 weeks
- Criterion Stain Free — Tris-HCl gels with unique trihalo compounds that allow rapid fluorescent detection of proteins with the Criterion Stain Free system, eliminating staining and destaining steps for faster results (see Chapter 5 for more details)

Gel Composition

| | |
|----------------|---|
| Gel buffer | 0.375 M Tris-HCl, pH 8.6 |
| Crosslinker | 2.6% C |
| Stacking gel | 4% T, 2.6% C |
| Storage buffer | 0.375 M Tris-HCl, pH 8.6 |
| Shelf life | ~12 weeks at 2–8°C; expiration date is printed on each cassette |

| Gel Percentage | Optimum Separation Range |
|----------------|--------------------------|
| 5% | 100–250 kD |
| 7.5% | 40–200 kD |
| 10% | 30–150 kD |
| 12.5% | 20–120 kD |
| 15% | 10–100 kD |
| 18% | 10–50 kD |
| 4–15% | 20–250 kD |
| 4–20% | 10–200 kD |
| 8–16% | 20–120 kD |
| 10–20% | 10–100 kD |
| 10.5–14% | 25–200 kD |



Migration charts for protein standards on Criterion Tris-HCl, Criterion TGX, and TGX Stain-Free gels.

3.2.3 Criterion XT Bis-Tris Gels


Criterion XT Bis-Tris gels are based on a Bis-Tris HCl buffer system (pH 6.4) that uses discontinuous chloride and MES or MOPS ion fronts to form moving boundaries that stack and separate denatured proteins by size. This chemistry of XT Bis-Tris gels allows maximum stability and consistent results with a shelf life of at least 12 months.

Running XT Bis-Tris gels with either XT MES or XT MOPS denaturing running buffer produces different migration patterns. A combination of these two running buffers and three XT Bis-Tris gels can generate up to six different migration patterns for small to midsize proteins.

Gel Composition

| | |
|----------------|---|
| Gel buffer | Bis-Tris HCl, pH 6.4 |
| Crosslinker | 5% C |
| Stacking gel | 4% T, 5% C |
| Storage buffer | Bis-Tris HCl, pH 6.4 |
| Shelf life | 12 months at ambient temperature; expiration date is printed on each cassette |

| Gel Percentage | Optimum Separation Range | |
|----------------|--------------------------|----------------|
| | XT MES Buffer | XT MOPS Buffer |
| 10% | 2.5–200 kD | 14–220 kD |
| 12% | 1–30 kD | 6–66 kD |
| 4–12% | 2.5–200 kD | 10–300 kD |



Migration charts for protein standards on Criterion XT Bis-Tris gels.

3.2.4 Criterion XT Tris-Acetate Gels

Criterion XT Tris-acetate gels are based on a Tris-acetate buffer system (pH 7.0). It uses discontinuous acetate and Tricine ion fronts to form moving boundaries that stack and separate large denatured proteins by molecular weight.

Gel Composition

| | |
|----------------|--|
| Gel buffer | Tris-acetate, pH 7.0 |
| Crosslinker | 3.8% C |
| Stacking gel | 4% T, 3.8% C |
| Storage buffer | Tris-acetate, pH 7.0 |
| Shelf life | 8 months at 2–8°C; expiration date is printed on each cassette |

Gel Percentage Optimum Separation range

| | |
|------|-----------|
| 7% | 36–200 kD |
| 3–8% | 40–400 kD |



Migration charts for protein standards on Criterion XT Tris-Acetate gels.

3.3 SDS-PAGE Buffers

Table 3.2. Recommended Criterion precast gels and buffers for SDS-PAGE.

| Gel Type | Sample Buffer | Running Buffer |
|---------------------------|--------------------------------------|--------------------------------|
| Criterion TGX | Laemmli (catalog #161-0737) | Tris/glycine/SDS |
| Criterion TGX Stain-Free | Optional: 2-mercaptoethanol | (catalog #161-0732) |
| Criterion Tris-HCl | (catalog #161-0710) or DTT | |
| Criterion Stain Free | (catalog #161-0611) | |
| Criterion XT Bis-Tris | XT sample buffer (catalog #161-0791) | XT MES (catalog #161-0789) |
| | Optional: XT reducing agent | XT MOPS (catalog #161-0788) |
| | (catalog #161-0792) | |
| Criterion XT Tris-acetate | | XT Tricine (catalog #161-0790) |

3.3.1 Running Buffers

See Appendix B for buffer formulations. Do not adjust pH.

| | |
|-----------------------------------|--|
| Tris/glycine/SDS (1x) (pH 8.3) | 25 mM Tris, 192 mM glycine, 0.1% SDS Dilute 100 ml 10x stock (catalog #161-0732) with 900 ml diH ₂ O |
| XT MES (pH 6.4) | Dilute 50 ml 20x stock (catalog #161-0789) with 950 ml diH ₂ O |
| XT MOPS (pH 6.9) | Dilute 50 ml 20x stock (catalog #161-0788) with 950 ml diH ₂ O |
| XT Tricine (pH 8.2) | Dilute 50 ml 20x stock (catalog #161-0790) with 950 ml diH ₂ O |

3.3.2 Sample Buffers

| | |
|---------|--|
| Laemmli | 62.5 mM Tris-HCl, pH 6.8, 2% SDS, 25% glycerol, 0.01% (catalog #161-0737) bromophenol blue, 5% β-mercaptoethanol or 100 mM DTT (added fresh) |
| XT | Use XT sample buffer (catalog #161-0791) and XT reducing agent (catalog #161-0792) |

3.4 Sample Preparation

1. Determine the appropriate concentration of sample to load (depends on the load volume and the detection method used; see Chapter 12 for approximate stain sensitivities).
2. Dilute the sample with at least an equivalent volume of sample buffer with added reducing agent. For nonreducing conditions, omit the reducing agent.

TGX and Tris-HCl Gels

4.75 μl Laemmli sample buffer

0.25 μl β-mercaptoethanol

5 μl sample

10 μl total volume

XT Gels (Bis-Tris and Tris-Acetate)

5 μl XT sample buffer

1 μl XT reducing agent

1–14 μl sample

Make up to 20 μl total volume with diH₂O

3. Heat the diluted sample at 90–95°C for 5 min, or at 70°C for 10 min.

3.5 Running Conditions

Run conditions and times are approximate. Conditions may vary depending on water and buffer conductivity, which vary from one lab setting to the next. Multiply current by the number of gels run.

Table 3.3. Running conditions for SDS-PAGE with Criterion gels in the Criterion cell. Do not run different gel formulations at the same time.

| | TGX | Tris-HCl | Bis-Tris | Tris-Acetate | |
|--|------------------|------------------|----------------|----------------|----------------|
| Running buffer | Tris/glycine/SDS | Tris/glycine/SDS | XT MOPS | XT MES | XT Tricine |
| Standard Conditions | | | | | |
| Power conditions | 200 V constant | 200 V constant | 200 V constant | 200 V constant | 150 V constant |
| Expected current (per gel) | | | | | |
| Initial | 55–80 mA | 90–120 mA | 165–175 mA | 185–200 mA | 170–180 mA |
| Final | 33–43 mA | 35–55 mA | 60–70 mA | 90–110 mA | 85–95 mA |
| Run time | 42–45 min | 50–55 min | 60 min | 45 min | 65 min |
| High Voltage (Rapid) Conditions | | | | | |
| Power conditions | 300 V constant | — | — | — | — |
| Expected current (per gel) | | | | | |
| Initial | 89–135 mA | — | — | — | — |
| Final | 66–99 mA | — | — | — | — |
| Run time | 20–26 min | — | — | — | — |

4

Native PAGE

4.1 Introduction

In native PAGE, proteins are prepared in nonreducing, nondenaturing sample buffer, which maintains native structure and mass-to-charge ratios. Separation is also performed in the absence of SDS and reducing agents. Though native PAGE uses the same moving boundary described for SDS-PAGE (see Section 3.1), protein mobility depends on a number of factors besides molecular weight, including the shape and charge of the protein. Protein-protein interactions may be retained during native PAGE, so some proteins may separate as multisubunit complexes. Consequently, native PAGE is not suitable for molecular weight determination.

The nonreducing and nondenaturing environment of native PAGE allows protein separation with retention of biological activity. Because native structure is retained, native PAGE can allow resolution of proteins with the same molecular weight.

4.2 Criterion™ Gel Selection Guide for Native PAGE

Table 4.1. Criterion precast gels for SDS-PAGE.

| Gel Formulation | Gels | Description |
|--------------------|----------------------------|--|
| Laemmli-like | Criterion™ TGX™ | Laemmli-like, extended shelf life gels |
| | Criterion™ TGX Stain-Free™ | Laemmli-like, extended shelf life gels with trihalo compounds for rapid fluorescence detection |
| Tris-HCl (Laemmli) | Criterion Tris-HCl | Tris-HCl Laemmli gels |
| | Criterion Stain Free™ | Tris-HCl Laemmli gels with trihalo compounds for rapid fluorescence detection |
| Tris-acetate | Criterion™ XT Tris-acetate | Based on a Tris-acetate buffer system (pH 7.0) |

4.2.1 Criterion TGX and Criterion TGX Stain-Free Gels

Criterion TGX (Tris-Glycine eXtended shelf life) gels are Tris-HCl, Laemmli-like gels with a proprietary modification that extends their shelf life to 12 months and enhances separation characteristics relative to conventional gel types. The TGX formulation yields Laemmli-like separation patterns with exceptionally straight lanes and sharp bands and has excellent staining quality and transfer efficiency.

These gels are run using native sample buffer and Tris/glycine running buffer.

Two types of TGX formulations are available:

- Criterion TGX — Laemmli-like gels with the TGX formulation
- Criterion TGX Stain-Free — Laemmli-like, extended shelf life gels that include unique trihalo compounds that allow rapid fluorescent detection of proteins with the Criterion Stain Free system, eliminating staining and destaining steps for faster results (see Chapter 5 for more details)

Gel Composition

| | |
|--------------|--|
| Crosslinker | 2.6% C |
| Stacking gel | 4% T, 2.6% C |
| Shelf life | ~12 months at 2–8°C; expiration date is printed on each cassette |

4.2.2 Criterion Tris-HCl and Criterion Stain Free Gels

These Tris-HCl Laemmli gels are run using native sample buffer and Tris/glycine running buffer.

- Criterion Tris-HCl — Tris-HCl, Laemmli-like formulation that offer a shelf life of 12 weeks
- Criterion Stain Free — Tris-HCl, Laemmli-like gels with unique trihalo compounds that allow rapid fluorescent detection of proteins with the Criterion Stain Free system, eliminating staining and destaining steps for faster results (see Chapter 5 for more details)

Gel Composition

| | |
|----------------|---|
| Gel buffer | 0.375 M Tris-HCl, pH 8.6 |
| Crosslinker | 2.6% C |
| Stacking gel | 4% T, 2.6% C |
| Storage buffer | 0.375 M Tris-HCl, pH 8.6 |
| Shelf life | ~12 weeks at 2–8°C; expiration date is printed on each cassette |

4.2.3 Criterion XT Tris-Acetate Gels

Criterion XT Tris-acetate gels can also be used to separate proteins by their charge-to-mass ratio (under native PAGE conditions). Separation by native PAGE with XT Tris-acetate gels uses discontinuous acetate and Tricine ion fronts to form moving boundaries to stack and separate proteins by both size and charge.

Gel Composition

| | |
|----------------|--|
| Gel buffer | Tris-acetate, pH 7.0 |
| Crosslinker | 3.8% C |
| Stacking gel | 4% T, 3.8% C |
| Storage buffer | Tris-acetate, pH 7.0 |
| Shelf life | 8 months at 2–8°C; expiration date is printed on each cassette |

4.3 Native PAGE Buffers

See Appendix B for buffer formulations. Do not adjust pH unless instructed to do so.

| | |
|---------------------|--|
| Running buffer (1x) | 25 mM Tris, 192 mM glycine Dilute 100 ml 10x stock (catalog #161-0734) with 900 ml diH ₂ O |
| Sample buffer | 62.5 mM Tris-HCl, pH 6.8, 25% glycerol, 0.01% bromophenol blue (catalog #161-0738) |

4.4 Sample Preparation

In the absence of SDS, the net charge of a polypeptide is determined by its amino acid composition and the pH of the sample buffer. Only polypeptides with a net negative charge migrate into Criterion gels under native conditions. Most polypeptides have an acidic or slightly basic pI (~3–8). These proteins can be separated using the following standard protocol:

1. Determine the protein concentration and load volume of your sample based on the detection method used (see Chapter 12 for approximate stain sensitivities).
2. Dilute the sample in twice the volume of native sample buffer (DO NOT HEAT SAMPLES).

For example, combine:

| |
|--|
| 5 μ l sample |
| <u>10 μl native sample buffer (catalog #161-0738)</u> |
| 15 μ l total volume |

Strongly basic proteins (pI >8.5) have a net positive charge and will not enter a Criterion gel under native conditions. To allow polypeptides with a net positive charge to migrate into a native Criterion gel, change the polarity of the electrodes by reversing the color-coded jacks when connecting to the power supply.

4.5 Running Conditions

Running conditions for native PAGE are similar to the standard running conditions used for SDS-PAGE (see Section 3.5). If high temperature is a concern, run native PAGE at lower voltage; at lower voltages, runs require more time to complete.

Table 4.1. Running conditions for native PAGE for Criterion gels in the Criterion cell.

| | Laemmli/Laemmli-like | Tris-Acetate |
|----------------------------|----------------------|----------------|
| Running buffer | Native | Native |
| Power conditions | 200 V constant | 200 V constant |
| Expected current (per gel) | | |
| Initial | 90–120 mA | 70–80 mA |
| Final | 35–55 mA | 25–35 mA |
| Run time | 55 min | 75 min |



Criterion Stain Free™ System

5.1 Introduction

The Criterion Stain Free system, which comprises the Gel Doc™ EZ imager, Image Lab™ software, and Criterion™ TGX Stain-Free™ and Criterion Stain Free precast gels, eliminates the time-consuming staining and destaining steps required by other protein detection methods. Criterion TGX (Tris-Glycine eXtended shelf life) Stain-Free gels include a proprietary modification that extends their shelf life to 12 months and enhances separation characteristics relative to conventional gel types. Criterion TGX Stain-Free and Criterion Stain Free gels also include unique trihalo compounds that allow rapid fluorescent detection of proteins with the Gel Doc EZ imager — without staining.


The trihalo compounds in the gels react with tryptophan residues in a UV-induced reaction to produce fluorescence, which can be easily detected by the Gel Doc EZ imager within gels or on low-fluorescence PVDF membranes. Activation of the trihalo compounds in the gels adds 58 Da moieties to available tryptophan residues and is required for protein visualization. Proteins that do not contain tryptophan residues cannot be detected using this system. The sensitivity of the Criterion Stain Free system is comparable to staining with Coomassie Brilliant Blue for proteins with a tryptophan content >1.5%; sensitivity superior to Coomassie staining is possible for proteins with a tryptophan content >3%.

Molecular weights of proteins are estimated by a regression method using Image Lab software. The software generates a standard curve using the molecular weight and relative mobility (R_f) of standard proteins (R_f = distance migrated by the protein/distance migrated by the dye front). The standard curve is then used to estimate the molecular weights of sample proteins.

Benefits of the Criterion Stain Free system include:

- Elimination of staining and destaining steps for faster results
- Automated gel imaging and analysis
- No background variability within a gel or between gels (as is often seen with standard Coomassie staining)
- Reduced organic waste by not requiring acetic acid and methanol for staining or destaining
- Visualization of transferred (blotted) proteins on low-fluorescence PVDF membranes

5.2 Criterion Stain Free Workflow



5.3 Electrophoresis with Criterion TGX Stain-Free Gels

Criterion TGX Stain-Free gels (and Criterion Stain Free gels) are made and packaged without SDS, so they can be used for both SDS and native PAGE applications. To perform electrophoresis with these gels, prepare the sample and running buffers, set up the Criterion cell, and perform the run as directed in Chapters 2–4.

Use unstained standards with Criterion TGX Stain-Free and Criterion Stain Free gels, as some prestained standards are not detected by Stain-Free technology. To monitor electrophoresis, use 10 µl of a 1:1 mixture of Precision Plus Protein™ unstained and Precision Plus Protein All Blue protein standards.

5.4 Using the Gel Doc EZ Imager

Image Criterion TGX Stain-Free and Criterion Stain Free gels and blots in the Gel Doc EZ imager. The imager activates the reaction between the proteins and trihalo compounds in the gel to enable visualization.

- Immediately place the gel in the tray of the imager; no fixation or rinsing steps are required. Prolonged rinsing may diminish image quality and lead to gel deformation
- If desired, stain the gel with any TGX-compatible stains after imaging. Certain stains, if used prior to imaging, eliminate detection capability

Refer to the Gel Doc EZ Stain-Free Sample Tray Instruction Manual (bulletin 10019634) for detailed instructions.

6

Peptide Analysis

6.1 Introduction

Criterion™ Tris-Tricine/peptide gels are optimized for separating peptides and proteins with molecular weight <10,000. Peptide-SDS complexes move more slowly through these gels, allowing the faster SDS micelles that normally interfere with peptide separations to separate completely from peptides. This enables resolution of distinct peptide bands.

6.2 Criterion Tris-Tricine/Peptide Gels

6.2.1 Gel Composition

| | |
|----------------|---|
| Gel buffer | 1.0 M Tris-HCl, pH 8.45 |
| Crosslinker | 2.6% C |
| Stacking gel | 4% T, 2.6% C |
| Storage buffer | 1.0 M Tris-HCl, pH 8.45 |
| Shelf life | ~12 weeks at 2–8°C; expiration date is printed on each cassette |

6.2.2 Gel Selection Guide

Criterion Tris-Tricine/peptide gels are available in either a single percentage or a linear gradient format.

| Gel Percentage | Optimum Separation Range |
|----------------|--------------------------|
| 16.5% | 1.5–30 kD |
| 10–20% | 1–40 kD |

6.3 Tris-Tricine/Peptide Buffers

See Appendix B for buffer formulations. Do not adjust pH unless instructed to do so.

| | |
|--------------------------------------|---|
| Running buffer (1x) | 100 mM Tris, 100 mM Tricine, 0.1% SDS Dilute 100 ml 10x stock (catalog #161-0744) with 900 ml diH ₂ O |
| Sample buffer (catalog #161-0739) | 200 mM Tris-HCl, pH 6.8, 2% SDS, 40% glycerol, 0.04% Coomassie Brilliant Blue G-250, 2% β-mercaptoethanol or 100 mM DTT (added fresh) |

6.4 Sample Preparation

1. Determine the appropriate concentration of sample to load (depends on the load volume and the detection method used; see Chapter 12 for approximate stain sensitivities).
2. Dilute the sample with at least an equivalent volume of sample buffer (catalog #161-0739) and reducing agent (β -mercaptoethanol, for example). Heat the diluted sample at 90–95°C for 5 min, or at 70°C for 10 min.

For example, combine:

5 μ l sample

4.75 μ l Tricine sample buffer (catalog #161-0739)

0.25 μ l β -mercaptoethanol (catalog #161-0710)

10 μ l total volume

6.5 Running Conditions

| | |
|------------------|----------------|
| Power conditions | 125 V constant |
| Starting current | 140–150 mA/gel |
| Final current | 60–70 mA/gel |
| Run time | 120 min |

7

Isoelectric Focusing (IEF)

7.1 Introduction

Isoelectric focusing (IEF) separates proteins by their net charge rather than molecular weight. Criterion™ IEF gels are cast with Bio-Rad's Bio-Lyte® ampholytes, amphoteric molecules that generate a pH gradient across the gels. Proteins migrate to their isoelectric point (pI), the pH at which the protein has no net charge. Criterion IEF gels contain no denaturing agents, so IEF is performed under native conditions.

7.2 Criterion IEF Gels

7.2.1 Gel Composition

| | |
|----------------|---|
| Gel buffer | 2% ampholyte, pH 3–10 or 5–8 |
| Crosslinker | 3.0% C |
| Stacking gel | None |
| Storage buffer | diH ₂ O |
| Shelf life | ~24 weeks at 2–8°C; expiration date is printed on each cassette |

7.2.2 Gel Selection Guide

| IEF gel | pH Range |
|---------|----------|
| 5–8 | 5–8.0 |
| 3–10 | 4–8.5 |

7.3 IEF Buffers

See Appendix B for buffer formulations. Do not adjust pH unless instructed to do so.

Running buffers:

| | |
|-------------------------|--|
| IEF cathode buffer (1x) | 20 mM lysine (free base), 20 mM arginine (free base) Dilute 100 ml 10x stock (catalog #161-0762) with 900 ml diH ₂ O |
| IEF anode buffer (1x) | 7 mM phosphoric acid Dilute 100 ml 10x stock (catalog #161-0761) with 900 ml diH ₂ O |
| Sample buffer | 50% glycerol |

7.4 Sample Preparation

1. Determine the appropriate concentration of sample to load (depends on the load volume and the detection method used).
2. Dilute the sample with at least an equivalent volume of sample buffer.

For example, combine:

| |
|-------------------------|
| 5 μ l sample |
| 5 μ l sample buffer |
| <hr/> |
| 10 μ l total volume |

7.5 Running Conditions

| | |
|-----------------------------|---|
| Power conditions (stepwise) | 100 V constant 60 min 250 V constant 60 min 500 V constant 30 min |
| Starting current | 5–25 mA/gel |
| Final current | 5–25 mA/gel |
| Run time | 150 min |



Protease Analysis by Zymogram PAGE

8.1 Introduction

Criterion™ zymogram gels are used to test for proteolytic activity. Gels are cast with gelatin or casein, which acts as a substrate for proteases that are separated in the gel under nonreducing conditions. Proteases are detected by first renaturing the enzymes and then allowing them to break down the substrate. Zymogram gels are stained with Coomassie Brilliant Blue R-250 stain, which stains the substrate while leaving clear areas around active proteases.

8.2 Criterion Zymogram Gels

8.2.1 Gel Composition

| | |
|----------------|--|
| Gel buffer | 0.375 M Tris-HCl, pH 8.6 |
| Crosslinker | 2.6% C |
| Stacking gel | 4% T, 2.6% C |
| Storage buffer | 0.375 M Tris-HCl, pH 8.6, 0.2% NaN ₃ |
| Shelf life | ~8 weeks at 2–8°C; expiration date is printed on each cassette |

8.2.2 Gel Selection Guide

| Zymogram Gel | Optimum Separation Range |
|-------------------|--------------------------|
| 10% with gelatin | 30–150 kD |
| 12.5% with casein | 20–120 kD |

8.3 Zymogram Buffers

See Appendix B for buffer formulations. Do not adjust pH unless instructed to do so.

| | |
|--------------------------------------|--|
| Running buffer (1x) | 25 mM Tris, 192 mM glycine, 0.1% SDS Dilute 100 ml 10x stock (catalog #161-0732) with 900 ml diH ₂ O |
| Sample buffer (catalog #161-0764) | 62.5 mM Tris-HCl, pH 6.8, 4% SDS, 25% glycerol, 0.01% Coomassie Brilliant Blue G-250 |

8.4 Sample Preparation

1. Determine the appropriate protein concentration of your sample based on the detection method and load volume used. (See Chapter 12 for approximate stain sensitivities.)
2. Dilute 1 part sample with 1 part sample buffer. Do not heat the samples.

8.5 Running Conditions

| | |
|------------------|----------------|
| Power conditions | 125 V constant |
| Starting current | 90–120 mA/gel |
| Final current | 35–55 mA/gel |
| Run time | 90 min |



Nondenaturing Nucleic Acid PAGE

9.1 Introduction

Criterion™ TBE gels are used to separate small double-stranded DNA (dsDNA) fragments, particularly PCR products. DNA molecules have nearly uniform mass-to-charge ratios, allowing nondenaturing nucleic acid PAGE to separate dsDNA by mass using a continuous TBE buffer system.

9.2 Criterion TBE Gels

9.2.1 Gel Composition

| | |
|----------------|---|
| Gel buffer | 89 mM Tris, 89 mM boric acid, 2 mM EDTA, pH 8.3 |
| Crosslinker | 3.3% C |
| Stacking gel | 4% T, 3.3% C |
| Storage buffer | 89 mM Tris, 89 mM boric acid, 2 mM EDTA |
| Shelf life | ~12 weeks at 2–8°C; expiration date is printed on each cassette |

9.2.2 Gel Selection Guide

| Gel Percentage | Optimum Separation Range |
|----------------|--------------------------|
| 5% | 200–2,000 bp |
| 10% | 50–1,500 bp |
| 15% | 20–1,000 bp |
| 4–20% | 10–2,000 bp |

9.3 Nondenaturing Nucleic Acid PAGE Buffers

See Appendix B for buffer formulations. Do not adjust pH unless directed to do so.

| | |
|---|---|
| Running buffer (1x) | 89 mM Tris, 89 mM boric acid, 2 mM EDTA Dilute 100 ml 10x stock (catalog #161-0733) with 900 ml diH ₂ O |
| Sample buffer (5x) (catalog #161-0767) | 50 mM Tris-HCl, pH 8.0, 5 mM EDTA, 25% glycerol, 0.2% bromophenol blue, 0.2% xylene cyanole FF |

9.4 Sample Preparation

Determine the DNA concentration of your sample based on the detection method used. (See Chapter 12 for approximate stain sensitivities.) Dilute 4 parts sample with 1 part sample buffer.

9.5 Running Conditions

Table 9.1. Running conditions for nondenaturing nucleic acid PAGE with Criterion gels in the Criterion cell.

| | 5% and 10% Gels | 15% and 4–20% Gels |
|----------------------------|-----------------|--------------------|
| Power conditions | 100 V constant | 150 V constant |
| Expected current (per gel) | | |
| Initial | 20–25 mA | 27–35 mA |
| Final | 14–18 mA | 20–35 mA |
| Run time | 90 min | 90 min |

10

Denaturing Nucleic Acid PAGE

10.1 Introduction

Criterion™ TBE-urea gels are used for separation of small RNA and single-stranded DNA (ssDNA) fragments. Applications include oligonucleotide analysis, RNase protection assays, and northern blotting.

10.2 Criterion TBE-Urea Gels

10.2.1 Gel Composition

| | |
|----------------|--|
| Gel buffer | 89 mM Tris, 89 mM boric acid, 2 mM EDTA, 7 M urea, pH 8.3 |
| Crosslinker | 3.3% C |
| Stacking gel | 4% T, 3.3% C |
| Storage buffer | 89 mM Tris, 89 mM boric acid, 2 mM EDTA, pH 8.3 |
| Shelf life | ~8 weeks at 2–8°C; expiration date is printed on each cassette |

10.2.2 Gel Selection Guide

| Gel Percentage | Optimum Separation Range |
|----------------|--------------------------|
| 5% | 50–1,000 nt |
| 10% | 25–300 nt |
| 15% | 10–50 nt |

10.3 TBE-Urea Buffers

See Appendix B for buffer formulations. Do not adjust pH unless directed to do so.

| | |
|---|---|
| Running buffer (1x) | 89 mM Tris, 89 mM boric acid, 2 mM EDTA Dilute 100 ml 10x stock (catalog #161-0733) with 900 ml diH ₂ O |
| Sample buffer (5x) (catalog #161-0768) | 89 mM Tris, 89 mM boric acid, 2 mM EDTA, pH 8.0, 12% Ficoll, 0.01% bromophenol blue, 0.02% xylene cyanole FF, 7 M urea |

10.4 Sample Preparation

Determine the desired ssDNA or RNA concentration for your sample based on the detection method used. Dilute 4 parts sample with 1 part sample buffer.

10.5 Running Conditions

Table 10.1. Running conditions for denaturing nucleic acid PAGE with Criterion gels in the Criterion cell.

| | 5% Gels | 10% Gels | 15% Gels |
|----------------------------|----------------|-----------------|-----------------|
| Power conditions | 200 V constant | 200 V constant | 200 V constant |
| Expected current (per gel) | | | |
| Initial | 40–45 mA | 20–33 mA | 18–22 mA |
| Final | 20–25 mA | 14–18 mA | 10–15 mA |
| Run time | 90 min | 90 min | 90 min |

11

2-D Electrophoresis

11.1 Introduction

Criterion™ precast gels are available for second-dimension PAGE in 2-D electrophoresis workflows. The IPG-well gels accommodate 11 cm IPG strips. Criterion™ Any kD™ gels are particularly well suited for 2-D electrophoresis applications.

The transition from first- to second-dimension gel electrophoresis involves:

- Equilibration of the resolved IPG strips in a reducing buffer containing SDS
- Placing the IPG strip on top of the second-dimension gel

11.2 Equilibration

Equilibration ensures that proteins in the IPG strips are coated with SDS and that cysteines are reduced and alkylated. Use the equilibration protocols (bulletin 4110009) and buffers in the ReadyPrep™ 2-D starter kit (catalog #163-2105), or other protocols and buffers used with Tris-HCl gels.

11.3 Agarose Overlay

Place the equilibrated IPG strip into the IPG well of the Criterion gel and overlay it with molten agarose to ensure good contact between the strip and gel.

- Criterion TGX™ and Tris-HCl gels: prepare 0.5% low-melt agarose (catalog #161-3111), 0.003% bromophenol blue (catalog #161-0404) in 1x Tris/glycine/SDS running buffer. (Alternatively, use ReadyPrep overlay agarose, catalog #163-2111)
 - Criterion™ XT gels: prepare 0.5% low-melt agarose (catalog #161-3111), 0.001% bromophenol blue (catalog #161-0404) in appropriate 1x XT running buffer
1. Following equilibration, place the IPG strip, gel side up, on the back plate of the Criterion gel, above the IPG well. The “+” and pH range on the IPG strip should be on the left.
 2. Using forceps, push the strip into the IPG well, taking care to not trap air bubbles under the strip. Push on the backing of the strip, not on the gel.
 3. Using a disposable pipet, add overlay agarose to the IPG well. Fill the well to the top of the inner plate. Dispense rapidly, as overlay agarose solidifies quickly. To avoid bubbles, tilt the cassette slightly to allow bubbles to escape. Push gently on the plastic backing of the strip to free any trapped bubbles.

11.4 Second-Dimension Electrophoresis

Place the cassettes into the Criterion cell and start the run using the run conditions for SDS-PAGE. Use the migration of the bromophenol blue in the overlay agarose to monitor the progress of the run.

12 Detection

12.1 SDS-PAGE and Native PAGE Detection

Following electrophoresis, either stain the gel or use the Criterion Stain Free™ system to visualize proteins in the gel.

- Refer to Table 12.1 for a comparison of total protein stains
- For Criterion™ TGX Stain-Free™ and Criterion Stain Free™ gels, immediately place the gel on the tray of the Gel Doc™ EZ imager; no additional fixation or rinsing steps are required. If desired, stain with any compatible stains (Table 12.1) following imaging. Some stains, if used prior to imaging, can impair imaging quality

Table 12.1. Total protein gel stains for use with Criterion gels.

| Stain | Sensitivity (Lower Limit) | Optimum Protein Load (µg/Band) | Advantages | Disadvantages | Imaging | Manual |
|--|---------------------------|--------------------------------|---|---|---|--------------------|
| Criterion Gels | | | | | | |
| Coomassie R-250 | 36–47 ng | ~0.5 | Laboratory standard | Requires methanol destaining | Photography with white light or transmission densitometry | Consult literature |
| Bio-Safe™ Coomassie | 8–28 ng | ~0.5 | Nonhazardous | | | 4307051 |
| Zinc stain ¹ | 6–12 ng | ~0.2 | High contrast, fast, reversible | Negative SDS-PAGE stain, must be photographed | | 4006082 |
| Silver Stain Plus™ kit | 0.6–1.2 ng | ~0.01 | Sensitive, robust, mass spectrometry compatible | Does not stain glycoproteins well | | LIT442 |
| Silver stain | 0.6–1.2 ng | ~0.01 | Stains complex proteins (glyco- or lipoproteins) | Not mass spectrometry compatible | | LIT34 |
| Dodeca™ silver stain kit | 0.5–1.2 ng | ~0.1 | Convenient staining for a large number of gels | | | 4110150 |
| Oriole™ fluorescent gel stain ¹ | ~2 ng | ~0.1 | High sensitivity, broad dynamic range, simple one-step protocol | | Fluorescence visualization with UV transillumination | 10017295 |

¹ Do not use zinc stain or Oriole fluorescent gel stain to stain native gels.

| Stain | Sensitivity (Lower Limit) | Optimum Protein Load (µg/Band) | Advantages | Disadvantages | Imaging | Manual |
|---|---------------------------|--------------------------------|---|---|--|----------|
| SYPRO Ruby protein gel stain | 1–10 ng | ~0.1 | Broad dynamic range | Requires laser- or LED-based imaging instrument for maximum sensitivity | Fluorescence visualization with UV, LED, or laser scanning | 4006173 |
| Flamingo™ fluorescent gel stain | 0.25–0.5 ng | ~0.02 | Broad dynamic range, mass spectrometry compatible | | | 10003321 |
| Criterion TGX Stain-Free and Criterion Stain Free Gels | | | | | | |
| Stain Free imaging | 2–28 ng | ~0.5 | Rapid (<5 min), compatible with blotting and mass spectrometry, simple protocol with no additional reagents | Requires tryptophan residues in proteins for detection | Fluorescence visualization using Criterion Stain Free imaging system | 10014472 |

12.2 Peptide Gel Staining

Peptides and small proteins are prone to diffusion and loss during staining. The following protocol includes a fixing step prior to staining to prevent sample loss and is suitable for detection of bands as low as 10–20 ng.

Fixative solution 40% methanol, 10% acetic acid

Stain solution 0.025% (w/v) Coomassie Blue G-250, 10% acetic acid

Destain solution 10% acetic acid

Place gels in fixative solution and equilibrate for 30 min. Stain gels with stain solution for 1 hr. Stain should be used only once; reuse may result in loss of sensitivity. Destain gels three times for 15 min or until the desired background is achieved. Some peptides may not be completely fixed and may diffuse out of the gels if fixing and staining times are greatly exceeded.

12.3 IEF Gel Staining

Samples on IEF gels can be detected using multiple methods. Use Table 12.2 as a guide to selecting an appropriate staining method.

Table 12.2. IEF gel detection methods.

| Method | Sensitivity (Lower Limit) | Optimum Protein Load (µg/Band) | Advantages | Disadvantages |
|--------------------------------|---------------------------|--------------------------------|---|---|
| Silver Stain Plus kit | 0.6–1.2 ng | ~0.01 | Sensitive, robust, mass spectrometry compatible | Requires TCA fixation |
| Silver stain | 0.6–1.2 ng | ~0.01 | Stains complex proteins (glyco- or lipoproteins) | Not mass spectrometry compatible |
| SYPRO Ruby protein gel stain | 1.0–10 ng | ~0.1 | Broad dynamic range, mass spectrometry compatible | Requires laser- or LED-based imaging instrument for maximum sensitivity |
| Flamingo fluorescent gel stain | 0.25–0.5 ng | ~0.2 | Broad dynamic range, mass spectrometry compatible | Requires laser- or LED-based imaging instrument for maximum sensitivity |

12.4 Zymogram Gel Staining

Prior to staining zymogram gels, sample proteases must first be renatured and allowed to break down the substrate contained in the gel. The following protocol provides basic guidelines for detection. Optimal results should be determined empirically.

| | |
|----------------------|--|
| Renaturing solution | 2.5% Triton X-100 |
| Development solution | 50 mM Tris, 200 mM NaCl, 5 mM CaCl ₂ (anhydrous), 0.02% Brij-35 Adjust to pH 7.5 |
| Staining solution | 40% methanol, 10% acetic acid, 0.5% Coomassie Blue R-250 |
| Destaining solution | 40% methanol, 10% acetic acid |

Place gels in renaturing solution for 30 min at room temperature. Incubate gels in development solution at 37°C for a minimum of 4 hr. Highest sensitivity is typically achieved with overnight incubation. Optimum conditions should be determined empirically. Stain gels with staining solution for at least 1 hr at room temperature. Destain until clear bands appear against the blue background (~30–60 min).

12.5 TBE Gel Staining

Use Table 12.3 as a guide to selecting an appropriate staining method.

Table 12.3. TBE gel detection methods.

| Method | Sensitivity (Lower Limit) | Advantages | Disadvantages |
|------------------|---------------------------|--------------------------------------|-------------------------------|
| Ethidium bromide | 50 ng | Classic fluorescent DNA stain | Carcinogenic |
| Silver stain | 1–2 ng | More sensitive than ethidium bromide | Requires multiple steps |
| SYBR® Green | 0.02–2 ng | High sensitivity | Multiple steps, –20°C storage |
| SYBR® Safe | 0.5 ng | Non-hazardous | Multiple steps |

12.6 TBE-Urea Gel Staining

Use Table 12.4 as a guide to selecting an appropriate staining method.

Table 12.4. TBE-urea gel detection methods.

| Method | Sensitivity (Lower Limit) | Advantages | Disadvantages |
|------------------|---------------------------|--------------------------------------|-------------------------|
| Ethidium bromide | 10 ng | Classic fluorescent DNA stain | Carcinogenic |
| Radiant™ Red | 10 ng | Fast, single-step protocol | RNA and ssDNA only |
| Silver stain | 1–2 ng | More sensitive than ethidium bromide | Requires multiple steps |

13 Blotting

13.1 Introduction

Western blotting is an electrophoretic technique used to move proteins from a gel onto a solid support such as a nitrocellulose or PVDF membrane. The membrane can be used for immunological or biochemical analyses or demonstration of protein-protein or protein-ligand interactions. Below are guidelines for western blotting of Criterion™ precast gels onto nitrocellulose or PVDF membranes using either wet or semi-dry transfer techniques.

Assess transfer efficiency using a total protein blot stain such as SYPRO Ruby stain (see Table 12.1). With Criterion™ TGX Stain-Free™ and Criterion Stain Free™ gels, transfer efficiency to low-fluorescence PVDF membranes may also be assessed using the Gel Doc™ EZ imager (see Chapter 5; activate the gel before blotting).

13.2 Transfer

13.2.1 Transfer Buffers

See Appendix B for buffer formulations. Do not adjust pH unless directed to do so.

Towbin buffer (1x) 25 mM Tris, 192 mM glycine, 20% (v/v) methanol (pH 8.3)
Dilute 100 ml 10x stock (catalog #161-0734) with 400 ml diH₂O.
Add 200 ml methanol, then adjust volume to 1 L with diH₂O.


Add SDS to 0.1% to promote transfer of high molecular weight proteins.

13.2.2 Wet Transfer Using the Criterion Blotter

1. Equilibrate the gels and membranes (for example, in Towbin buffer) for 15 min prior to blot assembly.
2. Assemble the transfer apparatus:
 - a) Fill the tank with transfer buffer to ~50% of the fill volume and place a magnetic stir bar inside the tank.
 - b) Place the ice block in the pocket in the back of the cell. Flip down the lever to hold the ice block. Alternatively, connect the optional cooling coil to an appropriate recirculating water chiller and place it in the grooves in the back of the tank.
3. Assemble the cassette:
 - a) Pour chilled transfer buffer into each compartment of the assembly tray, and then place the membrane (nitrocellulose, PVDF) into the front (small) compartment of the tray to soak.

Wet PVDF membranes in methanol before soaking in transfer buffer.

- b) Open the cassette and place it in the back (large) compartment of the tray so the red plate with the handle is vertical (anode) and the black plate (cathode) is horizontal and submerged in transfer buffer.
- c) Assemble the sandwich as shown, placing the gel closest to the black side of the cassette and the membrane closest to the red plate. Use a blot roller to remove air trapped between the layers of the blot assembly.




Assembly of the gel blot sandwich with the Criterion blotter.

4. Place the assembled cassette into the groove in the tank, aligning the red side of the card with the red electrode. Make sure the magnetic stirbar is free to move. Repeat steps 2–4 for a second blot.
5. Add the remaining transfer buffer to the fill level marked on the tank, place the tank on a stir plate, and begin stirring to maintain even buffer temperature and ion concentration during the transfer.
6. Connect the Criterion blotter to a PowerPac™ HC power supply and begin transfer.

For many proteins, excellent transfer efficiency is obtained in 30 min at a constant voltage of 100 V. For best results, optimize conditions for proteins of interest. Large proteins (>150 kD) may take 60 min, while smaller proteins (<30 kD) may transfer in 20 min. Refer to the Criterion Blotter Instruction Manual (bulletin 4006190) or the Protein Blotting Guide (bulletin 2895) for additional information.

13.2.3 Transfer Using the Trans-Blot® Turbo™ System

1. Open the transfer pack and assemble the components on the cassette in the order shown. For best results, use the roller to remove any air trapped between the layers. If using a single mini or midi sandwich, place the sandwich in the middle of the cassette bottom. With two mini gels, place the sandwiches on a midi stack with the foot of each gel facing the center of the stack.
2. Place the lid on the cassette and lock the lid in place by turning the knob clockwise, using the symbols on the lid as a guide. Slide the cassette into the appropriate bay of the Trans-Blot Turbo cell. Each cassette and bay can hold up to two mini gels or one midi gel (Table 13.1).
3. Start the transfer. With the cassette in the cell, press TURBO and select the gel type. Press A:RUN or B:RUN to begin the transfer. Press LIST to select a Bio-Rad optimized protocol (Table 13.2) or a user-defined protocol. Press NEW to create and run a new protocol.
4. When transfer completes, RUN COMPLETE appears. Pull the cassette straight out of the slot and unlock the lid. Disassemble the blotting sandwich.



Assembly of the gel blot sandwich with the Trans-Blot Turbo system.

Table 13.1. Placement of cassettes in the Trans-Blot Turbo cell.

| | Acceptable | | Unacceptable | |
|---------------|------------------------|---|-----------------------------|--------------------------------------|
| | Option 1 | Option 2 | Option 1 | Option 2 |
| Upper bay (A) | 1 mini gel -and/or- | 2 mini gels -or- 1 midi gel -and/or- | 1 mini gel -and- | 2 mini gels -or- 1 midi gel -and- |
| Lower bay (B) | 1 mini gel | 2 mini gels -or- 1 midi gel | 2 mini gels -or- 1 midi gel | 1 mini gel |

Table 13.2. Trans-Blot Turbo transfer protocols.

| Protocol Name | MW, kD | Time, Min | 1 Mini Gel | 2 Mini Gels or 1 Midi Gel |
|---------------|--------|-----------|----------------------------|----------------------------|
| | | | | |
| STANDARD SD | Any | 30 | Up to 1.0 A, 25 V constant | Up to 1.0 A, 25 V constant |
| 1.5 MM GEL | Any | 10 | 2.5 A constant, up to 25 V | 1.3 A constant, up to 25 V |
| HIGH MW | >150 | 10 | 2.5 A constant, up to 25 V | 1.3 A constant, up to 25 V |
| LOW MW | <30 | 5 | 2.5 A constant, up to 25 V | 1.3 A constant, up to 25 V |
| MIXED MW | 5–150 | 7 | 2.5 A constant, up to 25 V | 1.3 A constant, up to 25 V |
| 1 Mini TGX | 5–150 | 3 | 2.5 A constant, up to 25 V | N/A |


Refer to the Trans-Blot Turbo Instruction Manual (bulletin 10020688) for complete instructions or the Protein Blotting Guide (bulletin 2895) for additional information.

13.2.4 Semi-Dry Transfer Using the Trans-Blot® SD Cell

1. Equilibrate the gels and membranes (for example, in transfer buffer; see Appendix B for buffer recipes) for 20 min prior to blot assembly.
2. Assemble the blot for transfer using the Trans-Blot SD semi-dry transfer system.
3. Connect the Trans-Blot SD cell to a PowerPac HC power supply and begin transfer at 10–15 V.

For most proteins transferred from Criterion precast gels, optimum transfer efficiency is obtained in 30 min; smaller proteins (<30 kD) may transfer more quickly, while proteins >150 kD may show increased transfer efficiencies at up to 60 min. Run times longer than 60 min are NOT recommended for semi-dry transfers.

Refer to the Trans-Blot SD Instruction Manual (bulletin 1703940) or the Protein Blotting Guide (bulletin 2895) for additional information.



Assembly of the gel blot sandwich with the Trans-Blot SD cell.

13.3 Total Protein Blot Stains

Total protein staining of a membrane provides an image of the complete protein pattern, which is required for the full characterization of specific antigens detected in complex protein mixtures. Gels shrink during staining, so comparison of an immunologically probed membrane to a stained gel is not practical. Instead, the exact location of a specific antigen is determined by comparing two blotted membranes: one that has been probed with an antibody and the other stained for total protein.

Table 13.1. Total protein blot stains.

| Method | Sensitivity | Protein Load (µg/Band) | Advantages | Disadvantages | Imaging |
|--|--------------|------------------------|--|---|--|
| SYPRO Ruby protein blot stain | 2–8 ng | ~0.2 | Compatible with mass spectrometry, Edman-based sequencing, and standard immunological procedures | Multi-step protocol; requires UV, LED, or laser imaging for maximum sensitivity | Fluorescence visualization with UV, LED epi-illumination or laser scanning |
| Colloidal gold stain | 1 ng | ~0.1 | Highly sensitive, single-step protocol | Incompatible with nylon membranes | Photography with epi-illumination or reflectance densitometry |
| Anionic dyes (amido black, Coomassie R-250, Ponceau S, Fast Green FCF) | 100–1,000 ng | ~5.0 | Inexpensive, rapid | Low sensitivity | Photography with epi-illumination or reflectance densitometry |

To visualize total protein on blots using the Gel Doc EZ imager, refer to Section 5.4.

13.4 Immunodetection

After transfer, blots are ready for downstream processing. While all protein and antibody combinations are different and may require optimization, a general protocol for the immunodetection of a large number of protein and antibody combinations is provided (see Appendix B for buffer formulations).

1. Immediately after transfer, place the membrane into Tris-buffered saline with Tween 20 (TTBS) containing blocking agent (for example, 3% BSA, 5% nonfat dry milk, 1% casein, or 1% gelatin). Incubate at room temperature with agitation for 1 hr.
2. Dilute the primary antibody in blocking solution (suggested dilution is specified by the manufacturer). Incubate the blot in the primary antibody solution at room temperature and with agitation for 1 hr.
3. Wash the blot with TTBS as directed in the instructions for the detection method to be used (for example, 5 times, 5 min each at room temperature).
4. Dilute the secondary antibody into TTBS as specified by the manufacturer. Incubate the blot in the secondary antibody solution at room temperature with agitation for 1 hr.
5. Wash the blot with TTBS for 5 min at room temperature with agitation. Pour off the wash solution and repeat 5 times.
6. Follow the directions for the detection kit used to develop the blot. For the Immun-Star™ WesternC™ chemiluminescence kit (catalog #170-5070), mix 3 ml luminol/enhancer with 3 ml peroxide solution to make a 1x working solution for a 7 x 8.5 cm membrane. Incubate the membrane in the solution for 3–5 min. Prior to imaging, drain the excess substrate and place the membrane in a protective sleeve (such as plastic wrap) to prevent drying.

14

Troubleshooting

Table 14.1. Troubleshooting electrophoresis and detection with Criterion™ gels. For more troubleshooting tips, refer to the Criterion cell, Criterion blotter, and Trans-Blot® SD cell instruction manuals, or contact Technical Support.

| Problem | Cause | Solution |
|---|--|---|
| General Troubleshooting Tips | | |
| Current is zero or less than expected, and samples do not migrate into gel | Tape at bottom of cassette not removed | Remove tape |
| | Insufficient buffer in integral buffer chamber | Fill buffer chamber with running buffer |
| | Insufficient buffer in lower buffer chamber | Fill both halves of lower buffer tank with 400 ml running buffer when running two gels |
| | Electrical disconnection | Check electrodes and connections |
| Gels run faster than expected | Running buffer too concentrated or incorrect | Check buffer composition |
| | Gel temperature too high | Do not exceed recommended running conditions |
| Bands “smile” across gel: band pattern curves upward at both sides of the gel | Excessive heating of gel | Check buffer composition Do not exceed recommended running conditions |
| | Insufficient buffer | Fill both halves of lower buffer tank with 400 ml running buffer when running two gels |
| Bands “smile” or “frown” within gel lanes | Protein load too high | Load less protein |
| | Sample or buffer preparation issues | Minimize salts, detergents, and solvents in sample preparation and sample loading buffers |
| | Incorrect running conditions | Set correct voltage |
| Bands are skewed or distorted; lateral band spreading | Too much salt in samples | Remove salt from samples (dialysis, precipitation, or other method) |
| | Insufficient or wrong sample buffer | Check buffer composition and dilution instructions |
| | Sample precipitation | Selectively remove predominant proteins in sample Dilute sample in sample buffer |
| | Insoluble materials (for example, cell membranes) in samples | Centrifuge samples to remove particulates prior to sample loading |

| Problem | Cause | Solution |
|---|---|---|
| Artifactual bands at 60–70 kD | Skin keratin contamination | Clean all dishware, and wear gloves while handling and loading gels Filter all solutions (0.2–0.45 µm filter) |
| Poor resolution or fuzzy bands | Sample volume too high | If possible, load a more concentrated sample in a lower sample buffer volume |
| | Not enough running buffer in lower chamber | Add running buffer to fill line of lower buffer chamber |
| | Diffuse sample loading zone | Load sample with a syringe or gel loading pipet tip |
| | Sample diffusion during staining | Fix gel with 40% methanol, 10% acetic acid for 80 min prior to staining |
| | Incompatible sample components | Minimize salts, detergents, and solvents in sample preparation and sample loading buffers |
| | Expired gel | Use gels before expiration date printed on the cassette |
| Criterion™ TGX Stain-Free™ and Criterion Stain Free™ Gels | | |
| Low sensitivity for proteins | Low tryptophan content in proteins | After activation and imaging, stain gel with Bio-Safe™ Coomassie or similar to detect missing bands |
| Uneven sensitivity or fuzzy bands | Gel was soaked in water or buffer prior to activation and imaging | If possible, activate and image gel using the Gel Doc EZ imager immediately after electrophoresis |
| Bands are too light or missing from blot (membrane) | Proteins transferred through membrane | Use membrane with smaller pore size Decrease transfer time Decrease voltage |
| Standards not visible | Incorrect standards were used | Use unstained standards; some prestained standards are not detected by the Gel Doc EZ imager. To monitor electrophoresis, use a 1:1 mixture of unstained and prestained standards |
| Dye front at bottom of gels interfering with detection of proteins | Sample constituents present in gel interfering with imaging | Dilute sample in gel running buffer prior to loading Activate and image gel, rinse in fixation solution for 30 min, and repeat imaging |
| High background or low sensitivity on a blot of a Criterion TGX Stain-Free or Criterion Stain Free gel (imaged using the Gel Doc EZ imager) | Low tryptophan content in proteins | After activation and imaging, stain blot with a total protein blot stain to detect missing bands |
| | Membrane not low-fluorescence PVDF | Always use low-fluorescence PVDF membranes to blot Criterion TGX Stain-Free or Criterion Stain Free gels |
| | PVDF membrane is dry | Wet PVDF membrane briefly in methanol and wash in water for 1 min before imaging |



Quick Start Guide

The following instructions are for electrophoresis of Criterion™ precast gels using the Criterion system.

Prepare Buffers

Running buffer (1x) Dilute 100 ml 10x stock (catalog #161-0732) with 900 ml diH₂O.

Sample buffer Use Laemmli sample buffer (catalog #161-0737).

Prepare Gels and Assemble Electrophoresis Cell

Remove the comb and tape from the gels and assemble the electrophoresis cell.

Fill the inner and outer buffer chambers with running buffer.

Fill the lower buffer chamber with 400 ml running buffer (to the fill mark).

Prepare and Load Samples

| Component | Reducing | Nonreducing |
|--|----------|-------------|
| Sample | 5 µl | 5 µl |
| Laemmli sample buffer (catalog #161-0737) | 4.75 µl | 5 µl |
| β-Mercaptoethanol | 0.25 µl | — |
| Total volume | 10 µl | 10 µl |

Heat samples at 90–100°C for 5 min or at 70°C for 10 min.

Load the appropriate amount of sample on the gel.

Perform Electrophoresis

Connect the electrophoresis cell to the power supply and perform electrophoresis according to the conditions in Table A.1.

Table A.1. Running conditions for SDS-PAGE in the Criterion cell. Do not run different gel formulations at the same time.

| | TGX | Tris-HCl | Bis-Tris | | Tris-Acetate |
|--|------------------|------------------|----------------|----------------|----------------|
| Running buffer | Tris/glycine/SDS | Tris/glycine/SDS | XT MOPS | XT MES | XT Tricine |
| Standard Conditions | | | | | |
| Power conditions | 200 V constant | 200 V constant | 200 V constant | 200 V constant | 150 V constant |
| Expected current (per gel) | | | | | |
| Initial | 55–80 mA | 90–120 mA | 165–175 mA | 185–200 mA | 170–180 mA |
| Final | 33–43 mA | 35–55 mA | 60–70 mA | 90–110 mA | 85–95 mA |
| Run time | 42–45 min | 50–55 min | 60 min | 45 min | 65 min |
| High Voltage (Rapid) Conditions | | | | | |
| Power conditions | 300 V constant | — | — | — | — |
| Expected current (per gel) | | | | | |
| Initial | 89–135 mA | — | — | — | — |
| Final | 66–99 mA | — | — | — | — |
| Run time | 20–26 min | — | — | — | — |

Table A.2. Standard transfer (blotting) conditions.

| Method | Membrane | Time | Voltage |
|-------------------|---|-----------|--------------------|
| Tank blotting | Nitrocellulose (0.45 μ m) PVDF (0.2 μ m) | 20–60 min | 100 V ¹ |
| Semi-dry blotting | Nitrocellulose (0.45 μ m) PVDF (0.2 μ m) | 20–60 min | 25 V |

¹ For previously optimized protocols for proteins <30 kD, reduce transfer time by 20%.



Buffers

Running Buffers

10x SDS-PAGE (1 L)

(catalog #161-0732)

250 mM Tris, 1.92 M glycine, 1% SDS, pH 8.3

| | |
|--------------------|---------|
| Tris base | 30.3 g |
| Glycine | 144.1 g |
| SDS | 10 g |
| diH ₂ O | to 1 L |

Do not adjust the pH (~pH 8.3)

10x Native PAGE (1 L)

(catalog #161-0734)

250 mM Tris, 1.92 M glycine, pH 8.3

| | |
|--------------------|---------|
| Tris base | 30.3 g |
| Glycine | 144.1 g |
| diH ₂ O | to 1 L |

Do not adjust the pH (~pH 8.3)

10x Tris-Tricine (1 L)

(catalog #161-0744)

1 M Tris, 1 M Tricine, 1% SDS, pH 8.3

| | |
|--------------------|---------|
| Tris base | 121.1 g |
| Tricine | 179.2 g |
| SDS | 10 g |
| diH ₂ O | to 1 L |

Do not adjust the pH (~pH 8.3)

10x TBE (1 L)

(catalog #161-0741)

890 mM Tris, 890 mM boric acid, 20 mM EDTA

| | |
|--------------------|---------|
| Tris base | 107.8 g |
| Boric acid | 55.0 g |
| EDTA | 5.8 g |
| diH ₂ O | to 1 L |

Do not adjust the pH (~pH 8.3)

Sample Buffers

2x SDS-PAGE (Laemmli, 30 ml)
(catalog #161-0737)

62.5 mM Tris-HCl, pH 6.8, 2% SDS, 25% glycerol, 0.01% bromophenol blue, 5% β -mercaptoethanol (added fresh)

| | |
|------------------------|----------|
| 0.5 M Tris-HCl, pH 6.8 | 3.75 ml |
| 50% Glycerol | 15.0 ml |
| 1.0% Bromophenol blue | 0.3 ml |
| 10% SDS | 6.0 ml |
| diH ₂ O | to 30 ml |

Add β -mercaptoethanol (50 μ l to 950 μ l sample buffer) before use.

2x Native PAGE (30 ml)
(catalog #161-0738)

62.5 mM Tris-HCl, pH 6.8, 40% glycerol, 0.01% bromophenol blue

| | |
|------------------------|----------|
| 0.5 M Tris-HCl, pH 6.8 | 3.75 ml |
| 50% Glycerol | 24 ml |
| 1.0% Bromophenol blue | 0.3 ml |
| diH ₂ O | to 30 ml |

2x Tricine (30 ml)
(catalog #161-0739)

200 mM Tris-HCl, pH 6.8, 2% SDS, 40% glycerol, 0.04% Coomassie Brilliant Blue G-250, 2% β -mercaptoethanol (added fresh)

| | |
|------------------------|----------|
| 1.0 M Tris-HCl, pH 6.8 | 6.0 ml |
| 100% Glycerol | 12.0 ml |
| 10% SDS | 6.0 ml |
| Coomassie Blue G-250 | 12.0 mg |
| diH ₂ O | to 30 ml |

Add β -mercaptoethanol (20 μ l to 980 μ l sample buffer) before use.

5x Nucleic acid (10 ml)
(catalog #161-0767)

50 mM Tris-HCl, pH 8.0, 5 mM EDTA, 25% glycerol, 0.2% bromophenol blue, 0.2% xylene cyanole FF

| | |
|-----------------------|----------|
| Tris base | 78.8 mg |
| 50% Glycerol | 5 ml |
| EDTA | 14.6 mg |
| 1.0% Bromophenol blue | 2.0 ml |
| Xylene cyanole FF | 20.0 mg |
| diH ₂ O | to 10 ml |

TBE-urea (30 ml)
(catalog #161-0768)
Store at 4°C

89 mM Tris, 89 mM boric acid, 2 mM EDTA, pH 8.0, 12% Ficoll, 0.01% bromophenol blue, 0.02% xylene cyanole, 7 M urea

| | |
|--------------------|----------|
| Tris base | 0.32 g |
| Boric acid | 0.165 g |
| EDTA | 17.5 mg |
| Ficoll | 3.6 g |
| Bromophenol blue | 3 mg |
| Xylene cyanole FF | 6 mg |
| Urea | 12.6 g |
| diH ₂ O | to 30 ml |

Buffer Components

0.5 M Tris-HCl, pH 6.8 (1 L)
(catalog #161-0799)
Store at 4°C

| | |
|---------------------------|---------|
| Tris base | 60.6 g |
| diH ₂ O | ~900 ml |
| Adjust to pH 6.8 with HCl | |
| diH ₂ O | to 1 L |

10% SDS (250 ml)
(catalog #161-0416)

| | |
|--------------------|-----------|
| SDS | 25.0 g |
| diH ₂ O | to 250 ml |

1.0% Bromophenol blue (10 ml)
(10 g powder, catalog #161-0404)

| | |
|--------------------|----------|
| Bromophenol blue | 100 mg |
| diH ₂ O | to 10 ml |

Blotting Buffers

Towbin buffer (1 L)

25 mM Tris, 192 mM glycine, 20% methanol

| | | |
|-----------|--------------------|--------|
| Dissolve: | Tris base | 14.4 g |
| | Glycine | 3.03 g |
| | diH ₂ O | 500 ml |

| | | |
|-----------|--------------------|--------|
| Then add: | Methanol | 200 ml |
| | diH ₂ O | to 1 L |

Alternatively, use: 10x Tris/glycine (catalog #161-0734) 100 ml
Add 200 ml methanol and diH₂O to 1 L as above

Tris-buffered saline with Tween (TTBS, 1 L)

20 mM Tris, 500 mM NaCl, 0.05% Tween 20

| | |
|--------------------|--------|
| Tris base | 2.4 g |
| NaCl | 29.2 g |
| 10% Tween 20 | 5.0 ml |
| diH ₂ O | to 1 L |

Alternatively, use: 10x TBS (catalog #170-6435) 100 ml
10% Tween 20 (catalog #166-2404) 5 ml
diH₂O 895 ml



Related Literature

| Bulletin # | Title |
|-------------------|---|
| 4006183 | Criterion™ Cell Instruction Manual |
| 4006197 | Criterion™ Dodeca™ Cell Instruction Manual |
| 10014472 | Gel Doc™ EZ System Installation Guide |
| 10019634 | Stain-Free Sample Tray Instruction Manual |
| 4006190 | Criterion Blotter Instruction Manual |
| 4006066 | Trans-Blot® SD Semi-Dry Transfer Cell Quick Reference Guide |
| 10020688 | Trans-Blot® Turbo™ Blotting System Instruction Manual |
| 1703940 | Trans-Blot SD Semi-Dry Transfer Cell Instruction Manual |
| 6007 | Criterion™ TGX™ Precast Gels Product Information Sheet |
| 5974 | Criterion™ TGX Stain-Free™ Precast Gels Product Information Sheet |
| 2911 | Criterion™ XT Precast Gels Product Information Sheet |
| 2895 | Protein Blotting Guide |
| 6039 | Trans-Blot Turbo Transfer System Brochure |
| 3133 | Molecular Weight Determination by SDS-PAGE |
| 3144 | Using Precision Plus Protein™ Standards to Determine Molecular Weight |
| 1939 | Blotting Membrane Brochure |
| 2032 | Western Blotting Detection Reagents Brochure |
| 2317 | Ready-to-Run Buffers and Solutions Brochure |
| 2414 | The Little Book of Standards |



Ordering Information

Criterion TGX™ Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) | IPG+1 Well (11 cm IPG Strip) | Prep+2 Well (800 µl/well) |
|---------|----------------------------------|--------------------------------|--------------------------------|--|-------------------------------------|
| 7.5% | 567-1023 | 567-1024 | 567-1025 | 567-1021 | 567-1022 |
| 10% | 567-1033 | 567-1034 | 567-1035 | 567-1031 | 567-1032 |
| 12% | 567-1043 | 567-1044 | 567-1045 | 567-1041 | 567-1042 |
| 18% | 567-1073 | 567-1074 | 567-1075 | 567-1071 | 567-1072 |
| 4–15% | 567-1083 | 567-1084 | 567-1085 | 567-1081 | 567-1082 |
| 4–20% | 567-1093 | 567-1094 | 567-1095 | 567-1091 | 567-1092 |
| 8–16% | 567-1103 | 567-1104 | 567-1105 | 567-1101 | 567-1102 |
| 10–20% | 567-1113 | 567-1114 | 567-1115 | 567-1111 | 567-1112 |
| Any kD™ | 567-1123 | 567-1124 | 567-1125 | 567-1121 | 567-1192 |

Criterion™ TGX Stain-Free™ Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) | IPG+1 Well (11 cm IPG Strip) | Prep+2 Well (800 µl/well) |
|--------|----------------------------------|--------------------------------|--------------------------------|--|-------------------------------------|
| 7.5% | 567-8023 | 567-8024 | 567-8025 | 567-8021 | 567-8022 |
| 10% | 567-8033 | 567-8034 | 567-8035 | 567-8031 | 567-8032 |
| 12% | 567-8043 | 567-8044 | 567-8045 | 567-8041 | 567-8042 |
| 18% | 567-8073 | 567-8074 | 567-8075 | 567-8071 | 567-8072 |
| 4–15% | 567-8083 | 567-8084 | 567-8085 | 567-8081 | 567-8082 |
| 4–20% | 567-8093 | 567-8094 | 567-8095 | 567-8091 | 567-8092 |
| 8–16% | 567-8103 | 567-8104 | 567-8105 | 567-8101 | 567-8102 |
| 10–20% | 567-8113 | 567-8114 | 567-8115 | 567-8111 | 567-8112 |
| Any kD | 567-8123 | 567-8124 | 567-8125 | 567-8121 | 567-8122 |

Criterion™ Tris-HCl Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) | IPG+1 Well (11 cm IPG Strip) | Prep+2 Well (800 µl/well) |
|------------------------------|----------------------------------|--------------------------------|--------------------------------|--|-------------------------------------|
| Criterion Tris-HCl | | | | | |
| 5% | 345-0001 | 345-0002 | 345-0003 | — | — |
| 7.5% | 345-0005 | 345-0006 | 345-0007 | — | 345-0008 |
| 10% | 345-0009 | 345-0010 | 345-0011 | 345-0101 | 345-0012 |
| 12.5% | 345-0014 | 345-0015 | 345-0016 | 345-0102 | 345-0017 |
| 15% | 345-0019 | 345-0020 | 345-0021 | — | 345-0022 |
| 18% | 345-0023 | 345-0024 | 345-0025 | — | 345-0026 |
| 4–15% | 345-0027 | 345-0028 | 345-0029 | 345-0103 | 345-0030 |
| 4–20% | 345-0032 | 345-0033 | 345-0034 | 345-0104 | 345-0035 |
| 8–16% | 345-0037 | 345-0038 | 345-0039 | 345-0105 | 345-0040 |
| 10.5–14% | 345-9949 | 345-9950 | 345-9951 | 345-0106 | — |
| 10–20% | 345-0042 | 345-0043 | 345-0044 | 345-0107 | 345-0045 |
| Criterion Stain Free™ | | | | | |
| 10% | 345-1012 | 345-1018 | — | — | — |
| 4–20% | 345-0412 | 345-0418 | 345-0426 | — | — |
| 8–16% | 345-8162 | — | 345-8166 | 345-8161 | — |

Criterion™ XT Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) | IPG+1 Well (11 cm IPG Strip) | Prep Well (800 µl/well) |
|----------------------------------|----------------------------------|--------------------------------|--------------------------------|--|-----------------------------------|
| Criterion XT Bis-Tris | | | | | |
| 10% | 345-0111 | 345-0112 | 345-0113 | 345-0115 | — |
| 12% | 345-0117 | 345-0118 | 345-0119 | 345-0121 | 345-0120 |
| 4–12% | 345-0123 | 345-0124 | 345-0125 | 345-0127 | 345-0126 |
| Criterion XT Tris-Acetate | | | | | |
| 7% | 345-0135 | 345-0136 | 345-0137 | — | — |
| 3–8% | 345-0129 | 345-0130 | 345-0131 | 345-0133 | — |

Criterion Tris-Tricine Peptide Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) | Prep+2 Well (800 µl/well) |
|--------|----------------------------------|--------------------------------|--------------------------------|-------------------------------------|
| 16.5% | 345-0063 | 345-0064 | 345-0065 | 345-0066 |
| 10–20% | 345-0067 | 345-0068 | 345-0069 | — |

Criterion IEF Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) |
|---------|----------------------------------|--------------------------------|--------------------------------|
| pH 3–10 | 345-0071 | 345-0072 | 345-0073 |
| pH 5–8 | — | 345-0076 | — |

Criterion Zymogram Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) |
|------------------------|----------------------------------|--------------------------------|--------------------------------|
| 10% Zymogram, gelatin | 345-0079 | 345-0080 | 345-0081 |
| 12.5% Zymogram, casein | 345-0082 | 345-0083 | 345-0084 |

Criterion TBE Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) |
|--------|----------------------------------|--------------------------------|--------------------------------|
| 5% TBE | 345-0047 | 345-0048 | 345-0049 |
| 10% | 345-0051 | 345-0052 | 345-0053 |
| 15% | 345-0055 | 345-0056 | 345-0057 |
| 4–20% | 345-0059 | 345-0060 | 345-0061 |

Criterion TBE-Urea Gels

| | 12+2 Well (45 µl/well) | 18-Well (30 µl/well) | 26-Well (15 µl/well) |
|-------------|----------------------------------|--------------------------------|--------------------------------|
| 5% TBE-urea | — | 345-0086 | — |
| 10% | 345-0088 | 345-0089 | 345-0090 |
| 15% | 345-0091 | 345-0092 | 345-0093 |

Catalog # Description

Criterion Gel Accessories

| | |
|----------|---|
| 345-9901 | Criterion Empty Cassettes, 1.0 mm with 12+2 well comb, 10 |
| 345-9902 | Criterion Empty Cassettes, 1.0 mm with 18-well comb, 10 |
| 345-9903 | Criterion Empty Cassettes, 1.0 mm with 26-well comb, 10 |
| 345-9904 | Criterion Empty Cassettes, 1.0 mm with prep+2 well comb, 10 |
| 345-9906 | Criterion Empty Cassettes, 1.0 mm with IPG+1 well comb, 10 |
| 165-6006 | Criterion Sample Loading Guide, 12+2 well, 1 |
| 165-6007 | Criterion Sample Loading Guide, 18-well, 1 |
| 165-6008 | Criterion Sample Loading Guide, 26-well, 1 |

Equipment

| | |
|----------|---|
| 165-6001 | Criterion Cell, includes tank, lid with power cables, three sample loading guides |
| 165-4130 | Criterion™ Dodeca™ Cell |
| 170-4070 | Criterion Blotter with Plate Electrodes |
| 170-4071 | Criterion Blotter with Wire Electrodes |
| 170-4155 | Trans-Blot® Turbo™ Transfer Starter System |
| 170-3940 | Trans-Blot® SD Semi-Dry Electrophoretic Transfer Cell |
| 164-5050 | PowerPac™ Basic Power Supply |
| 164-5052 | PowerPac HC High Current Power Supply |
| 170-8270 | Gel Doc™ EZ Imaging System |
| 170-8274 | Stain-Free™ Sample Tray |

Catalog # Description

Protein Standards

| | |
|----------|--|
| 161-0363 | Precision Plus Protein™ Unstained Standards (10–250 kD), 1 ml, 100 applications |
| 161-0373 | Precision Plus Protein All Blue Prestained Standards (10–250 kD), 500 µl, 50 applications |
| 161-0374 | Precision Plus Protein Dual Color Standards (10–250 kD), 500 µl, 50 applications |
| 161-0375 | Precision Plus Protein™ Kaleidoscope™ Standards (10–250 kD), 500 µl, 50 applications |
| 161-0376 | Precision Plus Protein™ WesternC™ Standards (10–250 kD), 250 µl, 50 applications |
| 161-0377 | Precision Plus Protein Dual Xtra Standards (2.5–250 kD), 500 µl, 50 applications |
| 161-0385 | Precision Plus Protein WesternC Pack (10–250 kD), 50 applications each of standard and StrepTactin-HRP |
| 161-0317 | SDS-PAGE Standards, broad range, 200 µl |

Premixed Running Buffers

| | |
|----------|--|
| 161-0732 | 10x Tris/Glycine/SDS, 1 L |
| 161-0772 | 10x Tris/Glycine/SDS, 5 L |
| 161-0734 | 10x Tris/Glycine, 1 L |
| 161-0744 | 10x Tris/Tricine/SDS, 1 L |
| 161-0788 | XT MOPS Running Buffer, 20x, 500 ml |
| 161-0789 | XT MES Running Buffer, 20x, 500 ml |
| 161-0790 | XT Tricine Running Buffer, 20x, 500 ml |
| 161-0793 | XT MOPS Buffer Kit, includes 500 ml 20x XT MOPS running buffer, 10 ml 4x XT sample buffer, 1 ml 20x XT reducing agent |
| 161-0796 | XT MES Buffer Kit, includes 500 ml 20x XT MOPS running buffer, 10 ml 4x XT sample buffer, 1 ml 20 x XT reducing agent |
| 161-0797 | XT Tricine Buffer Kit, includes 500 ml 20x XT MOPS running buffer, 10 ml 4x XT sample buffer, 1 ml 20x XT reducing agent |
| 161-0761 | 10x IEF Anode Buffer, 250 ml |
| 161-0762 | 10x IEF Cathode Buffer, 250 ml |
| 161-0733 | 10x Tris/Boric Acid/EDTA, 1 L |
| 161-0770 | 10x Tris/Boric Acid/EDTA, 5 L |
| 161-0765 | Zymogram Renaturation Buffer, 125 ml |
| 161-0766 | Zymogram Development Buffer, 125 ml |

■ Criterion Precast Gels

| Catalog # | Description |
|--|---|
| Premixed Sample Buffers | |
| 161-0737 | Laemmli Sample Buffer, 30 ml ¹ |
| 161-0738 | Native Sample Buffer, 30 ml |
| 161-0791 | XT Sample Buffer, 4x, 10 ml |
| 161-0792 | XT Reducing Agent, 20x, 1 ml |
| 161-0739 | Tricine Sample Buffer, 30 ml |
| 161-0763 | IEF Sample Buffer, 30 ml |
| 161-0764 | Zymogram Sample Buffer, 30 ml |
| 161-0767 | Nucleic Acid Sample Buffer, 5x, 10 ml |
| 161-0768 | TBE-Urea Sample Buffer, 30 ml |
| Component Reagents | |
| 161-0719 | Tris, 1 kg |
| 161-0718 | Glycine, 1 kg |
| 161-0301 | SDS, 100 g |
| 161-0416 | SDS Solution, 10% (w/v), 250 ml |
| 166-2404 | 10% Tween 20, 5 ml |
| 170-6404 | Blotting-Grade Blocker, 300 g |
| 161-0710 | 2-Mercaptoethanol, 25 ml |
| 161-0611 | Dithiothreitol (DTT), 5 g |
| 161-0404 | Bromophenol Blue, 10 g |
| Total Protein Gel and Blot Stains | |
| 161-0786 | Bio-Safe™ Coomassie Stain, 1 L |
| 161-0400 | Coomassie Brilliant Blue R-250, 10 g |
| 161-0436 | Coomassie Blue R-250 Stain Solution, 1 L |
| 161-0438 | Coomassie Blue R-250 Destain Solution, 1 L |
| 161-0443 | Silver Stain Kit |
| 161-0449 | Silver Stain Plus™ Kit |
| 170-6527 | Colloidal Gold Total Protein Stain, 500 ml |
| 161-0440 | Zinc Stain and Destain Kit |
| 170-3127 | SYPRO Ruby Protein Blot Stain, 200 ml |
| 161-0491 | Flamingo™ Fluorescent Gel Stain (10x), 100 ml |
| 161-0496 | Oriole™ Fluorescent Protein Gel Stain (1x), 1 L |

¹ May require addition of 2-mercaptoethanol or DTT

| Catalog # | Description |
|--------------------------------------|--|
| Immunoblot Detection Reagents | |
| 161-0385 | Precision Plus Protein™ WesternC™ Pack |
| 170-5070 | Immun-Star™ WesternC™ Chemiluminescent Kit, 100 ml |
| 170-6431 | HRP Conjugate Substrate Kit, 4CN |
| 170-6535 | HRP Color Development Reagent, DAB, 5 g |
| 170-8238 | Amplified Opti-4CN™ Substrate Kit |
| 170-8235 | Opti-4CN Substrate Kit |
| 170-6432 | AP Conjugate Substrate Kit |
| 170-5012 | Immun-Star™ Substrate Pack |

Blotting Membranes

| | |
|----------|--|
| 162-0232 | 0.2 µm Nitrocellulose/Filter Paper Sandwich, 8.5 x 13.5 cm, 20 pack |
| 162-0233 | 0.2 µm Nitrocellulose/Filter Paper Sandwich, 8.5 x 13.5 cm, 50 pack |
| 162-0234 | 0.45 µm Nitrocellulose/Filter Paper Sandwich, 8.5 x 13.5 cm, 20 pack |
| 162-0235 | 0.45 µm Nitrocellulose/Filter Paper Sandwich, 8.5 x 13.5 cm, 50 pack |
| 162-0236 | Sequi-Blot™ PVDF/Filter Paper Sandwich, 8.5 x 13.5 cm, 20 pack |
| 162-0237 | Sequi-Blot PVDF/Filter Paper Sandwich, 8.5 x 13.5 cm, 50 pack |
| 170-4157 | Trans-Blot Turbo Midi PVDF Transfer Packs |
| 170-4159 | Trans-Blot Turbo Midi Nitrocellulose Transfer Packs |

For additional product sizes, please refer to the Bio-Rad catalog or website.



P4110001



**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site www.bio-rad.com **USA** 800 424 6723 **Australia** 61 2 9914 2800 **Austria** 01 877 89 01 **Belgium** 09 385 55 11 **Brazil** 55 31 3689 6600
Canada 905 364 3435 **China** 86 21 6169 8500 **Czech Republic** 420 241 430 532 **Denmark** 44 52 10 00 **Finland** 09 804 22 00
France 01 47 95 69 65 **Germany** 089 31 884 0 **Greece** 30 210 777 4396 **Hong Kong** 852 2789 3300 **Hungary** 36 1 459 6100 **India** 91 124 4029300
Israel 03 963 6050 **Italy** 39 02 216091 **Japan** 03 6361 7000 **Korea** 82 2 3473 4460 **Malaysia** 60 3 2117 5260 **Mexico** 52 555 488 7670
The Netherlands 0318 540666 **New Zealand** 64 9 415 2280 **Norway** 23 38 41 30 **Poland** 48 22 331 99 99 **Portugal** 351 21 472 7700
Russia 7 495 721 14 04 **Singapore** 65 6415 3170 **South Africa** 27 861 246 723 **Spain** 34 91 590 5200 **Sweden** 08 555 12700
Switzerland 061 717 95 55 **Taiwan** 886 2 2578 7189 **Thailand** 66 2 6518311 **United Kingdom** 020 8328 2000