

TaqMan™ SCID/SMA Assay

Catalog Numbers A47927, A47928, A47929

Pub. No. MAN0018922 Rev. C.0

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *TaqMan™ SCID/SMA Assay User Guide* (Pub. No. MAN0018921). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

This Quick Reference is intended as a benchtop reference for experienced users of the TaqMan™ SCID/SMA Assay. For detailed instructions and ordering information for additional products, see the *TaqMan™ SCID/SMA Assay User Guide* (Pub. No. MAN0018921).

Contents and storage

Table 1 TaqMan™ SCID/SMA Assay

| Cat. No. | Number of reactions | Storage |
|----------|---------------------|--|
| A47927 | 4,000 | <ul style="list-style-type: none"> • 2–8°C for up to 3 months • –30°C to –10°C for long-term storage |
| A47928 | 8,000 | |
| A47929 | 20,000 | |

Procedural guidelines

- Keep the PCR Reaction Mix on a chilled block during real-time PCR reaction setup.
- The PCR Reaction Mix can be formulated up to 48 hours before sample preparation and stored at 2–8°C. See “Prepare the real-time PCR reactions” on page 2.
- Ensure that the instrument is calibrated for each detector dye and passive reference dye, according to the instrument user guide. For more information, see the *TaqMan™ SCID/SMA Assay Calibration Guide* (Pub. No. MAN0019378).

Before you begin

DBS wash buffer and stabilizing buffer are stable at room temperature and can be prepared several days in advance. They do not need to be prepared fresh.

- Prepare DBS wash buffer—add 5 g of liquid Thesit™ reagent into 100 mL of luke warm 1X PBS to prepare a 5% stock, then perform serial dilution using luke warm 1X PBS to prepare 0.5% DBS wash buffer.

DBS wash buffer is 1X PBS and 0.5% (w/v) Thesit™ reagent.

Note:

- Alternatively, a higher or lower concentration and/or volume of stock can be prepared based on the throughput of your laboratory.
- Heat the Thesit™ reagent until liquid and the 1X PBS until luke warm to allow the Thesit™ reagent to dissolve readily.

- Prepare stabilizing buffer—add 0.5 mL of 10% Tween™-20 Surfact-Amps™ Detergent Solution to 100 mL of TE Buffer. Stabilizing buffer is modified TE buffer (10 mM Tris, 0.1 mM EDTA, pH 8) with 0.05% Tween™-20 detergent added.

Extract the DNA

1. Place a DBS punch (1.5 mm or 3.2 mm) from a blood card into the well of the 0.2 mL optical 96-well reaction plate.

Note: Use only one punch per well.

2. Wash the DBS punch.
 - a. Add the DBS wash buffer to each sample well.
 - 1.5 mm DBS punch: 70 μ L
 - 3.2 mm DBS punch: 150 μ L
 - b. Seal the plate with a clear adhesive film, then centrifuge at 2400 rpm for 1 minute to wet and submerge the punch.
 - c. Shake the plate for 30 minutes at 1500 rpm on a microplate shaker.
 - d. Centrifuge briefly to collect the contents at the bottom of the wells, then remove and discard the supernatant.
 - e. Add 150 μ L of water to each sample well, then remove and discard the water.
3. Add 5 μ L of DNA Extract All Reagents Kit Lysis Solution to each sample well.
4. Seal the plate, then centrifuge briefly to collect the contents at the bottom of the wells.
5. Incubate the reaction plate for 5 minutes at 95°C.
6. Centrifuge briefly to collect any condensed droplets at the bottom of the wells, then cool to room temperature before opening the seal.

Note: Optionally, place the plate on a chilled block to quickly cool the plate.
7. Add 35 μ L of stabilizing buffer to each sample well.

Note: The total lysate volume per 1.5 mm or 3.2 mm DBS punch is 40 μ L.
8. Seal the plate, vortex briefly, then centrifuge briefly to collect the contents at the bottom of the wells.
9. Transfer the lysate to a fresh optical reaction plate according to the following table.

| DBS punch | Lysate volume per well | Reaction plate |
|-----------|------------------------|----------------------------|
| 1.5 mm | 9 μ L | 96-well (0.1 mL or 0.2 mL) |
| 3.2 mm | 9 μ L | 96-well (0.1 mL or 0.2 mL) |
| 3.2 mm | 6.75 μ L | 384-well |

Prepare the real-time PCR reactions

1. Prepare the PCR Reaction Mix in an appropriately sized microcentrifuge tube according to the following table.

| Component | Volume per reaction ^[1] | |
|--|------------------------------------|-------------------------------|
| | 96-well plate | 384-well plate |
| 2X TaqPath™ ProAmp™ Multiplex Master Mix | 10 μ L | 7.5 μ L |
| 20X TaqMan™ SCID/SMA Assay | 1 μ L | 0.75 μ L |
| Total PCR Reaction Mix volume | 11 μL | 8.25 μL |

^[1] Add 10% overage for pipetting loss. Excess PCR Reaction Mix can be stored at 2–8°C for up to 48 hours.

2. Vortex to mix, then centrifuge briefly to collect the contents at the bottom of the tube.
3. Add the PCR Reaction Mix to each well containing the transferred lysate.
 - 96-well plate: 11 μ L
 - 384-well plate: 8.25 μ L
4. Seal the plate with an optical adhesive film, vortex to mix, then centrifuge briefly to collect the contents at the bottom of the wells.

Set up and run the real-time PCR instrument

See your instrument user guide for detailed instructions to program the thermal cycling conditions or to run the plate.

Note: The instrument must be configured with the block appropriate for the plate type.

1. Select or create dye detectors, then assign them to each well in the plate layout.

| Target | Reporter | Quencher |
|---------|----------|----------|
| SMN1 | FAM™ dye | MGB-NFQ |
| TREC | JUN™ dye | QSY™ |
| RNase P | VIC™ dye | QSY™ |

2. Set up the following conditions.
 - Cycling mode: Fast
 - Passive reference: MUSTANG PURPLE™ dye
3. Set the reaction volume.
 - Reaction volume for a 1.5 mm DBS punch in a 96-well plate: 20 µL
 - Reaction volume for a 3.2 mm DBS punch in a 96-well plate: 20 µL
 - Reaction volume for a 3.2 mm DBS punch in a 384-well plate: 15 µL
4. Set up the following thermal cycling protocol.

| Step | Temperature | Time | Cycles |
|--|-------------|------------|--------|
| Initial denaturation / Enzyme activation | 95°C | 20 seconds | 1 |
| Denature | 95°C | 1 second | 40 |
| Anneal / Extend | 60°C | 20 seconds | |

5. Load the plate into the real-time PCR instrument.
6. Start the run.

Analyze the results

See your instrument user guide for more detailed information about data analysis.

1. View the amplification plots.
2. Set the baseline and threshold values.
3. Calculate the C_t values on the instrument.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0018922

| Revision | Date | Description |
|----------|----------------|---|
| C.0 | 20 August 2020 | Updated wash steps. |
| B.0 | 10 August 2020 | <ul style="list-style-type: none">Added references to the <i>TaqMan™ SCID/SMA Assay User Guide</i> and the <i>TaqMan™ SCID/SMA Assay Calibration Guide</i>.Removed required materials not supplied.Updated quencher and disclaimer. |
| A.0 | 6 April 2020 | New document. |

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. Thesit is a trademark of Desitin Arzneimittel GmbH.