
TwisTorr 74 FS PCB 24 Vdc Controller

**Model
X3510-64050**

**Manuale di istruzioni
Bedienungshandbuch
Notice de mode d'emploi
Manual de instrucciones
Manual de instruções
Bedrijfshandleiding
Instruktionsbog
Bruksanvisning**

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Notices

© Agilent Technologies, Inc. 2018

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

Publication Number: 87-901-056-01

Edition

Edition 07/2018

Printed in ITALY

Agilent Technologies Italia S.p.A.

Vacuum Products Division

Via F.lli Varian, 54

10040 Leini (TO)

ITALY

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or

contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Trademarks

Windows and MS Windows are U.S. registered trademarks of Microsoft Corporation.

Safety Notices


CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

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Contents

Data Exchange Protocols	184
Example	194
Operation	199
Maintenance	201
Accessories and Spare Parts	202

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
Odlaganje odpadkov

16 Instructions for Use

General Information	164
Storage	165
Preparation prior to Installation	166
Installation	168
Use	169
Operating Procedures	169
Powering on the Controller	169
Starting the Pump	169
Stopping the Pump	169
Maintenance	170
Disposal	171

Original Instructions

General Information

This equipment is intended for professional use. Prior to using the equipment, the operator must carefully read this instruction manual and any other additional information provided by Agilent. Agilent declines all responsibility which may become as the result of the total or partial misinterpretation of the information provided in the documentation, by the improper use of the equipment by untrained personnel or by the use of the equipment by unauthorized personnel or not in compliance with specific national norms.

The TwisTorr 74 FS PCB controller is a microprocessor-controlled frequency converter made up of solid state components and capable of autodiagnostic and auto-protection functionality. It incorporates the circuitry needed for the automatic operation of TwisTorr 74 FS series pumps.

The TwisTorr 74 FS series pumps are driven by the controller (in a 10-step procedure) during the start-up phase that also checks the voltage and current according to the speed reached by the pump. Auxiliary connectors render the commands needed to start and stop the pump from remote, the signals that indicate the operating state of the pump and the locking commands (for pressure switches, water flow control switches, etc.) available. The controller must be provided with 24 Vdc.

In the next sections you will find all the information needed to guarantee operator safety during the use of the equipment. Detailed information is provided in the Appendix entitled "Technical Information".

The following conventions are adopted in this manual:

CAUTION!

Caution messages are provided before a procedure that could damage the equipment is not executed in the correct way.

WARNING!

Warning messages call the operator's attention to a specific procedure or operation that could cause serious personal injury if not executed in the correct way.

NOTE

Notes provide important information extracted from the text.

Storage

The controller must be transported or stored within the following environmental conditions:

- Temperature: -20 °C to +70 °C
- Relative humidity 0 to 95 % (non-condensing)

Preparation prior to Installation

The controller comes in a special protective packaging; if you notice any sign of damage that could have been caused during transportation, notify our local sales office immediately.


Figure 1 Controller Packaging

During the unpacking procedure, be careful to avoid dropping the controller or hitting it against anything. Do not discard the packaging material in the environment. The packaging material is completely recyclable and complies with the EEC environmental protection directive 85/399. Before installing the controller you need to prepare certain interconnection cables such as, for example, the power supply interconnection cables, the logic interconnection cables and the pump interconnection cable. The following figures provide the connection diagrams between the controller and pump and between the controller and power supply. For further information on these interconnections and for the technical characteristics of the transformer (available as option), refer to the Appendix entitled "Technical Information".


Figure 2 Controller – Pump Interconnection Cable

Installation

WARNING!

The controller develops internal voltages that could cause serious injury. Unplug the controller from the electrical wall outlet before beginning with any type of installation or maintenance procedure. In addition, the controller develops internal temperatures that could also cause serious injury. Adequately protect the controller in its final point of installation so as to avoid accidental contacts.

NOTE

When definitively installed in the system, the controller must be positioned so that the cooling air can circulate freely around the equipment. Do not install and/or use the controller in environments exposed to the atmosphere (rain, ice, and snow), to dust, to aggressive gases and in explosive areas where there is a high risk of fire.

The controller must be used within the following environmental conditions:

- Temperature: 0 °C to +40 °C
- Relative humidity: 0 to 95 % (non-condensing).

Use

This section provides all the major operating procedures. For further details and for the procedures that involve connections or specific options, refer to the section entitled "USE" in the Appendix "Technical Information".

Perform all electrical and pneumatic connections before using the controller, with reference to the manual of the connected pump.

WARNING!

Refer to the relevant pump manual to operate it in a safe way before starting.

Operating Procedures

Powering on the Controller

To power on the controller, simply provide electrical power to it.

Starting the Pump

To start the pump, leave pin 12 of connector J4 open so that the pump starts-up when 24 Vdc is provided to the controller.

Stopping the Pump

To stop the pump, short-circuit pins 12 and 15 of connector J4.

Maintenance

The TwisTorr 74 FS series controllers do not require any form of maintenance. Any intervention must be performed by authorized personnel.

In case of failure you can use the Agilent customer support service or the "Agilent advance exchange service", the latter consisting of replacing a faulty controller with a refurbished one.

WARNING!

Unplug the power cord from the electrical socket before performing any type of intervention on the controller.

If a controller needs to be scraped, proceed to do so in compliance with the specific national laws.

Disposal

Meaning of the "WEEE" logo found in labels. The following symbol is applied in accordance with the EC WEEE (Waste Electrical and Electronic Equipment) Directive. This symbol (**valid only in countries of the European Community**) indicates that the product it applies to must NOT be disposed of together with ordinary domestic or industrial waste but must be sent to a differentiated waste collection system. The end user is therefore invited to contact the supplier of the device, whether the Parent Company or a retailer, to initiate the collection and disposal process after checking the contractual terms and conditions of sale


For more information refer to:

<http://www.agilent.com/environment/product/index.shtml>

16 Instructions for Use

Disposal

17

Technical Information

General Description	174
TwisTorr 74 FS PCB 24 Vdc Controller Description	174
Controller Specifications	175
Controller Connector and Jumper	176
Jumper Configurations	177
Run-up Time Adjustment	177
Controller Outline	178
Installation	179
Data Exchange Protocols	184
RS-232 Protocol 1 Communication Description	184
Communication Protocol	185
Error Code Table	190
RS-232 Protocol 2 Communication Descriptions	190
Transmission Channel Characteristics	190
Message Structure	191
Example	194
Operation	199
Low Speed Activation/Deactivation	200
Maintenance	201
Accessories and Spare Parts	202

General Description

The TwisTorr 74 FS PCB controller is a microprocessor-controlled, solid-state, frequency converter with self-diagnosis and protection features.

TwisTorr 74 FS PCB 24 Vdc Controller Description

The controller is driven by a single chip microcomputer and is composed of a PCB that includes all the circuitry necessary for its operation. The microcomputer generates the variable output voltage according to the software and the gas load condition of the pump. Moreover, it manages signals from sensors, input/output connection information, and gives output for a fully automatic operation. The controller can be operated via remote signals through an RS-232 connection. The controller can be operated in local mode through suitable switches connected between the input pins of J4 connector and ground (see section J4 connector for a detailed description of the input signal present on J4 connector).

Controller Specifications

Tab. 1

Input:	
- Voltage	20.4 to 32.2 Vdc with 2 Vpp max ripple
- Power	128 W max.
Output:	
- Voltage	42 Vac nominal $\pm 10\%$, 3-phase
- Frequency	1167 Hz, $\pm 2\%$
- Power	70 W maximum
Fuse	T 4A
Operating temperature	0 °C to + 40 °C
Storage temperature	-20 °C to + 70 °C
Weight	0.5 Kg (1.1 lbs)

Controller Connector and Jumper

The following figure shows the controller and indicates its connectors and jumpers.


Figure 3 Connector and Jumper

Controller connector are named with prefix Px; in the following pages the connector counterpart are named as Jx to indicate the use of the appropriate connector housing and pins.

Jumper Configurations

The following figures show the configuration and meaning of the individual jumpers on the board.

Jumper W7 is used to select the data exchange rate, as shown in the figure.


Figure 4 W7 Jumper

Jumper W6 is used to set the Soft Start, as shown in the figure.

Figure 5 W6 Jumper

The W5 jumper is just for factory setting and must not be moved.

Run-up Time Adjustment

TR1 can be used to change the run-up time between 8 and 60 min. The controller is factory set to 8 min. If the gas load is large and the pump cannot reach the nominal speed in the set time, the pump is stopped and the controller gives Error 6 (Normal not reached at run-up time). In this case, the run-up time setting can be increased up to 60 min. by rotating TR1. The run-up time can be read through the serial communication at parameter 00700 (if protocol 1 is used) or at window 127 (if protocol 2 is used).

Controller Outline

The outline dimension for the controller is shown in the following figure.


Figure 6 Controller Outline

Installation

Inspect the controller for any shipping damage. The controller is designed to be installed into a rack. All the connections are fitted to J2, J3, J4 and J5 connectors. Should the controller be connected to a host computer via the RS-232 interface, a suitable cable must be prepared. In the following paragraphs are detailed the input/output signals.

NOTE

The PCB installed into the customer system must be positioned so that cold air (forced or natural convection) can flow through the PCB components.

J3 Connector


Figure 7 J3 – Power Connector

The signals of J3 connector are the following:

- **Pins 1/2** Power supply 24. Vdc (20 to 32 Vdc) (pin 2 positive).
- **Pins 4/3** K1 relay contacts. The relay switches when the output frequency overrides the set threshold (1117 Hz). The threshold overriding is indicated also on pin 6 of J4 connector.

J2 Connector


Figure 8 J2 – Connector to pump

The signals of J2 connector are the following:

- **Pin 1** 42 Vac 3-phase output to pump motor stator (phase T).
- **Pin 2** 42 Vac 3-phase output to. pump motor stator (phase S).
- **Pin 3** 42 Vac 3-phase output to pump motor stator (phase R).
- **Pins 5/4** Pump temperature sensor.

J4 Connector


Figure 9 J4 – Signal Connector

The signals of J4 connector are the following:

- **Pin 3** 24V 60mA open collector Output failure. When a fault condition is detected, the pump is stopped and this output is ON.
- **Pin 6** 24V 60mA open collector Output When the output frequency is less than 1117 Hz, this output is ON.
- **Pin 7** Analog output. A voltage from 0 to 10 Vdc proportional to the output pump supplied current (from 0 to 2.5 A).
- **Pin 8** Analog output. A voltage from 0 to 10 Vdc proportional to the output pump motor driving frequency (from 0 to 1167 Hz).

17 Technical Information

Installation

- **Pin 11** Error reset command. Should an operation fault be detected (Pin 3 of J4 is ON), the pump is stopped. To reset the fault you must connect pin 11 to ground (pin 15 of J4) for at least one sec.
- **Pin 12** ON/OFF pump. To stop the pump you must connect pin 12 to ground (pin 15 of J4); to start the pump you must remove the contact between pin 12 and ground.
- **Pin 15** Ground

The following figure shows a typical contact logic input connection and the related simplified circuit of the controller.


Figure 10 Typical Logic Input

The following figure shows a typical logic output connection, but any other device may be connected, and the related simplified circuit of the controller.


Figure 11 Typical Logic Output

J5 RS-232 Interface Connector


Figure 12 J2 – Interface Connector

The signal of J5, the RS-232 interface connector, are the following:

- **Pin 1** Ground.
- **Pin 2** RXD. Received data.
- **Pin 3** TXD. Transmitted data.
- **Pin 4** 24 Vdc Output (power supply to HHT, if necessary)

Data Exchange Protocols

This controller can manage two different serial communication protocols (1 or 2) on its RS-232 port. During the first communication frame after power on, the controller automatically recognizes the type of protocol used by the host to which it will interface.

NOTE

The identified protocol is kept active until the controller is powered off.

RS-232 Protocol 1 Communication Description

Transmission Channel Characteristics

Tab. 2

levels:	RS-232
baud rate:	9600/4800 programmable by W7 jumper on the board
character length:	8 bits
parity:	none
stop bit:	1 bit
protocol:	master (PC) / slave (converter)

Communication Protocol

The communication protocol is master-slave type, where:

- Host = Master
- Controller = slave

NOTE

To write data through the serial interface, the controller must be set to "REMOTE OFF" operation, otherwise only read data command will be available.

The exchanged frame format is composed of ASCII values. The frame composition is:

slave address			Parameter number					data field length	data						check sum			CR	

where:

- **Slave Address:** made of three characters containing the controller address requested by the master. In RS-232 version the address is "001";
- **Parameter Number:** made of five characters, it points out to the controller about which parameter the master wants to operate, and about the kind of operation (read/write). It consists of:

1	2	3	4	5
0= parameter reading 1 = parameter writing	always 0	000 to 299 = commands 300 to 699 = status 700 to 999 = parameter		

- **Data Field Length:** made of 2 characters, it points out the next field length in characters. In this protocol the possible combinations are:

02 = parameter reading request

06 = parameter, writing
- **Data:** variable length, it could be 2 or 6 characters according to the previous field content, it can be:
"=?": in the case of parameter reading request
"000000" or "111111": in the case of YES/NO command or parameter writing
a six numeric character string 0 to 9
pointing out the value of a configuration parameter or an analog measure
a six character string pointing out an error
- **Checksum:** made of 3 characters, it is the sum of the string ASCII values up to the first character of the checksum with a 8 bit cutting-off converted to decimal number.
- **CR:** made of 1 character < 0xD>

All the fields are aligned to the right with the "0" character used as pad character.

When the controller receives a frame with the correct address value (001 in this case) it answers as detailed in the following:

- for a parameter request: a frame having the same format described before, but the first character of the "Parameter Number" field set to 1.

Example for a low speed status request:

Master to Controller:
0010000202= ?097< 0xD >

Controller to Master:
0011000206111111016< 0xD >
where "111111" is for low speed status ON

- for a parameter writing: a frame including the parameter value that has been set out.

Example for a low speed status ON writing:

Master to Controller:
0011000206111111016< 0xD >

Controller to Master:
0011000206111111016< 0xD >

If the controller receives wrong frames it answers as following:

- in the case of frame with length greater than 40 characters:

			N	A	K	CR
Controller address						

- in case of frame with wrong checksum:

			N	A	K	CR
Controller address						

- in case of frame containing parameter number not present:

			1	0				0	6	N	O	-	D	E	F				
control. address				param. address												check- sum		CR	

- in case of writing values out of limits:

			1	0				0	6	-	R	A	N	G	E				
control. address				param. address												check- sum		CR	

- in the case of an attempt of writing a logic parameter with a value not equal to "000000" or "111111", or when the parameter is a read only one:

			1	0				0	6	-	L	0	G	I	C				
control. address						param. address										check- sum		CR	

- in case of an excessive time between two characters (greater than 1 second):

			N	A	K	CR
Controller address						

Information available on the serial interface are listed in the following table where the column "TYPE" points out a logic variable with "D" (allowed values "000000" or "111111") and a numeric variable with "N".

Variables with the first character of "Parameter Number" equal to X can be either written or read. Variables with that character equal to "0" can be only read, equal to "1" can be only written.

Tab. 3

PARAMETER NUMBER	TYPE	DESCRIPTION
X0002	D	Low speed ON/OFF
X0003	D	Pump ON/OFF
X0008	D	Remote ON/OFF
10009	D	Error acknowledgment: error status clearing
00302	D	K1 status
00303	N	Error code: error code displaying (see the following table)

PARAMETER NUMBER	TYPE	DESCRIPTION
00306	D	Normal YES/NO
00307	D	Starting YES/NO
00309	N	Output frequency (Hz)
00310	N	Pump current * 100 (A)
00311	N	Temperature (°C)
00312	N	Software version The controller sends "QExxxx" where "xxxx" is the CRC16 of EPROM
00313	N	Power (W)
00314	N	Pump life (hours)
00315	N	Number of cycles
00316	N	Last cycle time (minutes)
00317	N	Software version. The controller sends "PAxxxx" where "xxxx" is the CRC16 of EEPROM
00700	N	Run up time (minutes)
X0701	N	Threshold (Hz)
X0702	D	Soft-start YES/NO
X0703	D	Reduced output power (YES) Full output power (NO)

Error Code Table

Tab. 4

CODE	DESCRIPTION
00E000	No error
00E001	Output current > 15 A
00E002	Pump not connected
00E003	Output current > 1.5 A for more than 15 sec.
00E004	Bearing temperature > 60 °C
00E005	Heat exchanger temperature > 60 °C
00E006	Normal not reached within run up time
00E009	Vp voltage 8V lower than nominal for more than 5 sec.
00E010	Soft start ramp not ended within the expected time

RS-232 Protocol 2 Communication Descriptions

Transmission Channel Characteristics

Tab. 5

levels:	RS-232
baud rate:	9600/4800 programmable by W7 jumper on the board
character length:	8 bits
parity:	none
stop bit:	1 bit
protocol:	master (PC) / slave (converter)

Message Structure

The communication protocol is a MASTER/SLAVE type where:

- Host = MASTER
- Controller = SLAVE

The communication is performed in the following way:

1. the host (MASTER) send a MESSAGE to the controller (SLAVE);
2. the controller answer with a MESSAGE to the host.

The MESSAGE is a string with the following format:

<STX>+<ADDR>+<WIN>+<COM>+<DATA>+<ETX>+<CRC>

where:

NOTE

When a data is indicated between two quotes ('...') it means that the indicated data is the corresponding ASCII character.

-
- <STX> (Start of transmission) = one character containing 0x02
 - <ADDR> (Unit address) = one character containing 0x80 (for RS-232)
<ADDR> (Unit address) = one character containing 0x80 + device number (0 to 31) (for RS-485)
 - <WIN> (Window) = a string of 3 numeric characters indicating the window number (from '000' to '999'); for the meaning of each window see the relevant paragraph.
 - <COM> (Command) = one character containing 0x30 to read the window, 0x31 to write into the window.

- For host messages this field is present only in case of writing containing logic, numeric or alphanumeric data types according to the following table.
- For controller messages this field is always present containing the required data (logic, numeric or alphanumeric) in case of reading or the relevant protocol answer in case of writing acknowledge or protocol errors.

Tab. 6

	Length	Characters Permitted
Logic (L)	1	'0' = off '1' = on
Numeric (N)	6	'-', '.', '0' ... '9' right justified with '0'
Alphanumeric (A)	10	from blank to ' _ ' (ASCII)

- <ETX> (End of transmission) = one character containing 0x03
- <CRC> = two characters containing XOR of all characters subsequent to <STX> and including the <ETX> terminator. The value is hexadecimal coded and indicated by two ASCII character.

Only the addressed SLAVE will respond.

NOTE

Using the RS-485 interface, the message structure remains identical to the one used for the RS-232 interface, the only difference being that the value assigned to the ADDRESS <ADDR>

The controller can answers with the following response types:

Tab. 7 Response type

Type	Length	Value	Description
Logic	1 byte	-	Logic value
Numeric	6 bytes	-	Numeric value
Alphanumeric	10 bytes	-	Alphanumeric value
ACK	1 byte	(0x6)	The command execution has been successfully completed
NACK	1 byte	(0x15)	The command execu-tion has been failed
Unknown Window	1 byte	(0x32)	The specified window in the command is not a valid window
Data Type Error	1 byte	(0x33)	The data type specified in the command (Logic, Numeric or Alphanumeric) is not accorded with the specified Window
Out of Range	1 byte	(0x34)	The value expressed during a write command is out of the range value of the specified window
Win Disabled	1 byte	(0x35)	The specified window is Read Only or temporarily disabled (for example you can't write the Soft Start when the Pump is running)

Example

Command: START

Source: PC

Destination: Inverter

02	80	30	30	30	31	31	03	42	33
STX	ADDR	WINDOW			WR	ON	ETX	CRC	

Source: Inverter

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Command: STOP

Source: PC

Destination: Inverter

02	80	30	30	30	31	30	03	42	32
STX	ADDR	WINDOW			WR	OFF	ETX	CRC	

Source: Inverter

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Command: SOFT-START (ON)

Source: PC

Destination: Inverter

02	80	31	30	30	31	31	03	42	32
STX	ADDR	WINDOW			WR	ON	ETX	CRC	

Source: Inverter

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Command: SOFT-START (OFF)

Source: PC

Destination: Inverter

02	80	31	30	30	31	30	03	42	33
STX	ADDR	WINDOW			WR		OFF	ETX	CRC

Source: Inverter

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

17 Technical Information

Example

Command: CURRENT

Source: PC

Destination: Inverter

02	80	32	30	30	30	03	38	31
STX	ADDR	WINDOW			RD	ETX	CRC	

Source: Inverter

Destination: PC

02	80	32	30	30	30	30	30	30	2E	30	30	03	39	44
STX	ADD	WINDOW			RD	000.00						ETX	CRC	

Command: FREQUENCY

Source: PC

Destination: Inverter

02	80	32	30	33	30	03	38	32
STX	ADDR	WINDOW			RD	ETX	CRC	

Source: Inverter

Destination: PC

02	80	32	30	33	30	30	30	30	30	34	32	03	38	34
STX	ADD	WINDOW			RD	000042						ETX	CRC	

Command: ERR-CODE

Source: PC

Destination: Inverter

02	80	32	30	36	30	03	38	37
STX	ADDR	WINDOW			RD	ETX	CRC	

Source: Inverter

Destination: PC

02	80	32	30	36	30	30	30	30	30	30	30	03	38	37
STX	ADD	WINDOW			RD	000000						ETX	CRC	

Serial Communication Windows

Tab. 8

WIN	TYPE	R	W	Description
000	L	X	X	START/STOP (1= START ; 0= STOP)
008	L	X	X	REMOTE/SERIAL Configuration (1= Remote ; 0= Serial)
100	L	X	X	SOFT START YES/NO (1= YES ; 0= NO) Default = 0
101	N	X	X	SET POINT TYPE (0 =Frequency; 1=Current; 2=Time) Default = 0
102	N	X	X	SET POINT VALUE [Hz, Ma, SEC] Default = 1000 Hz
103	N	X	X	SET POINT MASK TIME [s] Default = 0
104	L	X	X	OUTPUT STATE (0= OFF if Val > threshold; 1 = ON if Val > threshold) Default = 0
105	N	X	X	HYSTERESIS IN THRESHOLD % [0-100] Default = 2
107	L	X	X	ACTIVE STOP (0=NO; 1=YES) Default = 1
108	N	X	X	BAUD RATE (3-4) [4800-9600] Default = 4
109	L		X	PUMP LIFE RESET [Write "1" to Reset]
111	N	X	X	SET ANALOG OUTPUT TYPE [0=Power; 1=Current] Default = 0
120	N	X	X	SET ROTATIONAL FREQUENCY [Hz] $1100 \text{ Hz} \leq F_{\text{imp}} \leq F_{\text{MAX}}$
121	N	X	X	MAX SETTABLE ROTATIONAL FREQUENCY [Hz] $F < 1167 \text{ Hz}$
127	N	X		RUN UP TIME (8 – 60 minute)
130	N	X		RAMP CURRENT [mA]
200	N	X		CURRENT [mA]
201	N	X		VOLTAGE [V]
202	N	X		POWER [W]
203	N	X		DRIVING FREQUENCY [Hz]

17 Technical Information

Example

WIN	TYPE	R	W	Description
204	N	X		PUMP TEMPERATURE [°C]
205	N	X		STATUS [0=stop; 1=interlock; 2=ramp; 3=regulation; 4=brake; 5=normal; 6=failure]
206	N	X		ERROR CODE:
				7 0 Too high load Short circuit Soft Start Not Ended Normal Not Reached No connection Pump overtemp Controller overtemp Power fail
211	N	X		CONTROLLER SENSOR TEMPERATURE READING [208= 25 °C - 128= 60 °C]
216	N	X		AMBIENT SENSOR TEMPERATURE READING
300	N	X		CYCLE TIME [min]
301	N	X		CYCLE NUMBER
302	N	X		PUMP LIFE [h]
319	A	X		Controller Model
320	A	X		Base Pump Model Number (8 characters)
321	A	X		Modified Standard Model Number (4 characters)
323	A	X		Controller Serial Number (5 characters)
325	A	X		Electrical Modification Level (10 characters)
400	A	X		CRC PROGRAM LISTING [QE7xxxx]
401	A	X		CRC BOOTLOADER [BL1xxxx]
402	A	X		CRC PARAMETER LISTING [PA7xxxx]
404	A	X		CRC FILE PARAMETER STRUCTURE
406	A	X		PROGRAM LISTING CODE & REVISION
407	A	X		PARAMETER LISTING CODE & REVISION

WIN = Window

R = Read

W = Write

T = Type

L = Logical

N = Numeric

A = Alphanumeric

Operation

Make all vacuum manifold and electrical connections and refer to Turbo pump instruction manual prior to operating the Turbo controller.

The controller operates completely automatically after the remote start command is given.

Switching On/Off the Pump

To switch on the pump is necessary to remove the connection between pin 12 of J4 connector and pin 15 (ground) of the same connector. To switch off the pump is necessary to connect pin 12 of J4 connector to pin 15 (ground) of the same connector.

Low Speed Activation/Deactivation

To activate the Low Speed it is necessary to connect pin 13 of J4 connector to pin 15 (ground) of the same connector.

To deactivate the Low Speed status it is necessary to disconnect pin 13 of J4 connector from pin 15 (ground) of the same connector.

The low speed status is available on pin 5 of the same connector.

The low speed frequency is 1100 Hz.

Error reset

If the controller identifies an error, the pump is switched off. To start again the pump the error cause has been removed and the error status has to be reset connecting pin 11 of J4 connector to pin 15 (ground) of the same connector for at least 1 second.

Maintenance

Replacement controllers are available on an advance exchange basis through Agilent. If necessary, information is provided to aid the operator in determining malfunctions and corrective steps to be taken.

WARNING!

Voltages developed in the unit are dangerous and may be fatal. Service must be performed by authorized personnel only.

Error Messages

The controller will self-diagnose some error causes. The controller shows the error occurred by means of a diagnostic LED located on the board, and through the serial communication port.

The LED blinks in a coded mode:

it flashes a number of time equal to the error code (see the following table) and then stays off for a while; this sequence has been repeated until the error status has been reset.

LED Error Code Table

Tab. 9

LED BLINKING NUMBER	DESCRIPTION
0	No error
1	Output current > 15 A
2	Not connected pump
3	Output current > 1.5 A for 15 sec.
4	Bearing temperature > 60 °C
5	Heat exchanger temperature > 60 °C
6	Normal not reached at run up time
9	Vp voltage 8V lower than nominal for 5 sec.
10	Soft start ramp not ended within the expected time

Accessories and Spare Parts

Tab. 10

DESCRIPTION	PART NUMBER
Pump Cable 1 m long with connector	9699869

Vacuum Products Division

Dear Customer,

Thank you for purchasing an Agilent vacuum product. At Agilent Vacuum Products Division we make every effort to ensure that you will be satisfied with the product and/or service you have purchased.

As part of our Continuous Improvement effort, we ask that you report to us any problem you may have had with the purchase or operation of our products. On the back side you find a Corrective Action request form that you may fill out in the first part and return to us.

This form is intended to supplement normal lines of communications and to resolve problems that existing systems are not addressing in an adequate or timely manner.

Upon receipt of your Corrective Action Request we will determine the Root Cause of the problem and take the necessary actions to eliminate it. You will be contacted by one of our employees who will review the problem with you and update you, with the second part of the same form, on our actions.

Your business is very important to us. Please, take the time and let us know how we can improve.

Sincerely,

Giampaolo LEVI

*Vice President and General Manager
Agilent Vacuum Products Division*

Note: Fax or mail the Customer Request for Action (see backside page) to Agilent Vacuum Products Division (Torino) – Quality Assurance or to your nearest Agilent representative for onward transmission to the same address.

Agilent Vacuum Products Division/Sales and Service Offices

United States

Agilent Technologies
121 Hartwell Avenue
Lexington, MA 02421 - USA
Ph. +1 781 861 7200
Fax: +1 781 860 5437
Toll-Free: +1 800 882 7426
vpl-customer-care@agilent.com

Netherlands

Agilent Technologies Netherlands B.V.
Customer Contact Center
Laan van Langerhuizen 1, toren A-8
1186 DS Amstelveen
Tel. +31 020 547 2600
Fax +31 020 654 5748
customer-care_netherlands@agilent.com

Belgium

Agilent Technologies Belgium S.A./N.V.
Customer Contact Center
Pegasus Park
De Kleetlaan 12A bus 12
B-1831 Diegem
Tel. +32 2 404 92 22
Fax +32 2 626 46 30
customer-care_belgium@agilent.com

Brazil

Agilent Technologies Brasil
Avenida Marcos Pentead de Ulhoa
Rodrigues, 939 - 6º andar
Castelo Branco Office Park
Torre Jacarandá - Tamboré
Barueri, São Paulo CEP: 06460-040
Toll free: 0800 728 1405

China

Agilent Technologies (China) Co. Ltd
No.3, Wang Jing Bei Lu, Chao Yang District
Beijing, 100102, China
Tel: +86 (0)10 64397888
Fax: +86 (0)10 64392765
Toll free: 800 820 6778 (ph); 400 820 6778 (cell.ph.)
vacuum.cnmarketing@agilent.com
vpc-customerservice@agilent.com

France

Agilent Technologies
Parc Technopolis - Z.A. de Courtaboeuf
3, avenue du Canada - CS 90263
91978 Les Ulis cedex, France
Tel: +33 (0) 1 64 53 61 15
Fax: +33 (0) 1 64 53 50 01
vpf.sales@agilent.com

Southeast Asia

Agilent Technologies Sales Sdn Bhd
Unit 201, Level 2 uptown 2,
2 Jalan SS21/37, Damansara Uptown
47400 Petaling Jaya, Selangor, Malaysia
Ph. +603 7712 6181
Fax: +603 7727 1239
Toll free: 1 800 880 805
vps-customerservice@agilent.com

India (Sales)

Agilent Technologies India Pvt. Ltd.
Unit Nos 110- 116, & Part of 101 & 109
First Floor, Splendor Forum,
Plot No.-3 , District Centre, Jasola
New Delhi-110025
Ph. +91 11 4623 7100
Fax: +91 4623 7105
Toll Free: 1 800 180 1517

Italy

Agilent Technologies Italia S.p.A.
Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel: +39 011 9979 111
Fax: +39 011 9979 350
Toll free: 00 800 234 234 00
vpt-customerservice@agilent.com

Japan

Agilent Technologies Japan, Ltd.
9-1 Takakura-cho Hachioji-city,
Tokyo, Japan
Tel.: +81 -3 -5232-1253
Fax: +81-120-565-154
Toll-Free: +81-120-477-111
jp-vvt-sales.pdl-ext@agilent.com

Singapore

Agilent Technologies Singapore Pte. Ltd.
1 Yishun Avenue 7,
Singapore 768923
Tel : (65) 6215 8045
Fax : (65) 6754 0574
Toll free: 1 800 2762622
vps-customerservice@agilent.com

Korea

Agilent Technologies Korea, Ltd.
Ilshin Building 4F
Yongsan-gu Hannam-daero
Seoul Korea 04418
Tel: +82 (0)2 2194 9449
Fax: +82 (0)2 2194 9853
Toll free: 080 222 2452
vpk-customerservice@agilent.com

UK and Ireland

Agilent Technologies LDA UK Limited
Lakeside Cheadle Royal Business Park
Cheadle, Cheshire SK8 3GR,
United Kingdom
Ph. +44 01865291570
Fax +44 01865291571
Toll free: 00 800 234 234 00
Toll free fax: 00 800 345 345 00
vpt-customer-care@agilent.com

India (Service)

Agilent Technologies India Pvt. Ltd.
C-Block, RMZ Centennial Plot Number-
8A, 8B, 8C, 8D,
Doddanakundi Industrial Area, ITPL
Road,
Mahadevapura Post, Bangalore- 560048
Ph. +91 80 40614000
Fax: +91 80 40148991

Taiwan

Agilent Technologies Taiwan Limited
No. 20 Gao-shuang Road, Ping-zhen Dist
Tao-Yuan City
32450 Taiwan, R.O.C.
Tel: +886 3 4959004
Toll free: 0800 018 768
vpw-customerservice@agilent.com

Germany and Austria

**Agilent Technologies
Sales & Services GmbH & Co. KG**
Lyoner Str. 20
60 528 Frankfurt am Main
GERMANY
Tel: +49 69 6773 43 2230
Fax: +49 69 6773 43 2250

Mexico

Agilent Technologies
Concepcion Beistegui No 109
Col Del Valle
C.P. 03100 – Mexico, D.F.
Tel.: +52 5 523 9465
Fax: +52 5 523 9472

Other Countries

Agilent Technologies Italia S.p.A.
Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel.: +39 011 997 9111
Fax: +39 011 997 9350
Toll-Free: 00 800 234 234 00
vpt-customer-care@agilent.com

Customer Support & Service

NORTH AMERICA:
Toll Free: 800 882 7426
vpl-ra@agilent.com
Lexington-service@agilent.com

EUROPE:

Toll Free: 00 800 234 234 00
vpt-customer-care@agilent.com

PACIFIC RIM:

please visit our website for individual
office information

<http://www.agilent.com/chem/vacuum>

**Worldwide Web Site, Catalog and Order
On-line:**

www.agilent.com/chem/vacuum
Representatives in most countries

© Agilent Technologies, Inc. 2016

Printed in ITALY

07/2018

Publication Number: 87-901-056-01

Agilent Technologies

CUSTOMER REQUEST FOR CORRECTIVE / PREVENTIVE / IMPROVEMENT ACTION

TO: AGILENT VACUUM PRODUCTS DIVISION TORINO – QUALITY ASSURANCE FAX

N°: XXXX-011-9979350

ADDRESS: AGILENT TECHNOLOGIES ITALIA S.p.A. – Vacuum Products Division –

Via F.lli Varian, 54 – 10040 Leini (TO) – Italy

E-MAIL: vpd-qualityassurance_pdl-ext@agilent.com

NAME	COMPANY	FUNCTION
ADDRESS:		
TEL. N° : FAX N° :		
E-MAIL:		
PROBLEM / SUGGESTION :		
REFERENCE INFORMATION (model n°, serial n°, ordering information, time to failure after installation, etc.):		
DATE		
CORRECTIVE ACTION PLAN / ACTUATION (by AGILENT VPD)		LOG N°

XXX = Code for dialing Italy from your country (es. 01139 from USA; 00139 from Japan, etc.)

Vacuum Products Division

Instructions for returning products

Dear Customer,

Please follow these instructions whenever one of our products needs to be returned.

Complete the attached **Request for Return form** and send it to Agilent Technologies (see below), taking particular care to include the completed **Health and Safety** declaration Section. No work can be started on your unit until we receive a completed copy of this form.

After evaluating the information, Agilent Technologies will provide you with a **Return Authorization (RA) number** via email or fax, as requested. Note: Depending on the type of return, a Purchase Order may be required at the time **the Request for Return is submitted**. We will quote any necessary services (evaluation, repair, special cleaning, eg).

Product preparation

- o Remove all accessories from the core product (e.g. inlet screens, vent valves).
- o Prior to shipment and if applicable for your product, drain any oils or other liquids, purge or flush all gasses, and wipe off any excess residue.
- o If ordering an Advance Exchange product, please use the packaging from the Advance Exchange to return the defective product.
- o Seal the product in a plastic bag, and package product carefully to avoid damage in transit. You are responsible for loss or damage in transit.
- o Include a copy of the Health and Safety Declaration in the shipping documentation on the outside of the shipping box of your returning product.
- o Clearly label package with RA number. Using the shipping label provided will ensure the proper address and RA number are on the package. Packages shipped to Agilent without a RA clearly written on the outside cannot be accepted and will be returned.
- o Return only products for which the RA was issued.

Shipping

- o Ship to the location specified on the printable label, which will be sent, along with the RA number, as soon as we have received all of the required information. Customer is responsible for freight charges on returning product.
- o Return shipments must comply with all applicable Shipping Regulations (IATA, DOT, ADR, etc.) and carrier requirements.

RETURN THE COMPLETED REQUEST FOR RETURN FORM TO YOUR NEAREST LOCATION:

EUROPE:

Fax: 00 39 011 9979 330
Fax Free: 00 800 345 345 00
Toll Free: 00 800 234 234 00

NORTH AMERICA:

Fax: 1 781 860 9252
Toll Free: 800 882 7426, Option 3

PACIFIC RIM:

Please visit our website for individual office information

vpt-customer@agilent.com

vpl-ra@agilent.com

<http://www.agilent.com>

TERMS AND CONDITIONS

Please read the terms and conditions below as they apply to all returns and are in addition to the Agilent Technologies Vacuum Product Division – Products and Services Terms of Sale.

- o Unless otherwise pre-negotiated, customer is responsible for the freight charges for the returning product. Return shipments must comply with all applicable **Shipping Regulations** (IATA, DOT, etc.) and carrier requirements.
- o Agilent Technologies is not responsible for returning customer provided packaging or containers.
- o Customers receiving an Advance Exchange product agree to return the defective, rebuildable part to Agilent Technologies **within 15 business days**. Failure to do so, or returning a non-rebuildable part (crashed), will result in an invoice for the non-returned/non-rebuildable part.
- o Returns for credit toward the purchase of new or refurbished Products are subject to prior Agilent approval and may incur a restocking fee. Please reference the original purchase order number.
- o Units returned for evaluation will be evaluated, and a quote for repair will be issued. If you choose to have the unit repaired, the cost of the evaluation will be deducted from the final repair pricing. A Purchase Order for the final repair price should be issued within 3 weeks of quotation date. Units without a Purchase Order for repair will be returned to the customer, and the evaluation fee will be invoiced.
- o Products returned that have not been drained from oil will be disposed.
- o A Special Cleaning fee will apply to all exposed products.
- o If requesting a calibration service, units must be functionally capable of being calibrated.

Customer information		
Company :		Contact Name:
Address:		Tel: Fax:
		Email:

Equipment			
Product description	Agilent PartNo	Agilent Serial No	Original Purchasing Reference
Failure description		Type of process (for which the equipment was used)	

Type of return	
☐ Non Billable ☐ Billable  New PO # (hard copy must be submitted with this form): _____	
☐ Exchange ☐ Repair ☐ Upgrade ☐ Demo ☐ Calibration ☐ Evaluation ☐ Return for Credit	

Health and safety		Substances (please refer to MSDS forms)			
The product has been exposed to the following substances: (by selecting 'YES' you MUST complete the table to the right)		* Agilent will not accept delivery of any product that is exposed to radioactive, biological, explosive substances or dioxins, PCB's without written evidence of decontamination.			
		Trade name	Chemical name	Chemical Symbol	CAS Number
Toxic	☐ YES ☐ NO				
Harmful	☐ YES ☐ NO				
Corrosive	☐ YES ☐ NO				
Reactive	☐ YES ☐ NO				
Flammable	☐ YES ☐ NO				
Explosive (*)	☐ YES ☐ NO				
Radioactive (*)	☐ YES ☐ NO				
Biological (*)	☐ YES ☐ NO				
Oxidizing	☐ YES ☐ NO				
Sensitizer	☐ YES ☐ NO				
Other dangerous substances	☐ YES ☐ NO				

Goods preparation	
If you have replied YES to one of the above questions. Has the product been purged? ☐ YES ☐ NO	
If yes, which cleaning agent/method:	
Has the product been drained from oil? ☐ YES ☐ NOT APPLICABLE	
I confirm to place this declaration on the outside of the shipping box. ☐	

I declare that the above information is true and complete to the best of my knowledge and belief. I understand and agree to the terms and conditions on page 2 of this document.	
Name:	Authorized Signature:
Position:	
Date:	
NOTE: If a product is received at Agilent which is contaminated with a toxic or hazardous material that was not disclosed, the customer will be held responsible for all costs incurred to ensure the safe handling of the product, and is liable for any harm or injury to Agilent employees as well as to any third party occurring as a result of exposure to toxic or hazardous materials present in the product.	

Agilent Vacuum Products Division/Sales and Service Offices

United States

Agilent Technologies

121 Hartwell Avenue
Lexington, MA 02421 - USA
Ph. +1 781 861 7200
Fax: +1 781 860 5437
Toll-Free: +1 800 882 7426
vpl-customer@agilent.com

Netherlands

Agilent Technologies Netherlands B.V.

Customer Contact Center
Laan van Langerhuizen 1, toren A-8
1186 DS Amstelveen
Tel. +31 020 547 2600
Fax +31 020 654 5748
customer@netherlands.agilent.com

Belgium

Agilent Technologies Belgium S.A./N.V.

Customer Contact Center
Pegasus Park
De Kleetlaan 12A bus 12
B-1831 Diegem
Tel. +32 2 404 92 22
Fax +32 2 626 46 30
customer@belgium.agilent.com

Brazil

Agilent Technologies Brasil

Avenida Marcos Penteado de Ulhoa
Rodrigues, 939 - 6° andar
Castelo Branco Office Park
Torre Jacarandá - Tamboré
Barueri, Sao Paulo CEP: 06460-040
Toll free: 0800 728 1405

China

Agilent Technologies (China) Co. Ltd

No.3, Wang Jing Bei Lu, Chao Yang District
Beijing, 100102, China
Tel: +86 (0)10 64397888
Fax: +86 (0)10 64392765
Toll free: 400 8206778 (mobile)
Toll free: 800 8206778 (landline)
vacuum.cnmarketing@agilent.com
vpc-customerservice@agilent.com

France

Agilent Technologies
Parc Technopolis - Z.A. de Courtaboeuf
3, avenue du Canada - CS 90263
91978 Les Ulis cedex, France
Tel: +33 (0) 1 64 53 61 15
Fax: +33 (0) 1 64 53 50 01
vpf.sales@agilent.com

Southeast Asia

Agilent Technologies Sales Sdn Bhd

Unit 201, Level 2 uptown 2,
2 Jalan SS21/37, Damansara Uptown
47400 Petaling Jaya, Selangor, Malaysia
Ph. +603 7712 6181
Fax: +603 7727 1239
Toll free: 1 800 880 805
vps-customerservice@agilent.com

India (Sales)

Agilent Technologies India Pvt. Ltd.

Unit Nos 110- 116, & Part of 101 & 109
First Floor, Splendor Forum,
Plot No.-3 , District Centre, Jasola
New Delhi-110025
Ph. +91 11 4623 7100
Fax: +91 4623 7105
Toll Free: 1 800 180 1517

Italy

Agilent Technologies Italia S.p.A.

Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel: +39 011 9979 111
Fax: +39 011 9979 350
Toll free: 00 800 234 234 00
vpt-customerservice@agilent.com

Japan

Agilent Technologies Japan, Ltd.

9-1 Takakura-cho Hachioji-city,
Tokyo, Japan
Tel.: +81- 3-5232-1253
Fax: +81-120-565-154
Toll-Free: +81-120-477-111
jp-vvt-sales.pdl-ext@agilent.com

Singapore

Agilent Technologies Singapore Pte. Ltd.

1 Yishun Avenue 7,
Singapore 768923
Tel : (65) 6215 8045
Fax : (65) 6754 0574
Toll free: 1 800 2762622
vps-customerservice@agilent.com

Korea

Agilent Technologies Korea, Ltd.

Ilshin Building 4F
Yongsan-gu Hannam-daero
Seoul Korea 04418
Tel: +82 (0)2 2194 9449
Fax: +82 (0)2 2194 9853
Toll free: 080 222 2452
vpc-customerservice@agilent.com

UK and Ireland

Agilent Technologies LDA UK Limited

Lakeside Cheadle Royal Business Park
Cheadle, Cheshire SK8 3GR,
United Kingdom
Ph. +44 01865291570
Fax +44 01865291571
Toll free: 00 800 234 234 00
Toll free fax: 00 800 345 345 00
vpt-customer@agilent.com

© Agilent Technologies, Inc. 2018
07/2018

India (Service)

Agilent Technologies India Pvt. Ltd.
C-Block, RMZ Centennial Plot Number- 8A, 8B, 8C,
8D,
Doddanakundi Industrial Area, ITPL Road,
Mahadevapura Post, Bangalore- 560048
Ph. +91 80 40614000
Fax: +91 80 40148991

Taiwan

Agilent Technologies Taiwan Limited

No. 20 Gao-shuang Road, Ping-zhen Dist
Tao-Yuan City
32450 Taiwan, R.O.C.
Tel: +886 3 4959004
Toll free: 0800 018 768
vpw-customerservice@agilent.com

Germany and Austria

Agilent Technologies Sales & Services GmbH & Co. KG

Lyoner Str. 20
60 528 Frankfurt am Main
GERMANY
Tel: +49 69 6773 43 2230
Fax: +49 69 6773 43 2250

Mexico

Agilent Technologies

Concepcion Beistegui No 109
Col Del Valle
C.P. 03100 – Mexico, D.F.
Tel.: +52 5 523 9465
Fax: +52 5 523 9472

Other Countries

Agilent Technologies Italia S.p.A.

Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel.: +39 011 997 9111
Fax: +39 011 997 9350
Toll-Free: 00 800 234 234 00
vpt-customer@agilent.com

Customer Support & Service

NORTH AMERICA:

Toll Free: 800 882 7426
vpl-ra@agilent.com
Lexington-service@agilent.com

EUROPE:

Toll Free: 00 800 234 234 00
vpt-customer@agilent.com

PACIFIC RIM:

please visit our website for individual office
information
<http://www.agilent.com/chem/vacuum>

**Worldwide Web Site, Catalog and Order On-
line:**

www.agilent.com/chem/vacuum
Representatives in most countries