Exotom-100

Instruction Manual

Manual No.: 15047001

Date of Release 30.07.2003

Exotom-100
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	19
Quick Reference Card	48

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2003.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Exotom-100

Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions on the cut-off wheels.
2. Use the supplied Allen key to remove the two red transport screws at the rear of the machine. Do not attempt to open the protection guard before the screws have been removed.
3. The machine must be placed on a safe and stable support. All safety functions of the machine must be in working order. The machine must be levelled by means of the adjustable legs provided.
4. The unit must be installed in compliance with local safety regulations.
5. Prior to lifting the unit by the built-in truck lifting point, ensure that the boom is properly secured with the locking pins provided. Before transport, secure the cutting arm with the locking system provided.
6. To achieve maximum safety and lifetime of the machine, use only original Struers consumables.
7. Use only intact cut-off wheels. The cut-off wheels must be approved for a spindle speed of min. 1950 rpm / 42 m/s. Do not use sawblades.
8. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
9. The workpiece must be securely clamped in the quick-clamping device or the like. Large or sharp workpieces must be handled in a safe way.
10. Struers recommend the use of an exhaust system as the cutting materials may emit harmful gasses or dust. See the instructions of the cut-off wheel.
11. The machine emits only moderate noise. However, the cutting process itself may emit noise, depending on the nature of the workpiece. In such cases, use of hearing protection is recommended.
12. The machine must be disconnected from the mains prior to any service.
13. Ensure that the cut-off wheel is secured before working on or around the cutting table.
14. Do not put hands through the rubber curtain during operation of the machine.

Exotom-100
Instruction Manual

- 15.** Protruding workpieces should be shielded or marked.
- 16.** Use of working gloves is recommended as workpieces may be both very hot and produce sharp edges.
- 17.** If any unusual noise is heard when the protection hood is operated, refrain from further use of the machine, and contact Struers service technician.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of Packing.....	3
Recirculation Cooling Unit.....	3
Unpacking and Placing Exotom-100.....	3
Mounting the Control Panel.....	4
Getting Acquainted with Exotom-100	5
Side view, left.....	6
Side view, right.....	6
Cooling unit compartment	7
Power Supply	8
Mounting the Cut-off wheel	8
Direction of the Cut-off Wheel	8
Connection to an External Exhaust System	8
Setting Up the Recirculation Cooling Unit	9
Direct Water Supply for the Cooling Unit.....	9
2. Operation	
Using the Controls.....	10
Control Panel of Exotom-100	10
Control Panel Functions.....	11
Display Types.....	11
Setting the Language	12
Reading the Cutting Display.....	13
Changing Cutting Mode and Cutting Parameters.....	14
Changing Cutting Mode	14
Changing Cutting Parameters.....	14
Reading the Motor Information	14
Sleep Mode	14
Changing the Cut-off Wheel.....	15
Clamping the Workpiece	15
Positioning the Cut-off Wheel	15
QuickPosition	16
Cutting Parameters	16
Feed	16
Force	16
Stop.....	16
AutoStop.....	16

Exotom-100
Instruction Manual

Setting the Cutting Parameters	17
Cutting on Exotom-100.....	17
Starting the Cutting.....	17
Fast advance	17
Stopping the Cutting (Manual Stop)	17
Re-starting cutting	17
Direct Cut.....	18
ExciCut.....	18
AxioCut (option).....	18
Additional cooling.....	18

1. Getting Started

Checking the Contents of Packing

The packing box contains the following items:

- 1 Exotom-100
- 1 Control Panel for Exotom-100
- 1 Allen key (5 mm) for mounting Control Panel
- 3 Screws (5 mm)
- 1 Fork spanner (30 mm) for cut-off wheel
- 1 Spanner (comb. 12/17 mm) for adjustment of the Control Panel
- 1 Drain hose, 2 m
- 1 Fitting for drain hose
- 2 Cabinet doors
- 1 Set of Instruction Manuals


Recirculation Cooling Unit

- 1 Recirculation cooling unit
- 1 Trolley for cooling unit


Unpacking and Placing Exotom-100

- Unscrew the nuts from the four transport bolts fixing the machine to the pallet.
- Lift the machine from the pallet by means of a forklift truck from the front, and place in a suitable location.
- Remove the safety-springs from the front crossbar, and remove bar.
- Use the supplied Allen key to remove the two red transport screws at the rear of the machine. Do not attempt to open the protection hood before the screws have been removed.
- Take out the loose parts (trolley, tank, drain hose etc.).
- Mount the cabinet doors.
- Turn the adjustable feet so that the machine stands firmly and is level.

NB!
Please remove the transportation support before use.


Mounting the Control Panel

- Unpack the Control Panel and mount on right side of unit using the two 5 mm screws and Allen key supplied (use holes ① and ② as indicated on drawing).
 - Mount the earth wire below the mounting plate using the third 5 mm screw ④.
 - Press cover towards cabinet wall until it snaps onto the mounting plate.
 - Connect Control Panel plug ⑥ to socket ⑤.
- The position of the Control Panel is controlled by a friction mounting. After mounting the control panel it is possible to adjust the mobility of the Control Panel this way:
- Adjust the main joint with the nut ③ using the combined 12/17 mm spanner supplied.
 - The joint just behind the Control Panel may be adjusted too, using the same spanner.

**Getting Acquainted with
Exotom-100**


Take a moment to familiarise yourself with the location and names of the Exotom-100 components.

- ① Cutting chamber
- ② Control Panel
- ③ Flexible water jets
- ④ Protection guards
- ⑤ Emergency stops
- ⑥ Flushing hose with flush nozzle
- ⑦ Cut-off wheel
- ⑧ Spindle lock button
- ⑨ Guard for cut-off wheel
- ⑩ Protection hood


Exotom-100
Instruction Manual

Side view, left


- ① Exhaust flange
- ② Water inlet
- ③ Water outlet

Side view, right

- ④ Name Plate
- ⑤ Type plate
- ⑥ Main Switch
- ⑦ Electrical connection box
- ⑧ Connection for electrical cable for power supply
- ⑨ Connection for external warning light
- ⑩ Fuse socket for external warning light
- ⑪ Connection for Control Panel

Cooling unit compartment

- ① Cooling unit drain pump on/off
- ② Filter drawer
- ③ Recirculation cooling unit
- ④ Trolley for cooling unit

Power Supply

Before connecting the machine, check that the mains voltage is correct by referring to the type plate.


Open the electric connection box and connect a 4-lead cable in the following way:

PE: earth
L1: phase
L2: phase
L3: phase

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Mounting the Cut-off wheel

- ① Nut
- ② Flange
- ③ Washer

- Press the knob for the spindle lock on the left-hand side of the cut-off wheel, while turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner (30mm).
- Remove the flange and the cut-off wheel.
- Mount the new cut-off wheel with a cardboard flange on each side.
- Mount the flanges and nut. Tighten carefully.

Direction of the Cut-off Wheel

To check that the cut-off wheel rotates in the direction indicated on the wheel-guard, do the following:

- Close the protection hood.
- Turn Main Switch to On.
- Press START .
- Check the rotational direction of the wheel.
- Press STOP .

If the direction of rotation is incorrect, turn mains off and switch two of the phases.

Connection to an External Exhaust System

Struers recommends the use of an exhaust system as workpieces may emit harmful gases when cut. The unit is prepared for connection to an exhaust system via an 80 mm fitting on the left hand side of the cabinet.

Recommended capacity for exhaust system: 350m³/h at 0mm water gauge.

Setting Up the Recirculation Cooling Unit

- Check that the nylon sieve in the filter drawer is properly placed and close the drawer.
- Mount the drain hose on the back of the machine.
- Place the recirculation cooling unit in the compartment.

IMPORTANT

The cooling unit must be positioned inside the compartment as indicated by the arrows on the edge of the cooling unit.

-
- Fill the cooling unit with 150 l of water (5-10 cm from the upper edge).
- Use additive enclosed with the machine and add to cooling water as specified on the bottle label and stir.

Direct Water Supply for the Cooling Unit

Exotom-100 may be directly connected to the water mains.


- Connect the water inlet fitting ($\frac{1}{2}$ " internal, $\frac{3}{4}$ " external thread).
- Open the external water supply.
- Fill the cooling unit with 150 l of water (5-10 cm from the upper edge).
- Shut off the water.
- Add additive in the quantity stated on the bottle label and stir.

IMPORTANT





When the cooling unit is pushed into the compartment, care should be taken that the flushing hose is not trapped behind it. To counter this, pull the flushing hose all the way out of its holder, then replace. The flushing hose, should then rest inside the cooling unit.

2. Operation

Using the Controls Control Panel of Exotom-100

Control Panel Functions

Name	Key	Function	Name	Key	Function
1 POWER	Power	Lights when the Main Switch has been turned to ON.	5 MENU		Push button to toggle between Cutting menu (daily use) and Configuration menu (basic settings).
2 START		Starts the machine and recirculation unit.	6 JOYSTICK		Move up- or downwards to position cut-off wheel.
3 STOP		Stops the machine and recirculation unit.	7 MULTI-FUNCTION KNOB		Push knob to select function. Turn knob to adjust settings.
4 CUTTING MODE		Push button to select desired cutting mode: Direct cutting, ExciCut, AxioCut/Step, AxioCut/Sweep.			

Display Types

The Cutting Display will be shown on the Control Panel when turning Exotom-100 on. This display is for everyday use. By pressing the MENU  button once, the CONFIGURATION menu will appear on the display. This menu will normally only be accessed during installation.

Setting the Language

Press MENU  button once to select CONFIGURATION Menu.

Turn knob to toggle between parameters in the CONFIGURATION Menu.


Push knob to select LANGUAGE. A pop-up menu appears.


From the pop-up menu, select your preferred language by pushing the knob.

Press MENU  button to move from CONFIGURATION Menu to Cutting Display.

Reading the Cutting Display

The Cutting Display offers three types of information:
 Cutting Mode
 Cutting Parameters, and
 Motor Information

- A** Cutting Mode.
- B** Cutting Parameters.
- C** Motor Information.

Cutting Mode

The top bar displays the selected Cutting Mode:
 Direct Cutting, ExciCut, AxioCut/Step or AxioCut/Sweep
 (AxioCut/Step and -/Sweep are optional).

Cutting Parameters and Motor Information

The large window on the Cutting Display, displays information about the Cutting Parameters (FEED, FORCE, STOP) as well as Motor information (TEMP, LOAD).

The Cutting Parameters (FEED, FORCE and STOP) can be set both prior to and during cutting. The set value is displayed below the column. Actual value is displayed above the bar graph.

The Motor Information columns TEMP and LOAD, inform about the condition of the motor during cutting. Read-out in relative value (%).


Changing Cutting Mode and Cutting Parameters

Changing Cutting Mode


Toggle between the four Cutting Modes, by pressing the CUTTING MODE button. The selected mode, Direct Cutting, ExciCut, AxioCut/Step or AxioCut/Sweep, appears on the top bar of the Cutting Display.

Changing Cutting Parameters

Toggle between the three Cutting Parameters by pushing the knob. The highlighted Cutting Parameter is selected.

Push knob to select Cutting Parameter.
Turn knob to adjust setting.

Actual values of the Cutting Parameters FEED, FORCE and STOP are displayed on top of the columns (A).
The actual position of the cut-off wheel (relative its starting position) is displayed graphically by the small icon  to the right of the STOP column.
Turn the knob to change the setting of the selected Cutting Parameter. The arrow on the left of the column will move to reflect the new setting (B).
By moving the pointer of the STOP column below the column, AUTO mode is selected.

Reading the Motor Information

The highlighted area of the bar graphs TEMP and LOAD, shows the status of the cutting motor:

TEMP. Temperature indicator of the cutting motor.

LOAD. Load indicator of the cutting motor.

Sleep Mode

To increase the lifetime of the display, the backlight is automatically switched off if Exotom-100 has not been used for 15 min. Push any key on the Control Panel to re-activate the backlight.

Changing the Cut-off Wheel

- Press the knob for the spindle lock on the left-hand side of the cut-off wheel, while turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner (30 mm).
- Remove the flange and the cut-off wheel.
- Mount the new cut-off wheel with a cardboard washer on each side.
- Mount the flange and nut. Tighten carefully and close the guard.

IMPORTANT

Conventional cut-off wheels based on Al_2O_3/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel.
For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

Clamping the Workpiece

- Clamp the workpiece with the clamping device of your choice e.g. a quick clamping device. Place the workpiece between the clamp and the back stop.
- Push the clamp towards the workpiece, and lock the quick clamping device with the locking handle.

Ensure that only one of the quick clamping devices is tight, the other device should only press lightly. Use support tools if the geometry of the workpiece makes support necessary.

IMPORTANT

When using the AxioCut cutting mode, the cut-off wheel cover may hit the jaw of the Quick Clamping Device, if same is mounted in the shown forward position; ie using the cross-directional T-slot closest to the front of the machine.

To avoid this situation, place the quick-clamping device as far back as possible, using the front-most T-slot.

Test possible infringement by starting the machine with the cut-off wheel completely clear of the clamping device.

Should the cut-off wheel hit the clamping device anyway, no danger is involved. Exotom-100 will stop automatically and display the following message: "AxioCut motor blocked".

Positioning the Cut-off Wheel


After clamping, test the correct position of the cut on the workpiece, by lowering the cut-off wheel into position, min. 1 mm above the workpiece:

- Push the joystick on the Control Panel downwards, to lower the cut-off wheel. The cut-off wheel stops when the joystick is released.

QuickPosition

Positioning of the cut-off wheel may be done, by simply pushing the joystick downwards until the cut-off wheel has contact with the workpiece. The force is automatically reduced on contact, preventing damage to the cut-off wheel and the workpiece. After contact with the workpiece, the cut-off wheel is automatically retracted 2 mm, to be ready for cutting.

Cutting Parameters

Feed

The feed speed can be set to values between 0.05-5.00 mm/sec (0.002-0.2 "/s).

Force

The maximum permitted force between cut-off wheel and workpiece can be set to values between 20-700 N (4-150 lbs). A built-in measuring cell constantly computes the force. If the force limit is reached, the feed speed will automatically be reduced to a value that allows the force to stay just below the set limit. As soon as the force drops below the set limit, the speed will be increased to the original setting.

NOTE

When cutting with a force below 50 N, beware of hysteresis (actual force applied may differ from selected cutting force).

Stop

There are two ways to set the stop position: AutoStop and Fixed Stop, explained under Advanced Operations.

AutoStop

When the AutoStop function is selected, the machine automatically stops when the workpiece has been cut through. For normal cutting, the AutoStop function is recommended.

Setting the Cutting Parameters

The Cutting Parameters, FEED, FORCE and STOP, are set in the Cutting Menu.

If in CONFIGURATION Menu, press MENU  once to select Cutting Menu.

Push knob to toggle between parameters in the Cutting Menu.

Turn knob to edit the value of the selected parameter.


Push knob to move to next parameter.

Cutting on Exotom-100

Starting the Cutting

- Position the cut-off wheel.
- Close the protection hood.
- Press START . The cut-off wheel starts rotating, the cooling water starts running and the cut-off wheel will slowly move down into the workpiece at the pre-set Feed Rate.

Fast advance

After START  has been pressed, the cut-off wheel may also be advanced towards the workpiece manually. To advance the cut-off wheel, push the joystick downwards. The cut-off wheel will advance towards the workpiece at maximum speed (5mm/s).

Stopping the Cutting (Manual Stop)

Exotom-100 automatically stops cutting when the workpiece is cut through.


- Press STOP  to manually interrupt the cutting process. The machine can be stopped at any time during operation by pressing STOP .

Re-starting cutting

Having interrupted cutting, the cut-off wheel has been automatically moved out of the workpiece.


- Press START  to resume cutting.

Direct Cut

Direct Cutting is the normal cutting mode. The cut-off wheel is moved into the workpiece in a slightly curved, vertical movement. Direct Cutting is intended for ordinary materials.

ExciCut

ExciCut cutting mode is ideal for cutting of very hard materials (HV >400). The oscillating movement of the cut-off wheel has three main advantages: less wear on the cut-off wheel, less risk of damage to workpiece and less risk of motor overheating. To select ExciCut, press the CUTTING MODE button, until ExciCut appears in the top bar of the display.


AxioCut (option)

The AxioCut option permits cutting of extra large workpieces (adding 150mm to max. depth). Two different cutting modes are provided, Step and Sweep. The Sweep-mode offers better cut-off wheel economy than the Step-mode. The latter method, however, is faster.

To select AxioCut/Step or AxioCut/Sweep, press the CUTTING MODE button, until the desired mode appears in the top bar of the display.

Additional cooling

Two flexible water jets are provided for additional cooling. Cooling water from the integrated water jets will flow away from the upper part of the cutting area, particularly when cutting hollow workpieces. To counter this, use the flexible water jets provided.

- Position the cooling jets to the left and right of the cutting area.
- Turn the valve on the jet to a position parallel to the hose to activate the cooling jet. The cooling water will flow as soon as the cutting starts.

Note

When using the flexible water jets, cooling water is diverted from the integrated water jets positioned over the cut-off wheel.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Configuration Menu	21
Display contrast.....	22
Language.....	22
Return position	22
Units.....	22
Operation mode	22
New Pass Code.....	23
Changing Operation Mode	24
Cutting Offset	25
Cutting Display.....	26
Cutting Modes.....	27
Direct Cut.....	27
ExciCut.....	27
AxioCut (option).....	27
AxioCut/Step.....	28
AxioCut/Sweep.....	28
Stop Settings	29
AutoStop.....	29
Fixed Stop	29
QuickPosition.....	30
OptiFeed.....	30
Cutting Long Workpieces	31
Clamping Irregular Workpieces.....	31
Connection to External Exhaust System.....	31
Other Safety Features	32
Optimising Cutting Results	33
2. Accessories	34
3. Consumables	
Cut-off Wheels.....	35
Other Consumables	35

4. Trouble-Shooting

Error Messages	39
Messages	39
Errors	39
Fatal Errors	39

5. Maintenance


Daily Service	42
Cleaning the Cutting Chamber	42
Weekly Service	42
Monthly Service	42
Replacing the Cooling Water	42
Maintenance of Recirculation Unit	43
Low Level Warning	43
Refilling the Cooling Unit	43
Emptying the Cooling Unit	43
Filter Drawer	44
Additive for Cooling Water	44
Maintenance of Cutting Table	45
Maintenance of Cut-off Wheels	45
Storage of Bakelite Bonded Al ₂ O ₃ Cut-off Wheels	45
Maintenance of Diamond and CBN Cut-off Wheels	45
Maintenance of Clamping Devices	45

6. Technical Data

Cutting Capacity	47
------------------------	----

1. Advanced Operations

Configuration Menu


Press MENU  button once to select CONFIGURATION Menu.

Turn knob to highlight different parameters in the CONFIGURATION Menu.

Push knob to edit the chosen parameter.

Press MENU  button to move from CONFIGURATION Menu to Cutting Display.

Display contrast

The contrast settings of the display can be adjusted to suit individual preferences (default value: 22, adjustment interval: 0-50).


Language

The language can be set to English (default), German, French, Spanish or Japanese.

Return position

After cutting or after pressing STOP , the return movement of the cut-off wheel can be set to three different functions:

Top: Exotom-100 automatically retracts the cut-off wheel to the top position.

Start: Exotom-100 automatically retracts the cut-off wheel to the original position of the cut-off wheel, at the time you pressed START  (default).

Stay: The cut-off wheel stays down.

Important

Use the Stay function for bakelite bonded diamond or CBN cut-off wheels, as retraction might destroy the rim of the cut-off wheel.

Units

The Feed, Force and Stop values in the display panel can be set to displayed in either mm/Newton (default) or inches/pounds force.

Operation mode

It is possible to select three different operation modes:

Configuration: Full functionality

Development: No access to parameters in CONFIGURATION menu, except Display contrast


Production: Access to START, STOP, Fixed Stop and movement of cut-off wheel, and to Display contrast in the CONFIGURATION menu

To select your own pass code, go to the CONFIGURATION menu. Select Operation mode to get access to CHANGE OPERATION MODE menu.

New Pass Code

To select your own pass code, go to the CONFIGURATION menu. Select Operation mode to get access to CHANGE OPERATION MODE menu.

 Push knob to select Pass code


 Turn knob until default pass code '176' appears.

 Push knob to select default pass code, and a new line (New pass code) appears in the CHANGE OPERATION MODE menu.


 Turn knob to move the cursor to New Pass Code.

 Push knob to select and turn knob to set your own three-digit pass code.


 Press Menu button to confirm your unique pass code and exit the CHANGE OPERATION MODE menu.

Changing Operation Mode

To change the operation mode, highlight the parameter in the CONFIGURATION menu.

 Push knob to move to Change Operation Mode menu.

 Enter your unique pass code by turning the knob.
Confirm by pushing knob.

 Push knob
A pop-up will appear.


 Select desired operation mode and push knob to confirm.

Press Cutting Mode button to move to CONFIGURATION menu.

Cutting Offset

Note:
 Cutting Offset requires installation of the AxioCut option.


The horizontal position (cutting offset) of the cut-off wheel can be set manually, prior to the cutting. In the ExciCut and Direct Cut cutting modes, the cut-off wheel may be moved 10 mm towards the front and 60 mm towards the back of the cutting chamber ('+10/-60', default value is 0). This allows for greater flexibility in terms of cutting difficult work pieces and clamping.


To change the position, select the parameter by pushing the knob. Turn knob to change the value. Push knob to confirm the desired value.

Once the new value has been confirmed, the cut-off wheel moves to the desired position - even if the protection hood is open.

If the cutting offset has been changed (i.e. not the default value of 0), the Cutting Display heading will remind the operator which offset has been selected.

Cutting Display

The Cutting Display appears on the Control Panel when Exotom-100 is turned on. The Cutting Display is designed as a single-layer user interface, i.e. without sub-menus. Whereas this menu is for everyday use, the CONFIGURATION menu (as described above) will normally only be used at installation.

If in CONFIGURATION menu, press MENU  to select Cutting Display.

Push knob to toggle between parameters in the Cutting Display. (FEED, FORCE, STOP).

Turn knob to edit the value of the selected parameter.

Push knob to move to the next parameter.

Cutting Modes

The Cutting menu has four cutting modes:

Direct Cut


ExciCut

AxioCut/Step (optional)

AxioCut/Sweep (optional)


To select Cutting Mode, press the CUTTING MODE button, until the desired mode appears in the top bar of the Cutting menu. The latest Cutting Mode used is saved, and will appear in the Cutting Display the next time Exotom-100 is switched on.

Direct Cut

Direct Cut is the normal cutting mode. The cut-off wheel is moved into the workpiece in a slightly curved, vertical movement, without any separate movement on the horizontal axis. Direct Cutting is intended for ordinary materials.

ExciCut

ExciCut cutting mode is ideal for fast cutting of very hard materials (HV >400). The oscillating movement of the cut-off wheel has three main advantages: less wear on cut-off wheel, less risk of damage to workpiece and less risk of motor overheating. To select ExciCut, press the CUTTING MODE button, until ExciCut appears in the top bar of the display.

AxioCut (option)


The AxioCut option permits cutting of extra large workpieces (adding 150 mm to max depth).


Two different modes, Step and Sweep, are provided; Step cutting offers faster cutting, whereas Sweep cutting offers better wheel economy.

Both methods are based on the Direct Cutting method, i.e. without oscillation of the cut-off wheel.

AxioCut/Step

AxioCut/Sweep

In Step-mode, the cut-off wheel enters the workpiece in three alternating, pre-programmed steps of 10mm. This method offers fast cutting of even very hard materials. To select AxioCut/Step, press the CUTTING MODE button, until the desired mode appears in the top bar of the display. Note that the two initial steps of a cycle, 1 & 2, are 5mm only. The cutting depth on steps 3 to 5 is 10mm. After completing step 5, steps 3 to step 5 are repeated until the workpiece has been cut through. The arrows on the drawing point towards the position of the cut-off wheel *after* cutting.

Sweep-mode offers better cut-off wheel economy than Step-mode, as a harder cut-off wheel can be used which in turn is more resistant to wear. The latter method, however, is faster.

AxioCut/Step is intended for cutting of large items. Attempting to cut small items in this mode may prove inefficient. Select alternative Cutting Mode.

In Sweep-mode, the cut-off wheel enters the workpiece at a pre-programmed feed speed. A Step column replaces the far-left column in the Cutting Display (Feed) when in Sweep-mode. FORCE is automatically set to maximum (700 N), to ensure utilisation of all available force, no matter which STEP size has been selected.


To select AxioCut/Sweep:

- Press the CUTTING MODE button, until the desired mode appears in the top bar of the display.
- Push knob to select STEP parameter.
- Turn knob to set STEP to the desired value (0.1-4.0 mm in steps of 0.1 mm).
- Push knob to select STOP parameter.
- Turn knob to set the desired STOP position.
- Press START ◊.

Note that the first step of a cycle moves the cut-off wheel from the back to the front of the workpiece without touching it. The arrows on the drawing point towards the position of the cut-off wheel *after* cutting/re-positioning. After completing step 6, steps 1 to 6 are repeated until the workpiece has been cut through.

Stop Settings

AutoStop

There are two ways to set the stop position: AutoStop or Fixed Stop.

When the AutoStop function is selected, the machine automatically stops when the workpiece has been cut through. For normal cutting, the AutoStop function is recommended. In the Cutting Menu, toggle between FORCE/FEED/STOP by pushing the knob, until STOP is highlighted. Turn knob to move pointer to setting AUTO at the bottom of the STOP column.

When the pointer in the STOP column is not on AUTO, Exotom-100 will only stop when it reaches the preset stop position or the STOP  key is pressed.

Note: AutoStop acts on sudden load changes. In some cases, eg using very low feed rate and/or low force, subtle load changes may not be detected by the machine, preventing AutoStop from working properly. This may particularly be the case when cutting soft materials, pipes or workpieces with changing cross-section. If you are facing problems repeatedly, the sensitivity of the AutoStop-function may be adjusted to fit your specific application. Please contact your Struers service engineer for details.

Fixed Stop


Fixed stop is used when a specific stop position is desired. When cutting pipes or other workpieces with changing cross sections, the cut-off wheel may retract before the workpiece has been cut. To overcome this, use the Fixed Stop function. Clamp the workpiece and position the cut-off wheel just above the workpiece.

This position is automatically set to 0 (zero) and indicated by the cut-off wheel icon  to the right of the Stop column. Accordingly, as soon as Start has been pressed, the actual position of the cut-off wheel, becomes a relative starting point (zero), from where the cutting depth is calculated.

Select the STOP parameter and set the desired stop position by turning the knob.

Exotom-100 will now stop when it reaches the pre-set stop position.

Remember to take the wear of the cut-off wheel into consideration.


QuickPosition

Positioning of the cut-off wheel may be done automatically, by simply pushing the joystick downwards until the cut-off wheel is in contact with the workpiece. The force is automatically reduced on contact, preventing damage to the cut-off wheel and the workpiece.

After contact with the workpiece, the cut-off wheel is automatically retracted 2mm, to be ready for cutting.

OptiFeed


During cutting, Exotom-100 continuously measures the load on the cutting arm.

The pre-set Feed and Force values are interpreted as maximum values. Throughout cutting, Exotom-100 will try to stay as close to these values as possible.

The factors that determine the load, are the shape and properties of the workpiece.

Whenever the set maximum Force limit is reached, Exotom-100 will reduce the Feed rate to be able to maintain the present force.

The figure below illustrates the increase in force as the cut-off wheel gets closer to the centre of a round workpiece.

Cutting Long Workpieces

Workpieces that are wider than the cutting chamber can be positioned so that one or both ends can protrude through the rubber curtains on both sides of the cutting chamber.

Clamping Irregular Workpieces

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the sample itself. Use the T-slots to mount the special clamping tools. Struers offers a kit of Clamping Tools (See Accessories).

To achieve faster cutting, position the workpiece so that the wheel will cut the smallest possible cross-section.

Connection to External Exhaust System

Struers recommends the use of an exhaust system as the workpieces may emit harmful gasses or dust during cutting. A joint for connection of a 80 mm dia. exhaust hose is located on the left side of Exotom-100.

- Remove the red cap.
- Mount an exhaust hose from your local exhaust system onto the flange.

Other Safety Features

To prevent accidental access to the cutting chamber via the cutting curtains, each opening is covered by a spring-loaded guard, which only opens outwards.

The protection hood has a safety switch to prevent the cut-off wheel from starting while the hood is open. Furthermore, a locking mechanism prevents the opening of the protection hood before the cut-off wheel has come to a complete stop.

The motors of Exotom-100 are protected against overload. Should the motors overheat and/or overload, the motors will disengage until a normal temperature has been obtained.

Exotom-100 is prepared for the connection of an external warning light. The connection socket is situated on the right hand side of the unit.

Optimising Cutting Results

The following table shows possible answers to a number of common questions:

Optimising the Cutting Results	
Question	Answer
How can I avoid discoloration or burning of the sample?	Use a lower Feed Speed.
	Change the cut-off wheel as the hardness of the present cut-off wheel may be inappropriate for the hardness of the sample. *)
How can I avoid burrs?	Use a softer cut-off wheel. *)
	Clamp the workpiece securely at the right hand clamping device. Tighten the left hand clamping device just enough to prevent the workpiece from moving when being cut.
How can I avoid the cut-off wheels wearing too quickly?	Use a lower Feed Speed, a different cutting mode or a harder cut-off wheel. *)
How can I achieve faster cutting?	Position the workpiece so as to cut the smallest possible cross-section. Use a high Feed Speed. If workpiece shape and properties permit, change to ExciCut or AxioCut/Step (optional) Cutting Modes

*) See Consumables Section, Cut-off Wheels

2. Accessories

Specification	Code
<p>AxioCut Kit for Exotom-100 Increases cutting capacity (depth) by 150 mm. Requires installation by Struers Service Technician. Gross weight 9 kg</p>	EXOMA
<p>Replacement Stainless Steel Bands For Cutting Table with 12 mm T-slots. Available in three sizes, each supplied in sets of two: 100 x 225 mm 100 x 300 mm 60 x 548 mm</p>	EXOBS EXOBM EXOBL
<p>Quick-clamping Device for 12 mm T-slots, Left For securing the workpiece. Complete with back stop. To be mounted on a cutting table and placed to the left of the cut-off wheel.</p>	EXOLE
<p>Quick-clamping Device for 12 mm T-slots, Right For securing the workpiece. Complete with back stop. To be mounted on a cutting table and placed to the right of the cut-off wheel.</p>	EXORI
<p>Clamping Tools for 12 mm T-slots For clamping irregularly shaped workpieces. Complete with clamps, supports and bolts.</p>	MAGOF
<p>Vertical Clamping System for 12 mm T-slots For clamping irregularly shaped workpieces. Complete with operating key and one flat clamping shoe.</p>	EXOVS
<p>Chain Spanner for 12 mm T-slots For clamping cylindrical or irregularly shaped workpieces. Complete with anchor block and operating key.</p>	EXOCS
<p>Swivel Shoes for Vertical Clamping System For clamping irregularly shaped workpieces. Set of four multi-shaped swivel shoes.</p>	TREVI
<p>Riser Block for 12 mm Vertical Clamping System For elevating Vertical Clamping System when clamping high workpieces</p>	EXOKS
<p>Arm Extension for 12 mm Vertical Clamping System For extending arm of Vertical Clamping System</p>	EXOAR
<p>Spare Nylon Filter for Recirculation Cooling System For insertion into filter-drawer in cooling unit compartment.</p>	EXOFI

3. Consumables

Cut-off Wheels

Application	Dimensions	Grinding agent	Code
For ExciCut and direct cutting of extremely hard steels (350-800 HV)	432x3x32mm	Al ₂ O ₃	101MA
For ExciCut and direct cutting of very hard and ductile ferrous metals. Strong, fibre-reinforced wheel (250-700 HV)	432x3x32mm	Al ₂ O ₃	202MA
For ExciCut and direct cutting of steels (< 500 HV). Extremely hard (white cast iron) or ductile (18/8 or St60) steel types are cut with ExciCut. Gives good wheel economy	432x3x32mm	Al ₂ O ₃	102MA
For ExciCut and direct cutting of steels (< 500 HV). Gives very good wheel economy.	432x3x32mm	Al ₂ O ₃	104MA
For ExciCut and direct cutting of non-ferrous metals (< 300 HV). Primarily for ExciCut.	432x3x32mm	SiC	106MA
For direct cutting of titanium and its alloys, as well as other non-ferrous metals (< 400 HV)	432x3x32mm	SiC	106MA
For direct cutting of sintered carbides and ceramics	305x1.8x32mm	Diamond	25EXO
	350x1.8x32mm	Diamond	26EXO

Other Consumables

Specification	Code
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE

4. Trouble-Shooting

Error	Explanation	Action
Machine Problems		
Water leaking.	Leak in recirculation water hose.	Check the hose and tighten the hose clamp.
	Water overflow in the recirculation water tank.	Remove the excess water in the tank.
Samples or cutting chamber corroded.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water, using the correct concentration. Check with a refractometer. Follow the instructions in the Maintenance section.
	The machine is left with closed protection hood.	Leave the protection hood open to let the cutting chamber dry.
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not in balance.	Adjust the screw underneath the clamping column. Use a 3mm Allen key.
	Clamping heart worn.	Call Struers service technician.

Exotom-100
Instruction Manual

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the sample.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the sample.	See Consumables section, Cut-off Wheels.
	Inadequate cooling.	Check that there is enough water in the recirculation cooling unit. Check the condition of the cooling tray.
	Feed Speed too high.	Reduce the Feed Speed.
Unwanted burrs.	Wheel too hard.	See Consumables section, Cut-off Wheels.
	Insufficient support of workpiece.	Add further support to the workpiece.
The cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
	Insufficient cooling water.	Refill tank with water. Remember to add Struers Additive.
The cut bends to one side.	Feed Speed is too high.	Reduce the Feed Speed.
The cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	Check that the centre-hole has the correct diameter. Check cardboard washer on both sides of the cut-off wheel. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only one of the quick-clamping devices is tight. The other device should only press lightly. Use support tools if the geometry of the workpiece makes support necessary.
	Wheel too hard.	See Consumables section, Cut-off Wheels.
	Feed Speed too high.	Reduce the Feed Speed.
	Inadequate cooling.	Check that there is enough water in the recirculation cooling unit. Check the cooling water hoses.
The cut-off wheel wears down too quickly.	Feed Speed too high.	Reduce the Feed Speed.
	Insufficient cooling.	Check that there is enough water in the recirculation cooling unit. Check the cooling water hoses.
	The cut-off wheel is too soft for the task.	See Consumables section, Cut-off Wheels.
	Exotom-100 vibrates (worn bearings).	Call Struers service technician.

Exotom-100
Instruction Manual

Error	Explanation	Action
The cut-off wheel does not cut through the sample.	Incorrect choice of cut-off wheel.	See Consumables section, Cut-off Wheels.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel in such a way to allow the cut to stay open.
	Incorrect choice of Cutting Mode. AxioCut/Step (option) is intended for large workpieces.	See section on Operation, Cutting Mode.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. NB! Always cut brittle workpieces very carefully.
The sample is corroded.	The sample is not resistant to water.	Use a neutral liquid as cooling fluid or cut without using cooling fluid at all. DO NOT USE AN INFLAMMABLE LIQUID
	The sample has been left in the cutting chamber for too long.	Leave the protection hood open, when you leave the machine.
	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water in the correct concentration. Check with a refractometer. See Maintenance section.
AutoStop does not stop the cutting action.	Feed Speed and/or Force is set too low to allow load change to be detected.	If problem persists, contact local Struers service engineer, for adjustment of AutoStop sensitivity.

Error Messages

Error messages are divided into three classes:
 Messages
 Errors
 Fatal Errors

Messages

Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.

Errors

Errors must be rectified before cutting can be continued.

Fatal Errors

In case of Fatal Errors, cutting cannot continue before an authorised technician has rectified the error. Turn off the unit at the Main Switch immediately. Do not attempt to operate unit before technician has rectified problem.

Message	Explanation	Action
Process in Progress.	Start is not possible as other process in progress e.g. emptying of recirculation tank.	
Process already stopping.	Appears if Stop is pressed repeatedly.	
Water level in tank too low; or pump filter clogged.	Water level in recirculation tank is too low. Alternatively, recirculation jumper filter clogged.	Refill circulation water if needed. Check condition of pump filter. Clean if needed.
Water level low; emptying finished.	Emptying of recirculation cooling unit completed.	Switch off recirculation pump to complete emptying process.
Flush nozzle not parked.	Start not possible. Start was pressed while flush nozzle was out of holder.	Park flush nozzle to continue.
Flushing was interrupted to avoid pump damage. Please park the flush nozzle to continue.	The pump automatically stops after 3 minutes to avoid pump damage.	Park flush nozzle to continue.
Emptying of tank was stopped to avoid pump damage.	The recirculation pump automatically stops after 8 minutes to avoid pump damage.	To re-start the pump, turn the pump handle horizontally, then re-activate the pump by turning the pump handle to the vertical position.
Valve handle in emptying tank pos.	Drain pump handle has not been moved to off position after use.	Re-position drain pump handle.
Do not operate handle now.	During cutting, an attempt is made to activate tank emptying function manually.	
Reference position not found.	Stop was pressed during search for reference point.	Press Start again to do new search. If message remains, switch unit off and on at Main Switch to start new search for reference point.
Cutting stopped by flange sensor.	Flange sensors next to the cut-off wheel flange will stop the movement of the cutting arm if it meets an obstacle.	Replace cut-off wheel or re-position workpiece.

Exotom-100
Instruction Manual

Message	Explanation	Action
Searching for reference position.	When unit is switched on, a reference position process is performed.	Wait until process has stopped.
Protection hood not closed.		Close hood properly.
Emergency stop is active.		Release emergency stop(s).
Out of Cutting Range.	Wheel at bottom position when Start is pressed.	Move cut-off wheel to higher start position.
Cutting finished.	Message after each cutting which has been completed normally. Protection hood can be opened.	
Process stopping.	Appears after Stop has been pressed once.	
Error Message	Explanation	Action
Main supply voltage too low.	Supply voltage insufficient.	Await normalisation of power grid.
Cut-off motor overloaded.	Wait until the motor has cooled down.	Press Start.
ExciCut motor overloaded.	Wait until the motor has cooled down.	Press Start.
Protection hood not locked.	Safety lock has failed to lock hood.	Re-close hood.
Cut-off motor blocked.	Cut-off motor may fail to move the wheel if caught in the workpiece.	Move cut-off wheel clear of workpiece before pressing Start.
Contact K3 not activated.	Contact fails to operate.	Re-start. If error remains, contact Struers service technician.
Contact K4 not activated.	Contact fails to operate.	Re-start machine using Start and Stop buttons. Alternatively, re-start using Main Switch If error remains, contact Struers service technician.
Step motor or encoder error.	Exotom-100 cannot complete reference point search.	Re-start using Main Switch. If error remains, contact Struers service technician.
Reference sensor not activated.	Search for reference point cannot be conducted.	Check whether cutting arm is blocked. Re-start using Main Switch. If error remains, contact Struers service technician.
Step motor may be blocked.	Cutting arm cannot move, most likely due to blockage of arm.	Remove blockage. If error remains, contact Struers service technician.
AxioCut motor blocked.	Cutting arm cannot move, most likely due to blockage of arm or AxioCut motor may fail to move the wheel if gets caught in the workpiece.	Move obstacle or move cut-off wheel clear of workpiece before pressing Start.
Recirculation pump overloaded.	Water hose blocked or too much swarf in coolant.	Replace coolant. Inspect pump.

Exotom-100
Instruction Manual

Fatal Error	Explanation	Action
Protection hood lock not released.		Re-start machine using Start and Stop buttons. Alternatively, re-start using Main Switch.
Cut-off motor will not stop.		Turn off unit at Main Switch. Contact Struers service technician.
Contactora K3 not deactivated.	Error in control system.	
Contactora K4 not deactivated.	Error in control system.	
Safety relay error during cutting.		
Force system uncalibrated.	Force-measuring system is not calibrated.	
No serial communication.	No contact between machine system and operating panel.	
Nominal motor current undefined.	Motor current has not been defined during set-up of unit.	
Program version mismatch.	Conflict between software in machine system and operating panel.	
AxioCut motor driver error.	When unit is switched on, AxioCut motor starts operating.	

5. Maintenance

Daily Service

- Clean the cutting chamber, especially the cutting table and the T-slots.

Cleaning the Cutting Chamber

To ensure a longer lifetime for your Exotom-100, Struers strongly recommends daily cleaning of the cutting chamber with the flushing hose:

- Lift the flushing hose from the holder on the rear wall. The recirculation pump is activated.
- Point the flushing hose towards the bottom of the cutting chamber.
- Press the lever on the flush nozzle to turn the water on. Adjust the flow using the lever.
- Clean the cutting chamber thoroughly.
- Turn off the water by releasing the lever on the flush nozzle.

IMPORTANT

Do **not** clean the lamp glass with alcohol.
Use a soft, damp cloth.

Weekly Service

- Refill the Recirculation Cooling Unit. The water level should be 5-10 cm below the upper edge of the tank.
- Clean the cutting chamber thoroughly.
- Clean the filter in the filter drawer.
- Check and if necessary empty and clean the sieve in the drawer of the Recirculation Cooling Unit.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

Monthly Service

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Fill the tank with 150 l of water. The water level should be 5-10 cm below the upper edge of the tank.
- Add Struers Additive according to the instructions.

Maintenance of Recirculation Unit

Low Level Warning

The Exotom-100 will automatically display a warning if the level of cooling water gets too low. However, it is recommended that the cooling water level is checked every week.

Refilling the Cooling Unit

- Refill the cooling unit with water using an external water hose or via direct connection to the water mains.
- Stop refilling when the water level is 5-10 cm from the upper edge.

NB!

Never refill by pouring the water into the cutting chamber, as it would be impossible to check the increasing level in the tank.

IMPORTANT


Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

Use only Struers Additive ADDUN in the Recirculation Cooling Unit.

ADDUN is specially selected for use with Struers cut-off machines.

Other additives may not be compatible with some of the components of the cut-off machine.

Emptying the Cooling Unit

When replacing the cooling water, the cooling unit is emptied using the recirculation pump. Turn the black pump handle in the cooling unit compartment to the vertical position to activate the pump.

By design, the pump cannot empty the tank completely, preventing cutting waste from entering the cooling system. Tilt the tank, to release the remaining water.


IMPORTANT

Remember to turn the pump handle horizontally (off position) after use.

Note

The recirculation pump automatically stops after 8 minutes to avoid pump damage. To re-start the pump, turn the pump handle horizontally, then re-activate the pump by turning the pump handle to the vertical position.

Filter Drawer

The filter drawer is situated on the left in the cooling unit compartment. A nylon sieve filters larger cutting waste particles from the cooling water. It is recommended to always have the filter drawer mounted during cutting.

Depending on use, the filter must be checked and emptied regularly, to ensure its effectiveness. Should the filter get blocked, an overflow function will prevent flooding, and direct the unfiltered cooling water into the tank.


Additive for Cooling Water

Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

The concentration of additive should always lie between 2.7 and 3.3%; the pH should be around 9. Replace the cooling water if pH is ≤ 8.0 .

To check the concentration of additive, use a refractometer (0-10 Brix, see diagram.). Concentration = $2 \times \text{Brix}$.

Add Struers Additive for Cooling Fluid if the concentration is lower than 2.7%. Add water if the concentration is higher than 3.3%.

Maintenance of Cutting Table

The stainless steel bands forming the cutting table should be replaced if they have become worn or damaged. The bands are available as spare parts.

To allow humidity to escape from the cutting table and chamber, it is recommended to leave the hood open when the machine is not in use.

Maintenance of Cut-off Wheels

*Storage of Bakelite Bonded
Al₂O₃ Cut-off Wheels*

This type of cut-off wheel is sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

*Maintenance of Diamond and
CBN Cut-off Wheels*

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.


6. Technical Data

Subject		Specification	
		Metric/International	US
CUTTING SPECIFICATIONS			
Workpiece Dimensions (Max)	<i>Workpiece inside cutting chamber:</i>		
	Height	250 mm	10"
	Width	650 mm	25.6"
	Depth	550 mm	20"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	140 mm	5.5"
	Depth	210 mm	8.3"
Cutting Capacity (Max) <i>Please refer to Cutting Diagram</i>	Max. workpiece diameter	160 mm	6.3"
	<i>Max. size of cut. In brackets, with AxioCut option</i>		
	Height	100 mm (100 mm)	3.9" (3.9")
	Depth	300 mm (450 mm)	11.8" (17.7")
PHYSICAL SPECIFICATIONS			
Cutting Motor	Cutting power*	10.5 kW	14 HP
	Maximum power*	15 kW	20 HP
	<small>*) Except 460-480 V / 60 Hz</small>		
Cut-off Wheel	Diameter x Thickness x Centre-hole	432 x 3 x 32 mm	17 x 0.12 x 1.26"
	Rotational speed (running idle)	1950 rpm	1950 rpm
Positioning & Feed	Positioning range (of cut-off wheel)	0 – 250 mm	0 – 9.8"
	Max. positioning speed	60 mm/s	2.4"/s
	Feed Speed range (adjustable in steps of)	0.05 – 5 mm/s (0.05mm/s)	0.002 – 0.2"/s (0.002"/s)
	Cutting force	50-700 N	10-150 lbf
Cutting table	Width	855 mm	33.6"
	Depth	550 mm	21.6"
	T-slots	12 mm	0.48"
Dimensions and Weight	Height (closed/open hood)	1756/2300 mm	69/90.5"
	Width (ex-/including control panel)	1050/1350 mm	41.3/53.1"
	Depth	1500 mm	59.1"
	Weight	820 kg	1810 lbs
Recirculation Cooling Unit	Tank volume	150 l	39.6 gallons
	Approx. flow	25 l/min	6.6 gallons/min

Exotom-100
Instruction Manual

Subject		Specification			
Electrical Data	Voltage / frequency:	Cutting Power*	Max. Power	Nom. Load	Max. Load
	3 x 200 V / 50 Hz	10.5 kW	15 kW	36 A	69 A
	3 x 200-210 V / 60 Hz	10.5 kW	15 kW	35 A	67 A
	3 x 220-230 V / 50 Hz	10.5 kW	15 kW	31 A	60 A
	3 x 220-240 V / 60 Hz	10.5 kW	15 kW	29 A	55 A
	3 x 380-415 V / 50 Hz	10.5 kW	15 kW	18 A	34 A
	3 x 380-415 V / 60 Hz	10.5 kW	15 kW	17 A	32 A
3 x 460-480 V / 60 Hz	12.6 kW	18 kW	18 A	34 A	
	*) at intermittent duty, S3 15%				
Environment	Noise level	Approx. 78 dB(A) running idle, at a distance of 1.0 m / 39.4" from the machine.			

Cutting Capacity

The graph shows the projected cutting capacity under the following conditions: A new cut-off wheel. The workpiece is laid directly on the cutting table, with overhang where appropriate. Vertical clamping is used. The actual cutting capacity depends on the sample material, cut-off wheel and clamping technique.

Quick Reference Card

Clamping the Workpiece

- Place the workpiece between the clamp and the back stop on the right hand cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Starting the Cutting

- Position the cut-off wheel.
- Close the protection hood.
- Press START \blacklozenge . The cut-off wheel starts rotating and the cooling water starts running.

Stopping the Cutting

Auto Stop

- In the Cutting Menu, toggle between FORCE/FEED/STOP pushing the knob, until STOP is highlighted.
- Turn knob to move pointer to setting AUTO at the bottom of the STOP column.

Manual Stop

- Press STOP \odot and the cutting process stops. The cut-off wheel stops rotating and the cooling water stops.

Setting the Cutting Parameters

- Press MENU \square once if CONFIGURATION MENU is selected.
- Push knob to toggle between parameters in the Cutting Menu.
- Turn knob to edit the value of the selected parameter.
- Push knob to move to next parameter.

Changing the Cut-off Wheel

- Press the knob for the spindle lock and turn the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner (30 mm).
- Remove the flange and the cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and nut. Tighten carefully and close the guard.

Cleaning the Cutting Chamber

- Point the flushing hose towards the bottom of the cutting chamber.
- Press the lever on the flush nozzle to turn the water on. Adjust the flow using the lever.
- Clean the cutting chamber thoroughly.
- Turn off the water by releasing the lever on the flush nozzle.

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
Date of Release 30.07.2003

Exotom-100
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2003.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents

Drawing

Exotom-100

Drawings

Exotom-100, complete.....	15040001E
Chassis, complete.....	15040003L
Pump, assembly.....	15040009C
Chassis, assembly.....	15040010K
Wire drum, assembly	15040013B
Counter weight, assembly.....	15040014
Cutting wheel guard.....	15040022D
Flushing Gun, complete	15040060
Shock absorber, complete.....	15040076A
Cutting table, complete	15040005G
Long-cut kit, complete.....	15040040C
Fixed rod	15040041A
Long-cut motor	15040044A
ExiCut motor	15040046D
High power electrical assembly	15040011D
Protection Guard, complete.....	15040016I
Lamp, complete.....	15040017
Top hat, assembly.....	15040036B
Cutting arm, assembly	15040021E
Tilting frame, assembly.....	15040025B
Cutting motor, complete.....	15040028D
Door for recirculation, RHS.....	15040031A
Door for recirculation, LHS.....	15040032A
Control panel, complete.....	15040050
Control panel assembly, front.....	15040055C
Water Panel, complete.....	15040065D
Feed Actuator, complete.....	15040070F
Stepper Motor, assembly	15040074B
Snap lock tube, complete.....	15040067
Quick clamping device for 12 mm T-slots, RHS compl.	15040081
Quick clamping device for 12 mm T-slots, LHS compl.	15040098
Back Stop for 12 mm T-slots, RHS complete.....	15040084
Back Stop for 12 mm T-slots, LHS complete.....	15040093

Exotom-100
Spare Parts and Diagrams

Diagrams

Exotom-100 Block Diagram	15043050
Exotom-100, Control voltage diagram (5 pages)	15043110A
Exotom-100, Main supply circuit	15043100A
Exotom-100, Control panel.....	15043200B
Exotom-100, Machine control board, A1 Power supply (8 pages)	15043210A
Exotom-100, LED circuit board	15043220
Exotom-100, Wiring diagram (5 pages)	15043450B
Primary connections on transformer AA-70951B	15043470A

Some of the drawings may contain position numbers
not used in connection with this manual.

Exotom-100
Spare Parts and Diagrams

Drawing	Spare Part list for Exotom-100		Cat no.
	Pos.		
		Rough Filter	15040397
		101MA Cut-off wheel ø432mm (5 pcs)	40009001
		104MA Cut-off wheel ø432mm (5 pcs)	40009004
		202MA Cut-off wheel ø432mm (5 pcs)	40009003
		ADDFI Additive 5L	49900003
		Fork Spanner 30mm	2GR00230
		Fork Spanner 13mm	2GR00213
		Pressure hose, water, white	2NU93020
15040001		Exotom-100, complete	
	150	Curtain, back wall	15040315
15040003		Chassis, complete	
	20	T1: Transformer 200-460V/ 600VA	2MT70951
	140	PCB A1, Machine control	15043393
	390	Fuse 4.00AT Glass (10 pcs)	RFU14200
	590	PD Eagle Belt W1792	2JT61792
	1090	Cardboard blotters 25 Pcs, ø32.5/120	381MP431
	1100	Loose Flange	R5040160
	1340	Shock Absorber with screws, complete	R5040076
15040009		Pump, assembly	
	20	Kit for impeller pump	15042901
15040010		Chassis, assembly	
	326	Snap Lock hose, complete	15040067
	390	Filter for impeller pump, complete	15040096
	400	Flushing gun, complete	15040060
	520	Perlator 3/8"	2YI20038
15040013		Wire drum, assembly	
	50-60	Bearing house, complete ø28	15042905
	100	Bearing Flange w. House PFT 25 FM	RBK40025
15040014		Counter weight, assembly	
	20	Glide guides for weight (4 pcs)	RBG00121
15040022		Cutting wheel guard, assembly	
	30	Sensor strap with Screws	R5040166
	50	Proximity Switch DJ2 (Namur)	2HQ00022

Exotom-100
Spare Parts and Diagrams

Drawing	Pos.		Cat no.
15040060		Flushing gun, complete	
	10	Flushing gun, complete	15040060
	30	Perlator for flushing gun	2YI26211
		Shower Brush	15042938
	40	Shower Hose	2NU99010
15040076		Shock Absorber, complete	
	40	Update shock absorber cartridge kit. Up to s/n 5040020.	RYS10016
	40	Replacement shock absorber cartridges. Use only from s/n 5040020.	15040319
	70-80	Bearing kit for Shock Absorber	15042906
15040005		Cutting table, complete	
	110	Timing Belt GT-36-5MR-09	2JT20450
	120	EXOBL Replacement steel bands for 12 mm T-slots 60x548 (Set of 2 pcs)	05046909
	150	EXOBM Replacement steel bands for 12 mm T-slots 100x300 (Set of 2 pcs)	05046908
	160	EXOBS Replacement steel bands for 12 mm T-slots 100x225 (Set of 2 pcs)	05046907
	340	Spring kit for side flap hinge	15042904
15040040		Long-cut kit, complete	
	100	M5, Long cut motor assembly	15040044
	120	Timing belt, GT 290 3MR09	2JT60270
15040041		Fixed rod	
	20	Vibration damper	2GS41260
	40	Joint bearing GA20	2BL20020
15040044		Long-cut Motor, M5	
	10	M5, Long cut motor assembly	15040044
15040046		Excicut Motor, M3	
	10	M3-Motor 1x230V/24W	15040167
15040011		High power electrical assembly	
F1 (3x200-240V/50-60Hz):	51	6.00 A Midget fuse 10x38 (3 pcs)	RFC10060
F1 (3x380-480V/50-60Hz):	51	4.00 A Midget fuse 10x38 (3 pcs)	RFC10040
F2 (3x200-240V/50-60Hz):	52	6.00 A Midget fuse 10x38 (3 pcs)	RFC10060
F2 (3x380-480V/50-60Hz):	52	4.00 A Midget fuse 10x38 (3 pcs)	RFC10040
F3 (3x200-240V/50-60Hz):	53	6.00 A Midget fuse 10x38 (3 pcs)	RFC10060
F3 (3x380-480V/50-60Hz):	53	4.00 A Midget fuse 10x38 (3 pcs)	RFC10040

Exotom-100
Spare Parts and Diagrams

Drawing


Pos.


Chain Spanner for 12 mm T-slots EXOCS

Cat no.

1	Chain Spanner without key	2GR26910
2	Spring for EXOCS, w. screw	2GR26023
3	Chain, 500 mm	2GR26925
4	Operating key for EXOCS (12 mm T-slots)	2GR26031
5	Anchor block	2GR26930

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


- Inputs:**
- 1x strain gauge for feed-force (analog)
 - 1x current ratio for main motor current (analog)
 - 1x potentiometer for γ-pos. sensing (analog)
 - 1x Pressure for main motor temp. (analog)
 - 1x incremental encoder for feed pos. and speed
 - 1x Namur sensor for flange in contact with workpiece
 - 1x Namur sensor for tank level
 - 1x Namur sensor for flush pistol in holder
 - 1x switch in rec. pump motor for overload sensing
 - 1x switch in excit. motor for overload sensing
 - 2x switch for emergency stop sensing
 - 2x switch for door lock and safety relay supervision
- Outputs:**
- 1x 24V/ac contactor for main motor (K1)
 - 3x 24V/ac contactors for star/delta motorstarter (K2-K4)
 - 1x 24V/ac contactor for rec. water pump (KS)
 - 4x 24V/ac contactors for excit. motor (K6)
 - 4x 24V/ac signal lamps
 - 1x 24V/ac door lock
 - 1x 24V/ac solenoid valve for flushing
 - 1x 5A stepper motor driver for feed
 - 1x 5A dc-motor driver for γ-cut


Synergi A/S Svangs Allé 176 DK-2600 Slagelse Denmark Mail: info@synergi.dk Web: www.synergi.dk	
EXOTOM-100, BLOCK DIAGRAM	
Scale	Sheet
150430050	1 of 1
Scale	Sheet
Monday, July 05, 1999	FTH/FTH

T1 600 VA

200-210V MODEL 460-480V MODEL 380-420V MODEL 220-240V MODEL

IEC742 Safety transformer

COLORCODES:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GR = GREY
 WH = WHITE

Rev.A: (FTH 04-10-99)
 XS7 added

Source: X35
 Vahnele Asia 1 76
 Vahnele Asia 1 76
 Denmark
 Relationship: 1-65-3670-3690

EXOTOM-100 CONTROL VOLTAGE DIAGRAM

Scale	1:1	Sheet	1	of	5
Scale	1:1	Sheet	1	of	5
Scale	1:1	Sheet	1	of	5

Tuesday, October 05, 1999
 FTH / FTH

15043110A

COLOR CODES:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE

EMERGENCY STOP

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac


From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac

From F7 (page 1) U4: 24Vac
 From X24 (page 1) U4: 0Vac
 From F6 (page 1) U3: 24Vac
 From X23 (page 1) U3: 0Vac


- Outputs:
- 1x 24Vac contactor for main motor (K1)
 - 3x 24Vac contactors for star/delta motorstarter (K2-K4)
 - 1x 24Vac contactor for rec. water pump (K5)
 - 1x 24Vac contactor for excit motor (K6)
 - 1x 24Vac solenoid valve for flushing (Y1)
 - 1x 24Vac solenoid valve for flushing
 - 1x 24Vac stepmotor driver for feed
 - 1x 5A dc-motor driver for y-cut

Rev. A: (FTH 04-10-99) XS7 added	STRUBER A/S Vejlevej 176 DK-8260, Hvidebø +45 3670 3900	EXOTOM-100, CONTROL VOLTAGE DIAGRAM	24 VAC OUTPUTS, STEP MOTOR AND DC MOTOR OUTPUT
Scale	Size	Sheet	Rev.
	A2	1	1
Scale	Size	Sheet	Rev.
Tuesday, October 05, 1999	FTH / FTH	2	5

15043110A

- Inputs:
- 1x strain gauge for feed-force (analog)
 - 1x current triad for main motor current (analog)
 - 1x potentiometer for y-pos. sensing (analog)
 - 1x potentiometer for x-pos. sensing (analog)
 - 1x incremental encoder for feed-pos. and speed
 - 1x Namur sensor for wheel speed
 - 1x Namur sensor for feed in top position
 - 1x Namur sensor for flange in contact with workpiece
 - 1x Namur sensor for flush pistol in holder
 - 1x Namur sensor for tank emptying with
 - 1x switch in recirculation tank for level sensing
 - 1x switch in rec. pump motor for overload sensing
 - 2x switch for emergency stop sensing
 - 3x switch for contactor supervision
 - 2x switch for door lock and safety relay supervision

COLOR CODES: BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE

Rev. A: (FTH 04-10-99) X57 added	Size: A2 Scale:	Sheet: 3 of 5
Supplier: S Drawing No.: 15043110A Drawing Date: 04-10-99 Drawing Location: Denmark Drawing Scale: 1:1	Project Code: EXOTOM-100, CONTROL VOLTAGE DIAGRAM	Rev. 1
Date: Tuesday, October 05, 1999		Scale: 1:1

15043110A

COLOR CODES:

- BK = BLACK
- BN = BROWN
- RD = RED
- OG = ORANGE
- YE = YELLOW
- GN = GREEN
- BU = BLUE
- VT = VIOLET
- GY = GREY
- WH = WHITE

S1, EMERGENCYSTOP


S2, EMERGENCYSTOP

To A1:X12-2 (page 3)









To A1:X12-1 (page 3)

To A1:X10-1 (page X)

To A1:X10-2 (page X)

Rev. (FTH 04-10-99) XS7 added		Source A/B Vejle, Aale 1.76 Vejle, Aale 1.76 Dinmark Halskøbenhavn - 46.300.9650
Size	Page Code	Page No
A2		15043110A
Scale		Sheet 4 of 5
Tuesday, October 05, 1999		FTH / FTH

COLOR CODES: BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE

Rev.A: (FTH 04-10-99) XS7 added		Source: A/S Vejlevej 179 DK-8000 Ålborg Danmark Illustration: L&S 99033659
Size	CAGE Code	EWG NO
A2		
Rev		15043110A
Sheet		5 of 5

Tuesday, October 06, 1999


FTH.FTH

CONTROL PANEL CIRCUIT


EXOTOM-100, CONTROL VOLTAGE DIAGRAM

All wires AWG 8

AWG16


For connection see other diagram

Must be short circuit protected with external fuses
Fmax: See note #1.

All wires AWG 10


VOLTAGE / FREQ. (from nameplate)	Note #1: (max. ext. fuse)	F1, F2, F3 (fuse size)	M1 (nom. power)	M1 S3, 10 min. (cutting power)	M2 (power)	Note #2: (M2-connection)
3x200V / 50Hz (CSA)	3x 80 AT	3x 6 AT	7.5 kW	10.5 kW	470 W	Delta
3x200V / 60Hz (CSA)	3x 80 AT	3x 6 AT	7.5 kW	13.5 kW	590 W	Delta
3x220-240V / 50Hz	3x 80 AT	3x 6 AT	7.5 kW	10.5 kW	480 W	Delta
3x220-240V / 60Hz	3x 80 AT	3x 6 AT	7.5 kW	10.5 kW	480 W	Delta
3x380V-415V / 50Hz	3x 50 AT	3x 4 AT	7.5 kW	10.5 kW	480 W	Star
3x380V-415V / 60Hz	3x 50 AT	3x 4 AT	7.5 kW	10.5 kW	480 W	Star
3x480V-480V / 60Hz	3x 50 AT	3x 4 AT	9.0 kW	12.6 kW	660 W	Star
3x480V / 60Hz (CSA)	3x 50 AT	3x 4 AT	9.0 kW	12.6 kW	660 W	Star


System A/E: 4.76
 Wire size corrected, CSA models added
 Date: 10/11/2009
 Drawing No.: 443_3070_3350

Rev: A: FTH-01-10-09
 Scale: 1:1
 Date: Monday, October 11, 1999



EXOTOM-100, MAIN SUPPLY CIRCUIT
 DWG NO: 15043100A
 Scale: 1:1

PAGE 2


PAGE 2


SCREW HOLES

Rev.A: Prototype core 481(05/09 PH)	STP/07/RS 2/5
	V. Højer, A. M. T. Rø
	Design: Rørdorpe
	Dokument: 15043210
	Abt.no: 43.3020.9500
EXOTOM-100 MACHINE CONTROL BOARD A1	
POWER SUPPLY	
Scale	1 A
DATE	15043210
Sheet	1 of 8
Tuesday, August 01, 2000	FTH/ETH

CPU AND MEMORY CIRCUIT


ANALOG/DIGITAL AND DIGITAL/ANALOG CONVERTERS


LINKS TO OTHER PAGES


RS-485 INTERFACE


RESET CIRCUIT


DECOUPLING

RevA: Prototype case (01/07/01/PT)	Sheet	4	of	8
Sheet No: 179 Dwg No: 15043210 Part No: 15043210	Scale	1:1	Rev	A
EXOTOM-100, MACHINE CONTROL BOARD 1 DC-MOTOR CIRCUIT	Size	A2	DWG NO	15043210
Rev: A Date: 01/07/01 Author: [Blank]	Scale	1:1	Rev	A

N14 IR2104 15V 100nF 332K	N15 IR2104 15V 100nF 332K	N13 TL084CN 15V 100nF 100K	N14 TL084CN 15V 100nF 100K	N15 TL084CN 15V 100nF 100K	C18 120uF 50V Low ESR	C19 100nF	C20 100nF	C21 100nF	C22 220nF	C23 100nF	C24 100nF	C25 100nF	C26 100nF	C27 100nF	C28 100nF	C29 100nF	C30 100nF	C31 100nF	C32 100nF	C33 100nF	C34 100nF	C35 100nF	C36 100nF	C37 100nF	C38 100nF	C39 100nF	C40 100nF	C41 100nF	C42 100nF	C43 100nF	C44 100nF	C45 100nF	C46 100nF	C47 100nF	C48 100nF	C49 100nF	C50 100nF	C51 100nF	C52 100nF	C53 100nF	C54 100nF	C55 100nF
---------------------------------------	---------------------------------------	--	--	--	-----------------------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------	--------------

Suresh A/S - Rev A, Prototype card, mail 06/07/99P114
 DK-2110, Redwood
 København - 4453, 30/03/9900.

EXOTOM-100 MACHINE CONTROL BOARD A1
STEP MOTOR CIRCUIT

Size: A2
 Scale: 1:1
 DWG NO: 15043210
 Rev: I A
 Sheet: 5 of 8

Tuesday, August 01, 2000

+30V
 +5V
 GND

C43 100n
 C44 100n
 C45 120pF/5V
 Low ESR
 C46 100n
 C47 100n
 C48 100n
 C49 100n
 C50 100n

+1.5V
 GND

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
R16 0R22/3W
 R17 0R22/3W
 R18 1K16
 R19 121R
 R20 100R
 R21 100R
 R22 100R
 R23 100R
 R24 100R
 R25 100R
 R26 100R
 R27 100R
 R28 100R
 R29 100R
 R30 100R
 R31 100R
 R32 100R
 R33 100R
 R34 100R
 R35 100R
 R36 100R
 R37 100R
 R38 100R
 R39 100R
 R40 100R
 R41 100R
 R42 100R
 R43 100R
 R44 100R
 R45 100R
 R46 100R
 R47 100R
 R48 100R
 R49 100R
 R50 100R
 R51 100R
 R52 100R
 R53 100R
 R54 100R
 R55 100R
 R56 100R
 R57 100R
 R58 100R
 R59 100R
 R60 100R
 R61 100R
 R62 100R
 R63 100R
 R64 100R
 R65 100R
 R66 100R
 R67 100R
 R68 100R
 R69 100R
 R70 100R
 R71 100R
 R72 100R
 R73 100R
 R74 100R
 R75 100R
 R76 100R
 R77 100R
 R78 100R
 R79 475R
 R80 121R
 R81 1K16
 R82 100R
 R83 100R
 R84 100R
 R85 100R
 R86 100R
 R87 100R
 R88 100R
 R89 100R
 R90 100R
 R91 100R
 R92 100R
 R93 100R
 R94 100R
 R95 100R
 R96 100R
 R97 100R
 R98 100R
 R99 100R
 R100 100R

V16 BC637
 V17 BC637
 V18 SUPRN06-18
 V19 SUPRN06-18
 V20 SUPRN06-18
 V21 SUPRN06-18
 V22 SUPRN06-18
 V23 SUPRN06-18
 V24 SUPRN06-18
 V25 SUPRN06-18
 V26 SUPRN06-18
 V27 BC637
 V28 SUPRN06-18
 V29 SUPRN06-18
 V30 SUPRN06-18
 V31 SUPRN06-18
 V32 SUPRN06-18
 V33 SUPRN06-18
 V34 SUPRN06-18
 V35 SUPRN06-18
 V36 SUPRN06-18


CL14 DSS306
 CL15 DSS306
 CL16 DSS306
 CL17 DSS306

INH1
 INH2
 INH3
 INH4
 SENSE1
 SENSE2
 SENSE3
 SENSE4

A
 B
 C
 D

Rev A: Prototype contr. board (05.07.00 FTH)	Drawn by: S. B. 1.76	EXCISE Code	DWG NO	Rev
	DK-2010 (Redrawn)	<Chipp Code>	15043210	1 A
	1481832000 - 4151_367038000	Scale	Sheet	6 of 6
EXOTOM-100 MACHINE CONTROL BOARD A1				
Transfery, August 01, 2000				
FTH / FTH				

page 2

Low level switch
in rec. bank

Overload switch
in rec.pump motor

Overload switch
in rec.pump motor

Emergency stop
switch (left side)


Emergency stop
switch (right side)

Status from
door lock switch


Z-POSITION
(INCR. IMPULSES)

Rev-A: Prototypa come del 05/07/99 FTH		Versioni: 05	
D:\K-2010\Problemi		V:\Materiale	
I:\MATERIALE\415_3670\3670		D:\K-2010\Problemi	
EXOTOM-100, MACHINE CONTROL BOARD.A1		I:\MATERIALE\415_3670\3670	
Sev	7.4	Stamp	7.4
Stamp	15/04/2010	Stamp	15/04/2010
Tuesday, August 01, 2000		FTH / FTH	
Sheet		7 of 8	


ANALOG INPUTS


INPUT FOR WHEEL SPEED DETECTION

VREF BUFFER


DECOUPLING

Rev A: Prototype come out 06.07.08 (PTH)		REV	
DESIGNED BY: N. V. P. (A18, 1.76)		DWG NO: 15043210	
DRAWN BY: N. V. P.		DATE: 06.07.08	
MANUFACTURED BY: 485_3102_3800		PTH/PTH	
Sheet	8	of	8

EXOTOM-100, MACHINE CONTROL BOARD A1

Tuesday, August 01, 2000


STRUERS A/S
 Valbøjs Allé 176
 DK-2610 Rødovre
 Denmark
 +45 3670 3500


Exotom-100, LED Circuit board

Size	CAGE Code	DWG NO	Rev
A4	<Cage Code>	15043220	
Scale		FTH / FTH	Sheet 1 of 1
Tuesday, August 01, 2000			

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







Herewith declares that

<i>Product Name:</i>	Exotom-100
<i>Type No:</i>	504
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN292-1:1991, EN292-2:1991, EN418:1992, EN954-1:1996, EN1050:1996, EN60204-1:1997.

EMC-Directive 89/336/EEC, 92/31/EEC and 93/68/EEC according to the following standard(s):
EN50081-1:1994, EN50082-1:1997.

Low Voltage Directive 73/23/EEC and 93/68/EEC according to the following standard(s):
EN60204-1:1997.

Supplementary Information The equipment also complies with the following standards:
IEC204-1, UL508:1988, UL508A:1993, NFPA70:1996, NFPA79:1997
and FCC part 15, subpart B.

The above has been declared according to the global method, module A

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Exotom-100

**Automatic cut-off
machine for cutting of
workpieces up to 160 mm
in diameter.
Designed for the
production environment**

A revolution in cutting technology

- **Extremely simple operation**
Just two operating keys
No programming required
- **Unprecedented capacity and accessibility**
- **Up to four different cutting modes**

Unique technologies and features:


- ExciCut** - easy through the hardest workpieces
- OptiFeed** - intelligent feed control
- AxioCut** - cutting mode module
- Variable cutting offset** - flexible wheel positioning
- Three operation modes** - block non-essential controls

Effortless operation even wearing working gloves

ExciCut -for cutting the hardest materials

The Exotom-100 comes standard with the unique ExciCut feature. The technology was first introduced by Struers a decade ago, and appears on the

Exotom-100 in a highly refined version. Through an oscillating movement of the cut-off wheel, the contact area between the wheel and the workpiece is kept at a minimum throughout the cutting process.

ExciCut is ideal for hardened materials and means faster cutting of standard materials. Since the contact area between sample and cut-off wheel constantly changes, cooling is greatly improved. In turn, this means reduced risk of damage to the workpiece, and reduced work in the processes following the cutting.




Direct Cut -the classic cutting mode

For smaller samples and some materials it might be an advantage to use the direct cutting mode where the cut-off wheel is fed through the workpiece in a linear movement.

AxioCut -for extra large workpieces

For cutting of larger workpieces, additional 150 mm cutting capacity is available with AxioCut. This new feature allows cutting of workpieces with a dimension of e.g. 100 x 450 mm (for


other dimensions, please refer to Cutting Diagram on the back of this brochure). The field-upgradable option includes two different cutting modes: Step and Sweep.

Step cutting offers very fast cutting, whereas Sweep cutting offers fast cutting combined with better wheel economy.

OptiFeed -Protection and speed

Cutting of workpieces with uneven cross-section or composition, represents a great challenge to the operator. The risk of damage to the workpiece or the cut-off wheel is high. Bearing this in mind, Struers has developed the new OptiFeed principle, based on a preset maximum force limit.


As the contact area (between workpiece and wheel, cf. below diagram) is getting larger, the


force increases to maintain the preset feed speed. If the preset force limit is reached, the feed speed is automatically reduced to avoid damage.

When the cut-off wheel has passed the centre of the workpiece, the contact area is getting smaller again, and with it the resulting force. In turn, the feed speed is automatically increased to the preset maximum.



Thanks to OptiFeed you will always cut through the workpiece at the fastest possible speed, within the set parameters - and without overloading the workpiece and the cut-off wheel.

Two flexible water jets for additional cooling. Ideal for hollow workpieces

*Single-level display:
Full control of cutting
process + motor information
in same view*

Controls of Unique Simplicity

Two simple keys combined with a large graphic display are incorporated on a swivelling control panel. With the multi-function knob all cutting parameters are easily controlled even wearing working gloves. The joystick is used for fast and precise positioning of the cut-off wheel. The side-mounted control panel is kept away from the machine for maximum safety and to avoid soiling of display and controls. Also the risk of slipping samples damaging the display is eliminated. The control panel can easily be moved if very long workpieces have to be placed in in the cutting chamber.

Large Cutting Chamber

The corrosion resistant cutting chamber is exceptionally spacious with its 855 x 550 mm cutting table, offering unrestricted access. Bulky workpieces are easily placed manually or using a forklift truck, since the cutting table is the highest point in the chamber. Maximum mounting flexibility is provided by multiple, bi-directional 12 mm T-slots. The powerful, built-in flush hose provides for easy cleaning of the cutting chamber. Extra large windows and fully, indirect lightning, give an excellent, non-glare view of the cutting chamber. A balanced mechanism ensures effortless opening and closing of the hood.

Variable Cutting Offset

The horizontal position (cutting offset) of the cut-off wheel can be set manually. In the ExciCut and Direct Cut cutting modes, the cut-off wheel may be moved 10 mm towards the front and 60 mm towards the

back of the cutting chamber. This allows for maximum flexibility when cutting difficult workpieces (e.g. turbine blades) and for clamping generally.

Three Operation Modes

Particularly in an industrial environment with many users, operator errors may cause damage to machine and samples. To overcome this, non-essential controls may be blocked partly or fully by using the pass code-protected operation modes.

High-capacity Cooling

The workpiece and cut-off wheel are effectively cooled via a large capacity (150 litres) cooling unit. For extreme cases, additional cooling may be added, using the two flexible water jets provided. This is particularly useful when cutting hollow workpieces. The cooling unit comes standard with an automatic emptying function and a

warning on low coolant level. A special filter tray is provided to collect cutting debris, greatly reducing the need for cleaning of the recirculation tank.

AutoStop

After the workpiece is securely clamped and Start has been pressed you can leave the machine unattended. OptiFeed makes sure that cutting is completely under control. With AutoStop, Exotom-100 automatically detects when the workpiece has been cut through, returns the wheel to its start position and stops all motors.

QuickPosition

Positioning the cut-off wheel can be sped up dramatically. Simply push the joystick downwards and the cut-off wheel will speed towards the workpiece. On contact with the workpiece, the cut-off wheel stops automatically without harming workpiece or wheel.

*Exceptionally spacious and accessible cutting chamber:
Even the largest workpieces are easily placed on the cutting table*

Rubber curtains mean trouble-free cutting of very long workpieces. Move the swivelling control panel for optimum access

After contact, the wheel is automatically retracted a little from the workpiece, to be ready for cutting. QuickPosition can be used both before and during cutting.

Cutting Long Workpieces

If the user wants to cut workpieces protruding the cutting chamber, it is possible to swivel the control panel aside, for easy access. The cutting chamber comes standard with side-mounted rubber curtains, for protruding workpieces.

Programmable Return Position

The return position of the cut-off wheel can be set to three different positions: Exotom-100 retracts the cut-off wheel to the top position, to the position where the cutting was started or it lets the cut-off wheel stay in the cut. In this way, it is very easy to get convenient access to the workpiece or an extremely quick exchange when cutting uniform workpieces.

Electronic Monitoring

More than a dozen electronic sensors monitor the complete cutting process and all machine functions to ensure flawless cutting. Never has sectioning on a large cut-off machine been this safe.

Safety Features

Exotom-100 complies with international safety standards, and safety features include emergency stops and safety sensors for all machine functions and motors. The cutting motor cannot be started before the protection hood is closed, and after Start, the hood cannot be opened before the machine has stopped. The windows of the protection hood are made of impact resistant material.



Cut-off wheels

A range of cut-off wheels for all ferrous and non-ferrous metals is available for Exotom-100.

Diamond cut-off wheels are available for cutting of ceramic and mineralogical materials and sintered carbides.

Application	Dimensions	Abrasive	Code
For ExciCut and direct cutting of extremely hard steels (350-800 HV)	432 x 3.0 x 32 mm	Al ₂ O ₃	101MA
For ExciCut and direct cutting of steels (< 500 HV). Extremely hard (white cast iron) or ductile (18/8 or St60) steel types are cut with oscillating cutting. Gives good wheel economy	432 x 3.0 x 32 mm	Al ₂ O ₃	102MA
For ExciCut and direct cutting of very hard and ductile ferrous metals. Strong, fibre-reinforced wheel (250-700 HV)	432 x 3.0 x 32 mm	Al ₂ O ₃	202MA
For ExciCut and direct cutting of steels (< 500 HV). Gives very good wheel economy	432 x 3.0 x 32 mm	Al ₂ O ₃	104MA
For ExciCut and direct cutting of non-ferrous metals (< 300 HV). Primarily for ExciCut	432 x 3.0 x 32 mm	SiC	106MA
For direct cutting of hardened steels	406 x 1.8 x 32 mm	CBN	108MA
For direct cutting of sintered carbides and ceramics	350 x 1.8 x 32 mm 305 x 1.8 x 32 mm	Diamond Diamond	26EXO 25EXO

Technical data

Subject		Metric/International	US
CUTTING SPECIFICATIONS			
Workpiece Dimensions (max.)	<i>Workpiece inside cutting chamber:</i>		
	Height	250 mm	10"
	Width	650 mm	25.6"
	Depth	550 mm	20"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	140 mm	5.5"
Depth	360 mm	14.2"	
Cutting Capacity*	Max. cutting diameter	160 mm	6.3"
	Rectangular sample:**		
	(In brackets with AxioCut option)		
	Height	100 mm (100 mm)	3.9" (3.9")
Depth	300 mm (450 mm)	11.8" (17.7")	
PHYSICAL SPECIFICATIONS			
Cutting Motor	Cutting power	10.5 kW	14 HP
	Maximum power	15-18 kW	20-24 HP
Cut-off Wheel	Diameter x Thickness x Centre-hole	432 x 3 x 32 mm	17 x 0.12 x 1.26"
	Rotational speed (running idle)	1950 rpm	1950 rpm
Positioning & Feed	Positioning range (of cut-off wheel)	0 – 250 mm	0 – 9.8"
	Max. positioning speed	50 mm/s	2.0"/s
	Feed Speed range	0.05 – 5 mm/s	0.002 – 0.2"/s
	(adjustable in steps of)	(0,05 mm/s)	(0.002"/s)
Cutting force	Max. 700 N	Max. 150 lbf	
Cutting table	Width	855 mm	33.6"
	Depth	550 mm	21.6"
	T-slots	12 mm	0.48"
Dimensions and Weight	Height (closed/open hood)	1756/2300 mm	69/90.5"
	Width (ex/including control panel)	1050/1350 mm	41.3"/53.1"
	Depth	1500 mm	59.1"
	Weight	820 kg	1810 lbs
Recirculation Cooling Unit	Tank volume	150 l	39.6 gallons
	Approx. flow	25 l/min	6.6 gallons/min
ENVIRONMENT			
Noise level	78 dB(A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.		
ELECTRICAL DATA			
Voltage/Frequency	200-240 V / Max. power 15 kW / Max. load 69 A		
	380-480 V / Max. power 18 kW / Max. load 34 A		

*Projected cutting capacity with new cut-off wheel. Actual cutting capacity depends on sample material, cut-off wheel and clamping technique

**Example of cutting capacity. Please refer to Cutting Diagram above

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

Specifications	Code
Exotom-100, AxioCut Automatic cut-off machine with OptiFeed, AxioCut, ExciCut and Direct Cut. 10.5 kW (14 HP) motor. For 432 mm cut-off wheels. With recirculation cooling unit (capacity 150 l) and cutting table with 12 mm T-slots. Quick-clamping devices (EXOLE + EXORI) and vertical clamping system (EXOVS) respectively are ordered separately.	EXOCM
Exotom-100 Automatic cut-off machine with OptiFeed, ExciCut and Direct Cut. 10.5 kW (14 HP) motor. For 432 mm cut-off wheels. With recirculation cooling unit (capacity 150 l) and cutting table with 12 mm T-slots. Quick-clamping devices (EXOLE + EXORI) and vertical clamping system (EXOVS) respectively are ordered separately. Field-upgradable to AxioCut by optional kit (EXOMA)	EXOCD
Accessories	Code
AxioCut Kit for Exotom-100 Increases cutting capacity (depth) by 150 mm. Requires installation by Struers Service Technician.	EXOMA
Replacement Stainless Steel Bands For Cutting Table with 12 mm T-slots. Available in three sizes, each supplied in sets of two: 100 x 225 mm 100 x 300 mm 60 x 548 mm	EXOBS EXOBM EXOBL
Quick-clamping Device for 12 mm T-slots, Left For securing the workpiece. Complete with backstop. To be mounted on a cutting table and placed to the left of the cut-off wheel.	EXOLE
Quick-clamping Device for 12 mm T-slots, Right For securing the workpiece. Complete with backstop. To be mounted on a cutting table and placed to the right of the cut-off wheel.	EXORI
Clamping Tools for 12 mm T-slots. For clamping irregularly shaped workpieces. Complete with clamps, supports and bolts.	MAGOF
Vertical Clamping System for 12 mm T-slots For clamping irregularly shaped workpieces. Clamping height up to 92 mm. Complete with operating key and one flat clamping shoe.	EXOVS
Swivel Shoes for Vertical Clamping System. For clamping irregularly shaped workpieces. Set of four multi-shaped swivel shoes.	TREVI
Riser Block for 12 mm Vertical Clamping System For elevating Vertical Clamping System when clamping high workpieces. To increase clamping height by 74 mm.	EXOKS
Arm Extension for 12 mm Vertical Clamping System For extending arm of Vertical Clamping System.	EXOAR
Chain Spanner For 12 mm T-slots. Complete with anchor block and operating key.	EXOCS
Spare Nylon Filter for Recirculation Cooling Unit For insertion into filter-drawer in cooling unit compartment.	EXOFI
Bandfilter connection kit	EXOCO

Clamping Tools

Struers offer a wide range of clamping tools for Exotom-100. Clamping requirements differ greatly, which is why Struers have chosen to supply Exotom-100 without any 'standard' clamping devices. Rather, the user can order exactly the clamping system matching his requirements. Please refer to separate brochure.

Large cylindrical and uneven samples (e.g. crankshafts) are easily clamped with the Struers chain spanner

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Struers Ltd.
Struers Inc.