

December 2014

artus[®] HSV-1/2 RG PCR Kit Handbook



24 (catalog no. 4500263)



96 (catalog no. 4500265)

Version 1



Qualitative in vitro diagnostics

For use with Rotor-Gene[®] Q instruments



4500263, 4500265



1060171EN



QIAGEN GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY

R2

MAT

1060171EN



Sample & Assay Technologies

QIAGEN Sample and Assay Technologies

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

QIAGEN sets standards in:

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit www.qiagen.com.

Contents

Kit Contents	4
Symbols	4
Storage	5
Intended Use	5
Product Use Limitations	6
Technical Assistance	6
Quality Control	6
Safety Information	7
Introduction	8
Principle	8
Pathogen information	8
Performance characteristics	9
Equipment and Reagents to Be Supplied by User	17
Important Notes	18
General precautions	18
DNA isolation	18
Internal control	19
Protocol	
■ PCR and Data Analysis	20
Troubleshooting Guide	30
References	33
Ordering Information	34

Kit Contents

artus HSV-1/2 RG PCR Kit		(24)	(96)
Catalog no.		4500263	4500265
Number of reactions		24	96
Blue	HSV-1/2 RG Master	2 x 300µl	8 x 300µl
Yellow	HSV-1/2 RG Mg-Sol* Mg-Sol	600 µl	600 µl
Red	HSV-1 RG PC [†] (100 cop/µl)	200 µl	200 µl
Brown	HSV-2 RG PC [†] (100 cop/µl)	200 µl	200 µl
Green	HSV-1/2 RG IC [‡] IC	1000 µl	2 x 1000 µl
White	Water (PCR grade)	1000 µl	1000 µl
Handbook 		1	1

* Magnesium solution.

[†] Positive control.

[‡] Internal control.

Symbols

 Σ <N> Contains reagents for <N> tests



Use by

IVD

In vitro diagnostic medical device

REF

Catalog number

LOT

Lot number

MAT


Material number

COMP

Components

CONT


Contains



Number



Global Trade Item Number



Temperature limitations



Legal manufacturer



Refer to information given in the handbook



Important note

Storage

The components of the *artus* HSV-1/2 RG PCR Kit should be stored at -15°C to -30°C and are stable until the expiration date stated on the label. Repeated thawing and freezing (>2 x) should be avoided, as this may reduce assay sensitivity. If the reagents are to be used only intermittently, they should be frozen in aliquots. Storage at $2-8^{\circ}\text{C}$ should not exceed a period of 5 hours.

Intended Use

The *artus* HSV-1/2 RG PCR Kit is a real-time based polymerase chain reaction (PCR) assay for the detection and discrimination of human herpes simplex virus 1 and 2 DNA on the Rotor-Gene Q instruments after fully automated purification of cerebrospinal fluid (CSF) samples from HSV infected individuals using the EZ1[®] DSP Virus Kit.



The *artus* HSV-1/2 RG PCR Kit may not be used with Rotor-Gene Q 2plex Instruments.

The *artus* HSV-1/2 RG PCR Kit is intended for use in conjunction with clinical presentation and other laboratory markers for disease prognosis.

Product Use Limitations

All reagents may exclusively be used in in vitro diagnostics.

The product is to be used by personnel specially instructed and trained in the in vitro diagnostics procedures (EN375) only.

Strict compliance with the user manual is required for optimal PCR results.

Attention should be paid to expiration dates printed on the box and labels of all components. Do not use expired components.

Technical Assistance

At QIAGEN, we pride ourselves on the quality and availability of our technical support. Our Technical Service Departments are staffed by experienced scientists with extensive practical and theoretical expertise in sample and assay technologies and the use of QIAGEN® products. If you have any questions or experience any difficulties regarding the *artus* HSV-1/2 RG PCR Kit or QIAGEN products in general, please do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information, please see our Technical Support Center at www.qiagen.com/Support or call one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit www.qiagen.com).

Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of *artus* HSV-1/2 RG PCR Kit is tested against predetermined specifications to ensure consistent product quality.

Warnings and Precautions

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.

Discard sample and assay waste according to your local safety regulations.

Introduction

The *artus* HSV-1/2 RG PCR Kit constitutes a ready-to-use system for the detection of HSV-1 and HSV-2 DNA using polymerase chain reaction (PCR) on Rotor-Gene Q instruments. The HSV-1/2 RG Master contains reagents and enzymes for the specific amplification of a 154 bp region of the HSV-1 and HSV-2 genomes, and for the direct detection of the specific amplicon in fluorescence channel Cycling Green (source 470 nm, detector 510 nm) and Cycling Orange (source 585 nm, detector 610 nm) of the Rotor-Gene Q instruments.

In addition, the *artus* HSV-1/2 RG PCR Kit contains a second heterologous amplification system to identify possible PCR inhibition. This is detected as an internal control (IC) in fluorescence channel Cycling Yellow (source 530 nm, detector 555 nm) of the Rotor-Gene Q instruments. The detection limit of the analytical HSV-1/2 RG PCR (see "Analytical sensitivity", page 9) is not reduced. External positive controls (HSV-1 RG PC and HSV-2 RG PC) are supplied.

Principle

Pathogen detection by the polymerase chain reaction (PCR) is based on the amplification of specific regions of the pathogen genome. In real-time PCR the amplified product is detected via fluorescent dyes. These are usually linked to oligonucleotide probes that bind specifically to the amplified product. Monitoring the fluorescence intensities during the PCR run (i.e., in real-time) allows the detection and quantitation of the accumulating product without having to re-open the reaction tubes after the PCR run.*

Pathogen information

Herpes simplex virus (HSV) is found in lesion fluids, saliva, cerebrospinal fluid (CSF), and vaginal secretions. It is transmitted primarily by direct contact with lesions and via sexual intercourse, as well as perinatally. Lesions on the skin and mucous membranes of the mouth and genitals characterize most HSV positive cases. HSV infection can be either primary (> 90 % of these cases are asymptomatic) or recurrent (secondary). Primary infection with HSV-1 can lead to, among others, gingivostomatitis, eczema herpeticum, keratoconjunctivitis and encephalitis; primary HSV-2 infection occurs as, among others, vulvovaginitis, meningitis and generalized herpes in newborns. The primary symptoms of a secondary infection are skin lesions in the nose, mouth and genital regions. Even more severe are the recurrent forms of keratoconjunctivitis and meningitis.

* Mackay, I.M. (2004) Real-time PCR in the microbiology laboratory. Clin. Microbiol. Infect. **10**, 190.

Performance characteristics

Analytical sensitivity

To determine the analytical sensitivity of the *artus* HSV-1/2 RG PCR Kit, a standard dilution series was set up from 10 to 0.001 copies/ μ l and analyzed on the Rotor-Gene Q/6000 in combination with the *artus* HSV-1/2 RG PCR Kit. Testing was carried out on 3 different days on 8 replicates. The results were determined by a probit analysis. The analytical detection limit of the *artus* HSV-1/2 RG PCR Kit in combination with the Rotor-Gene Q/6000 is consistently 0.12 copies/ μ l ($p = 0.05$) for HSV-1 and 0.16 copies/ μ l ($p=0.05$) for HSV-2. This means that there is a 95% probability that 0.12 copies/ μ l of HSV-1 DNA or 0.16 copies/ μ l of HSV-2 DNA will be detected. A graphical illustration of the probit analysis for HSV-1 is shown in Figure 1 below; the diagram of the probit analysis for HSV-2 is shown in Figure 2




Figure 1. Probit analysis: HSV-1 (Rotor-Gene Q/6000). Analytical sensitivity for HSV-1 of the *artus* HSV-1/2 RG PCR Kit on the Rotor-GeneQ/ 6000.




Figure 2. Probit analysis: HSV-2 (Rotor-Gene Q/6000). Analytical sensitivity for HSV-2 of the *artus* HSV-1/2 RG PCR Kit on the Rotor-Gene Q/6000.

Specificity

The specificity of the *artus* HSV-1/2 PCR Kit is first and foremost ensured by the selection of the primers and probes, as well as the selection of stringent reaction conditions. The primers and probes were checked for possible homologies to all published sequences in gene banks by sequence comparison analysis. The detectability of all relevant genotypes has thus been ensured by a database alignment and by a PCR run on Rotor-Gene instruments with the strains listed in Table 1.

Moreover, the specificity was validated with 30 different HSV-1 and HSV-2 negative CSF samples. These did not generate any signals with the HSV-1 and HSV-2 specific primers and probes, which are included in the HSV-1/2 RG Master.

A potential cross-reactivity of the *artus* HSV-1/2 RG PCR Kit was tested using the control group listed in Table 2. None of the tested pathogens has been reactive.

Table 1. Testing of the specificity of relevant genotypes

Virus	Strain	Source	HSV-1 (Cycling Green)	HSV-2 (Cycling Orange)	Internal control (Cycling Yellow)
HSV-1	HF	ATCC*	+	–	+
HSV-1	KOS	INSTAND†	+	–	+
HSV-1	MacIntyre	QCMD‡	+	–	+
HSV-2	HG-52	NCPV§	–	+	+
HSV-2	G	ATCC*	–	+	+
HSV-2	MS	QCMD‡	–	+	+

* ATCC American Type Culture Collection.

† INSTAND Society for Promotion of Quality Assurance in the Medical Laboratories.

‡ QCMD Quality Control for Molecular Diagnostics.

§ NCPV National Collection of Pathogenic Viruses.

Table 2. Testing the specificity of the kit with potentially cross-reactive pathogens

Control group	HSV-1 (Cycling Green)	HSV-2 (Cycling Orange)	Internal control (Cycling Yellow)
Human herpesvirus 3 (varicella-zoster virus)	–	–	+
Human herpesvirus 4 (Epstein-Barr virus)	–	–	+
Human herpesvirus 5 (cytomegalovirus)	–	–	+
Human herpesvirus 6A	–	–	+
Human herpesvirus 6B	–	–	+
Human herpesvirus 7	–	–	+
Human herpesvirus 8 (Kaposi's sarcoma-associated herpesvirus)	–	–	+
Hepatitis A virus	–	–	+
Hepatitis B virus	–	–	+
Hepatitis C virus	–	–	+
Human immunodeficiency virus (HIV)	–	–	+
Human T cell leukemia virus 1	–	–	+
Human T cell leukemia virus 2	–	–	+
Enterovirus	–	–	+
Parvovirus B19	–	–	+
West Nile virus	–	–	+

Precision

The precision data of the *artus* HSV-1/2 RG PCR Kit have been collected by means of Rotor-Gene instruments and allow the determination of the total variance of the assay. The total variance consists of the intra-assay variability (variability of multiple results of samples of the same concentration within one experiment), the inter-assay variability (variability of multiple results of the assay generated on different instruments of the same type by different operators within one laboratory) and the inter-batch variability (variability of multiple results of the assay using various batches). The data obtained were used to determine the standard deviation, the variance, and the coefficient of variation for the pathogen specific and the internal control PCR.

Precision data of the *artus* HSV-1/2 RG PCR have been collected using HSV-1 and HSV-2 DNA with the concentration of 10 copies/ μ l. Testing was performed with 8 replicates. The precision data were calculated on basis of the C_T values of the amplification curves (C_T : threshold cycle, see Table 3 and Table 4). Based on these results, the overall statistical spread of any given sample with the mentioned concentration is 1.82% (C_T) for HSV-1, 0.67% (C_T) for HSV-2, and 1.24% (C_T) and 1.58% (C_T) respectively for the detection of the internal control. These values are based on the totality of all single values of the determined variability.

Table 3. Precision data for HSV-1 on basis of the C_T values

	C_T value	Standard deviation	Coefficient of variation (%)
Intra-assay variability: HSV-1 10 copies/ μ l	30.46	0.25	0.81
Intra-assay variability: Internal control	25.29	0.08	0.3
Inter-assay variability: HSV-1 10 copies/ μ l	29.69	0.69	2.05
Inter-assay variability: Internal control	24.97	0.31	1.25
Inter-batch variability: HSV-1 10 copies/ μ l	29.95	0.40	1.35
Inter-batch variability: Internal control	24.90	0.30	1.20
Total variance: HSV-1 10 copies/ μ l	29.91	0.55	1.82
Total variance: Internal control	24.99	0.31	1.24

Table 4. Precision data for HSV-2 on basis of the C_T values

	C_T value	Standard deviation	Coefficient of variation (%)
Intra-assay variability: HSV-2 10 copies/ μ l	29.85	0.15	0.50
Intra-assay variability: Internal control	25.17	0.39	1.55
Inter-assay variability: HSV-2 10 copies/ μ l	29.92	0.15	0.49
Inter-assay variability: Internal control	25.11	0.41	1.63
Inter-batch variability: HSV-2 10 copies/ μ l	29.80	0.23	0.79
Inter-batch variability: Internal control	24.89	0.33	1.32
Total variance: HSV-2 10 copies/ μ l	29.88	0.20	0.67
Total variance: Internal control	25.07	0.40	1.58

Robustness

The verification of the robustness allows the determination of the total failure rate of the *artus* HSV-1/2 RG PCR Kit. To obtain very low virus titers of HSV-1 and HSV-2, 30 negative samples of CSF were spiked with 0.36 copies/ μ l elution volume of HSV-1 or with 0.48 copies/ μ l elution volume of HSV-2 DNA (threefold concentration of the analytical sensitivity limit). After extraction using the EZ1 DSP Virus Kit, these samples were analyzed with the *artus* HSV-1/2 RG PCR Kit. All 30 samples were assessed correctly as weak positives for each HSV type resulting in a failure rate of 0%. In addition, the robustness of the internal control was assessed by purification and analysis of 30 HSV-1 and HSV-2 negative CSF samples. No inhibition of the PCR was detected resulting in a total failure rate of 0%. Thus, the robustness of the *artus* HSV-1/2 RG PCR Kit is $\geq 99\%$.

Reproducibility

Reproducibility data permit a regular performance assessment of the *artus* HSV-1/2 RG PCR Kit as well as an efficiency comparison with other products. These data are obtained by the participation in established proficiency programs.

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- DNA isolation kit (see “DNA isolation”, page 18)
- Pipets (adjustable)*
- Sterile pipet tips with filters
- Vortex mixer*
- Benchtop centrifuge* with rotor for 2 ml reaction tubes
- Rotor-Gene Q or Rotor-Gene instrument*[†] with fluorescence channels for Cycling Green, Cycling Orange, and Cycling Yellow
- Rotor-Gene Q software version 1.7.94, and higher (Rotor-Gene 6000 software version 1.7.65, and higher)
- Strip Tubes and Caps, 0.1 ml, for use with 72-well rotor (cat. no. 981103 or 981106)
- Alternatively: PCR Tubes, 0.2 ml, for use with 36-well rotor (cat. no. 981005 or 981008)
- Cooling block (Loading Block 72 x 0.1 ml Tubes, cat. no. 9018901, or Loading Block 96 x 0.2 ml Tubes, cat. no. 9018905)

* Ensure that instruments have been checked and calibrated according to the manufacturer's recommendations.

[†] The *artus* HSV-1/2 RG PCR Kit may not be used with Rotor-Gene Q 2plex Instruments.

Important Notes

General precautions

The user should always pay attention to the following:

- Use sterile pipet tips with filters.
- Store and extract positive materials (specimens, positive controls, and amplicons) separately from all other reagents, and add them to the reaction mix in a spatially separated facility.
- Thaw all components thoroughly at room temperature (15–25°C) before starting an assay.
- When thawed, mix the components (by pipetting repeatedly up and down or by pulse vortexing) and centrifuge briefly.
- Work quickly and keep components on ice or in the cooling block (72/96-well loading block).

DNA isolation

The EZ1 DSP Virus Kit (QIAGEN, cat. no. 62724*) is validated for viral DNA purification from human CSF for use with the *artus* HSV-1/2 RG PCR Kit. Carry out the viral DNA purification according to the instructions in the *EZ1 DSP Virus Kit Handbook*.

❗ The *artus* HSV-1/2 RG PCR Kit should not be used with phenol-based isolation methods.

❗ The use of carrier RNA is critical for extraction efficiency and, consequently, for DNA yield. Add the appropriate amount of carrier RNA to each extraction following the instructions in the *EZ1 DSP Virus Kit Handbook*.

❗ The internal control of the *artus* HSV-1/2 RG PCR Kit can be used directly in the isolation procedure (see “Internal Control”, below).

* The EZ1 DSP Virus Kit is also available as CE-IVD-marked EASYartus® HSV-1/2 RG PCR Kits, combined with the *artus* HSV-1/2 RG PCR Kit (see page 34 for ordering information).

Internal control

An internal control (HSV-1/2 RG IC) is supplied. This allows the user both to control the DNA isolation procedure and to check for possible PCR inhibition. For this application, add the internal control to the isolation at a ratio of 0.1 μ l per 1 μ l elution volume. For example, using the EZ1 DSP Virus Kit, the DNA is eluted in 60 μ l Elution Buffer (AVE). Hence, 6 μ l of the internal control should be added initially.

 Do not add the internal control and the carrier RNA to the sample material directly.

The internal control can optionally be used exclusively to check for possible PCR inhibition. For this application, add the internal control directly to the mixture of HSV-1/2 RG Master and HSV-1/2 RG Mg-Sol, as described in step 2b of the protocol (page 21).

Protocol: PCR and Data Analysis

Important points before starting

- Before beginning the procedure, read “Important Notes”, page 18.
- Take time to familiarize yourself with the Rotor-Gene Q instrument before starting the protocol. See the instrument user manual.
- Make sure that the positive controls and one negative control (Water, PCR grade) are included per PCR run.

Things to do before starting

- Make sure that the cooling block (accessory of the Rotor-Gene Q instrument) is precooled to 2–8°C.
- Before each use, all reagents need to be thawed completely, mixed (by repeated up and down pipetting or by quick vortexing), and centrifuged briefly.

Procedure

1. Place the desired number of PCR tubes into the adapters of the cooling block.
2. If you are using the internal control to monitor the DNA isolation procedure and to check for possible PCR inhibition, follow step 2a. If you are using the internal control exclusively to check for PCR inhibition, follow step 2b.
Use the internal control according to step 2b for all positive and negative controls.
- 2a. The internal control has already been added to the isolation (see “Internal control”, page 18). In this case, prepare a master mix according to Table 5.

The reaction mix typically contains all of the components needed for PCR except the sample.

Table 5. Preparation of master mix (internal control used to monitor DNA isolation and check for PCR inhibition)

Number of samples	1	12
HSV-1/2 RG Master	25 μ l	300 μ l
HSV-1/2 RG Mg-Sol	5 μ l	60 μ l
HSV-1/2 RG IC	0 μ l	0 μ l
Total volume	30 μl	360 μl

- 2b. The internal control must be added directly to the mixture of HSV-1/2 RG Master and HSV-1/2 RG Mg-Sol. In this case, prepare a master mix according to Table 6.**

The reaction mix typically contains all of the components needed for PCR except the sample.

Table 6. Preparation of master mix (internal control used exclusively to check for PCR inhibition)

Number of samples	1	12
HSV-1/2 RG Master	25 μ l	300 μ l
HSV-1/2 RG Mg-Sol	5 μ l	60 μ l
HSV-1/2 RG IC	2 μ l	24 μ l
Total volume	32 μl*	384 μl*

* The volume increase caused by adding the internal control is neglected when preparing the PCR assay. The sensitivity of the detection system is not impaired.

- 3. Pipet 30 μ l of the master mix into each PCR tube. Then add 20 μ l of the eluted sample DNA (see Table 7), and mix well by pipetting repeatedly up and down. Correspondingly, 20 μ l of the HSV-1 RG PC and HSV-2 RG PC have to be used as positive controls and 20 μ l of water (Water, PCR grade) as a negative control.**

Table 7. Preparation of PCR assay

Number of samples	1	12
Master mix	30 μ l	30 μ l each
Sample	20 μ l	20 μ l each
Total volume	50 μl	50 μl each

4. Close the PCR tubes. Make sure that the locking ring (accessory of the Rotor-Gene Instrument) is placed on top of the rotor to prevent accidental opening of the tubes during the run.
5. For the detection of HSV-1 DNA or HSV-2 DNA, create a temperature profile according to the following steps.

Setting the general assay parameters	Figures 3, 4, 5
Initial activation of the hot-start enzyme	Figure 6
Amplification of the DNA	Figure 7
Adjusting the fluorescence channel sensitivity	Figure 8
Starting the run	Figure 9

All specifications refer to the Rotor-Gene Q software versions 1.7.94, and higher, Rotor-Gene 6000 software versions 1.7.65, and higher. Please find further information on programming Rotor-Gene instruments in the instrument user manual. In the illustrations these settings are framed in bold black. Illustrations are included for Rotor-Gene Q instruments.

6. First, open the “New Run Wizard” dialog box with the “Advanced” version (Figure 3). Check the “Locking Ring Attached” box and click “Next”.




Figure 3. The “New Run Wizard” dialog box.

7. Select 50 for the PCR reaction volume and click “Next” (Figure 4).




Figure 4. Setting the general assay parameters.

8. Click the “Edit Profile” button in the next “New Run Wizard” dialog box (Figure 5), and program the temperature profile as shown in Figures 6–7.




Figure 5. Editing the profile.




Figure 6. Initial activation of the hot-start enzyme.




Figure 7. Amplification of the DNA. Make sure to activate the touchdown function for 10 cycles in the Annealing step.

9. The detection range of the fluorescence channels has to be determined according to the fluorescence intensities in the PCR tubes. Click "Gain Optimisation" in the "New Run Wizard" dialog box (see , Step 2) to open the "Auto-Gain Optimisation Setup" dialog box (Figure 8). Set the calibration temperature to 65 to match the annealing temperature of the amplification program (Figure 7, Step 4b). Make sure that all three channels (Green, Orange, and Yellow) are selected for "Auto-Gain Optimisation". (Find channels in the drop down menu under "Channel Settings" and click "Add".) Click "Start" to begin the gain optimization. Click "Close" of the "Auto-Gain Optimisation Setup" dialog box when the gain calibration is completed.




Figure 8. Adjusting the fluorescence channel sensitivity.

10. The gain values determined by the channel calibration are saved automatically and are listed in the last menu window of the programming procedure (Figure 9). Click "Start Run".




Figure 9. Starting the run.

11. After the run is finished, analyze the data. The following results (11a, 11b, 11c, 11d, 11e, and 11f) are possible.

Examples of positive and negative PCR reactions are given in Figure 10, Figure 11, and Figure 12.

**11a. A signal is detected in fluorescence channel Cycling Green.
The result of the analysis is positive: the sample contains HSV-1 DNA.**

In this case, the detection of a signal in the Cycling Yellow channel is dispensable, since high initial concentrations of HSV-1 DNA (positive signal in the Cycling Green channel) can lead to a reduced or absent fluorescence signal of the internal control in the Cycling Yellow channel (competition).

**11b. In fluorescence channel Cycling Green no signal is detected. At the same time, a signal from the internal control appears in the Cycling Yellow channel.
In the sample no HSV-1 DNA is detectable. It can be considered negative.**

In the case of a negative HSV-1 PCR, the detected signal of the internal control rules out the possibility of PCR inhibition.

**11c. A signal is detected in fluorescence channel Cycling Orange.
The result of the analysis is positive: the sample contains HSV-2 DNA.**

In this case the detection of a signal in the Cycling Yellow channel is dispensable, since high initial concentrations of HSV-2 DNA (positive signal in the Cycling Orange channel) can lead to a reduced or absent fluorescence signal of the internal control in the Cycling Yellow channel (competition).

**11d. In fluorescence channel Cycling Orange no signal is detected. At the same time, a signal from the internal control appears in the Cycling Yellow channel.
In the sample no HSV-2 DNA is detectable. It can be considered HSV-2 negative.**

In the case of a negative HSV-2 PCR, the detected signal of the internal control rules out the possibility of PCR inhibition.

**11e. A signal is detected in the Cycling Green and Cycling Orange channels.
The result of the analysis is positive: the sample contains HSV-1 DNA as well as HSV-2 DNA.**

In this case, the detection of a signal in the Cycling Yellow channel is dispensable, since high initial concentrations of both HSV-1 and HSV-2 DNA (positive signal in the Cycling Green and Cycling Orange channels) can lead to a reduced or absent fluorescence signal of the internal control in the Cycling Yellow channel (competition).

**11f. No signal is detected in the Cycling Green, Cycling Orange or in the Cycling Yellow channels.
No result can be concluded.**

Information regarding error sources and their solution can be found in "Troubleshooting Guide", page 30.




Figure 10. Detection of the HSV-1 positive control (HSV-1 RG PC) in fluorescence channel Cycling Green. NTC: No template control (negative control).




Figure 11. Detection of the HSV-2 positive control (HSV-2 RG PC) in fluorescence channel Cycling Orange. NTC: No template control (negative control).




Figure 12. Detection of the internal control (IC) in fluorescence channel Cycling Yellow with simultaneous amplification of the positive controls (HSV-1 RG PC and HSV-2 RG PC). NTC: No template control (negative control).

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocols in this handbook or sample and assay technologies (for contact information, see back cover or visit www.qiagen.com).

Comments and suggestions

No signal with positive controls (HSV-1 RG PC and HSV-2 RG PC) in fluorescence channel Cycling Green or Cycling Orange

- | | |
|---|---|
| a) The selected fluorescence channel for PCR data analysis does not comply with the protocol |  For data analysis select the fluorescence channel Cycling Green and Cycling Orange for the analytical HSV-1/2 PCR and the fluorescence channel Cycling Yellow for the internal control PCR. |
| b) Incorrect programming of the temperature profile of the Rotor-Gene instrument |  Compare the temperature profile with the protocol. See "Protocol: PCR and Data Analysis", page 20. |
| c) Incorrect configuration of the PCR |  Check your work steps by means of the pipetting scheme, and repeat the PCR, if necessary. See "Protocol: PCR and Data Analysis", page 20. |
| d) The storage conditions for one or more kit components did not comply with the instructions given in "Storage" (page 5) |  Check the storage conditions and the expiration date (see the kit label) of the reagents and use a new kit, if necessary. |
| e) The <i>artus</i> HSV-1/2 RG PCR Kit has expired |  Check the storage conditions and the expiration date (see the kit label) of the reagents and use a new kit, if necessary. |

Comments and suggestions

Weak or no signal of the internal control of a negative CSF sample subjected to purification using the EZ1 DSP Virus Kit in fluorescence channel Cycling Yellow and simultaneous absence of a signal in channel Cycling Green or Cycling Orange

- | | |
|---|--|
| a) The PCR conditions do not comply with the protocol |  Check the PCR conditions (see above) and repeat the PCR with corrected settings, if necessary. |
| b) The PCR was inhibited |  Make sure that you use the recommended isolation method and closely follow the manufacturer's instructions. |
| c) DNA was lost during extraction |  If the internal control was added to the extraction, an absent signal of the internal control can indicate the loss of DNA during the extraction. Make sure that you use the recommended isolation method (see "DNA isolation", page 18) and closely follow the manufacturer's instructions. |
| d) The storage conditions for one or more kit components did not comply with the instructions given in "Storage" (page 5) |  Check the storage conditions and the expiration date (see the kit label) of the reagents and use a new kit, if necessary. |
| e) The <i>artus</i> HSV-1/2 RG PCR Kit has expired |  Check the storage conditions and the expiration date (see the kit label) of the reagents and use a new kit, if necessary. |

Comments and suggestions

Signals with the negative controls in fluorescence channel Cycling Green or Cycling Orange of the analytical PCR

- | | |
|---|--|
| a) Contamination occurred during preparation of the PCR | <ul style="list-style-type: none">① Repeat the PCR with new reagents in replicates.① If possible, close the PCR tubes directly after addition of the sample to be tested.① Make sure to pipet the positive controls last.① Make sure that work space and instruments are decontaminated at regular intervals. |
| b) Contamination occurred during extraction | <ul style="list-style-type: none">① Repeat the extraction and PCR of the sample to be tested using new reagents.① Make sure that work space and instruments are decontaminated at regular intervals. |

References

QIAGEN maintains a large, up-to-date online database of scientific publications utilizing QIAGEN products. Comprehensive search options allow you to find the articles you need, either by a simple keyword search or by specifying the application, research area, title, etc.

For a complete list of references, visit the QIAGEN Reference Database online at www.qiagen.com/RefDB/search.asp or contact QIAGEN Technical Services or your local distributor.

Ordering Information

Product	Contents	Cat. no.
<i>artus</i> HSV-1/2 RG PCR Kit (24)	For 24 reactions: Master, Mg Solution, 2 Positive Controls, Internal Control, Water (PCR grade)	4500263
<i>artus</i> HSV-1/2 RG PCR Kit (96)	For 96 reactions: Master, Mg Solution, 2 Positive Controls, Internal Control, Water (PCR grade)	4500265
EZ1 DSP Virus Kit — for purification of viral nucleic acids from human CSF for in vitro diagnostic purposes		
EZ1 DSP Virus Kit	For 48 viral nucleic acid preps: Prefilled Reagent Cartridges, Disposable Tip Holders, Disposable Filter-Tips, Sample Tubes, Elution Tubes, Buffers, Carrier RNA	62724
EASYartus HSV-1/2 RG PCR Kits — for fully CE-IVD-compliant integrated automated sample purification and pathogen detection		
EASYartus HSV-1/2 RG PCR Kit 1	For 48 viral nucleic acid preps and 24 assays: 1 x EZ1 DSP Virus Kit, 1 x <i>artus</i> HSV-1/2 RG PCR Kit (24)	EA10023
EASYartus HSV-1/2 RG PCR Kit 2	For 48 viral nucleic acid preps and 48 assays: 1 x EZ1 DSP Virus Kit, 2 x <i>artus</i> HSV-1/2 RG PCR Kit (24)	EA10024
Rotor-Gene Q and accessories		
Rotor-Gene Q 5plex HRM	Real-time PCR cycler and High Resolution Melt analyzer with 5 channels (green, yellow, orange, red, crimson) plus HRM channel, laptop computer, software, accessories, 1-year warranty on parts and labor	Inquire
Loading Block 72 x 0.1 ml Tubes	Aluminum block for manual reaction setup with a single-channel pipet in 72 x 0.1 ml tubes	9018901
Loading Block 96 x 0.2 ml Tubes	Aluminum block for manual reaction set-up in a standard 8 x 12 array using 96 x 0.2 ml tubes	9018905

Product	Contents	Cat. no.
Strip Tubes and Caps, 0.1 ml (250)	250 strips of 4 tubes and caps for 1000 reactions	981103
Strip Tubes and Caps, 0.1 ml (2500)	10 x 250 strips of 4 tubes and caps for 10,000 reactions	981106
PCR Tubes, 0.2 ml (1000)	1000 thin-walled tubes for 1000 reactions	981005
PCR Tubes, 0.2 ml (10000)	10 x 1000 thin-walled tubes for 1000 reactions	981008

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

This page intentionally left blank

This page intentionally left blank

This page intentionally left blank

The purchase of this product allows the purchaser to use it for the performance of diagnostic services for human in vitro diagnostics. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby.

Trademarks: QIAGEN®, *artus*®, EZ1®, EASYartus®, Rotor-Gene® (QIAGEN Group).

The *artus* HSV-1/2 RG PCR Kit and the EZ1 DSP Virus Kit are CE-marked diagnostic kits according to the European In Vitro Diagnostic Directive 98/79/EC. Not available in all countries.

Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the *artus* HSV-1/2 RG PCR Kit to the following terms:

1. The *artus* HSV-1/2 RG PCR Kit may be used solely in accordance with the *artus HSV-1/2 RG PCR Kit Handbook* and for use with components contained in the Kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *artus HSV-1/2 RG PCR Kit Handbook* and additional protocols available at www.qiagen.com.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see www.qiagen.com.

© 2009–2014 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ Orders 1-800-243-800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

Austria ■ Orders 0800-28-10-10 ■ Fax 0800-28-10-19 ■ Technical 0800-28-10-11

Belgium ■ Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

Brazil ■ Orders 55-11-5079-4000 ■ Fax 55-11-5079-4001 ■ Technical 0800-557779

China ■ Orders 800-988-0326 ■ Fax 800-988-0329 ■ Technical 800-988-0325

Denmark ■ Orders 80-885945 ■ Fax 80-885944 ■ Technical 80-885942

Finland ■ Orders 0800-914416 ■ Fax 0800-914415 ■ Technical 0800-914413

France ■ Orders 01-60-920-926 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930

Germany ■ Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

Hong Kong ■ Orders 800 933 965 ■ Fax 800 930 439 ■ Technical 800 930 425

Ireland ■ Orders 1800 555 049 ■ Fax 1800 555 048 ■ Technical 1800 555 061

Italy ■ Orders 800-789-544 ■ Fax 02-33430-4826 ■ Technical 800-787980

Japan ■ Telephone 03-6890-7290 ■ Fax 03-5547-0818 ■ Technical 03-6890-7300

Korea (South) ■ Orders 080-000-7146 ■ Fax 1544 7146 ■ Technical 1544 7145

Luxembourg ■ Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

Mexico ■ Orders 01-800-7742-639 ■ Fax 01-800-1122-330 ■ Technical 01-800-7742-436

The Netherlands ■ Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

Norway ■ Orders 800-18859 ■ Fax 800-18817 ■ Technical 800-18712

Singapore ■ Orders 1800-742-4362 ■ Fax +65-68548184 ■ Technical 1800-742-4368

Spain ■ Orders +34-91-630-7050 ■ Fax +34-91-630-5145 ■ Technical +34-91-630-7050

Sweden ■ Orders 020-790282 ■ Fax 020-790582 ■ Technical 020-798328

Switzerland ■ Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

UK ■ Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

