


Agilent OpenLab CDS

Quick Reference Sheet

Agilent
OpenLab



Part Number M8410-90061 Rev. D

OpenLab CDS – Basic Information

Terms	
ACQUISITION METHOD	Includes the instrument settings for an experiment. Edited in Acquisition .
AUDIT TRAIL	A representation of changes to a record that includes 1) identity of the user who made the change, 2) date/time of the change, 3) description of the change, and 4) reason for the change. Configure audit trails in OpenLab Control Panel .
LAYOUTS	Defines how and which information is displayed. Use preset Layouts or customize your personal layout. Changes in layouts are saved per user. Use Reset function to return to default layouts.
LINKED METHOD	Processing method that is assigned to an injection.
PROCESSING METHOD	Contains the information and parameters needed to process the data and generate results. Edited in Data Analysis .
PROJECT	Use to organize your data, e.g. by instrument, laboratory, or study. Apply individual access rights to each project
REPORT TEMPLATE	Defines the layout of a report. Edited in Data Analysis .
RESULT SET	Collection of raw data, methods, and injection list. Create custom result sets by combining single samples or sequence injections, and reprocess together.
SEQUENCE CREATION TEMPLATE	Defines a flexible pattern for creating new sequences. Useful for creating cyclic sequence. Edited in Acquisition .
SEQUENCE TABLE	Identifies the order of runs to be included in a sequence. Edited in Acquisition .

Help and Learning

GETTING STARTED



Access the interactive **Getting Started Modules** from the desktop by clicking the **OpenLab Help & Learning** icon.

The modules provide step-by-step instructions on the basic workflows and help you to familiarize yourself with OpenLab CDS.

CONTEXT SENSITIVE HELP


Use **F1** to get context sensitive help.

AGILENT COMMUNITY

<https://community.agilent.com>

Collaborate with others about applications, discuss Agilent products, and find in-depth documents and videos relevant to Agilent.

MENU COMMANDS



Find commands in the ribbon. **Data Analysis**: contextual ribbon tabs (highlighted in blue) offer commands for the selected window.

Prepare, Install, Configure

Prepare:

The OpenLab CDS Requirements guide contains information on Hardware and Software requirements, Firewall ports, OS configuration

Install, configure:

Click **setup.exe** on the USB stick to run the installation wizard.

OpenLab CDS – Quick Reference Sheet

 Access Instruments and Projects	 Acquire Data	 Analyze and Report Data
CONTROL PANEL Use the Desktop icon to launch OpenLab Control Panel. MANAGE USERS ROLE: SYSTEM ADMIN ■ Assign users to groups ■ Assign roles to groups ■ Use default roles or assign privileges to roles (use role Chemist for Acquisition and Data Analysis) ■ Edit security policy ■ Manage licenses PROJECTS ROLE: PROJECT ADMIN Create projects to organize your data: ■ Activate method audit trails ■ Apply project-specific options ■ Add custom parameters for samples and compounds Create shortcut for Data Analysis with a project. Launch Data Analysis for a project. INSTRUMENTS ROLE: INSTRUMENT ADMIN Create and configure instruments Monitor instruments Create shortcut for acquisition with an instrument. Launch Acquisition for an instrument.	CREATE ACQUISITION METHOD Define acquisition and instrument parameters for runs using this method. Download to instrument to use new parameters immediately. RUN A SINGLE SAMPLE ■ Enter sample parameters ■ Include acquisition method ■ Optional: Include processing method ■ Run RUN A SEQUENCE Create a new sequence, or Open and edit a previously saved sequence, or Import a CSV file (also possible via drag and drop), or From a result set, or Click Apply Template to apply a sequence creation template SEQUENCE Edit sequence table: Right-click and use Fill-down to edit rows, copy/paste entire rows, pin columns at the front Edit injections and sequence parameters: ■ For automated processing and reporting: enter existing processing method ■ Define name of the result set SEQUENCE CREATION TEMPLATE Create a new template in the Sequence Creation Template window: Use for cyclic sequences with varying numbers of samples. STATUS View Run Queue, Instrument status, Online signals LAUNCH DATA ANALYSIS Launch Data Analysis in review mode to show already completed injections of a running sequence: Click icon in the Run Queue.	REVIEW DATA ■ To load data, double-click node in the Data Selection view ■ Pin injections or result sets for comparison ■ Use Peak Explorer to view a sequence at a glance and identify trends or artifacts. ■ Create custom result sets based on existing injections (right-click nodes in the injection tree to add them, click Create new Result Set in the ribbon). INJECTIONS AND PROCESSING METHODS A processing method may already be linked (assigned) to the data. If not, create new master method Link injections to the master processing method A result set method (stored in result set folder) is generated on linking the master method to the injections. Reprocess (automatic after linking, or manually) EDIT PROCESSING METHOD Processing methods are divided in sections: e.g. ■ Integration Events for all or for specific signals ■ Compounds: Add compounds: Right-click a peak, or right-click the compound table. ■ Calibration: Set number of calibration levels in the General tab ■ Spectra: MS or UV reference spectra for confirmation or purity checks ■ Extraction > Spectrum: Background correction settings ■ Tools > Custom Calculation: Extend default calculations REPORT AND PRINT Import default templates into your project (start customizing templates from here): Data Selection > Import Default Templates Print results manually: In Data Analysis, edit the processing method Reports section, and link the method. Click Print all in the ribbon. Print processing method: Go to Data Processing > File ribbon tab > Print Method Print Individual PDF: Go to Reporting > Preview > Save as PDF Sequence summary report: submenu Save result set summary report