



PRODUCT INFORMATION

Long PCR Enzyme Mix

Pub. No. MAN0016323

Rev. Date 28 November 2016 (Rev. A.00)

#_

Lot _

Exp. _

Store at -20 °C

www.thermofisher.com

For Research Use Only. Not for use in diagnostic procedures.

CERTIFICATE OF ANALYSIS

Long PCR Enzyme Mix is functionally tested in PCR amplification of 47.4 kb fragment from lambda DNA.

Quality authorized by:



Jurgita Zilinskiene

CONTENTS	page
COMPONENTS	2
STORAGE.....	2
DESCRIPTION.....	2
APPLICATIONS	2
PROTOCOL	3
Reaction set-up	3
Recommended thermal cycling conditions	4
RECOMMENDATIONS FOR LONG PCR	5
Primer design	5
Template DNA	5
Magnesium concentration	5
dNTP concentration	5
DMSO.....	5
Cycling parameters	6
TROUBLESHOOTING	7

COMPONENTS

Component	#K0181 100 U	#K0182 500 U
Long PCR Enzyme Mix, 5 U/ μ L	20 μ L (100 U)	100 μ L (500 U)
10X Long PCR Buffer with 15 mM MgCl ₂	0.6 mL	2 \times 1.25 mL
10X Long PCR Buffer	0.6 mL	2 \times 1.25 mL
25 mM MgCl ₂ Solution	0.6 mL	2 \times 1.25 mL
Dimethylsulfoxide (DMSO)	0.3 mL	1 mL
Water, nuclease-free	1.25 mL	2 \times 1.25 mL

STORAGE

Store all components at -20 °C, except DMSO.
Store DMSO at room temperature.

DESCRIPTION

Long PCR Enzyme Mix is a unique blend of *Taq* DNA Polymerase and a thermostable DNA polymerase with proofreading activity. The two enzymes synergistically generate long PCR products with greater yield and three times higher fidelity than with *Taq* DNA Polymerase alone. The ratio of enzymes in the Long PCR Enzyme Mix is optimized for generation of very long amplicons: up to 47 kb with viral DNA and up to 21 kb with genomic DNA templates. The specially formulated Long PCR Buffer protects DNA from depurination and nicking during long thermal cycling. The products generated with the Long PCR Enzyme Mix are mostly 3'-dA tailed.

Long PCR Enzyme Mix is also used for efficient amplification of GC-rich DNA regions.

APPLICATIONS

- Long range PCR.
- PCR from GC-rich or difficult templates.
- Generation of PCR product for TA cloning.
- RT-PCR.

PROTOCOL

Reaction set-up

Gently vortex and briefly centrifuge all solutions after thawing. Set up the PCR reaction on ice. Setting reaction up at room temperature may result in primer degradation by 3'→5' exonuclease activity of the enzyme mix.

To prepare several parallel reactions and to minimize the possibility of pipetting errors, prepare a PCR master mix by adding water, buffer, dNTPs, and primers. Prepare enough master mix for the number of reactions and add one extra to compensate for pipetting errors. Aliquot the master mix into individual PCR tubes and add template DNA and enzyme mix.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Place a thin-walled PCR tube on ice and add the following components for each 50 μ L reaction:

Components	PCR <30 kb	PCR >30 kb and GC-rich PCR
10X Long PCR buffer with 15 mM MgCl ₂ *	5 μ L	5 μ L
dNTP Mix, 2 mM each (#R0241)	5 μ L (0.2 mM each)	5 μ L (0.2 mM each)
Forward primer	0.3-1 μ M	0.3-1 μ M
Reverse primer	0.3-1 μ M	0.3-1 μ M
Template DNA	10 pg – 1 μ g	1 ng – 1 μ g
DMSO	–	2 μ L (4%)
Long PCR Enzyme Mix	1.25-2.5 U	2.5 U
Water, nuclease-free	to 50 μ L	to 50 μ L

*If using 10X Long PCR Buffer (without MgCl₂), add following volumes of 25 mM MgCl₂ solution.

Final concentration of MgCl ₂ in 50 μ L reaction mix	1.0	1.25	1.5	1.75	2.0	2.5	3.0	4.0
Volume of 25 mM MgCl ₂ , μ L	2.0	2.5	3.0	3.5	4.0	5.0	6.0	8.0

3. Gently vortex and briefly centrifuge to collect all drops.
4. If the thermal cycler is not equipped with a heated lid, overlay the reaction mixture with a half volume of mineral oil.
5. Place the samples in a cycler and immediately start PCR.

Reaction volumes can be scaled up or down as long as the final concentrations of the reaction components remain the same.

Recommended thermal cycling conditions

Two-step cycling protocol

Step	Temperature, °C	Time	Number of cycles
Initial Denaturation	94	1-3 min	1
Denaturation	94-96	20 s	10
Annealing/ Extension	68	45-60 s/kb*	
Denaturation	94	20 s	15-25
Annealing/ Extension	68	45-60 s/kb+ x* s/cycle	
Final Extension	68	10 min	1

Three-step cycling protocol

Step	Temperature, °C	Time	Number of cycles
Initial Denaturation	94	1-3 min	1
Denaturation	94-96	20 s	10
Annealing	T _m -5	30 s	
Extension	68	45-60 s/kb*	
Denaturation	94	20 s	15-25
Annealing	T _m -5	30 s	
Extension	68	45-60 s/kb+x* s/cycle	
Final Extension	68	10 min	1

* Calculation of extension and auto-extension times:

PCR product length, kb	6	10	15	20	25	30	35	40	45
Extension time, min	4	7	10	14	17	20	24	27	30
Auto-extension per cycle, s	2	5	5	10	10	15	15	20	20

RECOMMENDATIONS FOR LONG PCR

Primer design

- PCR primers should be 27-36 nucleotides long.
- Optimal GC content of the primer is 40-60%. Ideally, C and G nucleotides should be distributed uniformly along the primer.
- Prefer one or two G or C at the 3'-end of the primer, but avoid placing more than three G or C nucleotides at the 3'-end to lower the risk of nonspecific priming.
- Avoid primer self-complementarities, complementarities between the primers and direct repeats in a primer to prevent hairpin formation and primer dimerization.
- Check for possible complementary sites between primers and template DNA.
- When designing degenerate primers, place at least 3 conservative nucleotides at the 3'-end.
- Differences in melting temperatures (T_m) of the two primers should not exceed 5 °C

Template DNA

Use 0.01-1 ng of plasmid and phage DNA, and 0.1-1 µg of genomic DNA in 50 µL of PCR. Higher amounts of template increase the risk of nonspecific PCR products. Lower amounts of template may reduce the accuracy of the amplification.

All routine DNA purification methods are suitable for template preparation e.g., Thermo Scientific GeneJET™ Genomic DNA Purification Kit (#K0721), GeneJET Plasmid Miniprep Kit (#K0502). Trace amounts of certain agents used for DNA purification, such as phenol, EDTA and proteinase K, may inhibit DNA polymerases. Ethanol precipitation and repeated washes of the DNA pellet with 70% ethanol normally remove trace contaminants from DNA samples. High quality and integrity of the template DNA is essential for reliable amplification of large DNA fragments. Extreme care must be taken during preparation and handling of template for long PCR. Nicked or damaged DNA can serve as a potential priming site resulting in high background. Avoid repeated freezing/thawing of template DNA.

Magnesium concentration

10X Long PCR Buffer with 15 mM $MgCl_2$ is optimal in most cases. If DNA samples contain EDTA or other chelators, increase $MgCl_2$ concentration in the reaction mixture accordingly. The 10X Long PCR buffer without $MgCl_2$ and 25 mM $MgCl_2$ solution are included for optimization, if necessary.

dNTP concentration

Final concentration of 0.2 mM for each dNTP is optimal in most cases.

Enzyme concentration

For PCR up to 20 kb, use 1-1.25 U of Long PCR Enzyme Mix per 50 µL reaction volume. For PCR of ≥ 20 kb, use up to 2.5 U per 50 µL.

DMSO

DMSO, a cosolvent, increases yields and improves reliability of the system for long PCR and PCR of complex targets. DMSO is used at a final concentration between 1 and 12%. For PCR fragments ≥ 30 kb, use 4% DMSO.

Cycling parameters

Initial denaturation step

It is essential to completely denature the template DNA at the beginning of PCR to ensure efficient utilization of the template during the first amplification cycle. If GC content of the template is 50% or less, an initial 1-3 min denaturation at 94 °C is sufficient. For GC-rich templates, this step may be prolonged to 5 min. Use shortest possible denaturation time. Exposure of DNA to high temperatures may cause depurination of single-stranded DNA, which leads to DNA truncation.

Denaturation step

Optimal denaturation time and temperature depend on PCR tubes and thermal cycler used. Thin-wall PCR tubes is the best choice. The optimal denaturation conditions for Applied Biosystems GeneAmp™ 9700 Thermocycler are 96 °C for 10 s, for Eppendorf Mastercycler™ – 95 °C for 15 s, for Perkin Elmer 480 – 94 °C for 20 s. When using other thermocyclers cycling conditions have to be determined experimentally.

Insufficient denaturation time or temperature may cause either diffuse smearing after electrophoresis or poor amplification efficiency. Excess denaturation time or temperature may result in nonspecific PCR products.

Annealing/Extension step

In most cases a two-step cycling protocol (denaturation followed by annealing/extension) is preferred over three-step cycling protocol (denaturation followed by annealing followed by extension). The two step protocol saves time and reduces exposure of enzymes and template to long cycling conditions. The three-step cycling protocol is used when annealing temperature of the primers is lower than 65 °C.

For annealing/extension at 68 °C 1 min/kb is recommended for amplicons up to 6 kb. For longer amplicons annealing/extension time, starting from 11th cycle, has to be prolonged by approximately 5 seconds to compensate for enzyme loss during cycling. See recommendations in cycling protocol p.4.

Too low extension temperature reduces specificity of amplification. Insufficient extension time results in no amplification product or short nonspecific products. Excessive extension time causes smeared bands after electrophoresis.

Cycle number

Optimal number of cycles is in a range of 25-35 and depends on the amount and complexity of template DNA and also on amplicon length. Too few cycles may not generate sufficient amounts of PCR product, while too many cycles may produce a diffuse smear in electrophoresis gel. As a guideline, use 30-35 cycles for single copy targets in genomic DNA and 25-30 cycles for plasmid or phage DNA.

Final extension step

After the last cycle, samples are usually incubated at 68-72 °C to fill-in the protruding ends of the PCR product.

TROUBLESHOOTING

Problem	Possible cause and solution
<p>Low or no PCR product</p>	<p>PCR component may be missing or degraded. Use a checklist when assembling reactions. Always perform a positive control to ensure that each component is functional.</p> <p>Poor template quality. Evaluate template integrity by electrophoresis on agarose gel. If necessary, repurify template using methods that minimize shearing and nicking. Isolate fresh template and resuspend it in TE buffer, pH 8.0 or sterile water. Also, try to increase enzyme amount in the reaction.</p> <p>Difficult template. If DNA is GC-rich or contains secondary structure, 1-12% dimethyl sulfoxide (DMSO) in the reaction mixture may help. Also, use higher enzyme concentration up to 2.5 U per 50 μL.</p> <p>Denaturation temperature not optimal. Optimize denaturation temperature by 1 $^{\circ}$C increments. Too high denaturation temperature can lead to degradation of the template, especially for long target sequences.</p> <p>Denaturation time not optimal. Optimize denaturation time by decreasing or increasing in 5 s increments. Too long denaturation time can lead to degradation of the template, especially for long target sequences.</p> <p>Insufficient number of cycles. Increase the number of cycles in 3-5 increments.</p> <p>Insufficient amount of template. Use more template or increase number of cycles.</p> <p>Suboptimal primer design. Check sequence information. Re-design the primers according to recommendations on p.5.</p> <p>Annealing temperature too high. Decrease the annealing temperature in 2-4 $^{\circ}$C increments.</p> <p>Insufficient extension time. Optimize by increasing the extension time in 1-minute increments. Make sure to include auto-extension per cycle starting from 11th cycle.</p> <p>Insufficient magnesium concentration. If DNA template contains EDTA or other chelators or RNA impurities, the MgCl₂ concentration should be increased in 0.2-0.3 mM increments.</p>

(continued)

Problem	Possible cause and solution
Band smearing	<p>Excess template. Reduce amount of template in reaction.</p> <p>Too many cycles. Reduce cycles in increments of 3-5.</p> <p>Extension time too long. Decrease elongation time in 1-2 min increments.</p> <p>Denaturation time not optimal. Optimize denaturation time by decreasing or increasing in 5 seconds increments. Too long denaturation time can lead to degradation of the template, especially for long target sequences.</p> <p>Denaturation temperature too low. Try increasing the denaturation temperature in 1 °C increments.</p> <p>Excess magnesium concentration. If dNTP concentration is lower than 0.2 mM of each dNTP, MgCl₂ concentration should be decreased accordingly. Make sure template DNA solution does not contain magnesium.</p> <p>Excess amount of enzyme mix. Be sure to add only 1.25-2.5 units of enzyme mix per 50 µL reaction and mix the reaction well after addition of the enzyme mix.</p> <p>Insufficient extension time Especially with longer templates, increase the elongation time in 1-2 min increments.</p> <p>Poor template quality Check template integrity by electrophoresis on agarose gel. If necessary, repurify the template using methods that minimize shearing and nicking. Isolate fresh template and resuspend it in TE buffer, pH 8 or sterile water.</p>
Non-specific PCR products	<p>Too much template When amplifying genomic DNA, the initial concentration of the template in the PCR reaction should not exceed 1 µg per 50 µL reaction volume.</p> <p>Suboptimal primer design Check sequence information. Re-design the primers according to recommendations on p.5.</p> <p>Too many cycles Reduce the number of cycles to eliminate nonspecific products.</p> <p>Annealing temperature too low Increase temperature in 2-3 °C increments.</p> <p>Excess magnesium concentration If dNTP concentration is lower than 0.2 mM, MgCl₂ concentration should be decreased accordingly. Make sure template DNA solution does not contain magnesium.</p>

NOTICE TO PURCHASER

Use of this product is covered by US Patent No. 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5781 Van Allen Way, Carlsbad, California 92008.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5781 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and *in vitro* use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Please refer to www.thermofisher.com for Material Safety Data Sheet of the product.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. GeneAmp is a trademark of Roche Diagnostics GmbH. MasterCycler is a trademark of Eppendorf AG. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.