Presented With Compliments From
www.DirectHomeMedical.com

simplyFLO

PHILIPS

RESPIRONICS

Table of Contents

Introduction	1
Intended Use.....	1
Contraindications.....	1
Device and Accessories Description.....	1
simplyFLO System	1
Warnings and Cautions	2
Symbol Key.....	5
How to Contact Philips Respironics.....	5
System Overview.....	6
Operating the simplyFLO Device.....	8
Alarm Indicators and Screen Symbols.....	11
Troubleshooting.....	14
Connecting the Optional Humidifier	15
Caring for your simplyFLO System	16
Cleaning the Device	16
Cleaning the Accessories.....	16
Device Storage	17
Preventive Maintenance.....	17
Device Disposal.....	17
Traveling With Your System	17
Specifications.....	18
Standards Compliance	19
Classification	19
EMC Information.....	20
Limited Warranty.....	22

Introduction

Your health care professional has determined that supplemental oxygen is of benefit to you and has prescribed a portable oxygen concentrator set at a specific flow setting to meet your needs. DO NOT change the flow settings without consulting with your health care professional. Please read and understand this entire manual before using the device.

Intended Use

The Philips Respironics simplyFLO transportable oxygen concentrator is for prescription use by patients requiring high concentrations of oxygen on a supplemental basis.

Contraindications

- This device is not intended to be life-sustaining or life-supporting.
- This device is not intended for newborn and infant use.

Device and Accessories Description

The simplyFLO transportable oxygen concentrator, its features, and its accessories are described in detail in this manual. Read and understand it completely before operating the device.

This manual applies to the AC Power Supply and AC Power Cord accessories.


This manual does not apply to the following accessory items sold separately:

- Standard Single Lumen Nasal Cannula
- Humidifier
- DC Power Supply
- Travel Case

simplyFLO System

The device comes with the following items. If any are missing, contact your equipment provider.

- This User Manual
- AC Power Supply and Cord

Warnings and Cautions

Caution: U.S. federal law restricts this device to sale by or on the order of a physician.

Warnings

A warning represents the possibility of harm to the operator or patient.

- The operator should read and understand this entire manual before using the device.
- The device is not intended for life support. Where the prescribing health care professional has determined that an interruption in the supply of oxygen, for any reason, may have serious consequences to the user, an alternate source of oxygen should be available for immediate use.
- Geriatric or any other patient unable to communicate discomfort, or hear or see the alarms while using this device, may require additional monitoring.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- Do not use the device in the presence of a flammable anaesthetic mixture in combination with oxygen or air, or in the presence of nitrous oxide.
- Do not use oil or grease on the concentrator or its components as these substances, when combined with oxygen, can greatly increase the potential for a fire hazard and personal injury.
- If you notice any of the following, discontinue use and contact your home care provider:
 - unexplained changes in the performance of this device
 - unusual or harsh sounds
 - dropped or mishandled device or the power supply
 - water spilled into the enclosure
 - broken enclosure
- Use only with Philips Respironics simplyFLO AC power supply.
- Use only approved simplyFLO accessories.

Warnings (continued)

- Repairs and adjustments must be performed by Philips Respironics authorized service personnel only. Unauthorized service could cause injury, invalidate the warranty, or result in costly damage.
- Periodically inspect electrical cords, cables, and the power supply for damage or signs of wear. Discontinue use and replace if damaged.
- To avoid electric shock, unplug the device before cleaning. DO NOT immerse the device in any fluids.
- Your home care provider is responsible for performing appropriate preventive maintenance at the intervals recommended by Philips Respironics.
- For proper operation, your device requires unobstructed ventilation. Always make sure any openings in the case are not obstructed by items which may impede ventilation. Do not place the device in a small closed space (such as a closet). The device should not be used adjacent to or stacked with other equipment. For more information, contact your home care provider.
- Do not use an extension cord.
- Device operation above or outside of the voltage, temperature, humidity and/or altitude values specified may decrease oxygen concentration levels.
- Never drop or insert any object into any opening.
- Be aware that the electrical cord and/or tubing could present a tripping or strangulation hazard.
- Use only power cords supplied by Philips Respironics for this device. Use of power cords not supplied by Philips Respironics may cause overheating or damage to the device and may result in increased emissions or decreased immunity of the equipment or system.
- The user should be aware that there is no backup power.
- Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual.
- Place the device in a location so as to avoid pollutants and fumes.
- Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment. See the EMC section of this manual for distances to observe between RF Generators and the simplyFLO device to avoid interference.

Cautions

A caution represents the possibility of damage to the equipment.

- Do not immerse the device or allow any liquid to enter the enclosure.

Note: Additional warnings, cautions, and notes are located throughout the manual.

Symbol Key

Consult accompanying instructions for use

Model number

No smoking

Serial number

No oil or grease

Type BF Applied Part

Keep away from open flames

IPX1

Drip proof equipment

Do not disassemble

AC Power

DC power

Class II (Double Insulated)

Caution: U.S. federal law restricts this device to sale by or on the order of a physician.

Note: Also see the section on alarm indicators and screen symbols.


How to Contact Philips Respironics

To have your device serviced, contact your home care provider. If you need to contact Philips Respironics directly, call the Customer Service department at 1-724-387-4000 or 1-800-345-6443 (US or Canada). You can also use the following address:

Respironics Inc.
1001 Murry Ridge Lane
Murrysville, PA 15668 USA


System Overview

Begin by becoming familiar with your simplyFLO device.


Item	Description	Function
1	Control Panel	Control switches and LCD display
2	Carrying Handle	Hand grip area for lifting and transporting the device
3	Patient Cannula Connector	Oxygen output and connection point for patient cannula
4	Power Input Connector	Connection point for Philips Respironics-provided external power supplies: AC line voltage
5	Air Exhaust Vent	Air outlet for enclosure ventilation fan
6	Air Intake Vent	Air inlet for enclosure ventilation fan
7	Data Port	Input-Output data communication with the device. Caution: For use by authorized service personnel only.

Control Panel

Item	Description	Function
1	Power	Turns the device On and Off
2	Alarm Silence	Turns audible alarms off
3	Display Screen	Shows information about the operating status of the device
4	Plus (+)/Minus (-)	+ Increases the displayed setting; - Decreases the displayed setting
5	Patient Cannula Connector	Oxygen output and connection point for patient cannula

Run Time Screen Description


Item	Symbol Description	Function
1	Operating Mode	For detailed descriptions, refer to the “Alarm Indicator and Screen Symbols” section.
2	Setting	
3	Audible Alarm Status	
4	Alarm Type	

Operating the simplyFLO Device

1. Select a location that allows the concentrator to draw in room air without being restricted. Make sure that the device is at least 15 to 30 cm away from walls, furniture, and especially curtains that could impede adequate airflow to the device. Do not place the device near any heat source.
2. After reading this entire manual, connect the output cable of the AC Power Supply to the simplyFLO power input connector.
3. Connect the AC Power Cord's connector into the AC power supply.
4. Connect the AC Power Cord's plug into an AC outlet installed to meet national or regional electrical code.


Note: The AC Power Supply is equipped with a green light, which illuminates when it is connected to an energized AC power source.

5. Connect a nasal cannula to the patient cannula connector on top of the device. Ensure that the cannula is routed to prevent it from being pinched or kinked to avoid a disruption of oxygen flow.


Note: A standard single lumen nasal cannula and tubing (not supplied) should be used to deliver oxygen from the device. The device works with cannula tubing up to 100 ft. (30.5 m).

6. To start the device, press the Power button. After the first press, the display screen and the 4 button backlights turn on and you will see a screen similar to the one shown below.

ITEM	DESCRIPTION
1	Push Power again to turn the device on.
2	Total number of operating hours for the device.
3	Software version in the device.
4	Indicates the allowable modes configuration of the device. Note: This may not display on all devices. This is normal.

- Press the Power button a second time to place the device in operating mode. This second press helps prevent unwanted starts or inadvertent changes in oxygen flow. The Introductory screen momentarily appears, shown at left below, followed by the final operating mode screen, shown at right.

Every time it starts up, the device starts producing oxygen to meet the specified setting. It may take up to twenty minutes to reach the specified purity levels for the setting. It is safe to begin breathing from the simplyFLO device immediately.

Note: When it is turned on, the device will automatically start at the flow rate setting used when the device was last turned off. As a precaution, each time you start the device, verify the flow setting.

- To set the device's flow setting, press either the + or – button.

To increase the flow setting, press the + button. To decrease the flow setting, press the – button.

The device can be set to .5, 1, 1.5, and 2 liters per minute (LPM) of continuous flow.

Warning: It is very important to set your device to your prescribed level of oxygen flow. Do not increase or decrease your flow rate from your prescribed level until you first consult with your health care professional.

- Be sure oxygen is flowing through the cannula. If it is not, refer to the Troubleshooting section of this manual.

simplyFLO User Manual

10. Position the nasal cannula on your face, as shown below, and breathe normally through your nose.

11. When you are ready to turn your simplyFLO device off, press the Power button. The screen below appears.

12. Press the Power button a second time and the simplyFLO device will turn off.
13. Rather than pressing the Power button a second time, when this screen is displayed, the backlight intensity may be increased or decreased by repeatedly pressing the Plus (+) or Minus (-) key until the desired backlight intensity is achieved. The device may then be powered off using steps 11 and 12 above.

Alarm Indicators and Screen Symbols

Visual, Audio Indicators	Description	What To Do
 A sequence of 3 long beeps + 1 short beep that repeats every 4 seconds.	Low Oxygen Concentration Alarm This alarm occurs when the device is delivering a lower concentration of oxygen than specified. This alarm can be silenced by pressing the Silence button. Note: This alarm occurs when the internal O₂ sensor reads <82% O₂. This alarm will turn off if the system is able to recover from the fault and achieve >84% O₂.	Change to another source of oxygen and contact your home care provider.
 A sequence of 3 short beeps + 3 long beeps that repeats every 4 seconds.	Technical Fault Alarm This alarm occurs when the device has a general malfunction and the device is no longer operating properly. The device may or may not shut down depending on the severity of the alarm. This alarm can be silenced if the  symbol appears on-screen. If not, the user must turn off the power to silence the alarm.	If this alarm occurs again, change to another source of oxygen and contact your home care provider.

Visual, Audio Indicators	Description	What To Do
 A sequence of 3 long beeps + 1 short beep that repeats every 4 seconds.	No Flow Alarm This alarm occurs when the device detects there is no oxygen flowing in the patient cannula. This alarm can be silenced by pressing the Silence button.	Check the cannula for kinks or other obstructions that are stopping oxygen flow through it.
	Alarm Symbol This symbol appears when an event occurs that generates an audible alarm.	Refer to the associated alarm indicator also being displayed.
	Alarm Silenced Symbol This symbol appears when you have pressed the alarm silence button on the control panel to silence an audible alarm.	Refer to the associated alarm indicator also being displayed.
	Attention This symbol appears when corrective action is required.	Ensure that the intake and exhaust vents are not blocked. Allow the device to warm up. If the alarm continues, contact your home care provider.

Visual, Audio Indicators	Description	What To Do
 1st press of the Power button = 1 beep. 2nd press of the Power button to turn device OFF = 3 beeps	On/Off Symbol This symbol appears when the device is powered on, or is ready to be turned on or off.	Press the Power button to turn the device on or off. Note: To prevent inadvertent starts and stops, the device requires two presses of the Power button to turn it on and off.
2.0 1 short beep	Flow Control Setting The large number on the right side of the screen shows the oxygen flow setting. (In this example, it is 2.) This number can range from .5 to 2 for Continuous Flow.	If this is your flow setting, no action is needed. If this is not your flow setting, press the + button to increase this number or press the – button to decrease this number.
Constant audible alarm	Loss of Power alarm The simplyFLO will sound a constant audible alarm if the device detects it has lost sufficient power to operate properly.	Verify that the device is plugged into an appropriate external power source. On the simplyFLO AC power supply, verify that the green LED is illuminated. If the device has in fact lost power, the alarm can be silenced by tilting the device 90 degrees to the left or right. If power cannot be restored to the device, change to another source of oxygen and contact your home care provider.

Troubleshooting

The table below lists common problems and actions you can take. If you are unable to resolve a problem, please contact your home care provider.

Problem	Possible Cause	What You Should Do
Device Won't Turn On	Power button not pressed twice.	Check that the power cord is plugged in. Press the power button twice.
Device Will Not Deliver Oxygen	simplyFLO is not turned on.	Press the power button twice.
	Cannula tubing is kinked or twisted.	Make sure the tubing is connected properly to the oxygen outlet port and that it is free of any obstruction.
	Device malfunction.	Contact your home care provider.
Oxygen Not At Full Concentration	Device is warming up.	Wait 10 minutes for the device to deliver oxygen at the prescribed concentration. If the condition persists, contact your home care provider.
Alarm Occurs	The device needs your attention.	See the Alarm Indicators and Screen Symbols section for information about specific alarms and what you should do.

Connecting the Optional Humidifier

Note: Application accessories shall include a means to reduce the propagation of fire.

Note: A bubble humidifier may be used with the simplyFLO.

Warnings

- Only use the humidifier connector tubing supplied with the SimplyGo/simplyFLO Humidifier Connecting Tube kit (REF 1101172).
- Do not overfill humidifier.
- Do not reverse the oxygen input and output connections. Water from the humidifier bottle will travel through the cannula back to the patient.

Note: Only connect the humidifier if prescribed.

To connect the humidifier:

1. Remove the cap from the humidifier bottle.
2. Fill the humidifier with distilled water to the level indicated by the manufacturer.
3. Replace the cap and make sure it is tightened securely.
4. Attach the humidifier bottle adapter tube to the humidifier bottle by turning the wing nut on the bottle counterclockwise until it is fastened securely.
5. Place the humidifier bottle in a secure location on a level surface, or use the simplyFLO Humidifier Bottle holder (REF 1076829).
6. Attach the oxygen tubing from the humidifier bottle assembly to the oxygen outlet connector on the simplyFLO.
7. Attach the cannula tubing to the humidifier bottle outlet.
8. After assembly is complete, make sure the oxygen is flowing through the cannula.

Caring for your simplyFLO System

Cleaning the Device

Warning: To avoid electrical shock, do not remove the simplyFLO covers. The covers should only be removed by authorized service personnel. Do not apply liquid directly to the covers. Do not use alcohol, solvents, polishes or any oily substances on the device, as they are flammable.

Caution: Do not allow liquids into any of the controls, the interior of the case, or the oxygen tubing connector. If this occurs, contact your home care provider for assistance.

The exterior cover of the device should be cleaned as needed:

1. Turn the device off and disconnect from the power source before cleaning.
2. Clean the device exterior using a damp cloth with a mild household cleaner and wipe it dry. Or, if you use medical disinfectants, be sure to follow the manufacturer's instructions. For disinfection, Philips Respironics recommends the following:
 - Discide Ultra Towelettes: Manufacture - Palmero, www.palmero-health.com (or equivalent)

Cleaning the Accessories

Follow the instructions below when cleaning the accessories.

AC Power Supply/Power Cord

The AC power supply and power cord should be cleaned as needed.

1. Disconnect cords from device before cleaning.
2. Clean the cords using a cloth with a mild household cleaner and wipe dry.

Cannula and Tubing

Clean and replace the cannula and tubing as instructed by the manufacturer and your home care provider.

Humidifier

Clean and replace both the humidifier bottle and humidifier tubing according to the manufacturer's and your home care provider's recommendations.

Device Storage

Store your device in a place where it will remain clean and dry.

Caution: Do not store your device or accessories in extreme temperatures, below -4° F (-20° C) or above 140° F (60° C).

Preventive Maintenance

This device does not require routine servicing. Please refer to the simplyFLO Service Manual for detailed maintenance instructions.

Device Disposal

Dispose of the device in accordance with local regulations or contact your home care provider.

Traveling With Your System

With your simplyFLO system and proper advance planning, you can enjoy traveling within your community and beyond. Before you leave, make sure you pack the following:

- AC power supply and connector cord
- Travel case (optional)
- Extra cannulas

Also, be sure to take the telephone numbers of your home care provider and physician in case of an emergency.

Specifications

Operating Conditions	Operational temperature: 41° F to 104° F (5° C to 40° C) Relative humidity: 15% to 95% Altitude: up to 10,000 ft. (3048 m.)
Storage and Transport Conditions	-4° F to 140° F (-20° C to 60° C) Relative humidity: up to 95%, non condensing

simplyFLO Oxygen Concentrator

Oxygen Concentration*	87%-96% at all settings
Continuous Flow Settings	Continuous Flow .5 = 500 ml/min; 1 = 1000 ml/min; 1.5 = 1500 ml/min; 2 = 2000 ml/min; +/- 15% or 150 ml/min, whichever is greater (3 minute running average) Note: Max recommended flow is 2 LPM (at nominal outlet pressure of 0 and 7 hPa)
User Interface	Push buttons, back-lit liquid crystal display (LCD)
Dimensions	11.5 in x 10 in x 6 in (29.2 cm x 25.4 cm x 15.2 cm)
Weight	<9.0 lbs (4.1 kg)
Sound Level	50 dBA typical at 2 LPM continuous
Audio Alarm	50 decibels (nominal) at 39 in. (1 m)
Outlet Pressure	6.4 psig maximum
*Based on an atmospheric pressure of 14.7 psia (101 kPa) at 70° F (21° C)	

AC Power Module

Type	Philips Respironics 1082661
Input	100 to 240 VAC, 50/60 Hz, 1.7 A max
Output	19 VDC, 7.9 A max

Standards Compliance

The device is designed to conform to the following standards

- IEC 60601-1, Medical Electrical Equipment, Part 1: General Requirement for Safety
- IEC 60601-1-2, 2nd edition, Medical Electrical Equipment, Part 1-2: General Requirement for Safety – Collateral Standard: Electromagnetic (EMC) Compatibility – Requirements and tests.
- ISO 8359, Oxygen Concentrators for Medical use – Safety Requirements.

Classification

The simplyFLO is classified as:

- IEC Class II Internally Powered Equipment
- Type BF Applied Part
- IPX1: Drip Proof Equipment
- Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Continuous Operation

EMC Information

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
simplyFLO is intended for use in the electromagnetic environment specified below. The user of simplyFLO should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	The Device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
simplyFLO is intended for use in the electromagnetic environment specified below. The user of simplyFLO should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV Contact ±8 kV Air	±6 kV Contact ±8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for Power Supply Lines ±1 kV for Input/Output Lines	±2 kV for Power Supply Lines ±1 kV for Input/Output Lines	Mains power quality should be that of a typical home or hospital environment.
Surge IEC 61000-4-5	±1 kV Line to Line ±2 kV Line to Ground	±1 kV Line to Line ±2 kV Line to Ground	Mains power quality should be that of a typical home or hospital environment.
Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	<5% U_T (>95% Dip in U_T) for 0.5 Cycle 40% U_T (60% Dip in U_T) for 5 Cycles 70% U_T (30% Dip in U_T) for 25 Cycles <5% U_T (>95% Dip in U_T) for 5 Seconds	<5% U_T (>95% Dip in U_T) for 0.5 Cycle 40% U_T (60% Dip in U_T) for 5 Cycles 70% U_T (30% Dip in U_T) for 25 Cycles <5% U_T (>95% Dip in U_T) for 5 Seconds	Mains power quality should be that of a typical home or hospital environment. If the user of the Device required continued operation during power mains interruptions, it is recommended that the Device be powered from an uninterruptible power supply or battery.
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical home or hospital environment.
Note: U_T is the A.C. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
simplyFLO is intended for use in the electromagnetic environment specified below. The user of simplyFLO should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the Device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ 150 kHz to 80 MHz $d = 0.35\sqrt{P}$ 80 MHz to 800 MHz $d = 0.7\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m ^c	
Note 1:	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.		
a	Field strength from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Device is used exceeds the applicable RF compliance level above, the Device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Device.		
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		
c	simplyFLO incorporates a robust design and conforms with Radiated RF Immunity Test up to 10 V/m between 80 MHz and 2.5 GHz while the test requirement is 3 V/m.		

Recommended Separation Distance Between Portable and Mobile RF Communication Equipment and the Device			
simplyFLO is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The user of simplyFLO can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and simplyFLO as recommended below, according to the maximum output power of the communications equipment.			
Rated Maximum Output Power of Transmitter (Watts)	Separation Distance According to Frequency of Transmitter (meters)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 0.35\sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7\sqrt{P}$
0.01	0.12	0.04	0.07
0.1	0.37	0.11	0.22
1	1.17	0.35	0.70
10	3.69	1.11	2.21
100	11.67	3.50	7.00
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1:	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.		

Limited Warranty

Respironics, Inc. ("Respironics") warrants the simplyFLO transportable oxygen concentrator system (the "Product") as set forth in the following paragraphs.

Limited Warranty

Respironics warrants that the Product, except as otherwise stated herein, is free from defects in materials and workmanship under normal and proper use and when correctly maintained in accordance with applicable instructions, for a period equal to the warranty period purchased with the device, or if not otherwise specified, for a period of two (2) years from the date of shipment.

Products Not Under Warranty

Accessories and replacement parts, including, but not limited to, filters and fuses, are not covered under this warranty.

Limitations

If any Product purchased from Respironics fails to conform to the warranties set forth herein during the Warranty Period, as determined by Respironics in its sole discretion, Respironics may discharge its warranty obligation by repairing or replacing the Product, in Respironics' sole discretion. This may be accomplished by installing new or remanufactured assemblies or components, or by other repairs deemed appropriate in the sole discretion of Respironics. The choice of repair or replacement by Respironics shall be the sole and exclusive remedy of the original purchaser. Respironics reserves the right, in its sole discretion, to refund the purchase price in lieu of repair or replacement of the Product. In no event shall Respironics' maximum liability under these warranties exceed the price paid to Respironics by the original purchaser for the Product.

Conditions

This warranty does not cover damage or injury whether to the Product or to personal property or persons caused by accident, misuse, abuse, negligence, failure to install in accordance with Respironics' installation instructions, failure to operate under conditions of normal use and in accordance with the terms of the operating manual and instructions, failure to maintain in accordance with the applicable service manuals, or alteration or any defects not related to materials or workmanship of the Product. This warranty does not cover damage which may occur in shipment. This warranty does not apply to any Product or individual part of a Product that may have been repaired or altered by anyone other than Respironics or an authorized Respironics service center. This warranty does not apply to any Product which is not purchased new.

Exclusions of Warranties

EXCEPT AS SET FORTH IN THIS LIMITED WARRANTY, RESPIRONICS MAKES NO WARRANTIES, EXPRESSED OR IMPLIED, STATUTORY OR OTHERWISE, REGARDING THE PRODUCT, ITS QUALITY OR PERFORMANCE. RESPIRONICS SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTY OF MERCHANTABILITY AND THE IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL RESPIRONICS BE LIABLE FOR THE COST OF PROCUREMENT OF SUBSTITUTE PRODUCTS OR FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL, PUNITIVE, EXEMPLARY OR INCIDENTAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, COMMERCIAL LOSS OR LOST REVENUES), FOR ANY CAUSE OF ACTION, WHETHER IN CONTRACT OR TORT, AND WHETHER OR NOT RESPIRONICS WAS AWARE OR SHOULD HAVE BEEN AWARE OF THE POSSIBILITY OF THESE DAMAGES. EXCEPT FOR THE OBLIGATIONS UNDER THIS LIMITED WARRANTY, RESPIRONICS SHALL NOT HAVE ANY OBLIGATION OR LIABILITY FOR ANY OTHER LOSS, DAMAGE OR INJURY RESULTING DIRECTLY OR INDIRECTLY FROM THE PRODUCT. PURCHASER'S SOLE AND EXCLUSIVE REMEDY FOR BREACH OF THE WARRANTIES SET FORTH IN THIS DOCUMENT SHALL BE AS PROVIDED IN THE PARAGRAPH DESCRIBING LIMITATIONS.

Purchaser is cautioned that no person or entity is authorized to make any warranties on behalf of Respironics and any such alleged warranties are hereby disclaimed by Respironics.

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