



Orbitrap Tribrid Series

Preinstallation Requirements Guide

80000-97024 Revision A • June 2018



© 2018 Thermo Fisher Scientific Inc. All rights reserved.

EASY-nLC, Fusion, ID-X, Lumos, and Orbitrap ID-X are trademarks; Unity is a registered service mark; and Dionex, Orbitrap Fusion, Orbitrap Fusion Lumos, Thermo Scientific, Tribrid, UltiMate, and Xcalibur are registered trademarks of Thermo Fisher Scientific Inc. in the United States. Optima is a registered trademark of Fisher Scientific Co.

The following are registered trademarks in the United States and other countries: Microsoft and Windows are registered trademarks of Microsoft Corporation. Teflon is a registered trademark of E.I. du Pont de Nemours & Co.

The following are registered trademarks in the United States and possibly other countries: Dranetz is a registered trademark of Dranetz Technologies, Inc. Oerlikon Leybold Vacuum is a registered trademark of OC Oerlikon Corporation AG. Parafilm is a registered trademark of Bemis Company, Inc. Powervar is a registered trademark of Powervar, Inc. SOGEVAC is a registered trademark of Oerlikon Leybold Vacuum. ShockWatch is a registered trademark of ShockWatch, Inc.

Triton is a trademark of Sigma-Aldrich Co.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history: Rev A, June 2018

For Research Use Only. Not for use in diagnostic procedures.

Orbitrap Tribrid Series Installation Request Form

Before completing this installation request form, read all of the *Orbitrap Tribrid Series Preinstallation Requirements Guide*. For laboratories outside the U.S., go to www.thermofisher.com, click **Contact Us**, select the **Instrument Support** option, type the product name, and request to be contacted by email to schedule your instrument installation. You can then send this completed and signed form to the office handling the installation. For laboratories in the U.S., send this completed and signed form to us.customer-support.analyze@thermofisher.com.

Laboratory setup

- 1. All laboratory remodeling is complete and complies with all relevant safety regulations.
- 2. The instrument is on site.
- 3. A principal operator will be on site during the installation/certification period.
- 4. Doorways, hallways, and so on are a minimum width of 1.4 m (54 in.).
- 5. Laboratory lighting is adequate.
- 6. Air conditioning is adequate for temperature, humidity, and particulate matter control.
- 7. Relative humidity is 40–80%, noncondensing.
- 8. The work area is free from magnetic disruption and electrostatic discharge.
- 9. A step stool is on site.
- 10. (Optional) The laboratory has Internet access.
- 11. Floor space is sufficient and flooring will support the load.

Power

- 12. Main power is installed and complies with local electrical codes.
- 13. Power is free from fluctuations due to slow changes in the average voltage or changes due to surges, sags, or transients.
- 14. Power outlets are of the correct configuration for the power cords. See [page 2](#).
Note NEMA type: _____
- 15. Voltage of power outlet has been measured.
Note measured voltage:
AC line-to-ground: _____ V
AC neutral-to-ground: _____ V
AC line-to-neutral: _____ V

- 16. Power outlets are available for testing and cleaning equipment.

Gas and exhaust

- 17. All required gases are on site, gas lines are installed, and appropriate gas regulators are available. For pressures, see [page 4](#).
Note gas types and actual purity levels:
Gas: Helium purity: _____
Gas: Nitrogen #1 purity: _____
Gas: Nitrogen #2 purity: _____
- 18. All gas lines are clean and have no leaks.
- 19. A suitable fume exhaust system is separate from the solvent waste and is within 2.4 m (8 ft) of the system. See [page 5](#) and [page 11](#).

System setup

- 20. Data system computer: The new computer shows no changes to ANY settings and has no additional software.
- 21. System setup provides for collecting solvent waste from the API source.
- 22. A new or recently cleaned HPLC system is available that produces pulse-free, continuous flow of 100–1000 $\mu\text{L}/\text{min}$.
- 23. Optima LC/MS-grade acetonitrile, isopropyl alcohol, methanol, and water are available for testing the instrument's performance.
- 24. Sufficient clearance is provided behind the instrument's table. See [page 10](#).

IMPORTANT Thermo Fisher Scientific reserves the right to invoice for the field service engineer's time if the installation requirements are not met by the installation date.

CONTINUED ON THE NEXT PAGE

For customized installations

Does your contract contain any special acceptance specifications? Yes No
If YES, attach full details of the specifications.

Does the system require additional equipment? Yes No
If YES, attach full details of the additional equipment.

I certify that the preinstallation requirements for the Orbitrap Fusion, Orbitrap Fusion Lumos, or Orbitrap ID-X are complete and accurate.

Signature _____ Date _____

Print name _____ Telephone _____

Email address _____

Principal instrument operator:

Print name _____ Telephone _____

Email address _____

Company _____ Telephone _____

Address _____

Address _____

City _____ State _____ Country _____

Sales order number _____

Note This form is intended to cover the essential components of your Orbitrap Tribrid Series installation. However, you must use the information in this guide and any additional information that your Thermo Fisher Scientific field service engineer provides to ensure the proper setup of your system. After receiving this form, the field service engineer contacts you to schedule the installation.

Contents

	Preface	vii
	Requesting Instrument Installation	viii
	Accessing Documentation	viii
	Viewing the Product Manuals	viii
	Viewing Online User Documentation	ix
	Providing Documentation Feedback	ix
	Special Notices, Symbols, and Cautions	ix
	Model Differences	x
	Contacting Us	xi
Chapter 1	Laboratory Requirements	1
	Electrical Requirements	1
	Environmental Requirements	3
	Gas Supply Requirements	4
	System Waste and Exhaust Requirements	5
Chapter 2	Installation Requirements	7
	Shipping Containers	7
	Electrical Ratings and Required Outlets	8
	Space and Load Requirements for the Workbenches	9
	Requirements for Workbench Clearances	10
	Recommended MS Layout	11
	Customer-Supplied Hardware	12
	Solvent Requirements and Recommendations	13
	Mobile Phase Requirements	13
	LC Additive Restrictions	13
	Recommended Solutions	14
Chapter 3	Instrument Shipments	15
	Inspecting Packages for Damage	15
	Reporting Damaged Packages	15
	Carrier to File Damage Claims	16

Chapter 4	Instrument Demonstration	17
Appendix A	Line Power Management	19
	Power Monitoring Devices	19
	Quality of Power	20
	Power Conditioning Devices	20
	Uninterruptible Power Supply	20

Preface

The *Orbitrap Tribrid Series Preinstallation Requirements Guide* is intended for the following Thermo Scientific™ mass spectrometers (MSs):

- Orbitrap Fusion™ (also known as Fusion™)
- Orbitrap Fusion Lumos™ (also known as Lumos™)
- Orbitrap ID-X™ (also known as ID-X™)

Before delivery and installation of the MS, read this guide carefully to ensure that your lab is fully prepared for the instrument. For a list of hardware differences for this MS series, see “[Model Differences](#)” on [page x](#).



CAUTION Operating an instrument or maintaining it outside the power and operating environment specifications described in this guide might cause failures of many types. The repair of such failures is specifically excluded from the Thermo Fisher Scientific standard warranty and service contract coverage.

Contents

- [Requesting Instrument Installation](#)
- [Accessing Documentation](#)
- [Providing Documentation Feedback](#)
- [Special Notices, Symbols, and Cautions](#)
- [Model Differences](#)
- [Contacting Us](#)

Requesting Instrument Installation

Before you schedule the instrument installation, ensure that the laboratory site meets the requirements in this guide.

❖ To schedule the instrument installation

1. Follow the procedure in [Inspecting Packages for Damage](#).

Do not open the shipping containers—the service engineer unpacks, inspects, and installs the system.

IMPORTANT If the instrument shipping container, ShockWatch™, or other indicators show any evidence of damage or mishandling during shipment, do NOT open the container. Follow the instructions in [Chapter 3, “Instrument Shipments.”](#)

2. If the instrument arrived safely, send the completed and signed “[Orbitrap Tribid Series Installation Request Form](#)” on [page iii](#) to your local office for Thermo Fisher Scientific San Jose products.

After receiving this form, the service engineer contacts you to schedule the instrument installation.

Accessing Documentation

The Orbitrap Tribid Series MS includes complete documentation.

- [Viewing the Product Manuals](#)
- [Viewing Online User Documentation](#)

For system requirements, refer to the release notes on the software DVD.

Viewing the Product Manuals

The Thermo Fisher Scientific service engineer installs the instrument control applications and the instrument manuals on the data system computer.

❖ To view the product manuals

From the Microsoft™ Windows™ taskbar, choose **Start > All Apps** (Windows 10) or **All Programs** (Windows 7) > **Thermo Instruments > model x.x**, and then open the applicable PDF file.

Viewing Online User Documentation

Visit the Thermo Fisher Scientific website for product manuals and more.

❖ To view user documentation from the Thermo Fisher Scientific website

1. Go to thermofisher.com.
2. Point to **Services & Support** and click **Manuals** on the left.
3. In the Refine Your Search box, search by the product name.
4. From the results list, click the title to open the document in your web browser, save it, or print it.

To return to the document list, click the browser **Back** button.

Providing Documentation Feedback

❖ To suggest changes to the documentation or to the Help

Complete a brief survey about this document by clicking the button below. Thank you in advance for your help.









Special Notices, Symbols, and Cautions

Make sure you understand the special notices, symbols, and caution labels in this guide. Most of the special notices and cautions appear in boxes; those pertaining to safety also have corresponding symbols. Some symbols are also marked on the instrument itself and can appear in color or in black and white. For complete definitions, see [Table 1](#).

Table 1. Notices, symbols, labels, and their meanings (Sheet 1 of 2)

Notice, symbol, or label	Meaning
IMPORTANT	Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the product.
Note	Highlights information of general interest.
Tip	Highlights helpful information that can make a task easier.

Table 1. Notices, symbols, labels, and their meanings (Sheet 2 of 2)




Notice, symbol, or label	Meaning
	Caution: Read the cautionary information associated with this task.
	Chemical hazard: Observe safe laboratory practices and procedures when handling chemicals. Only work with volatile chemicals under a fume or exhaust hood. Wear gloves and other protective equipment, as appropriate, when handling toxic, carcinogenic, mutagenic, corrosive, or irritant chemicals. Use approved containers and proper procedures to dispose of waste oil and when handling wetted parts of the instrument.
	Heavy object: The Orbitrap Tribrid Series MS, excluding its workbench, weighs over 227 kg (500 lb). Never try to detach and move the instrument from its workbench; you can suffer personal injury or damage the instrument.
	Risk of electric shock: This instrument uses voltages that can cause electric shock and personal injury. Before servicing the instrument, shut it down and disconnect it from line power. While operating the instrument, keep covers on.
	Risk of eye injury: Eye injury can occur from splattered chemicals, airborne particles, or sharp objects. Wear safety glasses when handling chemicals or servicing the instrument.
	Trip obstacle: Be aware of cords, hoses, or other objects located on the floor.

Model Differences

This table lists the required number of forepumps and the available options for the Orbitrap Tribrid Series MSs.

Instrument	Number of forepumps	Available options				
		EASY-IC ion source	EASY-ETD ion source	Advanced peak determination	1M Orbitrap resolution	UVPD laser module
Orbitrap Fusion Lumos	2	✓	✓	✓	✓	✓
Orbitrap Fusion	1	✓	✓	✓		
Orbitrap ID-X	1	✓				

Contacting Us

Contact	Email	Telephone	QR Code ^a
U.S. Technical Support	us.techsupport.analyze@thermofisher.com	(U.S.) 1 (800) 532-4752	
U.S. Customer Service and Sales	us.customer-support.analyze@thermofisher.com	(U.S.) 1 (800) 532-4752	
Global Support	<ul style="list-style-type: none"> ❖ To find global contact information or customize your request <ol style="list-style-type: none"> 1. Go to thermofisher.com. 2. Click Contact Us, select the country, and then select the type of support you need. 3. At the prompt, type the product name. 4. Use the phone number or complete the online form. ❖ To find product support, knowledge bases, and resources <p>Go to thermofisher.com/us/en/home/technical-resources.</p> ❖ To find product information <p>Go to thermofisher.com/us/en/home/brands/thermo-scientific.</p> 		

Note To provide feedback for this document, go to surveymonkey.com/s/PQM6P62 or send an email message to Technical Publications (techpubs-lcms@thermofisher.com).

^a You can use your smartphone to scan a QR Code, which opens your email application or browser.

Laboratory Requirements

Attention to the operating environment ensures continued high performance of the Thermo Scientific MS.

Contents

- [Electrical Requirements](#)
- [Environmental Requirements](#)
- [Gas Supply Requirements](#)
- [System Waste and Exhaust Requirements](#)

Electrical Requirements

Ensure that the laboratory site meets the requirements in [Table 2](#), including the wiring specification for the NEMA 6-15 or the 3-pole, CEE receptacle ([Figure 1](#)). For additional information, see [Appendix A](#), “[Line Power Management](#).”




CAUTION All devices connected between the power source and the instrument must be certified by recognized organizations for your country or territory (for example, UL, CSA, SEMKO, VDE, or TÜV). Such devices include the power cords, electrical outlets, circuit breakers, uninterruptible power supplies (UPSs), and so on.

To avoid an electric shock, always have a certified electrician install any new hardwired devices, such as electrical outlets.

Table 2. Electrical requirements

Parameter	Requirement
AC mains power system (nominal line power)	<ul style="list-style-type: none"> Data system: 120 or 230 Vac single phase, 50/60 Hz Instrument: 230 Vac single phase, 50/60 Hz <p>For optimal API source performance, the laboratory voltage must be at least 220 Vac. If it is not, order the Buck/Boost Transformer Kit (OPTON-01460) from Thermo Fisher Scientific San Jose. See also page 20.</p>
<p>Note The neutral-to-ground voltage must be 1 Vac or less.</p>	
Transient overvoltages	Category II
Earth ground	Hardwire all earth ground connections to the ground used for the main circuit breaker panel. Multiple external ground points can cause noise current to flow through the ground loop that is formed.
Overcurrent protection (circuit breakers)	Select suitably rated circuit breakers so that the connected equipment does not lose power by triggering a current overload condition. For added protection, install a surge protector at the input to the circuit breaker panel.
Electrical outlets (wall receptacles)	<p>Install outlets with a load rating that is suitable for the expected total current draw per outlet. See page 8.</p> <p>MS and forepump:</p> <ul style="list-style-type: none"> North America and other countries with the same configuration: NEMA 6-15 (250 Vac, 15 A) Other locations outside North America: 3-pole CEE (250 Vac, 16 A). The International Preinstallation Kit includes three CEE wall receptacles (Figure 1).
<p>IMPORTANT Connect the MS and each forepump to separate electrical outlets that have their own dedicated single-phase, circuit breakers rated 230 Vac, 15 A for North America and 230 Vac, 16 A for other countries and territories. For example, a dual-forepump system requires three separate, dedicated electrical outlets (one for the instrument and two for the forepumps). See page 11.</p>	
Power cord lengths	<ul style="list-style-type: none"> For the instrument and forepump: 2.5 m (8 ft) For the other system modules: 1.8 m (6 ft)

Figure 1. Pin-out for the 3-pole, CEE power receptacle (rated 230 Vac, 16 A)



Environmental Requirements

Ensure that the laboratory site meets these requirements.

Parameter	Requirement
Relative humidity	40–80%, noncondensing
Temperature	<ul style="list-style-type: none"> • Operating range: 15–27 °C (59–81 °F) • Optimum operating range: 18–21 °C (64–70 °F) • Ambient temperature fluctuations: Less than 1 °C or 2 °F over a one-hour period
Heat output (power)	<ul style="list-style-type: none"> • Data system^a: 889 W (3033 BTU/h) • Orbitrap Fusion and Orbitrap ID-X MS systems: 3800 W (12 966 BTU/h) • Orbitrap Fusion Lumos MS system: 5300 W (18 084 BTU/h) • (Optional) Thermo Scientific EASY-nLC™ 1200: 125 W (427 BTU/h) • (Optional) Thermo Scientific Dionex™ UltiMate™ 3000 RSLnano system: 265 W (904 BTU/h)

^a Mini-tower computer, monitor, and optional laser printer

Gas Supply Requirements

Ensure that the laboratory site meets these requirements.¹

Parameter	Requirement
Helium	<ul style="list-style-type: none">• Ultra-high-purity (UHP) (99.999%) with less than 1.0 ppm each of water, oxygen, and total hydrocarbons• 275 ±70 kPa (40 ±10 psi)• 9.22 L (0.325 ft³) typical daily consumption• Tubing provided: precleaned copper, 4.6 m (15 ft), 1/8 in. OD
Nitrogen #1	<ul style="list-style-type: none">• UHP (99.999%) with less than 3.0 ppm each of water, oxygen, and total hydrocarbons• 345 ±70 kPa (50 ±10 psi)• 46 L (1.62 ft³) typical daily consumption• Tubing provided: precleaned copper, 4.6 m (15 ft), 1/8 in. OD
Nitrogen #2	<ul style="list-style-type: none">• High-purity (99.5%)• 690 ±140 kPa (100 ±20 psi)• 11 500–26 700 L (406–943 ft³) typical daily consumption• Tubing provided: Teflon™ PFA, 4.6 m (15 ft), 1/4 in. OD

Note If you intend to use helium for sparging LC solvents, you must have a second tank and regulator.



CAUTION Contaminants in the laboratory gas lines can cause damage to the instrument. Make sure that all gas lines have been cleaned of all particulates and oils. You must use either copper (provided) or stainless steel for the UHP gas lines. Avoid using particulate filters, as they can be a source of contamination.

You are responsible for any damage to the instrument caused by contaminants introduced from your gas delivery system.

¹ If the location of your gas supplies exceeds the length of the provided tubing, you must provide suitable tubing for the gas connections. Consumption values are approximate when operating 24 hours and 7 days a week.

System Waste and Exhaust Requirements

In addition to meeting these requirements, the lab must have at least two fume exhaust systems.

Parameter	Requirement
#1: Solvent waste system	<ul style="list-style-type: none"> • API source solvent waste flow rate: maximum 30 L/min (64 ft³/h) • Solvent waste tubing: 1-3/8 in. OD, 1 in. ID, and 3 m (10 ft) long <p>Note Use a method of collecting the solvent waste that is compatible with the provided tubing and avoids pressure buildup in the API source.</p>
#2: Forepump fume exhaust system	<ul style="list-style-type: none"> • Forepump exhaust tubing: <ul style="list-style-type: none"> – (Orbitrap Fusion and Orbitrap ID-X MS systems) 0.75 in. OD, 0.5 in. ID, 3 m (10 ft) long – (Orbitrap Fusion Lumos MS system) 0.75 in. OD, 0.5 in. ID, 6 m (20 ft) long • Initial inrush flow rate: 3 L/min (6.4 ft³/hr) • Continuous flow rate: 1 L/min (2 ft³/hr)



CAUTION Do not run the forepump exhaust lines vertically near the forepump. Solvents and oils can condense in the line and flow back into the pump, causing pump damage and diminished pump capability. To maintain forepump integrity, route the exhaust tubing from the exhaust port down to the floor, not from the forepump vertically toward the ceiling. Run the hose at floor level for a minimum of 2 m (79 in.) before it reaches the external exhaust system.

The dedicated exhaust system must be actively vented and be able to accommodate the periodic purging of the accumulated solvents in the forepump. The frequency of the purging depends on the throughput of the instrument.

1 Laboratory Requirements

System Waste and Exhaust Requirements

Installation Requirements

Before the Thermo Fisher Scientific service engineer installs the Orbitrap Tribrid Series MS, read and understand these specific requirements and recommendations.

Contents

- Shipping Containers
- Electrical Ratings and Required Outlets
- Space and Load Requirements for the Workbenches
- Requirements for Workbench Clearances
- Recommended MS Layout
- Customer-Supplied Hardware
- Solvent Requirements and Recommendations

Shipping Containers

Table 3 lists the overall dimensions and weights of the shipping containers that are too large to carry by hand. Ensure that the width of all doorways and hallways is a minimum of 1.4 m (54 in.). Transport the container with the MS assembly (container #1) to the lab and do not open it until the service engineer is on site.

Table 3. Shipping container dimensions and weights^a

Container	Size	Weight	Content
1	175 × 119 × 160 cm (<i>l × w × h</i>) (69 × 47 × 63 in.) (<i>l × w × h</i>)	493 kg (1087 lb)	MS assembly
2	109 × 76 × 125 cm (<i>l × w × h</i>) (43 × 30 × 49 in.) (<i>l × w × h</i>)	68 kg (150 lb)	Syringe pump, divert/inject valve, and data system components
3	79 × 51 × 51 cm (<i>l × w × h</i>) (31 × 20 × 20 in.) (<i>l × w × h</i>)	66 kg (146 lb)	One forepump

^a The height and weight measurements include the shipping pallets.

Electrical Ratings and Required Outlets

This table lists electrical ratings for the individual modules and the number of required electrical outlets for the MS and data system. See [page 2](#) for the power cord lengths.

Note See “[Model Differences](#)” on [page x](#) for which instruments require two forepumps. Depending on whether the MS requires one or two forepumps, ensure that there are at least 10 or 11 outlets available.

Module	Voltage (Vac)	Current (A)	Required outlets ^a
MS system			
Orbitrap Tribrid Series MS	230	10	1
SOGEVAC™ SV65BI FC forepump ^b	230	8	1 or 2
Thermo Scientific divert/inject valve	110–220	1.5	1
Thermo Scientific syringe pump	110/220	0.2/0.1	1
Data system			
Computer, mini-tower	100–240	5.4	1
Monitor	100–240	1.5	1
Ethernet switch	100–240	< 1	1
(Optional) Laser printer	110 –or– 220	8.6 –or– 4.2	1
Optional devices^c			
High-intensity lamp (for instrument maintenance)	–	–	1
Laboratory stereoscope (for inspecting fused-silica parts)	–	–	1

^a Remember to add the number of outlets required for your LC system.

^b The values are for one forepump.

^c Refer to the equipment manual for the electrical ratings.

Space and Load Requirements for the Workbenches

You must provide workbenches for the LC and data system modules. Workbenches must have a load capacity of at least *twice* the combined weight of all expected devices (Table 4).

The workbench attached to the Orbitrap Tribrid Series MS measures approximately 132 cm wide, 74 cm high, and 78.8 cm deep (52 in. *w*, 29 in. *h*, and 31 in. *d*). It has a load capacity of 508 kg (1120 lb).



CAUTION Heavy object. The Orbitrap Tribrid Series MS, excluding its workbench, weighs over 227 kg (500 lb). Never try to detach and move the instrument from its workbench; you can suffer personal injury or damage the instrument. For additional information, contact your local Thermo Fisher Scientific service engineer.

Note Place the forepumps under the workbench, which was designed to minimize their noise without the need for sound-dampening boxes. Do not place the pumps on a shelf or other surface connected to the workbench.

Table 4. Space and load requirements for the system modules^a

Equipment	Width (<i>w</i>) cm (in.)	Height (<i>h</i>) cm (in.)	Depth (<i>d</i>) cm (in.)	Weight kg (lb)
LC/MS				
(Optional) EASY-nLC 1200 ^b	36 (14.2)	45 (17.7)	38 (15.1)	35 (77)
(Optional) UltiMate 3000 RSLnano system ^b	42 (16.5)	88 (34.6)	51 (20)	70 (154)
Orbitrap Tribrid Series MS ^c	127 (50)	70.31 (27.7)	76.7 (30)	245 (540)
Oerlikon Leybold Vacuum™ SOGEVAC SV65BI FC forepump (each)	32 (12.6)	26.4 (10.4)	48 (19)	52 (115)

Data system

Select a workbench that can hold the mini-tower computer, wide-screen monitor, keyboard, Ethernet switch, and optional laser printer.

^a Approximate values

^b These values exclude the solvent bottles and tubing.

^c These values exclude the workbench, API source, syringe pump, and modular valve.

Requirements for Workbench Clearances

Follow these clearance guidelines for the workbenches (Figure 2):


- Place the data system and MS workbenches close to each other to prevent strain on the interconnecting Ethernet communications cables.



CAUTION Safety and EMC regulations require the use of Category 5e shielded Ethernet communications cables, maximum 3 m (10 ft) long.

- For the LC and MS systems, allow for a minimum vertical clearance of 92 cm (36 in.) between the top of the system and any shelves above it.
- For the MS system, allow for these minimum horizontal clearances for proper air circulation and for the installed gas lines and tubing:
 - 25 cm (10 in.) behind the instrument
 - 61 cm (24 in.) on the right side of the instrument (between the MS and any solid barrier, such as a wall)

Figure 2. Top view (footprint) and recommended placement of the workbenches




CAUTION Trip hazard. Whenever possible, provide space under the workbench for the forepumps. If placed in front of the MS, the forepumps become a trip hazard.

Recommended MS Layout

Figure 3 shows an example of the recommended layout for the Orbitrap Fusion Lumos MS (dual-forepump system). Not all connections are shown; for installation instructions, refer to the Getting Connected Guide. For the other instruments, see “Model Differences” on page x for the required number of forepumps.

IMPORTANT The MS and forepumps must connect to separate electrical circuits.

Figure 3. Workbench for the Orbitrap Fusion Lumos MS



Customer-Supplied Hardware

Your purchase does not contain all parts, materials, or tools that are required for installation. To complete the installation, you must provide these additional parts.

IMPORTANT If the location of your gas supplies exceeds the length of the provided tubing, you must provide suitable tubing for the gas connections. See “[Gas Supply Requirements](#)” on page 4.

Item	Description
Communications cable	Connects from the Ready Out and Injection Hold pins (peripheral control connector) on the MS to a device not controlled by one of the Thermo Scientific MS applications, such as the Xcalibur™ data system.
Gas line fittings, UHP helium and UHP nitrogen	Connects the 1/8 in. ID copper or stainless steel tubings to the helium and nitrogen (#2) gas supplies.
Gas line fitting, HP nitrogen	Connects the 1/4 in. OD Teflon PFA tubing to the nitrogen (#1) gas supply.
LC system	IMPORTANT If your LC system is manufactured by another company, you must verify that it is suitable for use with the Orbitrap Tribrid Series MS. The output (start) signal from the external device must be <i>Normally Hi</i> (+5 Vdc) and momentarily go to <i>Low</i> . If you cannot configure the external device to go from <i>Normally Hi</i> to <i>Low</i> momentarily, you cannot use it with the Orbitrap Tribrid Series MS.
LC system, solvents	Used by the service engineer to calibrate the system during the initial setup. For details, refer to the appropriate LC manual.

Solvent Requirements and Recommendations

Before you start using the Orbitrap Tribrid Series MS, make sure that your lab has the necessary solvents and recommended solutions and equipment. Installation and maintenance of the MS requires LC/MS-grade or better methanol and water (see [Table 5](#)).

- [Mobile Phase Requirements](#)
- [LC Additive Restrictions](#)
- [Recommended Solutions](#)

Mobile Phase Requirements

Ensure the following:

- Use LC/MS-grade or better solvents and additives.
- Use dedicated glassware that you clean and store in the lab between uses.
- Clean solvent bottles before refilling them.
- Do not use plastic materials for preparing or storing solvents.
- Do not use Parafilm™ to seal mobile phase solution bottles.

LC Additive Restrictions

Note If your workflow requires the use of any of these reagents, contact Thermo Fisher Scientific Technical Support for advice.

Limit the type of reagents you use:

- Because it reacts with the PEEK material, do not use tetrahydrofuran (THF).
- Because these might cause corrosion of the API source, do not use the following:
 - Alkali-metal bases, such as sodium hydroxide (NaOH)
 - Inorganic acids, such as hydrogen chloride (HCl), phosphoric acid (H₃PO₄), or sulfuric acid (H₂SO₄)
- Because these reagents can suppress ionization, do not use the following:
 - Detergent solutions
 - Involatile buffer solutions, such as borate, citrate, or phosphate
 - Surfactants or other surface-active agents, such as polyethylene glycol (PEG), sodium dodecyl sulfate (SDS), or Triton™ X-100
 - Trifluoroacetic acid (TFA) at concentrations (volume/volume %) greater than 0.1

2 Installation Requirements

Solvent Requirements and Recommendations

Recommended Solutions

For a wide variety of LC/MS- and UHPLC/MS-grade solvents and consumables for purchase, visit www.fishersci.com.



CAUTION Avoid exposure to potentially harmful materials.

By law, producers and suppliers of chemical compounds are required to provide their customers with the most current health and safety information in the form of Material Safety Data Sheets (MSDSs) or Safety Data Sheets (SDSs). The MSDSs and SDSs must be freely available to lab personnel to examine at any time. These data sheets describe the chemicals and summarize information on the hazard and toxicity of specific chemical compounds. They also provide information on the proper handling of compounds, first aid for accidental exposure, and procedures to remedy spills or leaks.

Read the MSDS or SDS for each chemical you use. Store and handle all chemicals in accordance with standard safety procedures. Always wear protective gloves and safety glasses when you use solvents or corrosives. Also, contain waste streams, use proper ventilation, and dispose of all laboratory reagents according to the directions in the MSDS or SDS.

Table 5. Recommended solutions

Product	Grade	Size ^a	Part number
Solvents			
2-Propanol	Optima™ LC/MS	Amber glass, 1 L	A461-1
Acetonitrile	Optima UHPLC/MS	Amber glass, 1 L	A956-1
Methanol	Optima UHPLC/MS	Amber glass, 1 L	A458-1
Water	Optima UHPLC/MS	Amber glass, 1 L	W8-1
Blends			
0.1% Formic acid in acetonitrile	Optima LC/MS	Amber glass, 1 L	LS120-1
0.1% Formic acid in water	Optima LC/MS	Amber glass, 1 L	LS118-1
Additives			
Acetic acid (modifier)	Optima LC/MS	Ampule, 1 mL	A113-1AMP
Ammonium acetate (modifier)	Optima LC/MS	Amber glass, 50 g	A114-50
Ammonium formate (modifier)	Optima LC/MS	Amber glass, 50 g	A115-50
Formic acid (modifier)	Optima LC/MS	Ampule, 10 × 1 mL	A117-10X1AMP

^a Size for the stated part number

Instrument Shipments

Electronic equipment carriers that specialize in the handling and transport of delicate machinery ship the Orbitrap Tribrid Series MS to your site. When the instrument arrives, move it to a protected indoor location. If you have questions about moving your instrument, contact your local office for Thermo Scientific San Jose products (see “[Contacting Us](#)” on [page xi](#)).

On the occasion when shipments arrive that appear damaged, take these immediate actions.

Contents

- [Inspecting Packages for Damage](#)
- [Reporting Damaged Packages](#)
- [Carrier to File Damage Claims](#)

Inspecting Packages for Damage

When the shipment packages arrive, visually inspect them for any damage.

❖ To visually inspect for damage

Carefully inspect for obvious damage or evidence of rough handling.

If the instrument shipping container, ShockWatch, or other indicators show visible evidence of damage or mishandling, do NOT open the container.

Reporting Damaged Packages

Follow the next procedure, and then call your Thermo Fisher Scientific sales representative for further instructions.

IMPORTANT Freight insurance requires that you note obvious damage on the receiving documents. Thermo Fisher Scientific does not accept liability for damage if materials are received with obvious damage AND the damage is not recorded on the receiving documents.

❖ **To record damages on the receiving documents**

1. Read the information in [Carrier to File Damage Claims](#) to determine which parties might be responsible for filing a claim against the carrier.
2. On all copies of the receiving documents, note any apparent external damage and briefly describe the extent of the damage.
3. Have the driver sign or initial next to your comments to signify agreement with your observations.
4. Report the list of damages to your Thermo Fisher Scientific representative.

Carrier to File Damage Claims

If the instrument is damaged in transit, the shipment method determines the party who assumes the risk of damage and files a claim against the carrier—Thermo Fisher Scientific or the purchaser. The sales agreement or the sales quote identifies the shipment method for instruments sent from the San Jose, CA site.

[Table 6](#) lists the party who files the damage claim against the carrier for instruments damaged in transit based on the shipment method.

Table 6. Shipment methods for delivery from the San Jose, CA site to domestic and international destinations

Destination	Shipment method	Party responsible for filing a damage claim
Domestic (United States)	Destination or Origin—Thermo Fisher Scientific pays the carrier.	Thermo Fisher Scientific
	Origin—The purchaser pays the carrier.	Purchaser
International	Carriage Paid To (CPT) named destination ^a	Purchaser
	Carriage and Insurance Paid (CIP) to named destination ^b	Thermo Fisher Scientific

^a Unless specified differently, Thermo Fisher Scientific uses this shipment method for international shipments.

^b Under special circumstances, Thermo Fisher Scientific uses this shipment method for international shipments.

Instrument Demonstration

When your new Orbitrap Tribrid Series system is on site and ready for installation, the service engineer unpacks and installs it.

IMPORTANT Do not use the new system for sample analysis until the installation is complete and you have signed the Acceptance Form.

During the installation, the service engineer demonstrates the following:

- Basics of the instrument operation and routine maintenance
- Marketing specifications that are in effect when you purchased the instrument

Tip To receive maximum benefit from this on-site training opportunity, plan for the system operators to be available during the entire installation process.

4 Instrument Demonstration

Line Power Management

The quality of line power (ac mains power system) delivered to the MS can affect its performance and longevity. You are responsible for correcting any line power problems. Contact Thermo Fisher Scientific for assistance in monitoring the line voltage in your lab and in selecting a line conditioner.



CAUTION To support compliance and safety requirements, all devices connected between the power source and the MS must be certified by recognized organizations for your country or territory (for example, UL, CSA, SEMKO, VDE, or TÜV).

Such devices include the power cords, electrical outlets, circuit breakers, uninterruptible power supplies (UPSs), and so on.

Contents

- [Power Monitoring Devices](#)
- [Quality of Power](#)
- [Power Conditioning Devices](#)
- [Uninterruptible Power Supply](#)

Power Monitoring Devices

Several devices are available to monitor the quality of the line power. These devices provide a continuous record of line performance by analyzing and printing out data for the three most common voltage disturbances. Monitor the power line 24 hours a day for seven consecutive days. If inspection of the printout indicates disturbances, stop the test and take corrective action. Monitor the power again as previously described.

A power line disturbance analyzer detects and records most types of line power problems. The Dranetz™ system¹ is an example of a suitable analyzer. In some countries, you can rent power line analyzers from electrical equipment suppliers.

¹ Thermo Fisher Scientific does not endorse any power monitoring company, nor does it endorse products other than its own. Companies and products listed in this guide are given as examples only.

Quality of Power

Before the service engineer arrives to install your system, the line voltage must be stable and within the recommended specifications. The line voltage must be free of fluctuations due to slow changes in the average sags, surges, transients, or voltage. Establishing the quality of power supplied to the LC/MS system is very important for these reasons:

- Constant high line voltage, impulses, or surges in voltage can cause overheating and component failures.
- Constant low line voltage or sags in voltage can cause the equipment to function erratically or not at all.
- Transients—even a few microseconds in duration—can cause electronic devices to degrade or fail catastrophically, shortening the lifetime of the equipment.

Power Conditioning Devices

If the power regulation is good but the power line disturbance analyzer shows transient voltages, an isolation/noise-suppression transformer can resolve the problem. For both transient and regulation problems, consider the use of power conditioners to control these problems.



CAUTION Any conditioning device installed with the MS must be able to handle the potentially high currents that are drawn during the initial startup of the system. The maximum system inrush (start) current for one forepump (SOGEVAC SV65BI FC) is 12 A with an average duration of less than 1 s. Therefore, this initial energy demand from the ac power line is very low.

When the line voltage is free from voltage sags, surges, and impulses but is more than 10 percent outside of the voltage specifications, a buck/boost transformer can lower (buck 10 percent) or raise (boost 10 percent) the line voltage as appropriate for the rated voltage.



CAUTION Instruments installed in areas with 208 Vac power can experience voltage sags during high use periods that might place the line voltage below the operating parameters discussed in this section. In this case, protect the instrument by ordering a Buck/Boost Transformer Kit (P/N OPTON-01460) to ensure that power stays within the specified parameters.

Uninterruptible Power Supply

If your local area is susceptible to corrupted power or power disruptions, Thermo Fisher Scientific recommends installing an uninterruptible power supply (UPS) in the lab. For additional information, visit www.thermopowervar.com. For North American labs, you can order a Powervar™ UPS (5.2 kVA, 6.0 kVA, or 8.0 kVA) from [Unity Lab Services](#).