



Thermo Scientific
BioPharma Finder Software

Workflow driven informatics for complete biotherapeutic characterization

- Intact Protein Mass Analysis
- Subunit Mass Analysis
- Peptide Mapping

Thermo
SCIENTIFIC

Taking too much Time to Analyze your Data?

Simple software workflows guide you along the biocharacterization pathway all centrally managed with the Protein Sequence Manager feature, the integrated library for your sequences.

Struggle to Confidently Characterize your Intact Proteins?

Novel deconvolution algorithms ensure you determine the correct molecular weight, confidently identify glycoforms and accurately measure drug-antibody ratio (DAR).

Getting Lost in your Peptide Maps?

Combine in-depth knowledge of post-translational modifications (PTMs), disulfide bonds, and low level impurities with unique data visualization tools to make sure you never lose any data again.



Chances are you think you know your primary sequence, BioPharma Finder ensures you do.

Complete Protein Characterization...


Complete characterization is necessary for the reproducible and safe production of biologics. However this can be complex, from challenging sample preparation to lengthy data-processing routines.

From intact protein mass analysis to peptide mapping, we have created integrated software workflows delivering comprehensive interpretation and visualization allowing you to characterize your biologics with ease.

...trusted by the Experts

Orbitrap Power...



Thermo Scientific™ BioPharma Finder™ integrated software allows you to harness the power of the industry leading Thermo Scientific™ Orbitrap™ technology.



Deconvolve complex spectra

DAR Accuracy. Every time.

Why acquire all of your data and use only a small section to report a DAR? Sliding Window ensures you report the true DAR based upon all your data, not a subset.



Automatically annotate peptides

PTMs? See everything.

From powerful chromatographic shading, predicted vs. experimental fragmentation spectra to comprehensive sequence variant searching and *de novo* sequencing, you can be confident in your sequence.

...with
BioPharma Finder

Superior Deconvolution

Accurately deconvolute fully or partially resolved data using tested novel Xtract and ReSpect algorithms to provide confident molecular weights of proteins in both acidic and native conditions.




Historically it took weeks to evaluate data, ruling out many false positives. BioPharma Finder is a **unique solution** that has **efficient workflows** meaning that you don't need to spend hours fine tuning parameters to get the **right answer the first time**. For comparative quantitative modification analyses, we have gone from **two weeks** of analysis down to **one day** for multi-batch samples.

Guilong Charles Cheng,
Alexion Pharmaceutical Inc.






Intact Protein & Subunit Mass Analysis


The development of monoclonal antibodies (mAbs) and ADCs is challenging. There is now an easier path to take for the analysis of critical quality attributes in these challenging biologics – BioPharma Finder Software.



Sliding window is a novel algorithm for the interpretation of LC/MS data of intact proteins. Moving averaged spectrum across the chromatogram generates deconvoluted spectra, resulting in a single deconvoluted spectrum of all protein components. This is supported for both high resolution (Xtract) and low resolution (ReSpect) data.



Native separations can pose a challenge, especially when it comes to analyzing the resulting data. The novel sliding window algorithm ensures that all of your data is used, not just a slice, for deconvolution making sure you report the accurate DAR value.




Whether resolved or unresolved, your proteins are in safe hands with **superior deconvolution** algorithms (Xtract & ReSpect) to determine accurate information on various protein properties, such as intact molecular mass, glycosylation forms, amino acid sequence, posttranslational modifications, disulfide linkages as well as higher order structural information.

The software **provides confident deconvolved** molecular weight of proteins in both denaturing and native conditions and protein sequence matching **identifies N-linked glycoylations** and other common modifications via intact mass analysis, helping you to understand your biologics heterogeneity.




Simple, intuitive user interface



Deconvoluted high resolution Orbitrap MS spectra provides **low-ppm** mass accuracy molecular weights and relative quantitation of glycoform heterogeneity.

Are they really the same? When it comes to biosimilars you can't afford to miss a thing. Demonstrate equivalence or identify differences through unique comparison and visualization tools.




Easy, yet powerful mirror plotting



Peptide Mapping

Small changes in peptide sequence can result in changes to higher order protein structures, compromising biotherapeutic safety and efficacy. Complete peptide map coverage guarantees you see the whole picture.



Batch processing capabilities together with an interactive wizard for editing process methods allows easy and efficient data analysis with maximum throughput.

Thermo Scientific BioPharma Finder is workflow-driven, enabling beginner and expert users alike to obtain complete peptide coverage easily with high confidence.

Full MS, MS² experimental vs. predicted spectra and fragment coverage maps for each identification help to **increase confidence in your peptide sequence**. Compare extracted ion chromatograms for batch evaluation, with relative quantitation information for each sample.

Protein	Residue #	Modification	Category	Normalized Time Shift	Predicted Time Shift	Confidence	Recovery	% Abundance Sample	% Abundance Reference
1	48	Y48+125.8970	Unknown Modification	27.1%	N/A	100.0%	38.2%	0.0543%	0.0536%
2	157	N157(N157S)	Sequence Variant	1.3%	2.1%	100.0%	21.1%	0.6265%	0.0002%
3	166	D166+Isomerization	Modification	-3.0%	-1.0%	100.0%	13.2%	0.1201%	0.1285%
4	162	N162+Deamidation	Modification	1.1%	N/A	100.0%	7.8%	11.5107%	9.6140%
5	255	N255+Oxidation	Modification	-4.8%	-10.0%	100.0%	7.7%	81.0324%	75.4640%
6	318	N318+NH3 loss	Modification	5.8%	8.0%	100.0%	25.4%	2.1513%	2.1317%
7	384	W384+Double Oxidation	Modification	-8.2%	-8.0%	100.0%	14.9%	0.3724%	0.3042%

Interactive Modification Summary

Identify and localize glycoforms, deamidations, oxidations, other sequence variants and PTMs.

Comprehensive results table providing amino acid sequence confirmation with mass tolerance, modification, identification, retention time and confidence information. Real-time filtering to optimize your data, all exportable.



Informative peptide fragment map coverage.

Protein Coverage Map

Sequence Coverage Map

Proteins	Number of MS Peaks	MS Peak Area	Sequence Coverage
1:Light 1	314	23.3%	100.0%
2:Heavy 1	702	59.9%	100.0%



Color code for signal intensity
 >7.3e+006 >4.9e+005 >3.3e+004 >2.2e+003 >1.5e+002

Informative and interactive sequence coverage map, colored by peptide intensity.

Multiple interactive plot display options so you can see what you want, when you need it.

The ability to perform Sequence Variant Searching, disulfide link analysis and de novo sequencing with quantitation of post-translational modifications assures you won't miss anything.


“Being able to sort for sequence variants gives me the confidence that I'm not overlooking any.”
 Shannon Williams,
 Coherus Biosciences



Annotated peptide map



Powerful visualization tools such as chromatographic shading afford increased efficiency in data mining and interpretation.



Rapid User-friendly Sample Preparation

Whether you are analyzing intact proteins or peptides, we have sample preparation and chromatographic columns that allow you to get the most out of your samples.

Robust Separations

Purpose-built biocompatible systems yield precise retention times from sample to sample and increase confidence in peak assignments.



High Resolution Accurate Mass MS

Obtain excellent resolution, mass accuracy and dynamic range for the most challenging samples from the gold-standard in mass spectrometry.

Thermo Scientific BioPharma Finder


From novel deconvolution algorithms generating comprehensive results to easy to understand data visualization tools, BioPharma Finder helps you to take the right path. Every time.



Find out about BioPharma Finder Software
www.thermofisher.com/BioPharmaFinder

f facebook.com/ChromatographySolutions t twitter.com/ChromSolutions a analyteguru.com

Upgrade Your Maps! Scan the QR code or follow the URL to access the latest information on BioPharma Finder Software.



thermofisher.com

©2016 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand