| Parameters | | | | | | | | |
|-----------------------|----------------------------------|---|--|--------|--------|----------|---------|------------------|
| Classification | | Parameter | Measurement range | | | | | Measurement time |
| | | | Unit (A | A) | | Unit (B) | | (min.) |
| | Enzymes | ALP | 14 - 1183 | U/L | 0.23 - | 19.76 | μ Kat/L | 4 |
| | | v-AMYL | 100 - 2500 | U/L | 1.67 - | 41.75 | μ Kat/L | 3 |
| | | CPK | 10 - 2000 | U/L | 0.17 - | 33.40 | μ Kat/L | 4 |
| | | GGT | 10 - 1200 | U/L | 0.17 - | 20.04 | μ Kat/L | 5 |
| | | GOT/AST | 10 - 1000 | U/L | 0.17 - | 16.70 | μ Kat/L | 4 |
| | | GPT/ALT | 10 - 1000 | U/L | 0.17 - | 16.70 | μ Kat/L | 4 |
| | | v-LIP | 10 - 1000 | U/L | 0.17 - | 16.70 | μ Kat/L | 5 |
| | | LDH | 50 - 900 | U/L | 0.84 - | 15.03 | μ Kat/L | 2 |
| | General chemistry | ALB | 1.0 - 6.0 | g/dL | 10 - | 60 | g/L | 6 |
| | | BUN | 5.0 - 140.0 | mg/dL | 1.79 - | 49.98 | mmol/L | 4 |
| | | Ca | 4.0 - 16.0 | mg/dL | 1.00 - | 4.00 | mmol/L | 4 |
| E | | CRE | 0.2 - 24.0 | mg/dL | 18 - | 2122 | μ mol/L | 5 |
| | | GLU | 10 - 600 | mg/dL | 0.6 - | 33.3 | mmol/L | 6 |
| | | IP | 0.5 - 15.0 | mg/dL | 0.16 - | 4.84 | mmol/L | 5 |
| | | Mg | 0.2 - 7.0 | mg/dL | 0.08 - | 2.88 | mmol/L | 4.5 |
| | | NНз | 10 - 500 | μg/dL | 7 – | 357 | μ mol/L | 2 |
| | | TBIL | 0.2 - 30.0 | mg/dL | 3 - | 513 | μ mol/L | 6 |
| | | TCHO | 50 - 450 | mg/dL | 1.29 - | 11.64 | mmol/L | 6 |
| | | TCO ₂ | 5 - 40 | mmol/L | 5 - | 40 | mmol/L | 5 |
| | | TG | 10 - 500 | mg/dL | 0.11 - | 5.65 | mmol/L | 4 |
| | | TP | 2.0 - 11.0 | g/dL | 20 - | 110 | g/L | 6 |
| | | UA | 0.5 - 18.0 | mg/dL | 30 - | 1071 | μ mol/L | 4 |
| | Electrolytes (on only one slide) | Na | 75 – 250 | mEq/L | 75 – | 250 | mmol/L | |
| | | K | 1.0 - 14.0 | mEq/L | 1.0 - | 14.0 | mmol/L | 1 |
| | (on only one side) | CI | 50 - 175 | mEq/L | 50 - | 175 | mmol/L | |
| Immunological Test | | vc-CRP | 0.3 - 7.0 | mg/dL | 3 - | 70 | mg/L | 5 |
| Comprehensive S-panel | | TP, ALB, ALP, GLU, TBIL, IP, TCHO, GGT, GPT/ALT, Ca, CRE, BUN | | | | | | |
| Pre-Surgical S-panel | | TP, ALP, GLU, GPT/ALT, CRE, BUN | | | | | | |
| Calculated | Calculated parameters | | GLOB, ALB/GLOB, BUN/CRE, Na/K, Anion Gap | | | | | |

^{*} Unit (A) or (B) is available * There are parameters which may not be available in your area. For details please contact your local distributor. * v-LIP, v-AMYL and vc-CRP slides are designed for veterinary use.

Other slides are designed for use in analyzing human blood. Measurement results may vary depending on the animal species. For details, please contact your local distributors.

| Specifications of NX500V IC | | | | |
|-----------------------------|--|--|--|--|
| Measurement test | Colorimetry 23 tests Electrolytes 3 tests | | | |
| Throughput | Colorimetry 120 test/ hour Combined 128 test/ hour | | | |
| Number of sample rack | 1 | | | |
| Number of incubator cell | Colorimetry 12, Electrolytes 1 | | | |
| Measurement time | Colorimetry 2 to 6 minutes/ test Electrolytes 1 minute/ 3 tests (Na, K, Cl) | | | |
| Sample type | Plasma, Serum, Whole blood | | | |
| Sample volume | Colorimetry 10µL/ test Electrolytes 50µL/ 3 tests (Na, K, Cl) | | | |
| Data transmission to PC | USB 2.0 or RS-232C Serial D-Sub 9 pin - 9 pin cross cable | | | |

| Data print | Thermal printer |
|-------------------------|--------------------------------------|
| Electrical requirements | AC 100 - 240V, 50/ 60 Hz, 2.5 - 1.1A |
| Display | 5.7 inch color touch panel |
| Dimensions | 470 (W) x 360 (D) x 420 (H) mm |
| Weight | Approx. 24kg |
| Operating temperature | 15 to 32°C (59 to 89 °F) |
| Operating humidity | 30 to 80% RH (no vapor condensation) |
| | |

Option Barcode reader

Barcode reader is available as option item to read sample ID on sample tube.

DRI-CHEM NX500V IC (Product name: FUJI DRI-CHEM NX500V IC)

NX500V IC and AU10V Clinical Chemistry tests providing total in-house diagnostic testing





| Specifications of DRI-CHEM IMMUNO AU10V | | | | |
|---|--|--|--|--|
| Dimensions | 238 (W) x326 (D) x 310 (H) mm | | | |
| Weight | Approx. 7 kg | | | |
| Operating temperature | 15 to 30°C (59 to 86°F) | | | |
| Measurement time | Approx. 10 min (without dilution) | | | |
| Required sample volume | basic test: 100 μL dilution test: 20 μL | | | |
| Test parameter | v-PRG (Progesterone) v-T4 (Thyroxine (T4)) vc-TSH (Thyroid-Stimulating Hormone) v-COR (Cortisol) v-BA (Bile Acids) | | | |

The specifications and appearance of the present brochure may be changed without prior notification in order to improve the system. Please be sure to read the instruction manual carefully for proper use of the equipment. There are some different specifications between the veterinary analyzer and the human one. Please ask your local distributor for the availability.

Read more about our Medical Solutions: www.fujifilm.com



26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN http://www.fujifilm.com/products/medical/

Ref. No. NX500 VET E R2 (SK-18-07-F1079) ©2018 FUJIFILM Corporation



DRI-CHEM NX500V IC

For veterinary



18 tests in 12 minutes from a whole blood sample (15 chemistry and 3 electrolytes tests)



Professional results from effortless and efficient operation.

A maximum of 20 test slides in the slide cartridge will be automatically processed in a single run. Proceeding sample or test can be set up subsequently without waiting for the previous run to finish.

Proven quality from experience

Fujifilm has been developing analysis devices that work with dry chemistry for laboratory diagnostics since 1984. It has been the market leader in Japan ever since.



NEW Built-in centrifugal blood separator for direct use of whole blood samples.

The NX500V IC has a built-in centrifugal separator. After setting the whole blood sample and slides and pressing the START button, one has only to wait for the measurement results thereafter. The built-in centrifugal separator can cut the pre-treatment process of whole blood samples.











to external


centrifugal


Press the START button

Specially designed IC CUP for high quality plasma

1) IC CUP (Integrated Centrifuge Cup)

The blood cells, which have higher specific gravity than separating agent, are separated and trapped at the outer side of the separating agent after centrifugation. Only the plasma passes along the groove and gets accumulated at the bottom of the cup after the separation.





2) Straw


Function for efficient mixing of whole blood with heparin and dispensing of sample to the cup. A special straw comes with the cup.



Easy

■ Intuitive and easy operation on an LCD touch screen

FUJIFILM placed special emphasis in designing a simple, easy-to-use operation screen.





■ Microsample testing A very useful feature especially in small animals

Each test requires 10µL of sample. (CRP: 5µL/test, ISE: 50µL/test, dead volume is necessary using automatic measurement.) If there is less sample, measurement can be done by manual pipetting.



No calibration

The FUJI DRI-CHEM system uses dry type reagents which are very stable. Calibration is not necessary. Variability in the reagent lots are adjusted using the QC card.



■ Automatic dilution

Dilution, a time consuming process, is also automated. The only operation is to input the dilution ratio.

Easy maintenance

The analyzer is low maintenance. It does not require special skills or training to maintain.

Compact analyzer which requires only a small space



Fast

Quick measurement

The NX500V IC analyzer can quickly measure 23 different colorimetric parameters and 3 electrolytes. It has a 128 tests/hour throughput.



* In case of whole blood samples, it will take another 2 minutes for the centrifugation for each sample.

■ One pack panels (S-panels)

S-Panels are prepacked slides which can be easily mounted in the slide cartridge in just one step. This improves work efficiency and waste reduction.






■ Continuous measurement

As soon as the pipetting procedure for the first sample is completed, the slides for the next sample can already be set and the test started. This minimizes waiting time for your patients.

NX500V IC processing cycle in minutes for 18 clinical parameters per patient sample



^{*} In case of whole blood samples, it will take another 2 minutes for the centrifugation for each sample.

Height: 420 mm


Width: 470 mm


DRI-CHEM NX500V IC

► Accurate & Reliable

NX500V IC can measure lipase activity for canine

Lipase activity, which is useful for the screening of pancreatitis, can also be easily measured. This contributes to quick screening of pancreatitis onsite.







Citation: J. Vet. Med. Sci. 73(11): 1481-1483, 2011 Katsumi Ishioka J-vet vol. 291 2011.6 pp13-15, pp25-29

Precise results

The FUJI DRI-CHEM slide reagent has high reliability and stability brought by fine chemical technology cultivated through the long history of Fujifilm. The high quality of Fujifilm reagents is demonstrated by a comparison with the standard measurement methods.







50

Depth: 360 mm

00

50