



Pivotal Living™



Bluetooth™ Smart Scale

Welcome. Let's get started:



For more comprehensive instructions on how to setup, sync, and use your Pivotal Living™ Smart Scale, please visit:

www.PivotalLiving.com/Start

What's included:

- Pivotal Living™ Smart Scale
- AAA batteries (4)
- Surface caps (4)
- Instruction booklet

Measures:

- Weight (kg/lb)
- BMI (Body Mass Index)
- BMR (Basal Metabolic Rate)
- Body Fat %
- Lean body mass

Features and Functions:

- Tracks data for up to 16 users.
- Tempered glass (Supports up to 350 lbs)
- Connected user mode & guest mode



How it works:

For best results, place the Pivotal Living™ Smart Scale on a hard, flat surface to ensure safe and accurate measurements.



*CARPET NOT RECOMMENDED

For greater accuracy, stand barefoot on the Pivotal Living™ Smart Scale. The scale cannot measure body fat percentage if you are wearing shoes or socks.



Align your heel and balls of your feet to the BIA (Bioelectrical Impedance Analysis) sensors. Stand still on the scale until an accurate measurement is displayed.

NEED HELP? PivotalLiving.com/Support

Care and Maintenance:



- Please do not use any type of cleaning agent that would otherwise lead to damaging the tempered glass or electrode pads.
- Store the unit horizontally, away from direct sunlight and exposure to extreme temperature changes.

WARNING:

- **DO NOT USE** if you have a pacemaker or an electronic medical implant.
(The Pivotal Living Smart Scale uses a low-level electrical signal through the body, which may interfere with its operation)
- **DO NOT USE** if you are pregnant
- Not intended for use by small children

SAFETY & PRECAUTIONS:

- Do not step on the platform when wet
- Do not place the scale on slippery surfaces, such as a wet floor
- Avoid any excessive impact to the scale platform
- Never submerge the scale in water

NEED HELP? PivotalLiving.com/Support

FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation

Industry Canada

This device complies with Industry Canada licence-exempt RSS standard(s):

Operation is subject to the following Two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

The term "IC: " before the certification/registration number only signifies that the Industry Canada technical specifications were met.

This product meets the applicable Industry Canada technical specifications.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes : (1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement