1 Parts in the package ①

#	Denomination
1	NOVO IP unit
2	AC Adaptor
3	LAN cable*
4	SMILE unit
5	SMILE accessories (not shown)

* Depending on model/market

2 Connectors, buttons, LEDs etc. ②

#	Function
PS1	On/Off button (Power)
C1	IP (RJ45 Ethernet)
C2	DC (Micro USB)
B1	Alarm button (red)
B2	Reset button (green/yellow)*
B3	Extra button (yellow/green)*
M1	Microphone
L1	Network error
L2	AC/DC error
L3	Accumulator error
L4	Under Alarm button B1 (red)
L5	Under Reset button B2 (green/yellow)*
L6	Under Extra button B3 (yellow/green)*

* Depending on model/market

2.1 Sounds

These most common sounds used in NOVO IP are:

- OK – “Blipp”
- Disconnection – “Di-Du-Da”
- Error – “Blurp”
- Pre-call signal – “Melody”

3 Setting up the unit

1. Connect the network cable to the router/switch or equivalent and the network cable to the IP connector (C1).
2. Connect the AC adaptor to the wall socket.
3. Connect the AC adaptor’s Micro USB connector to the DC connector (C2).
4. Power on the unit by pressing PS1.
5. Test the Radio coverage by entering the Service Menu, see 4.2.
6. Make a test call to the alarm central to ensure proper operation.

3.1 Mount SMILE ③

SMILE is a portable alarm trigger that can be worn either in a necklace or with wristband. The supplied SMILE is already connected to NOVO IP and ready to use.

Picture ③ demonstrates how to mount the necklace (A) or wristband (B).

3.2 Send an alarm from NOVO IP

Press the Alarm button (B1). A pre-call signal is played whereafter NOVO IP starts dialling. During the pre-call signal it is possible to cancel the alarm by pressing the Reset button (B2).

4 Control Mode and Service Menu

NOVO IP is normally in Idle Mode, i.e. is not processing any user-initiated (e.g. user alarm) nor automatic event/action (e.g. test alarm). Control Mode and Service Menu are ways of easily configuring the unit.

4.1 Control Mode (CM) functions

Enter CM by shortly pressing button PS1 when unit is in Idle mode. When the unit is in CM all LEDs L4-6 are lit. The unit returns automatically to Idle Mode after after function execution or 5 seconds of inactivity.

#	Control Mode functions
1	Toggle Alarm button (B1) illumination
2	Change volume (conversation/signal)
3	Enter Service Menu

4.1.1 Toggle Alarm button (B1) illumination

Press and hold B1 for more than 3 seconds to toggle B1 illumination.

4.1.2 Change volume

Press and hold B2 for more than 3 seconds to toggle conversation level. Release B2 and press again (and hold) within 5 seconds to toggle signal level.

4.1.3 Enter Service Menu

Press B3 for more than 3 seconds to enter Service Menu.

4.2 Service Menu (SM) functions

Enter Service Menu by pressing the On/Off button (PS1) shortly (thus entering Control Mode) and then button B3 for more than 3 seconds. In Service Menu LED4 is lit and LED6 blinks. When entering each function, NOVO IP plays the function sound.

Enter the function by pressing B1, step to next by pressing B3. After finishing programming in a selected function NOVO IP automatically returns back to the SM. Step to next function by pressing B3 or exit by pressing B2.

The list below assumes you have already stepped into SM and have navigated to the requested function.

#	Service Menu functions
1	Add transmitter
2	Radio coverage mode
3	Remove transmitter
4	Change Home/Away/Passive status

4.2.1 Add a transmitter

Add a transmitter by pressing B3 to step to radio position 1 indicated with one short “beep”, press B3 again to step to position 2, now indicated with two short “beeps”. Step to desired position by pressing B3 and when NOVO IP beeps the number of times corresponding to the desired radio position. Add the transmitter by activating it (e.g. pressing the Alarm button on SMILE). A short beep confirms the addition and also the transmitter battery status.

- Short beep – action OK and transmitter battery OK
- Long beep – action OK and transmitter battery low
- Error sound – the transmitter is already present at other radio position.

Exit function by pressing B2.

4.2.2 Radio coverage

In Radio coverage it’s possible to test the radio connection between a connected Transmitter and NOVO IP. When in Radio coverage, activate the connected transmitter. (e.g. pressing the Alarm button on SMILE) and NOVO IP responds with a “blipp” if the radio coverage is OK. On SMILE the LED blinks green to acknowledge radio coverage is OK.

End Radio coverage test by pressing B2.

4.2.3 Remove transmitter

Remove a transmitter by pressing B3 to step to radio position 1 indicated with one short “beep”, press B3 again to step to position 2, now indicated with two short “beeps”. Step to desired position by pressing B3 and when NOVO IP beeps the number of times corresponding to the desired radio position. To remove a transmitter press B1 > 3 seconds at the desired radio. A short beep confirms the removal.

Exit to SM by pressing B2.

4.2.4 Change Home/Away/Passive On/off

This function is further explained in NOVO - Technical Handbook, NE41 15013-02.

Change status by pressing B1 > 3 seconds.

4.3 Error indications

If acoustical warning is activated an error is indicated with a blipp and/or voice message. It is repeated regularly until problem is solved or the B2 button is pressed to silence the warning sound. Errors are also indicated with the indication LEDs L1 to L3 which will indicate according to below table until the problem is fixed.

Indication	Possible reason	Suggested action
Green (fixed)	No IP communication	Check local network equipment (modem, switch, router etc.) Check Radio coverage
Red (blink)	Mains failure/ broken AC adaptor	Check AC adaptor and/or cable.
Yellow (fixed)	Accumulator failure/not charging	Check accumulator
Green-Red- Yellow (blink)	Radio interference	Check interfering radio sources. Move wireless (DECT) phones, wifi transmitters etc. at least 1 (one) meter away from NOVO IP.

5 Important

5.1 Use and maintenance

- Temperature range: +5 - +55°C
- Do not damage the unit or its parts. If damaged, immediately contact authorized personnel.
- Do not expose to direct sunlight.
- Keep away from dust, moist and dirt.
- Do not drop, knock, twist or shake the device.
- Do not warm up the device or use it near fire.

5.2 Cleaning

- All parts in the NOVO IP kit can be cleaned with a mild soap solution and a damp cloth. Dry with a dry cloth.
- Strong chemicals, grease and other harsh substances must not be used when cleaning or handling the parts in the NOVO IP kit.
- NOVO IP must be disconnected from the power socket before cleaning.
- After cleaning, control that the home care phone works properly by sending a test alarm to the alarm receiver/central.

5.3 Safety Notes

- Read instructions prior to use.
- Always test the system per instructions prior to use.
- Always check the function of the product after making adjustments.
- This product may not be suitable for all persons and should not be a substitute for the routine visual monitoring protocol by caregiver.
- Must not be used in situations where a delay in the arrival of appropriate medical care, could lead to a potentially life-threatening situation.
- Our units are NOT intended for any life support device, thus intending a device whose malfunction may result in damage to a life.
- Check the device regularly and replace when necessary.
- Do not integrate with other systems other than those specified in this document.
- Always keep the device dry. Exposure to excessive moist can cause malfunction.
- The product will not cause electromagnetic disturbances under normal working conditions.
- The product can be placed near other products or devices as long as mechanical vibration is not present.
- Remove batteries if the unit is to be out of use or stored for an extended period of time.

5.4 Disposal

At the end of the product’s use life, please dispose of it at appropriate collection points provided in your country. For disposal or recycling information, please contact your local authorities or the Electronic Industries Alliance (EIA, www.eiae.org). In the European Union, the bin label indicates that this product should not be disposed of with household waste. It should be deposited at an appropriate facility to enable recovery and recycling or returned to your distributor.

6 Technical data

Measures	175 x 130 x 68 (H x W x D)
Weight	250 g (incl. accumulator)
Frequency, RF	916 MHz
Communication	IP
Max transmission power	50.000 µV/m
Recommended AC-adaptor	
US	NE31 07006-21 5V _{DC} 1000mA

7 US Notes

FCC ID: ZAGLF1410104

Changes or modifications not expressly approved by the party responsible for compliance could void the user’s authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Ulrik Lundberg
Managing Director

8 Contact

NEAT Electronics AB
Varuvägen 2
SE-246 42 Löddeköpinge
Sweden

Phone: +46 (0)46 70 70 65
Fax: +46 (0)46 70 70 87

infosweden@neat-group.com
www.neat-group.com/se/en
© NEAT Electronics AB 2016

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