

Philips Sonicare AirFloss
Ultra
AirFloss Ultra - Trial

1 nozzle

AirFloss Ultra



HX8341/01

A triple burst of better gingival health

Efficacy meets compliance for infrequent flossers

Philips Sonicare AirFloss Ultra is clinically proven as effective as floss in improving gingival health.* In fact, up to 97% of inconsistent flossers who tried AirFloss Ultra showed improved gum health.*

Improves patients' oral health

- Up to 99.9% plaque removal***
- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings


The easiest way to effectively clean between teeth

- An easy way to improve cleaning between teeth
- Takes just 60 seconds: Point, Press, Clean!
- Fresh and clean with mouthwash

PHILIPS
sonicare


Highlights

Up to 99.9% plaque removal




Removes up to 99.9% of plaque from treated areas.**

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Better gum health



Philips Sonicare AirFloss Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Cleans in just 60 seconds



With the AirFloss Ultra, patients can clean their entire mouth in less than 60 seconds a day. Patients simply select their desired burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh and clean with mouthwash



Patients can use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities between



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Ultra helps prevent cavities from forming in the spaces between patients' teeth.

High performance nozzle



The new Airfloss Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Triple burst customization



AirFloss Ultra's powerful spray can be customized to your patient's liking with single, double, or triple bursts for each press of the activation button.

Specifications

Items included

- AirFloss Ultra handle: 1
- Nozzle: 1
- Charger: 1

Design and finishing

- Color: Gray

Cleaning performance

- Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction
- Performance: Removes up to 99.9% of plaque
- Health benefits: Cleans entire mouth in 60 sec, Improves gum health in 2 weeks
- Customizable settings: Auto-burst or manual burst, Single, double, triple burst
- For best results: Change nozzle every 6 months

Ease of use

- Handle: Slim ergonomic design
- Battery indicator: Handle LEDs show charge level
- Battery Life: Up to 2 weeks
- Nozzle system: Easily snaps on and off

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty



Issue date 2017-08-09

Version: 6.0.1

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* When used in conjunction with a manual toothbrush and antimicrobial rinse in patients with mild to moderate gingivitis. AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. For more information, please visit www.philips.com/airfloss/faq.