



Philips
BlueControl Advance

15 minutes treatment per plaque

Wearable blue light therapy
Clinically proven efficacy
Without chemicals or cortisone



PSK2103

Blue LED Light Therapy.

Changing the way you experience psoriasis therapy.

Meet the wearable Blue LED Light Therapy for treating mild to moderate plaque psoriasis. Philips BlueControl Advance* is gentle and UV-free. Discover its clinically proven effectiveness with our 90 day money-back guarantee.*

Get great results without disrupting your life

- Practical. Get great results without disrupting your day to day life.

Be kind to your skin while treating psoriasis.

- Gentle. Be kind to your skin while treating psoriasis.

Experience the change in psoriasis therapy

- Effective. Experience the change in psoriasis therapy.



PHILIPS


Highlights

Effective.



Successfully challenge your psoriasis with this innovative device. Simply treat your plaques for 15 minutes per day, for at least four weeks and you should notice a significant improvement.* Blue LED Light Therapy with Philips BlueControl Advance helps your skin to renew itself properly by slowing the accelerated production of skin cells.* Philips BlueControl Advance stimulates the natural relief process of the skin and therefore improves plaque symptoms like scaling and redness. Its effectiveness is clinically proven.*


Gentle.



Blue LED Light Therapy is UV-free and safe for your skin.* This gentle alternative offers effective relief

without any additional ointments or medication. The blue LED light of Philips BlueControl Advance induces natural relief processes in your skin. It helps control rapid cell proliferation and inflammation and helps your skin to regenerate at normal speed, providing relief from plaque psoriasis symptoms such as redness, scaling and thickness.* BlueControl Advance is UV-free and chemical-free - studies confirm that the blue LED light in BlueControl Advance is not harmful to your skin.*

Practical.



Thanks to the innovative design of the small, battery-powered device, you can benefit from clinically proven Blue LED Light Therapy while following your normal daily routines.* Just follow three easy steps: Attach BlueControl Advance to your arm or leg with the textile fixation strap. Press the button to start the treatment. Treat daily for 15 minutes per plaque. Experience Blue LED Light Therapy and get effective relief from plaque psoriasis symptoms such as redness, scaling and thickness without the use of chemical substances or the potential risk of UV radiation.*

Specifications

Accessories

- USB cable: Recharging
- Charging adapter: Recharging
- Fixation strap: Secure fixation
- Charger: With UK and EU plug
- Storage option: Washable carrying pouch

Technical specifications

- Frequency: 50 - 60 Hz
- Type of lamps: LED
- Voltage: 100-240 V
- Power consumption: 5 W
- Battery: Rechargeable
- Battery type: Lithium ION
- Operating time (full to empty): Up to 7 treatments

Safety and Regulations

- CE 0344 Certified
- UV-free: No UV or near UV radiation
- Optical output classification: Risk Group 1 (EN 60601-2-57)
- Medical device class IIa
- MDD 93/42/EEC
- Max. daily dose: 45 J/cm² per plaque
- Power density: 40 mW/cm² (average), 600 mW/cm² (peak)

Logistic data

- Country of Origin: The Netherlands

Clinical trials

- Efficacy shown in trials: Clinically proven

Portable

- Rechargeable battery
- Power adapters (EU and UK)
- Small, compact design

Dimensions and weight

- Product dimensions: W/H/L: 77 x 23 x 125 mm
- Product weight: 165 grams

Design

- Device colour: Pearlescent white
- Strap design: Black and blue

Easy to use

- Treatment timer

Performance

- Clinically proven
- Treatment time: 15 minutes per psoriasis plaque

Service

- 90-day money back guarantee
- Warranty: 2-year limited warranty



Issue date 2017-12-13

Version: 5.1.1

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* 15 minutes of daily treatment for at least 4 weeks should lead to effective relief of plaque psoriasis symptoms - Krings L, Liebmann J, M Born, Leverkus M, Felbert von V. Data on file, July 2017

* Blue light as used in Philips BlueControl is not toxic for the skin - Awakowicz P et al. Contributions to Plasma Physics. 2009; 49 (9): 641-647

* Blue light helps control rapid cell proliferation and inflammation - Liebmann J, Born M, Kolb-Bachofen MV. Journal of Investigative Dermatology. 2010; 130: 259-269. - Fischer Metal. Experimental Dermatology. 2013; 22: 554-563

* Effectiveness of Blue Light Therapy in reducing psoriasis symptoms is clinically proven in 3 different trials. - Weinstabl A et al. Dermatology. 2011; 223 (3): 251-9 - Pfaff S et al. Dermatology 2015; 231: 24-34 - Krings L, Liebmann J, M Born, Leverkus M, Felbert von V. Data on file, July 20

* Philips BlueControl Advance is an EU CE-marked class IIa medical device for the treatment of mild psoriasis vulgaris. It is not available outside the EU. Before altering any prescribed medication or using the device, patients should consult their physician and read the user manual. Oct 2017.