

MULTI-MODE FOOT/BODY THERAPY
MASSAGER
LA-439



Questions or concerns?

Contact the Liteaid Customer Service
Department by email at
customerservice@liteaid.com.

Or send us a contact form at www.liteaid.com:

Instruction Manual

LITEAID™

Introduction

- Thank you for your purchase of the MULTI-MODE FOOT/BODY THERAPY MASSAGER by Liteaid USA, Inc.
- Before using please read the instruction manual carefully for operating and safety instructions.
- Feet are one of the most sensitive and complicated parts of the Human body. The body's 12 pulsation points interact directly and indirectly with the feet.
- Stimulating the various points on the foot can help in improving blood circulation, metabolism the immune system. By stimulating the many points on the foot, muscles are relaxed and helps improve general health.
- The MULTI-MODE FOOT/BODY THERAPY MASSAGER incorporates electrical pulse technology with the Chinese medical tradition. The MULTI-MODE FOOT/BODY THERAPY MASSAGER allows the ability to have 99 levels of electromagnetic waves with 25 different modes. This wide range allows the use of the MULTI-MODE FOOT/BODY THERAPY MASSAGER for complete satisfaction.
- The breakthrough in technology combines electric pulse technology with Chinese medical knowledge provides 50 different electromagnetic impulses for the direct stimulation of pressure points on the feet. Wave intensity can be adjusted in 99 steps to meet the unique requirements of each user.

Caution: Please read all instructions carefully before operating.

- If you have concerns regarding your health, consult your doctor before using this product.
- Individuals with pacemakers should consult a doctor before use.

- Those with the following conditions should consult their doctors before using this product:
 - Pregnancy
 - Heart Disease
 - Malignant Tumor
 - Conditions requiring rest
- If at any time soreness, rash or an increase in pain occurs, discontinue use and consult your physician.
- A massage session longer than 15 minutes is not recommended.

ACCESSORIES:



Electrode Pad



Electrode Cable



Adapter



Remote Controller




Slimming Belt

Maintenance:

- Towel dry the unit with out any water
- Please do not sit or stand on the unit.

Instructions for Use

Foot Massage:

- Place your feet on the massage zone outlined on the device. Please do not wear socks or any other covering on your feet. Ensure that soles of your feet touch the massage zone fully.
- The red digital display screen indicates four functions: a) Time. b) Mode. c) Body. d) preset time



(🕒time:25 🗂mode:00 🧑body:00 🦶sole:00),
- When the mode is Auto. The 🗂 icon blinks and every minute changes to the next mode. Pressing the 🦶 icon allows the selection of 25 modes
- Pressing the Body function 🧑 adjusts the intensity
- Press 🕒 function to cycle through the various functions Time/Mode/Body/Sole. Mode cannot be adjusted when auto mode.

Body Massage:

- Ensure that the region you wish to massage is clean.. Connect the electrode wires to the output jacks for the device. Paste the pads onto the area you wish to massage. The pads need to have a firm contact with the skin for the device to work effectively
- Use the device control panel to select the modes and intensity of massage.

Product Diagram and Features

Parts identification



- | | |
|--------------------------------|-------------------------------------|
| ① Power switch | ⑧ Intensity display screen for body |
| ② Auto/Fixed button | ⑨ Intensity display screen for sole |
| ③ Function option dec(-)button | ⑩ Infrared area |
| ④ Function option add(+)button | ⑪ DC Jack |
| ⑤ Function choice button | ⑫ Frequency display screen |
| ⑥ Time display screen | ⑬ Jacks for electrode wire |
| ⑦ Mode display screen | |