

Lifeline® Tempo™ Watch

A Personal Response Timepiece



Lifeline Systems
111 Lawrence Street
Framingham, MA 01702-8156
Tel: 1-800-635-6156
www.lifelinesystems.com



USER MANUAL

Thank you for selecting the Lifeline® Tempo™ Watch. Your Lifeline Tempo Watch provides you with the independence and peace of mind of the Lifeline service in an attractive watch for daily use. Your waterproof watch should always be worn in the shower or tub and should only be removed when you're using your other Lifeline Personal Help Button.

USING YOUR LIFELINE TEMPO WATCH

Press the gold button on the face of your Lifeline Tempo Watch when you are within range of your Lifeline Communicator. The Lifeline Communicator will **beep** and dial the Lifeline Response Center. A Lifeline Personal Response Associate will come on the line and ask you if you need help.

Always wear your Lifeline Tempo Watch. It is water resistant to 100 feet, so it will never be damaged by contact with water (the stem of the watch must be pushed in). Remember, if you are not wearing your Lifeline Tempo Watch, you may not be able to contact Lifeline when help is needed.

WHEN YOU RECEIVE YOUR LIFELINE TEMPO WATCH


When you first receive your Lifeline Tempo Watch, you will need to set the time and conduct a signal range test to determine which areas of your home will be covered by the service.

TO SET THE TIME

Pull out the watch stem on the side of the watch.

Rotate the stem until the watch indicates the correct time.

Press the watch stem back to its original position. Please note that the watch may not retain its water resistance if the stem is left in the outward position.



TO CONDUCT A SIGNAL RANGE TEST

Attach the Lifeline Tempo Watch to your wrist.

Press and hold the Reset Bar on your Lifeline Communicator to start the test. The Communicator will **beep** and say “**Start range test**” or “**Start test.**”

Walk around the house with your watch. In each room, press the Help Button on the watch to see if the Communicator responds with a **beep**. The “beep” indicates that from that location in your home, the Lifeline Tempo Watch will alert the Lifeline Response Center that help is needed.

Continue moving around the house and be sure to check the bathrooms.

When you’ve completed the signal range test, press the blinking Reset Bar or Reset Button to end the range test.



CARING FOR YOUR LIFELINE TEMPO WATCH

CHANGING THE WATCH BAND

If the wristband or wrist strap on your Lifeline Tempo Watch becomes damaged, or you would like to replace it with another one, please contact a local retailer which carries and services watches. The wristband width is a standard size (18mm) so it will accommodate an off-the-shelf replacement.

SERVICING YOUR LIFELINE TEMPO WATCH

Your Lifeline Tempo Watch includes a limited one year warranty against defects in manufacturing (see next page) and can also be serviced or replaced by Lifeline after the warranty has expired.

Important: do not bring your Lifeline Tempo Watch to a commercial jeweler for servicing or changing of the battery. The watch should only be serviced by Lifeline.

If you have any questions about your Lifeline Tempo Watch, or need to return it for service, please contact Lifeline Subscriber Services at 1-800-635-6156.

LIFELINE WARRANTY

This Lifeline Tempo Watch (but not the crystal, band, or strap) is warranted to the owner for a period of **ONE YEAR** from the date of purchase against defects in manufacturing by Lifeline Systems.

What Lifeline Will Do:

If this watch develops such a defect within the one year period, it will be repaired (i.e. a new or thoroughly reconditioned and inspected module will be installed or replaced and a watch of equal value and similar appearance will be supplied) at our option, and you provide proof of purchase and date of purchase.

This Limited Warranty Does Not Cover:

Lifeline Systems will not repair defects relating to servicing not performed by Lifeline Systems. Lifeline Systems will not provide any warranty service if your watch shows evidence that it has been tampered with, misused, abused, or altered; for example: moisture damage sufficient to affect the proper function of the watch; damage to the case; or visible cracking of the crystal.

LIFELINE SYSTEMS IS NOT LIABLE FOR INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES. Some states do not allow the exclusion or limitation of incidental, special, or consequential, damages, so the above exclusion or limitation may not apply to you.

Your Rights Under Implied Warranties and State Law:

THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER EXPRESSED WARRANTIES, OBLIGATIONS OR LIABILITIES. LIFELINE LIMITS THE DURATION OF ANY WARRANTY IMPLIED BY STATE LAW, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, TO ONE YEAR FROM THE DATE OF ORIGINAL PURCHASE. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.