
# Smart Stretch PTX-0001



Ó

Buzzer

 USB Charging Port

#### Specifications

above

Model: PTX-0001
Connection mode: Bluetooth
Battery voltage: 3.7V
Battery capacity: 210mAh
Working time: 420 mins
Charging time: 60 mins
Charging voltage: 5V
Compatible Operating System:

## How to wear Smart Stretch



1. Place your foot on the pad




 Use the Velcro strips to wrap around your shin and press the Velcro strips together



 Turn on the electronic device, adjust your own stretching program and start to stretch

#### How to use the electronic device



 Hold down the power button on the right side of the device for a few seconds



 Hear the beep, the display shows S(Stretch), R(Rest), Rep(Repetition), and Set.



3. Hold down the button for a few second.



4. The time of Stretch will be flashing, press the + button to increase the time. Then press the 

button to proceed to next items.



 Hold down the button for a few seconds when all items are set.

### 4 simple steps to Smart Stretch



Play Store.







APP icon

1. Download and install Smart Stretch App by scanning


QR code for searching "Ptex" from App Store/ Google







3, Make sure the Bluetooth is on and Connect the Smart Stretch.



4, Select programme and starting.

#### **FCC Statement**

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- -- Reorient or relocate the receiving antenna.
- -- Increase the separation between the equipment and receiver.
- -- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- -- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

